ARCTIC SUN™
TEMPERATURE MANAGEMENT SYSTEM

MODEL 5000 OPERATOR'S MANUAL

Simply Advanced™

The 24/7 Helpline is intended to assist healthcare professionals with technical questions they may have regarding the use of the ARCTIC SUN™ Temperature Management System. While the Helpline is staffed by licensed critical care nurses, they are not able to provide medical or nursing advice or to prescribe treatment.

Languages

English	3
Français (French)	16
Deutsch (German)	30
Italiano (Italian)	44
Español (Spanish)	58
Nederlands (Dutch)	72
Português (Portuguese)	86
Português-Brasil (Portuguese)	100
Ελληνικά (Greek)	114
Dansk (Danish)	130
Svenska (Swedish)	144
Suomi (Finnish)	158
Norsk (Norwegian)	172
Polski (Polish)	186
Magyar (Hungarian)	200
Český (Czech)	214
Türkçe (Turkish)	228
Română (Romanian)	242
Slovenčina (Slovakian)	256
Русский (Russian)	270
中文简体 (Chinese - Simplified)	286
繁體中文 (Chinese - Traditional)	300
한국어 (Korean)	314
Hrvatski (Croatian)	328

Table of Contents

Chapter 1 – Getting Started	3
Indications for Use	3
Warnings and Cautions	3
Complications	4
System Setup	4
Chapter 2 – Patient Therapy	4
Place ARCTICGEL™ Pads	4
Connect ARCTICGEL™ Pads	4
Temperature Probe Placement	4
Patient Therapy Selection	5
Initiate Normothermia (Control Patient & Rewarm Patient)	5
Initiate Hypothermia (Cool Patient and Rewarm Patient)	5
End Therapy	5
Chapter 3 – Normothermia Settings	6
Normothermia Settings	6
Chapter 4 – Hypothermia Settings	6
Hypothermia Settings	6
Chapter 5 – Advanced Setup	6
Chapter 6 – Alarms and Alerts	6
Alarms	6
Main Safety Alarms	7
Non-Recoverable Alarms	7
Recoverable Alarms	7
Alerts	7
Chapter 7 – Maintenance and Service	7
Cleaning and Maintenance	7
Preventative Maintenance	7
Inspect Fluid Delivery Line	8
Software Update	8
Service	8
Calibration	8
Appendix A: Product Specifications	8
Technical Description	8
Environmental Conditions	8
Disposal	8
ARCTIC SUN™ Temperature Management System Specifications	9
Appendix B: Symbols	10
Appendix C: Electromagnetic Compatibility	13
Appendix D: Warranty	15

Chapter 1 – Getting Started

Indications for Use

The ARCTIC SUN™ Temperature Management System is a thermal regulating system, indicated for monitoring and controlling patient temperature in adult and pediatric patients of all ages.

Warnings and Cautions

Warnings

- Do not use the ARCTIC SUN™ Temperature Management System in the presence of flammable agents because an explosion and/or fire may result.
- Do not use high frequency surgical instruments or endocardial catheters while the ARCTIC SUN™ Temperature Management System is in use.
- There is a risk of electrical shock and hazardous moving parts. There are no user serviceable parts inside. Do not remove covers. Refer servicing to qualified personnel.
- Power cord has a hospital grade plug. Grounding reliability can only be achieved when connected to an equivalent receptacle marked "hospital use" or "hospital grade".
- When using the ARCTIC SUN™ Temperature Management System, note that all other thermal conductive systems, such as water blankets and water gels, in use while warming or cooling with the ARCTIC SUN™ Temperature Management System may actually alter or interfere with patient temperature control.
- Do not place ARCTICGEL™ Pads over transdermal medication patches as warming can increase drug delivery, resulting in possible harm to the patient.
- The ARCTIC SUN™ Temperature Management System is not intended for use in the operating room environment.
- Medivance supplies temperature simulators (fixed value resistors) for testing, training and demonstration purposes. Never use this device, or other method, to circumvent the normal patient temperature feedback control when the system is connected to the patient. Doing so exposes the patient to the hazards associated with severe hypo- or hyper-thermia.

Cautions

- This product is to be used by or under the supervision of trained, qualified medical personnel.
- Federal law (USA) restricts this device to sale, by or on the order of a physician.
- Use only sterile water. The use of other fluids will damage the ARCTIC SUN™ Temperature Management System.
- When moving the ARCTIC SUN™ Temperature Management System always use the handle to lift the controller over an obstacle to avoid over balancing.
- The patient's bed surface should be located between 30 and 60 inches (75 cm and 150 cm) above the floor to ensure proper flow and minimize risk of leaks.
- The clinician is responsible to determine the appropriateness of custom parameters. When the system is powered off, all changes to parameters will revert to the default unless the new settings have been saved as new defaults in the Advanced Setup screen. For small patients (≤ 30 kg) it is recommended to use the following settings: Water Temperature High Limit $\leq 40^{\circ}\text{C}$ (104°F); Water Temperature Low Limit $\geq 10^{\circ}\text{C}$ (50°F); Control Strategy = 2. It is recommended to use the Patient Temperature High and Patient Temperature Low alert settings.
- Manual Control is not recommended for patient temperature management. The operator is advised to use the automatic therapy modes (e.g. Control, Patient, Cooling, Rewarming) for automatic patient temperature monitoring and control.
- The ARCTIC SUN™ Temperature Management System will monitor and control patient core temperature based on the temperature probe attached to the system. The clinician is responsible for correctly placing the temperature probe and verifying the accuracy and placement of the patient probe at the start of the procedure.
- Medivance recommends measuring patient temperature from a second site to verify patient temperature. Medivance recommends the use of a second patient temperature probe connected to the ARCTIC SUN™ Temperature Management System Temperature 2 input as it provides continuous monitoring and safety alarm features. Alternatively, patient temperature may be verified periodically with independent instrumentation.
- The displayed temperature graph is for general information purposes only and is not intended to replace standard medical record documentation for use in therapy decisions.

- Patient temperature will not be controlled and alarms are not enabled in Stop Mode. Patient temperature may increase or decrease with the ARCTIC SUN™ Temperature Management System in Stop Mode.
- Carefully observe the system for air leaks before and during use. If the pads fail to prime or a significant continuous air leak is observed in the pad return line, check connections. If needed, replace the leaking pad. Leakage may result in lower flow rates and potentially decrease the performance of the system.
- The ARCTIC SUN™ Temperature Management System is for use only with the ARCTICGEL™ Pads.
- The ARCTICGEL™ Pads are only for use with the ARCTIC SUN™ Temperature Management Systems.
- The ARCTICGEL™ Pads are non-sterile for single patient use. Do not reprocess or sterilize. If used in a sterile environment, pads should be placed according to the physician's request, either prior to the sterile preparation or sterile draping. ARCTICGEL™ Pads should not be placed on a sterile field.
- Use pads immediately after opening. Do not store pads once the kit has been opened.
- Do not place ARCTICGEL™ Pads on skin that has signs of ulceration, burns, hives, or rash.
- While there are no known allergies to hydrogel materials, caution should be exercised with any patient who has a history of skin allergies or sensitivities.
- Do not allow circulating water to contaminate the sterile field when patient lines are disconnected.
- The water content of the hydrogel affects the pad's adhesion to the skin and conductivity, and therefore, the efficiency of controlling patient temperature. Periodically check that pads remain moist and adherent. Replace pads when the hydrogel no longer uniformly adheres to the skin. Replacing pads at least every 5 days is recommended.
- Do not puncture the ARCTICGEL™ Pads with sharp objects. Punctures will result in air entering the fluid pathway and may reduce performance.
- If accessible, examine the patient's skin under the ARCTICGEL™ Pads often, especially those at higher risk of skin injury. Skin injury may occur as a cumulative result of pressure, time and temperature. Possible skin injuries include bruising, tearing, skin ulcerations, blistering, and necrosis. Do not place bean bag or other firm positioning devices under the ARCTICGEL™ Pads. Do not place positioning devices under the pad manifolds or patient lines.
- The rate of temperature change and potentially the final achievable patient temperature is affected by many factors. Treatment application, monitoring and results are the responsibility of the attending physician. If the patient does not reach target temperature in a reasonable time or the patient is not able to be maintained at the target temperature, the skin may be exposed to low or high water temperatures for an extended period of time which may increase the risk for skin injury. Ensure that pad sizing/coverage and custom parameter settings are correct for the patient and treatment goals, refer to the ARCTICGEL™ pad instructions for use for the appropriate flow rate. For patient cooling, ensure environmental factors such as excessively hot rooms, heat lamps, and heated nebulizers are eliminated and patient shivering is controlled. Otherwise, consider increasing minimum water temperature, modifying target temperature to an attainable setting or discontinuing treatment. For patient warming, consider decreasing maximum water temperature, modifying target temperature to an attainable setting or discontinuing treatment.
- Due to underlying medical or physiological conditions, some patients are more susceptible to skin damage from pressure and heat or cold. Patients at risk include those with poor tissue perfusion or poor skin integrity due to diabetes, peripheral vascular disease, poor nutritional status, steroid use or high dose vasopressor therapy. If warranted, use pressure relieving or pressure reducing devices under the patient to protect from skin injury.
- Do not allow urine, antibacterial solutions or other agents to pool underneath the ARCTICGEL™ Pads. Urine and antibacterial agents can absorb into the pad hydrogel and cause chemical injury and loss of pad adhesion. Replace pads immediately if these fluids come into contact with the hydrogel.
- Do not place ARCTICGEL™ Pads over an electrosurgical grounding pad. The combination of heat sources may result in skin burns.
- If needed, place defibrillation pads between the ARCTICGEL™ Pads and the patient's skin.
- Carefully remove ARCTICGEL™ Pads from the patient's skin at the completion of use. Discard used ARCTICGEL™ Pads in accordance with hospital procedures for medical waste.
- The USB data port is to be used only with a standalone USB flash drive. Do not connect to another mains powered device during patient treatment.

- Users should not use cleaning or decontamination methods different from those recommended by the manufacturer without first checking with the manufacturer that the proposed methods will not damage the equipment. Do not use bleach (sodium hypochlorite) as it may damage the system.
- Medivance will not be responsible for patient safety or equipment performance if the procedures to operate, maintain, modify or service the Medivance ARCTIC SUN™ Temperature Management System are other than those specified by Medivance. Anyone performing the procedures must be appropriately trained and qualified.

Complications

Targeted Temperature Management can trigger pathophysiological side-effects on the body such as but not limited to: cardiac dysrhythmia, electrolyte and pH balance, metabolic changes, hemodynamic changes, blood-glucose balance, infection, shivering, and can affect the coagulation, respiratory, renal and neurological systems. The controlling of patient temperature should only be performed under the supervision of a qualified healthcare professional.


System Setup

Unpack

- Unpack the ARCTIC SUN™ Temperature Management System Control Module and accessories.
- Allow the control module to remain upright for at least 2 hours prior to completing the installation and setup procedure in order to allow the chiller oil to settle. Damage to the chiller compressor may result otherwise.

Connections

- Use only Medivance approved cables and accessories with the ARCTIC SUN™ Temperature Management System Control Module. Connect the Fluid Delivery Line, Patient Temp 1 cable, Patient Temp 2 cable (optional) and Fill Tube to the back of the control module.
- Plug the Power Cord into the wall outlet. Position ARCTIC SUN™ Temperature Management System so that access to the power cord is not restricted.

Power On

- Turn the power ON by activating the Power Switch.
- The control module will automatically go through a brief self-test of the independent safety alarm.
- A New User Training module option is available from the start up screen.
- When the self-test is complete, the **Patient Therapy Selection** screen will appear on the control panel.

Fill Reservoir

- Fill the reservoir with sterile water only.
- Four liters of water will be required to fill the reservoir at initial installation.
- Add one vial of ARCTIC SUN™ Temperature Management System Cleaning Solution to the sterile water.
- From the **Patient Therapy Selection** screen, press either the **Normothermia** or **Hypothermia** button, under the New Patient heading.
- From the **Hypothermia** or **Normothermia** therapy screen, press the **Fill Reservoir** button.
- The **Fill Reservoir** screen will appear. Follow the directions on the screen.

Functional Verification

Perform the following functional verification procedure after initial setup and installation of the control module.

- Power On the control module
- From the **Patient Therapy Selection** screen, press the **Hypothermia** button to display the **Hypothermia** therapy screen.
- From the **Hypothermia** therapy screen, press the **Manual Control** button to open the **Manual Control** window.
- Use the Up and Down arrows to set the **Manual Control** water target temperature to 40°C and the duration to 30 minutes.
- Press the **Start** button to initiate **Manual Control**. Allow at least 3 minutes for the system to stabilize.
- Monitor the flow rate and water temperature in the **System** status area on the **Hypothermia** therapy screen.
- Verify that the flow rate reaches at least 1.5 liters/minute.
- Verify that the water temperature increases to 30°C.
- Press the **Stop** button.
- Set the **Manual Control** water target temperature to 4°C and the duration to 30 minutes.
- Press the **Start** button to initiate **Manual Control**.
- Monitor the flow rate and water temperature in the **System** status area of the **Hypothermia** therapy screen. Verify that the water temperature drops to 6°C.
- Press the **Stop** button to stop **Manual Control**
- Press the **Cancel** button to close the **Manual Control** window
- Power Off the control module.

Chapter 2 – Patient Therapy

Place ARCTICGEL™ Pads

Read the Instructions for Use that accompany the ARCTICGEL™ Pads. Examine each pad for damage prior to placement.

Connect ARCTICGEL™ Pads

While holding the pad line tubing, insert the clear pad line connector into the Fluid Delivery Line manifold. Do not press or squeeze the wings when connecting. The connector will click into place. Check the surface of the device for mechanical damage prior to use.

Temperature Probe Placement

Patient temperature control with the ARCTIC SUN™ Temperature Management System requires patient temperature feedback provided by an indwelling patient temperature probe connected to the Patient Temperature 1 connector on the back of the Control Module. Any commercially-available Yellow Springs Instrument 400 Series (YSI 400) compatible patient temperature probe can be connected to the ARCTIC SUN™ Temperature Management System. Refer to the manufacturer's Instructions for Use for the specific indications and temperature probe placement.

Patient Therapy Selection

Use the Patient Therapy Selection screen to initiate a New Patient, Continue a Current Patient, or access the Advanced Setup screen.

New Patient - Normothermia

Select **Normothermia** if the therapy goal is to maintain a patient temperature at a pre-defined target temperature for an indefinite period of time. Press the **Normothermia** button to display the **Normothermia** therapy screen.

New Patient - Hypothermia

Select **Hypothermia** to reduce and maintain a patient temperature at a set target temperature for a defined period of time, then slowly re-warm the patient at a controlled re-warming rate. Press the **Hypothermia** button to display the **Hypothermia** therapy screen.

Additional Protocol Option

Two additional protocols (**Hypothermia** or **Normothermia**) may be visible on the Patient Therapy Selection screen.

Current Patient

The Continue Current Patient button and the date and time that the current therapy was paused will display on the Patient Therapy Selection screen if a patient therapy was paused within the past 6 hours.

Press the Continue Current Patient button to resume a paused patient therapy.

Initiate Normothermia (Control Patient & Rewarm Patient)

Normothermia therapy is initiated and managed, and patient temperature is automatically controlled to a set target temperature from the **Control Patient** window in the **Normothermia** therapy screen. The **Control Patient** window displays the patient target temperature and the duration since the initiation of normothermia therapy.

To initiate Normothermia therapy:

- 1) From the Patient Therapy Selection screen, press the **Normothermia** button to display the **Normothermia** therapy screen.
- 2) The default patient target temperature will display in the **Control Patient** window.
- 3) To modify the patient target temperature, press the **Adjust** button to display the **Control Patient-Adjust** window.
- 4) **Control Patient** to: Use the Up and Down arrows to set the desired patient target temperature to control the patient.
- 5) **Rewarm at a Rate of**: Use the Up and Down arrows on the right of the screen to set the rewarming rate.
- 6) Press the **Save** button to save the new settings and close the **Control Patient-Adjust** window.
- 7) Press **Start**, in the **Control Patient** window to initiate therapy. You will hear a tone and then a voice stating "Therapy Started". Additionally, the **Control Patient** window and the ARCTIC SUN™ Temperature Management System icon will blink, indicating that therapy is in progress.

Initiate Hypothermia (Cool Patient and Rewarm Patient)

Hypothermia therapy is initiated and managed, and patient temperature is automatically controlled to a set target temperature from the **Cool Patient** and **Rewarm Patient** windows in the **Hypothermia** therapy screen.

The **Cool Patient** window displays the cooling phase patient target temperature and the length of time remaining in the cooling phase of the **Hypothermia** therapy.

The **Rewarm Patient** window displays the rewarming phase patient target temperature and the length of time remaining in the rewarming phase of the **Hypothermia** therapy.

To initiate hypothermia therapy:

From the Patient Therapy Selection screen, press the **Hypothermia** button to display the **Hypothermia** therapy screen.

1. Cool Patient Settings

- The default patient target temperature and duration will display in the **Cool Patient** window.
- To modify the patient target temperature and duration, press the **Adjust** button to display the **Cool Patient-Adjust** window.
- **Cool Patient To**: Use the Up and Down arrows on the left side to set the desired patient target temperature to cool the patient

- **Cool Patient For**: Use the Up and Down arrows on the right side to set the cooling duration to cool the patient before rewarming begins.
- Press the **Save** button to save the new settings and close the **Cool Patient-Adjust** window

2. Rewarm Patient Settings

- The default patient target temperature and duration will display in the **Rewarm Patient** window.
- To change the rewarming phase patient target temperature and rewarming rate, press the **Adjust** button in the **Rewarm Patient** window to display the **Rewarm Patient-Adjust** screen. Use the Up and Down arrows on the left side to set the desired final patient target temperature.
- **Rewarm Patient To**: Use the Up and Down arrows on the right side to set the desired final patient target temperature.
- **Rewarm at a Rate of**: Use the Up and Down arrows in the center of the screen to set the rewarming rate.
- **Rewarm Patient From**: When cooling a patient, adjustment of the **Rewarm Patient From** setting on the left side of the screen is disabled and defaults to the **Cool Patient** target temperature.
- When rewarming a patient, the **Rewarm Patient From** adjustment is enabled and the value can be modified. The **Rewarm Patient From** setting is the temperature to which the system is currently controlling the patient. The **Rewarm Patient From** temperature will automatically increase as the rewarming process continues. This feature allows the rewarming procedure to be optimized by allowing complete control of the rewarming ramp.
- Using the **Rewarm Patient From** temperature, the **Rewarm Patient To** temperature and the rewarming rate settings, the system will calculate and display the rewarming duration and the date/time at which the patient will reach the final rewarming target temperature.
- Press the **Save** button to save the new settings and close the **Rewarm Patient-Adjust** window.

3. Initiate Patient Cooling

- Press **Start**, in the **Cool Patient** window to initiate therapy. You will hear a tone and then a voice stating "Therapy Started". Additionally, the **Cool Patient** window and the ARCTIC SUN™ Temperature Management System icon will blink, indicating that therapy is in progress.

4. Initiate Patient Rewarming

- Upon completion of the cooling phase, there are two options for initiation of patient rewarming, either Automatically or Manually, depending on the **Rewarming Begins** setting in **Hypothermia Settings**.
- If **Rewarming Begins** is set to **Automatically**, the rewarming process starts automatically when the **Cool Patient** therapy is complete and the duration reaches zero.
- If **Rewarming Begins** is set to **Manually**, the rewarming process starts when the **Start** button is pressed in the **Rewarm Patient** window. The cooling process will continue until the **Rewarm Patient Start** button is pressed. An **Alert** will occur when the **Cool Patient** duration reaches zero. When the **Rewarm Patient** duration timer reaches zero, the system will continue to control the patient to the target temperature until the **Stop** button is pressed. Once in **Normothermia**, the timer will reset and begin tracking the duration of **Normothermia** therapy.

End Therapy

- From the **Normothermia** therapy or **Hypothermia** therapy screen, press the **Stop** button to terminate water circulation to the pads.
- Press the **Empty Pads** button, and follow the instructions to purge the pads of water.
- Disconnect the pads from the Fluid Delivery Line.
- Slowly and carefully remove pads from the patient skin.
- Discard the used pads in accordance with hospital procedures for medical waste.
- Press the power switch **Off**.

If power is lost while the power switch is in the On position, an audible alert will be issued until it is switched Off. This alerts the user that the treatment may have been accidentally stopped.

Chapter 3 – Normothermia Settings

Normothermia Settings

Use the **Normothermia Settings** screen to view the current settings and modify the settings for the following parameters. To modify any parameter setting, press the **Adjust** button to the right of the parameter.

Normothermia Settings screen parameters:

Water Temperature Settings

- Pre-Condition Water
- Manual Control
- High Water Limit
- Low Water Limit

Patient Temperature Settings

- High Patient Alert
- Low Patient Alert
- Control Strategy

Display Settings

- Temperature Units
- Temperature Units Adjust
- Patient Temp 2

To access the Normothermia Settings screen:

- 1) Press **Adjust** on the **Control Patient** window.
- 2) Press the **More** button on the **Control Patient Adjust** window.
- 3) The **Normothermia Settings** screen will be displayed.
- 4) To save the new settings as the current patient therapy settings, press the **Close** button. For instructions on saving the settings as the system defaults, see **Advanced Setup**.

NOTE: Verify that the Patient Temperature 1 probe is correctly positioned and correctly connected to the system. If the patient temperature change is less than 0.15°C after the first hour of therapy, the system will generate **Alert 116 – Patient Temperature 1 Change Not Detected**. If Alert 116 is not acknowledged after 5 minutes, the system will generate **Alarm 117 – Patient Temperature 1 Change Not Detected**. Alarm 117 will stop therapy and an audible alarm will sound. After Alarm 117 is acknowledged, therapy will require a restart.

Chapter 4 – Hypothermia Settings

Hypothermia Settings

Use the **Hypothermia Settings** screen to view the current settings and modify the settings for the following parameters. To modify any parameter setting, press the **Adjust** button to the right of the parameter.

Hypothermia Settings screen parameters:

Therapy Settings

- Cooling Begins
- Rewarming Begins

Water Temperature Settings

- Pre-Condition Water
- Manual Control
- High Water Limit
- Low Water Limit

Patient Temperature Settings

- High Patient Alert
- Low Patient Alert
- Control Strategy

Display Settings

- Temperature Units
- Temperature Units Adjust
- Patient Temp 2

To access the Hypothermia Settings screen:

- 1) Press **Adjust** on the **Cool Patient** window or the **Rewarm Patient** window.
- 2) Press the **More** button on the **Cool Patient Adjust** window or **Rewarm Patient Adjust** window.
- 3) The **Hypothermia Settings** screen will be displayed.
- 4) To save the new settings as the current patient therapy settings, press the **Close** button. For instructions on saving the settings as the system defaults, see **Advanced Setup**.

NOTE: Verify that the Patient Temperature 1 probe is correctly positioned and correctly connected to the system. If the patient temperature change is less than 0.15°C after the first hour of therapy, the system will generate **Alert 116 – Patient Temperature 1 Change Not Detected**. If Alert 116 is not acknowledged after 5 minutes, the system will generate **Alarm 117 – Patient Temperature 1 Change Not Detected**. Alarm 117 will stop therapy and an audible alarm will sound. After Alarm 117 is acknowledged, therapy will require a restart.

Chapter 5 – Advanced Setup

Use the **Advanced Setup** screen to view the current settings and modify the settings for the following parameters. To modify any parameter setting, press the **Adjust** button to the right of the parameter.

Location / Time Settings

- Language
- Number Format
- Current Time
- Date Format
- Current Date

The following functions can be initiated from the Advanced Setup screen.

- Download Patient Data: The Patient Data for the last 10 (ten) cases are stored on the ARCTIC SUN™ Temperature Management System hard drive. This data is maintained when the ARCTIC SUN™ Temperature Management System is powered down, or in the event of a total loss of power.
- Calibration
- Total Drain
- Save All Settings As Default
- Upload custom file

Additionally, the following information can be viewed in the Advanced Setup screen.

- Software Versions
- Last Calibration Date
- Next Calibration Due

To access the Advanced Setup screen:

- 1) Press **Advanced Setup** button on the **Patient Therapy Selection** screen.
- 2) The **Advanced Setup** screen will be displayed.

To access the Additional Protocol Selection screen:

Refer to the ARCTIC SUN™ Temperature Management System Help screens for information regarding additional protocol setup.

Chapter 6 – Alarms and Alerts

The ARCTIC SUN™ Temperature Management System safety system continually monitors the state of the device and the patient, and issues alarms or alerts to notify the user of conditions that may interfere with patient safety or system performance.

There are two types of conditions: **Alarms** and **Alerts**.

An **Alarm** notifies the user that a condition that may potentially pose an unsafe situation with respect to the patient or the device. An **Alarm** is a High Priority condition that requires immediate operator response.

An **Alert** informs the user about patient and device status without interrupting the procedure. An **Alert** is a Medium Priority condition that requires prompt operator response.

Alarms

An **Alarm** is denoted by an audio signal that repeats every 10 seconds until the **Alarm** is cleared. The **Alarm** screen will appear that displays the alarm number, alarm title, a description of the problem or conditions that triggered the alarm, and solutions and instructions for troubleshooting and resolving the alarm condition. If certain **Alarm** conditions are not acknowledged by the operator within 2 minutes, a **Reminder** tone will sound. All **Alarm** settings are maintained in the event of a power interruption.

Main Safety Alarms

While there are multiple alarms and safety features in the ARCTIC SUN™ Temperature Management System, there are five main safety alarms that will place the device into Stop mode until the condition is addressed.

Alarm	Specification
High Patient Temperature	39.5°C (103.1°F) During system warming
Low Patient Temperature	31.0°C (87.8°F) During system cooling
High Water Temperature	42.5°C (108.5°F)
Low Water Temperature	3.5°C (38.3°F)
System Self-Test Failure	At device power ON

Consult the service manual for detailed Alarm description

Each time the ARCTIC SUN™ Temperature Management System is powered On, a system self test for the independent safety alarm is automatically run. This test simulates a "water high temperature" fault situation on both the primary and secondary water temperature sensors. Both the primary and secondary safety systems must respond to the fault and be verified by the opposing safety system. If either safety systems do not respond appropriately either an alarm 80 or 81 will be issued. Contact Customer Support.

Non-Recoverable Alarms

If an Alarm condition occurs that prevents proper use of the device or proper patient treatment (such as the five main safety alarms discussed above), the system is placed into Stop mode and will not allow therapy to continue. This type of Alarm is known as Non-Recoverable. If this situation occurs, cycle the device power (turn device Off then On). If the alarm recurs contact Customer Support.

Recoverable Alarms

Other Alarms that temporarily Stop the device until the user is able to correct the cause and clear the Alarm are classified as Recoverable. If the condition that initiated the alarm is not addressed and problem persists, the Alarm will recur.

If a Recoverable Alarm occurs:

- 1) When an alarm is issued the device is placed into **Stop** mode.
- 2) Read the displayed instructions.
- 3) Note the **Alarm** number.
- 4) Press the **Close** button to clear the alarm.
- 5) Follow the instructions to correct the alarm condition. Perform the actions in the order listed until the alarm condition is resolved.
- 6) Once you have cleared the alarm, press the Start button in the therapy window to restart therapy. You will hear a tone and a voice stating "Therapy Started". Additionally, the active therapy window and the ARCTIC SUN™ Temperature Management System icon will blink.
- 7) If the condition does not resolve, contact Customer Support.

Alerts

Alerts are denoted by an audio signal that repeats every 25 seconds. The Alert screen will appear that displays the alert number, alert title, a description of the problem that triggered the alert, and solutions and instructions for troubleshooting and resolving the alert condition.

If an Alert occurs:

- 1) Read the displayed instructions.
- 2) Note the Alert number.
- 3) Press the **Close** button to clear the alert.
- 4) Follow the instructions to correct the alert condition. Perform the actions in the order listed until the alarm condition is resolved. If the condition does not resolve, contact Customer Support.
- 5) Refer to the ARCTIC SUN™ Temperature Management System Help screens for additional information regarding alarms and alerts.

Chapter 7 – Maintenance and Service

Cleaning and Maintenance

Routine cleaning and preventive maintenance should be performed on the ARCTIC SUN™ Temperature Management System control module every 6 months at a minimum. This consists of cleaning the external surfaces, accessories and chiller condenser, inspecting the device, and replenishing the internal cleaning solution that suppresses microorganism growth in the water reservoir and hydraulic circuit. See the ARCTIC SUN™ Temperature Management System Service Manual for additional information.

External Surfaces

- Clean the exterior body of the control module, fluid delivery lines, power cords and temperature cables using a soft cloth and mild detergent or disinfectant according to hospital protocol.

Condenser

- A dirty chiller condenser will significantly reduce the cooling capacity of the control module.
- To clean the condenser, wipe the dust from the exterior grill using a soft cloth. Depending on the quality of your institution's air, periodically remove the back cover and vacuum or brush the condenser fins. At a minimum the condenser fins should be cleaned annually. Maintenance activities should be performed by qualified personnel.

Device Inspection

- Periodically inspect the external areas of the device for damaged, loose or missing parts, and frayed or twisted power cords and cables.
- Discontinue using the device displaying one or more of the above conditions until the problem is corrected and has been verified to be operating correctly.

Replenish Internal Cleaning Solution

Contact Customer Service to order internal cleaning solution. For information regarding safe handling please refer to <http://www.medivance.com/manuals> for the Cleaning Solution SDS.

To replenish the internal cleaning solution:

- 1) Drain the reservoir.
 - Turn control module power Off.
 - Attach the drain line to the two drain valves on the back of the control module. Place the end of the drain line into a container. The water will passively drain into the container.
- 2) Refill the reservoir.
 - From the **Hypothermia** therapy screen or the **Normothermia** therapy screen, press the **Fill Reservoir** button.
 - The **Fill Reservoir** screen will appear. Follow the directions on the screen.
 - Add one vial of ARCTIC SUN™ Temperature Management System cleaning solution to the first bottle of sterile water.
 - The filling process will automatically stop when the reservoir is full.
 - Continue to replace the bottles of sterile water until the filling process stops.
 - When the **Fill Reservoir** process is complete, the screen will close.

Do not use Cleaning Solution that has passed the use by date listed on the bottle.

The Cleaning Solution must be stored inside the UV resistant pouch provided.

Preventative Maintenance

Use of the ARCTIC SUN™ Temperature Management System in excess of 2,000 hours, without conducting preventative maintenance, may result in failure of certain system components and failure of the system to function as intended. To maintain system performance, the ARCTIC SUN™ Temperature Management System requires periodic service of key components.

For additional information, please refer to <http://www.medivance.com/manuals>, call 1-800-526-4455 or contact your local BARD representative.

Inspect Fluid Delivery Line

1. Power On the system
2. From the Patient Therapy Selection screen press the **Hypothermia** button to display the **Hypothermia** therapy screen.
3. From the **Hypothermia** therapy screen, press the **Manual Control** button to open the **Manual Control** window.
4. Set the **Manual Control** water target temperature to 28°C and the duration to 30 minutes.
5. Connect a shunt to a set of fluid delivery line ports.
6. Press the **Help** button and then press the **Help Index** button. Select the topic **Maintenance and Service** and sub topic **System Diagnostics** then press the **Display** button. Verify that inlet pressure is -7 ± 0.2.
7. Repeat on all valves. If inlet pressure is out of range, replace the two valves that the shunt is connected to.
8. Ensure that the shunt is removed before device is put back in service.

Software Update

Software updates will be provided on a flash drive. Installation of software updates will be performed via the USB port on the front of the control module.

The software update feature will automatically initiate if the control module detects the proper files on a flash drive inserted in the USB port at power On.

To install software update:

- 1) Insert flash drive into the USB port.
- 2) An image of a timer will display while the software update is being installed, and will disappear when the software installation process is completed.
- 3) After installation, the new software version will display in the **Software Version** field in **Advanced Setup**.

Service

Contact Customer Support for technical support and customer service instructions to enable appropriately qualified technical personnel to repair those parts of the equipment that Medivance considers repairable.

Calibration

See ARCTIC SUN™ Temperature Management System Service Manual for calibration requirements and instructions. Calibration is recommended after 2,000 hours of operation or 250 uses, whichever occurs first.

Appendix A: Product Specifications

Technical Description

The ARCTIC SUN™ Temperature Management System is a thermoregulatory device that monitors and controls patient temperature within a range of 32°C to 38.5°C (89.6°F to 101.3°F). The ARCTIC SUN™ Temperature Management System consists of the Control Module and disposable ARCTICGEL™ Pads.

A patient temperature probe connected to the Control Module provides patient temperature feedback to an internal control algorithm which automatically increases or decreases the circulating water temperature to achieve a pre-set patient target temperature determined by the clinician.

The ARCTIC SUN™ Temperature Management System pulls temperature-controlled water ranging between 4°C and 42°C (39.2°F and 107.6°F) through the ARCTICGEL™ Pads at approximately 0.7 liter per minute per pad. This results in heat exchange between the water and the patient.

The ARCTIC SUN™ Temperature Management System Control Module is a CLASS I mobile device (Type BF, IPX0 and Mode of Operation – Continuous) per classification scheme of IEC 60601-1.

The ARCTIC SUN™ Temperature Management System Control Module meets both the electromagnetic interference and susceptibility requirements of IEC 60601-1, and is compatible with other equipment that also conforms to that standard. There is no known failure mode in the ARCTIC SUN™ Temperature Management System Control Module associated with electromagnetic interference from other devices. See the ARCTIC SUN™ Temperature Management System Service Manual for the full declaration regarding electromagnetic compatibility.

Environmental Conditions

Temperature Range

Operating: 10°C to 27°C (50°F to 80°F)

Storage: -30°C to 50°C (-20°F to 120°F)

At operating temperatures higher than 27°C (80°F), the refrigeration system's cooling capacity and therefore its ability to cool a patient is compromised.

Humidity Range (relative humidity, non-condensing)

Operating: 5% to 70%

Storage: 5% to 95%

Atmospheric Pressure Range: 60 kPa to 110 kPa

Disposal

Upon end of life, dispose of in accordance with local WEEE regulations or contact your local BARD Supplier or Distributor to arrange for disposal.

ARCTIC SUN™ Temperature Management System Specifications



Parameter	Specification
Therapy Modes	Normothermia: Control Patient, Rewarm Patient Hypothermia: Cool Patient, Rewarm Patient
Heater Capacity	2500 BTU/hr / 750 Watts
Circulating Fluid	Sterile Water
Reservoir Capacity	3.5 liters
Patient Probe Type	YSI 400 Series compatible
Patient Temperature Inputs	Patient Temp 1: control, monitor, alarm Patient Temp 2: monitor, alarm
Patient Temperature Display Range	10°C to 44°C 50°F to 111.2°F 0.1°C/°F increments
Patient Temperature Measurement Accuracy	±0.4°C (10°C to 32°C) ±0.2°C (32°C to 38°C) ±0.4°C (38°C to 44°C) Includes ±0.1°C external probe
Responses of the PCLCS (Physiologic Closed-Loop Control System)	Settling Time: ~4.5 hrs Relative Overshoot: <0.5°C Command Overshoot: <0.5°C Response Time: Warming (max) 33°C to 37°C: ~6 hrs Cooling 37°C to 33°C: ~2 hrs Steady State Deviation: 0 Tracking Error: 0 Note: All values derived from testing in simulated use.
Patient Temperature Control Range	32°C to 38.5°C 89.6°F to 101.3°F 0.1°C/°F increments
Water Temperature Display Range	3°C to 45°C / 37.4°F to 113.0°F 0.1°C/°F increments
Water Temperature Control Range (Manual)	4°C to 42°C / 39.2°F to 107.6°F 1°C/°F increments
High Water Temperature Limit	36°C to 42°C / 96.8°F to 107.6°F 1°C/°F increments
Low Water Temperature Limit	4°C to 25°C / 39.2°F to 77°F 1°C/°F increments
Time to heat water from 20°C to 37°C	8 minutes (approximate)
Sound Pressure	Alarm Tone: 70dB to 80dB at 1 meter, repeats every 10 seconds Alert Tone: 63dB to 71dB at 1 meter, repeats every 25 seconds Reminder Tone: 65dB at 3 meters, 0.5 seconds on/20 seconds off
Mains Input	100-120VAC, 50-60Hz, 11A 220-240VAC, 50-60Hz, 5.5A
Leakage Current	<300 µA
Operating Relative Humidity Range	5% to 70% non-condensing
Storage Relative Humidity Range	5% to 95% non-condensing
Operating Temperature Range	10°C to 27°C / 50°F to 80°F
Storage Temperature Range	-30°C to 50°C / -20°F to 120°F
Atmospheric Pressure Range	60 kPa to 110 kPa
Dimensions	Height: 35 inches (89 cm) Width: 14 inches (36 cm) Depth: 18.5 inches (47 cm)
Weight	Empty: 43 kg / 95 lbs; Filled: 47 kg / 103 lbs

Appendix B: Symbols

The ARCTIC SUN™ Temperature Management System Control module bears the following symbols:

Symbol	Standard Reference	Standard Title	Symbol Title	Explanatory Text
	ISO 15223-1 Reference no. 5.1.1	Medical Devices-Symbols to be used with medical device labels, labelling and information to be supplied	Manufacturer	Indicates the medical device manufacturer.
	ISO 15223-1 Reference no. 5.1.2	Medical Devices-Symbols to be used with medical device labels, labelling and information to be supplied	Authorized representative in the European Community	Indicates the Authorized representative in the European Community.
	ISO 15223-1 Reference no. 5.1.3	Medical Devices-Symbols to be used with medical device labels, labelling and information to be supplied	Date of manufacture	Indicates the date when the medical device was manufactured.
	ISO 15223-1 Reference no. 5.1.4	Medical Devices-Symbols to be used with medical device labels, labelling and information to be supplied	Use-by date	Indicates the date after which the medical device is not to be used.
	ISO 15223-1 Reference no. 5.1.5	Medical Devices-Symbols to be used with medical device labels, labelling and information to be supplied	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1 Reference no. 5.1.6	Medical Devices-Symbols to be used with medical device labels, labelling and information to be supplied	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	ISO 15223-1 Reference no. 5.1.7	Medical Devices-Symbols to be used with medical device labels, labelling and information to be supplied	Serial number	Indicates the manufacturer's serial number so that a specific medical device can be identified.
	ISO 15223-1 Reference no. 5.3.7	Medical Devices-Symbols to be used with medical device labels, labelling and information to be supplied	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.
	ISO 15223-1 Reference no. 5.3.8	Medical Devices-Symbols to be used with medical device labels, labelling and information to be supplied	Humidity limitation	Indicates the range of humidity to which the medical device can be safely exposed.
	ISO 15223-1 Reference no. 5.4.2	Medical Devices-Symbols to be used with medical device labels, labelling and information to be supplied	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	IEC TR 60878 Reference ISO 7010-M002	Graphical Symbols for Electrical Equipment in Medical Practice	Refer to instruction manual/booklet	To signify that the instruction manual/booklet must be read.
	IEC TR 60878 Reference no. 6050	Graphical Symbols for Electrical Equipment in Medical Practice	Model number	To identify the model number or type number of a product. In the application of this symbol, the model number or type number of the product should be accompanied with this symbol.

Symbol	Standard Reference	Standard Title	Symbol Title	Explanatory Text
	IEC TR 60878 Reference no.5334	Graphical Symbols for Electrical Equipment in Medical Practice	Defibrillation-proof type BF applied part	To identify a defibrillation-proof type BF applied part complying with IEC 60601-1.
	IEC TR 60878 Reference no. 5041	Graphical Symbols for Electrical Equipment in Medical Practice	Caution, hot surface	To indicate that the marked item can be hot and should not be touched without taking care.
	IEC TR 60878 Reference ISO 7010-W012	Graphical Symbols for Electrical Equipment in Medical Practice	Warning; Electricity	To warn of electricity.
	IEC TR 60878 Reference 7010-P017	Graphical Symbols for Electrical Equipment in Medical Practice	No pushing	To prohibit pushing against an object.
	IEC TR 60878 Reference IEC 7010-W001	Graphical Symbols for Electrical Equipment in Medical Practice	General Warning Sign	To signify a general warning.
	IEC TR 60878 Reference no. 5007 & 5008	Graphical Symbols for Electrical Equipment in Medical Practice	Power On/Off	Indicates power on/off switch
	IEC TR 60878 Reference no. 5021	Graphical Symbols for Electrical Equipment in Medical Practice	Equipotentiality	To identify the terminals which, when connected together, bring the various parts of an equipment or of a system to the same potential, not necessarily being the earth (ground) potential, e.g. for local bonding.
	Directive 2002/96/EC (WEEE)	Marking of electrical and electronic equipment in accordance with Article 11(2) of Directive 2002/96/EC (WEEE)	Waste Stream Disposal Status	Do not dispose of electronic products in the general waste stream
	21CFR801.15	Labeling; Prescription devices	Prescription Use Only	Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.
	AAMI ES 60601-1 IEC 60601-1-8 IEC 60601-10 IEC 806012-35 CSA C22.2 No. 6060101	Medical electrical equipment-Part 1: General requirements for basic safety and essential performance	ETL Monogram	Per ETL Intertek, models of the ARCTIC SUN™ Temperature Management System that bear the ETL Monogram confirm to AAMI ES60601-1, IEC 60601-1-8, IEC 60601-10, IEC 80601-2-35 and are certified to CSA C22.2 No. 60601-1.

Symbol	Standard Reference	Standard Title	Symbol Title	Explanatory Text
	ISO 7000 Reference no. 3650	Graphical symbols for use on equipment	Universal Serial Bus (USB)	Indicates USB port
	N/A	N/A	N/A	Indicates the drain valve
	N/A	N/A	N/A	Indicates that only sterile water should be used when filling the ARCTIC SUN™ Temperature Management System Control Module.
	N/A	N/A	N/A	Identifies Patient Temperature 1, the patient temperature probe input for monitoring and control.
	N/A	N/A	N/A	Identifies Patient Temperature 2, the patient temperature probe input for monitoring.
	N/A	N/A	N/A	Identifies Patient Temperature Out, the patient temperature output to an external hospital monitor.
	N/A	N/A	N/A	Identifies mechanical hazard
	N/A	N/A	N/A	Unit quantity
	ISO 15223-1 Reference no. 5.3.9	Medical Devices-Symbols to be used with medical device labels, labelling and information to be supplied	Atmospheric Pressure Limit	Indicates the range of atmospheric pressure to which the medical device can be safely exposed.

Appendix C: Electromagnetic Compatibility

Medical electrical equipment needs special precautions regarding electromagnetic compatibility. Ensure that the ARCTIC SUN™ Temperature Management System is installed and used according to the electromagnetic compatibility information provided. The following are guidance and manufacturer's declarations regarding electromagnetic compatibility for the ARCTIC SUN™ Temperature Management System.

- The use of accessories or cables other than those specified or sold by Medivance (shown below) is not recommended. Use of unapproved accessories or cables may result in increased emissions or in decreased immunity of the ARCTIC SUN™ Temperature Management System.
- If the ARCTIC SUN™ Temperature Management System is used directly adjacent to or stacked with other equipment, the user should periodically observe the ARCTIC SUN™ Temperature Management System device to verify it operates normally in that environment.
- Portable and mobile RF communications equipment can affect Medical Electrical Equipment.

Accessories and cables approved by Medivance for use with the ARCTIC SUN™ Temperature Management System	Part #
Temperature In-Cable - Nellcor	735-02
Temperature In-Cable - BARD	735-03
Temperature In-Cable - Rusch	735-04
Temperature In-Cable - GE	735-05
Temperature In-Cable - Philips	735-06
Temperature Out-Cable - Nellcor	735-52
Temperature Out-Cable - BARD	735-53
Temperature Out-Cable - Rusch	735-54
Temperature Out-Cable - GE	735-55
Temperature Out-Cable - Philips	735-56
Power Cord, US, Canada, Mexico	733-00
Power Cord, Continental Europe	733-01
Power Cord, UK, Ireland	733-02
Power Cord, Australia, New Zealand	733-03
Power Cord, Mainland China	733-04
Power Cord, Brazil	733-05
Power Cord, Switzerland	733-07
Power Cord, South Africa	733-08
Transmission Interface Module (TIM) Kit	760-00
Transmission Interface Module (TIM)	761-00
RS232 Cord	762-00

1.1 EN/IEC 60601-1-2 Table 1		
Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The ARCTIC SUN™ Temperature Management System is intended for use in the electromagnetic environment specified below. The customer or the end user of the ARCTIC SUN™ Temperature Management System should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The ARCTIC SUN™ Temperature Management System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

1.1 EN/IEC 60601-1-2 Table 1 (cont.)		
Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
Emissions Test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Class A	The ARCTIC SUN™ Temperature Management System unit is suitable for use in all establishments other than domestic, establishments and those directly connected to the public low-voltage power supply network that supplies buildings for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

1.2 EN/IEC 60601-1-2 Table 2			
Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The ARCTIC SUN™ Temperature Management System unit is intended for use in the electromagnetic environment specified below. The customer or the end user of the ARCTIC SUN™ Temperature Management System unit should assure it is used only in such an environment.			
Immunity Test	IEC 60601 test level	Compliance Level	Intended Electromagnetic Environment
Electromagnetic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines ±1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential mode (line-line) ±2kV common mode (line-earth)	±1kV differential mode (line-line) ±2kV common mode (line-earth)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 seconds	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ARCTIC SUN™ Temperature Management System unit requires continued operation during power mains interruptions, it is recommended that the ARCTIC SUN™ Temperature Management System unit be powered from an uninterruptible power supply with sufficient capacity to run the unit for the maximum required time of interruption.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: UT is the a.c. mains voltage prior to application of the test level.

Appendix C: Electromagnetic Compatibility (continued)

1.3 EN/IEC 60601-1-2:2007 Sub-clause 5.2.2.2 Table 4:				1.4 EN/IEC 60601-1-2:2007 Sub-clause 5.2.2.2 Table 6:			
Guidance and Manufacturer's Declaration – Electromagnetic Immunity				Recommended separation distances between portable and mobile RF communications equipment and the ARCTIC SUN™ Temperature Management System unit			
Immunity Test	IEC 60601 test level	Compliance Level	Intended Electromagnetic Environment	RF communications equipment can effect medical electrical equipment. The ARCTIC SUN™ Temperature Management System unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ARCTIC SUN™ Temperature Management System unit can help prevent electromagnetic interference by maintaining a minimum distance between the portable and mobile RF communications equipment (transmitters) and the ARCTIC SUN™ Temperature Management System unit as recommended below, according to the maximum output power of the communications equipment.			
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80MHz	3Vrms 150kHz to 80MHz	Portable and mobile RF communications equipment should be used no closer to any part of the ARCTIC SUN™ Temperature Management System unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P} \quad 80\text{MHz to } 800\text{ MHz}$ $d = 2.3\sqrt{P} \quad 800\text{MHz to } 2.5\text{GHz}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended minimum separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 	150kHz to 80MHz $d = 1.2\sqrt{P}$	80MHz to 800MHz $d = 1.2\sqrt{P}$	800MHz to 2.5GHz $d = 2.3\sqrt{P}$	
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m 80MHz to 2.5GHz		0.01	0.12	0.12	0.23
				0.1	0.38	0.38	0.73
				1.0	1.2	1.2	2.3
				10	3.8	3.8	7.3
				100	12	12	23
				For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
NOTE 1 At 80MHz and 800MHz, the higher frequency range applies							
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from objects, structures and people.							
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ARCTIC SUN™ Temperature Management System unit is used exceeds the applicable RF compliance level above, the ARCTIC SUN™ Temperature Management System unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ARCTIC SUN™ Temperature Management System unit.							
^b Over the frequency range 150kHz to 80MHz, field strengths should be less than 3V/m.							

Appendix D: Warranty

Limited Warranty

BARD Medical Division, C. R. BARD, Inc. ("BARD") warrants to the original customer that each Arctic Sun™ Control Module ("Equipment") and ArcticGel Pad ("Disposable") will be free of defects in workmanship and materials for the period set forth in the labeling and if no such period is set forth, then one year from the date of purchase. If the Equipment or a Disposable proves to be so defective, such Equipment or Disposable may be repaired, replaced, refunded or credited, at BARD's option. A comprehensive post warranty service plan for Equipment is available for purchase. The warranty covers all parts and labor associated with defects in material and workmanship of the Equipment and Disposable. BARD will, at its discretion, determine if the Equipment or a Disposable is to be repaired on site, or at the BARD service center. If Equipment or a Disposable is to be returned for service, BARD will supply packaging materials and pay for ground shipping. However, it is the customer's responsibility to prepare and package the Equipment or Disposable for shipment at its own cost. Any expedited shipment request will be at the customer's expense. Any unauthorized Equipment or Disposable repair performed during the warranty period will void the warranty. All returns must be authorized in advance by BARD. The liability of BARD under this product warranty does not extend to any abuse, accidental damage, misuse, improper storage, alteration, further manufacture, packaging or processing, accidental damage or damage from misuse of Equipment, damage caused by using tap water rather than distilled water, routine maintenance, recalibration, or its repair by any person or entity not authorized by a BARD representative.

Disclaimers.

I. THE LIMITED WARRANTY PROVIDED ABOVE IS THE ONLY WARRANTY PROVIDED BY BARD AND IS IN LIEU OF ANY OTHER WARRANTIES, WHETHER EXPRESS, IMPLIED, OR STATUTORY, INCLUDING title, non-infringement, non-interference, interoperability, quality or condition, accuracy, completeness, merchantability, fitness for a particular purpose, or the absence of any defects, whether latent or patent.

II. CUSTOMER IS RESPONSIBLE FOR THE PERFORMANCE OF ITS RESEARCH AND THE CARE OF ITS PATIENTS, AND DETERMINING THE APPROPRIATENESS OF THE EQUIPMENT AND DISPOSABLES FOR ITS USE; CUSTOMER ACKNOWLEDGES THAT BARD is not responsible for the delivery of medical care or medical services to customer's patients or any other person. The Equipment and Disposables are a tool to be used by customer but do not replace professional skill or judgment. By providing the Equipment and Disposables to customer, neither BARD nor any employee of BARD is engaged in the practice of medicine. Customer is responsible for verifying the accuracy, completeness, and pertinence of any medical, patient, legal or other results, data or information entered in, received by, processed by, stored in, transmitted by, produced by, displayed by, or used in connection with the Equipment and Disposables. Customer assumes all risks and liabilities associated with the use of such information, whether such information is used alone or in combination with other information. BARD is not responsible for the performance, support, or any other aspect of customer's technology environment.

III. CUSTOMER ACKNOWLEDGES THAT COMPUTER AND TELECOMMUNICATIONS SYSTEMS ARE NOT FAULT-FREE AND OCCASIONAL PERIODS OF DOWNTIME OCCUR. BARD DOES NOT GUARANTEE THE USE OF THE EQUIPMENT AND ITS CONNECTIVITY FEATURES WILL BE UNINTERRUPTED, TIMELY, SECURE, OR ERROR-FREE OR THAT CONTENT LOSS WON'T OCCUR, NOR DOES BARD GUARANTEE ANY CONNECTION TO OR TRANSMISSION FROM THE CUSTOMER TECHNOLOGY ENVIRONMENT.

IV. Exclusions. The foregoing warranties will not apply to failure of any Equipment or Disposable caused by (i) customer's abuse, neglect or misuse or resulting from any failure to comply with the customer's responsibilities; (ii) malfunction or failure of any element of customer's technology environment or use other than as expressly authorized by BARD; (iii) customer's failure to maintain the physical environment for the Equipment (including normal maintenance) specified in the relevant documentation provided by BARD; (iv) malicious software not introduced by BARD; or (v) customer's failure to permit installation of any software updates or upgrades.

THE LIABILITY AND REMEDY STATED IN THIS LIMITED WARRANTY WILL BE THE SOLE LIABILITY OF BARD AND REMEDY AVAILABLE TO CUSTOMER WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE) OR OTHERWISE, AND BARD WILL NOT BE LIABLE TO CUSTOMER FOR ANY SPECIAL, INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF HANDLING OR USE OF BARD EQUIPMENT OR DISPOSABLES EVEN IF BARD HAS BEEN ADVISED OF THE POSSIBILITY OR LIKELIHOOD OF SUCH DAMAGES. IN NO EVENT WILL BARD'S LIABILITY UNDER THIS WARRANTY EXCEED THE PURCHASE PRICE PAID TO BARD BY CUSTOMER FOR SUCH EQUIPMENT AND DISPOSABLE.

Terms of Service

If Equipment availability is critical for patient treatment, it is the customer's responsibility to purchase back-up Equipment. Although BARD will attempt to promptly repair Equipment under warranty, the timeliness of repair is not guaranteed.

The customer is responsible for maintaining the Equipment according to the schedules and instructions in the documentation supplied with each system. BARD provides remote Technical Support from 8:00AM to 5:00PM Mountain Time and 24/7 emergency phone support. Contact customer Service for all service related requests. A detailed description of the problem or service required, the unit serial number, and contact information will be required to assist in providing efficient service of the unit. The customer must provide personnel to assist Technical Support with troubleshooting.

Loaned Equipment

If Equipment under warranty is returned for service, subject to availability, loaned Equipment may be available to the customer at no charge upon request for the duration of the service. The customer is responsible for setting up the loaned Equipment and to prepare and package the Equipment for return shipment according to the documentation. The customer is also responsible for the care and maintenance of the loaned Equipment and all accessories while the Equipment is in their possession. Any loss or damage will be the sole responsibility of the customer. Loaned Equipment must be returned within 7 days upon return of the repaired Equipment or rental charges will be applied at a rate of \$50 per day. Loaned equipment will be shipped ground at BARD's expense. Any expedited shipment request will be at the customer's expense.

Non-Warranty Service

Parts and service are available for a fee through customer Service for Equipment no longer under warranty. If requested, BARD can provide an estimate of the cost of factory repair. BARD will require a Purchase Order from the customer in order to initiate the repair service. If it is later determined the Equipment requires repair which exceeds the original estimate, BARD will contact the customer for authorization prior to proceeding with the repair.

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© 2021 BD. BD and the BD Logo are trademarks of Becton, Dickinson, and Company. All rights reserved.
BARD, ArcticGel, Arctic Sun, Medivance, and Simply Advanced are trademarks and/or registered trademarks of Becton, Dickinson, and Company.
Medivance, Inc.