
## Revisions


Rev from	Rev to	ECO #
0208	2010/07	5449-10

**Notes:**


1. BD Cat. Number 250160
2. Blank (Sheet) Size: Length: 8.5"      Width: 11"  
 Number of Pages: 24      Number of Sheets: 6  
 Page Size: Length 8.5"      Width 5.5"      Final Folded Size: N/A
3. Style (see illustrations below): # 5


#1


#2


#3


#4

#5

#6

#7

4. See Specification Control Number BALT0216446 for Material Information
5. Ink Colors: Printed two sides  Yes       No  
 No. of Colors: 1      PMS# 2755 Blue
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	<b>BD</b> Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 0216446JAA		Category and Description Package Insert, Directigen Neisseria Meningitidis Test	Sheet: 1 of 25 <hr/> Scale: N/A	A

# **BD Directigen™ *Neisseria meningitidis* Test**

For the detection of *Neisseria meningitidis* Groups A, C, Y and W135.

English:	pages	1 – 4	Italiano	pagine	11 – 14
Français:	pages	4 – 8	Español:	páginas	15 – 18
Deutsch:	Seiten	8 – 11	Tables:	pages	19 – 21
			Procedural Chart	page	22

 0216446JAA  
2010/07

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD representant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőtől. / Naudojamo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / **Kontakt din lokale BD-representant for mer informasjon.** / Aby uzyskać instrukcję użyj kownika, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / **Свържете се с местния представител на BD за инструкции.** / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндің жергілікті BD өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

## INTENDED USE

The **Directigen™ *N. meningitidis*** Groups A, C, Y and W135 Test is a presumptive latex slide agglutination test for the qualitative detection of antigens to *Neisseria meningitidis* Groups A, C, Y and W135 directly in cerebrospinal fluid (CSF), serum or urine. In addition, the test kit provides confirmation and serogrouping capabilities from suspected colonies of *N. meningitidis* Groups A/Y or CW135. Visible agglutination occurs when a sample containing any of these bacterial antigens is reacted with its respective antibody-coated latex beads.

## SUMMARY AND EXPLANATION

The diagnosis of bacteremia and meningitis, especially in young children, can be difficult. As many as 55% of children are seen by a physician and started on antibiotics before meningitis is detected.<sup>1</sup> Detection of microbial antigens in the CSF is a rapid and helpful method for diagnostic microbiology. It may be the single most important test in cases of partially treated meningitis since Gram stain and culture may be negative. Detection of a specific antigen is a clinically significant finding and a valuable aid when choosing antimicrobial therapy.<sup>2</sup>

*Haemophilus influenzae* type b, *Neisseria meningitidis*, and *Streptococcus pneumoniae* have been reported to be the three causative agents responsible for approximately 84% of cases of bacterial meningitis.<sup>3</sup>

Immunological methods for detecting characteristic exoantigens of pathogenic microorganisms in patient fluids (CSF, serum, urine) are typically faster than traditional methods such as culture. These techniques include counterimmunoelectrophoresis (CIE) and latex agglutination.<sup>4-7</sup> The latex agglutination procedure has been found to be more rapid and sensitive than CIE in the detection of purified antigen.<sup>8,9</sup>

## PRINCIPLES OF THE PROCEDURE

Specific antibodies are bound to the surface of latex beads. Latex particle aggregation becomes large enough to allow rapid visualization of positive agglutination in the presence of specific antigens. These specific soluble polysaccharide antigens accumulate in CSF, serum or urine as a result of infection by *N. meningitidis* Groups A, B, C, Y and W135, *H. influenzae* type b, *S. pneumoniae*, and Group B *Streptococcus*. These other antigens can be detected with the **Directigen Meningitis Combo Test Kit**, or other **Directigen Test Kits** (see "Availability").

## REAGENTS

**Directigen *N. meningitidis* Groups A, C, Y and W135:**

<b>Reagent 1</b>	(1.0 mL),	Anti- <i>N. meningitidis</i> Groups C and W135, Rabbit Antibody-Coated Latex Suspension, with 0.2% sodium azide (preservative),
<b>Reagent 2</b>	(1.0 mL),	Anti- <i>N. meningitidis</i> Groups A and Y, Rabbit Polyclonal Antibody-Coated Latex Suspension, with 0.2% sodium azide (preservative),
<b>Reagent A</b>	(0.5 mL),	Control Latex, Rabbit Immunoglobulin-Coated Latex Suspension, with 0.2% sodium azide (preservative),
<b>Control +</b>	(3.0 mL),	Polyvalent Positive Antigen Control, <i>N. meningitidis</i> Groups A/Y and CW135, <i>H. influenzae</i> type b, <i>S. pneumoniae</i> , and Group B <i>Streptococcus</i> Antigens, with 0.2% sodium azide (preservative),
<b>Control –</b>	(3.0 mL),	Negative Antigen Control, glycine buffered saline, with 0.2% sodium azide (preservative).

**Precautions:** For *in vitro* Diagnostic Use.

**Warning:** Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>10-13</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Prior to discarding, sterilize specimen containers and other contaminated materials by autoclaving.

**Reagents:** Do not use beyond the expiration date. Upon removal from the refrigerator, allow the reagents to warm to room temperature (15 – 30°C) before use.

To assure proper drop delivery, reagent dispensing bottles must be held vertically, dispensing one free-falling drop at a time.

**Warning:** Reagents contain sodium azide, which is very toxic by inhalation, in contact with skin and if swallowed. Contact with acids liberates very toxic gas. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

**Controls:** Do not use the kit if **Control +** and **Control –** do not yield appropriate results.

**Test Cards:** Cards must be flat for proper reactions. If necessary, flatten cards by bowing back in a direction opposite to that of the curl. Care should be taken not to finger-mark the test areas, since this may result in an oily deposit and improper test results. Use each card once and discard. Store cards in the original package in a dry area at room temperature. When spreading within confines of test areas, avoid scratching the card surface with the plastic stirrers. If the specimen does not spread to the outer perimeter of test area, use another test area of the card.

**Test Slide (Glass):** If the glass slide is used, disinfect the slide after each use, and wash before reuse. A phenolic disinfectant is suggested. Be sure all detergent and/or disinfectants are thoroughly rinsed from the slide before reuse.

**Rotation:** The recommended speed for mechanical rotation is  $100 \pm 2$  rpm, but rotation between 95 and 110 rpm does not significantly affect the results obtained. The rotator should circumscribe a circle approximately 2 cm in diameter in the horizontal plane. A moistened humidifying cover should be used to prevent drying of test specimens during rotation.

**Storage of Reagents:** Upon receipt, remove the carton containing reagents and refrigerate at  $2 - 8^{\circ}\text{C}$ . DO NOT FREEZE. Reagents should be recapped and returned to refrigeration when not in use, taking care not to mix color-coded caps.

#### SPECIMEN COLLECTION AND HANDLING

Refer to appropriate texts for details of specimen collection and handling procedures. Specimens should be tested as soon as possible; however, if the sample cannot be tested immediately, it should be stored at  $2 - 8^{\circ}\text{C}$  (for up to 48 h), or at  $-20^{\circ}\text{C}$ . Serum must be separated from whole blood prior to testing and storage.

#### SPECIMEN PRETREATMENT

The use of covered glass test tubes (borosilicate) in a  $100^{\circ}\text{C}$  (boiling) water bath is the most effective means of specimen treatment. Although glass test tubes in a heat block may be used, greater variation in the thorough heating of a sample may be noted as a result of variations in tube size and sample volume.

#### PROCEDURES

The testing area, all reagents, test specimens and test components should be at room temperature ( $15 - 30^{\circ}\text{C}$ ) when used.

**Materials Provided:** All materials as listed under "Reagents," work station and accessories.

**Materials Required But Not Provided:** Directigen Meningitis Combo Test Cards, Directigen Specimen Buffer, Rotator,  $100 \pm 2$  rpm, Humidifying cover, Micropipettor, 50  $\mu\text{L}$  delivery, Pipette tips. } (see "Availability")

Also required are the necessary equipment and labware used for preparation, storage and handling of CSF, serum and urine specimens.

#### Specimen Preparation (CSF):

1. Heat specimens for 3 min at  $100^{\circ}\text{C}$  (e.g., water bath or heat block) and allow to cool to room temperature before use.
2. For CSF specimens showing turbidity, centrifuge after heating for 10 min at  $1400 \times g$  prior to testing. The supernatant fluid is to be used as the test specimen.
3. Test specimens as described in "Test Procedure."

#### Specimen Preparation (Serum):

1. Dilute serum specimens of at least 0.6 mL 1:1 with Directigen Specimen Buffer and mix.
2. Heat specimens for 5 min at  $100^{\circ}\text{C}$  (e.g., water bath or heat block) and allow to cool to room temperature before use.
3. Using a wooden applicator stick, break up the protein "clot" formed, and vortex vigorously (approximately 5 s).
4. Centrifuge at a minimum of  $1400 \times g$  for 15 min.
5. Test supernatant fluid as described in "Test Procedure."

#### Specimen Preparation (Unconcentrated Urine):

1. Dilute urine specimens of at least 0.4 mL 1:1 with Directigen Specimen Buffer and mix.
2. Heat specimens for 5 min at  $100^{\circ}\text{C}$  (e.g., water bath or heat block) and allow to cool to room temperature before use.
3. Centrifuge at a minimum of  $1400 \times g$  for 10 min.
4. Test supernatant fluid as described in "Test Procedure."

#### Specimen Preparation (Concentrated Urine):

1. Urine samples that are turbid or have particulate material should be centrifuged at  $1400 \times g$  for 10 min before concentrating.
2. Urine samples may be concentrated 25-fold with a Minicon B-15 concentrator (Amicon Corporation, Danvers, Massachusetts).
3. Dilute at least 200  $\mu\text{L}$  of urine concentrate 1:1 with Directigen Specimen Buffer and mix.
4. Heat specimens for 5 min at  $100^{\circ}\text{C}$  (e.g., water bath or heat block) and allow to cool to room temperature before use.
5. Centrifuge at a minimum of  $1400 \times g$  for 10 min.
6. Test supernatant fluid as described in "Test Procedure."

#### Preparation for Confirmation of Colonies from Culture:

1. Locate suspected colonies on the agar surface from 18 – 24 h cultures that meet morphological and Gram stain characteristics of organisms that are appropriate for testing with Directigen meningitis latex reagents. **Note:** A Gram stain should be performed prior to testing to ensure that organisms are appropriate for testing with the Directigen meningitis latex reagents.
2. Pipette 0.5 mL (approximately 10 drops) of Control – reagent into a small glass test tube (10 x 75 mm or equivalent).
3. Select several (2 – 3) isolated colonies of similar morphology from the original or subculture plate using a sterile loop and suspend into the above tube to achieve a suspension equal to a McFarland #1 turbidity standard. **Note:** Over inoculating will yield an excessively heavy suspension which may cause erroneous results.
4. Heat the suspension for 3 min at  $100^{\circ}\text{C}$  (e.g., boiling water bath or heat block) and allow to cool to room temperature before testing.
5. Centrifuge at a minimum of  $1400 \times g$  for 10 min.
6. Test supernatant as described under "Test Procedure." **Note:** If atypical agglutination patterns are observed, refer to "Limitations of the Procedure."

**User Quality Control:** Include Control + and Control – testing with each batch of specimens tested as described in step 1, "Test Procedure."

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

### Test Procedure (CSF, serum, urine, concentrated urine, and colony confirmation):

Remove the reagents from refrigerated storage and slide the reagent tray into the slot in the kit Work Station.

For testing a specimen with Reagents 1 – 2, use the following procedure:

Refer to Procedural Chart Illustrations, page 22.

The Meningitis Combo test card or glass slide may be used to test samples. If using the Meningitis Combo card for more than one sample, label circles appropriately. Before use of the glass slide, thoroughly clean with a lint-free tissue.

1. Dispense one drop of **Control +** onto circles “+” in the NmC/W and the NmA/Y columns. Place one drop of **Control –** onto circles “-” in the NmC/W and the NmA/Y columns.
2. Micropipette 50 µL of test sample onto circle(s) of sample row S in the NmC/W and NmA/Y columns and circle A.
3. Holding the dispensing bottle by the cap, vigorously swing (without inverting) to thoroughly mix each reagent. Before uncapping the bottle, gently tap the base on the counter top to assure that no latex remains in the tip.
4. Dispense one drop of **Reagent A** onto circle A.
5. Dispense one drop of **Reagent 1** onto the circles in S, “+” and “-” in the NmC/W column. Dispense one drop of **Reagent 2** onto circles in S, “+” and “-” in the NmA/Y column.
6. Mix the samples and the latex reagents in each circle with a plastic stirrer, alternately using first one end of the stirrer and then the opposite end for the next circle. Discard the stirrer.
7. Place the test card or glass slide on a mechanical rotator and rotate at a speed of 100 ± 2 rpm for 10 min. Use a moistened humidifying cover to prevent evaporation.
8. Immediately at the end of 10 min, read the test results macroscopically under a high intensity incandescent light.

### INTERPRETATION OF TEST RESULTS

Record **Control +** and **Control –** test results first.

The **Control +** should yield strong agglutination within 10 min. The **Control –** should show no agglutination.

Agglutination in any of the circles containing **Reagent A** or **Control –** renders the reaction uninterpretable.

Record patient test results:

**Positive Test** – Should show agglutination. Any degree of agglutination present in one of the latex reagents indicates the presence of the corresponding antigen. *Agglutination in two or more latex reagents or the corresponding Directigen Reagent-A renders the reaction uninterpretable.*

**Negative Test** – Should show no agglutination.

### LIMITATIONS OF THE PROCEDURE

These latex agglutination tests are not intended as a substitute for bacterial culture. Confirmatory diagnosis of bacterial meningitis infection is only possible with appropriate culture procedures and should be routinely performed.

Samples with extremely low levels of antigen, for instance early in the course of infection, may yield negative results. Furthermore, samples with exceedingly high antigen concentrations may exhibit prozone effects producing inappropriately negative results. Although not extensively studied, prozone phenomena have only been observed in specimens seeded with extremely high antigen levels, and not in clinical specimens.

This assay is for the qualitative detection of antigens to *N. meningitidis* groups A, C, Y or W135 in CSF, serum or urine. In addition, the test kit can also be used with suspected colonies of *N. meningitidis* Groups A/Y or CW135. Performance with other specimen types has not been assessed. The **Directigen Neisseria meningitidis** test cannot be used with blood culture media.

*N. meningitidis* Groups C and W135 reagent is known to cross react in the presence of *Escherichia coli* K92. *Other cross reactions, uninterpretables and false positives may occur which are reduced or eliminated by proper specimen preparation (see “Specimen Pretreatment” and “Procedures-Specimen Preparation”).*

Untreated serum and urines may yield uninterpretable results with the latex reagents and should be treated as in “Procedures-Specimen Preparation.” CSF specimens or treated urine specimens showing turbidity should be centrifuged for 10 min at 1400 x g after heating, and prior to testing. The supernatant fluid is to be used as the test specimen.

Urine specimens (Processed according to “Procedures-Specimen Preparation”) yielding uninterpretable or suspected false positive results may be filtered using a MILLIPORE™ Millex-HA 0.45 µm filter (#SLHA-02505) and retested with the reactive test latex(es) and the corresponding control latex(es).

Colony confirmation samples that yield uninterpretable or atypical agglutination reactions may require a 1:10 dilution of the supernatant fluid with **Control –** reagent. Care should be taken to avoid bacterial suspensions greater than a McFarland #1 turbidity standard.

### EXPECTED VALUES

*Haemophilus influenzae* type b, *Neisseria meningitidis*, and *Streptococcus pneumoniae* have been reported to be the three causative agents responsible for approximately 84% of cases of bacterial meningitis.<sup>3</sup>

### PERFORMANCE CHARACTERISTICS

**Sensitivity:** **Directigen N. meningitidis** Groups A and Y, C and W135 latex reagents were compared with counter-immunoelectrophoresis for the ability to detect a partially purified *N. meningitidis* Group polysaccharide antigen in CSF, serum and urine. Commercial antisera specific for each *N. meningitidis* group were used in the CIE test. The **Directigen N. meningitidis** latex reagents were equal to or more sensitive than CIE in detecting *N. meningitidis* Group antigens.

A series of retrospective clinical evaluations of the *N. meningitidis* A and Y, C and W135 latex reagents were performed with specimens originally derived from patients with culture, CIE, or commercial latex reagent positive confirmed cases of *N. meningitidis* Groups A, C, and W135 infection. The data presented in **Table 1 (see page 19)** are an accumulation of these evaluations. In some cases CIE testing was not included in some of the evaluation.

**Specificity:** The specificity of the **Directigen N. meningitidis** Groups A and Y, C and W135 latex reagents was determined by testing retrospective and prospective CSF, serum and urine during a series of multiple site clinical evaluations (**Table 2, see page 20**). The latex reagents were tested against culture negative and nonindexed culture positive specimens. The data indicates specificity ranges of 97 – 100% depending on the latex reagent and specimen tested.

**Colony Confirmation:** Latex reagents were tested using suspensions of isolated colonies that met morphological characteristics of the organisms. Performance characteristics are listed in **Table 3 (see page 21)**.

## AVAILABILITY

Cat. No.	Description
250160	<b>Directigen™ Kits:</b> <i>N. meningitidis</i> Groups A, C, Y and W135 Test Kit, (30 patient tests, 60 controls).
255460	Group B Strep Test Kit, (30 patient tests, 60 controls).
252260	<i>H. influenzae</i> type b Test Kit, (30 patient tests, 60 controls).
255560	<i>N. meningitidis</i> Group B / <i>E. coli</i> K1 Tests Kit, (30 patient tests, 60 controls).
251960	<i>S. pneumoniae</i> Test Kit, (30 patient tests, 60 controls).
252360	Meningitis Combo Test Kit, (30 patient tests, 60 controls).
252480	<b>Directigen™ Meningitis Combo Test Card</b> (single use), Box 30.
256391	<b>Directigen™ Specimen Buffer</b> , 8 mL.
252350	<b>Directigen™ Meningitis Negative Control Reagent</b> , 9 mL, Box 4.
278051	<b>Macro Vue™ Card Test Rotator</b> with humidifying cover, 100 ± 2 rpm, automatic timer, friction drive, Model 51-II.

## REFERENCES

1. Rothrock, S.G., S.M. Green, J. Wren, D. Letaj, L. Daniel-Underwood, E. Pillar: Pediatric bacterial meningitis: is prior antibiotic therapy associated with unaltered clinical presentation?, *Annals Emerg. Med.* 21:146-152, 1992.
2. Fasola, E. and P. Ferriero: Laboratory diagnostic methods for central nervous system infections, *Neuro. Clinics N. Amer.* 3:279-290, 1992.
3. Morbidity and Mortality Weekly Report (CDC), Surveillance summary: bacterial meningitis and meningococemia, United States, 28:277-279, 1978.
4. Edwards, E.A.: Immunologic investigations of meningococcal disease, group-specific *Neisseria meningitidis* antigens present in the serum of patients with fulminant meningococemia., *J. Immunol.* 106:314-317, 1971.
5. Whittle, H.C., Egler, J.J., Tugwell, P. and Greenwood, B.M.: Rapid bacteriological diagnosis of pyogenic meningitis by latex agglutination., *Lancet* 2:619-621, 1974.
6. Ingram, D.L., Anderson, P. and Smith, D.H.: Counter-current immunoelectrophoresis in the diagnosis of disease caused by *Haemophilus influenzae* type b., *J. Ped.*, 81:1156-1159, 1972.
7. Newman, R.B., Stevens, R. W. and Gaafar, H.H.: Latex agglutination tests for the diagnosis of *Haemophilus influenzae* Meningitis., *J. Lab. Clin. Med.*, 76:107-113, 1970.
8. Data on file at BD Diagnostic Systems.
9. Tilton, R. C., Dias, F. and Ryan, R.W.: Comparative evaluation of three commercial products and counter-immunoelectrophoresis for the detection of antigens in cerebrospinal fluid, *J. Clin. Microbiol.*, 20:231-234, 1984.
10. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
11. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
12. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
13. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.

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Benex Limited