Revisions

SO 0191-5

Rev from	Rev to	ECO #
0109	201001	5254-10

Notes:

1. BD Cat. Number 441122
2. Blank (Sheet) Size: Length: 8.5" Width: 27.5"
 Number of Pages: 11 Number of Sheets: 1
 Page Size: Length 8.5" Width 5.5" Final Folded Size: 8.5" x 5.5"
3. Style (see illustrations below): # 4

4. See Specification Control Number TBD for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# Standard Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number:		Category and Description Package Insert, Vaginal Specimen Transport for the BD ProbeTec CT/GC Q ^X Amplified DNA Assays		Sheet: 1 of 12
61				Scale: N/A

BD Vaginal Specimen Transport for the **BD ProbeTec™ *Chlamydia trachomatis*/ *Neisseria gonorrhoeae* (CT/GC) Q^x** **Amplified DNA Assays**

English: pages 1 – 2 Français : pages 6 – 7
Dansk: side 2 – 3 Norsk: sider 7 – 8
Deutsch: Seiten 3 – 4 Svenska: sidan 8 – 9
Español: páginas 5 – 6

 0086

61
2010/01

Instrukcyo vám poskytnú miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinizi danişin. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкции свяжитесь с местным представителем компании BD. / Өзіндік жергілікті БД өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

The Vaginal Specimen Transport for the **BD ProbeTec** CT/GC Q^x Amplified DNA Assays offers a convenient way to collect and transport specimens from the patient to the laboratory. It is intended to be used for patient collection of vaginal swab specimens in a clinical setting according to the instructions provided. This transport system is for use with the **BD ProbeTec** CT/GC Q^x Amplified DNA Assays on the **BD Viper™** System in Extracted Mode. Patient-collected vaginal swab specimens are an option for screening women when a pelvic exam is not otherwise indicated. The Vaginal Specimen Transport for the **BD ProbeTec** CT/GC Q^x Amplified DNA Assays is not for home use.

SUMMARY AND EXPLANATION

The Vaginal Specimen Transport for the **BD ProbeTec** CT/GC Q^x Amplified DNA Assays is a sterile, self-contained system. The polyester-tipped swab on a plastic purple shaft facilitates obtaining the specimen and maintaining the organisms for nucleic acid extraction and amplification on the **BD Viper** System in Extracted Mode.

REAGENTS

Materials Provided: The Vaginal Specimen Transport for the **BD ProbeTec** CT/GC Q^x Amplified DNA Assays contains 100 units with one copy of patient collection instructions. Each unit contains one polyester-tipped swab in a transport sheath.

Storage Instructions: Store swabs at 5 – 25°C. Do not use beyond the expiration date. Swab sterility is guaranteed if primary swab container is intact.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

1. Pathogenic microorganisms including hepatitis viruses and Human Immunodeficiency Virus may be present in clinical specimens. "Standard Precautions"¹⁻⁴ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.
2. Optimal performance of the **BD ProbeTec** CT/GC Q^x Amplified DNA Assays requires proper specimen collection, handling and transport. Time and temperature conditions for storage must be maintained during transport.
3. Proper labeling should accompany each specimen to the laboratory.
4. Take care to avoid cross-contamination during the specimen handling steps. Specimens can contain high levels of organisms. If gloves come in contact with specimen, change gloves to avoid cross-contamination.
5. Specimens must be collected and processed before the expiration date of the Vaginal Specimen Transport.
6.  Vaginal Specimen Transport for the **BD ProbeTec** *Chlamydia trachomatis*/*Neisseria gonorrhoeae* (CT/GC) Q^x Amplified DNA Assays is for single use only; reuse may cause a risk of infection and/or inaccurate results.

SPECIMEN COLLECTION AND TRANSPORT

NOTE: Ensure that patients read the Patient Collection Instructions before providing them with a collection kit.

1. Wash hands with soap and water. Rinse and dry.
2. It is important to maintain a comfortable balance during the collection procedure.
3. Twist the cap to break the seal. Pull the cap with attached swab from the tube. Do not touch the soft tip or lay the swab down. If you touch or drop the swab tip or the swab is laid down, discard the swab and request a new vaginal swab.
4. Hold the swab by the cap with one hand so that the swab tip is pointing toward you.
5. With your other hand, gently spread the skin outside the vagina. Insert the tip of the swab into the vaginal opening. Point the tip toward your lower back and relax your muscles.

6. Gently slide the swab no more than two inches into the vagina. If the swab does not slide easily, gently rotate the swab as you push. **If it is still difficult, do not attempt to continue.** Make sure the swab touches the walls of the vagina so that moisture is absorbed by the swab.
7. Rotate the swab for 10 – 15 seconds.
8. Withdraw the swab without touching the skin. Place the swab in the tube and cap securely.
9. After collection, wash hands with soap and water, rinse, and dry.
10. Return tube with swab as instructed.

Specimen Storage and Transport:

Specimens collected with the Vaginal Specimen Transport must be stored and transported to the laboratory and/or test site at 2 – 30°C within 14 days of collection. Alternatively, specimens may be stored and transported to the laboratory frozen (-20°C) within 180 days of collection. Vaginal swab specimens must be transported to the laboratory in the swab tube provided.

For domestic and international shipments, specimens should be packaged and labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and etiologic agents/infectious substances. Time and temperature conditions for storage must be maintained during transport.

REFERENCES

1. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guide for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
3. U. S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/31/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.

Vaginal Specimen Transport for the BD ProbeTec *Chlamydia trachomatis*/ *Neisseria gonorrhoeae* (CT/GC) Q^x (Transportsæt til vaginalprøver til BD ProbeTec *Chlamydia trachomatis*/*Neisseria gonorrhoeae* [CT/GC] Q^x amplificerede DNA-analyser)

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Becton Dickinson France S.A.S.