BBL™ Collection and Transport Systems FALCON™ SWUBE™ Pinworm Paddle

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INTENDED USE

The **FALCON™ SWUBE™** Pinworm Paddle is for the collection and examination of the pinworm, *Enterobius vermicularis*. The paddle may also be used for the collection of fine particles, including pollen and spores.

PRINCIPLE OF THE PROCEDURE

The **SWUBE** Pinworm Paddle is a clear plastic paddle, coated on one side with an adhesive which collects and holds the specimen. The handle of the paddle is inserted through the cap of a 17 x 100 mm plastic tube. The paddle can be removed from the cap and used as a slide for direct microscopic examination.

Precautions: *in vitro* Diagnostic

Do not use if package is not intact.

SPECIMEN COLLECTION AND TRANSPORT

Observe established precautions against microbiological hazards throughout the procedures.

It must be assumed that all specimens collected might contain infectious organisms; therefore, all specimens should be handled with the necessary precautions in mind. Extreme CAUTION should be used in disposing of contaminated material. This product can be disposed of by autoclaving or incinerating.

Enterobius vermicularis, the causative agent of human enterobiasis, is a parasitic helminth that is widely distributed throughout the population. Because only five to ten percent of infected persons have demonstrable eggs in their stools, specimens must be collected from the perianal region. It is to this area that the gravid female migrates and lays large numbers of eggs while the host is resting.^{1,2} The sticky side of a **SWUBE** Pinworm Paddle is applied to the perianal area to collect the specimen. The paddle can be removed from the cap and examined directly under a low power microscope for the presence of eggs or adult female worms. Male worms are seen only upon purgation or at times with diarrhea.¹

Because of the migratory habits of the female pinworm, specimens are best obtained a few hours after the person has retired, between the hours of 9:00 p.m. and midnight, or in the morning immediately upon rising before bathing or bowel movement. Collection of three to six consecutive daily specimens is recommended.² Whenever possible, the laboratory should provide suitable collection containers and instructions for their proper use. Specimens should be obtained before therapy has begun. If therapy was initiated prior to collection of the specimen, this should be noted on forms accompanying the specimen.

Material should be collected from the site where the suspected causative agent is most likely to be found and in a clinical situation at the stage of infection most likely to yield a positive result. Specimens should be of sufficient quantity to permit complete examination.

Provisions must be made for prompt delivery to the laboratory. Specimen containers should be properly labeled with all pertinent information such as name of patient, time specimen taken, etc. The laboratory must be furnished sufficient clinical information to enable the technologist to select appropriate techniques; if at all possible, specimens should be transported directly to the laboratory.

For more detailed information, appropriate references should be consulted.²⁻⁴

PROCEDURE

Material Provided: **FALCON™ SWUBE™** Pinworm Paddle (See "Availability").

Materials Not Provided: Ancillary reagents and laboratory equipment as required.

Instructions: This method replaces collection by the cellulose tape slide and the petroleum jelly paraffin swab techniques.⁵

1. Hold the paddle by the cap and remove it from the tube.
2. Separate the buttocks and press the tacky surface against several areas of the perianal region.
3. Replace the paddle in the tube for transport to the laboratory; specimens should be refrigerated if examination is to be delayed for more than one day.
4. Remove the paddle from the cap and utilize it as a slide for direct low power microscopic examination of the specimen.

RESULTS

E. vermicularis eggs measure 50 – 60 µm in length and 20 – 30 µm in width and are elliptical-ovoid in shape with a thick, smooth colorless shell. Occasionally adult female worms may be observed. Female worms are transparent and measure 8 – 13 mm in length and 0.5 mm in diameter; they possess a ventrally curved posterior end, long, narrow, sharply pointed tail, and cephalic alae or wings at their anterior end.¹

AVAILABILITY

Cat. No. Description
220120 **FALCON™ SWUBE™** Pinworm Paddle, Case 200, Individually wrapped.

REFERENCES

1. Balows, A., W.J. Hausler, Jr., M. Ohashi and A. Turano (ed.). 1989. Laboratory diagnosis of infectious diseases, vol. 1. Springer-Verlag, New York.
2. Wentworth, B.B. (ed.). 1988. Diagnostic procedures for mycotic and parasitic infections, 7th ed. American Public Health Association, Inc., Washington, D.C.
3. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover and R.H. Tenover (ed.). 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.

4. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis.
5. Isenberg, H.D. (ed.). 1992. Clinical microbiology procedures handbook, vol. 2. American Society for Microbiology, Washington, D.C.

TECHNICAL INFORMATION: In the United States, telephone Technical Services, toll free (800) 638-8663.

© 1999 Becton Dickinson and Company
BBL, FALCON and SWUBE are trademarks
of Becton Dickinson and Company.

**BECTON
DICKINSON**

Becton Dickinson
Microbiology Systems
Becton Dickinson and Company
7 Loveton Circle
Sparks, Maryland 21152 USA

Becton Dickinson France S.A.
11 rue Aristide Bergès
38800 Le Pont de Claix, France

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