Revisions

Rev from	Rev to	ECO #
2009/09	2010/07	5388-10

Notes:

1. BD Cat. Number 441362
2. Blank (Sheet) Size: Length: 8.5" Width: 22"
 Number of Pages: 8 Number of Sheets: 1
 Page Size: Length 8.5" Width 5.5" Final Folded Size: 4.25" x 5.5"
3. Style (see illustrations below): # 4

4. See Specification Control Number 8085559 for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	<p style="font-size: small;">COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION</p>	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8085559		Category and Description Package Insert, BD ProbeTec Q ^X Urine Preservative Transport Kit	Sheet: 1 of 9 <hr/> Scale: N/A	A

BD Urine Preservative Transport for the BD ProbeTec™ Q^x Amplified DNA Assays

English: pages 1 – 2
Français : pages 2 – 3

Deutsch: Seiten 3 – 4
Español: páginas 4 – 5

 8085559
2010/07

Pokyny vám poskytnú miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндің жергілікті BD өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

The Urine Preservative Transport for the **BD ProbeTec Q^x Amplified DNA Assays (Q^x UPT)** is designed to preserve and transport *Chlamydia trachomatis* and *Neisseria gonorrhoeae* in male and female urine specimens from symptomatic and asymptomatic individuals prior to processing for analysis with the **BD ProbeTec CT Q^x** and **GC Q^x Amplified DNA Assays** on the **BD Viper™** System in Extracted Mode.

SUMMARY AND EXPLANATION

Testing for *C. trachomatis* and *N. gonorrhoeae* in screening programs and diagnostic situations is not always possible immediately after specimen collection. The Q^x UPT allows for an extended range of time and temperature conditions for specimen storage and transport.

REAGENTS

Materials Provided: Each Q^x UPT Kit contains 100 Q^x UPT tubes and 100 transfer pipettes. Each Q^x UPT contains 50 µL of **NAP Guard™** (≥ 742.5 mM K₂EDTA).

Materials Required But Not Provided: Sterile, plastic, preservative-free, urine collection cups and gloves.

Storage Instructions: Store kit at 2 – 33°C. Do not use past the expiration date.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

1. Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁻⁴ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.
2. The Q^x UPT is for use with male and female urine specimens to be tested using only the **BD ProbeTec CT Q^x** and **GC Q^x Amplified DNA Assays** on the **BD Viper System** in Extracted Mode. Do not use the Q^x UPT kit with other amplified DNA assays.
3. Optimal performance of the **BD ProbeTec CT Q^x** and **GC Q^x Amplified DNA Assays** requires proper specimen collection, handling and transport.
4. **NAP Guard** may be irritating to the eyes, skin and respiratory system. In case of contact with eyes, rinse opened eye immediately with plenty of water and seek medical advice if symptoms persist. After contact with skin, wash immediately with plenty of soap and water. If inhaled, seek medical attention in case of problems.
5. Take care to avoid cross-contamination during the specimen handling steps. Specimens can contain high levels of organisms. Ensure that specimen containers do not contact one another, and discard used materials without passing over open containers. If gloves come in contact with specimen, change gloves to avoid cross-contamination.
6. Specimens must be tested before the expiration date of the Q^x UPT.

SPECIMEN COLLECTION AND TRANSPORT

Urine Specimen Collection

1. The patient should not have urinated for at least 1 h prior to specimen collection.
2. Collect the specimen in a sterile, plastic, preservative-free specimen collection cup.
3. The patient should collect the first 20 – 60 mL of voided urine (the first part of the stream – NOT midstream) into a urine collection cup.
4. Label with patient identification and date/time collected.

Urine transfer to the Q^x UPT

NOTES:

Urine should be transferred from the collection cup to the Q^x UPT within 8 h of collection provided the urine has been stored at 2 – 30°C. Urine can be held for up to 24 h prior to transfer to the Q^x UPT provided that the urine has been stored at 2 – 8°C.

Wear clean gloves when handling the Q^x UPT and the urine specimen. If gloves come in contact with the specimen, immediately change gloves to prevent contamination of other specimens.

1. Remove the Q^x UPT and the transfer pipette from their packaging.

2. Label the Q^x UPT with the patient identification and date/time collected.
3. Hold the Q^x UPT tube upright and firmly tap the bottom of the tube on a flat surface to dislodge any large drops from inside the cap. Repeat if necessary.
4. Uncap the Q^x UPT and use the transfer pipette to transfer urine into the tube. The correct volume of urine has been added when the fluid level is between the purple lines on the fill window located on the Q^x UPT label. This volume corresponds to approximately 2.0 – 3.0 mL of urine. DO NOT overfill or under fill the tube.
5. Discard the transfer pipette. NOTE: The transfer pipette is intended for use with a single specimen.
6. Tighten the cap securely on the Q^x UPT.
7. Invert the Q^x UPT 3 – 4 times to ensure that the specimen and reagent are well mixed.

UPT Storage and Transport

Store and transport Q^x UPT urine specimens at 2 – 30°C and process within 30 days of collection. Specimens may be stored at -20°C for up to 180 days.

For domestic and international shipments, specimens should be packaged and labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and etiologic agents/infectious substances. Time and temperature conditions for storage must be maintained during transport.

REFERENCES

1. Clinical and Laboratory Standards Institute. 2005. Approved Guideline, M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, PA.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
3. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories. HHS Publication (CDC), 5th edition. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262, 17/10/2000, p. 021-0045.*

Urine Preservative Transport for the BD ProbeTec Q^x Amplified DNA Assays

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Benex Limited