Revisions

Rev from	Rev to	ECO #
0704	2010/07	5434-10

Notes:

- BD Cat. Number Various
- Blank (Sheet) Size : Length: 18" Width: 36"
 Number of Pages: 12 Number of Sheets: 1
 Page Size: Length 18" Width 6" Final Folded Size: 2.25" x 4"
- Style (see illustrations below): # 4

- See Specification Control No. 8810131 for Material Information.
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS #2755 Blue
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8810131JAA		Category and Description Package Insert, Febrile Antigens	Sheet: 1 of 13 Scale:	A

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INTENDED USE

Febrile Antigens are used in agglutination tests as an aid in the diagnosis of certain febrile diseases such as salmonellosis, brucellosis and rickettsial diseases. The patient's serum is tested directly for homologous antibodies by either a slide or tube agglutination test. These tests are qualitative and semi-quantitative. The rapid slide test is used primarily as a screening procedure especially useful when large numbers of sera must be examined. The tube test should be used to confirm positive results obtained by the slide test.

Febrile Control Antisera are used in agglutination tests in conjunction with Febrile Antigens as positive or negative controls.

SUMMARY AND EXPLANATION

Febrile antigen tests are serological applications of the classical Widal reaction devised for the diagnosis of typhoid fever and the Weil-Felix reactions where antigens prepared from a *Proteus* organism are used to detect related rickettsial antibodies.^{1,2} Serological diagnosis of patients suspected of having infectious diseases characterized by persistent fever is dependent upon demonstration of an agglutination reaction between the appropriate antigen and the patient's serum.

The natural response to the invasion by pathogenic organisms is the production of antibodies. This immune response is highly individualized and in addition to the host's physiological status and genetic capabilities, a number of other factors are involved with the production of antibodies to a particular stimulus. These include the antigenicity of the organism, the total amount introduced to the host and the route of introduction, and whether the host has had previous exposure to the organism. These factors will determine the rate of antibody formation, the amount of antibodies produced and their persistence in the circulatory system.

PRINCIPLES OF THE PROCEDURE

When the human body is invaded by any pathogenic microbiological agent, a variety of antibodies are formed. Among these antibodies are the agglutinins. Serum containing specific agglutinins in combination with homologous antigens under properly controlled conditions is capable of causing a visible agglutination. The degree to which agglutination occurs depends upon the concentration of the antigen, the number of antibodies present, the composition of the salt solution and the temperature.¹

Sera from normal patients may show positive agglutination with febrile antigens due to previous immunization, past infection or the presence of antibodies to related antigens. In general, the titers found in these instances will be lower and remain at a constant level. Titers detected as a result of active infection or recent immunization with an organism containing homologous antigens will usually be higher and tend to rise over a period of time.

It is, therefore, necessary to evaluate two or more serum samples taken at 3- to 5-day intervals after the onset of the disease. A progressive increase in titer is the prime evidence of recent infection or immunization.³

REAGENTS

Febrile Antigen Set Contains Febrile Antigens:

Salmonella O Group D (9-12) (Typhoid O)
Salmonella Flagellar Group a
Salmonella Flagellar Group b
Salmonella Flagellar Group d (Typhoid H)
Brucella abortus
Proteus OX19

Control Antisera:

Salmonella Somatic Polyvalent Antiserum (A, B, D)
Salmonella Flagellar Polyvalent Antiserum (a, b, d)
Brucella Positive Control Antiserum (AMS)
Proteus Polyvalent Antiserum (OXK, OX2, OX19)
Febrile Antigen Negative Control

Febrile Antigens are nonviable bacterial cells in 0.5% phenolized saline with added crystal violet, brilliant green and chemical stabilizers. The Salmonella Flagellar Antigens (Cat. No. 240785, 240834, 240835) and Francisella Tularensis Antigen (Cat. No. 241050) contain 0.5% formaldehyde in place of the phenol.

Febrile Control Antisera are preparations of rabbit antiserum with 30% - 50% glycerin as preservative.

Febrile Antigen Negative Control is Bovine Albumin (2.5%) with 50% glycerin as preservative.

Warnings and Precautions

For *in vitro* Diagnostic Use.

The Packaging of This Product Contains Dry Natural Rubber.

Warning: May cause sensitization by skin contact. Avoid contact with skin. Wear suitable gloves. If swallowed, seek medical advice immediately and show this container or label.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁴⁻⁷ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Prior to discarding, sterilize specimen containers and other contaminated materials by autoclaving.

Storage Instructions: On receipt, store at 2 to 8°C. Febrile Antigens are ready for use as supplied. Do not open until ready to use. Control sera are ready for use as supplied. Protect from evaporation and contamination. The expiration date on each vial applies to product in intact container stored as directed.

SPECIMEN COLLECTION

Using aseptic techniques, collect 5 – 10 mL of blood in a blood tube without anticoagulant (i.e. Vacutainer™ Brand Glass Evacuated Tube). This type of tube is recommended because the collection, clotting and centrifugation of the specimen can all be done without transferring the blood. Allow the blood to clot for 20 to 30 min. Loosen the clot from the side of the container with a clean glass rod or applicator stick, stopper the tube and centrifuge. Draw off the serum and refrigerate at 2 to 8°C until testing, or immediately freeze the serum at -15°C or below if the test is delayed more than 4 h.

PROCEDURE

Materials Provided: Depending upon which product is ordered, one or more of the products listed under "Availability" is provided.

Materials Required But Not Provided:

- 3" x 6" glass slide (circled)
- Pipettes
- Toothpicks, glass rods or other applicator sticks
- Rotary shaker
- 0.85% saline
- Water bath (37 ± 1°C, 49 ± 1°C)
- Refrigerator (2 – 8°C)
- Tubes

Test Procedure

NOTE: ALL MATERIALS AND EQUIPMENT MUST BE AT ROOM TEMPERATURE AT THE TIME OF TEST PERFORMANCE.

Directions for use should be followed carefully.

Slide Test

- Dispense serum (patient or control) onto slides. Use 3" x 6" (circled). Dispense serum samples into separate circles as follows:

Circle	Serum
1	0.08 mL
2	0.04 mL
3	0.02 mL
4	0.01 mL
5	0.005 mL

- Invert the bottle several times to mix antigen.
- Add one drop of antigen to each circle using the dropper supplied with the vial.
- Using a new toothpick or applicator for each circle, mix antigen and serum.
- Rock the slide for 1 min near a light source and observe for agglutination.
- Repeat steps 1 through 5 for each antigen.

Tube Test

- Prepare dilutions of patient and control sera.

Patient Sera

- Add 0.5 mL of 0.85% saline to each of 8 tubes.
- Place 0.8 mL of saline in a separate tube and add 0.2 mL of serum. Mix well with a separate pipette and add 0.5 mL from this tube to the first tube of 8 tubes prepared in step "a".
- Using a clean pipette for each patient serum, mix and transfer 0.5 mL from tube 1 to tube 2 and serial dilute through tube 7. Discard 0.5 mL from tube 7. Tube 8 is the antigen control.

Note: Reading References may be set up as described in step 3 of RESULTS, "Reading References."

Control Sera

Prepare serial dilutions of the control sera in 0.5 mL amounts in tubes in the following manner:

- Place 8 tubes in a rack for each serum to be tested.
 - Pipette 0.9 mL of 0.85% saline into the first tube of each row and 0.5 mL into each of the remaining tubes.
 - Add 0.1 mL of the serum to the first tube containing 0.9 mL of saline solution.
 - Mix well with a pipette and transfer 0.5 mL of the first tube to the second tube. Mix thoroughly.
 - Continue carrying the 0.5 mL of the serum dilution through tube 7. Discard 0.5 mL from tube 7 after mixing thoroughly. Tube 8 is the antigen control tube.
- Prepare Antigen Dilutions
 - Dilute Brucella, Salmonella, Proteus and Francisella antigens to one part antigen in forty-nine parts 0.85% saline (1:50), e.g., 0.25 mL antigen and 12.25 mL saline.
 - Mix Antigen and Serum
Add 0.5 mL of diluted antigen to each tube containing patient sera or control sera, as well as to the antigen control tube. Shake tubes for a period of 10 sec. The resultant dilutions are 1:20 through 1:1280, respectively.

4. Incubate.

Antigen	Incubation
Salmonella O Group A Salmonella O Group B	18 to 24 h in a 48 to 50°C water bath
Salmonella O Group D Salmonella Flagellar a Salmonella Flagellar b Salmonella Flagellar d	16 to 18 h in a 48 to 50°C water bath
Brucella abortus Brucella melitensis Brucella suis	48 h in a 37°C water bath
Proteus OX19 Proteus OX2 Proteus OXK	Two h in a 37°C water bath followed by overnight refrigeration at 2 to 8°C
Francisella tularensis	20 h in a 37°C water bath

User Quality Control

At time of use, apply both positive and negative antisera controls to check performance of the antigens, techniques and methodology.

The use of Francisella tularensis Antiserum (Cat. No. 240939) can aid in the control of both positive results with Francisella tularensis Antigen and cross reactions, especially with Brucella abortus Antigen.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Slide test

1. Record Agglutination

Degree of Agglutination	
100 %	background clear to slightly hazy 4+
75 %	background slightly cloudy 3+
50 %	background moderately cloudy 2+
25 %	background cloudy 1+
None	Negative

2. Determination of Antibody Titer

The titer is the reciprocal of the highest dilution producing 2+ agglutination.

Serum, mL	Correlation Dilution
0.08	1:20
0.04	1:40
0.02	1:80
0.01	1:160
0.005	1:320

Tube Test

1. Record Agglutination

Very gently agitate each tube and record agglutination as follows:

100 %	Clear liquid above agglutinins	4+
75 %	Slightly cloudy	3+
50 %	Moderately cloudy	2+
25 %	Cloudy	1+
Negative	Very cloudy with no agglutinins	

2. Measure Titer

The titer is the reciprocal of the highest final dilution producing a 2+ agglutination.

Tube Number	Dilution
1	1:20
2	1:40
3	1:80
4	1:160
5	1:320
6	1:640
7	1:1280

3. Reading References

The following Reading References table will be helpful in determining how the liquid portion of the tube test should appear for differing degrees of agglutination. This is for the liquid portion only.

Antigen, mL	Saline, mL	Liquid Reading
0.5 mL	0.5 mL	Negative
0.25 mL	0.75 mL	2+
0.125 mL	0.875 mL	3+

NOTE: If auto-agglutination of the antigen is suspected, prepare the "Reading References" mixtures exactly as indicated above. Agglutinated material in any or all of the tubes indicates that the antigen is unstable and should be discarded.

As in other agglutination tests, no single test or determination can be considered diagnostic. Many variables are involved (See "Limitations of the Procedure"). However, a fourfold rise in antibody titer demonstrable in paired acute and convalescent phase sera collected at approximately five-day intervals is strongly suggestive of a specific presumptive diagnosis. The chart below gives an approximate indication of the significance of serum titers.³

Disease	Febrile Antigen	Serum Agglutinins		
		Appear	Maximum	Titer and Significance
Typhoid Fever	Salmonella O Group D	7 to 10 days	3 to 5 weeks	1:80 [†] (in early stages) = suspicious 1:160 [†] and rising = strongly suggestive
	Salmonella Flagellar d	later	later	1:40* = suspicious 1:160* = strongly suggestive
Paratyphoid Fever and other Salmonella Infections	Salmonella Flagellar a Salmonella O Group A Salmonella O Group B	Those characterized by prolonged fever and typhoid-like symptoms present antibodies of titers similar to above; lower titers may be more significant depending on the prevalence of a particular <i>Salmonella</i> species.		
Tularemia	Francisella tularensis	7 to 14 days	4 to 8 weeks	1:160 = strongly suggestive
Typhus Fevers	**Proteus OX19	7 to 10 days	by 14th day	1:40 to 1:80 early = suspicious 1:160 = strongly suggestive
Rocky Mountain Spotted Fever	**Proteus OX19	7 to 10 days	by 14th day	peak titers usually not above 1:160 to 1:320
Brucellosis	Brucella abortus	2 to 3 weeks	3 to 5 weeks	1:80 to 1:160 = strongly suggestive

[†]May be higher in vaccinated individuals

*Much higher in vaccinated individuals

**The Weil-Felix reactions that occur in rickettsial disease are shown in the chart below:

Disease	Proteus OX19	Proteus OX2	Proteus OXK
Epidemic typhus	4+	1+	0
Murine typhus	4+	1+	0
Scrub Typhus	0	0	3+
Spotted fever group	4+ or 1+	1+ or 3+	0
Q fever	0	0	0
Rickettsialpox	0	0	0

LIMITATIONS OF THE PROCEDURE

Prozone reactions are possible. Care should be taken to observe agglutination at the highest dilutions if no clumping is seen in lower dilutions.

Patients occasionally fail to develop any serum agglutinins.

In certain geographic regions and occupations, typhoid fever, *Salmonella* and *Brucella* are endemic and a high level of natural agglutinins may be present.

Cross reactions may occur between Brucella and Francisella antigens and antisera. Therefore, parallel tests should be run with those antigens. Generally a higher titer is obtained with the homologous antigen.

Cross reactions, previous vaccinations, anamnestic responses,⁸ antibiotic therapy, other diseases known or unknown, prozones, and autoagglutinins, as well as other factors, may affect these tests.

Consider the stage of disease when the specimen is collected. It is possible that peak titers occur only during convalescence. During infection, agglutinins to Salmonella O antigens usually appear earlier and disappear sooner than agglutinins to the H antigens.

Weil-Felix reactions may vary widely from case to case of spotted fever and therefore may be of little help in either detecting the disease or differentiating it from murine typhus.⁹

For more complete information, consult appropriate references.⁹⁻¹¹

These products are not intended as a substitute for culture. An appropriate attempt should be made to recover and identify the etiologic organism.

Only the product catalog numbers listed in this insert are to be used together. Do not use with other product catalog numbers.

PERFORMANCE CHARACTERISTICS

1. The positive control antisera should produce 2+ or greater agglutination at 1:80 in the slide or tube test when tested with the corresponding homologous febrile antigens.

Febrile Antigen	Homologous Control Antiserum
Salmonella O Group A	Salmonella Somatic Polyvalent
Salmonella O Group B	
Salmonella O Group D	
Salmonella O Flagellar a Salmonella O Flagellar b Salmonella O Flagellar d	Salmonella Flagellar Polyvalent
Brucella abortus Brucella melitensis Brucella suis	Brucella Positive Control
Proteus OX19 Proteus OX2 Proteus OXK	Proteus Polyvalent Control
Francisella tularensis	Francisella tularensis

- The negative control should show no agglutination with any of the Febrile Antigens.

Sensitivity of the Febrile Antigens is determined by demonstrating appropriate reactivity in both the slide and tube tests, as defined in the "Results" section, using the respective homologous Febrile Control Antisera. Specificity of the Febrile Antigens is determined by demonstrating non-reactivity against non-related (heterologous) Febrile Control Antisera.

AVAILABILITY

Cat. No.	Description
240786	Febrile Antigen Set, one set of 11 vials.

Febrile Antigens

240731	Salmonella O Group A Antigen (1-2-12), one 5 mL vial.
240732	Salmonella O Group B Antigen (1-4-5-12), one 5 mL vial.
240734	Salmonella O Group D Antigen (9-12) (Typhoid O), one 5 mL vial.
240834	Salmonella Flagellar a Antigen, one 5 mL vial.
240835	Salmonella Flagellar b Antigen, one 5 mL vial.
240785	Salmonella Flagellar d Antigen (Typhoid H), one 5 mL vial.
241049	Brucella abortus Antigen, one 5 mL vial.
240943	Brucella melitensis Antigen, one 5 mL vial.
240944	Brucella suis Antigen, one 5 mL vial.
240782	Proteus OX19 Antigen, one 5 mL vial.
240783	Proteus OX2 Antigen, one 5 mL vial.
240784	Proteus OXK Antigen, one 5 mL vial.
241050	Francisella tularensis Antigen, one 5 mL vial.

Febrile Control Antisera

240941	Salmonella Somatic Polyvalent Antiserum (A,B,D), one 5 mL vial.
240942	Salmonella Flagellar Polyvalent Antiserum (a,b,d), one 5 mL vial.
240934	Brucella Positive Control Antiserum (AMS), one 5 mL vial.
240940	Proteus Polyvalent Antiserum (OXK, OX2,OX19), one 5 mL vial.
240939	Francisella tularensis Antiserum, one 5 mL vial.
240937	Febrile Antigen Negative Control, one 5 mL vial.

REFERENCES

- Davidsohn, I, and J.B. Henry (ed.). 1969. Todd-Sanford clinical diagnosis by laboratory methods, 14th ed. W.B. Saunders Company, Philadelphia.
- Schubert, J.H., L. Holdeman, and D.S. Martin. 1954. J. Lab. Clin. Med. 44:194.
- Sonnenwirth, A.C. 1970. Gradwohl's clinical laboratory methods and diagnosis, 7th ed., p. 1482. The C.V. Mosby Co., St. Louis.
- National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
- U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
- Foome, J., and H.R. Morgan. 1954. Amer. J. Med. Soc. 228:520.
- Ormsbee, R.A. 1985. Rickettsiae, p. 845-855. In E.H. Lennette, A. Balows, W. J. Hausler, Jr., and H.J. Shadomy (ed.), Manual of clinical microbiology, 4th ed. American Society for Microbiology, Washington, D.C.
- Kolmer, J.A., E.H. Spaulding, and H.W. Robinson. 1951. Approved laboratory technic, 5th ed. Appleton Century-Crofts, New York.
- Sack, R.B. 1986. Serologic tests for the diagnosis of enterobacterial infections, p. 359-362. In N.R. Rose, H. Friedman, and J.L. Fahey (ed.), Manual of clinical laboratory immunology, 3rd ed. American Society for Microbiology, Washington, D.C.

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Benex Limited