
Revisions

SO 0191-5

Rev from	Rev to	ECO #
0408	0608	4823-08

Notes:

1. BD Cat. Number Various
2. Blank (Sheet) Size: Length: 25.5" Width: 16.5"
 Number of Pages: 18 Number of Sheets: 1
 Page Size: Length 8.5" Width 5.5" Final Folded Size: 4.25" x 5.5"
3. Style (see illustrations below): # 7

4. See Specification Control Number N/A for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# Standard Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	<p style="font-size: small;">COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION</p>	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8820191JAA		Category and Description Package Insert, Gram Stain Kits and Reagents	Sheet: 1 of 19 <hr/> Scale: N/A	A

BD Gram Stain Kits and Reagents

English: pages 1 – 3 Italiano: pagine 10 – 12
Français: pages 4 – 6 Español: páginas 13 – 15
Deutsch: Seiten 7 – 9

8820191JAA
2008/06

Pokyny vám poskytnú miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Інструкції зискате у місцевого зástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva.

Gram Stain Kit		1	Cat. No. 212539
Gram Crystal Violet	For the differential staining of bacteria.	1 x 250 mL	
Gram Iodine (Stabilized)		1 x 250 mL	
Gram Decolorizer		1 x 250 mL	
Gram Safranin		1 x 250 mL	
Gram Stain Kit		1	212524
Gram Crystal Violet	For the differential staining of bacteria.	1 x 250 mL	
Gram Iodine (Unstabilized)		1 x 250 mL	
Gram Decolorizer		1 x 250 mL	
Gram Safranin		1 x 250 mL	
Gram Crystal Violet	For staining microorganisms by the differential Gram method.	4 x 250 mL 1 x 3.8 L	212525 212526
Gram Iodine (Stabilized)	For staining microorganisms by the differential Gram method.	4 x 250 mL 1 x 3.8 L	212542 212543
Gram Iodine (Unstabilized)	For staining microorganisms by the differential Gram method.	4 x 250 mL 1 x 3.8 L	212529 212530
Gram Decolorizer	For staining microorganisms by the differential Gram method.	4 x 250 mL 1 x 3.8 L	212527 212528
Gram Safranin	For staining microorganisms by the differential Gram method.	4 x 250 mL 1 x 3.8 L	212531 212532
Gram Basic Fuchsin	For staining microorganisms by the differential Gram method.	4 x 250 mL 1 x 3.8 L	212544 212545

INTENDED USE

Gram Stain Kits and Reagents are used to stain microorganisms from cultures or specimens by the differential Gram method.

SUMMARY AND EXPLANATION

The Gram stain was devised in 1884 by Christian Gram in an attempt to differentiate bacterial cells from infected tissue. Although Gram observed what is now called the "Gram reaction," he did not recognize the taxonomic value of his technique.

The Gram stain is now used to differentiate intact, morphologically similar bacteria into two groups based on cell color after staining. In addition, cell form, size and structural details are evident. Such preliminary information provides important clues to the type of organism(s) present and the further techniques required to characterize them.

Because inorganic iodine is rapidly oxidized and loses its effectiveness as a mordant,¹ the Gram Stain Kit (Cat. No. 1212539) differs from Gram's original formulation by offering a more stable organic iodine complex, L-polyvinylpyrrolidone-iodine.

PRINCIPLES OF THE PROCEDURE

The Gram stain procedure² consists of:

Staining a fixed smear with crystal violet.

Applying iodine as a mordant.

Decolorizing the primary stain with alcohol/acetone; and, counterstaining with safranin or basic fuchsin.

A crystal violet-iodine complex forms in the protoplast (not the cell wall) of all organisms stained by this procedure. Organisms able to retain this dye complex after decolorization are classified as gram-positive while those that can be decolorized and counterstained are classified as gram-negative.

Upon disruption or removal of the cell wall, the protoplast of gram-positive (as well as gram-negative) cells can be decolorized and the gram-positive attribute lost. Thus, the mechanism of the Gram stain appears to be related to the presence of an intact cell wall able to act as a barrier to decolorization of the primary stain.

Generally, the cell wall is nonselectively permeable. It is theorized that during the Gram stain procedure, the cell wall of gram-positive cells is dehydrated by the alcohol in the decolorizer and loses permeability, thus retaining the primary

stain. However, the cell wall of the gram-negative cells has a higher lipid content and becomes more permeable when treated with alcohol, resulting in loss of the primary stain.

The molecular basis for the Gram stain has not yet been determined.

REAGENTS

Approximate Formula* Per Liter

Gram Crystal Violet

PRIMARY STAIN

Crystal Violet	3.0 g
Isopropanol	50.0 mL
Ethanol/Methanol	50.0 mL
Distilled Water	900.0 mL

Gram Iodine

MORDANT

(Working solution prepared from Gram Diluent and Gram Iodine 100X)

Iodine Crystals	3.3 g
Potassium Iodide	6.6 g
Distilled Water	1.0 L

Stabilized Gram Iodine

MORDANT

Polyvinylpyrrolidone-Iodine Complex	100.0 g
Potassium Iodide	19.0 g
Distilled Water	1.0 L

Gram Decolorizer

DECOLORIZER

Acetone	250.0 mL
Isopropanol	750.0 mL

Gram Safranin

COUNTERSTAIN

Safranin O Powder (pure dye)	4.0 g
Ethanol/Methanol	200.0 mL
Distilled Water	800.0 mL

Gram Basic Fuchsin

COUNTERSTAIN

Basic Fuchsin	0.08 g
Phenol	2.6 g
Isopropyl Alcohol	4.5 mL
Distilled Water	993.0 mL

*Adjusted and/or supplemented as required to meet performance criteria.

Warnings and Precautions: For *in vitro* Diagnostic Use.

Over time, a fine precipitate may develop in Gram Basic Fuchsin. Product performance will not be affected.

Gram Crystal Violet: Harmful by inhalation and if swallowed. Irritating to eyes, respiratory system and skin. Possible risk of irreversible effects. Keep container tightly closed. Keep away from sources of ignition - No Smoking. Do not breathe gas / fumes / vapor / spray. Wear suitable protective clothing and gloves.

Gram Iodine 100X: Harmful by inhalation, in contact with skin and if swallowed. Keep container tightly closed. Do not breathe gas / fumes / vapor / spray. Avoid contact with skin and eyes. Wear suitable protective clothing.

Stabilized Gram Iodine: Harmful in contact with skin. Causes burns. Irritating to eyes, respiratory system and skin. Keep container tightly closed. Do not breathe gas / fumes / vapor / spray. Avoid contact with skin and eyes. Wear suitable protective clothing.

Gram Decolorizer: Irritating to eyes, respiratory system and skin. Keep container tightly closed and in a well ventilated place. Keep away from sources of ignition - No smoking. Do not breathe gas / fumes / vapor / spray. Avoid contact with skin and eyes.

Gram Safranin: Harmful by inhalation and if swallowed. Irritating to eyes, respiratory system and skin. Keep container tightly closed. Keep away from sources of ignition - No Smoking. Wear suitable protective clothing and gloves.

Gram Basic Fuchsin: May cause cancer. May cause heritable genetic damage. Irritating to eyes. Keep container tightly closed. Do not breathe dust. Avoid contact with skin and eyes. Wear suitable protective clothing.

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Gram Iodine 100X: Immediately remove all contaminated clothing.

After contact with skin, wash immediately with plenty of water.

If inhaled, remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Seek medical advice.

If swallowed seek medical advice immediately and show the container or label.

Storage: On receipt, store at 15 – 30°C. The expiration date is for product in unopened bottles stored as directed. Do not open until ready to use.

Use the traditional Gram Iodine working solution within 3 months of preparation, not exceeding the Expiry of either component.

Product Deterioration: Unstabilized Gram Iodine when reconstituted may cause variability in the Gram stain when sufficient iodine is no longer available in solution. Protect the iodine solution from undue exposure to air, light and heat, to ensure that the solution is providing proper mordant activity.

SPECIMEN COLLECTION AND PREPARATION

Apply the test specimen to a clean glass slide in a manner that will yield a thin, uniform smear. Emulsify colonies from an 18 – 24 h culture in saline to obtain the proper density.

Allow the smear to air dry.

Fix the smear to the slide using one of the following techniques:

- Heat fix by passing the slide through a low flame 2 – 3 times. Cool the slide to room temperature before staining.
NOTE: Do not overheat the slide; excessive heating will cause atypical staining.
- Methanol fix the slide by flooding with absolute methanol for 1 – 2 min and rinse with tap water before staining.³
NOTE: For proper fixation, store absolute methanol in a brown screw-capped bottle and replenish the working supply every two weeks.

PROCEDURE

Reagent Preparation

Prepare the traditional Gram Iodine working solution by adding an entire 2.5 mL ampule of Gram Iodine 100X to 250 mL Gram Diluent or an entire 40 mL vial of Gram Iodine 100X to 3.8 L of Gram Diluent; mix thoroughly.

Materials Provided: Gram Crystal Violet, Gram Iodine or Stabilized Gram Iodine, Gram Decolorizer and Gram Safranin or Gram Basic Fuchsin.

Materials Required But Not Provided: Microscope slides, bunsen burner or methanol, bacteriological loop, swabs, blotting paper, microscope with oil immersion lens and Gram slide.

Test Procedure:

1. Flood the fixed smear with primary stain (Gram Crystal Violet) and stain for 1 min.
2. Remove the primary stain by gently washing with cold tap water.
3. Flood the slide with mordant (Gram Iodine or Stabilized Gram Iodine) and retain on the slide for 1 min.
4. Remove the mordant by gently washing with tap water.
5. Decolorize (Gram Decolorizer) until solvent running from the slide is colorless (3 – 60 sec).
6. Wash the slide gently in cold tap water.
7. Flood the slide with counterstain (either Gram Safranin or Gram Basic Fuchsin) and stain for 30 – 60 sec.
8. Wash the slide with cold tap water.
9. Blot with blotting paper or paper towel or allow to air dry.
10. Examine the smear under an oil immersion lens.

User Quality Control

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

Run controls using BBL™ Gram Slide (Cat. No. 231401) or 18 – 24 h cultures of known gram-positive and gram-negative microorganisms. The following test strains are recommended:

Organism	ATCC™	Expected Results
<i>Staphylococcus aureus</i>	25923	gram-positive cocci
<i>Escherichia coli</i>	25922	gram-negative rods

LIMITATIONS OF THE PROCEDURE

The Gram stain provides preliminary identification information only and is not a substitute for cultural studies of the specimen. Gram stain results must be confirmed with additional procedures such as direct antigen tests and culture on media.

Prior treatment with antibacterial drugs may cause gram-positive organisms from a specimen to appear gram-negative.

Use of an 18 – 24 h culture is advisable for best results since fresh cells have a greater affinity than old cells for most dyes. This is particularly true of many spore formers, which are strongly gram-positive when examined in fresh cultures but which later become gram-variable or gram-negative.

The Gram stain reaction is altered by physical disruption of the bacterial cell wall or protoplast. The cell walls of gram-positive bacteria interpose a barrier which prevents leaching of the dye complex from the cytoplasm. Cell walls of gram-negative bacteria contain lipids soluble in organic solvents, which are then free to decolorize the cytoplasm. Therefore, a microorganism that is physically disrupted by excess heating will not react to Gram staining as expected. "Careful adherence to procedure and interpretive criteria is required for accurate results. Accuracy is highly dependent on the training and skill of the microbiologist."²

Gram stain results, including organism morphology, can be affected by the age of the isolate, bacteria containing autolytic enzyme systems, cultures transferred from antibiotic-containing media, as well as specimens collected from patients on antibiotics.⁴ "Background material and artifacts can also interfere with interpretation. Precipitated gram-positive stain generally appears as irregular coccoid shapes or as asters resembling fungal hyphae."⁴

EXPECTED RESULTS AND PERFORMANCE CHARACTERISTICS¹⁻⁴

Reaction	Using Gram Safranin	Using Gram Basic Fuchsin
Gram-positive	Purple-black cells	Bright purple to purple-black cells
Gram-negative	Pink to red cells	Bright pink to fuchsia cells

REFERENCES

1. Magee, C.M., G. Rodenheaver, M.T. Edgerton, and R.F. Edlich. 1975. A more reliable Gram staining technique for diagnosis of surgical infections. *Am. J. Surg.* 130:341-346.
2. Kruczak-Filipov, P., and R.G. Shively. 1994. Gram stain procedures, p. 1.5-1-1.5.18. *In* H.D. Isenberg (ed.), *Clinical microbiology procedures handbook*, vol.1. American Society for Microbiology, Washington, D.C.
3. Mangels, J.I., M.E. Cox, and L.H. Lindley. 1984. Methanol fixation. An alternative to heat-fixation of smear. *Diag. Microbiol. Infect. Dis.* 2:129-137.
4. Chapin, K. 1995. Clinical microscopy, p. 33-51, *In* P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover and R.H. Tenover (ed.), *Manual of clinical microbiology*, 6th ed., American Society for Microbiology, Washington, D.C.

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ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo and BBL are trademarks of Becton, Dickinson and Company. ©2008 BD.
BENEX Limited