## Revisions


SO 0191-5

Rev from	Rev to	ECO #
2010/06	2011/03	5727-11


**Notes:**

1. BD Cat. Number Various
2. Blank (Sheet) Size : Length: 25.5"    Width: 22"  
 Number of Pages: 24    Number of Sheets: 1  
 Page Size: Length 8.5"    Width 5.5"    Final Folded Size: 4.25" x 5.5"
3. Style (see illustrations below): # 7


Accordion folds first, two horizontal parallel folds to form a size of 8 1/2" x 5 1/2", and then final horizontal fold to form a final size of 4 1/4" x 5 1/2" Showing English header.


#1


#2


#3


#4

#5

#6

#7

4. See Specification Control Number N/A for Material Information
5. Ink Colors: Printed two sides  Yes     No  
 No. of Colors: 1    PMS# Standard Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	<small>COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION</small>	<b>Becton, Dickinson and Company</b> 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8820201JAA		Category and Description Package Insert, TB Stain Kits and Reagents	Sheet: 1 of 25 <hr/> Scale: N/A	A

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TB Stain Kit K		1	<b>Cat. No. 212522</b>
TB Carbofuchsin KF	For staining mycobacteria by the Kinyoun (cold)	1 x 250 mL	
TB Decolorizer	acid-fast procedure.	1 x 250 mL	
TB Brilliant Green K		1 x 250 mL	
TB Stain Kit ZN		1	<b>212520</b>
TB Carbofuchsin ZN	For staining mycobacteria by the Ziehl-Neelsen	1 x 250 mL	
TB Decolorizer	(hot) acid-fast procedure.	1 x 250 mL	
TB Methylene Blue		1 x 250 mL	
TB Fluorescent Stain Kit M		1	<b>212519</b>
TB Auramine M	For staining mycobacteria by the Morse, Blair,	1 x 250 mL	
TB Decolorizer TM	Weiser and Sproat fluorescent procedure.	1 x 250 mL	
TB Potassium Permanganate		1 x 250 mL	
TB Fluorescent Stain Kit T		1	<b>212521</b>
TB Auramine-Rhodamine T	For staining mycobacteria by the Truant, Brett	1 x 250 mL	
TB Decolorizer TM	and Thomas fluorescent procedure.	1 x 250 mL	
TB Potassium Permanganate		1 x 250 mL	
TB Auramine M	For staining mycobacteria by the Morse, Blair,	4 x 250 mL	<b>212514</b>
	Weiser and Sproat fluorescent procedure.		
TB Auramine-Rhodamine T	For staining mycobacteria by the Truant, Brett	4 x 250 mL	<b>212515</b>
	and Thomas fluorescent procedure.		
TB Brilliant Green K	For staining mycobacteria by the Kinyoun (cold)	4 x 250 mL	<b>212523</b>
	and Ziehl-Neelsen (hot) acid-fast procedures.		
TB Carbofuchsin KF	For staining mycobacteria by the Kinyoun (cold)	4 x 250 mL	<b>212518</b>
	acid-fast procedure.		
TB Carbofuchsin ZN	For staining mycobacteria by the Ziehl-Neelsen	4 x 250 mL	<b>212511</b>
	(hot) acid-fast procedure.		
TB Decolorizer	For staining mycobacteria by the Kinyoun (cold)	4 x 250 mL	<b>212517</b>
	and Ziehl-Neelsen (hot) acid-fast procedures.		
TB Decolorizer TM	For staining mycobacteria by the Truant, Brett	4 x 250 mL	<b>212512</b>
	and Thomas and the Morse, Blair, Weiser and		
	Sproat fluorescent procedures.		
TB Methylene Blue	For staining mycobacteria by the Kinyoun (cold)	4 x 250 mL	<b>212516</b>
	and Ziehl-Neelsen (hot) acid-fast procedures.		
TB Potassium Permanganate	For staining mycobacteria by the Truant, Brett	4 x 250 mL	<b>212513</b>
	and Thomas and the Morse, Blair, Weiser and		
	Sproat fluorescent procedures.		

## INTENDED USE

TB Stain Kits are used to stain smears prepared from specimens and cultures suspected of containing mycobacteria for early presumptive diagnosis of mycobacterial infection and characterization of bacterial isolates.

## SUMMARY AND EXPLANATION

The microscopic staining technique is one of the earliest methods devised for detecting the tubercle bacillus and it remains a standard procedure.<sup>1</sup> The unique acid-fast characteristic of mycobacteria makes the staining technique valuable in early presumptive diagnosis as well as in providing information about the number of acid-fast bacilli present.

Fluorescent microscopy offers numerous advantages over classic methods for detecting mycobacteria because of its speed and simplicity, the ease of examining the slide and the reliability and superiority of the method.<sup>2</sup>

The TB Stain Kit K uses the Kinyoun (cold) acid-fast procedure.<sup>3</sup>

The TB Stain Kit ZN uses the Ziehl-Neelsen (hot) acid-fast procedure.<sup>4</sup>

The TB Fluorescent Stain Kit M uses the auramine O acid-fast fluorescent procedure.<sup>5</sup>

The TB Fluorescent Stain Kit T uses the auramine-rhodamine acid-fast fluorescent procedure described by Truant, Brett and Thomas.<sup>6</sup>

## PRINCIPLES OF THE PROCEDURE

A major characteristic of mycobacteria is that they are "acid-fast," i.e., once stained with the aniline dye, basic fuchsin, they are difficult to decolorize, and retain the red color even when treated with a mixture of acid and alcohol.

In the Ziehl-Neelsen stain, methylene blue is the counterstain. In the Kinyoun stain, brilliant green provides the necessary contrast.

Fluorescent dyes possess the property of emitting visible radiation when stimulated by shorter wavelength ultra violet light, i.e., a mercury vapor lamp or a strong, blue appropriately filtered light source.<sup>7</sup>

Auramine O and rhodamine B are reported to bind to mycolic acids in the mycobacterial cell wall.<sup>8</sup> Unbound dye is removed by the acid-alcohol decolorizer solution. Potassium permanganate serves as a counterstain to suppress nonspecific background fluorescence.

## REAGENTS

Approximate Formula\* Per Liter

### TB Stain Kit K

#### TB Carbofuchsin KF

Basic Fuchsin	15.0 g
Phenol, USP	45.0 g
Isopropanol	200.0 mL
Ethanol	50.0 mL
Distilled Water	750.0 mL

#### TB Decolorizer

Hydrochloric Acid	30.0 mL
Ethanol/Methanol	970.0 mL

#### TB Brilliant Green K

Brilliant Green	2.0 g
Sodium Hydroxide	0.02 g
Distilled Water	1000.0 mL

### TB Stain Kit ZN

#### TB Carbofuchsin ZN

Basic Fuchsin	1.7 g
Phenol, USP	50.0 g
Isopropanol	95.0 mL
Distilled Water	905.0 mL

#### TB Decolorizer

Hydrochloric Acid	30.0 mL
Ethanol/Methanol	970.0 mL

#### TB Methylene Blue

Methylene Blue, USP	2.4 g
Ethanol/Methanol	300.0 mL
Distilled Water	700.0 mL

### TB Fluorescent Stain Kit M

#### TB Auramine M

Auramine O	2.0 g
Phenol, USP	4.0 g
Glycerine, USP	100.0 mL
Isopropanol	250.0 mL
Distilled Water	650.0 mL

#### TB Decolorizer TM

Hydrochloric Acid	5.0 mL
Isopropanol	700.0 mL
Distilled Water	300.0 mL

#### TB Potassium Permanganate

Potassium Permanganate	5.0 g
Distilled Water	1000.0 mL

### TB Fluorescent Stain Kit T

#### TB Auramine-Rhodamine T

Auramine O	12.0 g
Rhodamine B	6.0 g
Phenol, USP	80.0 g
Glycerine, USP	600.0 mL
Isopropanol	140.0 mL
Distilled Water	260.0 mL

#### TB Decolorizer TM

Hydrochloric Acid	5.0 mL
Isopropanol	700.0 mL
Distilled Water	300.0 mL

#### TB Potassium Permanganate

Potassium Permanganate	5.0 g
Distilled Water	1000.0 mL

\*Adjusted and/or supplemented as required to meet performance criteria.

## Warnings and Precautions:

For *in vitro* Diagnostic Use.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>9-12</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Biosafety Level 2 practices and procedures, containment equipment and facilities are required for non-aerosol-producing manipulations of clinical specimens such as preparation of acid-fast smears. All aerosol-generating activities must be conducted in a Class I or II biological safety cabinet. Biosafety Level 3 practices, containment equipment and facilities are required for laboratory activities in the propagation and manipulation of cultures of *M. tuberculosis* and *M. bovis*. Animal studies also require special procedures.<sup>11</sup>

These reagents are harmful or fatal if swallowed and can cause eye irritation or blindness. In the event of eye contact, flush thoroughly with an eye wash system or tap water for at least 15 min and seek medical advice. The dye solutions and decolorizer solution are highly flammable and vapors may be harmful; use away from open flame. Avoid contact with skin, eye or mucous membranes. Use all reagents in a well ventilated area, avoiding prolonged inhalation, skin contact and ingestion.

**TB Carbofuchsin ZN:** Toxic in contact with skin and if swallowed. Causes burns. Keep container tightly closed. Avoid contact with skin and eyes. Wear suitable protective clothing.

**TB Decolorizer TM:** Causes burns. Keep container tightly closed and in a well ventilated place. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves and eye/face protection. In case of accident or if you feel unwell, seek medical advice immediately.

**TB Potassium Permanganate:** Irritating to eyes, respiratory system and skin. Keep container tightly closed. Avoid contact with skin and eyes. Wear suitable protective clothing.

**TB Auramine M:** Flammable. Irritating to eyes, respiratory system and skin. After contact with skin, wash immediately with plenty of water. Keep container tightly closed. Keep away from sources of ignition - No Smoking. Wear suitable protective clothing.

**TB Auramine-Rhodamine T:** Flammable. Toxic in contact with skin and if swallowed. Causes burns. Possible risks of irreversible effects. Keep container tightly closed. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves and eye/face protection. In case of accident or if you feel unwell, seek medical advice immediately.

**TB Decolorizer:** Harmful by inhalation and if swallowed. Keep container tightly closed and in a well ventilated place. Keep away from sources of ignition - No smoking. Wear suitable protective clothing.

**TB Methylene Blue:** Harmful by inhalation and if swallowed. Irritating to eyes, respiratory system and skin. Possible risks of irreversible effects. Keep container tightly closed. Keep away from sources of ignition - No Smoking. Wear suitable protective clothing and gloves.

**TB Carbofuchsin KF:** Harmful in contact with skin and if swallowed. Irritating to eyes and skin. Possible risks of irreversible effects. Keep container tightly closed. Keep away from sources of ignition - No Smoking. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing and gloves.

**Storage:** Store TB Stain Kits and Reagents at 15 - 30°C. Reagents that have been removed from the packing carton should be stored in the dark. The expiration date applies to the product in its intact container when stored as directed.

## SPECIMEN COLLECTION AND PREPARATION

Apply a thin smear of the specimen directly on a clean microscope slide. Sputum (direct specimen or concentrated sputum), spinal or other body fluid (catheterized urine), homogenized tissue and impression slides are satisfactory specimens for smears.<sup>1</sup>

Allow smear to dry.

Fix the smear to the slide by passing the slide through a low flame 2 - 3 times, avoiding excessive heat, or use an electric slide warmer at 65 to 75°C for at least 2 h.<sup>13</sup>

## PROCEDURE

**Materials Provided\*:** TB Stain Kit K, TB Stain Kit ZN, TB Quick Stain Kit, TB Fluorescent Stain Kit M, TB Fluorescent Stain Kit T.

*\*Individual reagents available separately.*

**Materials Required But Not Provided:** Microscope slides - new or cleaned in acid dichromate solution, staining rack, filter paper and microscope with oil immersion lens or fluorescent microscope.

## Test Procedure

### KINYOUN STAIN

TB Stain Kit K

#### Instructions:

1. Place slides on a staining rack and flood with TB Carbofuchsin KF for 4 min. Do not heat.
2. Wash gently in running water.
3. Decolorize with TB Decolorizer for 3 - 5 s.
4. Wash gently in running water.
5. Counterstain with either TB Brilliant Green K or TB Methylene Blue (available separately) for 30 s.
6. Wash gently in running water.
7. Air dry. If using TB Methylene Blue, dry over gentle heat.

## ZIEHL-NEELSEN STAIN

### TB Stain Kit ZN

#### Instructions:

1. Place slides on a staining rack and place a piece of filter paper, larger than the size of the smear, on each slide. Flood with TB Carbofuchsin ZN. Heat gently to steaming and allow to steam for 5 min.
2. Remove the filter paper. Wash gently in running water.
3. Decolorize with TB Decolorizer with two changes of reagent for 1 - 2 min until no more red color appears in washing.
4. Wash gently in running water.
5. Counterstain with either TB Methylene Blue or TB Brilliant Green K for 30 s.
6. Wash gently in running water.
7. Dry over gentle heat.

## MORSE STAIN

### TB Fluorescent Stain Kit M

#### Instructions:

1. Place slides on a staining rack and flood with TB Auramine M for 15 min.
2. Wash gently in running water.
3. Decolorize with TB Decolorizer TM for 30 - 60 s.
4. Wash slides gently in running water.
5. Counterstain with TB Potassium Permanganate for 2 min.
6. Wash gently in running water.
7. Air dry.
8. Examine under a microscope fitted with filter sets as described by Chapin and Murray.<sup>14</sup>

## TRUANT STAIN

### TB Fluorescent Stain Kit T

#### Instructions:

1. Place slides on a staining rack and flood with TB Auramine-Rhodamine T that has been thoroughly shaken prior to use. Leave undisturbed for 20 - 25 min at room temperature.
2. Wash gently in running water.
3. Decolorize with TB Decolorizer TM for 2 - 3 min.
4. Wash gently in running water.
5. Counterstain with TB Potassium Permanganate for 4 - 5 min.
6. Wash gently in running water.
7. Blot lightly. Dry in air or very gently over a flame.
8. Examine under a microscope fitted as described by Truant et. al.,<sup>6</sup> or Chapin and Murray.<sup>14</sup>

**User Quality Control:** It is recommended that controls be run each day of use using BBL™ AFB Slide (Cat. No. 231391) or known acid-fast and nonacid-fast microorganisms.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

Organism	Expected Results
<i>Mycobacterium tuberculosis</i> H37Ra ATCC™ 25177	Positive
<i>Bacillus subtilis</i> ATCC 6633	Negative

## LIMITATIONS OF THE PROCEDURE

A positive staining reaction provides presumptive evidence of the presence of mycobacteria in the specimen. A negative staining reaction does not necessarily indicate that the specimen will be culturally negative for mycobacteria. For positive identification of mycobacteria, other methods must be employed (cultural, molecular, etc.).

## EXPECTED RESULTS AND PERFORMANCE CHARACTERISTICS<sup>1-8,14</sup>

	Acid-fast Bacilli	Nonacid-fast Bacilli
TB Stain Kit K	dark red to pink	green
TB Stain Kit ZN	dark pink to red	blue
<b>Mycobacteria</b>		
TB Fluorescent Stain Kit M	bright yellow-green fluorescence	
TB Fluorescent Stain Kit T	reddish-orange fluorescence	

## REFERENCES

1. Nolte, F.S. and B. Metchock. 1995. *Mycobacterium*, p. 400 – 436. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Towner and R.H. Tenover (ed.), Manual of clinical microbiology, 6<sup>th</sup> ed. American Society for Microbiology, Washington, D.C.
2. Towner, F.C., J.T. Crawford, R.E. Huebner, L.J. Geiter, C.R. Horsburgh, Jr., and R.C. Good. 1993. The resurgence of tuberculosis: is your laboratory ready? *J. Clin. Microbiol.* 31: 767 – 770.
3. Ebersole, L.L. 1994. Acid-fast stain procedures, p.3.5.1. – 3.5.11. In H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
4. Murray, R.G.E., R.N. Doetsch and C.F. Robinow. 1994. Determinative and cytological light microscopy, p. 21 – 41. In P. Gerhardt, R.G.E. Murray, W.A. Wood and N.R. Krieg (ed.), Methods for general and molecular bacteriology, American Society for Microbiology, Washington, D.C.
5. Heifets, L.B. and R.C. Good. 1994. Current laboratory methods for the diagnosis of tuberculosis, p. 85 – 110. In B.R. Bloom (ed.), Tuberculosis: pathogenesis, protection, and control, American Society for Microbiology, Washington, D.C.
6. Truant, J.P., W.A. Brett and W. Thomas, Jr. 1962. Fluorescence microscopy of tubercle bacilli stained with auramine and rhodamine. *Henry Ford Hosp. Med. Bull.* 10: 287 – 296.
7. Cernoch, P.L., R.K. Enns, M.A. Saubolle and R.J. Wallace, Jr. 1994. Cumitech 16A, Laboratory diagnosis of the mycobacterioses. Coord. ed. A.S. Weissfeld. American Society for Microbiology, Washington, D.C.
8. Baron, E.J. and S.M. Finegold. 1990. Bailey & Scott's diagnostic microbiology, 8<sup>th</sup> ed. The C.V. Mosby Company, Baltimore, MD.
9. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
10. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17: 53-80.
11. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
12. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021-0045.
13. Kent, P.T. and G.P. Kubica, 1985. Public Health Mycobacteriology. A guide for the level III laboratory, p. 58, U.S. Dept. of Health and Human Services, Centers for Disease Control, Atlanta Ga.
14. Chapin, K.C. and P.R. Murray. 1999. Stains, p. 1674 – 1686. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Towner and R.H. Tenover (ed.), Manual of clinical microbiology, 7<sup>th</sup> ed. American Society for Microbiology, Washington, D.C.

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ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo and BBL are trademarks of Becton, Dickinson and Company © 2011 BD.
Benex Limited