Revisions

SO 0191-5

Rev from	Rev to	ECO #
0908	2010-01	5254-10

Notes:

1. BD Cat. Number 440461
2. Blank (Sheet) Size : Length: 8.5" Width: 27.5"
 Number of Pages: 10 Number of Sheets: 1
 Page Size: Length 8.5" Width 5.5" Final Folded Size: 8.5" x 5.5"
3. Style (see illustrations below): #4

4. See Specification Control Number N/A for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS # Standard Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	<small>COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION</small>	BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: L008319		Category and Description Package Insert BD ProbeTec CT/GC Urethral Dry Transport Kit	Sheet: 1 of 11 Scale:	A

 **BD ProbeTec™ ET *Chlamydia trachomatis* /
Neisseria gonorrhoeae (CT/GC) Amplified DNA
Assay Male Urethral Specimen Collection and
DRY TRANSPORT Kit**

 L008319
0086 2010/01

English:	pages	1 – 2	Italiano:	pagine	5 – 6
Dansk:	side	2 – 3	Español:	páginas	6 – 7
Français :	pages	3 – 4	Svenska:	sidan	7 – 8
Deutsch:	seiten	4 – 5			

Europe: Medical Device
93/42/EEC

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinizi danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндік жергілікті BD өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

The **BD ProbeTec™ ET *Chlamydia trachomatis*/Neisseria gonorrhoeae** (CT/GC) Amplified DNA Assay Male Urethral Specimen Collection and DRY TRANSPORT Kit offers a convenient way to collect and transport specimens from the patient to the laboratory. It is designed for the collection of male urethral specimens. This transport system is for use with the **BD ProbeTec ET *Chlamydia trachomatis* and *Neisseria gonorrhoeae*** amplified DNA assays on the **BD ProbeTec ET System**.

SUMMARY AND EXPLANATION

The **BD ProbeTec ET *Chlamydia trachomatis*/Neisseria gonorrhoeae** (CT/GC) Amplified DNA Assay Male Urethral Specimen Collection and DRY TRANSPORT Kit is a sterile, self-contained, disposable system. A flexible aluminum wire and small rayon-tipped swab facilitate obtaining the specimen and maintaining the organisms for nucleic acid amplification on the **BD ProbeTec ET System**. A plastic handle housing the aluminum wire and secured in the cap provides a slide action for extension of the swab if a longer swab is required.

REAGENTS

Materials Provided: The **BD ProbeTec ET *Chlamydia trachomatis*/Neisseria gonorrhoeae** (CT/GC) Amplified DNA Assay Male Urethral Specimen Collection and DRY TRANSPORT Kit contains 100 units. Each unit contains one rayon-tipped swab on a flexible and extendable aluminum wire in a tube.

Materials Required But Not Provided: Gloves.

Storage Instructions: Store swabs at 15 – 30°C. Do not use past the expiration date. Product sterility guaranteed if primary container is intact.

Precautions: For *in vitro* Diagnostic Use

Warnings:

1. Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁻⁴ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.
2. Optimal performance of the **BD ProbeTec ET** amplified DNA assays requires proper specimen collection, handling and transport.
3. Use only collection and transport systems labeled for use with the **BD ProbeTec ET *Chlamydia trachomatis* and *Neisseria gonorrhoeae*** Amplified DNA Assays.
4. Proper labeling should accompany each specimen to the laboratory.
5.  **BD ProbeTec ET *Chlamydia trachomatis*/Neisseria gonorrhoeae** (CT/GC) Amplified DNA Assay Male Urethral Specimen Collection and DRY TRANSPORT Kit is for single use only; reuse may cause a risk of infection and/or inaccurate results.

SPECIMEN COLLECTION AND TRANSPORT

Male Urethral Swab Specimen Collection

1. Insert the male urethral specimen collection swab 2 – 4 cm into the urethra and rotate for 3 – 5 s.
2. Withdraw the swab and place the cap/swab into the transport tube. Make sure the cap is tightly secured to the tube.
3. Label the tube with patient information and date/time collected.

Swab Storage and Transport

The male urethral specimen collection swab must be stored and transported to the laboratory and/or test site at 2 – 27°C within 4 – 6 days of collection. Storage up to 4 days has been validated with clinical specimens; storage

up to 6 days has been demonstrated with seeded specimens. In addition, storage up to 30 days at 2 – 8°C has been demonstrated with seeded specimens.

NOTE: If specimens cannot be transported directly to the testing laboratory under ambient temperatures (15 – 27°C) and must be shipped, an insulated container with ice should be used with either an overnight or 2-day delivery vendor.

For domestic and international shipments, specimens should be labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and etiologic agents/infectious substances. Time and temperature conditions for storage must be maintained during transport.

REFERENCES

1. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, PA.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
3. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.

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Becton Dickinson France S.A.S.