## Revisions

SO 0191-5

Rev from	Rev to	ECO #
2009/08	2010/04	5232-09

**Notes:**

1. BD Cat. Number 256030
2. Blank (Sheet) Size: Length: N/A      Width: N/A  
 Number of Pages: N/A      Number of Sheets: N/A  
 Page Size: Length N/A      Width N/A      Final Folded Size: N/A
3. Style (see illustrations below): N/A

4. See Specification Control Number BALTL009004 for Material Information
5. Ink Colors: Printed two sides  Yes       No  
 No. of Colors: N/A
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	<p style="font-size: small;">COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION</p> <p style="font-weight: bold; font-size: large;">BD</p> <p>Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA</p>	<p style="font-size: x-large; font-weight: bold;">A</p>	
Proofer	Date			
Checked By	Date		<p>Category and Description</p> <p style="font-weight: bold;">Package Insert, Directigen EZ RSV</p>	
Part Number: L009004JAA		<p>Sheet: 1 of 65</p> <hr/> <p>Scale:      N/A</p>		

For the direct detection of  
Respiratory Syncytial Virus (RSV)

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# BD Directigen™ EZ RSV

For the direct detection of Respiratory Syncytial Virus (RSV)

L009004JAA  
2010/04

## INTENDED USE

The **Directigen™ EZ RSV** test is a rapid chromatographic immunoassay for the direct and qualitative detection of Respiratory Syncytial Virus (RSV) antigen in nasopharyngeal washes, nasopharyngeal aspirates, nasopharyngeal swabs and nasopharyngeal swab/washes from patients suspected of having a viral respiratory infection. This test is intended for *in vitro* diagnostic use to aid in the diagnosis of Respiratory Syncytial Virus (RSV) infections in neonatal and pediatric patients under the age of 20. It is recommended that negative test results be confirmed by cell culture.

## SUMMARY AND EXPLANATION

Viral respiratory tract infections are responsible for widespread disease. Respiratory syncytial virus (RSV) is a leading cause of lower respiratory tract infections in young children during the cold season.<sup>1</sup> RSV has also been implicated in severe respiratory infections in the elderly and immunocompromised.<sup>2,3</sup> More recently, RSV has been identified as causing 20% of “influenza-like” illness in people 15–44 years of age.<sup>4</sup>

Traditional methods for detection of respiratory viruses have included cell culture and direct fluorescent antibody (DFA).<sup>5-7</sup> More recently, PCR has shown clinical utility for the detection of respiratory viruses including RSV.<sup>8</sup> Enzyme immunoassay (EIA) and rapid manual systems are available for specific viruses such as Influenza A/B and RSV. Rapid tests allow for a quick diagnosis so that patients may be appropriately isolated and treated to prevent the nosocomial spread of infections to fellow patients with compromised cardiac, respiratory or immune functions.<sup>9</sup> In addition, rapid tests assist with the selection of appropriate antiviral therapy.

The most common specimens collected for RSV testing include nasopharyngeal washes, nasopharyngeal aspirates and nasopharyngeal swabs. Nasopharyngeal washes and aspirates have been shown to be superior to swabs and are the specimens of choice.<sup>10</sup>

The **Directigen EZ RSV** antigen detection test is a chromatographic assay to detect RSV antigens extracted from various specimens of symptomatic patients. The speed and workflow of the **Directigen EZ RSV** test make it applicable as a “STAT” RSV antigen detection test, providing rapid, relevant information to assist with antiviral intervention and other clinical or support decisions.

## PRINCIPLES OF THE PROCEDURE

The **Directigen EZ RSV** test is a chromatographic assay to qualitatively detect RSV antigen in samples extracted from respiratory specimens. When extracted specimens are added to the test device, RSV A and/or B antigens bind to the antibody-colloidal gold conjugate in the test strip forming an antigen-antibody complex. This complex migrates across the test strip to the reaction area and is captured by the line of RSV antibody on the membrane. Excess conjugate binds to a second line consisting of inactivated RSV antigen that serves as a functional control. A positive result is indicated by the appearance of two reddish-purple lines in the read window, one line next to the Test “T” and the other next to the Control “C”. The absence of a reddish-purple line next to the “T” and the presence of a reddish-purple line next to the “C” indicate a negative result. The test is considered uninterpretable if no visible reddish-purple line is present next to the “C”.

## REAGENTS

The following are included in the **Directigen EZ RSV** kit.

<b>BD RSV Devices</b>	30	Each device containing a control line of inactivated RSV antigen and a test line of monoclonal RSV antibody on the membrane.
Extraction Reagent E	4.0 mL	Detergent with 0.2% sodium azide (preservative).
Control + Swab	1 ea.	Positive Control, detergent-treated RSV infected HEP-2 cells with < 0.1% sodium azide (preservative).
Control – Swab	1 ea.	Negative Control with < 0.1% sodium azide (preservative).
<b>DispensTube™ Tubes</b>	30	Tubes for specimen extraction and sample delivery into device.
<b>DispensTube Tips</b>	30	Tips to filter sample when delivered into device.

**Materials Required But Not Provided:** Pipette (capable of delivering 250 µL), timer, vortex and Transport Media (see Specimen Collection and Handling).

**Warnings and Precautions:**

For *in vitro* Diagnostic Use.

1. The RSV Control (+) swab and the control line on the **BD** RSV device have been prepared from RSV-infected tissue culture cells which have been inactivated by detergent treatment and sonication then subsequently tested by bioassay procedures.
2. Do not use the kit if the Controls do not yield appropriate results.
3. Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>11-14</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.
4. Do not use kit components beyond the expiration date.
5. Do NOT mix reagents from different kit lot numbers. Do not reuse the device.
6. Reagents contain sodium azide, which is harmful by inhalation, in contact with skin and if swallowed. Contact with acids liberates very toxic gas. If there is contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

**Storage and Handling:** Kits may be stored at 2–30°C. DO NOT FREEZE. Reagents and **BD** RSV devices must be at room temperature (15–30°C) when used for testing.

**SPECIMEN COLLECTION AND PREPARATION**

**Specimen Transport and Storage:** Transport fresh specimens to the laboratory as rapidly as possible in a suitable liquid transport system. Process specimens as soon as possible after collection. If necessary, specimens may be stored at 2–8°C for up to 72 h or at –20°C for up to seven days after collection.

It is essential that the correct specimen collection and preparation methods be followed. Do not centrifuge specimens prior to use with the **Directigen** EZ RSV Test, as the removal of cellular material will adversely affect the sensitivity of the test.

**Transport Media:** The following transport media have been tested and found to be compatible with the **Directigen** EZ RSV Test:

Normal Saline	<b>Trypticase</b> <sup>™</sup> Soy Broth + 0.5% gelatin*
Phosphate Buffered Saline (PBS)	<b>Trypticase</b> Soy Broth + 0.5% BSA
PBS + 0.5% gelatin*	Earle's Minimal Essential Medium (EMEM)
PBS + 0.5% Bovine Serum Albumin (BSA)	EMEM + 0.5% BSA
Veal Infusion Broth (VIB)	EMEM + 1% BSA
VIB + 0.5% BSA	EMEM + 0.5% lactalbumin hydrolysate
Hank's Balanced Salt Solution	EMEM + 1.0% lactalbumin hydrolysate
M4 Media*	Modified Stuart's (liquid) <b>CultureSwab</b> <sup>™™</sup>
M4-RT Media*	Amies (liquid) <b>CultureSwab</b>
M5 Media	Starplex Multitrans <sup>™™</sup> *
Bartels ViraTrans <sup>™</sup> Media	<b>BD</b> <sup>™</sup> Universal Viral Transport Medium*
Sucrose Phosphate (2-SP)	<b>BD</b> <sup>™</sup> ESwab

\*Medium contains gelatin.

**Specimen Collection and Preparation:** Nasopharyngeal washes and aspirates have been shown to be superior to nasopharyngeal swabs and are the specimens of choice.<sup>9,10</sup> Acceptable specimens for testing with the **Directigen** EZ RSV test include nasopharyngeal washes, nasopharyngeal aspirates, nasopharyngeal swabs and nasopharyngeal swab/washes.

**NOTE:** For nasopharyngeal swabs (NPS), polyester or rayon-tipped swabs with an aluminum wire are recommended. Flocked swabs have been analytically tested and have been found to be compatible with the **Directigen** EZ RSV test.<sup>15</sup> Calcium alginate-tipped swabs are not acceptable for collection of viral specimens.<sup>16</sup>

**Procedure for Nasopharyngeal Washes and Nasopharyngeal Swab/Washes:**

1. Sample volumes of 2–3 mL are recommended.
2. Excessive wash or lavage volumes should be avoided as they may result in decreased test sensitivity.
3. Process specimen as described in "Test Procedure."

**Procedure for Nasopharyngeal Swabs:**

1. Swab specimens may be added to 350 µL–2 mL of transport medium.

2. Vortex the swab and transport medium.
3. Remove as much liquid from the swab as possible.  
**NOTE:** Alternatively swab may be left in tube.
4. Dispose the swab into appropriate container.
5. Process specimen as described in "Test Procedure."

**Procedure for Nasopharyngeal Aspirates:**

1. Specimens may be dispersed in 1–3 mL of transport medium or saline prior to processing.
2. Process specimen as described in "Test Procedure."

**PROCEDURE**

**Test Procedure**

**NOTES:**

- Reagents, specimens and **BD** RSV devices must be at room temperature (15–30°C) when used for testing.
- Thoroughly mix all specimens prior to removal of an aliquot for extraction. **Do not centrifuge specimens** prior to use with the **Directigen EZ** RSV test, as the removal of cellular material will adversely affect the sensitivity of the test.
- To assure proper delivery, **DispensTube** tubes and reagent bottles must be held vertically (approximately one inch from the **BD** RSV device sample well or **DispensTube**), while gently dispensing one drop at a time, in quick succession.


1. Remove a **BD** RSV device from its foil pouch immediately before use.
2. Label a **BD** RSV device and a **DispensTube** tube for the control or specimen to be tested.
3. Place the labeled **DispensTube** in the designated area of the workstation or rack.
4. Gently mix Extraction Reagent E by inversion. Dispense 3 drops into the **DispensTube**. Hold reagent bottle vertically (approximately one inch from the **DispensTube**) while dispensing drops.

5. Process specimens and controls as directed below:
  - a. For specimens:
    1. Vortex or thoroughly mix specimen. **Do not centrifuge.**
    2. Pipette 250 µL of specimen into **DispensTube** (containing Reagent E).

- b. For Controls:
  1. Add 250 µL of normal saline, to the **DispensTube** tube (containing Reagent E).
  2. Insert control swab and express by rotating swab 6-8 times while pinching tube.
  3. Remove control swab while pinching tube to remove excess fluid from the swab tip.


6. Insert a **DispensTube** Tip into each **DispensTube**.  
**NOTE: Do not use tips from other Directigen products.**
7. Vortex or mix thoroughly.

8. Invert the **DispensTube** and holding the tube on the upper half, away from the tip, gently squeeze three (3) drops of the extracted specimen into the appropriately labeled **BD RSV** sample well.

To assure proper delivery, **DispensTube** tubes must be held vertically (approximately one inch from the **BD RSV** device sample well), while gently dispensing one drop at a time, in quick succession.

**NOTE: Squeezing the tube close to the tip may result in ejection of the tip and leakage of the contents from the tube.**

9. Read results at 15 min or up to 60 min. Positive results may be reported as soon as 5 min as long as the test and control line are visible.

**NOTE: If using a gelatin-containing transport medium, results must be read at 15 min. See "Limitations of the Procedure."**

10. Read in a well-lighted area and record the test result.

#### Quality Control:

Quality control requirements must be performed in accordance with local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to CLSI (formerly NCCLS) EP12-A<sup>17</sup> and 42 CFR 493.1202 (c)<sup>18</sup> for guidance on appropriate Quality Control practices.

Each **BD RSV** device contains built-in controls. Built-in control features include:

**Internal Control:** The appearance of a reddish-purple control line provides an internal antigen control (positive internal control) that serves as a reagent check for the conjugate and capture antibody as well as ensuring sufficient capillary flow has occurred. The absence of this line indicates an uninterpretable test.

**Test Membrane Control (negative internal control):** The unreactive membrane surface surrounding the positive internal control and test lines contrasts with a positive reaction and therefore serves as a background reference for interpreting reaction color.

#### External Positive and Negative Controls:

Positive and Negative controls are also supplied with each kit. These controls are provided as a means of additional quality control to demonstrate a positive or negative reaction. The external controls should be run as a quality control procedure for each new lot or shipment received. The formation of the reddish-purple line in the read window next to the "T" when the Control (+) is employed further indicates that the RSV antigen binding property of the membrane is functional. Do not report results if either the Control (+) or the Control (-) do not give appropriate results.

If expected control results are not obtained, do not report patient results. Contact your local BD representative or Technical Services for assistance.

#### INTERPRETATION OF RESULTS

**Positive Test (antigen present)** – A visible reddish-purple line appears in the read window next to the Test "T", and a reddish-purple line next to the Control "C". This indicates RSV antigen was detectable in the specimen. The background area should be white to light pink.

+

**Negative Test (no antigen detected)** – No reddish-purple line is visible next to the Test "T". This indicates RSV antigen was not detectable in the specimen. A reddish-purple line next to the Control "C" indicates proper performance of test procedure and reagents. The background area should be white to light pink.

-

uninterpretable

**Uninterpretable Test** – If a reddish-purple line next to the Control "C" is not visible or the background color interferes with interpretation of the test or control lines, the test is uninterpretable. An uninterpretable test must be repeated, a new specimen obtained and retested, or the specimen sent to the clinical laboratory for culture isolation.

## REPORTING OF RESULTS

- Positive** Positive for the presence of RSV antigen. A positive result may occur in the absence of viable virus.
- Negative** Negative for the presence of RSV antigen. Infection due to RSV cannot be ruled-out since the antigen present in the sample may be below the detection limit of the test. Culture confirmation of negative samples is recommended.

## LIMITATIONS OF THE PROCEDURE

1. The **Directigen EZ RSV** test is capable of detecting both viable and non-viable RSV particles. The **Directigen EZ RSV** test performance depends on antigen load and may not correlate with tissue culture performed on the same specimen. The etiology of respiratory infection caused by microorganisms other than Respiratory Syncytial Virus cannot be established with this test.
2. Inadequate specimen collection, improper sample handling/transport or low levels of virus shedding may yield a false-negative result. Accordingly, a negative test result does not eliminate the possibility of an RSV infection. As with all diagnostic procedures, the results obtained with the **Directigen EZ RSV** test should be used in conjunction with other clinical information available to the physician.
3. In chromatographic immunoassays, weakly visible lines are more likely to be false positives than are strongly visible lines. As with all diagnostic procedures, the results obtained with the **Directigen EZ RSV** test should be used in conjunction with other clinical information available to the physician.
4. The validity of the **Directigen EZ RSV** test has not been proven for identification/confirmation of tissue culture isolates and should not be used in this capacity.
5. Contents of the transport medium being used should be carefully reviewed to determine whether or not it contains gelatin. If using a gelatin-containing transport medium, results must be read at 15 min. Extending read time beyond 15 min may result in the appearance of faint lines in the test (T) position of the read window.
6. Amies Gel **CultureSwab™** Plus medium is not compatible with the **Directigen EZ RSV** test. Storing specimens in this medium over time may cause a potential false negative result.
7. Cross-reactivity of this assay with human metapneumovirus has not been studied.
8. Monoclonal antibodies may not detect all antigenic variants or new strains of RSV.
9. The **Directigen EZ RSV** test has not been evaluated with a sufficient number of lower respiratory specimens, lower nasal swabs, nose/throat swabs and pharyngeal swabs to estimate performance characteristics with these specimen types.

## EXPECTED VALUES

The rate of positivity observed in RSV testing will vary depending on the method of specimen collection, handling/transport system employed, detection method utilized, time of year, age of the patient, geographic location and most importantly, local disease prevalence. The prevalence observed with culture during the clinical study ranged from 4.2–46.5%.

## PERFORMANCE CHARACTERISTICS

Performance characteristics for the **Directigen EZ RSV** test were established in a multi-center study conducted at seven clinical centers during the 2001–2002 RSV season. The clinical centers were located in Canada and geographically diverse areas within the United States.

All specimens were collected and transported to the laboratory according to each laboratory's procedure. Each specimen was evaluated with the **Directigen EZ RSV** test, cell culture and the laboratory's routine method (i.e., DFA, EIA). Any remaining specimen was archived at or below -20°C.

For cell culture, a portion of the specimen was inoculated into cell lines appropriate for respiratory viral culture according to each laboratory's procedure. Cells were examined for the appearance of cytopathic effects (CPE). Infected cells were confirmed for RSV by immunofluorescent staining. Specimens negative for CPE at fourteen days were stained for negative confirmation by species-specific immunofluorescent staining.

A total of 1176 specimens, consisting of nasopharyngeal wash, nasopharyngeal aspirate, nasopharyngeal swab and nasopharyngeal swab/wash, from patients suspected of having RSV were evaluated with the **Directigen EZ RSV** test.

**Clinical Performance:** For all specimens evaluated, the overall sensitivity and specificity of the **Directigen EZ RSV** test for RSV when compared to culture were 80% and 91%, respectively. The overall uninterpretable rate was 0%. The performance characteristics of the **Directigen EZ RSV**

as compared to cell culture for each specimen type are shown in Table 1. The performance characteristics as compared to cell culture by clinical site are shown in Table 2.

**Table 1: Summary of the Performance of the Directigen EZ RSV (EZ) Test Compared to Culture for all Specimen Types**

Specimen Type	n	Culture/EZ				Sensitivity (%) (95% CI)	Specificity (%) (95% CI)
		+/+	-/+ <sup>a</sup>	+/-	-/-		
Nasopharyngeal Wash (NPW)	348	129	29	19	171	87.2 (80.7–92.1)	85.5 (79.8–90.1)
Nasopharyngeal Aspirate (NPA)	401	90	24	27	260	76.9 (68.2–84.2)	91.6 (87.7–94.5)
Nasopharyngeal Swab/Wash (NPSW)	160	28	6	11	115	71.8 (55.1–85.0)	95.0 (89.7–98.2)
Nasopharyngeal Swab (NPS)	267	20	20	10	217	66.7 (47.2–82.7)	91.6 (87.3–94.8)

(+) = RSV positive (-) = RSV negative

<sup>a</sup> The PCR method used to test the culture negative, **Directigen EZ RSV** positive specimens is not an FDA-cleared method.

There were 79 specimens that were culture negative, **Directigen EZ RSV** test positive. PCR testing was performed on 73 of the 79 specimens; a total of 46 of the 73 specimens were positive by PCR.

**Table 2: Summary of the Performance of the Directigen EZ RSV (EZ) Test Compared to Culture by Clinical Site**

Site	Specimen Type	n	Culture/EZ				Sensitivity (%) (95% CI)	Specificity (%) (95% CI)
			+/+	-/+	+/-	-/-		
Site 1	50% NPW 45% NPA 5% NPS	38	11	2	0	25	100 (71.5–100)	92.6 (83.8–99.4)
Site 2	100% NPW	173	56	18	15	84	78.9 (67.6–87.7)	82.4 (73.6–89.2)
Site 3	59% NPA 37% NPW 4% NPS	195	79	24	13	79	85.9 (77.1–92.3)	76.7 (67.3–84.5)
Site 4	100% NPS	207	4	17	5	181	44.4 (13.7–78.8)	91.4 (86.6–94.9)
Site 5	73% NPSW 24% NPS 2% NPW 1% NPA	219	46	8	15	150	75.4 (62.7–85.5)	94.9 (90.3–97.8)
Site 6	100% NPW	77	31	0	3	43	91.2 (76.3–98.1)	100 (91.8–100)
Site 7	100% NPA	267	40	10	16	201	71.4 (57.8–82.7)	95.3 (91.5–97.7)

(+) = RSV positive (-) = RSV negative

**Directigen EZ RSV** positive test results were distributed over a range of weakly visible to strongly visible lines. Of the 70 **Directigen EZ RSV** positive specimens with weakly visible lines, 39 were culture positive and a total of 58 were positive by one or more methods (i.e., culture, DFA, EIA, or PCR<sup>a</sup>). Of the 170 **Directigen EZ RSV** positive specimens with visible lines, 130 were culture positive while a total of 166 were positive by one or more methods (i.e., culture, DFA, EIA, or PCR<sup>a</sup>). Of the 106 **Directigen EZ RSV** positive specimens with strongly visible lines, 98 were culture positive while a total of 106 were positive by one or more methods (i.e., culture, DFA, EIA, or PCR<sup>a</sup>).

<sup>a</sup> The PCR method used to test the culture negative, **Directigen EZ RSV** positive specimens is not an FDA-cleared method.

**Reproducibility:** The reproducibility of the **Directigen EZ RSV** test was evaluated at three sites. The reproducibility panel was composed of 12 simulated RSV samples and included four replicates of a low positive level (near the Limit of Detection [LOD]) and four replicates of a moderate positive level of RSV antigen seeded into the negative control diluent. The panel also included four simulated RSV negative samples that were the unseeded negative control diluent. The overall reproducibility for the **Directigen RSV EZ** test was 99.1%.

**Analytical Studies**

**Analytical Sensitivity (Limit of Detection)**

The limit of detection (LOD) for the **Directigen EZ RSV** test was established for a total of five RSV strains; two RSV A and three RSV B strains.

Type	RSV Viral Strain	LOD (TCID <sub>50</sub> )*
A	RSV (Long)	3.95 x 10 <sup>3</sup>
A	RSV (A-2)	7.9 x 10 <sup>3</sup>
B	RSV (9320)	4.05 x 10 <sup>2</sup>
B	RSV (Washington)	7.03 x 10 <sup>3</sup>
B	RSV (Wild-type)	5.56 x 10 <sup>2</sup>

\* TCID<sub>50</sub> = Tissue Culture Infectious Dose at which 50% of cells are lysed.

**Analytical Specificity**

The **Directigen EZ RSV** test was evaluated using a total of 99 microorganisms (58 bacteria, two yeasts and 39 viruses). Bacteria and yeasts were tested at concentrations of 10<sup>8</sup> CFU/mL. *Mycoplasma orale* was tested at >10<sup>7</sup> CFU/mL. *Mycoplasma pneumoniae* was tested at >10<sup>6</sup> CCU/mL (CCU – color changing units). *Mycobacterium tuberculosis* was tested at 7.0 x 10<sup>6</sup> CFU/mL. Viral isolates were tested at titers between 10<sup>3.5</sup> to 10<sup>9.5</sup> TCID<sub>50</sub>/0.20 mL. Influenza viruses were tested at titers of 10<sup>3.5</sup> to 10<sup>9.5</sup> CEID<sub>50</sub>\*\*/0.20 mL. None of the microorganisms listed gave a positive result in the **Directigen EZ RSV** test.

\*\* CEID<sub>50</sub> = Chick Embryo Infectious Dose at 50% of chick embryo infected.

Bacteria and Yeast Panel		
<i>Acinetobacter baumannii</i> (calcoaceticus)	<i>Klebsiella pneumoniae</i>	<i>Proteus mirabilis</i>
<i>Actinobacillus suis</i>	<i>Lactobacillus casei</i>	<i>Proteus vulgaris</i>
<i>Bacteriodes fragilis</i>	<i>Lactobacillus fermentum</i>	<i>Pseudomonas aeruginosa</i>
<i>Bordetella pertussis</i>	<i>Lactobacillus plantarum</i>	<i>Salmonella choleraesuis</i> subsp <i>minnesota</i>
<i>Candida albicans</i>	<i>Legionella pneumophila</i>	<i>Serratia marcescens</i>
<i>Candida glabrata</i>	<i>Listeria monocytogenes</i>	<i>Staphylococcus aureus</i>
<i>Cardiobacterium hominis</i>	<i>Moraxella catarrhalis</i>	<i>Staphylococcus aureus</i> - Cowan
<i>Chlamydia psittaci</i>	<i>Mycobacterium avium</i>	<i>Staphylococcus epidermidis</i>
<i>Chlamydia trachomatis</i> LGVII	<i>Mycobacterium intracellulare</i>	<i>Streptococcus bovis</i> II Group D
<i>Corynebacterium diphtheriae</i>	<i>Mycobacterium tuberculosis</i>	<i>Streptococcus mutans</i>
<i>Eikenella corrodens</i>	<i>Mycoplasma orale</i>	<i>Streptococcus oralis</i>
<i>Enterococcus faecalis</i>	<i>Mycoplasma pneumoniae</i>	<i>Streptococcus pneumoniae</i>
<i>Enterococcus gallinarum</i>	<i>Neisseria gonorrhoeae</i>	<i>Streptococcus pyogenes</i> Group A
<i>Escherichia coli</i>	<i>Neisseria meningitidis</i>	<i>Streptococcus sanguis</i>
<i>Fusobacterium nucleatum</i>	<i>Neisseria mucosa</i>	<i>Streptococcus sp.</i> Group B
<i>Gardnerella vaginalis</i>	<i>Neisseria sicca</i>	<i>Streptococcus sp.</i> Group C
<i>Haemophilus aphrophilus</i>	<i>Neisseria subflava</i>	<i>Streptococcus sp.</i> Group F
<i>Haemophilus influenzae</i>	<i>Peptostreptococcus anaerobius</i>	<i>Streptococcus sp.</i> Group G
<i>Haemophilus parainfluenzae</i>	<i>Porphyromonas asaccharolyticus</i>	<i>Veillonella parvula</i>
<i>Haemophilus paraphrophilus</i>	<i>Prevotella oralis</i>	
<i>Kingella kingae</i>		

---

### Viral Panel

---

Adenovirus, Type 3	Influenza B/GL/1739/54
Adenovirus, Type 5	Influenza B/Hong Kong/5/72
Adenovirus, Type 7	Influenza B/Lee/40
Adenovirus, Type 10	Influenza B/Allen/45
Adenovirus, Type 18	Influenza B/Maryland/1/59
Coronavirus	Influenza B/Taiwan/2/62
Coxsackievirus Type A9 (Griggs)	Influenza C/Taylor/1233/47
Coxsackievirus Type A9 (P.B. Bozek)	Measles virus
Coxsackievirus Type B5	Mumps virus
Coxsackievirus Type B6	Parainfluenza Type 1
Coxsackievirus Type A21	Parainfluenza Type 2
Cytomegalovirus	Parainfluenza Type 3
Echovirus Type 2	Rhinovirus Type 1A
Echovirus Type 3	Rhinovirus Type 2
Echovirus Type 6	Rhinovirus Type 13
Echovirus Type 11	Rhinovirus Type 15
Herpes Simplex Virus Type 1	Rhinovirus Type 16
Herpes Simplex Virus Type 2	Rhinovirus Type 37
Influenza A (H1N1)	Varicella-Zoster Virus
Influenza A (H3N2)	

---

#### Interfering Substances

Various substances were tested in the **Directigen** EZ RSV test including whole blood at 2%, over the counter (OTC) drugs and prescription drugs. No interference was noted in the assay for any substance at the levels tested.

#### Non-Prescription Medications (OTC)

Three OTC mouthwashes at 25%	Diphenhydramine HCl at 5 mg/mL
Three OTC throat drops at 25%	Pseudoephedrine HCl at 20 mg/mL
Three OTC nasal sprays at 10%	Guaiaicol Glyceryl Ether at 20 mg/mL
4-Acetamidophenol at 10 mg/mL	Ibuprofen at 10 mg/mL
Acetylsalicylic acid at 20 mg/mL	Oxymetazoline at 0.05 mg/mL
Chlorpheniramine maleate at 5 mg/mL	Phenylephrine at 1 mg/mL
Dextromethorphan at 10 mg/mL	

#### Prescription Medications

Zanamavir at 1 mg/mL	Ribavirin at 500 ng/mL
Amantadine at 500 ng/mL	Synagis at 0.1 mg/mL
Rimantadine at 500 ng/mL	Albuterol Sulfate at 0.083 mg/mL

#### AVAILABILITY

Cat. No.	Description
256030	<b>Directigen™</b> EZ RSV 30 Test Kit
220115	<b>BBL™ CultureSwab™</b> Sterile Single Swab, Pkg. of 100
220131	<b>BBL™ CultureSwab™</b> Liquid Amies, Flexible Aluminum Wire, Pkg. of 50
220134	<b>BBL™ CultureSwab™</b> Liquid Stuart, Flexible Aluminum Wire, Pkg. of 50
220220	<b>BD™</b> Universal Viral Transport 3 mL Vial, Ctn. of 50
220244	<b>BD™</b> Universal Viral Transport 1 mL Vial, Ctn. of 50
220526	<b>BD™</b> Universal Viral Transport Kit, 1 mL Vial and a sterile nylon flocced flexible minitip swab, Ctn. of 50
220531	<b>BD™</b> Universal Viral Transport Kit, 3 mL Vial and a sterile nylon flocced flexible minitip swab, Ctn. of 50
220532	<b>BD™</b> ESwab Collection Kit, 1 mL Liquid Amies and a sterile nylon flocced flexible minitip swab, Ctn. of 50
256034	<b>Directigen™</b> EZ RSV Control Kit, Box of 10

## REFERENCES

1. CDC. 1999. Update: respiratory syncytial virus activity—United States, 1998-1999 season. *Morbidity and Mortality Weekly Report*. 48:1104-1115.
2. Falsey, A.R. and E. E. Walsh. 2000. Respiratory syncytial virus infection in adults. *Clinical Microbiology Review*. 13: 371-384.
3. Englund, J.A., P. A. Peidra, A. Jewell, K. Patel, B. B. Baxter and E. Whimbey. 1996. Rapid diagnosis of respiratory syncytial virus infections in immunocompromised adults. *J. Clin. Microbiol.* 34:1649-1653.
4. Crowcroft, N.S., F. Cutts and M.C. Zambon. 1999. Respiratory syncytial virus: an underestimated cause of respiratory infection, with prospects for a vaccine. *Communicable Disease Public Health*. 2: 234-241.
5. Hughes, J.D., D. A. Mann and V. V. Hamparian. 1988. Detection of respiratory syncytial virus in clinical specimens by viral culture, direct and indirect immunofluorescence, and enzyme immunoassay. *J. Clin. Microbiol.* 226:588-591.
6. Ahluwalia, G.J., J. Embree, P. McNicol, B. Law and G. W. Hammond. 1987. Comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. *J. Clin. Microbiol.* 25: 763-767.
7. Waner, J.L., N. J. Whitehurst, S. J. Todd, H. Shalaby and L. V. Wall. 1995. Comparison of Directigen RSV with viral isolation and direct immunofluorescence for the identification of respiratory syncytial virus. *J. Clin. Microbiol.* 25: 763-767.
8. Kehl, S.C., K. J. Henrickson, W. Hua and F. Jiang. 2001. Evaluation of the hexaplex assay for detection of respiratory viruses in children. *J. Clin. Microbiol.* 39:1696-1701.
9. Barenfanger, J., C. Drake, N. Leon, T. Mueller and T. Troutt. 2000. Clinical and financial benefits of rapid detection of respiratory viruses: an outcomes study. *J. Clin. Microbiol.* 38:2824-2828.
10. Heikkinen, T., J. Marttila, A. Salmi, and O. Ruuskanen. 2002. Nasal Swab versus Nasopharyngeal Aspirate for Isolation of Respiratory Viruses. *J. Clin. Microbiol.* 40:4337-4339.
11. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections. CLSI, Wayne, PA.
12. Garner, J. S. 1996. Hospital Infection Control Practices Advisory Committee, U. S. Department of Health and Human Services, Center for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
13. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
14. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262, 17/10/2000*, p. 0021-0045.
15. Data on file, BD Diagnostics.
16. Smith, T.F. 2000. Specimen requirements; selection, collection, transport, and processing *In* S. Spector, R.L. Hodinka, and S.A. Young (eds.), *Clinical virology manual*, 3rd ed., ASM Press, Washington, D.C.
17. Clinical and Laboratory Standards Institute (formerly NCCLS). 2002. Approved guideline: EP12-A. User protocol for evaluation of qualitative test performance, CLSI, Wayne, PA.
18. Federal Register. 2001. Quality control for tests of moderate complexity (including the subcategory), high complexity, or any combination of these tests. Code of Federal Regulations 42(3):493.1202.

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©2010 BD.
ViraTrans is a trademark of Trinity Biotech, PLC.
Multitrans is a trademark of Starplex Scientific, Inc.
BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company.
Benex Limited