### Revisions

Rev from	Rev to	ECO #
0503	0604	2912-04

**Notes:**

1. BD Cat. Number 254004
2. Blank (Sheet) Size : Length: 8 1/2"    Width: 11"  
 Number of Pages: 44    Number of Sheets: 11  
 Page Size: Length 8 1/2"    Width 5 1/2"                      Final Folded Size: 8 1/2" x 5 1/2"
3. Style (see illustrations below): # 5

4. See Specification Control Number M765441JAA for Material Information
5. Ink Colors: Printed two sides     Yes     No  
 No. of Colors: 1                      PMS# 327 Teal
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: M765441JAA		Category and Description Package Insert, Culturette Toxin CD	Sheet: 1 of 45 <hr/> Scale:    N/A	A

**BD Culturette™ Toxin CD**

**English:** pages 1 - 8  
**Français :** pages 9 - 16  
**Deutsch:** Seiten 17 - 24  
**Italiano:** pagine 25 - 32  
**Español:** páginas 33 - 40

CLIA COMPLEXITY: HIGH  
 CDC IDENTIFIER CODES  
 ANALYTE: 1022  
 TEST SYSTEM: 07486, 07487

M765441JAA  
 2004/06

See symbol glossary at end of insert. / Viz popis symbolů na konci příbalového letáku. / Se symbolglossaret i slutningen af indlægssedlen. / Zie lijst met symbolen aan het einde van de bijsluiting. / Vaadake sūmbolite seletust infolehe lõpus. / Katso pakkausselosteen lopussa olevaa kuvamerkkien sanastoa. / Voir le glossaire des symboles à la fin de la notice. / Siehe Symbol-Erklärungen am Ende der Packungsbeilage. / Δείτε το γλωσσάριο των συμβόλων στο τέλος του ένθετου. / A jelmagyarázat a használati utasítás végén található. / Vedere il glossario dei simboli alla fine del foglio illustrativo. / Zr. informacino lapelio pabaigoje pateikiama simbolių glosarijų. / Se i symbolforklaringen på slutten av produktvedlegget. / Zobacz objaśnienie symboli na końcu ulotki. / Consulte o glossário de símbolos no fim do folheto informativo. / Pozri slovník symbolov na konci letáka. / Consulte el glosario de símbolos al final del prospecto. / Se symbolförteckningen vid slutet av bipacksedeln.

Pokyny vám poskytné místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar.

#### INTENDED USE

The **Culturette™** Toxin CD test is a rapid 70-min enzyme immunoassay (EIA) for the detection of *Clostridium difficile* toxin A (enterotoxin) in human stool. This assay is intended as an aid in the diagnosis of *C. difficile* associated disease.

#### SUMMARY AND EXPLANATION

*C. difficile* is an important cause of antibiotic-associated diarrhea, which in its most serious form can result in the clinical syndrome of pseudomembranous colitis and significant mortality. Although *C. difficile* may be a part of the normal bacterial intestinal flora, it may become an opportunistic pathogen following the patient's treatment with antibiotics and subsequent alteration of the normal intestinal flora. Under the proper conditions, toxin-producing strains of *C. difficile* produce two toxins – toxin A, a tissue-damaging enterotoxin, and toxin B, an *in vitro* cytotoxin.<sup>1</sup> At the present time, the literature states that both toxin A and toxin B are produced at the same time.<sup>2</sup> The clinical symptoms associated with the disease are thought to be mainly due to toxin A, and to date, there is not convincing evidence that toxin B has any important biologic activity in naturally-occurring disease.<sup>3</sup>

The most common clinical diagnostic aids for *C. difficile* antibiotic-associated colitis have been cell culture cytotoxicity (CTA) and latex agglutination (LA) assays.<sup>4</sup> The CTA assay detects toxin B through the cytopathic effect on cell culture. While very sensitive, this assay requires a minimum of 2 days to complete. Latex agglutination detects the antigens of *C. difficile* rather than the specific toxins, but is regarded as a valuable rapid assay in defining an etiologic role for *C. difficile* in patients with diarrhea.<sup>5</sup> Neither of these assays detect toxin A, the putative cause of disease.

#### PRINCIPLES OF THE PROCEDURE

The **Culturette** Toxin CD test utilizes an anti-toxin A capture antibody coated on microwells. The patient sample and an enzyme-conjugated anti-toxin A detector antibody are added to the microwells and incubated at 35 to 37°C for 1 h. If toxin A is present, a reactive anti-toxin enzyme complex develops. After washing to remove unbound conjugate, a substrate and a chromogen are added and incubated at room temperature for 10 min. Color develops in the presence of bound enzyme. The reactions are then read visually, or spectrophotometrically if desired.

#### REAGENTS

**Culturette** Toxin CD Kit:

<b>Reagent 1</b>	(20.0 mL)	Sample Buffer, protein stabilizer, with 0.02% thimerosal (preservative).
<b>Reagent 2</b>	(10.0 mL)	Conjugate Reagent, polyclonal (rabbit) antibody specific for <i>C. difficile</i> toxin A conjugated to horseradish peroxidase in a buffered protein stabilizer, with 0.01% thimerosal (preservative).
<b>Reagent 3</b>	(225.0 mL)	Wash Reagent, nonreactive detergent in deionized water.
<b>Reagent 4</b>	(10.0 mL)	H <sub>2</sub> O <sub>2</sub> Reagent, buffered solution containing hydrogen peroxide.
<b>Reagent 5</b>	(10.0 mL)	TMB Reagent, solution containing tetramethylbenzidine.
<b>Reagent 6</b>	(10.0 mL)	Stop Reagent, 1.0 N sulfuric acid. <i>Caution:</i> Avoid skin contact; if contact occurs, wash immediately with water.
<b>Control -</b>	(1.9 mL)	Negative Control, contains protein stabilizers and 0.02% thimerosal (preservative).
<b>Control +</b>	(1.9 mL)	Low Positive Control, contains inactivated <i>C. difficile</i> toxin A, protein stabilizers, and 0.01% thimerosal (preservative).
<b>Control ++</b>	(1.9 mL)	High Positive Control, contains inactivated <i>C. difficile</i> toxin A, protein stabilizers, and 0.01% thimerosal (preservative).

Microwell strips coated with polyclonal rabbit antibodies specific for *C. difficile* toxin A.

**Precautions:** For *in vitro* Diagnostic Use.

After review by the Centers for Disease Control and Prevention (CDC), and the Food and Drug Administration (FDA) under CLIA '88, this product has been identified as high complexity. The CDC Analyte Identifier Code is 1022; the CDC Test System Identifier Code is 07486 (visual) and 07487 (spectrophotometric).

**Reagents:** Do not use beyond the expiration date. Upon removal from the refrigerator, allow reagents to warm to room temperature before use. Do NOT mix reagents from different kit lot numbers.

Avoid prolonged exposure of **Reagents 4 and 5** to strong light.

To assure proper drop size delivery, hold the dispensing bottle vertically above the microwell, dispensing one free-falling drop at a time.

Avoid skin contact with **Reagent 6** (1 N H<sub>2</sub>SO<sub>4</sub>). Flush with water immediately if contact occurs.

**Controls:** Do not use the kit if positive and negative controls do not yield appropriate results. Positive controls are made with inactivated toxin A antigen; however, they should be handled as potentially biohazardous material.

**Transfer Pipets:** Single use; do not reuse.

**Storage of Reagents:** Upon receipt, refrigerate reagents at 2 to 8°C. DO NOT FREEZE. Reagents should be recapped and returned to refrigeration when not in use, taking care not to mix color-coded caps.

**Reagent 3 (Wash)** may be stored at room temperature (with original cap) or refrigerated at 2 to 8°C. **Allow 90 min to warm to room temperature following refrigeration.**

**Microwells:** Do not reuse. Perform microwell washing exactly as described under "Performance of Test." Inadequate washing may result in elevated background reading or random color development in toxin A negative samples.

The microwells are contained in a resealable foil envelope designed to protect strips from moisture. The envelope should be opened as described under "Assay Procedure." Keep microwells in the resealable foil envelope until the envelope reaches room temperature. Immediately replace all unused microwells into the resealable foil envelope containing the desiccant, reseal and return to refrigeration.

#### SPECIMEN COLLECTION AND HANDLING

Collect stool specimens into a clean air-tight container with no preservative. Testing of specimens should be performed as soon as possible upon receipt in the laboratory; however, storage for up to 48 h at 2 to 8°C is permissible. If specimens are to be tested after 48 h, specimens should be frozen at -70°C immediately upon receipt in the laboratory. Although the majority of positive specimens will show little or no reduction in toxin A detected after one freeze-thaw cycle,<sup>6</sup> *fresh specimens are preferred.*

Allow stool specimens to warm to room temperature and mix as thoroughly as possible prior to use.

**WARNING:** Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>7-10</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

In keeping with good laboratory practice, wear disposable gloves throughout the assay and wash hands thoroughly afterwards.

Dispose of all materials used in performing the test by autoclaving for 60 min at 121°C or by treatment with a 0.05% solution of sodium hypochlorite (1:100 dilution of household bleach) for 30 min. *Do not autoclave materials containing sodium hypochlorite.* **WARNING:** Liquid waste containing sulfuric acid must be neutralized (e.g., with 1 N sodium bicarbonate) before addition of sodium hypochlorite.

Decontamination of the reusable microwell strip holder between assays can be accomplished by treatment with 0.05% sodium hypochlorite solution as described above.

#### PROCEDURE

**Materials Provided:** All materials as listed under "Reagents," Work Station, plastic Transfer Pipets, microwell Strip Holder, strip sealer, Sample Identification Template and accessories.

**Materials Required But Not Provided:** A 35 to 37°C incubator, plastic or glass test tubes (12 x 75 mm) for sample dilution, vortex mixer, EIA microwell reader capable of reading absorbance at 450 nm or 450/630 nm (optional). Also required are the necessary equipment and labware used for preparation, storage and handling of clinical specimens.


**Performance of Test:** Review "Precautions" and "Specimen Collection and Handling" prior to performing procedures. The testing area, reagents, test specimens and test components should be at room temperature when used.

Assure reagents are in the proper designated wells of the work station. Remove the bottle cap from **Reagent 3 (Wash)**, and replace with the squirt nozzle assembly provided in the kit. After the completion of daily tests, the squirt nozzle assembly on the Wash bottle should be replaced with the original bottle cap.

Gently mix by inverting all reagents prior to use. Avoid foaming.

#### SPECIMEN PREPARATION:

1. Add 200 µL of **Reagent 1** (Sample Buffer) to a clean 12 x 75 mm tube using a calibrated micropipettor or a transfer pipet (up to the second mark from the tip of the pipet).


2. Mix stool sample as thoroughly as possible.

- a. *Liquid, loose or semi-solid stool*: Using a transfer pipet, transfer 100  $\mu$ L (up to the first mark) of sample into the tube containing **Reagent 1**. Using the same transfer pipet, mix by withdrawing and expelling the tube contents 2 – 3 times.

While keeping the transfer pipet in the tube, vortex the tube at moderate speed for 15 sec.

- b. *Solid stool*: Using a wood spatula, transfer approximately 200 mg (5 – 6 mm diameter, "pea" size) of stool specimen into the tube containing **Reagent 1**. Emulsify the sample and discard the wood spatula into a biohazard receptacle.


Vortex tube at moderate speed for 15 sec and place a transfer pipet into the tube.


**Note:** Stool samples, once diluted in **Reagent 1**, should be assayed as soon as possible.

**Preparation of Microwells:**

1. Open the foil envelope containing the microwells by tearing or cutting at the notches between the heat seal area and the resealable band and remove the contents.


2. Remove all of the microwell strips from the holder using the row of buttons at the front of the workstation.


3. Break off the necessary number of microwells (one per patient, plus three controls per batch) and position them in the strip holder.
4. Return all unused microwells to the foil envelope, **reseal**, and immediately return to refrigeration.
5. Identify the location of each well using the workstation Sample Identification Template.

**ASSAY PROCEDURE:**

1. Add 2 free-falling drops of **Reagent 2** (Conjugate) to all microwells.


2. Using the transfer pipet or a calibrated micropipettor, draw up 100  $\mu$ L of diluted stool (up to the first mark) and allow sample to run slowly down the side of the designated microwell. Discard the transfer pipet or micropipettor tip into a biohazard receptacle.


3. Using a separate transfer pipet or micropipettor tip for each control, draw up 100  $\mu$ L of **Control ++**, **Control +** and **Control -** and allow control to run slowly down the side of the designated microwells.

4. Mix the microwells by gently tapping the plate on a firm surface for 10 sec.


5. Trim the microwell strip sealer to size and press firmly atop the microwells.


6. Incubate the plate in a 35 to 37°C incubator for 60 min.

**Note:** Ensure that proper incubator temperature is maintained throughout the 1 h incubation period.


7. Gently remove the microwell strip sealer and wash the microwells with **Reagent 3** (Wash) as follows:

- DECANT by inverting (or aspirate) microwells into a biohazard receptacle.
- FIRMLY SLAP** the inverted plate on clean paper towels placed on a solid surface.
- FILL each microwell with **Reagent 3** (Wash), decant (or aspirate), and firmly slap the inverted plate on clean paper towels.
- REPEAT step "c" 4 additional times.
- After the final wash, vigorously slap the inverted plates on clean paper towels with sufficient force to remove as much excess **Reagent 3** as possible *without dislodging the strips*. Do not allow microwells to completely dry at any time.


- Clean the underbottom of all microwells with a lint-free tissue.
- Add 2 free-falling drops of **Reagent 4** (H<sub>2</sub>O<sub>2</sub>) to each microwell, followed by 2 free-falling drops of **Reagent 5** (TMB).
- Mix the microwells by gently tapping the plate for 10 sec and incubate at room temperature for 10 min.

**READING RESULTS:**

- If visual determinations are to be performed, read blue color within 5 min (after the 10 min incubation period in step 10). If yellow color is preferred for visual reading, add two drops **Reagent 6** (Stop) to each microwell immediately after the 10 min incubation period in step 10, gently tap for 10 sec and read yellow color within 10 min after adding **Reagent 6**.
- If EIA microwell reader determinations are to be performed, clean the underbottom of all microwells, zero the EIA reader on air (without microwells), and read absorbance at 450 nm or 450/630 nm within 10 min after adding **Reagent 6**.

3. After completion of tests, used microwell strips may be removed from the strip holder by pressing the holder against the row of buttons at the front of the workstation. Used microwells should be disposed of as described under "Specimen Collection and Handling."

**Quality Control:** The **Control +** (Low Positive) and **Control -** should be included with each batch of specimens to verify performance of the reagents and proper procedural technique. A **Control ++** (High Positive) is also provided as an optional means of verifying two levels of positive reactions. This is of particular value when results are read visually as an added assurance of procedural integrity.

Patient results should not be reported if positive and negative controls do not yield appropriate results as described in the table opposite.

	VISUAL	450 nm	450/630 nm
Control ++	Strong Color (Blue or Yellow)	>1.500	>1.455
Control +	Moderate Color (Blue or Yellow)	0.300 – 1.000	0.250 – 0.900
Control -	Colorless	<0.100	<0.070

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

#### INTERPRETATION OF RESULTS

The **Culturette** Toxin CD test is designed to enable visual reading of a blue color after step 10 ("Assay Procedure") or alternatively, a yellow color following the addition of **Reagent 6**. The test can also be read by an EIA Microwell Reader.

- Blue Visual Reading:**  
Negative = Colorless  
Positive = Blue color of **any** intensity
- Yellow Visual Reading (Reagent 6 addition):**  
Negative = Colorless  
Positive = Yellow color of **any** intensity
- EIA Microwell Reader, Single Wavelength (450 nm):**  
Negative =  $OD_{450} < 0.1$   
Indeterminate =  $OD_{450} \geq 0.1$  but  $< 0.15$   
Positive =  $OD_{450} \geq 0.15$
- EIA Microwell Reader, Dual Wavelength (450/630 nm):**  
Negative =  $OD_{450/630} < 0.070$   
Indeterminate =  $OD_{450/630} \geq 0.070$  but  $< 0.100$   
Positive =  $OD_{450/630} \geq 0.100$

When using an EIA microwell reader to interpret reactions, indeterminate results should be repeated. If the results are still indeterminate, a fresh patient specimen should be obtained and tested. When reading visually, reactions should be recorded as positive if color of any intensity develops. Colorless reactions are recorded as negative.

**Note:** Extremely strong positive reactions may produce a dark brown precipitate within 5 – 10 min after the addition of **Reagent 6**, which does not interfere with interpretation of results.

#### LIMITATIONS OF THE PROCEDURE

The **Culturette** Toxin CD test detects the presence of toxin A in human stool. Since the level of toxin A has not been shown to have a definitive correlation with the presence or severity of disease, the assay results should be interpreted by a physician in conjunction with clinical and other laboratory findings.

Some isolates of *Clostridium sordellii* have been shown to produce a hemorrhagic toxin (HT) which has similar biologic, physicochemical and immunochemical properties as toxin A. The HT may cross react in tests for toxin A.<sup>4</sup> The *C. sordellii* HT strains have not been detected in patients with antibiotic-associated diarrhea and colitis. Infants and cystic fibrosis patients may have *C. difficile* toxin present in their stool without clinical significance.<sup>11,12</sup>

As with other brands of tests, specimen handling is important for the maintenance of toxin A titers. If testing is delayed, freezing of samples at -70°C is recommended (see "Specimen Collection and Handling").

#### EXPECTED VALUES

*Clostridium difficile* is an opportunistic pathogen that exerts its toxigenic effects when the intestinal tract has been compromised in some manner, such as by antibiotic therapy. Therefore, patients with recent antibiotic therapy or those in chronic care situations are most often infected. Up to 1% of healthy adults may have positive stool toxin tests on culture.<sup>13</sup> The rate of nosocomial infections may vary with the institution, section, housekeeping practices and patient population.<sup>13,14</sup> The **Culturette** Toxin CD test has given positive rates of 3.8% to 18.5% in symptomatic patients.<sup>15</sup> Positive rates above or below these values may be found.

#### PERFORMANCE CHARACTERISTICS

The performance of the **Culturette** Toxin CD test was determined in prospective evaluations conducted at four major independent medical centers located in geographically distinct areas of the United States. A total of 905 fresh stool specimens from symptomatic patients with suspected *C. difficile*-associated disease that were specifically submitted for *C. difficile* cytotoxin B testing were utilized in this evaluation. The mean prevalence of *C. difficile* in patient samples tested was 13%, ranging from 4% to 20%. Cytotoxin B test results were used as the primary reference method. Discordant specimens were retested in duplicate, and if still unresolved, the **Culturette** Toxin CD test result was compared to a consensus result. The consensus result was determined by reviewing the clinical diagnosis, toxigenic culture result and other testing information when available. Results of the study are presented in Table 1 where the **Culturette** Toxin CD test is interpreted by an EIA microwell reader and in Tables 2 and 3 where interpreted by blue and yellow visual color, respectively. The **Culturette** Toxin CD test delivered performance levels of 87 – 92% sensitivity and 95 – 98% specificity depending on the reading method and reference utilized. Indeterminate results are not included in the calculations for sensitivity and specificity.

**Table 1**  
**Comparison of Culturette Toxin CD Test Results Using an EIA Microwell Reader Versus Cytotoxin Test Results and Cytotoxin Results Plus Resolution of Discordants**

Culturette Toxin CD	Cytotoxin Result		Cytotoxin Plus Resolution of Discordants	
	Positive	Negative	Positive	Negative
Positive	99	24	110	12
Negative	17	732	9	760
Indeterminate	5	19	2	3

\*Sensitivity 85% **92%**  
 \*Specificity 97% **98%**  
 Positive Predictive Value 80% **89%**  
 Negative Predictive Value 98% **99%**  
 Indeterminate 3% **<1%**  
 Overall Accuracy 93% **97%**  
 \*\*Total Number 896 **896**

**Table 2**  
**Comparison of Culturette Toxin CD Test Results Using a Visual Blue Color Endpoint Versus Cytotoxin Test Results**

Culturette Toxin CD	Cytotoxin Result	
	Positive	Negative
Positive	103	21
Negative	15	727

\*Sensitivity 87%  
 \*Specificity 97%  
 Positive Predictive Value 83%  
 Negative Predictive Value 98%  
 Overall Accuracy 95%  
 \*\*\*Total Number 871

**Table 3**  
**Comparison of Culturette Toxin CD Test Results Using a Visual Yellow Color Endpoint Versus Cytotoxin Test Results**

Culturette Toxin CD	Cytotoxin Result	
	Positive	Negative
Positive	104	37
Negative	14	710

\*Sensitivity 88%  
 \*Specificity 95%  
 Positive Predictive Value 74%  
 Negative Predictive Value 98%  
 Overall Accuracy 93%  
 \*\*\*Total Number 871

\*Weighted averages based on the number of positive or negative specimens tested at each trial site in relation to the number of specimens tested.

\*\*Spectrophotometric readings were not recorded for 9 specimens.

\*\*\*Visual endpoint results were not recorded for 34 specimens.

The reproducibility of the assay was tested using three control samples containing high, low and zero amounts of toxin A. The concordance of single and dual wavelength readings for positive and negative was 100%. The intra- and inter-assay precision in liquid stool samples is presented in Table 4.

**Table 4: Culturette Toxin CD Assay Precision**

Sample	Mean Single	Mean Dual	%CV Intraassay Single	%CV Interassay Single	%CV Intraassay Dual	%CV Interassay Dual
Control -	0.0440	0.0075	4.34	4.77	10.6	18.67
Control +	0.6692	0.6285	7.6	10.55	8.06	11.06
Control ++	2.38	2.32	9.23	13.92	9.18	13.95

A breakdown of the stool specimen consistency encountered from the symptomatic patients in this study is presented in Table 5. The information is presented from a total specimen perspective and for positive specimens only.

**Table 5: Stool Specimen Consistency**

	Total Population N/%	Positive Specimens N/%
Liquid	100 (11%)	6 (5%)
Loose	371 (41%)	67 (55%)
Semi-Solid	263 (29%)	35 (29%)
Solid	150 (17%)	11 (9%)
Not Classified	21 (2%)	3 (2%)

**Table 6**

Microorganism	Absorbance at 450 nm	
	In Negative Stool	In Positive Stool
<i>Aeromonas hydrophilia</i> ATCC™ 7966	0.054 (-)	0.980 (+)
<i>Bacillus cereus</i> ATCC 11778	0.053 (-)	0.951 (+)
<i>Bacillus subtilis</i> ATCC 33608	0.062 (-)	0.993 (+)
<i>Bacteroides fragilis</i> ATCC 25285	0.054 (-)	0.902 (+)
<i>Candida albicans</i> ATCC 10231	0.053 (-)	0.920 (+)
<i>Clostridium botulinum</i> ATCC 17786	0.056 (-)	0.966 (+)
<i>Clostridium butyricum</i> ATCC 8260	0.057 (-)	0.904 (+)
<i>Clostridium histolyticum</i> ATCC 19401	0.055 (-)	0.959 (+)
<i>Clostridium innocuum</i> ATCC 14501	0.054 (-)	0.990 (+)
<i>Clostridium novyi</i> ATCC 19402	0.054 (-)	0.832 (+)
<i>Clostridium perfringens</i> ATCC 13124	0.055 (-)	0.984 (+)
<i>Clostridium septicum</i> ATCC 12464	0.055 (-)	1.006 (+)
<i>Clostridium sordellii</i> VPI 9048	0.280 (+)	1.124 (+)
<i>Clostridium sporogenes</i> ATCC 3584	0.057 (-)	0.858 (+)
<i>Clostridium subterminale</i> ATCC 29748	0.058 (-)	0.891 (+)
<i>Clostridium tetani</i> ATCC 19406	0.057 (-)	0.839 (+)
<i>Enterococcus faecalis</i> ATCC 29212	0.053 (-)	0.903 (+)
<i>Escherichia coli</i> ATCC 35218	0.058 (-)	0.884 (+)
<i>Escherichia coli</i> ATCC 43889	0.057 (-)	0.297 (+)
<i>Escherichia coli</i> ATCC 43894	0.056 (-)	0.312 (+)
<i>Escherichia coli</i> ATCC 43895	0.058 (-)	0.281 (+)
<i>Peptostreptococcus anaerobius</i> ATCC 27337	0.061 (-)	0.808 (+)
<i>Pseudomonas aeruginosa</i> ATCC 14207	0.056 (-)	0.835 (+)
<i>Salmonella choleraesuis</i> ATCC 14028	0.051 (-)	0.838 (+)
<i>Shigella dysenteriae</i> ATCC 29027	0.054 (-)	0.942 (+)
<i>Shigella flexnerii</i> ATCC 12661	0.053 (-)	0.941 (+)
<i>Shigella sonnei</i> ATCC 25931	0.056 (-)	0.969 (+)
<i>Staphylococcus aureus</i> ATCC 12598	0.051 (-)	1.069 (+)
<i>Vibrio cholerae</i> ATCC 11623	0.052 (-)	0.944 (+)
<i>Vibrio parahaemolyticus</i> ATCC 17802	0.053 (-)	0.799 (+)
<i>Yersinia enterocolitica</i> ATCC 23715	0.056 (-)	1.016 (+)

**Cross Reactivity:** The **Culturette** Toxin CD test reacts with all known toxigenic reference strains of *C. difficile*.<sup>15</sup> The test did not react with the nontoxigenic strain of VPI 11186. Cross reactivity was demonstrated to be negative with the bacteria shown in Table 6. Testing was performed by taking *C. difficile* toxin A negative and positive stools and adding the various organisms to final concentrations of 10<sup>7</sup> to 10<sup>9</sup> CFU/mL stool. As expected, the only organism shown to cross react in the test is a highly toxigenic isolate of *Clostridium sordellii* VPI 9048. This isolate elaborates high levels of hemorrhagic and lethal toxins which have been shown to be immunologically and biologically similar to *C. difficile* toxins A and B. The *Staphylococcus aureus* Cowan strain ATCC 12598, which produces protein A, did not show cross reactivity. Also, *Escherichia coli* ATCC 43889, 43894, and 43895 which produce Shiga-like toxins (SLT) did not show any cross reactivity.

#### AVAILABILITY

Cat. No.	Description
254004	<b>Culturette</b> <sup>™</sup> Toxin CD, 96 test kit.
264001	<b>Culturette</b> <sup>™</sup> <b>CDT</b> <sup>™</sup> Rapid <i>C. difficile</i> latex test, 25 test kit.
273310	Pipette Tips, box of 1000.

#### REFERENCES

- Sullivan, N.M., S. Pellett and T.D. Wilkins. 1982. Purification and characterization of toxin A and B of *Clostridium difficile*. *Infect. Immun.* 35:1032-1040.
- Lyerly, D.M., N.M. Sullivan and T.D. Wilkins. 1983. Enzyme-linked immunosorbent assay for *Clostridium difficile* toxin A. *J. Clin. Microbiol.* 17:72-78.
- Bartlett, J.G. 1990. *Clostridium difficile*: clinical considerations. *Rev. Infect. Dis.* 12:S243-S251.
- Lyerly, D.M., H.C. Krivan and T.D. Wilkins. 1988. *Clostridium difficile*: its disease and toxins. *Clin. Microbiol.* 1:1-18.
- Baron, E.J. 1989. Assessment of currently-available laboratory tests for *Clostridium difficile*-associated diarrhea. *Clin. Microbiol. Newsl.* 11:118-120.
- Lyerly, D.M., D.E. Lockwood, S.H. Richardson and T.D. Wilkins. 1982. Biological activities of toxins A and B of *Clostridium difficile*. *Infect. Immun.* 35:1147-1150.
- National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
- U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021-0045.
- COOPERSTOCK, M. 1988. *CLOSTRIDIUM DIFFICILE* IN INFANTS AND CHILDREN. *IN* R.O. ROLFE, AND S.M. FINEGOLD (EDS.). *CLOSTRIDIUM DIFFICILE: ITS ROLE IN INTESTINAL DISEASE*. ACADEMIC PRESS, INC., SAN DIEGO, CA., P. 45-64.
- Peach, S.L., S.P. Borriello, H. Gayla, F.E. Barclay, and A.R. Welch. 1968. Asymptomatic carriage of *Clostridium difficile* in patients with cystic fibrosis. *J. Clin. Pathol.* 39:1013-1018.
- Delmee, M., B. Vandercam, V. Avesani and J.L. Michaux. 1987. Epidemiology and prevention of *Clostridium difficile* infections in a leukemia unit. *Eur. J. Clin. Microbiol.* 6:623-627.
- Gerding, D.M. 1989. Disease associated with *Clostridium difficile* infection. *Ann. Intern. Med.* 110:255-257.
- Data on file at BD Diagnostics.

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ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, CDT and Culturette are trademarks of Becton, Dickinson and Company © 2004 BD.
BENEX Limited