Revisions

SO 0191-5

Rev from	Rev to	ECO #
0603	2010/07	5441-10

Notes:

1. BD Cat. Number 211737, 211738, 221001, 221031, 221041, 221051, 221061
2. Blank (Sheet) Size: Length: 9" Width: 32"
 Number of Pages: 16 Number of Sheets: 1
 Page Size: Length 9" Width 4" Final Folded Size: 2 1/4" x 4"
3. Style (see illustrations below): # 4

4. See Specification Control Number S1160 for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# 2755 (BD Blue)
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	<small>COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION</small>	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: S1160JAA		Category and Description Package Insert, Difco QC Antigen Shigella	Sheet: 1 of 17 <hr/> Scale: N/A	A

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INTENDED USE

Difco™ QC Antigens Shigella are used in the quality control testing of Difco Shigella Antisera Poly by slide agglutination tests.

SUMMARY AND EXPLANATION

The use of Difco Shigella Antisera Poly in the serologic identification of *Shigella* requires the use of quality control test suspensions to verify that the antisera are performing as expected. Most laboratories are required to test antisera with positive and negative controls prior to use.^{1,2} Difco QC Antigens Shigella are chemically stabilized and inactivated suspensions of known strains of the genus *Shigella* designed for use as positive controls for testing the efficacy of the *Shigella* grouping antisera employed in routine laboratory procedures.

Difco QC Antigens Shigella may also be used as negative controls by using a heterologous antigen (possessing no common antigen) with a given test antiserum. However, cross-reactivity may occur. Consult appropriate references for further details on cross-reactivity.¹

The chart below includes the Difco QC Antigens Shigella recommended for the homologous (positive) control antigen for Difco Shigella Antisera Poly. The homologous control antigen has in common certain identifying antigen(s) with the antiserum.

Difco Antiserum	Difco QC Antigen Homologous Control
Shigella Antiserum Poly Group A	Shigella Group A
Shigella Antiserum Poly Group A ₁	Shigella Group A ₁
Shigella Antiserum Poly Group B	Shigella Group B
Shigella Antiserum Poly Group C	Shigella Group C
Shigella Antiserum Poly Group C ₁	Shigella Group C ₁
Shigella Antiserum Poly Group C ₂	Shigella Group C ₂
Shigella Antiserum Poly Group D	Shigella Group D
Alkalescens-Dispar Group 1	Alkalescens-Dispar Group 1

PRINCIPLES OF THE PROCEDURE

Serologic procedures that confirm the identification of an organism are usually agglutination reactions. Agglutination reactions may be either homologous or heterologous. Homologous reactions occur between a microorganism (antigen) and the corresponding antibody. These reactions occur rapidly and are strong. Heterologous reactions occur when a microorganism (antigen) reacts with an antibody produced in response to some other species or serotype. These reactions occur slowly and are weak.

Heterologous reactions may be unexpected and unpredictable and may lead to confusion in serologic identification. Therefore, only strongly positive homologous agglutination reactions should be regarded as significant.

REAGENTS

Difco QC Antigen	Organism Used for Antigen Preparation
Shigella Group A	<i>Shigella dysenteriae</i> type 1
Shigella Group A ₁	<i>Shigella dysenteriae</i> type 8
Shigella Group B	<i>Shigella flexneri</i> type 6
Shigella Group C	<i>Shigella boydii</i> type 3
Shigella Group C ₁	<i>Shigella boydii</i> type 8
Shigella Group C ₂	<i>Shigella boydii</i> type 15
Shigella Group D	<i>Shigella sonnei</i>

When used as described (see PROCEDURE) each vial of antigen is sufficient for 20 slide tests.

Difco QC Antigens Shigella contain 0.5% formaldehyde as a preservative.

Warnings and Precautions

For *in vitro* Diagnostic Use.

The Packaging of This Product Contains Dry Natural Rubber.

Observe aseptic technique and established precautions against microbiological hazards throughout all procedures. After use, specimens, containers, slides, tubes and other contaminated materials must be sterilized by autoclaving. Directions for use should be followed carefully.

Difco QC Antigens Shigella are not to be used for immunization of man or animals.

Storage

Store **Difco QC Antigens Shigella** at 2 – 8°C.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

Expiration date applies to product in its intact container when stored as directed. Do not use if the product is discolored or shows other signs of deterioration.

Difco QC Antigens Shigella are stable through the expiration date on the label when stored as described. Exposure to temperatures below 2°C can result in autoagglutination. Antigens must be smooth uniform suspensions; examine antigen vials for agglutination before use. Suspensions with agglutination are not usable and should be discarded.

PROCEDURE

Materials Provided: **Difco QC Antigens Shigella**

Materials Required But Not Provided: Agglutination slides, Applicator sticks, Sterile 0.85% saline, **Difco Febrile Negative Control** (optional).

Reagent Preparation

Difco QC Antigens Shigella are ready for use. Shake well before use to suspend the organisms.

Before using **Difco QC Antigens Shigella**, examine the **Difco Shigella Antiserum Poly** to determine that the antisera meet all specifications for the product (see product information for **Difco Shigella Antisera Poly**).

Equilibrate all materials to room temperature prior to performing the tests.

Ensure that all glassware and pipettes are clean and free of residues such as detergents.

Test Procedure

1. Positive control – Dispense 1 drop (approximately 35 µL) of the **Difco Shigella Antiserum Poly** to be tested on an agglutination slide. After shaking well, add 1 drop of the appropriate **Difco QC Antigen Shigella** and mix thoroughly.
2. Negative control – Dispense 1 drop of sterile 0.85% saline or **Difco Febrile Negative Control** (optional) on the agglutination slide. After shaking well, add 1 drop of the appropriate **Difco QC Antigen Shigella** and mix thoroughly.
3. Rotate the slides for 1 min and read for agglutination.

User Quality Control

At the time of use, apply both homologous and heterologous controls to check performance of the antigen, techniques and methodology.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

Interpretation of Results

Read and record results as follows:

4+	100% agglutination (background clear to slightly hazy)
3+	75% agglutination (background slightly cloudy)
2+	50% agglutination (background moderately cloudy)
1+	25% agglutination (background cloudy)
-	No agglutination

LIMITATIONS OF THE PROCEDURE

1. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False positive reactions can occur.
2. Exposure to temperatures below 2°C can cause autoagglutination. Antigens must be smooth, uniform suspensions. Examine antigen vials for

agglutination before use. Suspensions with agglutination are not usable and should be discarded.

3. Allow the **Difco** QC Antigens, the antisera and all equipment used to be at room temperature at the time of testing. The test reagents, if cold may cause false-negative reactions.
4. Shake the antigen well before use to suspend the organisms.
5. Rough reactions can occur when using the **Difco** QC Antigen Shigella and 0.85% saline. If this occurs, repeat the negative control test using **Difco** Febrile Negative Control.

EXPECTED RESULTS AND PERFORMANCE CHARACTERISTICS

1. The positive control should show 3+ or greater agglutination.
2. The negative control should show no agglutination.

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Difco is a trademark of Difco Laboratories, Inc., a subsidiary of Becton, Dickinson and Company. BD and BD Logo are trademarks of Becton, Dickinson and Company. © 2010 BD.
Benex Limited