

BD GeneOhm™ MRSA ACP

Nasal Specimen Collection Protocol


one



Use the recommended swab to collect the nasal specimen.


The swab can be moistened with two drops (about 50 µL) of sterile physiological serum (saline) or used dry.

two




Carefully insert the swab into the patient's nostril (the swab tip should be inserted up to 2.5 cm (1 inch) from the edge of the nares). Roll the swab 5 times along the mucosa inside the nostrils.

three



Insert the same swab into the second nostril and repeat sampling as in the preceding step.

four



Return swab to its container and send to the laboratory immediately.



US:
BD Diagnostics
11085 N. Torrey Pines Road, Suite 210
La Jolla, California 92037
www.bd.com/geneohm
888-GENE0HM

Belgium:
BD Diagnostics
Erembodegem- Dorp 86
9230 Erembodegem
(+32) 53 720 211

UK:
BD Diagnostics
The Danby Building
Edmund Halley Road
Oxford Science Park
Oxford OX4 4DQ
(+44) 1865 78 16 66

Germany:
BD Diagnostics
Tuillastr. 8-12
69126 Heidelberg
(+49) 6221 3050

France:
BD Diagnostics
11 rue Aristide Berges
38800 Le Pont de Claix
(+33) 476 68 36 36

Not for Lab Use · See package insert for test protocol.

BD, BD logo and all other trademarks are property of Becton, Dickinson and Company. Copyright © 2009 BD
Patent protection pending