BD GeneOhm™ Cdiff Assay

REF 441400
REF 441401

48 Tests
200 Tests

TABLE OF CONTENTS

ENGLISH	4-15
INTENDED USE	4
SUMMARY AND EXPLANATION OF THE TEST	4
PRINCIPLE OF THE PROCEDURE	4
REAGENTS	5
PRECAUTIONS	6
MATERIALS PROVIDED	6
STORAGE, HANDLING AND STABILITY	7
COLLECTED SPECIMENS.....	7
REAGENTS.....	7
MATERIALS REQUIRED BUT NOT PROVIDED	7
INSTRUCTIONS FOR USE	8
SPECIMEN COLLECTION	8
SPECIMEN PREPARATION	8
BD GENE ^{OHM} ™ CDIFF ASSAY PROCEDURE	9
QUALITY CONTROL	10
POSITIVE AND NEGATIVE CONTROLS	10
SPECIMEN PROCESSING CONTROLS.....	10
CULTURING OF CLINICAL SPECIMENS	10
INTERPRETATION OF RESULTS	10
INVALID ASSAY RUN	11
UNRESOLVED SPECIMEN	11
SPECIMEN NOT DETERMINED DUE TO I-CORE® MODULE FAILURE	11
LIMITATIONS OF THE PROCEDURE	11
INTERFERING SUBSTANCES	12
PERFORMANCE CHARACTERISTICS	12
ANALYTICAL SENSITIVITY	14
ANALYTICAL SPECIFICITY.....	14
REPRODUCIBILITY.....	14
PRECISION.....	145
FRANÇAIS	16-27
APPLICATION	16
RÉSUMÉ ET EXPLICATION DU TEST	16
PRINCIPE DE LA PROCÉDURE	16
RÉACTIFS	17
PRÉCAUTIONS	17
MATÉRIEL FOURNI	18

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English

Intended Use

The BD GeneOhm™ Cdiff Assay is a rapid *in vitro* diagnostic test for the direct, qualitative detection of *C. difficile* toxin B gene (*tcdB*) in human liquid or soft stool specimens from patients suspected of having *Clostridium difficile*-associated disease (CDAD). The test, based on real-time PCR, is intended for use as an aid in diagnosis of CDAD. The test is performed directly on the specimen, utilizing polymerase chain reaction (PCR) for the amplification of specific targets and fluorogenic target-specific hybridization probes for the detection of the amplified DNA.

Summary and Explanation of the Test

A liquid or soft stool specimen is collected and transported to the laboratory. A sterile dry swab is dipped into the liquid or soft stool material and processed. For testing, the swab is eluted in sample buffer and the specimen is lysed. An aliquot of the lysate is added to PCR reagents which contain the *tcdB* specific primers used to amplify the genetic target of *Clostridium difficile*, if present. The assay also includes an internal control (IC) to detect PCR inhibited specimens and to confirm the integrity of assay reagents. Amplified targets are detected with hybridization probes labelled with quenched fluorophores (molecular beacons). The amplification, detection and interpretation of the signals are done automatically by the Cepheid SmartCycler® software. The entire procedure takes about 75 to 90 minutes, depending on the number of specimens processed.

C. difficile is the etiologic agent of several diseases. Epidemiology studies demonstrate that an increased number of *C. difficile* associated outbreaks have been reported worldwide¹, some with increased mortality and morbidity². Disease symptoms range from mild nuisance diarrhea to severe colitis and even to bowel perforation and death by dehydration. This pathogen is the major cause of antibiotic-associated diarrhea (AAD) and pseudomembranous colitis³. The predisposing risk factors are numerous and include age, length and number of hospital stays, invasive medical procedures, immunosuppressive treatments, and chronic pathologies (diabetes, cardiovascular syndromes or AIDS)⁴. These factors explain why some wards (intensive care, surgery, long-term care) in hospitals or institutions are more affected than others. The bacterium *C. difficile*, as well as other intestinal flora bacteria, are killed by antimicrobial treatment but not the *C. difficile* spores, which are insensitive to the majority of antimicrobial agents. When antimicrobial concentration is lower than the minimum inhibitory concentration, spores germinate into viable bacteria which produce toxins⁵.

The diagnosis of toxigenic *C. difficile* is usually done by tissue culture cytotoxicity assay and/or by *C. difficile* culture identification and/or by enzyme immunoassay (EIA). The tissue culture cytotoxicity assay and *C. difficile* culture identification are labor and time consuming and results are obtained within 24 to 96 hours. Toxin EIA assays display low sensitivity and glutamate dehydrogenase antigen enzyme immunoassays (GDH EIA) display low specificity. Development of sensitive and specific molecular amplification techniques now allows the detection of only a few copies of bacterial DNA in clinical samples, and the discrimination of specific species⁶. In addition, rapid PCR technology can achieve this in about 1 hour. The combination of these characteristics may allow the prompt targeted treatment of CDAD patients and thus a potential improved patient outcome and reduced recovery time.

Principle of the Procedure

Following specimen lysis, amplification of the *tcdB* target, if present, occurs. Amplification of the internal control (IC), a DNA fragment of 333 base pairs (bp) including a 277 bp sequence not found in *C. difficile*, will also take place unless there are PCR inhibitory substances.

The amplified DNA target is detected with a molecular beacon, a hairpin-forming single-stranded oligonucleotide labelled at one end with a quencher and at the other end with a fluorescent reporter dye (fluorophore). In the absence of target, the fluorescence is quenched. In the presence of target, the hairpin structure opens upon beacon/target hybridization, resulting in emission of fluorescence. For the detection of *tcdB* amplicons, the molecular beacon contains the fluorophore FAM at the 5' end and the non-fluorescent quencher DABCYL at the opposite 3' end of the oligonucleotide. For the detection of the IC amplicons, the molecular beacon contains the fluorophore TET at the 5' end and the quencher moiety DABCYL at the 3' end. Each beacon-target hybrid fluoresces at a wavelength characteristic of the fluorophore used in the particular molecular beacon. The amount of fluorescence at any given cycle, or following cycling, depends on the amount of specific amplicons present at that time. The SmartCycler® software simultaneously monitors the fluorescence emitted by each molecular beacon, interprets all data, and provides a final result at the end of the cycling program (see Section "Interpretation of Results").

Reagents

BD GeneOhm™ Cdiff Assay	48 Tests	200 Tests
<i>Sample Buffer</i>	60 X 1 mL	240 X 1 mL
Tris-EDTA buffer		
<i>Lysis tube</i>	50 tubes	200 tubes
Glass beads		
<i>Master Mix</i>	8 tubes	34 tubes
< 0.0005% DNA polymerase complex		
< 0.001% Internal Control (non-infectious DNA containing <i>tcdB</i> primer binding sequences with a unique sequence for probe hybridization)		
< 0.002% Primers		
< 0.002% Molecular probes		
< 0.05% Nucleotide mix (dATP, dCTP, dGTP, dTTP)		
Bovine serum albumin		
Carbohydrate		
MgCl ₂		
< 0.001% Non-infectious genomic DNA of <i>Escherichia coli</i> (ATCC 25922)		
<i>Control DNA</i>	8 tubes	34 tubes
Tris-EDTA buffer		
Carbohydrate		
< 0.001% Non-infectious genomic DNA of <i>C. difficile</i> bearing <i>tcdB</i> gene (ATCC 43255)		
<i>Diluent</i>	8 X 700 µL	34 X 700 µL
Tris-HCl buffer		
MgCl ₂		
(NH ₄) ₂ SO ₄		
KCl		

Precautions

- This test is for *in vitro* diagnostic use only.
- Do not use the kit if the outer carton safety seal is broken.
- Do not use reagents if their protective pouches are open or torn upon arrival.
- Close protective pouches of Master Mix and Control DNA quickly with the zip seal after each use.
- Do not remove desiccant from Master Mix and Control DNA pouches.
- Do not use reagents if desiccant is not present inside Master Mix and Control DNA pouches.
- Reagents are not interchangeable between lots.
- Never pool reagents from different tubes even if they are from the same lot.
- Do not use the reagents after their expiration date.
- Do not interchange caps among reagents as contamination may occur and compromise test results.
- Avoid microbial and deoxyribonuclease (DNase) contamination of reagents when removing aliquots from tubes. The use of sterile DNase-free disposable filter-blocked or positive displacement pipettor tips is recommended.
- To avoid contamination of the environment with *tcdB* gene amplicons, do not open the reaction tubes post-amplification.
- Use a new pipettor tip for each specimen or reagent.
- Performing the assay outside of the recommended time ranges can produce invalid results. Assays not completed within specified time ranges should be repeated.
- Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organizations.
- In cases where open-tube PCR tests are conducted in the same general area by the laboratory, separated and segregated working areas should be used for specimen preparation and amplification/detection activities. Supplies and equipment should be dedicated to each area and should not be moved from one area to another. Gloves must always be worn and must be changed before going from one area to another. Gloves must be changed before manipulating lyophilized reagents.
- Always handle specimens as if they are infectious and in accordance with safe laboratory procedures such as those described in *Biosafety in Microbiological and Biomedical Laboratories*⁷ and in the CLSI Document M29⁸.
- Wear protective clothing and disposable gloves while handling kit reagents. Wash hands thoroughly after performing the test.
- Do not pipet by mouth.
- Do not smoke, drink, or eat in areas where specimens or kit reagents are being handled.
- Dispose of unused reagents and waste in accordance with local, state and/or federal regulations.
- Users should complete proficiency testing and training prior to performing the BD GeneOhm™ Cdiff Assay according to local, state and/or federal regulations or accrediting organizations.

Materials Provided

- Sample Buffer
- Lysis tube
- Master Mix
- Control DNA
- Diluent
- SmartCycler® reaction tubes, 25 µL
- Specimen identification labels

Storage, Handling and Stability

Collected Specimens

- Specimens should be kept between 2 °C and 25 °C during transport. Protect against freezing or exposure to excessive heat.
- Specimens can be stored up to 5 days at 2-8 °C before testing. Specimens can be kept at room temperature (15-25 °C) up to 48 hours before testing. Specimens can be tested after one freeze and thaw cycle.

Reagents

Note: Storage conditions must follow the specifications written on each pouch.

Kit Component		Master Mix, Diluent and Control DNA (white, black strip and red strip labels, respectively)	Lysis Tubes (yellow cap)	Sample Buffer (blue cap)
Sealed pouch	Temperature	2-8 °C	2-25 °C	2-25 °C
	Stability	Expiration date	Expiration date	Expiration date
Opened pouch ^A	Temperature	2-8 °C	2-25 °C	2-25 °C ^B
	Stability	1 month ^C	Expiration date	2 months ^C

^A Once the original seal on the pouch is broken, carefully close the pouch with the zip seal after each use and store at appropriate temperature.

^B Although these reagents can be stored at room temperature, they should be kept with their accompanying reagents of the same lot between 2-8 °C.

^C Provided the bags are properly closed with the zip seal after each use.

Kit Components outside of their protective pouch		Master Mix and Control DNA (white and red strip labels)	
Tubes containing unreconstituted reagents	Temperature	15-25 °C	
	Stability	2 hours	
Tubes containing reconstituted reagents ^A	Original container	Temperature	2-8 °C
		Stability	3 hours
	SmartCycler [®] tube	Temperature	2-8 °C
		Stability	1 hour

^A Discard unused tubes after expiration of the indicated stability.

Materials Required but not Provided

- Dry sterile swab
- Dry sterile container for the collection of liquid or soft stool specimens
- Vortex Genie 2 (Scientific Industries Inc.) with 1.5 mL microtube holder or equivalent; for processing multiple samples, adapter with multiple holding sites can be used
- Micropipettors (accurate range between 1-10 µL, 10-100 µL and 100-1000 µL)
- Sterile DNase-free filter-blocked or positive displacement micropipettor tips
- DNase free microcentrifuge tubes
- Scissors (optional)
- Gauze
- Disposable gloves, powderless
- Microcentrifuge for low speed centrifugation
- Dry heating block for 1.5 mL tubes or water bath
- Ice or cooling block for 1.5 mL tubes
- Stopwatch or timer
- SmartCycler[®] starter system with Dx Software (processing block, user manual⁹, accessory kit, and desktop computer) (Cepheid, Sunnyvale, CA, USA)

Instructions for Use

Specimen Collection

In order to obtain an adequate specimen, the procedure for specimen collection must be followed closely.

Liquid stool or soft stool specimen

Using a dry sterile container, liquid stool or soft stool specimens are collected according to the following procedure:

1. **Transfer liquid or soft stool (but not urine) into the container. Avoid mixing toilet paper, water, or soap with the sample.**
2. **Label the container.**
3. **Ship the container to the laboratory according to hospital standard operating procedures.**
4. **Refer to the section entitled Storage, Handling and Stability – Collected Specimens for storage and handling.**

Specimen Preparation

NOTE: One (1) sample buffer tube (**blue cap**) and one (1) lysis tube (**yellow cap**) are required for each specimen to be tested. An additional sample buffer tube (**blue cap**) is also required to dilute specimens. This sample buffer tube will also be used for the BD GeneOhm™ Cdiff Assay procedure. Remove the required number of tubes from their protective pouches, remove the excess air and close the pouches quickly with the zip seal.

1. **Prior to starting specimen preparation for the BD GeneOhm™ Cdiff Assay, vortex the specimen at high speed for 15 seconds and dip a sterile dry swab into the fecal material for testing. Remove excess stool.**
2. **Place the swab in a sample buffer tube (blue cap).**
Identify the sample buffer tube on the cap and/or the tube label.
3. **Break the swab stem and close the sample buffer tube tightly.**
Hold the swab by the stem near the rim of the tube (use gauze to minimize risks of contamination). Lift the swab a few millimeters from the bottom and bend the stem against the edge of the tube to break it. Alternative method: use clean scissors to cut the stem. Make sure the cap will close tightly.
4. **Vortex the sample buffer tube (containing the swab) at high speed for one (1) minute.**
For processing multiple samples, adapters with multiple holding sites can be used.
5. **Add 40 µL from the additional sample buffer tube (blue cap) to the lysis tube (yellow cap) in order to dilute the specimen.**
6. **Transfer 10 µL of cell suspension (from Step 4) to the lysis tube (yellow cap) already containing 40 µL of sample buffer (blue cap) (from Step 5).**
7. **Vortex the lysis tube for five (5) minutes at high speed.**
For processing multiple samples, adapters with multiple holding sites can be used.
8. **Centrifuge the lysis tube briefly (quick spin).**
Centrifuge at low speed for two (2) to five (5) seconds to bring the solid contents at the bottom of the tube.
9. **Heat the lysis tube at 95 ± 2 °C between five (5) and seven (7) minutes.**
Use a dry heating block for 1.5 mL tubes or a water bath.
10. **Place the lysis tube on ice or on a cooling block.**
Lysates are stable up to four (4) hours at 2-8 °C.

BD GeneOhm™ Cdiff Assay Procedure

NOTE: One (1) reconstituted Master Mix tube (**white label**) will yield enough reagents to run **eight (8)** reactions. Allow one SmartCycler® tube per specimen to be tested and two (2) additional SmartCycler® tubes for the **positive and the negative controls**. One (1) **positive** and one (1) **negative** control must be included in each BD GeneOhm™ Cdiff Assay run. One (1) Control DNA (**red strip label**) is required per assay run. One (1) diluent tube (**black strip label**) is required for the reconstitution of up to three (3) Master Mix tubes. Remove the required number of tubes from their protective pouches, **remove the excess air, and close the pouches quickly with the zip seal**.

Prepare only enough SmartCycler® tubes to fill available I-CORE® modules on the SmartCycler® instrument.

- 1. Place the required number of Master Mix tubes on ice or on a cooling block for 1.5 mL tubes.**
- 2. Add 225 µL of diluent (black strip label) to each Master Mix tube.**
Insert the micropipettor tip through the septum of the cap of the Master Mix tube. Do not insert the tip too deeply into the cap. Deliver the diluent. Discard the unused diluent afterward.
- 3. Vortex the tube(s) for 5-10 seconds.**
- 4. Place the tube(s) on ice or on a cooling block for 1.5 mL tubes until ready to use.**
Place a Control DNA tube (red strip label) on ice or on a cooling block for 1.5 mL tubes.
- 5. Add 225 µL of sample buffer (blue cap) to the Control DNA tube.**
Use the additional sample buffer tube (blue cap) from the specimen preparation (Step 4). Insert the micropipettor tip through the septum of the cap of the Control DNA tube. Do not insert the tip too deeply into the cap. Deliver the sample buffer.
- 6. Vortex the tube for 5-10 seconds.**
Place the tube on ice or on a cooling block for 1.5 mL tubes until ready to use.
- 7. Place the required number of SmartCycler® tubes on the SmartCycler® cooling block.**
Allow one (1) SmartCycler® tube per specimen and two (2) more SmartCycler® tubes for the controls. Avoid touching the optical detection windows at the bottom edges of the tube and the lower diamond-shaped area.

THE FOLLOWING STEPS MUST BE PERFORMED WITHIN A ONE (1) HOUR TIME FRAME:

- 8. Add 25 µL of reconstituted Master Mix to the SmartCycler® tubes (appropriate pipetting technique is required to ensure proper transfer of the solution).**
Remove the cap before pipetting the reagent. Deliver the liquid into the reservoir (upper part) of the SmartCycler® tubes. Identify the SmartCycler® tubes on the cap. Specimen identification labels can be used (provided with the kit). Discard the unused Master Mix.
- 9. Add 3.0 µL of each lysed specimen to a different SmartCycler® tube (previously filled with reconstituted Master Mix); close the tubes.**
Be careful not to aspirate beads when pipetting into the lysis tube. After addition of the specimen, pipet up and down 2-3 times in the reservoir to ensure transfer of the complete volume. Use a new micropipettor tip for each specimen.
- 10. Add 3.0 µL of the reconstituted Control DNA to the next to last SmartCycler® tube (Positive Control); close the tube.**
After addition of the DNA, pipet up and down 2-3 times in the reservoir to ensure transfer of the complete volume. Identify as the positive control. Discard the unused control DNA.
- 11. Add 3.0 µL of sample buffer (blue cap) to the last SmartCycler® tube (Negative Control); close the tube.**
Use the sample buffer tube from Step 5. This will monitor PCR contamination that might occur during the manipulation of the specimens. Identify as the negative control. Discard the unused sample buffer.
- 12. Centrifuge all reaction tubes for 5-10 seconds.**
Use the specially adapted microcentrifuge provided with the SmartCycler® instrument.
- 13. Keep the tubes at 2-8 °C on the SmartCycler® cooling block before loading on the instrument.**
The remaining lysates should be frozen at -20 ± 5 °C for later use, if necessary.
- 14. Create a run with the BD GeneOhm™ Cdiff Assay protocol.**
Refer to the SmartCycler® Dx Software Operator Manual⁹ if needed. It is recommended that identification parameters for the specimens be entered before starting the run.
- 15. Insert each reaction tube into an I-CORE® module of the SmartCycler® and close the I-CORE® lid.**
Place the positive and negative controls at their appropriate positions (see the section entitled "Quality control"). Press all the tubes firmly down into place.
- 16. Start the run.**

Quality Control

Positive and Negative Controls

Quality control procedures are designed to monitor assay performance. The positive control is intended to monitor substantial reagent failure. The negative control is used to detect reagent or environmental contamination (or carry-over) by either DNA containing *tcdB* genes or *tcdB* amplicons. Positive and negative controls are assay controls (run controls). An invalid control invalidates the run. Finally, an internal control incorporated into each reaction mixture is intended to monitor the reagent integrity and PCR inhibition in each specimen.

One positive control and one negative control must be included in each assay run on the SmartCycler®. The software automatically assigns the position of the controls on the instrument (refer to the SmartCycler® Dx Software Operator Manual⁹).

Specimen Processing Controls

Additional control strains may be tested according to guidelines or requirements of local, state and/or federal regulations or accreditation organizations. A reference toxigenic *C. difficile* strain bearing the *tcdB* gene (e.g. American Type Culture Collection, ATCC 43255 or a well characterized *C. difficile* clinical isolate identified to be carrying the *tcdB* gene) may be used as a positive specimen processing control while a culture of a non toxigenic *C. difficile* strain (e.g. ATCC 700057) may be used as a negative specimen processing control.

Colonies are isolated after 18 to 24 hours of anaerobic incubation on 5% sheep blood agar. Resuspend colonies in saline to a turbidity of 0.5 McFarland (~1.5 X 10⁷ CFU/mL). Dilute with saline to obtain a suspension of ~10⁶ CFU/mL. Dip a dry swab into the bacterial suspension, press out the excess fluid. Process and test as a clinical specimen (refer to the sections entitled "Specimen Preparation" and "BD GeneOhm™ Cdiff Assay Procedure"), including controls. All specimens and controls should yield valid results (no invalid positive or negative control; no failed internal control; and no incorrect specimen processing control results – when specimen processing controls are performed). In the event of an incorrect specimen processing control result, it is recommended that a new aliquot be obtained from the original stool specimen, and that this specimen be retested along with new controls before reporting results.

For general QC guidance, the user may wish to refer to CLSI MM3¹⁰ and C24¹¹.

Culturing of Clinical Specimens

To perform species identification directly from stools, clinical specimens can be cultured using the hospital procedures.

Interpretation of Results

The decision algorithm for the BD GeneOhm™ Cdiff Assay is embedded in the SmartCycler® DX software. The interpretation of assay results is done according to the following criteria:

Sample Type	Instrument-Reported Assay Result ⁹	Instrument-Reported IC Result ⁹	User Interpretation of Results
Clinical Specimen	NEG	PASS	No <i>tcdB</i> gene DNA detected
	POS	NA	<i>tcdB</i> gene DNA detected
	Unresolved	FAIL	Unresolved – inhibitory specimen or reagent failure
	ND	ND	Not determined due to I-CORE® Module failure (with Warning or Error Codes ⁹)
Positive Control	Valid	NA	Valid Positive Control; valid run when Negative Control is also valid.
	Invalid	NA	Invalid Positive Control; invalid run. ^A Assay results are invalid and must not be reported.
Negative Control	Valid	PASS	Valid Negative Control; valid run when Positive Control is also valid.
	Invalid	FAIL	Invalid Negative Control; invalid run. ^A Assay results are invalid and must not be reported.

IC - Internal Control; NA – not applicable; ND – not determined

^A Invalid assay run or instrument error codes or warnings are flagged on-screen and on reports. Before reporting *C. difficile* results, always verify that the assay run is valid.

Invalid Assay Run

Using frozen lysate(s), prepare new reaction tubes for all clinical specimens within that assay run along with new control tubes.

Unresolved Specimen

Repeat testing with the corresponding frozen specimen lysate. The effect of a freeze-thaw cycle has been shown to reduce the effect of PCR inhibitory substances.

Specimen Not Determined Due to I-CORE® Module Failure

Repeat testing with the corresponding frozen specimen lysate. For the interpretation of warning or error code messages, refer to the SmartCycler® Dx Software Operator Manual⁹.

Limitations of the procedure

- The performance characteristics of this assay have not been established with automated real-time PCR instruments other than the SmartCycler® instrument.
- This test is for use only with liquid or soft stools; performance characteristics of other clinical specimen types have not been established.
- Negative test results may also occur from improper specimen collection, handling or storage, presence of inhibitors, technical error, sample mix-up or because the number of organisms in the specimen is below the analytical sensitivity of the test. Careful compliance to the instructions given in this insert and in the SmartCycler® Dx Software Operator Manual⁹ is necessary to avoid erroneous results. Use of this test should be limited to personnel trained on the procedure and on the use of the SmartCycler®.
- The BD GeneOhm™ Cdiff Assay may generate unresolved or invalid results due to an invalid control; therefore, retesting of the lysate kept between -20 ± 5 °C is required and will lead to a delay in obtaining results.
- A positive test result does not necessarily indicate the presence of viable organisms. It does however indicate the presence of the *tcdB* gene and allows for presumptive detection of the *C. difficile* toxigenic organism. The BD GeneOhm™ Cdiff Assay cannot be used for species identification as it does not contain primers and probes specific to *C. difficile*.
- Mutations or polymorphisms in primer or probe binding regions may affect detection of *C. difficile tcdB* gene variants, resulting in a false negative result with the BD GeneOhm™ Cdiff Assay.
- Although there is no need for reagent preparation and the main technical operation is pipetting, good laboratory techniques are essential for the proper performance of this assay. Due to the high analytical sensitivity of this test, extreme care should be taken to preserve the purity of all reagents, especially in cases where multiple aliquots are taken from a tube.
- Variant toxigenic *C. difficile* without the *tcdB* gene or with a non-functional Toxin B protein are very rare¹²⁻¹⁵. The BD GeneOhm™ Cdiff Assay targets the *tcdB* gene and it is unknown whether it would detect Toxin A+/ Toxin B- variant strains.
- The BD GeneOhm™ Cdiff Assay has not been evaluated with the following viruses which may be present in stool specimens: the RNA viruses Enterovirus, Astrovirus, Norovirus, Rotavirus and Hepatitis A virus; the DNA virus Adenovirus. However, sequence homology between the BD GeneOhm Cdiff Assay primers and target probe relative to these viruses is very low.
- BD GeneOhm™ Cdiff Assay primers and target probe were designed to discriminate the *C. difficile* toxin B gene from the *C. sordellii* lethal toxin gene. However, cross reactivity may occur due to the high sequence homology between *C. difficile* Toxin B and *C. sordellii* lethal toxin. Due to the rare occurrence and limited availability of *C. sordellii*, only one strain was available for cross reactivity testing with the assay. No cross reactivity was observed.

Interfering Substances

Twenty-six (26) biological and chemical substances occasionally used or found in perianal, rectal and/or stool specimens were evaluated for interference with the BD GeneOhm™ Cdiff Assay. Potentially interfering substances include, but are not limited to, blood and mucus. The presence of excessive blood may inhibit PCR and may give unresolved results. The remaining twenty-four (24) substances illustrated in the table below showed no detectable interference with the BD GeneOhm™ Cdiff Assay.

Endogenous and Commercial Exogenous Substances Tested with the BD GeneOhm™ Cdiff Assay

Substance	Result	Substance	Result
Anusol ^{MC} Plus *	NI**	Monistat™ Derm Miconazole nitrate cream USP 2% (McNeil)	NI
Atlas Ihle's Paste * Zinc oxide 25 % w/w paste (Laboratoire Atlas Inc.)	NI	Palmitic acid * Fresh solution from powder (LabMat)	NI
Barium sulfate Fresh solution from powder form (LabMat)	NI	Preparation H [®] with Bio-Dyne [®] * Cream (Wyeth)	NI
Exact™ Hydrocortisone acetate * Cream USP 0.5 % (Taro Pharmaceuticals Inc.)	NI	Preparation H [®] with Bio-Dyne [®] * Ointment (Wyeth)	NI
Exact™ stomach relief Bismuth subsalicylate liquid (Perrigo [®])	NI	Rougier Neo-Laryngobis *Suppositories (Rougier Pharma)	NI
Fecal fat	NI	SAB-Dimenhydrinate [®] * Suppositories (SABEX [®])	NI
Fresh control [®] Moist towelettes pH 5,5 (Blue Skin)	NI	Steric acid * Fresh solution from powder (LabMat)	NI
Gyne Moistrin [®] Vaginal moisturizing gel (Schering)	NI	Trojan [®] latex condoms (with nonoxynol-9) Spermicidal lubricant (Church & Dwight Co., Inc.)	NI
Imodium AD [®] * Loperamide hydrochloride oral solution (McNeil)	NI	Tucks ^{MC} personal cleansing pads Moist, soft cloth pads (Pfizer)	NI
Kaopectate [®] Oral attapulgite suspension (Pharmacia & Upjohn)	NI	Vagisil [®] Anti-itch cream (Combe Incorporated)	NI
K-Y [®] Jelly (Johnson & Johnson Inc.)	NI	Vancomycin Liquid (MP Biomedicals, LLC)	NI
Metronidazole Fresh solution from powder form (Acros Organics)	NI	Vaseline™ * White petroleum jelly U.S.P. (Lever Pond's)	NI

* Substance tested with two strains of *C. difficile* (Tox 0 and Tox VIII)

** NI: No detectable interference with the BD GeneOhm™ Cdiff Assay

Performance Characteristics

Performance characteristics of the BD GeneOhm™ Cdiff Assay were determined in a multi-site prospective investigational study. Four (4) medical centers, two (2) in Canada and two (2) in the United States, participated in the study. To be enrolled in the study, specimens had to be from individuals for whom *Clostridium difficile* testing was indicated and/or ordered, according to institutional policies. The Reference Cytotoxicity Assay was performed using a tissue culture Cytotoxicity assay on liquid or soft stool specimens within 48 hours of collection. The procedure was performed according to the Manufacturer's Instructions for Use.

A total of 1108 specimens were tested with both the Reference Assay described above and the BD GeneOhm™ Cdiff Assay, producing 1090 reportable results. The first dataset includes 835 fresh specimens tested at three (3) of the four (4) clinical sites (Table 1). In comparison to the Reference Assay, the BD GeneOhm™ Cdiff Assay identified 93.8% of the *C. difficile* positive specimens and 95.5% of the negative specimens (Table 2). For the population tested this resulted in a Negative Predictive Value (NPV) of 99.1% and a Positive Predictive Value (PPV) of 67.3%. Testing at the fourth clinical site revealed that the Reference Cytotoxicity Assay was not reporting accurate results. Due to the high number of inaccurate reference assay results, samples were retested from aliquots of the original stool specimens which had been frozen after the initial testing. These frozen aliquots were tested with both the Reference Assay and the BD GeneOhm™ Cdiff Assay. The second dataset includes results from 255 frozen stool specimens available for analysis (Table 3). In comparison to the Reference Assay, the BD GeneOhm™ Cdiff Assay identified 100% of the *C. difficile* positive specimens and 97.7% of the negative specimens in the frozen dataset; resulting in a NPV of 99.2% and PPV of 81.5% (Table 4).

Out of 852 fresh specimens tested with the BD GeneOhm™ Cdiff Assay, 39 were initially reported as unresolved (4.6%). Upon repeat testing from the frozen lysates, 22 were resolved and 17 remained unresolved (2.0%) (Table 5). Out of 256 frozen specimens tested with the BD GeneOhm™ Cdiff Assay, only one (1) specimen (0.4%) was initially reported unresolved. The specimen remained unresolved upon repeat testing from the frozen lysate (0.4%) (Table 6). One (1) run was reported invalid due to Run Control failure (0.6%). The run was reported valid upon repeat testing of the specimen lysates (Table 7).

Table 1: Fresh Stool Results Obtained with the BD GeneOhm™ Cdiff Assay in Comparison with the Reference Assay

		Reference Cytotoxicity Assay		
		+	-	
BD GeneOhm™ Cdiff Assay	+	76	34 [†]	110
	-	5 [‡]	720	725
		81	754	835

[†] Cytotoxicity Assay on isolated strains was positive for 21 out of the 34 samples, verifying the presence of toxigenic *C. difficile*. For the remaining 13 samples, standard PCR with alternative primers followed by bi-directional sequencing revealed that 11 out of the 13 samples contained the expected *tcdB* gene.

[‡] For two (2) of the five (5) false negative specimens, *C. difficile* was recovered by culture, and only one (1) of these two (2) was reported as toxigenic. Of the remaining three (3) false negative PCR specimens, no *C. difficile* was recovered by culture.

Table 2: Performance Obtained with Fresh Stools using the BD GeneOhm™ Cdiff Assay in Comparison with the Reference Method

Clinical Sites	Prevalence	Sensitivity with 95% CI*	Specificity with 95% CI*
Site 1	11.0% (40/365)	90.5% (38/42) (77.4% - 97.3%)	95.7% (309/323) (92.8% - 97.6%)
Site 2	6.7% (16/240)	94.4% (17/18) (72.7% - 99.9%)	96.4% (240/249) (93.2% - 98.3%)
Site 3	11.1% (18/162)	100% (21/21) (83.9% - 100%)	94.0% (171/182) (89.4% - 96.9%)
Overall	9.6% (74/767)	93.8% (76/81) (86.2% - 98.0%)	95.5% (720/754) (93.8% - 96.9%)

* CI: Confidence Intervals

Table 3: Frozen Stool Results Obtained with the BD GeneOhm™ Cdiff Assay in Comparison with the Reference Assay

		Reference Cytotoxicity Assay		
		+	-	
BD GeneOhm™ Cdiff Assay	+	34	5	39
	-	0	216	216
		34	221	255

Table 4: Performance Obtained with Frozen Stools using the BD GeneOhm™ Cdiff Assay in Comparison with the Reference Method

Clinical Site	Prevalence	Sensitivity with 95% CI*	Specificity with 95% CI*
Site 4	12.7% (34/267)	100.0% (34/34) (89.7% - 100%)	97.7% (216/221) (94.8% - 99.3%)

* CI: Confidence Intervals

Table 5: Fresh Stool Unresolved Rates

Clinical Sites	Initial unresolved rate with 95% CI*	Unresolved rate after repeat with 95% CI*
Site 1	0.8% (3/367) (0.2% - 2.4%)	0.5% (2/367) (0.1% - 2.0%)
Site 2	6.6% (18/273) (4.0% - 10.2%)	2.2% (6/273) (0.8% - 4.7%)
Site 3	8.5% (18/212) (5.1% - 13.1%)	4.2% (9/212) (2.0% - 7.9%)
Overall	4.6% (39/852) (3.3% - 6.2%)	2.0% (17/852) (1.2% - 3.2%)

* CI: Confidence Intervals

Table 6: Frozen Stool Unresolved Rates

Clinical Site	Initial unresolved rate with 95% CI*	Unresolved rate after repeat with 95% CI*
Site 4	0.4% (1/256) (0.0% - 2.2%)	0.4% (1/256) (0.0% - 2.2%)

* CI: Confidence Intervals

Table 7: Overall Invalid Run Rates

Site	Invalid Run Rates with 95% CI*	
Site 1	2.6% (1/38)	(0.1% - 13.8%)
Site 2	0.0% (0/41)	(0.0% - 8.6%)
Site 3	0.0% (0/58)	(0.0% - 6.2%)
Site 4	0.0% (0/23)	(0.0% - 14.8%)
Overall	0.6% (1/160)	(0.0% - 3.4%)

* CI: Confidence Intervals

Analytical Specificity

Genomic DNA from one non toxigenic *C. difficile* strain, two strains of Toxinotype XI lacking *tcdB* gene¹⁷ and 29 other-*Clostridium* strains (including one strain of *C. sordellii*), along with 99 closely related organisms and other pathogenic and commensal flora found in the intestine and stools (representing 96 species) were tested. All strains were tested at a concentration of approximately 1X10⁸ CFU/mL or 1X10⁸ target copies/mL. None of these species tested positive with the BD GeneOhm™ Cdiff Assay.

Analytical Sensitivity

Quantitated culture and purified genomic DNA diluted in BD GeneOhm™ Cdiff Assay sample buffer were tested in five (5) replicates. The LOD was defined as the lowest concentration, in DNA copy number per reaction and CFU per reaction, at which five replicates out of five were found positive.

The analytical sensitivity (limit of detection or LOD) of the BD GeneOhm™ Cdiff Assay was determined with one strain of Toxinotype 0 *Clostridium difficile* carrying the *tcdB* gene (ATCC 43255).

The BD GeneOhm™ Cdiff Assay LOD is 10 DNA copies per reaction. The LOD in Colony Forming Units (CFU) is established at 4 CFU per reaction.

The analytical sensitivity in CFU per reaction was confirmed with a second Toxinotype 0 (ATCC 9689) and with Toxinotypes IIIa (SE844¹⁶), V (SE881¹⁶), VII (57267¹⁶) and VIII (1470¹⁶) *Clostridium difficile* toxigenic strains.

In addition to strains used for LOD determination, one hundred (100) other toxigenic *C. difficile* strains (including 17 other Toxinotypes), representing 21 countries, from well-characterized clinical isolates or public collections were evaluated using the BD GeneOhm™ Cdiff Assay. *C. difficile* strains were tested at a concentration of approximately 6.7 DNA copies/μL or 1 CFU/μL. The assay correctly identified all 100 *C. difficile* strains carrying the *tcdB* gene.

Reproducibility

The reproducibility panel consisted of three (3) simulated specimen categories where each tube contained 100 μL of simulated bowel flora; the two positive panel members were also inoculated with *C. difficile* (ATCC 43255). Additionally, two (2) Specimen Processing Controls (ATCC 9689 and ATCC 25922) and, two (2) Run Controls (Positive and Negative) were included. The specimens were tested in triplicate per panel run, on five (5) distinct days (consecutive or not), wherein each day two (2) panels were tested, one for each of two (2) technologists, at three (3) clinical sites with one (1) lot of reagents. One (1) of these clinical sites participated in the extended study where two (2) additional lots of reagents were tested.

The overall percent agreement for the low positive *C. difficile* specimen category is 96.7%; the moderate positive *C. difficile* specimen category is 100% and the negative specimen category is 100% for the Site-to-Site Reproducibility (Table 8).

The overall percent agreement for the low positive *C. difficile* specimen category is 100%; the moderate positive *C. difficile* specimen category is 97.8% and the negative specimen category is 100% for the Lot-to-Lot Reproducibility (Table 9).

Cycle threshold (Ct), an internal criteria used to determine a final assay result, was selected as an additional means of assessing assay reproducibility. Overall mean Ct values with variance components (SD and %CV) are shown in Tables 8 and 9.

An additional reproducibility study was performed, in accordance with the original reproducibility study protocol, to assess high negative specimens below the BD GeneOhm™ Cdiff Assay limit of detection (LOD). A sample containing simulated bowel flora was inoculated with *C. difficile* (ATCC 43255) at a concentration equivalent to the assay LOD. 100-fold and 10-fold dilutions of this sample were prepared, respectively, to obtain the two (2) high negative panel members. Overall percent agreement for negative test results and overall mean Ct values with variance components (SD and %CV) are shown in Table 10. As expected, the more dilute panel member (100-fold below the LOD) containing lower levels of target, demonstrates a higher percent agreement for negative test results than the less dilute panel member (10-fold below the LOD) which contains higher levels of target. Although high negative panel members are below the analytical LOD of the assay, positive test results may still be observed due to the presence of target in these specimens.

Table 8: Site-To-Site Reproducibility Study Results using One Lot

Category	SITE						Overall Percent Agreement		Ct Values		
	Site 1		Site 2		Site 3						
	Percent Agreement	Percent Agreement	Percent Agreement	Percent Agreement	Percent Agreement	Percent Agreement	Overall Mean	SD	%CV		
NEG	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%	36.2 [†]	0.3 [†]	0.8% [†]
LOW POS	28/30	93.3%	29/30	96.7%	30/30	100.0%	87/90	96.7%	38.8	0.9	2.3%
MOD POS	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%	38.3	1.0	2.7%

[†]Data represent values from the internal control.

Table 9: Lot-To-Lot Reproducibility Study Results using Three Lots

Category	LOT						Overall Percent Agreement		Ct Values		
	Lot 1		Lot 2		Lot 3						
	Percent Agreement	Percent Agreement	Percent Agreement	Percent Agreement	Percent Agreement	Percent Agreement	Overall Mean	SD	%CV		
NEG	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%	36.1 [†]	0.3 [†]	0.8% [†]
LOW POS	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%	38.6	1.0	2.5%
MOD POS	29/30	96.7%	29/30	96.7%	30/30	100.0%	88/90	97.8%	37.8	1.1	2.8%

[†]Data represent values from the internal control.

Table 10: Additional Reproducibility Study Using a High Negative Sample Panel

High Negative Panel Member	Site 1		Site 2		Site 3		Overall Percent Agreement*		Ct Values		
	Percent Agreement*	Percent Agreement*	Percent Agreement*	Percent Agreement*	Percent Agreement*	Percent Agreement*					
1:100 dilution	25/30	83.3%	21/30	70.0%	26/30	86.7%	72/90	80.0%	41.3	0.9	2.1%
1:10 dilution	11/30	36.7%	5/30	16.7%	5/30	16.7%	21/90	23.3%	40.2	1.4	3.4%

*Percent agreement for a negative result.

Precision

Within-laboratory precision was evaluated for the BD GeneOhm Cdiff Assay at one (1) site. The study was performed over 12 days, with two (2) runs per day and two (2) sample replicates per run. Samples included simulated specimens representing low and moderate positive *C. difficile* as well as negative *C. difficile*. One (1) out of 24 runs was excluded due to failure of the positive control (PC). One (1) moderate positive sample produced an unresolved result. All remaining samples and controls produced reportable results for a total of 46 replicates. Precision study results for low and moderate positive samples demonstrated agreement for (46/46) and (45/46) replicates, respectively; negative sample results demonstrated agreement for (46/46) replicates.

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1. Cloud J and Kelly CP., Update on *Clostridium difficile* associated disease. *Curr Opin Gastroenterol* 2007 Jan; 23 (1):4-9.
2. Warny M, *et al.* Toxin production by an emerging strain of *Clostridium difficile* associated with outbreaks of severe disease in North America and Europe. *Lancet*. 2005 Sep 24-30; 366(9491):1079-84.
3. Bartlett, John G. MD, *Clostridium difficile*: Old and New Observations. *J of Clin Gastroent* 41 Supplement 1:S24-S29, May/June 2007.
4. Hurley BW and Nguyen CC. The spectrum of pseudomembranous enterocolitis and antibiotic-associated diarrhea. *Arch Intern Med*. 2002 Oct 28; 162 (19):2177-84.
5. Polgreen PM, *et al.* An outbreak of severe *Clostridium difficile*-associated disease possibly related to inappropriate antimicrobial therapy for community-acquired pneumonia. *Infect Control Hosp Epidemiol*. 2007 Feb; 28 (2):212-4.
6. Nagy E, *et al.* The place of molecular genetic methods in the diagnostics of human pathogenic anaerobic bacteria. A minireview. *Acta Microbiol Immunol Hung*. 2006 Jun; 53 (2):183-94.
7. Centers for Disease Control and Prevention. Biosafety in microbiological and biomedical laboratories. Richmond JY and McKinney RW (eds) (1993). HHS Publication number (CDC) 93-8395.
8. Clinical and Laboratory Standards Institute. Protection of laboratory workers from occupationally acquired infections; Approved Guideline. Document M29 (Refer to the latest edition). Clinical Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2006.
9. SmartCycler® Dx Software Operator Manual D4198 Rev. C, Cepheid, Sunnyvale, CA, USA.
10. Clinical and Laboratory Standards Institute. Molecular Diagnostic Methods for Infectious Diseases; Approved Guideline, document MM3-A2 (Refer to the latest edition). Clinical Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2006.
11. Clinical and Laboratory Standards Institute Statistical Quality Control for Quantitative Measurements: Principles and Definitions; Approved Guideline. Document C24 (Refer to the latest edition). Clinical Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2006.
12. Cohen SH *et al.* (1998) Isolation of a Toxin B-deficient mutant strain of *Clostridium difficile* in a case of recurrent *C. difficile* – Associated Diarrhea, *Clin. Infect. Dis*. 26: 410-412.
13. Cohen SH *et al.* (2000) Analysis of the pathogenicity locus in *Clostridium difficile* strains, *J. Infect. Dis*. 181: 659-63.
14. Maccannell D. *et al.* (2006) Characterization of a novel, TcdB-deficient, NPA1 variant strain of *Clostridium difficile*, 46th Annual ICAAC, San Francisco, Sept. 2006.
15. McFarland, L.V., *et al.*, Implications of the changing face of *Clostridium difficile* disease for health care practitioners. *Am J Infect Control*, 2007. 35(4): p. 237-53.
16. Rupnik M. *et al.* A novel toxinotyping scheme and correlation of toxinotypes with serogroups of *Clostridium difficile* isolates, *J Clin Microbiol*, 1998, vol 36 no 8 : 2240-2247 .
17. Rupnik M. *et al.* Comparison of toxinotyping and PCR ribotyping of *Clostridium difficile* strains and description of novel toxinotypes, *Microbiology* 2001, 147, 439-447.

This kit is sold under license from the Public Health Research Institute of the City of New York, Inc. and may be used under the PHRI patent rights only for human *in vitro* clinical diagnostics.

The purchase of this product allows the purchaser to use it for amplification and detection of nucleic acid sequences for providing human *in vitro* diagnostics. No general patent or other license of any kind other than this specific right of use from purchase is granted hereby.


This product is subject to an agreement between Molecular Probes, Inc. and GeneOhm Sciences, Inc., and the manufacture, use, sale or import of this product may be subject to one or more of U.S. patents, pending applications and corresponding foreign equivalents, owned by Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corporation). The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product, for use in assays for nucleic acid detection for the purpose of identifying microorganisms in human diagnostics. The buyer cannot use this product or its components for manufacturing or for therapeutic or prophylactic use, or sell or otherwise transfer this product or its components to any third party, or use for any use other than use in human diagnostics. For information on purchasing a license to this product for purposes other than use in human diagnostics, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, USA, Tel: (541) 465-8300. Fax: (541) 335-0504.

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
Australian Representative:
Becton Dickinson Pty Ltd.
4 Research Park Drive,
Macquarie University Research Park,
North Ryde,
NSW 2113 Australia