**BD GeneOhm™ Lysis Kit**  
**Trousse de lyse BD GeneOhm™**  
**BD GeneOhm™-Lyseatz**  
**Equipo BD GeneOhm™ *Lysis***  
**Kit di lisi BD GeneOhm™**

**REF 441243**

**IVD**

## Table of Contents / Table des matières / Inhaltstabelle / Índice / Indice

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## English

### Intended use

BD GeneOhm™ Lysis kit is for the rapid lysis of cells and spores from different sample types.

### Reagents

#### BD GeneOhm™ Lysis Kit

**100 reactions**

#### *Sample Buffer*

**120 X 1 mL**

Tris-EDTA buffer

#### *Lysis tube*

**100 tubes**

Glass beads

### Precautions

- Do not use the kit if the outer carton safety seal is broken.
- Do not use reagents if the protective pouches are open or torn upon arrival.
- Close protective pouches with the zip seal after each use.
- Do not use the reagents after their expiration date.
- Do not interchange caps among reagents as contamination may occur and compromise test results.
- The use of sterile disposable filter-blocked or positive displacement pipettor tips is recommended.
- Use a new tip for each sample or reagent.
- Always handle specimens as if infectious in accordance with safe laboratory procedures such as those described in *Biosafety in Microbiological and Biomedical Laboratories*<sup>1</sup> and in the CLSI Document M29<sup>2</sup>.
- Wear protective clothing and disposable gloves while handling kit reagents. Wash hands thoroughly after use.
- Do not pipet by mouth.
- Do not smoke, drink, or eat in areas where specimens are being handled.
- Dispose of unused reagents and waste in accordance with country, federal, provincial, state and local regulations.

### Materials provided

- Sample Buffer
- Lysis tube

### Storage, handling and stability

Kit Component		Sample Buffer (blue cap)	Lysis Tube (yellow cap)
Sealed pouch	Temperature	2-25 °C	2-25 °C
	Stability	Expiration date	Expiration date
Opened pouch	Temperature	2-25 °C	2-25 °C
	Stability	2 months <sup>1</sup>	Expiration date

<sup>1</sup> Provided that the pouch is properly closed with the zip seal after each use.

## Materials required but not provided

- **Vortex** Genie 2 (Fisher) with 1.5 mL microtube holder or equivalent; for processing multiple samples, adapter with multiple holding sites can be used
- **Micropipettors** (accurate range between 10-100 µL and 100-1000 µL)
- **Sterile** filter-blocked or positive displacement **micropipettor tips**
- **Sterile transfer pipets** (optional)
- **Sterile tubes** (optional)
- **Sterile Saline (0.85% NaCl)** (optional)
- **Scissors** (optional)
- **Gauzes** (optional)
- Disposable **gloves**
- **Microcentrifuge** for high (must reach 14 000 x g) (optional) and low speed centrifugation
- **Dry heating block** for 1.5 mL tubes or water bath
- Ice or **cooling block** for 1.5 mL tubes
- Stopwatch or **timer**

## Instructions for use

**Note:** One Lysis tube (**yellow cap**) and at least one Sample Buffer tube (**blue cap**) are required for each sample to be prepared depending on the procedure. Remove the required number of tubes from their protective pouch, and close the pouches with the zip seal.

(Preparation time approx. 15 min)

### A- Specimen preparation for lysis

#### 1) The following procedure is suggested for isolated colonies:

- i. Resuspend isolated colonies in sterile saline (0.85% NaCl) to a turbidity of 0.5 McFarland.
- ii. Proceed with point C if specimen dilution is needed. If not, proceed with point D.

#### 2) The following procedure is suggested for sample collected on a swab:

- i. Place the collection device (swab) in a Sample Buffer tube (blue cap).

Identify the Sample Buffer tube on the cap and/or the tube label.

- ii. Break the swab stem and close the tube tightly.

Hold the swab by the stem near the rim of the tube (use gauze to minimize risks of contamination). Lift the swab a few millimeters from the bottom and push it against the edge of the tube to break it. Alternative method: use clean scissors to cut the stem. Make sure the cap will close tightly.

- iii. Vortex at high speed for one minute.

Time may vary with the type of sample. For processing multiple samples, an adapter with multiple holding sites can be used.

- iv. Proceed with point B to concentrate the specimen. Proceed with point C to dilute the specimen.
- v. Proceed with point E (Lysis Method) if step iv is not required.

#### 3) The following procedure is suggested for other types of sample:

- i. Proceed with point B for specimen concentration, if needed.
- ii. Proceed with point C for specimen dilution, if needed.
- iii. Proceed with point D for specimen washing, if needed.
- iv. Proceed directly with point E (Lysis Method) if the type of sample requires no treatment prior to lysis.

**B- Concentration Method**

- 1) **Transfer** a maximum of 1000 µL of sample to the Lysis tube, **close tightly** and **centrifuge**, at high speed (**between 14 000 x g and 21 000 x g**) for **5 minutes** at room temperature.
- 2) **Remove the supernatant and discard it.**  
Use a sterile transfer pipet; take care not to touch the pellet. Use a new transfer pipet for each sample.
- 3) **If necessary, repeat steps 1 and 2 as many times as needed. If not, proceed with step 4.**
- 4) **Add 50 µL to 100 µL of sample buffer to the Lysis tube; close tightly. If a washing step is necessary, proceed with point D.**
- 5) **Continue with step 2 of the Lysis Method (E).**

**C- Dilution Method**

- 1) **Use the sample buffer to dilute a concentrated sample.**
- 2) **Continue with Step 1 of the Lysis Method (E).**

**D- Washing Method**

- 1) **Centrifuge** the Lysis tube containing a maximum of 1000 µL of sample, at high speed (**between 14 000 x g and 21 000 x g**) for **5 minutes** at room temperature.
- 2) **Remove the supernatant and discard it.**  
Use a sterile transfer pipet; take care not to touch the pellet. Use a new transfer pipet for each sample.
- 3) **Add 50 µL to 100 µL of sample buffer to the lysis tube; close tightly.**  
Use a new pipettor tip for each sample.
- 4) **Repeat steps 1 to 3 as many times as needed.**
- 5) **Continue with step 2 of the Lysis Method (E).**

**E- Lysis Method**

- 1) **Add 50 µL to 100 µL of sample to the Lysis tube (yellow cap); close tightly.**  
Use appropriate micropipettor and a new pipettor tip for each sample. Identify the Lysis tube on the cap and/or the tube label. Do not add more than 100 µL to the Lysis tube.
- 2) **Vortex at high speed for 5 minutes.**  
Time may vary with the type of sample to be lysed. For processing multiple samples, an adapter with multiple holding sites can be used.
- 3) **Centrifuge the Lysis tube at low speed for 2 to 5 seconds, to bring the content at the bottom of the tube.**
- 4) **Heat at 95 ± 2 °C for two (2) minutes.**  
Use a dry heating block or a water bath.
- 5) **Keep the Lysis tube on ice or on a cooling block.**

**Notes:**

- 1- **Dilution of lysate:** To dilute a lysate, add the appropriate volume of sample buffer to the Lysis tube and store on ice or on a cooling block.
- 2- **Storage of lysate:** The lysate could be frozen at -20 ± 5 °C for later use, if necessary. Freeze-thaw cycle may also eliminate PCR inhibitory substances.

## Results

Table 1. Expected yields

Species	Type of Sample	Amount (CFU <sup>1</sup> /mL)	Lysis efficiency (%)
<i>Streptococcus agalactiae</i>	Vaginal/anal swab	Between 10 <sup>5</sup> and 10 <sup>6</sup>	98.8
<i>Staphylococcus aureus</i>	Bacterial suspension		>99
<i>Bacillus anthracis</i> (vegetative cells)	Nasal swab		>99
<i>Bacillus anthracis</i> (Spores)	Silica powder (0.24 mg)		>99
<i>Bacillus subtilis</i> (Spores)	Spore suspension		98.8

<sup>1</sup> CFU : Colony Forming Unit

Lysis achieved with the BD GeneOhm™ Lysis kit results in the release of intracellular components as demonstrated by the amplification of genomic nucleic acids present in resulting lysates<sup>3, 4, 5, 6, 7, 8, 9</sup>.

Table 2. CFU detected by PCR reaction

Species	Type of sample	Calculation				
		Initial Concentration	Volume added to the Lysis tube (µL)	Dilution of lysate	Volume added to reaction tube (µL)	Final concentration
<i>S. saprophyticus</i> <sup>3</sup>	<i>Urine</i>	1.4x10 <sup>3</sup> CFU/mL	50	N/Ap	1	≤ 1.4 CFU/rx
<i>C. difficile</i> <sup>4</sup>	<i>Feces</i>	10 <sup>7</sup> CFU/g of feces	75 (liquid stool)	N/Ap	1.5	10 <sup>4</sup> CFU/g of feces
<i>M. smegmatis</i>	<i>Suspension</i>	2.4 x10 <sup>5</sup> CFU/mL	100	1/2500	1	9.6 CFU/rx
<i>S. epidermidis</i>	<i>Suspension</i>	2.4 x10 <sup>8</sup> CFU/mL	100	1/10000	1	24 CFU/rx
<i>E. coli</i> <sup>8</sup>	<i>Feces</i>	10 <sup>8</sup> CFU/g of feces	100 (liquid stool)	N/Ap	1.5	10 <sup>5</sup> CFU/g of feces
<i>E. coli</i>	<i>Suspension</i>	10 <sup>8</sup> CFU/mL	100	1/100000	1	1 CFU/rx
<i>E. faecium</i>	<i>Suspension</i>	8 x10 <sup>8</sup> CFU/mL	100	1/100000	1	≤ 8 CFU/rx
<i>C. albicans</i>	<i>Suspension</i>	10 <sup>8</sup> CFU/mL	10	1/10000	1	≤ 10 CFU/rx

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
rx = reazione

## References / Références / Referenzen / Referencias / Riferimenti

- 1 Centers for disease Control and prevention. Biosafety in Microbiological and Biomedical Laboratories. Richmond JY and McKinney RW (eds) (1993). HHS Publication number (CDC) 93-8395.
- 2 Clinical and Laboratory Standards Institute. Protection of laboratory workers from occupationally acquired infections; Approved Guidelines – Second edition. Document M29 (Refer to the latest edition).
- 3 Martineau, F. *et al.* Development of a rapid PCR Assay Specific for *Staphylococcus saprophyticus* and Application to Direct Detection from Urine Samples. Journal of Clinical Microbiology. Sept 2000; 38(9): 3280-3284.
- 4 Bélanger, S D. *et al.* Rapid Detection of *Clostridium difficile* in Feces by Real-Time PCR. Journal of Clinical Microbiology. Feb. 2003; 41(2): 730-734.
- 5 BD Diagnostics, Sainte-Foy, Québec, Canada. Package insert of **BD GeneOhm™ MRSA Assay (cat. # 441242** for 200 tests kit and **441244** for 48 tests kit).
- 6 BD Diagnostics, Sainte-Foy, Québec, Canada. Package insert of **BD GeneOhm™ Strep B Assay (cat. # 441240)**.
- 7 Comparison of nine commercial kits for rapid nucleic acid extraction from microbial cultures. (M. Dion, C. Ménard, F.J. Picard *et al.* Abstr. 99th Annual Meeting of ASM, Chicago, abstr. C-481, 1999).
- 8 Bélanger, S D. *et al.* Rapid Detection of Shiga Toxin-Producing Bacteria in Feces by Multiplex PCR with Molecular Beacon on the SmartCycler. Journal of Clinical Microbiology. Apr. 2002; 40(4) : 1436-1440.
- 9 Martineau, F. *et al.* Species-Specific and ubiquitous DNA-Based Assays for Rapid Identification of *Staphylococcus epidermidis* Journal of Clinical Microbiology. Dec. 1996; 34(12): 2888-2893.

## Index of symbols / Table des symboles / Symbolindex / Índice de símbolos / Indice dei simboli

Symbol / Symbole / Símbolo / Simbolo	Meaning / Signification / Bedeutung / Significado / Significato
REF	Catalog number (US) / Catalogue number (UK) / Référence du Catalogue / Katalognummer / Número de Catálogo / Numero di catalogo
	In Vitro Diagnostic use / Aux fins de Diagnostic In Vitro / In-vitro-Diagnostika / Para fines de Diagnóstico in Vitro / Uso diagnostico In Vitro
	Manufacturer / Fabricant / Hersteller / Fabricante / Produttore
	Authorized european representative / Représentant européen autorisé / Bevollmächtigter in der Europäischen Gemeinschaft / Representante europeo autorizado / Rappresentante europeo autorizzato
	Contains sufficient for "n" tests / Contenu suffisant pour « n » tests / Inhalt ausreichend für « n » Prüfungen / Contiene cantidad suficiente para « n » pruebas / Il contenuto è sufficiente per "n" test
	Batch code / Code de lot / Chargenbezeichnung / Código de lote / Codice del lotto
	Use by / Utiliser avant / Verwendbar bis / Fecha de caducidad / Data di scadenza
	Temperature limitation / Limites de température / Temperaturbegrenzung / Límites de temperatura / Limiti di temperatura
	Reseal bag after use / Refermer le sac après utilisation / Beutel nach Gebrauch wieder verschließen / Volver a cerrar la bolsa después de la utilización / Richiudere la sacca dopo l'uso
	Consult instructions for use / Se référer aux instructions d'utilisation / Gebrauchsanweisung beachten / Ver las instrucciones de utilización / Consultare le istruzioni per l'uso

Customer Service/Service à la Clientèle/Kundendienst/ Servicio de Atención al Cliente/ Servizio Assistenza:  
1.888.436.3646

BENEX Limited, Bay K 1a/d,  
Shannon Industrial, Estate,  
Shannon, Country Clare,  
Ireland