
Revisions

BALTSO0191 Version 11.0 Template 4
Inserts

Rev from	Rev to	JOB #
01	02	9509-17

Notes:

- BD Catalog Number: 212536, 212537
- Blank (Sheet) Size: Length: 16 7/8 " Width: 8.5"
- Number of Pages: 12 Number of Sheets: 1
- Page Size: Length: 8.5" Width: 5.5" Final Folded Size: 4.25" x 5.5"
- Ink Colors: No. of Colors: 2 PMS#: 032 Red, Standard Black
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 7
One vertical fold and two horizontal folds to form a final size of 4 1/4" x 5 1/2" showing English header

- Vendor Printed ☒ Online/In House Printed ☐ Web ☐
- See specification control no. N/A for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Nancy Carlsen at 3:39 pm, Sep 19, 2017	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY By Natalie Morio at 7:17 am, Sep 20, 2017		
Checked By	THIRD EYE BY By Kim Stone at 8:23 am, Sep 20, 2017		
Part Number:	8820211JAA	Category and Description Package Insert, Acridine Orange Stain	Sheet: 1 of 2 Scale: N/A A

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BD BBL™ Acridine Orange Stain	For detecting microorganisms in direct smears by the fluorescent staining technique		Cat. No.
			212536
		1 x 250 mL	212536
		4 x 250 mL	212537

INTENDED USE

BD BBL Acridine Orange Stain is recommended for use in the fluorescent microscopic detection of microorganisms in direct smears prepared from clinical and non-clinical materials. It is particularly useful in the rapid screening of normally sterile specimens, such as cerebrospinal fluid, where few organisms may be present and in the rapid examination of blood smears or smears containing proteinaceous material where differentiation of organisms from background material may be more difficult.

SUMMARY AND EXPLANATION

Fluorochromatic staining of microorganisms using acridine orange was first described by Strugger and Hilbrich in 1942. It has since been widely used in the examination of soil and water for microbial content. In 1975, Jones and Simon evaluated epifluorescent methods used in direct counts of aquatic bacteria and determined that acridine orange provided the best estimate of the bacterial population in lake, river and seawater samples.¹ Acridine orange direct count (AODC) methodology has been used in the enumeration of landfill bacteria.^{2,3} Heidelberg et al. used AODC in a study of seasonal changes in marine bacterial populations and concluded that the acridine orange stain compared favorably to fluorescent oligonucleotide direct counting (FODC) procedures.⁴ Direct epifluorescent filter technique (DEFT) using acridine orange is specified in methods for the microbial examination of food and water.^{5,6,7,8}

Acridine orange has also been used in clinical applications and its use in highlighting bacteria in blood cultures has become widely accepted. In 1980, McCarthy and Senne compared acridine orange staining with blind subcultures for the detection of positive blood cultures.⁹ Their results showed acridine orange staining to be a rapid, inexpensive alternative to blind subcultures. They also reported that the acridine orange stain appeared to be more sensitive than the Gram stain for detecting microorganisms and was able to detect bacteria in concentrations of approximately 1 x 10⁴ CFU (colony-forming units)/mL. Lauer, Reller and Mirret compared acridine orange with the Gram stain for detecting microorganisms in cerebrospinal fluid and other clinical materials.¹⁰ Their results were in agreement with those reported by McCarthy and Senne and showed acridine orange to be a simple, rapid staining procedure which was more sensitive than the Gram stain in detecting microorganisms in clinical materials.

Acridine orange has also been used for the detection of *Trichomonas vaginalis* in vaginal smears,¹¹ diagnosis of malaria^{12,13} and mycoplasmas.¹⁴

PRINCIPLES OF THE PROCEDURE

Acridine orange is a fluorochromatic dye which binds to nucleic acids of bacteria and other cells.¹⁵ Under UV light, acridine orange stains RNA and single-stranded DNA orange; double-stranded DNA appears green.

When buffered at pH 3.5–4.0, acridine orange differentially stains microorganisms from cellular materials. Bacteria and fungi uniformly stain bright orange, whereas human epithelial and inflammatory cells and background debris stain pale green to yellow. Nuclei of activated leukocytes stain yellow, orange or red due to increased RNA production resulting from activation. Erythrocytes either do not stain or appear pale green.

Due to this differential staining property, acridine orange-stained smears prepared from clinical materials may be rapidly screened using fluorescent microscopy at 100X to 400X magnification for the presence of microorganisms fluorescing bright orange against a black or pale green to yellow background.

REAGENTS

BD BBL Acridine Orange Stain

Approximate Formula*

Acridine Orange 0.1 g
Acetate Buffer, 0.5 M 1,000 mL

*Adjusted and/or supplemented as required to meet performance criteria.

Warnings and Precautions: For *in vitro* Diagnostic Use.

Follow proper, established laboratory procedures in handling and disposing of infectious materials.

BD BBL Acridine Orange Stain:

WARNING

H315+H320 Causes skin and eye irritation. **H335** May cause respiratory irritation.

P261 Avoid breathing dust/fume/gas/mist/vapors/spray. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P264** Wash thoroughly after handling. **P271** Use only outdoors or in a well-ventilated area. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P321** Specific treatment (see on this label). **P304+P340** IF INHALED: Remove person to fresh air and keep comfortable for breathing. **P312** Call a POISON CENTER or doctor/physician if you feel unwell. **P332+P313** If skin irritation occurs: Get medical advice/attention. **P337+P313** If eye irritation persists: Get medical advice/attention. **P302+P352** IF ON SKIN: Wash with plenty of soap and water. **P362+P364** Take off contaminated clothing and wash it before reuse. **P405** Store locked up. **P403+P233** Store in a well-ventilated place. Keep container tightly closed. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage Instructions: Store at 15–30 °C. The expiration date applies to the product in its intact container when stored as directed.

Product Deterioration: Do not use if there is evidence of a precipitate or the solution shows other signs of deterioration.

SPECIMEN COLLECTION AND HANDLING

Specimens should be collected in sterile containers or with sterile swabs and transported immediately to the laboratory in accordance with recommended guidelines.¹⁶

PROCEDURE

Materials Provided: BD BBL Acridine Orange Stain.

Materials Required But Not Provided: Fluorescent microscope suitable for use with acridine orange, glass microscope slides and methanol.

Preparation, Staining and Examination of Smears

1. Prepare a smear of the specimen to be stained on a clean glass slide.
2. Allow to air dry.
3. Fix smear with 50% or 100% methanol for 1–2 min.
4. Drain excess methanol and allow smear to dry.
5. Flood slide with acridine orange stain for 2 min.
6. Rinse thoroughly with tap water and allow to dry.
7. Smears may be initially examined at 100X to 400X magnification using a fluorescent microscope. Findings should be confirmed by examination at 1000X with an oil immersion objective.

USER QUALITY CONTROL

1. Examine the acridine orange staining solution for color and clarity. The solution should be clear, orange and without evidence of a precipitate.
2. Determine the pH of the solution. The pH should be 3.5–4.0.
3. Check the performance of the stain using 4–6 h Tryptic Soy Broth with 5% sheep blood cultures of the organisms indicated below. Prepare smears, one culture per slide, and proceed as described under Preparation, Staining and Examination of Smears.

Organisms	Bacteria	Background
<i>Escherichia coli</i> ATCC® 25922	Bright orange	Pale green erythrocytes and yellow, yellow-green or orange leukocytes against a black field.
<i>Enterococcus faecalis</i> ATCC 33186	Bright orange	Green, yellow, orange or red staining debris may be observed.

LIMITATIONS OF THE PROCEDURE

Acridine orange staining provides presumptive information on the presence and identification of microorganisms which may be present on the specimen. Since microorganisms seen in smears, including non-viables, may arise from external sources, i.e., specimen collection devices, slides, or water used for rinsing, all positive smears should be confirmed by culture.

Approximately 10⁴ CFU/mL are required for detection by this method.

Acridine orange staining does not distinguish between gram-positive and gram-negative organisms. The Gram reaction may be determined by Gram staining directly over the acridine orange after removal of the immersion oil.

Nuclei or granules from disintegrated activated leukocytes may resemble cocci at lower magnifications, i.e., 100X to 400X. They may be distinguished on the basis of morphology at higher magnifications, i.e., 1000X.

Certain types of debris may fluoresce in acridine orange stained smears. This debris may be distinguished from microorganisms on the basis of morphology when viewed at higher magnifications.

EXPECTED RESULTS AND PERFORMANCE CHARACTERISTICS

Bacteria and fungi stain bright orange. The background appears black to yellow green. Human epithelial and inflammatory cells and tissue debris stain pale green to yellow. Activated leukocytes will stain yellow, orange or red depending upon the level of activation and amount of RNA produced, whereas erythrocytes either do not stain or stain pale green.

REFERENCES

1. Jones, J.G. and B.M. Simon. 1975. An investigation of errors in direct counts of aquatic bacteria by epifluorescence microscopy, with reference to a new method for dyeing membrane filters. J. Appl. Bacteriol. 39: 317-329.

2. Barlaz, M.A. 1997. Microbial studies of landfills and anaerobic refuse decomposition, p. 541-557. In C.J. Hurst, G.R. Knudsen, M.J. McInerney, L.D. Stetzenbach, and M.V. Walter (ed.), Manual of environmental microbiology, American Society for Microbiology, Washington, D.C.

3. Palmisano, A.C., D.A. Maurusick, and B.S. Schwab. 1993. Enumeration of fermentative and hydrolytic micro-organisms from three sanitary landfills. J. Gen. Microbiol. 139:387-391.

4. Heidelberg, J.F., K.B. Heidelberg, and R.R. Colwell. 2002. Seasonality of Chesapeake Bay bacterioplankton species. Appl. Environ. Microbiol. 68:5488-5497.

5. Splittstoesser, D.F. 1992. Direct microscopic count, p. 97-104. In C.V. Vanderzant and D.F. Splittstoesser (ed.), Compendium of methods for the microbiological examination of foods, 3rd ed., American Public Health Association, Washington, D.C.

6. Packard, V.S., Jr., S. Tatini, R. Fugua, J. Heady, and C. Gilman. 1992. Direct microscopic methods for bacteria or somatic cells, p. 309-325. In R.T. Marshall (ed.), Standard methods for the examination of dairy products, 16th ed. American Public Health Association, Washington, D.C.

7. Duffy, G., Kilbride, B., Fitzmaurice, J., Sheridan, J.J. 2001. Routine diagnostic tests for food-borne pathogens. The National Food Centre, Dublin.

8. Eaton, A.D., L.S. Clesceri, and A.E. Greenberg (ed.). 1995. Standard methods for the examination of water and wastewater, 19th ed. American Public Health Association, Washington, D.C.

9. McCarthy, L.R. and J.E. Senne. 1980. Evaluation of acridine orange stain for detection of microorganisms in blood cultures. J. Clin. Microbiol. 11:281-285.

10. Lauer, B.A., L.B. Reller, and S. Mirrett. 1981. Comparison of acridine orange and Gram stains for detection of microorganisms in cerebrospinal fluid and other clinical specimens. J. Clin. Microbiol. 14:201-205.

11. Greenwood, J.R., and K. Kirk-Hillaire. 1981. Evaluation of acridine orange stain for detection of *Trichomonas vaginalis* in vaginal specimens. J. Clin. Microbiol. 14:699.

12. Keiser, J., J. Utzinger, Z. Premji, Y. Yamagata, and B.H. Singer. 2002. Acridine orange for malaria diagnosis: its diagnostic performance, its promotion and implementation in Tanzania, and the implications for malaria control. Ann. Trop. Med. Parasitol. 96:643-654.

13. Bosch, I., C. Bracho, and H.A. Perez. 1996. Diagnosis of malaria by acridine orange fluorescent microscopy in an endemic area of Venezuela. Mem. Inst. Oswaldo Cruz. 91:83-86.

14. Rosendal, S. and A. Valdivieso-Garcia. 1981. Enumeration of mycoplasmas after acridine orange staining. Appl. Environ. Microbiol. 41:1000-1002.

15. Kasten, F.H. 1967. Cytochemical studies with acridine orange and the influence of dye contaminants in the staining of nucleic acids. Internat. Rev. Cytol. 21:141-202.

16. Shea, Y.R. 1994. Specimen collection and transport, p. 1.1.1-1.1.30. In H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1, American Society for Microbiology, Washington, D.C.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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R_x Only
This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Australian Sponsor:
Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

ATCC is a trademark of the American Type Culture Collection.
© 2017 BD. BD, the BD Logo and BBL are trademarks of Becton, Dickinson and Company.