
Rev from	Rev to	JOB #
01	02	7278-15

NOTES:

- BD Catalog Number: 260640
- Blank (Sheet) Size: Length: 11" Width: 8.5" (Cover: Horiz.: 21.75"; Vert.: 8.5")
- Number of Pages: 20 Number of Sheets: 4 (Including 4-panel, 8-page cover)
- Page Size: Length: 8.5" Width: 5.5" Final Folded Size: 5.5" x 8.5"
- Ink Colors: No. of Colors: 2 PMS#: 300 Blue; Black; (Note to printer: 100% magenta copy, color labeled Non-printing, in slug area of file does not print.)
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 6

- See specification control no. VS 8801351JAA for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Nancy Carlsen at 4:01 pm, Mar 12, 2015	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY By Helaine Jones at 5:16 pm, Mar 12, 2015		
Checked By	THIRD EYE BY By Kim Stone at 8:23 am, Mar 13, 2015		
Part Number:	8801351JAA	Category and Description Package Insert, BBL Sensi-Disc Designer Dispenser	Sheet: 1 of 21 Scale: N/A A

BD BBL™ Sensi-Disc™ Designer Dispenser

**Designer Dispenser BD BBL Sensi-Disc
BD BBL Sensi-Disc Designer-Dispenser
BD BBL Sensi-Disc Designer-Dispenser Dispensador
BD BBL Sensi-Disc Designer**

REF 260640
8801351JAA(02)
2015-03

**one unit, 12 place, complete
(with protective base and cover and one desiccant canister)**

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ANATOMY OF THE BBL™ SENSI-DISC™ DESIGNER DISPENSER

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ANATOMÍA DEL DISPENSADOR DE DISCOS BBL SENSI-DISC DESIGNER

Lock latch with finger grips permits easy movement of locking mechanism.

Loquet de verrouillage avec poignée permettant des mouvements aisés du mécanisme de verrouillage / Verschlussriegel mit Griffflächen erleichtert die Betätigung des Verriegelungsmechanismus / Dispositivo di bloccaggio, facile da afferrare e da spostare / Llave conmutadora diseñada anatómicamente para facilitar el accionamiento del mecanismo de cierre.

Hardened steel dispensing plate does not bend as easily when empty cartridge is encountered.


Plaque de distribution en acier durci ne se plie pas aussi facilement quand la cartouche est vide / Die Abgabeplatte aus Hartstahl ist bei leeren Kartuschen biegefest / Piastra di distribuzione, di acciaio indurito. Non si piega facilmente se incontra cartucce vuote / La placa dispensadora de acero templado no se dobla fácilmente cuando se encuentra con cartuchos vacíos.

Nonbreakable cartridges fit directly into dispenser; no attachments required.

Cartouches incassables se logent directement dans le distributeur, pas de fixation requise / Die bruch sicheren Kartuschen werden direkt in den Dispenser eingelegt; keine Haltevorrichtung erforderlich / Cartucce infrangibili, entrano direttamente e perfettamente nel distributore / Los cartuchos irrompibles calzan directamente en el dispensador; no se requiere ningún accesorio.

Specially designed metal rods press discs into uniform contact with surface of media. This eliminates manual tamping of each disc with forceps.

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Tamper latch prevents contamination of metal tamper rod when cartridge position is empty.

Loquet d'apporteur empêche la contamination du barreau applicateur métallique quand la cartouche est vide / Eine Stempelsperre verhindert die Kontamination der Metallstäbe bei leeren Kartuschenkanälen / Sistema di bloccaggio dell'applicatore – per impedire la contaminazione dell'asticella metallica quando la cartuccia non è inserita / El mecanismo de cierre de los apisonadores impide la contaminación de la varilla metálica de éstos cuando el cartucho queda vacío.

INTENDED USE

The **Sensi-Disc™ Designer Dispenser** accommodates up to twelve cartridges of **Sensi-Disc** susceptibility test discs. It is used to deliver discs onto 150 mm-style Petri dishes when performing disc susceptibility testing, especially by the standardized method.¹⁻⁴

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PRINCIPLES OF THE PROCEDURE

Discs are dispensed when downward pressure is applied to the operating knob. This unit is designed to provide 33 mm distance between disc centers.

No adjusting ring is required to control the height of the dispenser base above the agar surface. The discs are placed in a consistently uniform pattern as long as the depth of the agar is 3 to 5 mm thick and the height of the dish is less than 26 mm. The unit will automatically tamp the discs onto the agar surface and, therefore, eliminates the need for this manual step.

Each dispenser is provided complete with interlocking base and cover to protect the unit from moisture when not in use.

It is recommended that a desiccant be placed in the dispenser base to help maintain a low-moisture environment. Each cartridge is packed with desiccant capsules that should be replaced each time the cartridges are replaced, or more often, as required. Each dispenser is packed with a desiccant canister that should be regenerated when the silica gel, as viewed through the monitor window on the top of the canister, has turned pink. To regenerate, heat the canister in a vented oven at 300°F for approximately 3 h, or until the deep blue color returns. Upon cooling, the canister may be reused and will be capable of re-absorbing moisture to its full capacity. Replacement desiccant canisters can also be ordered.

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Warnings and Precautions:

For in vitro Diagnostic Use.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, prepared tubes, specimen containers and other contaminated materials must be sterilized by autoclaving before discarding.

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OPERATING PROCEDURE

This dispenser has been carefully manufactured and thoroughly tested. Trouble-free performance is assured by understanding and practicing the correct operating procedures.

The following instructions describe the operation of the dispenser and precautions to be observed.

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
Preparation

1. Remove the dispenser from its storage case and place on work surface.
2. Slide the black button to the "Lock" position. The dispenser will not operate if in the "Unlock" position.
3. Lift the dispenser about an inch and set it down sharply on the work surface. Do this with the cartridges removed at the start of the day's usage or whenever the tampers have become displaced by turning dispenser upside down. This action assures that the tampers are returned to their operating positions.

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


Dry Cycling

1. Lift the dispenser and place it over an empty Petri dish.
2. Depress the knob firmly until the sound of the tampers is heard. (If the dispenser will not dry cycle, tap the side of the unit against the side of the work surface and repeat the procedure.)
3. Allow the knob to return to the upper position.
4. Repeat steps 1, 2, and 3 several times to get the feel of the dispenser.

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
Inserting Cartridges

1. Slide the black button to the "Unlock" position.
2. Insert the cartridges into the numbered holes and press down lightly.
3. Slide the black button to the "Lock" position. If resistance is encountered, the cartridges are not properly seated. To correct, steps 1 and 2 should be repeated.

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
Cartridges

The last disc in each cartridge is marked with an "X". A green plastic slug has been incorporated in each cartridge after the "X" marked disc. This prevents the dispenser from operating when an empty cartridge is in place. Watch for the appearance of discs marked "X" and immediately insert full cartridges.

When inserting cartridges into the dispenser, make certain that the first disc is centered in the end of the cartridge.

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Dispenser

The dispenser must always complete a full cycle before it can be operated again. However, if greater than normal resistance is encountered in depressing the knob, STOP and immediately begin the "CLEARING PROCEDURE."

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
INSTRUCTIONS FOR DISPENSING DISCS

1. Lift dispenser and place over agar plate.
2. Depress knob **firmly** until the sound of the tampers is heard.
Note: Partial depression of the knob may cause discs to jam in the dispenser. (See "CLEARING PROCEDURE.")
3. Allow knob to return to the upper position and lift dispenser from agar plate.
4. Inspect the agar plate for the "X" on any of the dispensed discs.
5. Replace any empty cartridges.

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CLEARING PROCEDURE

1. Reversing the cartridge, slide the black button to the "Unlock" position. Remove all the cartridges and reinsert into the dispenser upside down.

Place the black button in the "Lock" position.

Carefully inspect the end of each cartridge for the presence of a green plastic slug or deformed discs. The green slug indicates that the cartridge is empty and must be replaced. Deformed discs must be manually removed and discarded.

2. Dry Cycling the Dispenser

With the cartridges upside down in the dispenser, dry cycle the dispenser several times to assure that it is operating freely. This completes the "CLEARING PROCEDURE."

Note: Occasionally, a disc will be delivered on edge and embedded in the agar; immediately remove the upright disc from the medium with sterile forceps and lay it flat on the plate centered over the point of initial contact. Horizontally embedded discs may occasionally occur. This could affect the accurate measurement of the zone area. Therefore, it is recommended that the test be repeated before accepting a minimally intermediate or minimally susceptible result.

If a disc is embedded in the medium, decontamination of the tamper is recommended. If no disc is delivered, or if a disc is delivered on edge, the tamper may be contaminated; to determine if tamper contamination has occurred, carefully examine the medium around the disc for tamper marks. If any marks are noted, it must be assumed that the tamper has become contaminated and the tamper should be decontaminated by following the "CLEANING AND DECONTAMINATION PROCEDURE."

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CLEANING AND DECONTAMINATION PROCEDURE

Each time that cartridges are changed, or if a tamper becomes contaminated during use by touching the inoculated agar surface, clean and decontaminate the tamper(s) by the following procedure:

1. Remove all **Sensi-Disc** cartridges from the unit and replace the lock button in the "Lock" position.
2. Place a **Sensi-Disc** cartridge upside down in the dispenser in the position that is contaminated; or, for periodic cleaning, place cartridges upside down in all positions.
3. Fill one empty 150 mm style Petri dish with 65 to 70 mL of a 3% solution of *Lysol.™
4. Fill one empty 150 mm style Petri dish with 65 to 70 mL of an 85% solution of isopropyl alcohol.
5. Fill two empty 150 mm style Petri dishes with 65 to 70 mL of sterile distilled water.
6. Place the dispenser over the dish containing Lysol solution.
7. Push the operating knob down until the tamper(s) are released. Hold the knob in the depressed condition for 30 seconds.
8. Repeat steps 6 and 7 with the isopropyl alcohol solution.
9. Repeat steps 6 and 7 with each of the dishes of distilled water. This will help ensure that the cleaning agents have been completely removed.

The tampers have now been cleaned and decontaminated. Any remaining alcohol and water left on the tampers will evaporate. The unit may be stored, or if testing has not been completed, dry the tampers by releasing them several times (depressing the knob fully) while the unit is placed over a sterile towel, gauze pad, swab, or other sterile absorbent material. The tampers can also be allowed to air dry, but this will require additional time, depending on the temperature and relative humidity in the area in which the dispenser is being stored.

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CARE OF THE DISPENSER

1. If desired, the dispenser and protective base and cover may be wiped clean with water, isopropyl alcohol, or a 3% Lysol solution. *Do not use alkaline cleaning agents, organic solvents, or abrasive cleaners.* After cleaning, always air dry the unit thoroughly before inserting cartridges.
2. Avoid the use of either an excessively cold or warm dispenser. Always allow it to reach room temperature before using. Return loaded dispenser to refrigerator as quickly as possible after use. Do not store in freezer.

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TROUBLE SHOOTING GUIDE

Read the insert before operating the BBL™ Sensi-Disc™ Dispenser.

If one or more discs fail to dispense:

1. Ensure the operating knob is UP.
2. Ensure the lock button is on "Unlock."
3. Examine for empty or damaged cartridges.
4. Examine for loose discs.

To clear loose or jammed discs:

1. Remove cartridges.
2. Set lock button on "Lock."
3. Push knob several times.
4. Examine bottom for trapped discs; if present, place inverted cartridges in appropriate positions and operate dispenser over a sterile surface.

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SERVICE AND MAINTENANCE INFORMATION

Should any malfunction occur that cannot be corrected by the "CLEARING PROCEDURE," "CLEANING AND DECONTAMINATION PROCEDURE," or "TROUBLE SHOOTING GUIDE," in the United States, telephone Technical Services for instructions, (800) 638-8663.

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WARRANTY INFORMATION

Becton Dickinson extends a two-year warranty for the **BBL™ Sensi-Disc™** Dispenser from date of receipt. This warranty entitles the holder to repair or replacement at no charge if found to be defective as a result of normal wear. To obtain repair or replacement within the terms of this warranty, in the United States please telephone Technical Services for instructions (800) 638-8663.

This warranty is null and void if the **BBL™ Sensi-Disc™** Dispenser has been disassembled or subjected to misuse or abuse. Becton Dickinson makes no warranty expressed or implied except as stated above.

To activate this warranty, please complete and mail the Business Reply Card within 10 days of receipt of your **BBL™ Sensi-Disc™** Dispenser. If the Business Reply Card is not returned, the warranty will expire two years after the date of manufacture.

Complete this section for your records. Serial Number: _____ Purchase Date: _____

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REFERENCES

1. Bauer, A.W., W.M.M. Kirby, J.C. Sherris, and M. Turck. 1966. Antibiotic susceptibility testing by a standardized single disk method. *Am. J. Clin. Pathol.* 45: 493-496.
2. Ryan, K.J., F.D. Schoenkecht, and W.M.M. Kirby. 1970. Disc sensitivity testing. *Hospital Practice.* 5: 91-100.
3. Food and Drug Administration. Standardized disc susceptibility test. Section 460.1, Title 21 of the Code of Federal Regulations.
4. National Committee for Clinical Laboratory Standards. 1997. Approved Standard M2-A6, Performance standards for antimicrobial disk susceptibility tests, 6th ed. National Committee for Clinical Laboratory Standards, Wayne. Pa.

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Technical Information: In the United States, contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Lysol is a trademark of Reckitt and Colman, Inc.

BD, BD Logo, and all other trademarks are property of Becton, Dickinson and Company. ©2015 BD.

**BBL™ Sensi-Disc™ Designer Dispenser
Certificate of Performance**

Becton, Dickinson and Company
extends this certificate to

Name of Hospital/Laboratory

for the period of two (2) years commencing _____

Month Day Year

This certificate entitles the holder to the repair or replacement at no charge, for a period of two years, of the **BBL™ Sensi-Disc™** Designer Dispenser if found to be defective as a result of normal wear.

To obtain repair or replacement within the terms of this warranty, in the United States, telephone Becton Dickinson Technical Services at 800-638-8663. During the call, arrangements will be made for repair/ replacement of the dispenser. When the dispenser is returned, a return address label (complete with Building, Room, Person Responsible, Purchase Order Number, if needed, etc.) must be included with the dispenser. The serial number in the center of the bottom of the dispenser should be recorded below.

This warranty is null and void if the **BBL™ Sensi-Disc™** Designer Dispenser has been disassembled or subjected to misuse or abuse. The company makes no warranty expressed or implied except as stated above.


To activate this Warranty the enclosed card must be completed and mailed within 10 days of receipt of your **BBL™ Sensi-Disc™** Designer Dispenser. Complete and retain this certificate for your records.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152


Serial No: _____

TRIM

PERF

NO POSTAGE
NECESSARY
IF MAILED
IN THE
UNITED STATES

BUSINESS REPLY MAIL

FIRST CLASS MAIL PERMIT NO. 25 SPARKS, MD 21152

POSTAGE WILL BE PAID BY ADDRESSEE

BD Diagnostic Systems
7 Loveton Circle
Sparks, MD 21152 USA

BUSINESS REPLY CARD

The **BBL™ Sensi-Disc™** Designer Dispenser has been purchased by the undersigned who hereby requests a Warranty be placed on file for 2 years.

Signature _____

Hospital/Laboratory _____

Address _____

City, State, Zip _____

Telephone Number _____

Distributor _____

Date of Receipt _____

Serial No. _____

Designer Dispenser BBL Sensi-Disc / BBL Sensi-Disc Designer-Dispenser / BBL Sensi-Disc Designer Dispenser / Dispensador BBL Sensi-Disc Designer

CERTIFICAT DE GARANTIE / GARANTIEBESCHEINIGUNG / CERTIFICATO DI PERFORMANCE / CERTIFICADO DE DESEMPEÑO

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