**BD**

Veritor™ System

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

Step 6B: Starting Walk Away mode:

- Turn the Analyzer on by pressing the blue power button once.
- When the display window reads: “INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE,”
 - **Double-click** the blue power button

Step 7B: Adding the specimen

- When the display window reads, “ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY”
 - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System RSV device sample well).
 - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System RSV device.

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
- “DO NOT DISTURB TEST IN PROGRESS” appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- A countdown timer in the display window shows the remaining analysis time.

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ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes when AC power adapter is connected).

Step 10C: Remove the test device

- Pull the device out. The display will show “INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE” to indicate the Analyzer is ready to perform another test.

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
The display window will briefly show “SCAN CONFIG BARCODE.” This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration instructions. Ignore this message and postpone this process when an assay is awaiting analysis.

- When the display window reads: “INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE”
 - Double-click** the blue power button.

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
- When the display window reads: “ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY”:
 - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System RSV device sample well).
 - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System RSV device.
- NOTE: Squeezing the tube close to the tip may cause leakage.**

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
- “DO NOT DISTURB TEST IN PROGRESS” appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- A countdown timer in the display window shows the remaining analysis time.
- Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.**
- When analysis is complete, a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed. **If a printer is not connected, record the result before removing the assay device.**

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Step 10D: Remove the test device

- Pull the device out. The display will show INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE to indicate the Analyzer is ready to perform another test. Note that the Analyzer returns to **Analyze Now** mode at the conclusion of each read sequence.
- If the Veritor Plus Analyzer is connected to an LIS, a steady ENVELOPE symbol will appear to indicate that results are awaiting transmission. In the event that a network connection is not detected while the ENVELOPE symbol is still displayed, the Analyzer will queue all untransmitted results and attempt to transmit them when reconnected. If it is powered off during this time, it will attempt to transmit as soon as power is restored, and connection is re-established. A flashing envelope indicates that data are in the process of being transmitted.

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Each BD Veritor System RSV device contains both positive and negative internal/procedural controls:

1. The internal positive control validates the immunological integrity of the device, proper reagent function, and assures correct test procedure.
2. The membrane area surrounding test lines functions as a background check on the assay device.

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Change History

Revision	Date	Change Summary
(11)	2020-04	Changed name of BD Synapsys™ Microbiology Informatics Solution to BD Synapsys™ Informatics Solution. Additional result documentation capabilities made possible with the implementation of the BD Synapsys Informatics Solution, and with the addition of the BD Veritor InfoScan module and BD Veritor Plus Connect (replacing the BD Veritor System Reader and the BD InfoSync Module). In sections C and D, revised Steps 10C and 10D regarding the use of the BD Veritor Plus Analyzer with the BD Veritor InfoScan. Added access information to obtain the document from bd.com/e-labeling .
(12)	2020-05	Corrected typographical error - no content affected.

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US Customers only: For symbol glossary, refer to bd.com/symbols-glossary

Manufacturer / Производител / Yürobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Tootja / Fabricant / Proizvođač / Gyártó / Fabriccante / Аткарушы / 제조업체 / Gamintojas / Ražotājs / Tilvirker / Producent / Producător / Производител / Vyrobca / Proizvodač / Tillverkare / Üretici / Виробник / 生产厂商

Use by / Използвайте до / Spotřebujte do / Brug før / Verwendbar bis / Χρήση έως / Usar antes de / Kasutada enne / Date de péremption / 사용 기한 / Upotřebiti do / Felhasználhatóság dátuma / Usare entro / Дейін пайдалануға / Naudokite iki / Izljetiti līdz / Houdbaar tot / Bruks for / Stosować do / Prazo de validade / A se utiliza până la / Исползовать до / Použite do / Upotřebiti do / Använd före / Son kullanna tarhi / Використати до / 使用截止日期

YYYY-MM-DD / YYYY-MM (MM = end of month)
 YTTT-MM-DD / TTTT-MM (MM = kraj na meseca)
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
 AAAA-MM-DD / AAAA-MM (MM = slutning af måned)
 JJJJ-MM-TT / JJJJ-MM (MM = Monatsende)
 EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα)
 AAAA-MM-DD / AAAA-MM (MM = fin del mes)
 AAAA-KK-PP / AAAA-KK (KK = kuų lõpp)
 AAAA-MM-JJ / AAAA-MM (MM = fin du mois)
 GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)
 EEEE-HH-NN / EEEE-HH (HH = lõnar utolsó napja)
 AAAA-MM-GG / AAAA-MM (MM = fine mese)
 JOJOJO-AA-KK / JOJOJO-AA / (AA = айдың соңы)
 YYYY-MM-DD/YYYY-MM(MM = 월말)
 MMMM-MM-DD / MMMM-MM (MM = mēnesio pabaiga)
 GGGG-MM-DD/GGGG-MM (MM = mēneša beigas)
 JJJJ-MM-DD / JJJJ-MM (MM = einde maand)
 AAAA-MM-DD / AAAA-MM (MM = sluffen av måneden)
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
 AAAA-MM-DD / AAAA-MM (MM = fin do mês)
 AAAA-LL-ZZ / AAAA-LL (LL = sfârşitul lunii)
 TTTT-MM-DD / TTTT-MM (MM = конец месяца)
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
 GGGG-MM-DD / GGGG-MM (MM = kraj meseca)
 AAAA-MM-DD / AAAA-MM (MM = slutet av månaden)
 YYYY-AA-GG / YYYY-AA (AA = аяын сопу)
 PPPP-MM-DD / PPPP-MM (MM = кінець місяця)
 YYYY-MM-DD / YYYY-MM (MM = 月末)

Catalog number / Каталоген номер / Katalogové číslo / Katalognummer / Αριθμός καταλόγου / Número de catálogo / Katalognummer / Numéro catalogue / Kataloški broj / Katalogusszám / Numero di catalogo / Каталог номери / 카탈로그 번호 / Katalog / numeris / Kataloga numurs / Catalogous number / Numer katalogowy / Număr de catalog / Номер по каталогу / Katalogové číslo / Kataloški broj / Katalog numarası / Номер за каталогом / 目录号

Authorized Representative in the European Community / Оторизиран представител в Европейската общност / Autorizovaný zástupce pro Evropském společenství / Autoriseret repræsentant i De Europæiske Fællesskaber / Autorisierter Vertreter in der Europäischen Gemeinschaft / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Reprezentante autorizado en la Comunidad Europea / Volitatud esindaja Euroopa Nõukogu / Répresentant autorisé pour la Communauté européenne / Autorizirani predstavnik u Europskoj uniji / Meghatalmazott képviselő az Európai Közösségekben / Reprezentante autorizzato nella Comunità Europea / Европа қауымдастығында уәкілетті екіл / 유권 공동체의 위임 대표 / Igallotais atstovas Europos Bendrijos / Pilnvarotais pārstāvis Eiropas Kopienā / Bevoegde vertegenwoordiger in de Europese Gemeenschap / Autoriseret repræsentant i EU / Autoryzowane przedstawicielstwo we Wspólnocie Europejskiej / Reprezentante autorizado na Comunidade Europeia / Reprezentantul autorizat pentru Comunitatea Europeană / Уполномоченный представитель в Европейском сообществе / Autorizovaný zástupca v Európskom spoločenstve / Autorizovano predstavništvo u Evropskoj uniji / Autoriserad representant i Europeiska gemenskapen / Avrupa Topuluğu Yetkili Temsilcisi / Уполноважений представник у країнах ЄС / 歐洲共同體授權代表

In Vitro Diagnostic Medical Device / Медицински уред за диагностика ин витро / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicensk anordning / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / Dispositivo médico para diagnóstico in vitro / In vitro diagnostika meditsiinaparatuur / Dispositif médical de diagnostic in vitro / Medicinska pomagala za In Vitro Dijagnostiku / In vitro diagnostikai orvosi eszköz / Dispositivo medicale per diagnostica in vitro / Жасанды жағдайда жұпиретін медициналық диагностика аспабы / In Vitro Diagnostik 의료 기기 / In vitro diagnostikos prietaisais / Medicinas ierīces, ko lieto in vitro diagnostikā / Medisch hulpmiddel voor in-vitro diagnostiek / In vitro diagnostisk medinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Dispositiv medical pentru diagnostic in vitro / Медицинский прибор для диагностики ин витро / Medicinska pomôcka na diagnostiku in vitro / Medicinski uređaj za in vitro dijagnostiku / Medicinteknisk produkt för in vitro-diagnostik / In Vitro Diyagnostik Tibbi Cihaz / Медицинский пристрій для діагностики ин витро / 体外诊断医疗设备

Temperature limitation / Температурни ограничения / Teplotní omezení / Temperaturbegrænsning / Temperaturbegrenzung / Περιορισμοί θερμοκρασίας / Limitation de température / Temperaturi pirang / Limites de température / Dozvoljena temperatura / Hömörséklet határ / Limiti di temperatura / Температураны шектеу / 온도 제한 / Laikymo temperatūra / Temperaturus ierobežojumi / Temperaturlimit / Temperaturbegrænsning / Ograniczenie temperatury / Limites de temperatura / Limite de temperatură / Ограничение температуры / Ohraničenie teploty / Ograničenje temperature / Temperaturgräns / Sıcaklık sınırlaması / Обмеження температури / 温度限制

Batch Code (Lot) / Код на партидата / Kód (číslo) šarže / Batch-code (lot) / Batch-Code (Charge) / Κωδικός παρτίδας (παρτίδα) / Código de lote (lote) / Partii kod / Numéro de lot / Lot (kod) / Tétel száma (Lot) / Codice batch (lotto) / Τοπταιμα κωды / 배치 코드 (로트) / Partijos numeris (LOT) / Partijas kods (laidiens) / Lot nummer / Batch-code (parti) / Kod partii (seria) / Código do lote / Cod de serie (Lot) / Kod partii (not) / Kód série (šarža) / Kod serije / Partinummer (Lot) / Parti Kodu (Lot) / Код партії / 批号 (亚批)

Contains sufficient for <n> tests / Содержит е достатъчно за <n> теста / Dostatečné množství pro <n> testů / Indeholder tilstrækkelig til <n> tests / Ausreichend für <n> Tests / Περιέχει επαρκή ποσότητα για <n> εξετάσεις / Contenuto sufficiente para <n> pruebas / Küllaldane <n> testide jaoks / Contenu suffisant pour <n> tests / Sadržaj za <n> testova / <n> testezhez elegendő / Contenido suficiente per <n> test / <n> testteeri piisn жеткілікті / <n> 테스트가 충분히 포함됨 / Pakankamas kiekis atlikti <n> testų / Satur pietiekami <n> pārbaudām / Inhoud voldoende voor "n" testen / Inneholder tilstrækkelig til <n> tester / Zawiera ilość wystarczającą do <n> testów / Conteúdo suficiente para <n> testes / Conținut suficient pentru <n> teste / Достаточно для <n> тестов(a) / Obsah vystačí na <n> testov / Sadržaj dovoljan za <n> testova / Innehåller tillräckligt för <n> analyser / <n> test için yeterli miktarda içerir / Вистачить для аналізів: <n> / 足够进行 <n> 次检测

Consult Instructions for Use / Направете справка в инструкциите за употреба / Prostudujte pokyny k použití / Se brugsanvisningen / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar las instrucciones de uso / Lugeda kasutusjuhendi / Consulter la notice d'emploi / Koristi upute za upotrebu / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Пайдалану нұсқаулығымен танысып алыңыз / 사용 지침 참조 / Skaitykite naudojimo instrukcijas / Skatīt lietošanas pamācību / Raadpleeg de gebruiksaanwijzing / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consultar as instruções de utilização / Consultați instrucțiunile de utilizare / См. руководство по эксплуатации / Pozri Pokyny na používanie / Pogledajte uputstvo za upotrebu / Se bruksanvisningen / Kullanim Talimatları'na başvurun / Див. інструкції з використання / 请参阅使用说明

Do not reuse / Не използвайте отново / Nepoužívejte opakovaně / Ikke til genbrug / Nicht wiederverwenden / Μην επαναχρησιμοποιείτε / No reutilizar / Keine kasutada korduvalt / Ne pas réutiliser / Ne koristiti ponovo / Egyezsr használatos / Non riutilizzare / Пайдаланбаңыз / 재사용 금지 / Tik vienkartiniam naudojimui / Nelietot atkārtoti / Nie opnieuw gebruiken / Kun til engangsbruk / Nie stosować powtórnie / Não reutilize / Nu refolosiți / Не использовать повторно / Nepoužívať opakované / Ne upotrebljavajte ponovo / Får ej återanvändas / Tekrar kullannayin / Не використовувати повторно / 请勿重复使用

Serial number / Серийн номер / Sériové číslo / Seriennummer / Seriennummer / Σειριακός αριθμός / N° de serie / Seerianumber / Numéro de série / Serijski broj / Sorozatszám / Numero di serie / Топтамалық нөмірі / Serijos numeris / Serijas numurs / Serie nummer / Numer serijny / Número de série / Număr de serie / Серийный номер / Seri numarası / Номер серії / 序列号

For IVD Performance evaluation only / Само за оценка качеството на работа на работа на IVD / Pouze pro vyhodnocení výkonu IVD / Kun til evaluering af IVD ydelse / Nur für IVD-Leistungsbewertungszwecke / Móvo για αξιολόγηση απόδοσης IVD / Sólo para la evaluación del rendimiento en diagnóstico in vitro / Ainult IVD seadme hindamiseks / Réserve à l'évaluation des performances IVD / Samo u znanstvene svrhe za In Vitro Dijagnostiku / Kizárólag in vitro diagnosztikához / Solo per valutazione delle prestazioni IVD / Жасанды жағдайда «пробирка ішінде» диагностикада тек жұмысты бағалау үшін / IVD 성능 평가에 대해서만 사용 / Tik IVD prietaisų veikimo charakteristikoms tikrinti / Vienigti IVD darības novērtēšanai / Uitsluitend voor doeltreffendheidsonderzoek / Kun for evaluering af IVD-ydelse / Tylko do oceny wydajności IVD / Uso exclusivo para avaliação de IVD / Numai pentru evaluarea performanței IVD / Только для оценки качества диагностики ин витро / Určené iba na diagnostiku in vitro / Samo za procenu učinka u in vitro dijagnostici / Endast för utvärdering av diagnostisk användning in vitro / Yalnızca IVD Performans değerlendirmesi için / Тільки для оцінювання якості діагностики ин витро / 仅限 IVD 性能评估

For US: "For Investigational Use Only"

Lower limit of temperature / Долен лимит на температурата / Dolní hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Κατώτερο όριο θερμοκρασίας / Limite inferior de temperatura / Alumine temperatuuripiir / Limite inférieure de température / Najnižia dozvoljena temperatura / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Температураның төменгі рұқсат шегі / 하한 온도 / Žemiausia laikymo temperatūra / Temperatūras zemākā robeža / Laagste temperatuurlimiet / Nedre temperaturgrænse / Dolna granica temperatury / Limite minimo de temperatura / Limită minimă de temperatură / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sıcaklık alt sınırı / Мінімальна температура / 温度下限

Control / Контроль / Kontrola / Kontrol / Kontrolle / Μάρτυρας / Kontroll / Contrôle / Controllo / Бақылау / 컨트롤 / Kontrolé / Kontrolle / Controle / Controllo / Контроль / kontroll / Контроль / 对照

Positive control / Положительный контроль / Pozitívny kontrol / Positiv kontrol / Positive Kontrolle / Ετικός μάρτυρας / Control positivo / Positive kontrol / Contrôle positif / Pozitívna kontrola / Pozitív kontroll / Controllo positivo / Оң бақылау / 양성 컨트롤 / Teigiamia kontrolė / Pozitívna kontrola / Positive controle / Kontrola dodatnia / Controllo positivo / Control pozitiv / Положительный контроль / Pozitív kontrol / Позитивный контроль / 阳性对照试剂

Negative control / Отрицательный контроль / Negatívny kontrola / Negativ kontrol / Negative Kontrolle / Αρνητικός μάρτυρας / Control negativo / Negative kontrol / Contrôle négatif / Negatívna kontrola / Negatív kontroll / Controllo negativo / Негативтік бақылау / 음성 컨트롤 / Neigiamia kontrolė / Negatívna kontrola / Negative controle / Kontrola ujemna / Controllo negativo / Control negativ / Отрицательный контроль / Negatív kontrol / Негативный контроль / 阴性对照试剂

	Method of sterilization: ethylene oxide / Μέθοδος στεριλίστασης: εθιλενοξείδιο / Způsob sterilizace: etylenoxid / Steriliseringsmetode: ethylenoxid / Sterilisationsmethode: Ethylenoxid / Μέθοδος αποστείρωσης: αιθυλενοξείδιο / Método de esterilización: óxido de etileno / Sterilisieringsmetode: etüleenoksiid / Méthode de stérilisation: oxyde d'éthylène / Metoda sterilizacije: etilen oksid / Sterilizálás módszere: etilén-oxid / Metodo di sterilizzazione: ossido di etilene / Sterilizacija adici – этилен тотыгы / 소독 방법: 에틸렌옥사이드 / Sterilizavimo būdas: etileno oksidas / Sterilizācijas metode: etilēnoksīds / Gesteriliseerd met behulp van ethyleenoxide / Steriliseringsmetode: etylenoksid / Metoda sterylizacji: tlenek etylu / Método de esterilização: óxido de etileno / Metoda de sterilizare: oxid de etilenă / Metod sterilizacije: этиленоксид / Metoda sterilizacije: etylenoxid / Metoda sterilizasyon yöntemi: etilen oksit / Method sterilizacji: этиленоксид / 灭菌方法: 环氧乙烷
	Method of sterilization: irradiation / Μέθοδος στεριλίστασης: ακτινοβολία / Způsob sterilizace: záření / Steriliseringsmetode: bestraling / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Sterilisieringsmetode: kirkus / Méthode de stérilisation: irradiation / Metoda sterilizacije: zračenje / Sterilizálás módszere: besugárzás / Metodo di sterilizzazione: irradiazione / Sterilizacija adici – сәуле түсірі / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizācijas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringsmetode: bestraling / Metoda sterylizacji: napromienianie / Método de esterilização: irradiação / Metoda de sterilizare: iradiere / Metod sterilizacije: облучение / Metoda sterilizacije: ožiarovanie / Metoda sterilizacije: ozračevanje / Steriliseringsmetode: stråling / Sterilizasyon yöntemi: ırradyasyon / Method sterilizacji: opromienienie / 灭菌方法: 辐射
	Biological Risks / Биологични рискове / Biologická rizika / Biologisk fare / Biogefährdung / Βιολογικοί κίνδυνοι / Riesgos biológicos / Biologiskie ryzyko / Biologiskie risici / Biologisk risiko / Biologisk risico / Zagrożenia biologiczne / Perigo biológico / Riscuri biologice / Биологическая опасность / Biologické riziko / Biološki rizici / Biologisk risk / Biyolojik Riskler / Біологічна небезпека / 生物学风险
	Caution, consult accompanying documents / Внимание, направте справка в придружаващите документи / Pozor! Prostudujte si přílohu dokumentaci! / Forsigtig, se ledsagende dokumenter / Achtung, Begleitdokumente beachten / Προσοχή, συμβουλευτείτε τα συνοδούσα έγγραφα / Precaución, consultar la documentación adjunta / Ettevaatusti! Lugea kaasnevat dokumentatsiooni / Attention, consulter les documents joints / Upozorenje, koristi prateću dokumentaciju / Fügeyen! Olvassa el a mellékelt tájékoztatót / Attenzione: consultare la documentazione allegata / Абайлаңыз, тиісті құжаттармен танысыңыз / 주의, 동봉된 설명서 참조 / Demoso, žiūrėkite priedamus dokumentus / Piesardzība, skatīt pavaddokumentus / Voorzichtig, raadpleeg bijgevoegde documenten / Forsiktig, se vedlagte dokumentasjon / Należy zapoznać się z dołączonymi dokumentami / Cuidado, consulte a documentação fornecida / Atentie, consultați documentele însoțitoare / Внимание: см. прилагаемую документацию / Vytřaďte, poriz správné dokumenty / Pažnja! Pogledajte priložena dokumenta / Obs! Se medföljande dokumentation / Dikkat, birlirte verilen belgeler bağvurun / Увага: див. супутню документацию / 小心, 请参阅附带文档。
	Upper limit of temperature / Горен лимит на температурата / Horní hranice teploty / Øvre temperaturgrense / Temperaturobergrenze / Ανώτερο όριο θερμοκρασίας / Limite superior de temperatura / Ulemine temperatuuripiiri / Limite supérieure de température / Gornja dozvoljena temperatura / Felső hőmérsékleti határ / Limite superiore di temperatura / Температураның рұқсат етілген жоғарғы шегі / 상한 온도 / Aukščiausia laikymo temperatūra / Augšējā temperatūras robeža / Hoogste temperatuurlimiet / Øvre temperaturgrense / Gorna granica temperatury / Limite máximo de temperatura / Limitā maxīmā de temperatūrā / Верхний предел температуры / Horná hranica teploty / Gornja granica temperature / Øvre temperaturgräns / Sıcaklık üst sınırı / Максимальна температура / 温度上限
	Keep dry / Пазете сухо / Skladujte v suchém prostředí / Opbevaras tørt / Trocklegen / Φυλάξτε το στεγνό / Mantener seco / Hoida kuivas / Conserver au sec / Držati na suhom / Száraz helyen tartandó / Tenere all'asciutto / Кудрак күйнде уста / 건조 상태 유지 / Laikyti sausai / Uzglabāt saus / Droog houden / Hoides tørt / Przechowywać w stanie suchym / Manter seco / A se feri de umezeală / Не допускать попадания влаги / Uchovávajte v suchu / Držite na suvom mestu / Förvaras tørt / Kuru bir şekilde muhafaza edin / Бергеті від вологості / 请保持干燥
	Collection time / Време на събиране / Čas odběru / Opsamlingsstidspunkt / Entnahmeuhrzeit / Ώρα συλλογής / Hora de recogida / Koğumisaeg / Heure de prélèvement / Sati prikupljanja / Mintavétel időpontja / Ora di raccolta / Жинау уакыты / 수집 시간 / Paėmimo laikas / Savākšanas laiks / Verzameltijd / Tid prøvetaking / Godzina pobrania / Hora de coleta / Ora colectării / Время сбора / Doba odberu / Vreme prikupljanja / Uppsamlingsstid / Toplama zamanı / Час забору / 采集时间
	Peel / Обелете / Otefete zde / Åbn / Abziehen / Αποκολλήστε / Desprendre / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Ystígní kabátyn alyp tasta / 뗀기 / Plésti čia / Atfílmēt / Schillen / Trekk av / Oderwac / Destacar / Se dezlipze / Отклеить / Odtrhnite / Oljuštiti / Dra isår / Ayırma / Відклеїти / 撕下
	Perforation / Перфорация / Perforace / Perforering / Διάτρηση / Perforación / Perforatsioon / Perforacija / Perforálás / Perforazione / Tescik tescy / 찢히선 / Perforacja / Perforăcia / Perforatie / Perforacja / Perfuração / Perforare / Перфорация / Perforácia / Perforasyon / Перфорация / 穿孔
	Do not use if package damaged / Не използвайте, ако опаковката е повредена / Nepoužívejte, je-li obal poškozený / Må ikke anvendes hvis emballagen er beskadiget / Inhal beschädigter Packungnicht verwenden / Μην χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά. / No usar si el paquete está dañado / Mitte kasutada, kui pakend on kahjustatud / Ne pas l'utiliser si l'emballage est endommagé / Ne koristiti ako je oštećeno pakiranje / Ne használnia, ha a csomagolás sérült / Non usare se la confezione è danneggiata / Егер пакет бузылган болса, пайдаланба / 패키지가 손상된 경우 사용 금지 / Jei pakuotė pažeista, nenaudoti / Nelietot, ja ierakojums bojāts / Niet gebruiken indien de verpakking beschadigd is / Må ikke brukes hvis pakke er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Nào usar se a embalagem estiver danificada / A nu se folosi dacă pachetul este deteriorat / Не использовать при повреждении упаковки / Nepoužívať, ak je obal poškodený / Ne koristite ako je pakovanje oštećeno / Använd ej om förpackningen är skadad / Ambalaj hasar görmüşse kullanmayın / Не використовувати за пошкодженої упаковки / 如果包装破损, 请勿使用
	Keep away from heat / Пазете от топлина / Nevystavujte přílišnému teplu / Må ikke udsættes for varme / Vor Wärme schützen / Κρατήστε το μακριά από τη θερμότητα / Mantener alejado de fuentes de calor / Hoida eemal valgusest / Protéger de la chaleur / Držati dalje od izvora topline / Oveja a melegtől / Tenere lontano dal calore / Салкын жерде сакта / 열을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargatí no karstuma / Beschermen tegen warmte / Må ikke utsettes for varme / Przechowywać z dala od źródła ciepła / Manter ao abrigo do calor / A se feri de căldură / He harpevat / Uchovávajte mimo zdroja tepla / Držite dalje od toplote / Får ej utsättas för värme / Isidan uzak tutun / Бергеті від тепла / 请远离热源
	Cut / Срежете / Odstřihněte / Kutt / Schneiden / Κόψτε / Cortar / Lõigata / Découper / Reži / Vája / Tagliare / Keciñiz / 잘라내기 / Kirpi / Nogriezet / Knippen / Kutt / Odciać / Cortar / Decupati / Orzepsat / Odstřihnite / Iseći / Klipp / Kesme / Розрізати / 剪下
	Collection date / Дата на събиране / Datum odběru / Opsamlingsdato / Entnahmedatum / Ημερομηνία συλλογής / Fecha de recogida / Koğumiskupavēn / Date de prélèvement / Dani prikupljanja / Mintavétel dátuma / Data di raccolta / Жинаған тізбекүні / 수집 날짜 / Paėmimo data / Savākšanas datums / Verzameldatum / Dato prøvetaking / Data pobrania / Data de coleta / Data colectării / Дата сбора / Datum odberu / Datum prikupljanja / Uppsamlingsdatum / Toplama tarihi / Дата забору / 采集日期

µL/test / µL/test / µL/Test / µL/εξέταση / µL/prueba / µL/teszt / µL/테스트 / мкл/тест / µL/tyrimas / µL/pārbaude / µL/teste / мкл/анализ / µL/检测

Keep away from light / Пазете от светлина / Nevystavujte světlu / Må ikke udsættes for lys / Vor Licht schützen / Κρατήστε το μακριά από το φως / Mantener alejado de la luz / Hoida eemal valgusest / Conserver à l'abri de la lumière / Drżać dalej od światła / Fény nem érheti / Tenere al riparo dalla luce / Қараңғыланған жерде ұста / 빛을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no gaismas / Niet blootstellen aan zonlicht / Må ikke utsettes for lys / Przechowywać z dala od źródeł światła / Manter ao abrigo da luz / Ferijti de lumină / Хранить в темноте / Uchovávať mimo dosahu svetla / Držite dalej od svetlosti / Får ej utsättas för ljus / Išiklanti užak tutun / Берегти від дії світла / 请远离光线

Hydrogen gas generated / Образуван е водород газ / Možnost úniku plyného vodíku / Frembringer hydrogengas / Wasserstoffgas erzeugt / Δημιουργία αερίου υδρογόνου / Producción de gas de hidrógeno / Vesinikgaasi tekitatud / Produit de l'hydrogène gazeux / Sadržaji hidrogen vodik / Hidrogeen gázt fejeleszt / Produzione di gas idrogeno / Гастектес сытери пайда болды / 수소 가스 생성됨 / Išiskiria vandenilio dujas / Rodas ūdenradis / Waterstofgas gegenereerd / Hydrogengass generert / Powoduje powstawanie wodoru / Produção de gás de hidrogénio / Generare gaz de hidrogen / Выделение водорода / Vyrobené použitím vodíka / Oslobada se vodonik / Genererad vätegas / Açıda çıkan hidrojen gazı / Реакція з виділенням водню / 会产生氢气

Patient ID number / ИД номер на пациента / ID pacienta / Patientens ID-nummer / Patienten-ID / Αριθμός αναγνώρισης ασθενούς / Número de ID del paciente / Patsiendi ID / No d'identification du patient / Identifikacijski broj pacijenta / Beteg azonosító száma / Numero ID paziente / Пациенттің идентификациялық нөмірі / 환자 ID 번호 / Pacienta identifikavimo numeris / Pacienta ID numurs / Identificatienummer van de patiënt / Patientsens ID-nummer / Numer ID pacienta / Número da ID do doente / Număr ID pacient / Идентификационный номер пациента / Identifikačné číslo pacienta / ID broj pacijenta / Patientnummer / Hasta kimlik numarası / Идентификатор пациента / 患者标识号

Fragile, Handle with Care / Чупливо, Роботете с необходимото внимание. / Křehké. Při manipulaci postupujte opatrně. / Forsigtig, kan gå i stykker. / Zerbrechlich, vorsichtig handhaben. / Εύθραστο. Χειριστείτε με προσοχή. / Frágil. Manipular con cuidado. / Őrn, kásitsege ettevaatlikult. / Fragile. Manipuler avec précaution. / Lomljivo, rukujte pažljivo. / Törékeny! Óvatosan kezelendő. / Fragile, maneggiare con cura. / Сығыш, абайлап пайдаланыңыз. / 조심 깨지기 쉬운 처리 / Trapu, elkitės atsargiai. / Trausis; rūkotiės uzmanigi / Breekbaar, voorzichtig behandelen. / Ømtålig, håndter forsigtig. / Krucha zawartość, przenosić ostrożnie. / Frágil, Manuseie com Cuidado. / Frágil, manipulați cu atenție. / Хрупкое! Обращаться с осторожностью. / Křehké, vyžaduje sa opatrná manipulácia. / Lomljivo - rukujte pažljivo. / Bräckligt. Hantera försiktigt. / Kolay Kırılır, Dikkatli Taşın. / Тендітна, звертатися з обережністю / 易碎，小心轻放

R_x Only

This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

bd.com/e-labeling

Europe, CH, GB, NO:		+800 135 79 135
International:		+31 20 794 7071
AR	+800 135 79 135	LT 8800 30728
AU	+800 135 79 135	MT +31 20 796 5693
BR	0800 591 1055	NZ +800 135 79 135
CA	+1 855 805 8539	RO 0800 895 084
CO	+800 135 79 135	RU +800 135 79 135
EE	0800 0100567	SG 800 101 3366
GR	08000 161 22015 7799	SK 0800 606 287
HR	0800 804 804	TR 00800 142 064 866
IL	+800 135 79 135	US +1 855 236 0910
IS	800 8996	UY +800 135 79 135
LI	+31 20 796 5692	VN 122 80297

Becton, Dickinson and Company
7 Loveton Circle
Sparks, Maryland 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

BD, the BD logo, Synapsys, and Veritor are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2020 BD. All rights reserved.

Insert ends here

DO NOT PRINT THIS PAGE

CLSI document should not be included in print file to vendor.

NOTES:

The CLSI document is to be updated at the same time as the package insert and should be released as a separate PDF document to the web.

The package insert should also be released as a separate PDF document to the web.

The package insert and CLSI document should be released as one document to SAP.

Final art provided to Buyer/print vendor should not include this page or the CLSI portion of this document.

Proofreader: please proof the CLSI document along with the package insert.

Procedure* BD Veritor™ System For Rapid Detection of Respiratory Syncytial Virus (RSV)

Laboratory kit configured for testing liquid nasopharyngeal wash, aspirate and swab in transport media samples.

Prepared by	Date Adopted	Supersedes Procedure #

Review Date	Revision Date	Signature

Distributed to	# of Copies	Distributed to	# of Copies

***Any modifications to this document are the sole responsibility of the facility. This “Sample Procedure” is not intended as a substitute for your facility procedure manual, instrument manual, or reagent labeling/package insert. This “Sample Procedure” is intended as a model for use by your facility to meet the needs of your laboratory.**

For Rapid Detection of Respiratory Syncytial Virus (RSV)

For in vitro diagnostic use only.

INTENDED USE

The BD Veritor™ System for Rapid Detection of Respiratory Syncytial Virus (RSV) is a chromatographic immunoassay with an instrumented read for the direct and qualitative detection of RSV fusion protein from nasopharyngeal wash, aspirate and swab in transport media samples from patients suspected of having a viral respiratory infection. This test is intended for *in vitro* diagnostic use to aid in the diagnosis of RSV infections in infants and pediatric patients under the age of 20 years. Negative results do not preclude RSV infection and should not be used as the sole basis for treatment or for other management decisions. A negative test is presumptive. It is recommended that negative test results be confirmed by viral cell culture or an alternative method, such as a FDA-cleared molecular assay. Outside the U.S. a negative test is presumptive and it is recommended that these results be confirmed by viral culture or a molecular assay cleared for diagnostic use in the country of use. FDA has not cleared this device outside the U.S. The test is intended for professional and laboratory use. It is to be used in conjunction with the BD Veritor System Instrument.

SUMMARY AND EXPLANATION

Viral respiratory tract infections are responsible for widespread disease. Respiratory syncytial virus (RSV) is a leading cause of lower respiratory tract infections (LRI) in young children in both the developed and developing worlds. Worldwide, it is estimated that RSV is responsible for greater than 30 million cases of LRI in children under 5 years of age each year.^{1,2} RSV has also been implicated in severe respiratory infections in the elderly and immunocompromised.^{3,4} RSV has been identified as causing 20% of "influenza-like" illness in people 15–44 years of age and is responsible for greater than 17,000 deaths per year in the United States, almost 80% of which occur in adults over age 65.^{5,6}

Diagnostic methods for detection of respiratory viruses include viral cell culture, direct fluorescent antibody (DFA), rapid immunoassays, and nucleic acid amplification assays such as the polymerase chain reaction (PCR).^{7,8} Each has been demonstrated to have clinical utility for the detection of respiratory viruses including RSV. Rapid immunoassays available for specific viruses such as influenza A/B and RSV allow a quick diagnosis so that patients may be appropriately isolated and treated to prevent the nosocomial spread of infections to fellow patients with compromised cardiac, respiratory or immune functions.⁹ In addition, rapid tests assist with the selection of appropriate antiviral therapy. The most common specimen types collected for RSV testing include nasopharyngeal washes, nasopharyngeal aspirates, nasal swabs and nasopharyngeal swabs.

The BD Veritor System for Rapid Detection of RSV (also referred to as the BD Veritor System and BD Veritor System RSV) is a digital immunoassay (DIA) to qualitatively detect RSV fusion protein from various specimens of symptomatic patients with a time to result of 10 minutes.

All BD Veritor System RSV test devices are interpreted by a BD Veritor System Instrument, either a BD Veritor Reader or BD Veritor Plus Analyzer (the "Analyzer"). When using an Analyzer, procedures to evaluate test devices depend on the workflow configuration chosen. In Analyze Now mode, the instrument evaluates assay devices after manual timing of their development. In **Walk Away** mode, devices are inserted immediately after application of the specimen, and timing of assay development and analysis is automated. Connection of the Analyzer to a printer is possible if desired. Additional result documentation capabilities are possible with the implementation of the BD Synapsys™ Informatics Solution, and with the addition of the BD Veritor InfoScan module and BD Veritor Plus Connect. Please refer to the Analyzer *Instructions for Use* for details on these features and contact BD technical support for more information.

PRINCIPLES OF THE PROCEDURE

The BD Veritor System for Rapid Detection of RSV is a qualitative, digital immunoassay assay for the detection of RSV fusion protein in samples processed from respiratory specimens. When specimens are processed and added to the test device, RSV antigen binds to anti-RSV antibodies conjugated to detector particles in the RSV test strip. The antigen-conjugate complex migrates across the test strip to the reaction area and is captured by the line of RSV antibody on the membrane. A positive result for RSV is determined by the BD Veritor System Instrument (purchased separately) when antigen-conjugate is deposited at the Test "T" position and the Control "C" position on the BD Veritor System RSV assay device. The instrument analyzes and corrects for non-specific binding and detects positives not recognized by the unaided eye to provide an objective digital result.

REAGENTS

The following components are included in the BD Veritor System for Rapid Detection of RSV kit:

BD Veritor System RSV Devices	30 devices	Foil pouched device containing one reactive strip. Each strip has one test line of monoclonal antibody specific to RSV viral antigen and murine monoclonal control line antibodies.
RV Reagent C	30 tubes with 100 µL reagent	Detergent with <0.1% sodium azide (preservative).
300 µL Pipette	30 each	Transfer pipette
RSV Positive Control Swab	1 each	RSV Positive Control Swab, RSV antigen (noninfectious cell lysate) with <0.1% sodium azide (preservative)
RSV Negative Control Swab	1 each	RSV Negative Control Swab, (detergent-treated non-infected cells) with <0.1% sodium azide (preservative)

Materials Required But Not Provided: BD Veritor™ Plus Analyzer (Cat. No. 256066), Timer, Tube Rack for specimen testing.

Optional Equipment: BD Veritor™ InfoScan Module (Cat. No. 256068), USB Printer Cable for BD Veritor™ Analyzer (Cat. No. 443907), Epson Printer model TM-T20 II, BD Veritor Plus Connect (contact BD Technical Services for details).

Warnings and Precautions:

Warning

H302 Harmful if swallowed. **H402** Harmful to aquatic life. **H412** Harmful to aquatic life with long lasting effects.

P273 Avoid release to the environment. **P264** Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P330** Rinse mouth. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

1. For *in vitro* Diagnostic Use.
2. Test results are not meant to be visually determined. **All test results must be determined using the BD Veritor System Instrument.**
3. The RSV Positive Control Swab and the positive control line on the BD Veritor System for Rapid Detection of RSV device have been prepared from RSV-infected tissue culture cells which have been inactivated by detergent treatment and sonication then subsequently tested by bioassay procedures.
4. Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁰⁻¹³ and institutional guidelines should be followed in handling, storing and disposing of all specimens and all items contaminated with blood and other body fluids.
5. Dispose of used BD Veritor System test devices as biohazardous waste in accordance with federal, state and local requirements.
6. Reagents contain sodium azide, which is harmful if inhaled, swallowed or exposed to skin. Contact with acids produces very toxic gas. If there is contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.
7. Do not use kit components beyond the expiration date.
8. Do not reuse the BD Veritor System test device.
9. Do not use the kit if the Control RSV Positive Swab and Control RSV Negative Swab do not yield appropriate results.
10. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
11. To avoid erroneous results, specimens must be processed as indicated in the assay procedure section. The addition of excess sample may give invalid test results.
12. Proper specimen collection, storage and transport are critical to the performance of this test.
13. Specific training or guidance is recommended if operators are not experienced with specimen collection and handling procedures.

Storage and Handling: Kits may be stored at 2–30 °C. DO NOT FREEZE. **Reagents and devices must be at room temperature (15–30 °C) when used for testing.**

SPECIMEN COLLECTION AND HANDLING

Specimen Collection and Preparation: Acceptable specimens for testing with the BD Veritor System for Rapid Detection of RSV include nasopharyngeal (NP) washes, aspirates and swab specimens in transport media. It is essential that correct specimen collection and preparation methods be followed. Specimens obtained early in the course of the illness will contain the highest viral titers.

Inadequate specimen collection, improper specimen handling and/or transport may yield a false negative result; therefore, training in specimen collection is highly recommended due to the importance of specimen quality for generating accurate test results.

Specimen Transport Media: The following transport media have been tested and found to be compatible using moderate positive samples with the BD Veritor System for Rapid Detection of RSV:

- Modified Amies Medium (liquid) ESwab, Liquid Stuart Medium, Amies, Bartel ViraTrans™, BD Universal Viral Transport, Hank's Balanced Salt solution, M4, M4-RT, M5, M6, Normal Saline, Phosphate Buffered Saline.

Samples in these transport media can be stored at 2–8 °C for up to 72 hours.

Other transport media may be utilized if an appropriate validation exercise is performed.

Specimen Transport and Storage:

Freshly collected specimens should be processed within 1 hour. If necessary, specimens may be stored at 2–8 °C for up to 72 hours. It is essential that correct specimen collection and preparation methods be followed. Do not centrifuge specimens prior to use, as the removal of cellular material may adversely affect test sensitivity.

Procedure for Nasopharyngeal Washes/Aspirates:

- For NP washes/aspirates, sample volumes of 1 to 3 mL are recommended. If transport medium is used, minimal dilution of specimens is recommended.
- Excessive wash volumes should be avoided as they may result in decreased test sensitivity.
- Process specimen as described in "Test Procedure".

Procedure for Nasopharyngeal Swabs in Transport Media:

- For NP swabs in transport media, a minimal volume of transport media (1 mL) is recommended.
- Process specimen as described in "Test Procedure".

TEST PROCEDURE


NOTES: Reagents, specimens and devices must be at room temperature (15–30 °C) for testing. Thoroughly mix all specimens prior to removal of an aliquot for processing. Do not centrifuge specimens.

The following steps assume that Users of a BD Veritor Plus Analyzer have chosen and set all configuration options, and that the Analyzer is ready to use. To choose or change these settings, see the BD Veritor Plus Analyzer Instructions for Use, section 4.7. A printer is not necessary to display results. However, if your facility has chosen to connect the BD Veritor Plus Analyzer to a printer, check that the Analyzer is plugged into a power source, paper supply is adequate and any necessary network connections are enabled before testing.

For each patient specimen or control swab:

- Step 1:** Remove one RV Reagent C tube/tip and one BD Veritor System RSV device from its foil pouch immediately before testing.
- Step 2:** Label one BD Veritor System device and one RV Reagent C tube for each specimen and control to be tested.
- Step 3:** Place the labeled RV Reagent C tube(s) in the designated area of the tube rack.


All reagents and specimens must be at room temperature prior to processing

Step 4: Process the Sample or Control:


a. For NP washes, aspirates and swab specimens in transport media:

1. Vortex or thoroughly mix specimen. Do not centrifuge.
2. Remove and discard the cap from the RV Reagent C tube corresponding to the sample to be tested.
3. Using the transfer pipette, transfer 300 µL of the specimen into the RV Reagent C tube. Discard pipette after use.


b. For kit swab controls:

1. Remove and discard the cap from the RV Reagent C tube corresponding to the sample to be tested.
2. Using the transfer pipette add 300 µL of distilled or deionized water to the RV Reagent C tube.
3. Insert the control swab into the tube and vigorously plunge the swab up and down in the fluid for a minimum of 15 seconds.
4. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.

Step 5:

5. Press the attached tip firmly onto the RV Reagent C tube containing the processed specimen or control (threading/twisting not required).
6. Vortex or mix thoroughly by swirling or flicking the bottom of the tube.


NOTE: Do not use tips from any other product, including other products from BD or other manufacturers.

After step 5, choose from the model and workflow option below before continuing to step 6:				
	BD Veritor Reader or Analyzer in Analyze Now mode	BD Veritor Plus Analyzer in Walk Away mode	BD Veritor Plus Analyzer with the BD Veritor InfoScan module In Analyze Now mode—or— Walk Away mode	
Instructions in section:	A	B	C	D

Step 6A: Adding the specimen

- Invert the RV Reagent C tube and hold the tube vertically (approximately one inch above the labeled BD Veritor System RSV device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well of a labeled BD Veritor System RSV device.

NOTE: Squeezing the tube too close to the tip may cause leakage.

Step 7A: Timing development


- After adding the sample, allow the test to run for **10 minutes** before inserting into the BD Veritor Instrument.

NOTE: If running test under laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.

Step 8A: Using the BD Veritor Instrument

- During incubation time, turn the BD Veritor Instrument on by pressing the power button once.
- Insert assay device when the 10-minute assay development time is complete.
- The status of the assay analysis process appears in the display window. Follow the on-screen prompts to complete the procedure.

Do not touch the Instrument or remove the test device

Step 9A: Record the Result


- When analysis is complete, the test result appears in the display window

ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).

To use Walk Away mode – connect the AC power adapter to the Analyzer and a power source


Step 6B: Starting Walk Away mode:

- Turn the Analyzer on by pressing the blue power button once.
- When the display window reads: "INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE,"
 - **Double-click** the blue power button

Step 7B: Adding the specimen

- When the display window reads, "ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY"
 - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System RSV device sample well).
 - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System RSV device.


NOTE: Squeezing the tube close to the tip may cause leakage.

Step 8B: Start the development and reading sequence

- Immediately insert the test device into the slot on the right side of the Analyzer.

The test device must remain horizontal to prevent spilling the specimen out of the sample well.

- "DO NOT DISTURB TEST IN PROGRESS" appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- A countdown timer in the display window shows the remaining analysis time.

Do not touch the Analyzer or remove the test device during this time. Doing so will abort the assay analysis.


Step 9B: Record the Result

- When analysis is complete, the test result appears in the display window

ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 60 minutes (when AC power adapter is connected).

Step 6C: Adding the specimen

- Invert the tube, holding it vertically (approximately one inch above the BD Veritor System RSV device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well of a labeled BD Veritor System RSV device.

NOTE: Squeezing the tube close to the tip may cause leakage.

Step 7C: Timing development

- Allow the test to develop for **10 minutes**.
- If running the test in a laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.



Step 8C: Using the Analyzer

During incubation time, turn the BD Veritor Plus Analyzer on by pressing the blue button once.


The display window briefly shows “SCAN CONFIG BARCODE.” This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration instructions. Ignore this message and postpone this process when an assay is awaiting analysis.

- When assay development time is complete and the Analyzer display window reads: “INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE”
 - Insert the BD Veritor System RSV device into the BD Veritor Plus Analyzer.

Step 9C: Using the Bar Code scanner

- Follow the prompts on the display window to complete any required barcode scans of:
 - OPERATOR ID
 - SPECIMEN ID and/or
 - KIT LOT NUMBER
 In accordance with site requirements and Analyzer settings.

- Prompts for each scanning step appear in the display window for only 30 seconds. Failure to complete scans during that time will cause the Analyzer to default to the beginning of step 8C. To restart this step, remove and reinsert the test device to initiate a new sequence.
- Move the barcode slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.
- The Analyzer can record the Kit Lot Number in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user. BD recommends against the use of expired materials.

After required scans are completed, the Analyzer displays a countdown timer and test analysis begins.

- **Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.**
- When analysis is complete a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed.

If a printer is not connected, record the result before removing the assay device.

ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).

Step 10C: Remove the test device

- Pull the device out. The display will show “INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE” to indicate the Analyzer is ready to perform another test.
- If the Veritor Plus Analyzer is connected to an LIS, a steady ENVELOPE symbol will appear to indicate that results are awaiting transmission. In the event that a network connection is not detected while the ENVELOPE symbol is still displayed, the Analyzer will queue all untransmitted results and attempt to transmit them when reconnected. If it is powered off during this time, it will attempt to transmit as soon as power is restored, and connection is re-established. A flashing envelope indicates that data are in the process of being transmitted.


To use Walk Away mode – connect the AC power adapter to the Analyzer and a power source

Step 6D: Starting Walk Away mode

- Turn the Analyzer on by pressing the blue power button once.


The display window will briefly show "SCAN CONFIG BARCODE." This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration instructions. Ignore this message and postpone this process when an assay is awaiting analysis.

- When the display window reads: "INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE"
 - Double-click** the blue power button.

Step 7D: Using the Bar Code scanner


- Follow the prompts on the display window to complete any required barcode scans of:
 - OPERATOR ID
 - SPECIMEN ID and/or
 - KIT LOT NUMBER
 In accordance with site requirements and Analyzer settings.

- Prompts for each scanning step appear in the display window for only 30 seconds. Failure to complete scans during that time will cause the Analyzer to default to the beginning of step 6D. To restart this step, double-click the power button.
- Move the barcode slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.
- The Analyzer can record the Kit Lot Number in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user. BD recommends against the use of expired materials.

Step 8D: Add the specimen to the test device

- When the display window reads: "ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY":
 - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System RSV device sample well).
 - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System RSV device.
- NOTE: Squeezing the tube close to the tip may cause leakage.**

Step 9D: Start the development and reading sequence


- Immediately insert the test device into the slot on the right side of the Analyzer.

The test device must remain horizontal to prevent spilling the specimen out of the sample well.

- "DO NOT DISTURB TEST IN PROGRESS" appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- A countdown timer in the display window shows the remaining analysis time.

Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.

- When analysis is complete, a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed. **If a printer is not connected, record the result before removing the assay device.**

ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 60 minutes (when AC power adapter is connected).

Step 10D: Remove the test device

- Pull the device out. The display will show INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE to indicate the Analyzer is ready to perform another test. Note that the Analyzer returns to **Analyze Now** mode at the conclusion of each read sequence.

If the Veritor Plus Analyzer is connected to an LIS, a steady ENVELOPE symbol will appear to indicate that results are awaiting transmission. In the event that a network connection is not detected while the ENVELOPE symbol is still displayed, the Analyzer will queue all untransmitted results and attempt to transmit them when reconnected. If it is powered off during this time, it will attempt to transmit as soon as power is restored, and connection is re-established. A flashing envelope indicates that data are in the process of being transmitted.

OPTIONAL TEST PROCEDURE: Use this procedure to test for both RSV and INFLUENZA A+B using a single NP wash, aspirate or swab specimen in transport medium.

Note: The BD Veritor™ System for Rapid Detection of Flu A+B (Cat. No. 256041) is required for this procedure in addition to the BD Veritor™ System for Rapid Detection of RSV (Cat. No. 256042).

IMPORTANT NOTE: This optional procedure allows for use of the remaining processed sample from Step 5 above to test additionally for Flu A+B. The sample to be tested must be from a patient less than 20 years of age as indicated in the **Indications for Use** on page 2 of this Package Insert. **THE PROCESSED SAMPLE SHOULD BE TESTED WITHIN 15 MINUTES.**

Allow specimens that have been stored at 2–8 °C to warm to room temperature prior to use. Thoroughly mix all specimens prior to removal of an aliquot for processing. Do not centrifuge specimens.

1. Collect specimen from the patient and follow Steps 1–5 of the test procedure above to prepare the sample for testing.
2. Using the sample from Step 5, continue the test procedure using the test device for Flu A+B, and the same workflow configuration used to obtain the RSV result.
3. Refer to the product insert for BD Veritor System for Rapid Detection of Flu A+B, (Cat. No. 256041) for the test procedure and full description of the BD Veritor Flu test. Follow the Instructions in the insert and the Instrument on screen prompts to complete the test procedure and obtain results. Refer to the product insert for the BD Veritor System Flu A+B Kit for result interpretation.

INTERPRETATION OF RESULTS

The BD Veritor System Instrument (purchased separately) must be used for interpretation of all test results. Operators should not attempt to interpret assay results directly from the test strip contained within the BD Veritor System RSV assay device.

Display	Interpretation
RSV: +	Positive Test for RSV (RSV antigen present)
RSV: -	Negative Test for RSV (no antigen detected)
CONTROL INVALID	Test Invalid. Repeat the Test.

Invalid Test – If the test is invalid, the BD Veritor System Instrument will display “CONTROL INVALID” and the test or control must then be repeated. If the “CONTROL INVALID” reading recurs, contact BD Technical Support.

REPORTING OF RESULTS

- Positive Test** Positive for the presence of RSV antigen. A positive result may occur in the absence of viable virus.
- Negative Test** Negative for the presence of RSV antigen. Infection due to RSV cannot be ruled out because the antigen present in the sample may be below the detection limit of the test. In the U.S., a negative test is presumptive and it is recommended that these results be confirmed by viral cell culture or an FDA-cleared RSV molecular assay.
- Control Invalid** Do not report results. Repeat the test.

QUALITY CONTROL:

To utilize the Analyzer’s QC documentation capability, specimen barcode scanning must be enabled on an Analyzer equipped with the BD Veritor InfoScan module. Please refer to the Analyzer *Instructions for Use*, section 4, to choose or change this configuration.

- Each BD Veritor System RSV device contains both positive and negative internal/procedural controls:
1. The internal positive control validates the immunological integrity of the device, proper reagent function, and assures correct test procedure.
 2. The membrane area surrounding test lines functions as a background check on the assay device.

The BD Veritor System Instrument evaluates the positive and negative internal/procedural controls after insertion of each test device. The BD Veritor System Instrument prompts the operator if a quality issue occurs during assay analysis. Failure of the internal/procedural controls will generate an invalid test result. **NOTE:** The internal controls do not assess proper sample collection technique.

External Positive and Negative Controls:

RSV positive and RSV negative control swabs are supplied with each kit. These controls provide additional quality control material to assess that the test reagents and the BD Veritor System Instrument perform as expected. Prepare kit control swabs and test using the same procedure (either **Analyze Now** or **Walk Away** mode) as used for patient specimen swabs. When using the barcode scanning feature to document QC procedures, scan the barcode on the control swab packaging when prompted for a Specimen ID.

Your laboratory’s standard Quality Control procedures and applicable local, state and/or federal regulations or accreditation requirements dictate the performance of external quality control procedures.

BD recommends external controls be run once for:

- each new kit lot,
- each new operator,
- each new shipment of test kits
- as required by internal quality control procedures and in accordance with local, state and federal regulations or accreditation requirements.

LIMITATIONS OF THE PROCEDURE

- Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.
- The contents of this kit are to be used for the qualitative detection of RSV antigens from NP wash, aspirate and swab in transport media specimens.
- The BD Veritor System for Rapid Detection of RSV is capable of detecting both viable and non-viable RSV particles. The BD Veritor System for Rapid Detection of RSV performance depends on antigen load and may not correlate with other diagnostic methods performed on the same specimen.
- Results from the BD Veritor System for Rapid Detection of RSV test should be correlated with the clinical history, epidemiological data and other data available to the clinician evaluating the patient.
- A false-negative test result may occur if the level of viral antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly; therefore, a negative test result does not eliminate the possibility of RSV infection, and should be confirmed by viral cell culture or an FDA-cleared RSV molecular assay.
- Positive test results do not rule out co-infections with other pathogens.
- Negative test results are not intended to rule in other non-RSV viral or bacterial infections.
- Positive and negative predictive values are highly dependent on prevalence rates. Positive test results are more likely to represent false positive results during periods of little/no RSV activity when disease prevalence is low. False negative test results are more likely during peak RSV activity when prevalence of disease is high.
- This device has been evaluated for use with human specimen material only.
- Monoclonal antibodies may fail to detect, or detect with less sensitivity, RSV viruses that have undergone minor amino acid changes in the target epitope region.
- The performance of this test has not been evaluated for use in patients without signs and symptoms of respiratory infection.
- The validity of the BD Veritor System for Rapid Detection of RSV test has not been proven for identification/confirmation of tissue culture isolates and should not be used in this capacity.
- Therapeutic anti-RSV monoclonal antibodies may interfere with the BD Veritor System for Rapid Detection of RSV.
- Performance characteristics have not been established for use with patients older than 20 years of age and for immunocompromised patients.

EXPECTED VALUES

The rate of positivity observed in RSV testing will vary depending on the method of specimen collection, handling/transport system employed, detection method utilized, time of year, age of the patient, geographic location and most importantly, local disease prevalence. In the 2011/2012 clinical trial, the overall prevalence of RSV as determined by viral cell culture for the nasopharyngeal swabs (NPS) in transport media was 24.5% (range of 5.6% to 31.8%). The overall prevalence of RSV as determined by viral cell culture for the nasopharyngeal washes and aspirates (NPWA) was 37.7% (range of 10.5% to 49.6%).

PERFORMANCE CHARACTERISTICS

Explanation of Terms

P: Positive

N: Negative

C.I.: Confidence Interval

Clinical Performance:

Performance characteristics for the BD Veritor System for Rapid Detection of RSV test were established in multi-center clinical studies conducted at five U.S. trial sites during the 2011–2012 respiratory season. A total of 1,174 prospectively collected specimens received in the laboratory with an order for respiratory virus testing were enrolled in the study, of which, 26 were noncompliant with the study protocol and one was noncompliant on the viral cell culture reference testing level. Removal of these specimens yields a total of 1,147 specimens. One additional specimen had a final undetermined viral cell culture reference result which could not be verified. Removal of this specimen results in a total of 1,146 specimens. A total of 1,146 were evaluated using the BD Veritor System for Rapid Detection of RSV test and viral cell culture. The prospective specimens consisted of 440 NPWA and 706 NPS in transport media from symptomatic patients. 44.3% of the samples were from females and 55.7% from males. 80% of patients were 2 years and under.

The performance of the BD Veritor System for Rapid Detection of RSV test was compared to an FDA cleared D³ Duet™ DFA on R-Mix cell culture and is presented in the following tables.

Table 1. Summary of the performance of the BD Veritor System for Rapid Detection of RSV Test compared to viral cell culture by specimen type, all sites.

		Viral Cell Culture		
Specimen Type	BD Veritor RSV	P	N	Total
NPS	P	153	9*	162
	N	20	524	544
	Total	173	533	706
Reference Method: Viral Cell Culture Sensitivity: 88.4% (95% CI: 82.8–92.4%) Specificity: 98.3% (95% CI: 96.8–99.1%)				
NPWA	P	152	15**	167
	N	14	259	273
	Total	166	274	440
Reference Method: Viral Cell Culture Sensitivity: 91.6% (95% CI: 86.3–94.9%) Specificity: 94.5% (95% CI: 91.2–96.7%)				

*Of the 9 BD Veritor RSV Positive, Viral Cell Culture negative specimens, 6 were positive by FDA cleared Prodesse ProFlu+ molecular assay.

**Of the 15 BD Veritor RSV Positive, Viral Cell Culture negative specimens, 8 were positive by FDA cleared Prodesse ProFlu+ molecular assay.

Reproducibility

The reproducibility of the BD Veritor System for Rapid Detection of RSV test was evaluated at three clinical laboratory sites. The reproducibility panel was composed of 12 simulated RSV samples. These included moderate positive samples, low positive samples (near the assay limit of detection), high negative samples (i.e., containing very low concentrations of virus) and negative samples. The panel was tested by two operators at each site for five consecutive days. The results are summarized below.

Reproducibility Results – Percent of RSV Positives				
Sample	Site 1	Site 2	Site 3	Total
High negative RSV	0.0% (0/30) (95% CI: 0–11.3%)	3.3% (1/30) (95% CI: 0.6–16.7%)	3.3% (1/30) (95% CI: 0.6–16.7%)	2.2% (2/90) (95% CI: 0.6–7.7%)
Low positive RSV	93.3% (28/30) (95% CI: 78.7–98.2%)	76.7% (23/30) (95% CI: 59.1–88.2%)	93.3% (28/30) (95% CI: 78.7–98.2%)	87.8% (79/90) (95% CI: 79.4–93%)
Moderate positive RSV	100% (30/30) (95% CI: 88.6–100%)	100% (30/30) (95% CI: 88.6–100%)	100% (30/30) (95% CI: 88.6–100%)	100% (90/90) (95% CI: 95.9–100%)
Negative	0.0% (0/30) (95% CI: 0–11.3%)	0.0% (0/30) (95% CI: 0–11.3%)	0.0% (0/30) (95% CI: 0–11.3%)	0.0% (0/90) (95% CI: 0–4.1%)

Analytical Studies

Analytical Sensitivity (Limit of Detection)

The limit of detection (LOD) for the BD Veritor System for Rapid Detection of RSV test was established for the following RSV strains The LOD for each strain represents the lowest concentration producing a positivity rate of ≥95% based on testing 60 to 80 replicates.

Viral Strain	Calculated LOD (TCID ₅₀ /mL)	No. Positive / Total	% Positive
VR-26 (Long Subgroup A)	1.43 x 10 ⁵	57/60	95.0
VR-955 (9320 subgroup B)	3.98 x 10 ⁴	57/60	95.0
VR-1540 (A-2)	1.94 x 10 ³	59/60	98.3
VR-1580 (Washington subgroup B)	1.08 x 10 ⁴	58/60	96.7
VR-1400 (Wild Type subgroup B)	2.96 x 10 ³	76/80	95.0

TCID₅₀/mL = 50% Tissue Culture Infectious Dose

Analytical Specificity (Cross Reactivity)

The BD Veritor System for Rapid Detection of RSV test was evaluated with bacteria and yeast at a target concentration of approximately 10^6 CFU/mL (CFU – Colony Forming Units) with the exception of *Fusobacterium nucleatum* which was tested at 1.5×10^6 . The viruses were evaluated at concentrations of 10^3 TCID₅₀/mL or greater. Of the microorganisms tested, none showed cross-reactivity in the RSV test.

<i>Bacteroides fragilis</i>	<i>Neisseria mucosa</i>	<i>Veillonella parvula</i>
<i>Bordetella pertussis</i>	<i>Neisseria</i> sp. (<i>Neisseria perflauis</i>)	Adenovirus, type 1
<i>Candida albicans</i>	<i>Neisseria subflava</i>	Adenovirus, type 7
<i>Chlamydia pneumoniae</i>	<i>Peptostreptococcus anaerobius</i>	Cytomegalovirus
<i>Corynebacterium diphtherium</i>	<i>Porphyromonas asaccharolyticus</i>	Enterovirus
<i>Escherichia coli</i>	<i>Prevotella oralis</i>	HSV Type 1
<i>Fusobacterium nucleatum</i>	<i>Propionibacterium acnes</i>	Human Coronavirus OC43
<i>Haemophilus influenzae</i>	<i>Proteus mirabilis</i>	Human metapneumovirus (HMPV-27 A2)
<i>Haemophilus parainfluenzae</i>	<i>Pseudomonas aeruginosa</i>	Human Parainfluenza
<i>Kingella kingae</i>	<i>Serratia marcescens</i>	Influenza A/Brisbane/10/2007 H3N2
<i>Klebsiella pneumoniae</i>	<i>Staphylococcus aureus</i>	Influenza A/California/7/2009 H1N1
<i>Lactobacillus</i> sp.	<i>Staphylococcus epidermidis</i>	Influenza A/Victoria/3/75 H3N2
<i>Legionella</i> sp.	<i>Streptococcus mutans</i>	Influenza B/Brisbane/60/2008
<i>Moraxella catarrhalis</i>	<i>Streptococcus pneumoniae</i>	Influenza B/Florida/4/2006
<i>Mycobacterium tuberculosis</i>	<i>Streptococcus pyogenes</i>	Influenza B/Lee/40
<i>Mycoplasma pneumoniae</i>	<i>Streptococcus</i> sp. Group C	Measles virus
<i>Neisseria gonorrhoeae</i>	<i>Streptococcus</i> sp. Group G	Mumps virus
<i>Neisseria meningitidis</i>	<i>Streptococcus salivarius</i>	Rhinovirus

Interfering Substances

Various substances were evaluated with the BD Veritor System for Rapid Detection of RSV test. These substances included whole blood (2%) and various medications. No interference was noted with this assay for any of the substances at the concentrations tested.

Substance	Concentration	Substance	Concentration
4-Acetamidophenol	10 mg/mL	Ibuprofen	10 mg/mL
Acetylsalicylic acid	20 mg/mL	Loratidine	100 ng/mL
Albuterol	0.083 mg/mL	Menthol Throat Lozenges	10 mg/mL
Amantadine Hydrochloride	500 ng/mL	Mometasone	500 ng/mL
Ayr Saline Nasal Gel	10 mg/mL	Mupirocin	500 ng/mL
Beclomethasone	500 ng/mL	Oseltamivir	500 ng/mL
Budesonide	500 ng/mL	Oxymetazoline	0.05 mg/mL
Chlorpheniramine maleate	5 mg/mL	Phenylephrine	1 mg/mL
Dexamethasone	10 mg/mL	Pseudoephedrine HCl	20 mg/mL
Dextromethorphan	10 mg/mL	Purified Mucin Protein	1 mg/mL
Diphenhydramine HCl	5 mg/mL	Ribavirin	500 ng/mL
Fexofenadine	500 ng/mL	Rimantadine	500 ng/mL
FluMist™	1%	Synagis	4 µg/mL
Flunisolide	500 ng/mL	Tobramycin	500 ng/mL
Fluticasone	500 ng/mL	Triamcinolone	500 ng/mL
Four OTC nasal sprays	10%	Two OTC mouthwashes	5%
Four OTC throat drops	12.5%	Whole Blood	2%
Guaiacol Glyceryl Ether	20 mg/mL	Zanamivir	1 mg/mL
Homeopathic Allergy Medicine	10 mg/mL		

AVAILABILITY

Cat. No.	Description
256041	BD Veritor™ System for Rapid Detection of Flu A+B, 30 tests
256042	BD Veritor™ System for Rapid Detection of RSV, 30 tests
256061	BD Veritor™ System RSV Control Swab Set, 10 pairs of swabs
256066	BD Veritor™ Plus Analyzer
256068	BD Veritor™ System InfoScan Module
443907	USB Printer Cable for BD Veritor™ Analyzer

To network a BD Veritor Plus Analyzer to an LIS, contact BD Technical Services for details.

REFERENCES

1. Hall, C.B, Weinberg, G.A., Iwane, M.K., et al. 2009. The Burden of Respiratory Syncytial Virus Infection in Young Children. *N Engl J Med* 2009; **360**:588–598.
2. Nair, H., Nokes, D.J., Gessner, B.D., et al., 2010. Global burden of acute lower respiratory infections due to respiratory syncytial virus in young children: a systematic review and meta-analysis. *Lancet* **375** (9725), 1545–1555.
3. Falsey, A.R. and E. E. Walsh. 2000. Respiratory syncytial virus infection in adults. *Clin. Microbiol. Rev.* **13**: 371–384.
4. Murata Y. and A. R. Falsey. 2006. RSV Infection in Elderly Adults, In: Patricia Cane, Editor(s), *Perspectives in Medical Virology*, Elsevier, Volume 14, Pages 163–182.
5. Crowcroft, N.S., F. Cutts and M.C. Zambon. 1999. Respiratory syncytial virus: an underestimated cause of respiratory infection, with prospects for a vaccine. *Commun Dis Public Health.* **2**: 234–241.
6. Thompson WW, Shay DK, Weintraub E, et al. 2003. Mortality associated with influenza and respiratory syncytial virus in the United States. *JAMA.* 2003. **289**(2):179–186.
7. Henrickson KJ and C.B.Hall, Diagnostic assays for respiratory syncytial virus disease. 2007. *Pediatr Infect Dis J.* Nov; **26**(11 Suppl):S36–40.
8. Therese Popow-Kraupp, T and J. H. Aberle. 2011. Diagnosis of Respiratory Syncytial Virus Infection. *Open Microbiol J.*; **5**:128–134. Published online 2011 December 30.
9. Barenfanger, J., C. Drake, N. Leon, T. Mueller and T. Troutt. 2000. Clinical and financial benefits of rapid detection of respiratory viruses: an outcomes study. *J. Clin. Microbiol.* **38**:2824–2828.
10. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed., CLSI, Wayne, PA.
11. Garner, J. S. 1996. Hospital Infection Control Practices Advisory Committee, U. S. Department of Health and Human Services, Center for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* **17**:53–80.
12. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
13. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021–0045.

Technical Information: In the United States, contact BD Technical Service and Support at 1.800.638.8663 or bd.com.