

## A Unique Design For greater comfort and confidence in use

BD Flashback needle  
gives you instantaneous  
confirmation of the needle  
placement, adding greater  
comfort in use.

- Avoids needle probing and the multiplication of attempts to draw blood from a patient.
- Ensures easy training of new staff or inexperienced personnel, and full utilization of closed-system collection.

Now, the technique of closed-system blood collection is easier for everyone.






## Reliable results start with reliable samples

### Closed blood collection system:

- Better specimen quality, less needs for redraws
- Less manipulations
- Less exposure to blood
- Faster procedures

### BD Vacutainer® Flashback needle

Ref.	NPC	Gauge	Length	Colour Coding	Packaging
301746	KFK361	21G - 0.8 MM	1" - 25 MM		50/box 1000/case
301747	KFK362	22G - 0.6 mm	1" - 25 mm		50/box 1000/case
301748	KFK364	22G - 0.6 mm	1,25" - 32 mm		50/box 1000/case

BD Vacutainer® Flashback needle ideally complements our range of blood collection devices and accessories. Ask your BD representative for more information.




### BD Diagnostics, Preanalytical Systems

East Europe, Middle East & Africa  
Arif & Bintok Building  
1st Floor, Zabeel Road  
Karama, PO BOX 52279  
Dubai, United Arab Emirates  
Tel: +971 4 3379525  
Fax: +971 4 3379551  
ema\_pas@europa.bd.com  
[www.bd.com](http://www.bd.com)




# BD Flashback Needle


Evacuated blood collection  
made easy




Twist and remove white shield




Screw on holder.  
Prepare venipuncture site.  
Twist and pull needle shield off.



Perform venipuncture. Observe  
blood flashback, confirming  
needle placement in the vein.



Introduce tube into holder.  
Remove tube from holder when  
blood flow ceases. If more sam-  
ples are needed insert addi-  
tional tubes. Remove last full  
tube from holder before removal of  
needle from vein.



Discard the needle and holder  
together into a sharps con-  
tainer.

#### CAUTION:

All biologic samples should be handled carefully. Exposure to such samples may transmit HBV (Hepatitis), HIV (AIDS) or other infectious diseases. All samples and collection materials should be disposed of in accordance with the policies and procedures of your facility.

Customer Service

Tel: +44 (0) 1865 781699

