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Referanslar

1. Stankovic AK, DiLauri E. Quality improvements in the preanalytical phase: Focus on urine specimen workflow. *Clin Lab Med.* 2008;28(2):339-350.
2. Howanitz PJ, Saladino AJ, Dale, JC. Timeliness of urinalysis: A college of American pathologists Q-probes study of 346 small hospitals. *Arch Pathol Lab Med.* 1997;121(7):667-672.
3. American Diabetes Association. Standards of medical care in diabetes—2010. *Diabetes Care.* 2010;33(suppl 1):s11-s61.
4. American Diabetes Association Website. <http://www.diabetes.org/diabetes-basics/diabetes-statistics>. Accessed October 29, 2010.
5. Levey AS, Coresh J, Balk E, et al. National Kidney Foundation practice guidelines for chronic kidney disease: evaluation, classification, and stratification. *Ann Intern Med.* 2003;139(2):137-147.
6. National Kidney Foundation Website. <http://www.kidney.org/kidneydisease/ckd>. Accessed October 29, 2010.
7. Warren JW, Abrutyn E, Hebel R, et al. Guidelines for antimicrobial treatment of uncomplicated acute bacterial cystitis and acute pyelonephritis in women. *Clin Infect Dis.* 1999;29(4):745-758.
8. Mehnert-Kay SA. Diagnosis and management of uncomplicated urinary tract infections. *Am Fam Physician.* 2005;72(3):451-456.
9. National Institute of Diabetes and Digestive and Kidney Diseases Website. <http://kidney.niddk.nih.gov/kudiseases/pubs/kustats>. Accessed October 29, 2010.
10. Hooton TM, Bradley SF, Cardenas DD, et al. Diagnosis, prevention, and treatment of catheter-associated urinary tract infection in adults: 2009 international clinical practice guidelines from the Infectious Diseases Society of America. *Clin Infect Dis.* 2010;50(5):625-663.
11. Jacobsen SM, Stickler DJ, Mobley HLT, Shirtliff ME. Complicated catheter-associated urinary tract infections due to *Escherichia coli* and *Proteus mirabilis*. *Clin Microbiol Rev.* 2008;21(1):26-59.
12. Bekeris LG, Jones BA, Walsh MK, Wagar, EA. Urine culture contamination: A college of American pathologists Q-probes study of 127 laboratories. *Arch Pathol Lab Med.* 2008;132(6):913-917.
13. CLSI; CLSI GP-16 A3. Urinalysis: Approved Guideline. 3rd Ed., February 2010.
14. Basic Laboratory Procedures in Clinical Bacteriology. 2nd Ed., World Health Organization. 2003.
15. Levinson, W. Growth. In Review of Microbiology and Immunology. 10th ed. San Francisco, CA: McGraw-Hill Companies, Inc; 2008:15.
16. Data on file. Becton, Dickinson and Company: Franklin Lakes, NJ; 2010.
17. Weinstein MP. Clinical evaluation of a urine transport kit with lyophilized preservative for culture, urinalysis, and sediment microscopy. *Diagn Microbiol Infect Dis.* 1985;3(6):501-508.
18. Cronin M. Automated urinalysis technology improves efficiency and patient care. *MLO Med Labs Obs.* 2008;40(10):30,32.
19. Amacher K. Advanced urinalysis technology and lean management help a hospital lab improve productivity. *MLO Med Lab Obs.* 2006;38:33-34.
20. Simerville JA, Maxted WC, Pahira JJ. Urinalysis: a comprehensive review. *Am Fam Physician.* 2005;71(6):1153-1162.
21. Costelloe C, Metcalfe C, Lovering A, Mant D, Hay AD. Effect of antibiotic prescribing in primary care on antimicrobial resistance in individual patients: Systematic review and meta-analysis. *BMJ.* 2010;340:c2096.
22. Crowe D. CEO voicing caution despite Iris's urinalysis breakthrough. Los Angeles Business Journal. September 12, 2005.
23. Valenstein P and Meier F. Urine culture contamination: A college of American Pathologists Q-Probes study of contaminated urine cultures in 906 institutions. *Arch Pathol Lab Med.* 1998;122(2):123-129.

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