SUREDETECT™

Slide Preparation Buffer 10X

REF 090-11008-10

 500 mL

IVD

2-8°C

PRESENTATION

500 mL of Slide Preparation Buffer, 10x Concentrate

INTENDED USE

This product, after dilution, is to be used on cytology and histology specimens, which have been mounted on glass for target retrieval prior to immunocytochemistry (ICC) or immunohistochemistry (IHC) procedures. Cytology specimens should be liquid-based cytology specimens prepared on glass slides. Histology specimens should be formalin-fixed paraffin-embedded tissue sections mounted on glass slides.

PRINCIPLES OF PROCEDURE

Target retrieval for cytology specimens is achieved by immersing the prepared slide in the diluted Slide Preparation Buffer for a minimum of 1 hour to a maximum of 72 hours at room temperature. Histology slides, which are immersed in the diluted Slide Preparation Buffer, do not require incubation at room temperature. The slides, while still in the Slide Preparation Buffer, are then heated at 95°C for 15 minutes using the appropriate equipment. After cooling down, slides are ready for ICC or IHC staining procedures. Target retrieval is important in ICC and IHC procedures as a means of improving the sensitivity and specificity of primary antibodies.

REAGENT PREPARATION

Prior to use, this reagent must be diluted 1:10 with deionized (diH₂O) or distilled water (one part Slide Preparation Buffer to nine parts water).

WARNINGS AND PRECAUTIONS

Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. Refer to the Material Safety Data Sheet (MSDS) for additional information.

Take reasonable precautions when handling reagents. Solution is HOT during use and may cause burns to exposed areas of eyes and skin.

PROCEDURE


Recommended Cytology Procedure: Black and Decker Vegetable Steamer

1. Fill Coplin jar or other suitable container with a sufficient quantity of the diluted Slide Preparation Buffer.
2. Immerse cytology slides and incubate for a minimum of 1 hour to a maximum of 72 hours at room temperature.
3. Place 1 liter of room temperature deionized water in the Base of the steamer to the Hi Fill mark.
4. Place Steaming Bowl onto Base.
5. Place Coplin jar, containing slides, into the Steaming Bowl.
6. Place lid on Steaming Bowl and set timer for 45 minutes (30 minutes heating time, 15 minutes at 95°C).
7. Remove Lid and transfer Coplin jar, containing slides, to the bench. Allow slides to cool at room temperature for 20 minutes.
8. Rinse the slides with diH₂O and transfer to an appropriate buffer (TBS, PBS, etc.).
9. Proceed with staining.

Recommended Histology Procedure: Black and Decker Vegetable Steamer

1. Fill Coplin jar or other suitable container with a sufficient quantity of the diluted Slide Preparation Buffer.
2. Deparaffinize and rehydrate tissue sections.
3. Immerse tissue sections in the diluted Slide Preparation Buffer.
4. Place 1 liter of room temperature deionized water in the base of the steamer to the Hi Fill mark.
5. Place Steaming Bowl onto Base.
6. Place Coplin jar, containing slides, into the Steaming Bowl.
7. Place lid on Steaming Bowl and set timer for 45 minutes (30 minutes heating time, 15 minutes at 95°C).
8. Remove Lid and transfer Coplin jar, containing slides, to the bench. Allow slides to cool at room temperature for 20 minutes.
9. Rinse the slides with diH₂O and transfer to an appropriate buffer (TBS, PBS, etc.).
10. Proceed with staining.

GLOSSARY OF SYMBOLS

	Catalog number
	For <i>in vitro</i> diagnostic use
	Consult instructions for use
	Contains 500 mL
	Caution, consult accompanying document
	Storage Temperature Limitations
	Batch Code
	Use by YYYY-MM-DD OR YYYY-MM
	Manufacturer

TECHNICAL SUPPORT

USA
Telephone: 1-877-822-7771
Fax: 1-336-290-8333
E-mail: tripathtechsupport@bd.com

Europe
Telephone: +32 (0)53 720 673
Fax: +32 (0)53 720 678

TriPath Imaging, Inc.
780 Plantation Drive
Burlington, NC 27215 USA

EC REP

Benex Limited
Rineanna House
Shannon Free Zone
Shannon, County Clare
Ireland

SureDetect is a trademark of TriPath Imaging, Inc.

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SureDetect är ett varumärke som tillhör TriPath Imaging, Inc.