# SUREDETECT™

## SiHa Cell Controls

REF 005-11012-10

8mL / 75 tests

IVD

2°C → 8°C

### INTENDED USE

The SureDetect™ SiHa Cell Control may be used in immunoassays that require an epithelial-based cell line as an assay staining control. This product is to be used in slide-based assays to allow the user to verify the staining procedure.

### SUMMARY AND EXPLANATION

The SureDetect™ SiHa Cell Control was derived from a patient with a grade II squamous cell carcinoma of the cervix (ATCC Number: HTB-35™).<sup>1-5</sup> The cell line is reported to contain an integrated human papillomavirus type 16 genome (HPV-16, 1 to 2 copies per cell). The product is designed to be used in slide-based assays to allow the user to verify the staining procedure. Light microscopy is used to detect the staining reaction.

### PRINCIPLES OF PROCEDURE

The SiHa Cell Control is packaged in a dropper-tip bottle and designed to deliver a consistent number of cells per drop. The cells have been gently fixed in alcohol-based preservative fluid. Three (3) drops per slide can be processed in a similar manner as clinical specimens in immunocytochemistry assays.

### REAGENT PROVIDED

8 mL of SureDetect™ SiHa Cell Control supplied in a buffered ethanol-based preservative fluid.

SiHa Cell Control is supplied ready to use.

### WARNINGS AND PRECAUTIONS

The product contains material of human origin and should be handled as a potential carrier and transmitter of disease and disposed of with proper precautions. Never pipette by mouth.

The cell line has NOT been screened for Hepatitis B, human immunodeficiency viruses or other adventitious agents.

Wear appropriate Personal Protective Equipment to avoid contact with skin and mucous membranes. If product comes in contact with sensitive areas, wash with copious amounts of water. Refer to the Material Safety Data Sheet for additional product safety information.

Universal Precautions should be followed when handling this product.

The product contains an aqueous solution of denatured ethanol. The mixture contains small amounts of methanol and isopropanol. Do not ingest.

Avoid microbial contamination of product to avoid non-specific staining.

Dispose of product according to local, state and federal regulations.

### STORAGE INSTRUCTIONS

Store at 2-8°C upon receipt and after each use. Do not freeze.

Do not use the product beyond the expiration date for the prescribed storage method.

There are no definitive signs to indicate instability of this product. Contact TriPath Technical Services if there is an indication of product instability.

### PROCEDURE

1. Allow SiHa Cell Control bottle to come to room temperature.
2. Gently invert the SiHa Cell Control bottle several times and apply 3 drops to the appropriate region of a control slide. The use of a hydrophobic barrier may aid in the correct localization of the SiHa cells on the slide.
3. Allow the cells to settle for up to 10 minutes. Do not allow the SiHa cells to dry out at any time.
4. Place the control slide into the appropriate pretreatment buffer and process as per standard test samples.

### REFERENCES

1. Baker CC, et al. Structural and transcriptional analysis of human papillomavirus type 16 sequences in cervical carcinoma cell lines. *J. Virol.* 61: 962-971, 1987.
2. Pater MM, Pater A. Human papillomavirus types 16 and 18 sequences in carcinoma cell lines of the cervix. *Virology* 145: 313-318, 1985.
3. Yee C, et al. Presence and expression of human papillomavirus sequences in human cervical carcinoma cell lines. *Am. J. Pathol.* 119: 361-366, 1985.
4. Friedl F, et al. Studies on a new human cell line (SiHa) derived from carcinoma of uterus. I. Its establishment and morphology. *Proc. Soc. Exp. Biol. Med.* 135: 543-545, 1970.
5. Scheffner M, et al. The state of the p53 and retinoblastoma genes in human cervical carcinoma cell lines. *Proc. Natl. Acad. Sci. USA* 88: 5523-5527, 1991.

# SUREDETECT™

SiHa Cell Controls

## GLOSSARY OF SYMBOLS

	Catalog number
	For <i>in vitro</i> diagnostic use
	Consult instructions for use
	Contains / Sufficient for
	Caution, consult accompanying documents
	Biological risk
	Storage Temperature Limitations
	Batch Code
	Use by YYYY-MM-DD OR YYYY-MM
	Manufacturer

## TECHNICAL SUPPORT

USA

Telephone: 1-877-822-7771

Fax: 1-336-290-8333

E-mail: [tripathtechsupport@bd.com](mailto:tripathtechsupport@bd.com)

Europe

Telephone: +32 (0)53 720 673

Fax: +32 (0)53 720 678

TriPath Imaging, Inc.  
780 Plantation Drive  
Burlington, NC 27215 USA

**EC REP**

Benex Limited  
Rineanna House  
Shannon Free Zone  
Shannon, County Clare  
Ireland

SureDetect™ is a trademark of TriPath Imaging, Inc.

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
### REFERENCIAS

1. Baker CC, et al. Structural and transcriptional analysis of human papillomavirus type 16 sequences in cervical carcinoma cell lines. *J. Virol.* 61: 962-971, 1987.
2. Pater MM, Pater A. Human papillomavirus types 16 and 18 sequences in carcinoma cell lines of the cervix. *Virology* 145: 313-318, 1985.
3. Yee C, et al. Presence and expression of human papillomavirus sequences in human cervical carcinoma cell lines. *Am. J. Pathol.* 119: 361-366, 1985.
4. Friedl F, et al. Studies on a new human cell line (SiHa) derived from carcinoma of uterus. I. Its establishment and morphology. *Proc. Soc. Exp. Biol. Med.* 135: 543-545, 1970.
5. Scheffner M, et al. The state of the p53 and retinoblastoma genes in human cervical carcinoma cell lines. *Proc. Natl. Acad. Sci. EE. UU.* 88: 5523-5527, 1991.

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1. Baker CC, et al. Structural and transcriptional analysis of human papillomavirus type 16 sequences in cervical carcinoma cell lines. J. Virol. 61: 962-971, 1987.
2. Pater MM, Pater A. Human papillomavirus types 16 and 18 sequences in carcinoma cell lines of the cervix. Virology 145: 313-318, 1985.
3. Yee C, et al. Presence and expression of human papillomavirus sequences in human cervical carcinoma cell lines. Am. J. Pathol. 119: 361-366, 1985.
4. Friedl F, et al. Studies on a new human cell line (SiHa) derived from carcinoma of uterus. I. Its establishment and morphology. Proc. Soc. Exp. Biol. Med. 135: 543-545, 1970.
5. Scheffner M, et al. The state of the p53 and retinoblastoma genes in human cervical carcinoma cell lines. Proc. Natl. Acad. Sci. USA 88: 5523-5527, 1991.

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SureDetect™ är ett varumärke som tillhör TriPath Imaging, Inc.