## MANUAL METHOD

### A cell preparation process for gynecologic applications

REF 490529

REF 490635

#### INTENDED USE

##### For In Vitro Diagnostic Use

The Manual Method is a method for producing liquid-based cell preparations (LBPs). The Manual Method is intended as a replacement for the conventional Pap smear preparation method for use in cervical cancer screening.

SUREPATH® Preservative Fluid is an appropriate collection and transportation medium for gynecologic specimens tested with BD ProbeTec™ *Chlamydia trachomatis* (CT) Q<sup>+</sup> and *Neisseria gonorrhoeae* (GC) Q<sup>+</sup> Amplified DNA Assays. Refer to the assay package inserts for instructions on using SUREPATH® Preservative Fluid to prepare specimens for use with these assays.

#### SUMMARY AND EXPLANATION

Cervical cytology screening by the Papanicolaou (Pap) method involves the microscopic examination of cell samples that have been taken primarily from the ecto- and endocervix, smeared on glass slides and stained using the Pap procedure.<sup>1,2,3</sup> Cervical cytology screening with the Pap smear has decreased the mortality rates of invasive cervical carcinoma by 50 to 70 percent.<sup>4</sup> Because cervical cytology is a screening test, abnormal findings must be confirmed histologically.

Specimen collection and preparation are extremely important to accuracy in Pap smears. Randomization or uniform sub-sampling is essential for complete accuracy. The conventional Pap smear technique does not provide for mixing of the sample prior to slide preparation. Due to the entanglement of cells in mucus on the sampling device, the cells actually transferred to the slide may not be representative of the total population collected. The cells are transferred to the slide in relation to where they happen to be on the sampling device. Many cells are left on the device.<sup>5</sup>

The non-homogeneity of a typical cervical sample can make conventional smears difficult to prepare, screen and interpret. Large areas of the conventional slide are often covered with debris, inflammatory cells and sheets of epithelial cells that can obscure valuable diagnostic material. In addition, if the smear is not fixed immediately after preparation, cellular morphology can be distorted as the smear dries (air-drying artifact).

The Manual Method is a method for converting a liquid suspension of a cervical sample into a consistently stained, homogeneous SUREPATH® slide while maintaining diagnostic cell clusters.<sup>6,7,8,9</sup> The process includes cell preservation, randomization, enrichment of diagnostic material, pipetting, and sedimentation to create a cellular preparation. The result of the preparation process is a SUREPATH® slide for use in routine cytology screening and categorization as defined by the Bethesda System.<sup>10</sup>

#### PRINCIPLES OF THE PROCEDURE

The Manual Method is a procedure for the preparation of LBPs of cervical cells. Gynecologic specimens are collected by qualified medical personnel using broom-like devices (e.g. Cervex-Brush®) or combination plastic spatula and endocervical brush devices (e.g. Cytobrush® Plus GT and Pap-Perfect® spatula, MedScand (USA), Inc.) with detachable heads. The head of the brush is removed from the handle and placed into a vial of SUREPATH® Preservative Fluid. The vial is capped, labeled and sent with appropriate paperwork to the laboratory for processing.

In the laboratory, the preserved sample is mixed by vortexing, and transferred into a tube containing PREPSTAIN® Density Reagent. An enrichment step, consisting of centrifugal sedimentation through PREPSTAIN® Density Reagent, partially removes non-diagnostic debris and excess inflammatory cells from the sample. After centrifugation, the tube containing the enriched cellular component is reconstituted with D.I. water and the cellular material is re-suspended with a pipettor, using an aspirate / dispense sequence. The sample material is then transferred to a PREPSTAIN® Settling Chamber mounted on a SUREPATH® PreCoat slide. Gravity sedimentation occurs during a short incubation. Excess material is decanted. The SUREPATH® PreCoat slide is stained, cleared and coverslipped, with the cells presented in a circle, 13 mm in diameter. The SUREPATH® Slide is examined by trained cytotechnologists and pathologists in conjunction with other relevant patient information.

#### LIMITATIONS

- Gynecologic specimens for preparation using the Manual Method should be collected using a broom-type device according to the standard collection procedure provided by the manufacturer.
- The production and evaluation of TRIPATH Imaging®, Inc liquid-based preparations should be performed only by personnel who have been trained by TRIPATH or others authorized by TRIPATH to provide such training.
- Proper performance of the device requires the use of only those supplies supported by TRIPATH Imaging®, Inc, or recommended by TRIPATH Imaging®, Inc. Used supplies should be disposed of properly in accordance with institutional and governmental regulations.
- All supplies are for single use only and cannot be reused.
- A volume of 8.0 ± 0.5 mL of the sample collected in the SUREPATH® Preservative Fluid Collection Vial is required for the SUREPATH® LBC Test process.

#### WARNINGS

- SUREPATH® Preservative Fluid contains a dilute solution of denatured ethanol and is not intended for human consumption. The mixture contains small amounts of methanol and isopropanol which can be harmful and cause blindness if ingested.
- PREPSTAIN® Density Reagent contains sodium azide. Sodium azide may react with lead or copper plumbing to

form highly explosive metal azides. On disposal, flush with a large volume of water to prevent buildup of azide. For further information, refer to the Manual Guide issued by the Centers for Disease Control<sup>11</sup>.

#### PRECAUTIONS

- It is expected that good laboratory practices will be followed and that all procedures for use of the Manual Method will be strictly observed.
- All reagents are stable until stated expiration dates, provided recommended storage conditions are followed and maintained.
- Microbial contamination of reagents may give incorrect results.
- Substitution of other than SUREPATH® PreCoat slides may result in less than optimal results.
- Avoid splashing or generating aerosols. Use appropriate hand, eye and clothing protection.
- SUREPATH® Preservative Fluid is bactericidal and has been tested against: *Escherichia coli*, *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Candida albicans*, *Mycobacterium tuberculosis*, and *Aspergillus niger*. However, universal precautions for safe handling of biological fluids should be practiced at all times.

#### OPTIONAL ALIQUOT REMOVAL

Sufficient volume is available in the SUREPATH® Preservative Fluid Collection Vial to allow removal of up to 0.5 mL of homogeneous mixture of cells and fluid for ancillary testing, prior to the SUREPATH® Pap Test while still allowing sufficient volume for Pap testing.

While there is no evidence that removal of an aliquot from the SUREPATH® Preservative Fluid Collection Vial affects the quality of the specimen for cytology testing, rare instances of misallocation of pertinent diagnostic material may occur during this process. Healthcare providers may need to acquire a new specimen if the results do not correlate with the clinical history of the patient. Furthermore, cytology addresses different clinical questions than sexually transmitted disease (STD) testing; therefore, aliquot removal may not be suitable for all clinical situations. If necessary, a separate specimen may be collected for STD testing rather than taking an aliquot from the SUREPATH® Preservative Fluid Collection Vial.

Aliquot removal from low-cellularity specimens may leave insufficient material in the SUREPATH® Preservative Fluid Collection Vial for preparation of a satisfactory SUREPATH® Pap test.

Aliquot must be removed prior to processing the SUREPATH® Pap Test. Only one aliquot may be removed from the SUREPATH® Preservative Fluid Collection Vial prior to performing the PapTest, regardless of the volume of the aliquot.

## Procedure

1. In order to ensure a homogenous mixture, the SUREPATH® Preservative Fluid Collection Vial must be vortexed for 10-20 seconds and the 0.5 mL aliquot must be removed within one minute of vortexing.
2. A polypropylene aerosol barrier pipette tip that is sized appropriately for the volume being withdrawn must be used for aliquot removal. *Note:* Serological pipettes should not be used. Good laboratory practices must be followed to avoid introducing contaminants into the SUREPATH® Preservative Fluid collection vial or the aliquot. Aliquot removal should be performed in an appropriate location outside an area where amplification is performed.
3. Visually check the aliquot material in the pipette for evidence of large gross particulates or semi-solids. Evidence of such material encountered while withdrawing the aliquot material should prompt return of all the material to the specimen vial and disqualify the specimen for ancillary testing prior to performing the Pap test.
4. For instructions on processing the aliquot using the BD ProbeTec™ CT Q<sup>x</sup> and GC Q<sup>x</sup> Amplified DNA Assays, refer to the assay Package Inserts provided by the manufacturer.

## MATERIALS REQUIRED

### Materials Provided

- SUREPATH® Preservative Fluid Collection Vial
- PREPSTAIN® Density Reagent
- PREPSTAIN® Settling Chambers
- SUREPATH® PreCoat slides
- Centrifuge Tubes
- PREPSTAIN® Syringing Pipettes
- Aspirator Tips

### Materials Required But Not Provided

- Centrifuge
- Slide Racks
- Easy Aspirator (optional)
- Processing Tray (optional)
- Broom type sampling device or endocervical brush/plastic spatula with detachable head(s)
- Vortex Mixer
- Precision Pipettes with Disposable Tips
- Deionized Water (pH 7.5 to 8.5)
- Isopropanol and Reagent Grade Alcohol
- Staining Reagents
- Clearing Agent, Mounting Media, Coverslips

## STORAGE

- SUREPATH® Preservative Fluid without cytologic samples may be stored at room temperature (15° to 30° C) for up to

36 months from date of manufacture.

- The storage limit for SUREPATH® Preservative Fluid with cytologic samples is 6 months at refrigerated temperatures (2° to 10° C) or 4 weeks at room temperature (15° to 30° C).
- SUREPATH® Preservative Fluid containing cytologic sample intended for use with the BD ProbeTec™ CT Q<sup>x</sup> and GC Q<sup>x</sup> Amplified DNA Assays can be stored and transported for up to 30 days at 2° – 30° C prior to transfer to the Liquid-Based Cytology Specimen (LBC) Dilution Tubes for the BD ProbeTec™ Q<sup>x</sup> Amplified DNA Assays.

## PROCEDURES

1. After the sample is collected using a Rovers Cervex-Brush® or equivalent sampling device, the brush head is rinsed directly into the fluid, removed from the handle and dropped into a SUREPATH® Preservative Fluid vial. The vial is then tightly capped, labeled, and sent to the laboratory.
2. When sample vials have been accessioned into the lab, place each vial into the processing tray with a labeled centrifuge tube pre-filled with 4 ml of PREPSTAIN® Density Reagent and a labeled SUREPATH® PreCoat slide. PREPSTAIN® Density Reagent must be added to the tube before the sample is added or performance will be reduced.
3. Vigorously vortex each sample vial for 10 – 20 seconds. (Sufficient volume is available in the SUREPATH® Preservative Fluid Collection Vial to allow removal of up to 0.5 mL of homogenous mixture of cells and fluid for ancillary testing, while still allowing sufficient volume for Pap testing. The aliquoting may be performed after this vortexing step in the SUREPATH® LBC Test process.)
  - Use the PREPMATE™ Automated Accessory and PREPSTAIN® Syringing Pipettes to transfer the sample. See the PREPMATE™ Operators Manual for instructions. **Or**
  - Remove vial cap. Hold a SUREPATH® Preservative Vial in one hand, gently push a PREPSTAIN® Syringing Pipettes into the vial until it stops, pointing the end of the PREPSTAIN® Syringing Pipettes away from your face. Invert the vial / PREPSTAIN® Syringing Pipettes assembly into the appropriately numbered PREPSTAIN® Density Reagent tube. Allow all of the sample solution to completely drain from the PREPSTAIN® Syringing Pipettes before proceeding. **Or**
  - Remove vial cap. Slowly pour approximately 8 ml of specimen onto PREPSTAIN® Density Reagent.
4. Place the tubes into the centrifuge racks according to the placement sequence diagram in the Procedure Manual (each tube rack has a capacity of 12 tubes). Placement sequence is critical and must be balanced.
5. Balance the centrifuge tubes by adding Preservative Fluid if necessary, and centrifuge specimens for 2 minutes at 200 x g.
6. If the Easy Aspirator is used, turn on the Easy Aspirator

system and adjust pressure to 9 ± 2 in Hg. Place clean tips onto the Easy Aspirator aspirator.

7. Remove centrifuge tube racks from the centrifuge. Slowly lower the Easy Aspirator tips (or disposable transfer pipettes) into the centrifuge tubes to aspirate supernatant. The aspirate device should touch the top of the tubes at completion. Rinse Aspirator Tips between samples with water.
8. Centrifuge the tubes for 10 minutes at 800 x g to concentrate the diagnostic component into a cell pellet at the bottom of the tube.
9. Remove the tube rack from the centrifuge. Gently and quickly decant into the sink. Keeping rack inverted, blot the tubes carefully on absorbent paper making sure that the cell pellet stays in the tube.
10. Label a SUREPATH® PreCoat slide with each specimen number, being careful not to touch the surface of the SUREPATH® PreCoat slide. Place slides into Slide Rack and lock a PREPSTAIN® Settling Chamber onto each slide. The position of each numbered SUREPATH® PreCoat slide on the platter must match to the position of the corresponding centrifuge tube.
11. Add 4 ml of Deionized Water (pH 7.5-8.0) to each specimen tube and mix well by vortexing.
12. Working with one sample tube at a time, vortex tube and immediately transfer 800 µl of cell suspension into the correspondingly numbered PREPSTAIN® Settling Chamber/SUREPATH® PreCoat slide. Repeat for each sample.
13. Allow 10 minutes for full sedimentation to occur. After the sedimentation, gently invert the slide platter(s) over the sink to decant the remaining fluid and blot the excess liquid on absorbent paper.
14. Rinse each PREPSTAIN® Settling Chamber with 500 µl of denatured ethanol and decant. Repeat alcohol rinse and decant the remaining fluid and blot the excess liquid on absorbent paper, allowing them to remain inverted for at least 1 minute.
15. Remove the PREPSTAIN® Settling Chamber.
16. Stain and coverslip SUREPATH® Slides.

## RESULTS AND INTERPRETATION

- All diagnostic criteria currently utilized in cytology laboratories for conventional Pap smears are applicable to TRIPATH Imaging®, Inc. liquid-based preparations.
- Any abnormal or questionable screening observations should be referred to a pathologist for review and diagnosis. Any cellular morphologic changes are significant and should be noted.

## BIBLIOGRAPHY

1. Papanicolaou GN: A New Procedure for Staining Vaginal Smears. Science 1942; 95: 438-439.
2. King A, Clay K, Felmar EG, Heustis DG, Karns RM, Krahl P, Tench WD: The Papanicolaou Smear. West J Med 1992; 156: 202-204.
3. Mandelblatt J, Gopaul I, Wistreich M: Gynecological Care of Elderly Women: Another Look at Papanicolaou Smear Testing. J Am Med Assoc 1986; 256: 367-371.
4. Koss LG: The Papanicolaou Test for Cervical Cancer Detection: A Triumph and a Tragedy. J Am Med Assoc 1989; 261: 737-743.
5. Hutchinson ML, Isenstein LM, Goodman A, Hurley AA, Douglass KL, Mui KK, Patten FW, Zahniser DJ: Homogeneous Sampling Accounts for the Increased Diagnostic Accuracy Using the ThinPrep® Processor. Am J Clin Pathol 1994; 101: 215-219.
6. McGoogan E, Reith A: Would Monolayers Provide More Representative Samples and Improved Preparations for Cervical Screening? Overview and Evaluation of Systems Available. Acta Cytol 1996; 40: 107-119.
7. Bishop JW: Comparison of the CytoRich System with Conventional Cervical Cytology: Preliminary Data on 2,032 Cases from a Clinical Trial Site. Acta Cytol 1997; 41: 15-23.
8. Geyer JW, Hancock F, Carrico C, Kirkpatrick M: Preliminary Evaluation of CytoRich: An Improved Automated Cytology Preparation. Diagn Cytopathol 1993; 9: 417-422.
9. Grohs HK, Zahniser DJ, Geyer JW: Standardization of Specimen Preparation Through Mono/Thin-Layer Technology in Automated Cervical Cancer Screening. Edited by HK Grohs, OAN Husain. New York, Igaku-shoin, 1994, pp 176-185.
10. Kurman RJ, Solomon D: The Bethesda System for Reporting Cervical/Vaginal Cytologic Diagnoses. Springer-Verlag, 1994.
11. Manual Guide--Safety Management No. CDC-22, Decontamination of Laboratory Sink Drains to Remove Azide Salts, Center for Disease Control, Atlanta, Georgia, April 30, 1976.

## Technical Support

USA

Telephone: 1-877-822-7771

Fax: 1-336-290-8333

Europe

Telephone: +32 (0)53 720 673

Fax: +32 (0)53 720 678

TriPath Imaging, Inc.

780 Plantation Drive

Burlington, NC 27215 USA

Benex Limited

Rineanna House

Shannon Free Zone

Shannon, County Clare

Ireland

Australian Representative:

Becton Dickinson Pty Ltd.

4 Research Park Drive,

Macquarie University Research Park,

North Ryde,

NSW 2113 Australia

©2011 TRIPATH IMAGING, INC.

ALL RIGHTS RESERVED.

BD, BD Logo, and BD ProbeTec are trademarks of  
Becton, Dickerson and Company. © 2011 BD

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#### **BIBLIOGRAPHIE**

1. Papanicolaou GN: A New Procedure for Staining Vaginal Smears. Science 1942; 95: 438-439.
2. King A, Clay K, Felmar EG, Heustis DG, Karns RM, Krahl P, Tench WD: The Papanicolaou Smear. West J Med 1992; 156: 202-204.
3. Mandelblatt J, Gopaul I, Wistreich M: Gynecological Care of Elderly Women: Another Look at Papanicolaou Smear Testing. J Am Med Assoc 1986; 256: 367-371.
4. Koss LG: The Papanicolaou Test for Cervical Cancer Detection: A Triumph and a Tragedy. J Am Med Assoc 1989; 261: 737-743.
5. Hutchinson ML, Isenstein LM, Goodman A, Hurley AA, Douglass KL, Mui KK, Patten FW, Zahniser DJ: Homogeneous Sampling Accounts for the Increased Diagnostic Accuracy Using the ThinPrep® Processor. Am J Clin Pathol 1994; 101: 215-219.
6. McGoogan E, Reith A: Would Monolayers Provide More Representative Samples and Improved Preparations for Cervical Screening? Overview and Evaluation of Systems Available. Acta Cytol 1996; 40: 107-119.
7. Bishop JW: Comparison of the CytoRich System with Conventional Cervical Cytology: Preliminary Data on 2,032 Cases from a Clinical Trial Site. Acta Cytol 1997; 41: 15-23.
8. Geyer JW, Hancock F, Carrico C, Kirkpatrick M: Preliminary Evaluation of CytoRich: An Improved Automated Cytology Preparation. Diagn Cytopathol 1993; 9: 417-422.
9. Grohs HK, Zahniser DJ, Geyer JW: Standardization of Specimen Preparation Through Mono/Thin-Layer Technology in Automated Cervical Cancer Screening. Edited by HK Grohs, OAN Husain. New York, Igaku-shoin, 1994, pp 176-185.

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1. Papanicolaou GN: A New Procedure for Staining Vaginal Smears. *Science* 1942; 95: 438-439.
2. King A, Clay K, Felmar EG, Heustis DG, Karns RM, Krahl P, Tench WD: The Papanicolaou Smear. *West J Med* 1992; 156: 202-204.
3. Mandelblatt J, Gopaul I, Wistreich M: Gynecological Care of Elderly Women: Another Look at Papanicolaou Smear Testing. *J Am Med Assoc* 1986; 256: 367-371.
4. Koss LG: The Papanicolaou Test for Cervical Cancer Detection: A Triumph and a Tragedy. *J Am Med Assoc* 1989; 261: 737-743.
5. Hutchinson ML, Isenstein LM, Goodman A, Hurley AA, Douglass KL, Mui KK, Patten FW, Zahniser DJ: Homogeneous Sampling Accounts for the Increased Diagnostic Accuracy Using the ThinPrep® Processor. *Am J Clin Pathol* 1994; 101: 215-219.
6. McGoogan E, Reith A: Would Monolayers Provide More Representative Samples and Improved Preparations for Cervical Screening? Overview and Evaluation of Systems Available. *Acta Cytol* 1996; 40: 107-119.
7. Bishop JW: Comparison of the CytoRich System with Conventional Cervical Cytology: Preliminary Data on 2,032 Cases from a Clinical Trial Site. *Acta Cytol* 1997; 41: 15-23.
8. Geyer JW, Hancock F, Carrico C, Kirkpatrick M: Preliminary Evaluation of CytoRich: An Improved Automated Cytology Preparation. *Diagn Cytopathol* 1993; 9: 417-422.
9. Grohs HK, Zahniser DJ, Geyer JW: Standardization of Specimen Preparation Through Mono/Thin-Layer Technology in Automated Cervical Cancer Screening. Edited by HK Grohs, OAN Husain. New York, Igaku-shoin, 1994, pp 176-185.
10. Kurman RJ, Solomon D: The Bethesda System for Reporting Cervical/Vaginal Cytologic Diagnoses. Springer-Verlag, 1994.
11. Manual Guide--Safety Management No. CDC-22, Decontamination of Laboratory Sink Drains to Remove Azide Salts, Center for Disease Control, Atlanta, Georgia, April 30, 1976.

#### Technical Support

USA  
Telephone: 1-877-822-7771  
Fax: 1-336-290-8333

Europe  
Telephone: +32 (0)53 720 673  
Fax: +32 (0)53 720 678

TriPath Imaging, Inc.  
780 Plantation Drive  
Burlington, NC 27215 USA

Benex Limited  
Rineanna House  
Shannon Free Zone  
Shannon, County Clare  
Ireland

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1. Papanicolaou GN: A New Procedure for Staining Vaginal Smears. Science 1942; 95: 438-439.
2. King A, Clay K, Felmar EG, Heustis DG, Karns RM, Krahl P, Tench WD: The Papanicolaou Smear. West J Med 1992; 156: 202-204.
3. Mandelblatt J, Gopaul I, Wistreich M: Gynecological Care of Elderly Women: Another Look at Papanicolaou Smear Testing. J Am Med Assoc 1986; 256: 367-371.
4. Koss LG: The Papanicolaou Test for Cervical Cancer Detection: A Triumph and a Tragedy. J Am Med Assoc 1989; 261: 737-743.
5. Hutchinson ML, Isenstein LM, Goodman A, Hurley AA, Douglass KL, Mui KK, Patten FW, Zahniser DJ: Homogeneous Sampling Accounts for the Increased Diagnostic Accuracy Using the ThinPrep<sup>®</sup> Processor. Am J Clin Pathol 1994; 101: 215-219.
6. McGoogan E, Reith A: Would Monolayers Provide More Representative Samples and Improved Preparations for Cervical Screening? Overview and Evaluation of Systems Available. Acta Cytol 1996; 40: 107-119.
7. Bishop JW: Comparison of the CytoRich System with Conventional Cervical Cytology: Preliminary Data on 2,032 Cases from a Clinical Trial Site. Acta Cytol 1997; 41: 15-23.
8. Geyer JW, Hancock F, Carrico C, Kirkpatrick M: Preliminary Evaluation of CytoRich: An Improved Automated Cytology Preparation. Diagn Cytopathol 1993; 9: 417-422.
9. Grohs HK, Zahniser DJ, Geyer JW: Standardization of Specimen Preparation Through Mono/Thin-Layer Technology in Automated Cervical Cancer Screening. Edited by HK Grohs, OAN Husain. New York, Igaku-shoin, 1994, pp 176-185.
10. Kurman RJ, Solomon D: The Bethesda System for Reporting Cervical/Vaginal Cytologic Diagnoses. Springer-Verlag, 1994.
11. Manual Guide--Safety Management No. CDC-22, Decontamination of Laboratory Sink Drains to Remove Azide Salts, Center for Disease Control, Atlanta, Georgia, April 30, 1976.

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#### BIBLIOGRAFIE

1. Papanicolaou GN: A New Procedure for Staining Vaginal Smears. Science 1942; 95: 438-439.
2. King A, Clay K, Felmar EG, Heustis DG, Karns RM, Krahl P, Tench WD: The Papanicolaou Smear. West J Med 1992; 156: 202-204.
3. Mandelblatt J, Gopaul I, Wistreich M: Gynecological Care of Elderly Women: Another Look at Papanicolaou Smear Testing. J Am Med Assoc 1986; 256: 367-371.
4. Koss LG: The Papanicolaou Test for Cervical Cancer Detection: A Triumph and a Tragedy. J Am Med Assoc 1989; 261: 737-743.
5. Hutchinson ML, Isenstein LM, Goodman A, Hurley AA, Douglass KL, Mui KK, Patten FW, Zahniser DJ: Homogeneous Sampling Accounts for the Increased Diagnostic Accuracy Using the ThinPrep® Processor. Am J Clin Pathol 1994; 101: 215-219.
6. McGoogan E, Reith A: Would Monolayers Provide More Representative Samples and Improved Preparations for Cervical Screening? Overview and Evaluation of Systems Available. Acta Cytol 1996; 40: 107-119.
7. Bishop JW: Comparison of the CytoRich System with Conventional Cervical Cytology: Preliminary Data on 2,032 Cases from a Clinical Trial Site. Acta Cytol 1997; 41: 15-23.
8. Geyer JW, Hancock F, Carrico C, Kirkpatrick M: Preliminary Evaluation of CytoRich: An Improved Automated Cytology Preparation. Diagn Cytopathol 1993; 9: 417-422.
9. Grohs HK, Zahniser DJ, Geyer JW: Standardization of Specimen Preparation Through Mono/Thin-Layer Technology in Automated Cervical Cancer Screening. Edited by HK Grohs, OAN Husain. New York, Igaku-shoin, 1994, pp 176-185.
10. Kurman RJ, Solomon D: The Bethesda System for Reporting Cervical/Vaginal Cytologic Diagnoses. Springer-Verlag, 1994.
11. Manual Guide-Safety Management No. CDC-22, Decontamination of Laboratory Sink Drains to Remove Azide Salts, Center for Disease Control, Atlanta, Georgia, April 30, 1976.

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