BD Microtainer® MAP

Microtube for Automated Process with K2 EDTA

For Determinations of Hematology Parameters and Blood Lead*

REF 363706

Code	Language	Page	Code	Language	Page
EN	English	1 - 2	KK	Қазақша	30 - 31
BG	Български	3 - 4	KO	한국어	32
CS	Česky	5 - 6	LT	Lietuvių k.	33 - 34
DA	Dansk	7 - 8	LV	Latviski	35 - 36
DE	Deutsch	9 - 10	NL	Nederlands	37 - 38
EL	Ελληνικά	11 - 12	NO	Norsk	39 - 40
ES	Español	13 - 14	PL	Polski	41 - 42
ET	Eesti	15 - 16	PT	Português	43 - 44
FI	Suomi	17 - 18	RO	Română	45 - 46
FR	Français	19 - 20	RU	Русский	47 - 48
HR	Hrvatski	21 - 22	SK	Slovensky	49 - 50
HU	Magyar	23 - 24	SL	Slovenski jezik	51 - 52
IS	Íslenska	25 - 26	SV	Svenska	53 - 54
IT	Italiano	27 - 28	TR	Türkçe	55 - 56
JA	日本語	29	UK	Українська	57 - 58

Symbol and Mark Key 59

INTENDED USE:

BD Microtainer® MAP Microtube for Automated Process with K2 EDTA is used to collect, anticoagulate, transport and store skin puncture blood specimens for measurement of the following hematology parameters:

WBC, RBC, Hgb, HCT, MCV, MCH, MCHC, Platelets, 5-part WBC differential (Neutrophils, Lymphocytes, Monocytes, Eosinophils, Basophils), Reticulocytes and Whole Blood Lead testing.

SUMMARY AND EXPLANATION:

BD Microtainer® MAP consists of a specially designed 13x75 mm plastic tube with penetrable color-coded cap for automated processes. The tube has an integrated collector and markings to indicate fill volumes of 250 µL, 375 µL, and 500 µL. The lavender color cap indicates K2 EDTA, which is spray-dried onto the interior surface of the tube.

Reagent

BD Microtainer® MAP tube contains 1.0 mg of K2 EDTA, which is sufficient to anticoagulate 250 µL to 500 µL of skin puncture blood.

EDTA functions as an anticoagulant by chelating mono and divalent cations. Calcium is chelated by EDTA, which prevents its activity in the coagulation process. EDTA also reduces platelet activation, but does not eliminate it completely.

WARNINGS AND PRECAUTIONS:

1. Practice Standard Precautions, using gloves and appropriate apparel for protection from exposure to bloodborne pathogens.
2. Handle all biologic samples and blood collection products in accordance with policies and procedures of your facility. Obtain appropriate medical attention in the event of any exposure to biologic samples (e.g., through a puncture injury) as samples may transmit HBV (viral hepatitis), HIV (AIDS), or other infectious diseases. Discard all products contaminated with blood in biohazard containers approved for their disposal.
3. Do not use tubes if foreign matter is present.
4. Do not use the tubes after expiration date.

Storage

Store tubes at 4 - 40°C (39 to 104°F). Do not use tubes after their expiration dates. Avoid exposure to direct sunlight.

BLOOD COLLECTION AND HANDLING:

This entire package insert should be read and understood before performing skin puncture.

Blood collected in the BD Microtainer® MAP with K2 EDTA is stable at room temperature for up to 12 hours and will provide accurate and clinically useful test results for the hematology parameters listed above.

*Contains less than 1 ng lead per tube when evaluated using Inductively Coupled Plasma Mass Spectroscopy.

PROCEDURE:

Materials Provided:

BD Microtainer® MAP Microtube for Automated Process with K2 EDTA

Materials Required But Not Provided:

1. Gloves for the person performing the collections.
2. Alcohol swabs (lead-free if lead testing is performed).
3. Soap and water (for lead testing).
4. Dry, clean gauze pads.
5. BD Lancet appropriate for site and volume of blood required.
6. Sharps and contaminated materials disposal container.

Optional Materials:

1. Warming device.
2. Adhesive bandage. Avoid use of bandage with patients likely to place fingers or feet in their mouths, as ingestion/aspiration may occur.

Directions for Use:

1. Gloves must be worn during skin puncture blood collection.
2. Order of draw: Fill other BD Microtainer® Tubes as required following directions provided in their respective product circulars. NOTE: If capillary blood gas test is to be done, the capillary tubes used to collect samples must be filled before other samples are collected. Otherwise, the K2 EDTA samples for hematology testing should be the first tube collected.
3. Select puncture site, warm as appropriate.
4. Cleanse with alcohol and allow site to air-dry. Do not dry by wiping, as disinfection occurs during air-drying.
5. For lead testing, wash area with soap and water and dry thoroughly to minimize surface contamination of skin by environmental lead. Cleanse with lead-free alcohol swab. Institutional protocols for collection of specimens for blood lead should be followed.
6. Remove cap from BD MAP and place on a convenient lead-free surface.
7. Puncture skin with the appropriate BD Lancet, following instructions supplied.
8. Dispose of used lancet in an approved biohazard sharps container.
9. Wipe away first drop of blood with gauze. Hold BD MAP at an angle from surface of puncture site. Touch integrated collector end to drop of blood. After collecting 2 or 3 drops, blood will freely flow down the interior tube wall. **CAUTION:** "Milking" or "Scraping" of skin puncture site may cause hemolysis and adversely affect test result accuracy.
10. Fill BD MAP between 250 µL and 500 µL to achieve proper blood to additive ratio. Over or under filling of tube may result in clotting and/or erroneous test results.
11. A minimum blood volume of 375 µL is required for automated processing on most hematology systems. Refer to your instrument manufacturer's manual for minimum volume requirements.
12. Replace cap onto BD MAP. Push cap down to completely cover collar for a secure snap fit. Immediately gently invert 8 times after collection.
13. Reminder: Dispose of used BD Lancet into an approved biohazard sharps container. Dispose of any contaminated materials into appropriate container.

14. The specimen should be mixed thoroughly just before analysis in accordance with instrument manufacturers' instructions.
15. It is not necessary to remove the penetrable cap from the tube for automated processing. This minimizes potential exposure to bloodborne pathogens.

LIMITATIONS:

1. Users should determine the stability of the specimen as specimen stability can vary depending upon the reagents, system, instrumentation, storage conditions, and patient disease state.
2. Blood collected by skin puncture should only be used for lead screening. Elevated blood lead results obtained on capillary specimens should be considered presumptive and must be confirmed using venous blood. When lead testing is performed, follow the manufacturer's recommended procedure.
3. For lead, other analytes or alternate test methods, users must establish that test results meet their institution's criteria for clinically acceptable values. The user should consult their instrument manufacturer for expected specimen stability.
4. Some patients have shown sensitivity to EDTA, resulting in significant and irreversible platelet clumping. This clumping may also artificially elevate WBC counts.
5. It is the laboratory's ultimate responsibility to determine reference intervals for all analytes based upon the tubes used for sample collection by that laboratory. The clinical laboratory should establish or verify its reference ranges if changing specimen collection tube types and sizes. Changes could potentially affect analytical results from patient samples.

Technical Services

U. S. customers, please call BD Global Technical Services at 1.800.631.0174.

Canadian customers, please call 1.866.979.9408.

Customers outside the U.S., please contact your local BD sales consultant.

All customers: www.bd.com/vacutainer

REFERENCES:

1. CLSI. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition. CLSI document H04-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
2. Preventing Lead Poisoning in Young Children. A Statement by the Centers for Disease Control – October 1991, U.S. Department of Health and Human Services. p.41.

 Becton, Dickinson and Company, Franklin Lakes, NJ USA
 Becton, Dickinson and Company, Belliver Industrial Estate, Plymouth. PL6 7BP. UK.
BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. ©2010 BD
Patent Pending. Mfg by BD, USA 4/2010
40149-03

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REFERENCES:

- CLSI. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition. CLSI document H04-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- Preventing Lead Poisoning in Young Children. A Statement by the Centers for Disease Control – October 1991, U.S. Department of Health and Human Services. p.41.

Becton, Dickinson and Company, Franklin Lakes, NJ USA

EC REP Becton, Dickinson and Company, Belliver Industrial Estate, Plymouth. PL6 7BP. UK. BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. ©2010 BD Patent Pending.

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REFERENCES:

1. CLSI. Procedures and Devices for the Collection of Diagnostic
Capillary Blood Specimens; Approved Standard-Sixth Edition. CLSI
document H04-A6. Wayne, PA: Clinical and Laboratory Standards
Institute; 2008.
2. Preventing Lead Poisoning in Young Children. A Statement by the
Centers for Disease Control – October 1991, U.S. Department
of Health and Human Services. p.41.

Becton, Dickinson and Company, Franklin Lakes, NJ USA
 Becton, Dickinson and Company, Belliver Industrial Estate, Plymouth. PL6 7BP. UK.
BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. ©2010 BD
Patent Pending. Mfg by BD, USA 4/2010
40149-03

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 Becton, Dickinson and Company, Franklin Lakes, NJ USA
 Becton, Dickinson and Company, Belliver Industrial Estate, Plymouth. PL6 7BP. UK.
 BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. ©2010 BD
 Patent Pending. Mfg by BD, USA 4/2010
 40149-03

Symbol and Mark Key

EN For reference material, www.bd.com/vacutainer
BG За материали за справка, вижте <http://www.bd.com/vacutainer>
CS Referenční materiál najdete na www.bd.com/vacutainer
DA Referencemateriale findes på www.bd.com/vacutainer
DE Als Referenzmaterial, www.bd.com/vacutainer
EL Για υλικό αναφοράς, www.bd.com/vacutainer
ES Para consultar material de referencia, www.bd.com/vacutainer
ET Viitematerjal: www.bd.com/vacutainer
FI Viitemateriaalia on sivustolla www.bd.com/vacutainer
FR Pour matériel de référence, www.bd.com/vacutainer
HR Za referentni materijal pogledajte www.bd.com/vacutainer
HU Referencianyagok: www.bd.com/vacutainer
IS Sjá nánari upplýsingar, www.bd.com/vacutainer
IT Per il materiale di riferimento: www.bd.com/vacutainer
JA 参考資料: www.bd.com/vacutainer
KK Сілтемелік материалдар үшін осы вебсайтты қараңыз www.bd.com/vacutainer
KO 참고 자료 사이트: www.bd.com/vacutainer
LT Daugiau informacijos pateikta www.bd.com/vacutainer
LV Materiāli atsaucēi atrodami šeit: www.bd.com/vacutainer
NL Voor referentiemateriaal, www.bd.com/vacutainer
NO For dokumentasjon, gå til, www.bd.com/vacutainer
PL Materiały referencyjne znaleźć można pod adresem www.bd.com/vacutainer
PT Para material de referência, www.bd.com/vacutainer
RO Pentru materiale de referință, www.bd.com/vacutainer
RU Справочные материалы см. на Web-сайте www.bd.com/vacutainer
SK Viac informácií nájdete na webovej stránke www.bd.com/vacutainer
SL Referenčno gradivo najdete pod www.bd.com/vacutainer
SV För referensmaterial, se www.bd.com/vacutainer
TR Referans materyeli için, www.bd.com/vacutainer
UK За довідковими матеріалами звертайтеся на сайт: www.bd.com/vacutainer

 **EN** Do Not Use If Package Is Damaged
BG Да не се използва, ако опаковката е нарушена
CS Nepoužívejte, je-li obal poškozen
DA Må ikke anvendes, hvis pakningen er beskadiget
DE Bei beschädigter Verpackung nicht verwenden
EL Μην χρησιμοποιείτε το προϊόν αν η συσκευασία είναι κατεστραμμένη
ES No utilizar si el envase está dañado
ET Mitte kasutada, kui pakend on kahjustatud
FI Ei saa käyttää, jos pakkaus on vahingoittunut
FR Ne pas utiliser si l'emballage est endommagé
HR Ne koristite ako je ambalaža oštećena
HU Ne használja fel amennyiben a csomag sérült
IS Notið ekki ef pakkningin er skemmd
IT Non utilizzare se la confezione è danneggiata
JA パッケージが破損している場合は使用しないこと
KK Қорабы бұзылған жағдайда бұйымды қолданбаңыз
KO 포장이 손상된 경우 사용하지 마십시오
LT Negalima naudoti, jei pakuotė pažeista
LV Neizmantot, ja iepakojums ir bojāts
NL Niet gebruiken indien de verpakking beschadigd is
NO Må ikke brukes hvis pakken er skadet
PL Nie stosować, jeżeli opakowanie zostało uszkodzone
PT Não utilizar se a embalagem estiver danificada
RO Nu utilizați dacă ambalajul este deteriorat
RU Не использовать, если упаковка повреждена
SK Nepoužívajte, ak je obal poškodený
SL Ne uporabljajte, če je ovojnina poškodovana
SV Får ej användas om förpackningen är skadad
TR Ambalaj Hasarlıysa Kullanmayın
UK Не використовувати, якщо упаковка пошкоджена

EN Recyclable	KK Қайтадан өнделіне алынады
BG Може да се подлага на рециклиране	KO 재활용
CS Recyklovatelné	LT Perdirbamas
DA Genanvendelig	LV Derīgs otrreizējai pārstrādei
DE Recyclen	NL Recycleerbaar
EL Ανακυκλώσιμο	NO Gjenvinnbar
ES Reciclable	PL Podlega utylizacji
ET Ümbertöödeldav.	PT Reciclável
FI Kierrätettävä	RO Reciclabil
FR Recyclable	RU Пригодно для переработки
HR Pogodno za recikliranje	SK Recyklovateľné
HU Visszaforgatható	SL Primerno za reciklažo
IS Endurvinnanlegt	SV Återanvändbar
IT Riciclabile	TR Geri Dönüşümlü
JA リサイクル可能	UK Підлягає вторинній переробці

Symbol and Mark Key

	Do Not Reuse
	Keep Away From Sunlight
	Recyclable
	Use By
	Batch Code
	Temperature Limitation
	Non-Sterile
	Manufacturer
REF	Catalogue Number
	Authorized Representative
	In Vitro Diagnostic Medical Device
	Do Not Use if Package is Damaged