Bard Access Systems, Inc.

605 North 5600 West
Salt Lake City, Utah 84116 USA
801-595-0700

Clinical Information / Ligne directe d'informations cliniques / Hotline für klinische Informationen / Linea diretta per le informazioni cliniche / Teléfono de información clínica / Telefonische hulplijn voor klinische informatie / Linha de informações clínicas / Γραμμή βοήθειας κλινικών πληροφοριών / Hotline for klinisk information / Telefonnummer för klinisk information / Kliiniset tiedustelut numerosta / klinisk informasjon / Klinikai tájékoztató forródrót / klinické informace / Informacje kliniczne / Klinik Bilgi Acil Hattı:

800-443-3385

Ordering Information / Informations relatives aux commandes / Bestellinformationen / Per ordinare / Información para pedidos / Bestelinformatie / Informação sobre encomendas / Πληροφορίες παραγγελίας / Bestillingsinformation / Beställningsnummer / Tilaustiedot / Bestillingsinformasjon / Rendeléssel kapcsolatos információ / Informace k objednávkám / Informacje nt. zamawiania / Sipariş Bilgisi:

800-545-0890

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NEW IMPORTANT INFORMATION:

Recommended Flushing/Maintenance Procedure(s)

The catheter should be maintained in accordance with standard hospital protocols. Recommended catheter flushing/maintenance is as follows:

- Flush the catheter after every use, or at least weekly when not in use. Use a 10 ml or larger syringe.
- Flush the catheter with a minimum of 10 ml of 0.9% sodium chloride, using a “pulse” or “stop/start” technique.
 - Use of heparinized saline to lock each lumen of the catheter is optional.
- Disconnect the syringe and attach a sterile end cap to the catheter hub and tighten securely.
Caution: Always remove needles or syringes slowly while injecting the last 0.5 ml of saline.
- Prior to blood sampling when infusing TPN, follow routine maintenance procedure except use 20 ml saline and flush to clear TPN from the catheter.
- If resistance is met when flushing, no further attempts should be made. Further flushing could result in catheter rupture with possible embolization. Refer to institution protocol for clearing occluded catheters.

NOTE: When injecting or infusing medications that are incompatible, you should always flush the catheter with a minimum of 10 ml saline before and after each medication.

NOTE: When maintained in accordance with these instructions, the **PowerPICC SOLO*** catheter does not require the use of heparinized saline to lock the catheter lumens. However, use of heparinized saline will not adversely effect the catheter and may be necessary based on patient status or use of alternate flushing and locking techniques.

Caution: Use aseptic techniques whenever the catheter lumen is opened or connected to other devices.

Caution: The **PowerPICC SOLO*** catheter is designed for use with needleless injection caps or “direct-to-hub” connection technique. Always apply a sterile end cap on the catheter hub to prevent contamination when not in use. **Use of a needle longer than 1.6 cm may cause damage to the valve.**

Warning: Alcohol should not be used to lock, soak or declot polyurethane PICCs because alcohol is known to degrade polyurethane catheters over time with repeated and prolonged exposure.

Product Description

A family of peripherally inserted central catheters made from specially formulated and processed medical grade materials. Each **PowerPICC SOLO*** catheter has a kink resistant, reverse tapered design. Catheters are packaged in a tray with accessories for reliable long (greater than 30 days) or short (less than 30 days) term vascular access.

Sterilized by ethylene oxide. Do not re-sterilize.

PowerPICC SOLO* Catheter Valve Function

The **PowerPICC SOLO*** catheter valve controls the flow of fluids to provide clamp-free infusion therapy. Positive pressure into the catheter (gravity, pump, syringe) will open the valve, allowing fluid infusion. When negative pressure (aspiration) is applied, the valve opens allowing for the withdrawal of blood into a syringe.

Indications

The **PowerPICC SOLO*** catheter is indicated for short or long term peripheral access to the central venous system for intravenous therapy, power injection of contrast media, and allows for central venous pressure monitoring. For blood sampling, infusion, or therapy use a 4 French or larger catheter. The maximum recommended infusion rate is 5 ml/ sec for power injection of contrast media. For central venous pressure monitoring, it is recommended that a catheter lumen of 20 gauge or larger be used.

Contraindications

- The device is contraindicated whenever:
- The presence of device-related infection, bacteremia, or septicemia is known or suspected.
 - The patient's body size is insufficient to accommodate the size of the implanted device.
 - The patient is known or is suspected to be allergic to materials contained in the device.
 - Past irradiation of prospective insertion site.
 - Previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site.
 - Local tissue factors will prevent proper device stabilization and/or access.

Chloraprep* Solution One-Step Applicator Contraindications:

- Do not use in children less than 2 months of age because of the potential for excessive skin irritation and increased drug absorption.
- Do not use on patients with known allergies to chlorhexidine gluconate or isopropyl alcohol.
- Do not use for lumbar puncture or allow contact with meninges.
- Do not use on open skin wounds or as a general skin cleanser.

Warnings

General Warnings

- When using alcohol or alcohol containing antiseptics with polyurethane PICCs, care should be taken to avoid prolonged or excessive contact. Solutions should be allowed to completely dry before applying an occlusive dressing. Chlorhexidine gluconate and/or povidone iodine are the suggested antiseptics to use.
- Alcohol should not be used to lock, soak or declot polyurethane PICCs because alcohol is known to degrade polyurethane catheters over time with repeated and prolonged exposure.
- Do not use the catheter if there is any evidence of mechanical damage or leaking. Damage to the catheter may lead to rupture, fragmentation, possible embolism, and surgical removal.
- If signs of extravasation exist, discontinue injections. Begin appropriate medical intervention immediately.
- Do not wipe the catheter with acetone based solutions or polyethylene glycol containing ointments. These can damage the polyurethane material if used over time.
- Intended for Single Patient Use. DO NOT REUSE. The **PowerPICC SOLO*** is a single use device and should never be re-implanted. Reuse carries with it the attendant concern of cross-infection regardless of the cleaning or sterilization method. Re-sterilization of incompletely cleaned devices may not be effective. Any device that has been contaminated by blood must not be reused or re-sterilized.
- After use, this product may be a potential biohazard. Handle and discard in accordance with accepted medical practice and applicable local, state and federal laws and regulations.
- The fluid level in the catheter will drop if the catheter connector is held above the level of the patient's heart and opened to air. To help prevent a drop in the fluid level (allowing air entry) while changing injection caps, hold the connector below the level of the patient's heart before removing the injection cap.
- CVP Monitoring should always be used in conjunction with other patient assessment metrics when evaluating cardiac function.

Placement Warnings

- If the artery is entered, withdraw the needle and apply manual pressure for several minutes.
- Place a finger over the orifice of the sheath to minimize blood loss and risk of air aspiration. The risk of air embolism is reduced by performing this part of the procedure with the patient performing the Valsalva maneuver until the catheter is inserted into the sheath.
- This is not a right atrium catheter. Avoid positioning the catheter tip in the right atrium. Placement or migration of the catheter tip into the right atrium may cause cardiac arrhythmia, myocardial erosion or cardiac tamponade. The risk of these complications may be more likely in neonatal patients.

Power Injection Warnings

- Exceeding the maximum flow rate of 5 ml/sec, or the maximum pressure of power injectors of 300 psi (2068 kPa), may result in catheter failure and/or catheter tip displacement.
- Failure to ensure patency of the catheter prior to power injection studies may result in catheter failure.
- Failure to warm contrast media to body temperature prior to power injection may result in catheter failure.
- Use of lumens not marked “Power Injectable” for power injection of contrast media may cause failure of the catheter.
- Power injector machine pressure limiting feature may not prevent over-pressurization of an occluded catheter, which may cause catheter failure.
- PowerPICC SOLO*** catheter indication for power injection of contrast media implies the catheter's ability to withstand the procedure, but does not imply appropriateness of the procedure for a particular patient. A suitably trained clinician is responsible for evaluating the health status of a patient as it pertains to a power injection procedure.

Chloraprep* Solution One-Step Applicator Warnings

- Flammable, keep away from fire or flame.
- Do not use with electrocautery procedures.
- For external use only.
- When using this product keep out of eyes, ears, and mouth. May cause serious or permanent injury if permitted to enter and remain. If contact occurs, rinse with cold water right away and contact a physician.
- Stop use and ask a doctor if irritation, sensitization, or allergic reaction occurs. These may be signs of a serious condition.
- Keep out of reach of children. If swallowed, get medical help or contact a poison control center right away.

Precautions

General Precautions

- Sterilized by ethylene oxide. Do not re-sterilize.
- Carefully read and follow all instructions prior to use.
- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- Only qualified health care practitioners should insert, manipulate and remove these devices.
- Follow Universal Precautions when inserting and maintaining the catheter.
- Follow all contraindications, warnings, cautions, precautions and instructions for all infusates, including contrast media, as specified by their manufacturer.
- Precautions are intended to help avoid catheter damage and/or patient injury.
- To minimize the risk of catheter breakage and embolization, the catheter must be secured in place.
- Always remove needles or syringes slowly while injecting the last 0.5 ml of saline.
- Acetone and tincture of iodine should not be used.
- Use aseptic techniques whenever the catheter lumen is opened or connected to other devices.
- Examine the package carefully before opening to confirm its integrity and that the expiration date has not passed. The device is supplied in a sterile package and is non-pyrogenic. Do not use if package is damaged, opened or the expiration date has passed.
- Inspect kit for inclusion of all components.
- Flush the catheter with sterile normal saline prior to use. Catheter stylet must be wetted prior to stylet repositioning or withdrawal.
- Accessories and components used in conjunction with this device should incorporate luer lock connections.
- DO NOT USE A SYRINGE SMALLER THAN 10 ml. Prolonged infusion pressure greater than 25 psi (172 kPa) may damage blood vessels or viscus.

Precautions Related to Device Placement Procedure

- The **PowerPICC SOLO*** catheter features a reverse-taper catheter design. Placement of larger catheters at or below antecubital fossa may result in an increased incidence of phlebitis. Placement of the **PowerPICC SOLO*** catheter above antecubital fossa is recommended.
- Avoid placement or securement of the catheter where kinking may occur, to minimize stress on the catheter, patency problems or patient discomfort.
- Do not place suture around the catheter. Sutures may damage the catheter or compromise catheter patency.
- Do not cut stylet.
- Do not advance the guidewire past the axilla without fluoroscopic guidance or other tip locating methods.
- If the guidewire must be withdrawn while the needle is inserted, remove both the needle and wire as a unit to prevent the needle from damaging or shearing the guidewire.
- Never use force to remove the stylet. Resistance can damage the catheter. If resistance or bunching of the catheter is observed, stop stylet withdrawal and allow the catheter to return to normal shape. Withdraw both the catheter and stylet together approximately 2 cm and reattempt stylet removal. Repeat this procedure until the stylet is easily removed. Once the stylet is out, advance the catheter into the desired position.
- Avoid accidental device contact with sharp instruments and mechanical damage to the catheter material. Use only smooth-edged atraumatic clamps or forceps.
- Avoid perforating, tearing or fracturing the catheter when using a guidewire.
- Do not use the catheter if there is any evidence of mechanical damage or leaking.
- Avoid sharp or acute angles during implantation which could compromise the patency of the catheter lumen.
- The **PowerPICC SOLO*** catheter is designed for use with needleless injection caps or “direct-to-hub” connection technique. Always apply a sterile end cap on the catheter hub to prevent contamination when not in use. **Use of a needle longer than 1.6 cm may cause damage to the valve.**

Possible Complications

The potential exists for serious complications including the following:

- | | | |
|---|--|--|
| <ul style="list-style-type: none">Air EmbolismBleedingBrachial Plexus InjuryCardiac ArrhythmiaCardiac TamponadeCatheter Erosion Through the SkinCatheter EmbolismCatheter OcclusionCatheter Related SepsisEndocarditis | <ul style="list-style-type: none">Exit Site InfectionExit Site NecrosisExtravasationFibrin Sheath FormationHematomaIntolerance Reaction to Implanted DeviceLaceration of Vessels or ViscusMyocardial ErosionPerforation of Vessels or Viscus | <ul style="list-style-type: none">PhlebitisSpontaneous Catheter Tip Malposition or RetractionThromboembolismVenous ThrombosisVessel ErosionRisks Normally Associated with Local or General Anesthesia, Surgery, and Post-Operative Recovery |
|---|--|--|

Insertion Instructions


1. Identify the Vein and Insertion Site

- Apply a tourniquet above the anticipated insertion site.
- Select and mark the vein based on patient assessment. Recommended veins are basilic, cephalic and median cubital veins.
- Caution:** The **PowerPICC SOLO*** catheter features a reverse-taper catheter design. Placement of larger catheters at or below antecubital fossa may result in an increased incidence of phlebitis. Placement of the **PowerPICC SOLO*** catheter above antecubital fossa is recommended.
- Caution:** Avoid placement or securement of the catheter where kinking may occur, to minimize stress on the catheter, patency problems or patient discomfort.
- Release tourniquet.

2. Patient Position / Catheter Measurement

- Position the arm at a 90° angle.
- For central placement, the recommended target tip location is in the lower 1/3 of the Superior Vena Cava (SVC). Measure from the planned insertion site to the right clavicular head, then down to the third intercostal space.
- Note:** The external measurement can never exactly duplicate the internal venous anatomy.

3. Skin Preparation


- Don prep gloves, head cover and mask.
- Apply underdrape.
- Prepare the site with the Chloraprep* Solution One-Step Applicator or according to institutional policy using sterile technique.
- Pinch the wings of the Chloraprep* Solution One-Step Applicator to break the ampule and release the antiseptic. **Do not touch the sponge.**
- Wet the sponge by repeatedly pressing and releasing the sponge against the treatment area until fluid is visible on the skin.
- Use repeated back-and-forth strokes of the sponge for approximately 30 seconds. Completely wet the treatment area with antiseptic. Allow the area to dry for approximately 30 seconds. **Do not blot or wipe away.**
- Maximum treatment area for one applicator is approximately 130 cm² (10 cm x 13 cm). Discard the applicator after a single use.
- When alcohol is used as a skin prep, it must be allowed to completely air dry before proceeding with insertion.
- Remove and discard gloves.

4. Sterile Field Preparation

- Apply the tourniquet above the intended insertion site to distend the vessel.
- Don sterile gloves and gown.
- Apply fenestrated and body drapes & complete sterile field preparation.

5. Preflush the Catheter


- Attach prefilled syringe to the luer attachment on the T-Lock extension set.
- Flush catheter lumen(s) until saline drips out of distal end of the catheter.
- The syringe may be left attached during procedure.

6. Modification of Catheter Length

Note: Catheters can be cut to length if a different length is desired due to patient size and desired point of insertion according to hospital protocol. Catheter depth markings are in centimeters.

- Measure the distance from the zero mark to the desired tip location.
- Disconnect the T-Lock from the stylet funnel.
- Withdraw the entire T-Lock connector/stylet assembly, as one unit.
- Retract the stylet to well behind the point the catheter is to be cut.
- Using a sterile scalpel or scissors, carefully cut the catheter according to institutional policy if necessary.
Caution: Do not cut stylet.
- Inspect cut surface to assure there is no loose material.
- Re-advance the T-lock connector/stylet assembly locking the connector to the stylet funnel. Assure stylet tip is intact.
- Gently retract the stylet through the locked T-lock connector until the stylet tip is contained inside the catheter.
- Assure proper alignment of the stylet to the distal end of the trimmed catheter.

7. Perform Venipuncture

- Anesthetize with local anesthesia as required.
- Insert the safety introducer needle into the desired vein.
Alternate Technique: The safety IV catheter may be used as an alternate to the safety introducer needle. Remove the needle from the catheter after the vein is accessed.
Warning: If the artery is entered, withdraw the needle or safety IV catheter and apply manual pressure for several minutes.
- Release tourniquet.
- Remove the guidewire tip protector from the guidewire hoop and insert the flexible end of the guidewire into the introducer needle or catheter and into the vein. Advance the guidewire to the desired depth.
Caution: Do not advance the guidewire past the axilla without fluoroscopic guidance or other tip locating methods.
- Gently withdraw and remove the safety introducer needle or catheter, while holding the guidewire in position.
Caution: If the guidewire must be withdrawn while the needle is inserted, remove both the needle and wire as a unit to prevent the needle from damaging or shearing the guidewire.
- Advance the small sheath and dilator together as a unit over the guidewire, using a slight rotational motion. If necessary, a small incision may be made adjacent to the guidewire to facilitate insertion of the sheath and dilator. Verify institutional guide lines concerning the use of a safety scalpel prior to making incision.
- Withdraw the dilator and guidewire, leaving the small sheath in place.
Warning: Place a finger over the orifice of the sheath to minimize blood loss and risk of air aspiration. The risk of air embolism is reduced by performing this part of the procedure with the patient performing the Valsalva maneuver until the catheter is inserted into the sheath.

8. Insert and Advance the Catheter

- Insert the catheter into the introducer sheath.
- Advance the catheter slowly.
- Complete catheter advancement into the desired position.

9. Complete Catheter Insertion

- Continue to advance the catheter. For central placement, when the tip has advanced to the shoulder, have the patient turn head (chin on shoulder) toward the insertion side to prevent possible insertion into the jugular vein. The **PowerPICC SOLO*** catheter features a reverse-taper catheter design. Placement of larger catheters at or below antecubital fossa may result in an increased incidence of phlebitis. Placement of the **PowerPICC SOLO*** catheter above antecubital fossa is recommended.
Note: Resistance may be felt approximately 7cm distal of catheter hub when introducing the catheter into the sheath due to an increase in outer diameter (O.D.) The introducer may be partially split, but not removed to facilitate insertion of the catheter past this point if necessary.
Note: PICCs should be positioned with the catheter tip in the lower 1/3 of the SVC. Verify correct catheter tip position using radiography or appropriate technology.
Warning: This is not a right atrium catheter. Avoid positioning the catheter tip in the right atrium. Placement or migration of the catheter tip into the right atrium may cause cardiac arrhythmia, myocardial erosion or cardiac tamponade. The risk of these complications may be more likely in neonatal patients.

10. Retract and Remove the Introducer Sheath

- Stabilize the catheter position by applying pressure to the vein distal to the introducer sheath.
- Withdraw the introducer sheath from the vein and away from the site.
- Split the introducer sheath and peel it away from the catheter.

11. Remove the Stylet, Stylet Funnel, and T-Lock Assembly

- Disconnect the T-Lock and stylet funnel from the catheter luer connector.
- Stabilize the catheter position by applying light pressure to the vein distal to the insertion site.
- Slowly remove the T-Lock, stylet funnel, and stylet, as a unit. Do not remove stylet through T-Lock.
Caution: Never use force to remove the stylet. Resistance can damage the catheter.
If resistance or bunching of the catheter is observed, stop stylet withdrawal and allow the catheter to return to normal shape. Withdraw both the catheter and stylet together approximately 2 cm and reattempt stylet removal. Repeat this procedure until the stylet is easily removed. Once the stylet is out, advance the catheter into the desired position.

12. Aspirate and Flush

- Attach primed extension set and/or saline-filled syringe.
- Aspirate for adequate blood return and flush catheter with 10ml normal saline to ensure patency.
Caution: The **PowerPICC SOLO*** catheter is designed for use with needleless injection caps or “direct-to-hub” connection technique. Always apply a sterile end cap on the catheter hub to prevent contamination when not in use. Use of a needle longer than 1.6 cm may cause damage to the valve.
Caution: Always remove needles or syringes slowly while injecting the last 0.5 ml of saline.
- Cap catheter.
Warning: The fluid level in the catheter will drop if the catheter connector is held above the level of the patient’s heart and opened to air. To help prevent a drop in the fluid level (allowing air entry) while changing injection caps, hold the connector below the level of the patient’s heart before removing the injection cap.

13. Securing the PowerPICC SOLO* Catheter

The StatLock® catheter stabilization device is included in **PowerPICC SOLO*** catheter kits. Please refer to Instructions For Use on the proper use and removal. The StatLock® catheter stabilization device should be monitored daily and replaced at least every seven days.

Caution: To minimize the risk of catheter breakage and embolization, the catheter must be secured in place.




Warning: When using alcohol or alcohol containing antiseptics with polyurethane PICCs, care should be taken to avoid prolonged or excessive contact. Solutions should be allowed to completely dry before applying an occlusive dressing. Chlorhexidine gluconate and/or povidone iodine are the suggested antiseptics to use.

Warning: Alcohol should not be used to lock, soak or de clot polyurethane PICCs because alcohol is known to degrade polyurethane catheters over time with repeated and prolonged exposure.

Warning: Do not wipe the catheter with acetone based solutions or polyethylene glycol containing ointments. These can damage the polyurethane material if used over time.

The StatLock® Stabilization Device Procedure

Single Lumen





1. Secure catheter with the StatLock® stabilization device.

2. Cover site and the StatLock® stabilization device with transparent dressing.

3. Place anchor tape sticky side up, under hub. Wedge tape between hub and wings.

4. Chevron anchor tape on top of transparent dressing.

Dual Lumen





1. Secure catheter with the StatLock® stabilization device.

2. Cover site and the StatLock® stabilization device with transparent dressing.

3. Place 1st anchor tape sticky side up, under one extension leg. Wedge tape between hub and wings. Chevron anchor tape on top of transparent dressing.

4. Place 2nd anchor tape sticky side up under hub. Wedge tape between hub and wings. Chevron anchor tape on top of transparent dressing.

Triple Lumen

1. Secure catheter with the StatLock® stabilization device.





2. Cover site and the StatLock® stabilization device with transparent dressing.

3. Place 1st anchor tape sticky side up, under one extension leg. Wedge tape between hub and wings. Chevron anchor tape on top of transparent dressing.

4. Place 2nd and 3rd anchor tapes sticky side up under remaining hubs. Wedge tape between hubs and wings. Chevron anchor tape on top of transparent dressing.

Tape Strip Stabilization Procedure

Single Lumen





1. Place 1st anchor tape over wings or junction.

2. Cover site and 1st anchor tape with transparent dressing up to hub, but not over hub.

3. Place 2nd anchor tape sticky side up under hub and close to transparent dressing. Wedge tape between hub and wings.

4. Chevron 2nd anchor tape on top of transparent dressing and place 3rd anchor tape over hub

Dual Lumen





1. Place 1st anchor tape over wings or bifurcation.

2. Cover site and 1st anchor tape with transparent dressing up to hub, but not over hub.

3. Place 2nd anchor tape sticky side up under hub and close to transparent dressing. Wedge tape between hub and wings. Anchor only one hub of dual lumen catheter.

4. Chevron 2nd anchor tape on top of transparent dressing and place 3rd anchor tape over hub

Triple Lumen

1. Place 1st anchor tape over wings or trifurcation.

2. Cover site and 1st anchor tape with transparent dressing up to hub, but not over hub.

3. Place 2nd anchor tape sticky side up under hub and close to transparent dressing. Wedge tape between hub and wings. Anchor only one hub of the triple lumen catheter. Chevron anchor tape on top of transparent dressing.

4. Place 2nd and 3rd anchor tapes sticky side up under remaining hubs. Wedge tape between hubs and wings. Chevron anchor tape on top of transparent dressing and place third anchor tape over hub.

14. Verify Placement

- PICCs should be positioned with the catheter tip in the lower 1/3 of the SVC. Verify correct catheter tip position using radiography or appropriate technology.

15. Power Injection Procedure

Warning: **PowerPICC SOLO*** catheter indication for power injection of contrast media implies the catheter’s ability to with stand the procedure, but does not imply appropriateness of the procedure for a particular patient. A suitably trained clinician is responsible for evaluating the health status of a patient as it pertains to a power injection procedure.

- Remove the injection/needleless cap from the **PowerPICC SOLO*** catheter.
- Attach a 10 ml or larger syringe filled with sterile normal saline.
- Aspirate for adequate blood return and vigorously flush the catheter with the full 10 ml of sterile normal saline.
Warning: Failure to ensure patency of the catheter prior to power injection studies may result in catheter failure.
- Detach syringe.
- Attach the power injection device to the **PowerPICC SOLO*** catheter per manufacturer’s recommendations.
- Contrast media should be warmed to body temperature prior to power injection.
Warning: Failure to warm contrast media to body temperature prior to power injection may result in catheter failure.
- Use only lumens marked “Power Injectable” for power injection of contrast media.
Warning: Use of lumens not marked “Power Injectable” for power injection of contrast media may cause failure of the catheter.
- Complete power injection study taking care not to exceed the flow rate limits. Do not exceed the maximum flow rate of 5ml/sec.
Warning: Exceeding the maximum flow rate of 5 ml/sec, or the maximum pressure of power injectors of 300 psi (2068 kPa), may result in catheter failure and/or catheter tip displacement.
Warning: Power injector machine pressure limiting feature may not prevent over-pressurization of an occluded catheter, which may cause catheter failure.
- Disconnect the power injection device.
- Replace the injection/needleless cap on the **PowerPICC SOLO*** catheter.
- Flush the **PowerPICC SOLO*** catheter with 10 ml of sterile normal saline, using a 10 ml or larger syringe. Use of heparinized saline to lock each lumen of the catheter is optional.

- The **PowerPICC SOLO*** catheter testing included 10 power injection cycles.

16. Suggested Catheter Maintenance

- Dressing Changes**
 - Assess the dressing in the first 24 hours for accumulation of blood, fluid or moisture beneath the dressing. During all dressing changes, assess the external length of the catheter to determine if migration of the catheter has occurred. Periodically confirm catheter placement, tip location, patency and security of dressing.
 - Maintain according to hospital protocol. Avoid using acetone based solutions, or ointment. These substances are known to degrade polyurethane.
 - Chlorhexidine gluconate is the suggested antiseptic to use. 2% Chlorhexidine gluconate /70% isopropyl alcohol swab sticks may be used for dressing changes. Povidone-iodine may also be used as an antiseptic.
 - Allow all cleaning agents / antiseptics to dry completely before applying dressing.
Caution: Acetone and tincture of iodine should not be used.
Warning: When using alcohol or alcohol containing antiseptics with polyurethane PICCs, care should be taken to avoid prolonged or excessive contact. Solutions should be allowed to completely dry before applying an occlusive dressing. Chlorhexidine gluconate and/or povidone iodine are the suggested antiseptics to use.
- Flushing**
 - Flush the catheter after every use, or at least weekly when not in use. Use a 10 ml or larger syringe.
 - Flush the catheter with a minimum of 10 ml of 0.9% sodium chloride, using a “pulse” or “stop/start” technique. Use of heparinized saline to lock each lumen of the catheter is optional.
 - Disconnect the syringe and attach a sterile end cap to the catheter hub and tighten securely.
Caution: Always remove needles or syringes slowly while injecting the last 0.5 ml of saline.
 - Prior to blood sampling when infusing TPN, follow routine maintenance procedure except use 20 ml saline and flush to clear TPN from the catheter.
 - If resistance is met when flushing, no further attempts should be made. Further flushing could result in catheter rupture with possible embolization. Refer to institution protocol for clearing occluded catheters.
NOTE: When injecting or infusing medications that are incompatible, you should always flush the catheter with a minimum of 10 ml saline before and after the medication.
NOTE: When maintained in accordance with these instructions, the **PowerPICC SOLO*** catheter does not require the use of heparinized saline to lock the catheter lumens. However, use of heparinized saline will not adversely effect the catheter and may be necessary based on patient status or use of alternate flushing and locking techniques.
Caution: Use aseptic techniques whenever the catheter lumen is opened or connected to other devices.
Caution: The **PowerPICC SOLO*** catheter is designed for use with needleless injection caps or “direct-to-hub” connection technique. Always apply a sterile end cap on the catheter hub to prevent contamination when not in use.

Use of a needle longer than 1.6 cm may cause damage to the valve.

Warning: Alcohol should not be used to lock, soak or de clot polyurethane PICCs because alcohol is known to degrade polyurethane catheters over time with repeated and prolonged exposure.

C. **Occluded or Partially Occluded Catheter**

Catheters that present resistance to flushing and aspiration may be partially or completely occluded. Do not flush against resistance. If the lumen will neither flush nor aspirate and it has been determined that the catheter is occluded with blood, a de clotting procedure per institution protocol may be appropriate.

D. **When Cleaning the Exit Site**

- Warning:** Do not wipe the catheter with acetone based solutions or polyethylene glycol containing ointments. These can damage the polyurethane material if used over time.
- Maintain according to hospital protocol. Avoid using acetone based solutions, or ointment. These substances are known to degrade polyurethane.
 - Use chlorhexidine gluconate or povidone iodine to clean the exit site around the catheter.
 - Allow all cleaning agents/antiseptics to dry completely before applying dressing.

17. Central Venous Pressure Monitoring

- Prior to conducting central venous pressure monitoring:
- Ensure proper positioning of the catheter tip.
- Flush catheter vigorously with sterile normal saline.
- Ensure the pressure transducer is at the level of the right atrium.
- It is recommended that a continuous infusion of saline (3 ml/hr) is maintained through the catheter while measuring CVP to improve accuracy of CVP results.
- Use your institution's protocols for central venous pressure monitoring procedures.

Warning: CVP Monitoring should always be used in conjunction with other patient assessment metrics when evaluating cardiac function.

18. Catheter Removal

- Remove dressing and StatLock® catheter stabilization device or tape securement strips.
- Grasp catheter near insertion site.
- Remove slowly. Do not use excessive force.
- If resistance is felt, stop removal. Apply warm compress and wait 20-30 minutes.
- Resume removal procedure.
- Examine catheter tip to determine that the entire catheter has been removed.

Bard Access Systems, Inc. warrants to the original purchaser that this product will be free from defects in material and workmanship for a period of one (1) year from the date of purchase. If this product proves to be so defective, purchaser may return same to Bard Access Systems, Inc. for repair or replacement, at Bard Access Systems, Inc. option. All returns must be authorized in advance in accordance with Bard Access Systems, Inc. Returned Goods Policy found in its then current Price List. The liability of Bard Access Systems, Inc. under this limited product warranty does not extend to any abuse or misuse of this product or its repair by anyone other than an authorized Bard Access Systems, Inc. representative.

THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED (INCLUDING, WITHOUT LIMITATION ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.THE LIABILITY AND REMEDY STATED IN THIS LIMITED PRODUCT WARRANTY WILL BE THE SOLE LIABILITY OF BARD ACCESS SYSTEMS AND REMEDY AVAILABLE TO PURCHASER FOR THIS PRODUCT, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE) OR OTHERWISE, AND BARD ACCESS SYSTEMS WILL NOT BE LIABLE TO PURCHASERS FOR ANY INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF ITS HANDLING OR USE.

An issued or revision date for these instructions is included for the user's information. In the event two years have elapsed between this date and product use, the user should contact Bard Access Systems, Inc. to see if additional product information is available.

Revised date: April 2009

*Bard, PowerPICC, PowerPICC SOLO, "The Power of Purple", Power. Saline. One., "The Simplicity of Saline", "One Valve. One PICC.", StatLock, Sherlock, and the color purple are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

ChloraPrep is a trademark and/or registered trademark of Cardinal Health, Inc or an affiliate.

© Copyright 2009 C. R. Bard, Inc. All rights reserved.

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