


Clearsight Jr Cuff Technology Overview

How does it work?


Volume clamp method

- The essence is to dynamically provide equal pressures on either side of the wall of the artery by clamping the artery to a certain constant volume
- 1000 times each second the cuff pressure is adjusted to keep the diameter of the finger arteries constant
- Continuous recording of the cuff pressure results in a real-time finger pressure waveform¹




The Physiological method – Physiological Calibration

- The Physiological method is the real-time expert system that determines the proper arterial 'unloaded' volume, i.e. no pressure gradient across the arterial wall
- Automatic, periodic adjustments are essential to track the unloaded volume clamp setpoint when smooth muscle tone changes (e.g. during vasoconstriction)
- Calibration interval starts at 10 beats, but it increases to every 70 beats as stability increases
- Physiological interval >30 beats is considered reliable²




Radial pressure reconstruction

- Clinical standard for noninvasive BP is radial level
- The ClearSight Jr cuff reconstructs the radial arterial pressure waveform from the finger arterial pressure waveform
- The reconstruction algorithm is based on a clinical database³



Cardiac output calculation

- Stroke volume is calculated by an algorithm based on a pulse contour method using:
 - Mean arterial pressure (MAP) and heart rate (HR)
 - A physiological model to calculate afterload individualized by age, gender, height and weight
- Cardiac output results from stroke volume times heart rate and is updated every beat⁴



References

- 1 Peñáz J. Photoelectric measurement of blood pressure, volume and flow in the finger. 1973; Dresden 1973, p. 104
- 2 Wesseling KH, Wit B, Hoeven GMA, Goudoever J, Settels JJ. Physiological, calibrating finger vascular physiology for Finapres. Homeostasis. 1995;36:67–82.
- 3 Gizdulich P, Prentza A, Wesseling KH. Models of brachial to finger pulse wave distortion and pressure decrement. Cardiovasc Res. 1997;33:698–705. doi: 10.1016/S0008-6363(97)00003-5
- 4 Truijen J, van Lieshout JJ, Wesseling WA, Westerhof BE. Noninvasive continuous hemodynamic monitoring. J Clin Monit Comput. 2012 Jun 14.

Validation of monitoring blood pressure and cardiac output with noninvasive finger cuff technology versus traditional methods

Blood pressure validation studies

BP can be measured reliably according to AAMI standard, bias <5 and SD < 8 mmHg

Vs. Noninvasive upper arm cuff	
Kouz, K. et al. Anesthesiology (2023) ⁵	242 non-cardiac surgery patients
Maheshwari K. et al Anesthesiology (2018) ⁶	320 non-cardiac surgery patients
Meidert, A.S. et al. European Journal of Anaesthesiology (2017) ⁷	160 orthopedic patients with a history of chronic hypertension
Chen, G. et al. Journal of Clinical Monitoring and Computing (2012) ⁸	25 ASA I-III patients undergoing abdominal or orthopedic surgery
Vs. Invasive radial line	
Kho, et al. Journal of Clinical Anesthesia (2022) ⁹	33 patients non-cardiac surgery patients
Shu Y. Lu et al. Journal of Cardiothoracic and Vascular Anesthesia (2021) ¹⁰	20 transcatheter aortic valve replacement patients

Cardiac output validation studies

Measurement of relative changes in CO performs comparably with invasive measurements

Vs. Noninvasive upper arm cuff	
Bubenek-Turconi et al – Anesthesia-Analgesia 2013 ¹¹	28 cardiac surgery patients
Bogert et al – Anaesthesia 2010 ¹²	25 post CABG patients
Stover et al – BMC Anesthesiology 2009 ¹³	10 severely ill ICU patients
Vs. Transpulmonary thermodilution	
Broch et al – Anaesthesia 2012 ¹⁴	40 cardiac surgery patients
Hofhuizen et al. – J of Critical Care 2014 ¹⁵	20 post-cardiac patients
Vs. Trans-thoracic echo-Doppler	
van der Spoel et al. – J Clin Anesth 2012 ¹⁶	40 ASA 1-2 patients
Vs. Esophageal echo-Doppler	
Chen et al. J Clin Anesth 2012 ¹⁷	25 ASA 1-3 patients

References (continued)

- Kouz, K. et al. (2023) 'Continuous finger-cuff versus intermittent oscillometric arterial pressure monitoring and hypotension during induction of anesthesia and noncardiac surgery: The detect Randomized Trial', *Anesthesiology*, 139(3), pp. 298–308.
- Maheshwari K, Khanna S, Bajracharya GR, Makarova N, Riter Q, Raza S, et al. A randomized trial of continuous noninvasive blood pressure monitoring during noncardiac surgery. *Anesth Analg*. 2018;127(2):424-431.
- Meidert, A.S. et al. (2017) 'The impact of continuous non-invasive arterial blood pressure monitoring on blood pressure stability during general anaesthesia in orthopaedic patients', *European Journal of Anaesthesiology*, 34(11), pp. 716–722.
- Chen, G. et al. (2012) 'Impact of noninvasive and beat-to-beat arterial pressure monitoring on intraoperative hemodynamic management', *Journal of Clinical Monitoring and Computing*, 26(2), pp. 133–140.
- Kho E, van der Ster BJP, van der Ven WH, Vlaar APJ, Immink RV, Veelo DP. Clinical agreement of a novel algorithm to estimate radial artery blood pressure from the non-invasive finger blood pressure. *J Clin Anesth*. 2022 Dec;83:110976. doi: 10.1016/j.jclinane.2022.110976. Epub 2022 Sep 26. PMID: 36174389.
- Shu Y. Lu, Adam A. Dalia, Continuous Noninvasive Arterial Pressure Monitoring for Transcatheter Aortic Valve Replacement, *Journal of Cardiothoracic and Vascular Anesthesia*, Volume 35, Issue 7, 2021, Pages 2026-2033, <https://doi.org/10.1053/j.jvca.2021.01.012>.
- Bubenek-Turconi SI, Craciun M, Miclea I, Perel A. Noninvasive Continuous Cardiac Output by the Nexfin Before and After Preload-Modifying Maneuvers: A Comparison with Intermittent Thermodilution Cardiac Output. *Anesth Analg*. 2013 Aug;117(2):366-72.
- Bogert LW, Wesseling KH, Schraa O, Van Lieshout EJ, De Mol BA, Van GJ, Westerhof BE, Van Lieshout JJ. Pulse contour cardiac output derived from non-invasive arterial pressure in cardiovascular disease. *Anaesthesia*. 2010;65:1119–25.
- Stover JF, Stocker R, Lenherr R, Neff TA, Cottini SR, Zoller B, Bécher M. Noninvasive cardiac output and blood pressure monitoring cannot replace an invasive monitoring system in critically ill patients. *BMC Anesthesiol*. 2009 Oct 12;9:6. doi: 10.1186/1471-2253-9-6.
- Broch O, Renner J, Gruenewald M, Meybohm P, Schoettler J, Caliebe A, Steinfath M, Malbrain M, Bein B. A comparison of the Nexfin and transcardiopulmonary thermodilution to estimate cardiac output during coronary artery surgery. *Anaesthesia* 2012;67:377–83.
- Hofhuizen C, Lansdorp B, van der Hoeven JG, Scheffer GJ, Lemson J. Validation of noninvasive pulse contour cardiac output using finger arterial pressure in cardiac surgery patients requiring fluid therapy. *Journal of Critical Care* 2014 Feb;29(1):161-5.
- van der Spoel AG, Voogel AJ, Folkers A, Boer C, Bouwman RA. Comparison of noninvasive continuous arterial waveform analysis (Nexfin) with transthoracic Doppler echocardiography for monitoring of cardiac output. *J Clin Anesth*. 2012 Jun;24(4):304-9. doi: 10.1016/j.jclinane.2011.09.008.
- Chen G, Meng L, Alexander B, Tran NP, Kain ZN, Cannesson M. Comparison of noninvasive cardiac output measurements using the Nexfin monitoring device and the esophageal Doppler. *J Clin Anesth*. 2012 Jun;24(4):275-83. doi: 10.1016/j.jclinane.2011.08.014.

Medical device for professional use. See instructions for use. CE marked medical device 

BD, Eysins, Vaud, Switzerland

bd.com



BD, the BD Logo, ClearSight Jr, and PhysioCal are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2026 BD. All rights reserved. BD-162554 v2.0 EU