

BD Product Security - Patch Whitelist / Blacklist

BD Biosciences Product Name: BD FACSLyric™ Cell Analyzers (IVD and RUO)

April 2020

Critical and/or Security Patches: IE Scripting Engine patch was released with February 2020 security updates, SMBv3 Remote Code Execution patch was released with March 2020 security updates and the Adobe Type Manager Remote Code Execution patch was released with the April 2020 security updates. The April 2020 security rollup/cumulative update was released in April 2020.

Microsoft® and Third-Party Patches

BD has identified patches from Microsoft that have been indicated as critical or security related for the IE Scripting Engine, Server Message Block and Adobe Type Manager vulnerabilities. Microsoft has also released a security update for the month of April. The BD products listed above have been tested with security patches from the March and April rollup/cumulative releases installed. The patches in these updates were not found to adversely affect operation of the listed BD products.

Note: There were no updates released in March or April for Windows 7. Windows 7 support officially ended as of January 14, 2020.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s). For products that are supported remotely by BD Service, patches will be applied according to customers' service agreements.



Patch Name	Patch Description	Patch ID	Notes
Feb 11, 2020 – KB4532691 (CVE-2020-0674)	<p>This security update includes quality improvements. Key changes include:</p> <p>Security updates to Microsoft Edge, Internet Explorer, Microsoft Graphics Component, Windows Input and Composition, Windows Media, Windows Shell, the Microsoft Scripting Engine, Windows Fundamentals, Windows Management, Windows Cryptography, Windows Virtualization, Windows Hyper-V, Windows Core Networking, Windows Peripherals, Windows Network Security and Containers, Windows Storage and Filesystems, and Windows Server.</p>	KB4532691	1,2
Mar 12, 2020 – KB4551762 (CVE-2020-0796)	This security update resolves a Microsoft Server Message Block 3.1.1 protocol issue that provides shared access to files and printers.	KB4551762	1,3
Apr 14, 2020 – KB4550930 (CVE-2020-0938, CVE-2020-1020)	<p>This security update includes quality improvements. Key changes include:</p> <p>Security updates to the Microsoft Scripting Engine, Windows Kernel, Windows App Platform and Frameworks, Microsoft Graphics Component, Windows Media, Windows Fundamentals, Windows Core Networking, Windows Update Stack, and the Microsoft JET Database Engine.</p>	KB4550930	1,4
April 2020 Monthly Cumulative Updates	Security and quality updates for Microsoft Windows 10 Pro (SAC) v1909 and Enterprise (LTSC) v1809	KB4549951 (Professional) KB4549949 (Enterprise)	1,5

Notes

1. Patches may require the system to reboot.
2. For further information from Microsoft regarding the IE Scripting Engine vulnerabilities, please refer to the following link –
<https://portal.msrc.microsoft.com/en-US/security-guidance/advisory/CVE-2020-0674>
3. For further information from Microsoft regarding the SMBv3 Remote Code vulnerabilities, please refer to the following link –
<https://portal.msrc.microsoft.com/en-US/security-guidance/advisory/CVE-2020-0796>
4. For further information from Microsoft regarding the ATM Remote Code vulnerabilities, please refer to the following links –
<https://portal.msrc.microsoft.com/en-US/security-guidance/advisory/CVE-2020-0938>
<https://portal.msrc.microsoft.com/en-US/security-guidance/advisory/CVE-2020-1020>
5. For further information from Microsoft regarding CVEs covered by the April Security updates, please refer to the following link –
<https://portal.msrc.microsoft.com/en-us/security-guidance/releasenotedetail/2020-Apr>

BD FACSLyric™ Flow Cytometers are Class 1 Laser Products.

The BD FACSLyric™ Flow Cytometer is for In Vitro Diagnostic Use with BD FACSuite™ Clinical Application for up to six colors. The BD FACSLyric™ Flow Cytometer is for Research Use Only with BD FACSuite™ Application for up to 12 colors. Not for use in diagnostic or therapeutic procedures.

23-23052-00

BD, the BD Logo and FACSLyric are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2020 BD. All rights reserved.

