




BD Nexiva™

Closed IV Catheter System – Dual Port


Points to practice



- 1 Prepare the connector**
- Connect a flush syringe to the BD Q-Syte™ needle-free connector using a straight-on approach and prime the device.



- 2 Prepare the catheter**
- Secure the vent plug and BD Q-Syte™ needle-free connector.
 - The clamp should not be engaged.
 - Twist to remove the needle cover.



- Holding as shown, pull back approximately 1/8" on the finger grips (Fig. 2a).
- Push the finger grips back to their original position so the stabilization platform and finger grips are snugly together (Fig. 2b).




- 3 Insertion and flashback**
- Hold the system as shown and access the vessel (Fig. 3).
 - Initial blood return is along the catheter, then up the extension tube. Look at the catheter for initial blood return.
 - Lower and advance the catheter and needle unit slightly to ensure the catheter tip is within the vessel.




- 4 Advancement**
- Place the index finger pad behind the push-tab and push the catheter off the needle into the vessel (Fig. 4).

Tip: Do not pull back on the needle during advancement.



- 5 Needle removal**
- Stabilize the catheter system, by holding down the extension set, and pull back on the finger grips until the push-tab component releases from the stabilization platform (Fig. 5).
 - Discard the shielded needle into a puncture-resistant, leak-proof sharps container.

Tip: Do not hold onto the push-tab component of the device as this will prevent the release of the needle shield.



- 6 Stabilization and preparation for use**
- Apply a transparent dressing to cover the septum, allowing maximum use of the extension tubing (Fig. 6).
 - Engage the clamp.
 - Remove the vent plug and attach the primed BD Q-Syte™ device (Fig. 6).



- 7 Flushing**
- Release the clamp and flush (Fig. 7).



- 8 Disconnection**
- Clamp before disconnecting to minimize the reflux of blood.
 - Hold the BD Q-Syte™ device during disconnection and use a straight-off approach.

NOTE: Single lumen. Compatible meds only.

Consult the product insert for complete instructions, warnings and cautions.

BD–Canada, Mississauga, ON, L5N 0B3

bd.com

BD, the BD Logo, Nexiva, and Q-Syte are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved, BD-61526(04/22)

