

You must first download this PIVC Excellence PDF for interactive features to work



Peripheral IV Catheter (PIVC) Excellence




User guide/navigational tips

The downloaded version of this **PIVC Excellence PDF** is interactive.

That means that searching for the information you need has never been easier!

You can use key word searches to target topics for quick reference (ex. Phlebitis – Infusate – Hemodilution – Securement – Infiltration – Extravasation – Irritant – Vesicant – Venipuncture – Flush – Occlusion – Dislodgement – CABSIs – Infection).



You can also jump to scroll through the table of contents and jump into each section for additional information about the topics there.

2	Points to Practice	9
	Expected infusion therapy outcomes	10
	Therapy considerations	11
	Site preparation	14
	Device selection	19

The index of each section will allow you to hone your search by digging into individual topics that are covered in the section.

Expected infusion therapy outcomes	10
Therapy considerations	11
Site preparation	14
Device selection	19

Expected infusion therapy outcomes

The ability to achieve the expected outcomes of infusion therapy are influenced by a variety of factors, many of which are shown below.

Infusion Therapy Process Map

- 1. Successful completion of PIVC-BSI therapy
- 2. Minimize therapy-related complications
- 3. Minimize the number of venipunctures
- 4. Minimize supply and labor costs
- 5. Patient satisfaction
- 6. Reduce workload of nurses and health care workers

BD

User guide/navigational tips (cont'd)

Navigate to other sections of the **PIVC Excellence PDF** at any time through the table of contents by opening this menu in the lower left-hand corner of the page you're on.

Simply close the menu when you're done by selecting the arrow in the upper right-hand corner of the pop-up.

Expected infusion therapy outcomes

The ability to achieve the expected outcomes of infusion therapy are influenced by a variety of factors, many of which are shown below.

Infusion Therapy Process Map

Expected Outcomes

1. Successful completion of prescribed therapy
2. Minimize therapy-related complications
3. Minimize the number of venipunctures
4. Minimize supply and labor costs
5. Patient satisfaction
6. Reduced needlestick injuries and blood exposure to healthcare workers

PIVC complications

Introduction to complications	29
Phlebitis	30
CABSI	32
Occlusion	34
Infiltration/extravasation	35
Dislodgement	37

Expected infusion therapy outcomes

The ability to achieve the expected outcomes of infusion therapy are influenced by a variety of factors, many of which are shown below.

Infusion Therapy Process Map

Expected Outcomes

1. Successful completion of prescribed therapy
2. Minimize therapy-related complications
3. Minimize the number of venipunctures
4. Minimize supply and labor costs
5. Patient satisfaction
6. Reduced needlestick injuries and blood exposure to healthcare workers

Points to practice



TABLE OF CONTENTS

1

Introduction to PIVC excellence

Challenges of vascular access

The burden of PIVC insertion failure

6

7

8

2

Points to practice

Expected infusion therapy outcomes

Therapy considerations

Site selection

Device selection

9

10

11

14

19

3

Insertion considerations

Vein distention and stabilization

Venipuncture

Securement

20

21

22

23

4	PIVC management	24
	Flushing	25
	Dressing management	27
5	PIVC complications	28
	Introduction to complications	29
	Phlebitis	30
	CABSI	32
	Occlusion	34
	Infiltration/extravasation	35
	Dislodgement	37
6	References	39
	References for tabs 1–5	40



Introduction to PIVC excellence

Challenges of vascular access

7

The burden of PIVC insertion failure

8

Challenges of vascular access

Peripheral intravenous catheters (PIVCs) are used in one of the most frequently performed invasive procedures for:¹



Central vascular access devices (CVADs) are used for necessary clinical interventions and are indicated for:²



More than 400 million PIVCs are sold in the US each year³

Up to 90% of hospitalized patients receive a PIVC¹

Contributing factors to complications...



Variation in policy, practice and training⁷
 Inappropriate device selection, placement, care and maintenance^{7,16}
 Variation in patients (>50% may be difficult intravenous access)⁸

Up to **50%** of catheters require replacement before the completion of therapy¹

Up to **25%** of central lines need to be removed prior to completion of therapy in pediatric patients⁴

More than 4 million CVADs are sold in the US each year³

Up to 1 in 4 patients who acquire a CLABSI will die⁵



Up to 28% of central line-associated blood stream infections (CLABSIs) are antimicrobial resistant in acute care settings⁶

Occlusions occur in up to **36%** of long-term central venous catheters⁹

72% of CLABSIs occur more than 5 days after¹⁰

On average, up to 82 patients acquire a CLABSI everyday¹¹

Vascular access complications cause substantial burdens

Economic



CLABSIs are estimated to cost the U.S. healthcare system \$45.814 per occurrence or ~\$1.9 billion every year.¹²

Unnecessary PIV restarts can cost a 200 bed hospital more than \$980,000 annually.¹

Patient



Multiple insertion attempts increase pain to the patient.¹³

Clinical



CLABSI patients have been associated with 2.27x greater risk of mortality than non-CLABSI patients.¹⁴

Treatment of CRBSIs can extend a patient's length of stay.¹²



Quality assurance and performance improvement plans may include multidisciplinary functions and an integrated multimodal approach such as:

- Education and training of healthcare workers¹⁵
- Products that support and align to industry best practices¹⁵
- Surveillance and feedback¹⁶

Care and maintenance of catheters should be the focus of performance improvement and quality assurance.¹⁶



The burden of PIVC insertion failure

Risk factors

contributing to PIVC insertion failure in adult and pediatric populations include:



Patient
Age, skin/vein characteristics



Procedural
Insertion site and catheter gauge



Clinician/Operator expertise¹⁷




Repeated PIVC insertion attempts can negatively impact **three key areas** of the healthcare system: clinical, economic and patient experience.

Three areas of impact



1 Clinical Impact

First-attempt insertion failure occurs up to:



Attempts per successful IV insertion occurs up to:

1.6 adults¹⁹ **2.1 pediatrics²⁰**

Repeated attempts can lead to:

- Vessel trauma²¹
- Venous depletion²¹
- Increased risk of catheter failure²¹

Less experienced/confident nurses and clinicians can:

- Be fearful of the risks of puncturing arteries or nerves²²
- Lead to first time insertion failure²³

2 Economic impact

First-attempt insertion failure can lead to increased costs.

\$28-\$35 per adult insertion²¹
\$43-\$47 per pediatric insertion²⁰

Failure to insert a PIVC on the first attempt increases procedural time.²⁴

3 Patient impact

Adult patients and caregivers of pediatric patients report that inadequate insertion competency impacted patient well-being, especially with multiple insertion attempts.²⁵

Each IV insertion can result in:

- Heightened anxiety²⁶⁻²⁹
- A sense of vulnerability²⁶
- Great pain²⁷⁻³⁰



Reducing PIVC insertion attempts through proper strategies and techniques is paramount to PIVC success rates and patient care.






Points to practice

Expected infusion therapy outcomes	10
Therapy considerations	11
Site selection	14
Device selection	19

Expected infusion therapy outcomes

The ability to achieve the expected outcomes of infusion therapy are influenced by a variety of factors, many of which are shown below.



Not intended to be comprehensive.
BD and the BD Logo are trademarks of Becton, Dickinson and Company. © 2020 BD. All rights reserved. BD-2432 (04/23)



Points to practice

Therapy considerations

Infusate considerations

Knowledge regarding medications and solutions to be administered through the PIVC, their purpose, properties, duration and potential complications can lead to informed decisions and influence expected outcomes.



When considering the therapies to be administered through a PIVC, clinicians should ask:

- How many medications and solutions are prescribed?
- Are they compatible with one another?
- What are the chemical properties of the infusates (pH and osmolarity)?
- How often will they be administered? How quickly do they need to be administered?
- Do they require any special equipment (e.g., power injector, filter, etc.)?



Check out our new BD Infusate Consideration Companion App



Your guide to Infusates



Download on the App Store

Scan this with your iPhone or iPad to download the free app from the Apple App Store



GET IT ON Google Play

Scan this with your Android phone or tablet to download the free app from the Google Play Store



Points to practice


Infusate considerations: pH and osmolarity

Consider the pH of an infusion therapy:

- Normal serum pH falls between 7.35–7.45.³⁵
- pH of infusates can vary widely between very acidic and very basic.³⁵
- The farther the pH of an infusate is from normal serum pH, the greater the risk of complications when that infusate is administered through a PIVC.³¹

Consider the osmolarity of an infusion therapy:


- The concentration of solutes in solution is written as mOsm/L.
- Normal serum osmolarity (use osmolarity when describing fluids within the body) is 280–295 mOsm/L.³⁵
 - A perfect isotonic solution will be 280–295 mOsm/L.
 - A near isotonic solution is between 240 and 450 mOsm/L.
- The farther the osmolarity of an infusate is from normal serum osmolarity, the greater the risk of complications when that infusate is administered through a PIVC.³²



Irritants and vesicants

Many drugs and solutions can be damaging to the venous endothelium irrespective of their osmolarity or pH.

Irritant drugs and solutions can cause itching, phlebitis, or reaction along the vessel.³³



Vesicant drugs can cause blistering, tissue sloughing, or necrosis when they escape from the normal vascular pathway.³⁴





Site selection

Anatomy of the vein

The vein is made up of three layers of tissues and fibers:

- The **tunica intima**, or inner layer, has a **smooth, elastic endothelial lining, allowing blood to easily flow** through the vein.
- The **tunica media**, or middle layer, is comprised of **smooth muscle and elastic tissue, allowing the vein to get wider or narrower** as blood passes through.
- The **tunica adventitia**, or outer layer, which is made up of **fibrous connective tissue, provides structure to the vein** and helps shape it.

The **valves** prevent the backward flow of blood and are formed by the tunica intima. They can prevent catheter advancement and may prevent blood withdrawal. Damage to the valves increases risk of thrombus formation.





Hemodilution


Hemodilution refers to the ability of blood to flow around the catheter at a sufficient rate to disperse infusates quickly and to prevent venous stasis. The larger the catheter diameter or the smaller the vein diameter, the less blood will be able to flow around it. This can cause complications due to the increased risk for coagulation and clot formation as well as the potential for caustic infusates to stay in one place for longer and cause irritation or degradation of the vessel wall.




You can help prevent these complications by selecting a vessel large enough to accommodate the device required for therapy, and by choosing a device with a small enough gauge to allow appropriate hemodilution in the target vein, as shown in the second image below.³⁵



Axillary vein
16 mm
150–350 mL/min



Upper arm basilic vein
8 mm
90–150 mL/min




Lower forearm cephalic vein
6 mm
0–40 mL/min



Peripheral vasculature³⁶

The veins of the hand and arm are generally used for peripheral IV catheter insertion. The forearm has increased likelihood of lasting the full length of therapy. Review the following commonly used insertion sites.

- 1 Cephalic vein:** The cephalic vein is a superficial vein that begins in the hand and ascends the outer side of the arm. Use caution at the radial wrist and at or above the antecubital fossa to avoid the risk of nerve damage.
- 2 Basilic vein:** The basilic vein begins in the hand and ascends the inner side of the arm. It is superficial up to halfway through the upper arm, where it becomes a deep vein. Use caution at or above the antecubital fossa to avoid the risk of nerve damage.
- 3 Median cubital vein:** The median cubital vein is often used for infusion therapies because it is easily accessible. **However**, it is also located at an area of high flexion, so stabilization is especially important. Use caution at or above the antecubital fossa to avoid the risk of nerve damage.
- 4 Hand:** Although easily accessible, the metacarpal, cephalic and basilic veins of the hand are associated with high rates of failure over time. Consider these veins for therapies less than 24 hours.



Insertion site considerations



Optimal vein conditions for an insertion site include the following qualities:³⁹

- Soft, straight and elastic
- Springy, easily palpated
- No visible valves
- Supported by intact, elastic skin
- Easily stabilized



Avoid the following:⁴⁰

- Dominant extremity
- Areas of flexion, including antecubital space and wrist
- Compromised areas
- Sites below recent venipunctures
- Areas where restraints are being used
- Lower extremities




Vascular visualization technology

There are various devices that use sound or light waves to allow for the location and identification of blood vessels and guide device insertion.⁴¹

1 Visible light devices: These devices function by providing transillumination of the peripheral veins. These devices aid in locating superficial veins in neonates. However, their usefulness in infants, older children and adults is limited due to the thickness of subcutaneous tissue and size of the arm circumference.⁴²

2 Near infrared (nIR) light technology: Near infrared works by shining a frequency of light onto the patient that is absorbed by blood but reflected by surrounding tissue. The reflected light forms an image, either on the patient's skin or on a screen, that shows the path of fluid-filled vessels. nIR light technology can aid in locating superficial veins and helping assess their pathway, including to identify bifurcations and valves.⁴³

3 Ultrasound: Ultrasound technology transmits sound waves through tissue. These waves bounce off of tissue in different ways depending on its density, and the image the reflected waves creates can be used to visualize anatomical structures under the skin. Ultrasound provides real-time guidance for catheter placement, unlike visible light devices and near infrared technology, which only identify the venous pathway.⁴⁴



Device selection

Size considerations

Different catheters work best for different infusion therapies. Choosing the right catheter accomplishes the following:⁴⁵



Has the **greatest** likelihood of **remaining intact** throughout the anticipated **length of therapy**



Accommodates the **requirements** of the therapy



Is the **least invasive** option possible



Utilizes the **fewest number of catheters** over the **course of the therapy**



Meets a **benefit versus risk** assessment

Catheter gauge	INS recommendation for use ⁴⁶
20–24 gauge	Most infusion therapies*
22–26 gauge	Neonates, pediatric patients, older adults and patients with limited venous options
large (>20)–gauge catheters	Consider when rapid fluid replacement is required
Fenestrated catheters	Consider for a contrast-based radiographic study

*Catheters larger than 20-gauge are more likely to cause phlebitis

Follow your facility policy for specific procedures regarding catheter size requirements.



Points to practice



Insertion considerations

Vein distention and stabilization

21

Venipuncture

22

Securement

23

Vein distention and stabilization

Vein distention:



Remember the following when applying a tourniquet:⁴⁷

- Apply 5–6 inches above the intended insertion site.
- It should be tight enough to impede venous flow while maintaining arterial circulation.
- Apply loosely or consider not using with patients at risk for bleeding, who have fragile skin/veins or have compromised circulation.

Additional methods to promote venous distention include the following:

- 1 Position the extremity lower than the heart for a minute or two.
- 2 Have the patient open and close their fist repeatedly.
- 3 Lightly stroke the vein downward or lightly tap on it.
- 4 Apply a controlled heat source to the extremity.

Vein stabilization:

- Anchors the vein
- Minimizes rolling
- Straightens curvy veins
- Facilitates catheter advancement

Stabilize by gently pulling the skin down and away from the insertion site. Maintain vein stabilization while advancing the catheter. Avoid contamination of the insertion site by re-palpating or when applying skin traction.






Venipuncture considerations

- 1** **Blood return** indicates the **needle tip has entered the vein**, but the **catheter tip may not be in the vein**.
- 2** Once **blood return is seen**, **lower the insertion angle** and **advance the needle and catheter as a unit**, slightly. This **helps ensure that the catheter tip is also within the vein** and **helps prevent puncturing through the back wall of the vein**.
- 3** **Puncture of the back wall** of the vein (also known as **through puncture** or **backwalling**) occurs when the **catheter/needle unit is inserted too far** or at **too great of an angle**.



Tips for success:

- Use only one catheter for each venipuncture attempt.
- Make no more than two insertion attempts per clinician.
- Release the tourniquet after the catheter is advanced off of the needle.



Securement


Proper securement is essential to the device's longevity. Effective securement provides the following benefits:⁴⁸

- Preserves the integrity of the PIVC.
- Minimizes catheter movement at the hub.
- Helps prevent catheter dislodgement and loss of access.



Points to practice:

- The insertion site should be readily visible for routine assessment and monitoring.⁴⁹
- The securement method should not impede vascular circulation.⁴⁹
- Avoid the use of rolled bandages, with or without elastic properties, as the primary securement method.⁵⁰
- Secure extensions, tubings and other add-on devices that could get caught and dislodge the PIVC.





PIVC management

Flushing

25

Dressing management

27

Flushing

Appropriate catheter flushing is important to decreasing complications and maintaining the patency of PIVCs.⁵¹

When should a catheter be flushed?

- 1 Prior to each infusion to assess catheter function.
- 2 After each infusion to clear the infused medication.

Refer to your institutional policy and manufacturer instructions for use for more information.

Tips for catheter flushing:⁴²



Use prefilled **single-use saline syringes** with preservative-free **0.9% sodium chloride** to reduce contamination risk and save time.



Use a minimum flush volume of **2x** the internal volume of the catheter system. **Larger volumes** may also remove more **fibrin, drug precipitate**, and **debris** from the lumen.



Consider using a **pulsatile technique**: Multiple short **1 mL boluses** with **brief pauses** may more effectively **remove solid deposits** from the lumen.




Use **positive pressure** techniques to **minimize blood reflux** upon disconnect.




Flushing

Fluid displacement

Fluid displacement is the movement of fluid upon disconnection when no clamping is performed.



Blood reflux can play a significant role in thrombotic catheter occlusions.⁵⁴




“Understanding NFC design and the appropriate clamping and disconnecting sequence for the various device types is necessary to ensure the best outcomes and prevent catheter-related complications.⁵³”



Dressing management

It's essential to secure the catheter in a way that will minimize catheter movement. If the catheter is moving constantly, it can cause mechanical phlebitis or create a risk for dislodgement. Consider the potential impacts for each of the following examples:

- 1 The label is covering the insertion site which could interfere with site assessment.
- 2 The extension tubing crosses the wrist joint, increasing the risk of dislodgement and impeding the patient's use of their hand.
- 3 The catheter and extension tubing are well secured. The needle-free connector can be easily accessed, and the insertion site is visible.





PIVC complications

Introduction to complications	29
Phlebitis	30
CABSI	32
Occlusion	34
Infiltration/extravasation	35
Dislodgement	37

Intro to complications

This tab will introduce common complications that may affect PIVCs:

- **Phlebitis**
- **CABSI**
- **Occlusion**
- **Infiltration/extravasation**
- **Dislodgement**

As you read through each complication, you'll see links allowing you to quickly find and review the relevant best practices.



Most vascular access complications can be prevented, and it starts by learning to recognize them early.



Dislodgment



Occlusion



Infection



Phlebitis



Infiltration/extravasation




Phlebitis

Phlebitis is inflammation of the intimal lining of the vessel.⁵⁵ This irritation can have several causes:⁵⁶

- 1 Mechanical:** manipulation during insertion, improper catheter stabilization, catheter size and tip location.
- 2 Chemical:** irritating solutions, inadequate hemodilution, skin antiseptics not allowed to dry prior to insertion.
- 3 Infectious:** breaks in aseptic technique, poor hand washing, inadequate cleaning for connectors, non-occlusive dressing.



- Select appropriate device for intended therapy.⁵⁷
- Manage vein/catheter ratio to allow adequate hemodilution.⁵⁷
- Stabilize catheter.⁵⁷
- Practice proper aseptic technique during insertion, access and maintenance.⁵⁷



Phlebitis (cont'd)

Identification and progression⁵⁸

Phlebitis follows a grading system which is used to assess the severity of the symptoms.

Grade 1 phlebitis presents as erythema at the insertion site with or without pain.

Grade 2 phlebitis may look similar to grade 1 phlebitis, but consists of pain at the access site with erythema and/or edema.

Grade 3 phlebitis consists of pain at the access site with erythema, streak formation, or a palpable venous cord.

Grade 4 phlebitis includes pain at the access site, erythema, streak formation, a palpable venous cord greater than 1 inch in length, and/or purulent drainage.

Intervention⁵⁹

When you identify phlebitis:



Remove the catheter.



Elevate the extremity.



Apply a warm compress to the site.



Provide analgesics as needed.



Consider other pharmacologic interventions such as anti-inflammatory agents.

Grade 1–2



Grade 3



Grade 4



CABSI

Catheter-associated bloodstream infection (CABSI) is the presence and growth of a pathogenic microorganism in the bloodstream originating from the catheter.⁶⁰

Extraluminal:⁶¹


- Seeding from a distant site
- Contaminated device
- Ingress of bacteria at the insertion site

Intraluminal:⁶¹

- Caused by the patient or healthcare worker
- Contaminated infusate (fluids infused through the catheter), caused by the manufacturer, pharmacy or other healthcare workers



- Continue education⁶²
- Practice proper aseptic technique during insertion, access and maintenance⁶²
- Disinfect skin properly⁶²
- Select site properly⁶²
- Remove VAD promptly when no longer needed⁶²



Identification

Common **signs** and **symptoms of infection** include:⁶³

- Fever
- Increased white blood cell count
- Redness around the insertion site
- Drainage from the insertion site
- Pain or tenderness associated with the catheter

Intervention⁶⁴




Remove a PIVC when:



Patient **develops symptoms of complication** and **failure** such as **infection** (e.g., erythema extending at least 1 cm from the insertion site, induration, exudate, fever with no other obvious source of infection).



Patient **reports any pain or tenderness** associated with the **catheter**.



Occlusion

Occlusion is the obstruction of the catheter lumen, preventing or limiting the ability to flush and/or administer solutions or withdraw blood.⁶⁵

Occlusion may be caused by:⁶⁶

- 1 Improper **flushing** and **locking** procedures
- 2 Infusion of **incompatible** medications
- 3 Mechanical causes, such as having the **clamp engaged** or **tubing kinked**



- Use vigorous pulsatile flush technique⁶⁷
- Disallow reflux of blood or fluid after flush⁶⁷
- Use adequate flush volumes⁶⁷
- Flush before or after each use or at least once every 24 hours⁶⁷

Identification

Common **signs** and **symptoms** of occlusion include:⁶⁸

- Frequent pump occlusion alarm
- Swelling/leaking from the insertion site
- Inability to withdraw blood or sluggish blood return
- Sluggish flow: resistance or inability to infuse or flush the catheter

Intervention⁶⁹

If occlusion is suspected, start by straightening any kinks in the catheter tubing and ensuring that the clamp is not engaged. If occlusion persists, remove catheter.



Infiltration and extravasation

Infiltration: Nonvesicant solution or medication is inadvertently administered into surrounding tissue.⁷⁰

Extravasation: Vesicant solution or medication inadvertently infiltrates into surrounding tissue.⁷¹

Risk factors for infiltration/extravasation:⁷²

Inappropriate insertion site	Irritants/vesicants get infused
Inadequate stabilization	Difficult access history
Medications that alter pain sensation	Diseases that produce change in vasculature
Use of deep veins with insufficient catheter length	Inability of patient to report symptoms

Also note that **early identification is essential** to limiting the **potential impact of infiltration** and **extravasation**. Monitor the **insertion site closely** if **risk factors** for **infiltration** or **extravasation** are present.



Select an appropriate vein for insertion^{73, 74}



Use ultrasound-guidance for deep and/or difficult to visualize veins^{73, 74}



Select appropriate catheter (based on type, gauge, length, etc.)^{73, 74}



Stabilize catheter adequately^{73, 74}



Assess site routinely to recognize infiltration early^{73, 74}



Infiltration and extravasation (cont'd)

Identification

The following are symptoms of infiltration or extravasation:⁷⁵

- Swelling of the extremity at or around the insertion site.
- Blanching and taut skin.
- Patient complaint of pain, burning, tightness or discomfort.

Intervention⁷⁶



Immediately stop the infusion upon identification of infiltration/extravasation injury and initiate appropriate interventions.



Aspirate for blood return from the peripheral catheter.



Do not flush the VAD, as this will inject additional medication into the tissue.



Disconnect the administration set and aspirate from the catheter, even though a very small amount of fluid may be retrieved.



Remove the peripheral catheter.



Avoid application of pressure.



Elevate the extremity.



Dislodgement

Dislodgement occurs when the catheter becomes dislodged from the vein. Dislodgement may be **incomplete** (meaning that the catheter tip is no longer in the vessel but is still in the patient) or **complete** (meaning that the catheter falls out entirely).⁷⁷

Dislodgement may be caused by:^{78, 79}

- 1 **Patient manipulation:** body habitus or arm movement causing VAD dislodgement
- 2 **Loosening:** due to inadequate stabilization and lack of proper securement
- 3 **Poor site selection**




- Avoid areas of flexion while selecting insertion sites.^{80, 81}
- Use a securement device ISD, SASS, TA or ASD, in addition to primary dressing to stabilize and secure VADs.^{80, 81}
- Use ISD, SSAS, TA or ASD for peripherally inserted catheters (PICCs) as an alternative to sutures.^{80, 81}

ISD - integrated securement device

SASS - subcutaneous anchor securement system

TA - tissue adhesion

ASD - adhesive securement device




Dislodgement (cont'd)

Identification⁸²

The following are signs and symptoms of dislodgement:

- Inability to continue infusion
- Increased external catheter length since last assessment
- Leakage at the insertion site
- If catheter moves out of the vein but not the skin, it can also result in:
 - Edema
 - Pain
 - Changes in skin color
 - Progress to blistering ulceration



Intervention⁸³



Remove the peripheral catheter.



Never re-advance a vascular access device that has been pulled out back into the vein as it could introduce bacteria.





References

References for tabs 1–5

40



References for tabs 1–5

Tab 1: Introduction to PIVC excellence

1. Helm RE, Klausner JD, Klemperer JD, Flint LM, Huang E. Accepted but unacceptable: peripheral IV catheter failure. *J Infus Nurs.* 2015;38(3):189–203.
2. Cheung E, Baerlocher MO, Asch M, Myers A. Venous access: a practical review for 2009. *Can Fam Physician.* 2009;55(5):494.
3. Tian M. Millenium Research Group, Inc. Vascular Access Devices: US Market Analysis. Toronto, Canada: Decision Resources Group, 2019. (Internal Data)
4. Ullman AJ, Marsh N, Mihala G, Cooke M, Rickard CM. Complications of central venous access devices: a systematic review. *Pediatrics.* 2015;136(5):e1331–e1344. doi: 10.1542/peds.2015-1507
5. Centers for Disease Control and Prevention. Vital signs: central line-associated blood stream infections. March 1, 2011. Accessed April 18, 2023. <https://www.cdc.gov/vitalsigns/pdf/2011-03-vitalsigns.pdf>.
6. Weiner LM, Fridkin SK, Aponte-Torres Z, et al. Vital signs: preventing antibiotic-resistant infections in hospitals — United States, 2014. *MMWR Morbidity and Mortal Weekly Report.* 2016;65(9):235–241. doi: 10.15585/mmwr.mm6509e1
7. Nickel B. Hiding in plain sight: peripheral intravenous catheter infections. *Crit Care Nurse.* 2020;40(5):57-66. doi: 10.4037/ccn2020439
8. Armenteros-Yeguas V, Garate-Echenique L, Tomas-Lopez MA, et al. Prevalence of dicult venous access and associated risk factors in highly complex hospitalized patients. *J Clin Nurs.* 2017;26(23-24):4267-4275. doi: 10.1111/jocn.13750
9. Baskin JL, Pui CH, Reiss U, et al. Management of occlusion and thrombosis associated with long-term indwelling central venous catheters. *Lancet.* 2009;374(9684):159–169.
10. Davis J. Central-line-associated bloodstream infection: comprehensive, data-driven prevention. *Pennsylvania Patient Safety Advisory.* 2011;8(3):100–5.
11. Centers for Disease Control. Bloodstream infection event (central line-associated bloodstream infection and non-central line associated bloodstream infection). January 2023. Accessed on April 18, 2023. https://www.cdc.gov/nhsn/pdfs/pscmanual/4psc_clabscurrent.pdf.
12. Zimlichman E, Henderson D, Tamir O, et al. Health care-associated infections: a meta-analysis of costs and financial impact on the US health care system. *JAMA Intern Med.* 2013;173(22):2039–2046.
13. Fields JM, Piela NE, Ku BS. Association between multiple IV attempts and perceived pain levels in the emergency department. *J Vasc Access.* 2014;15(6):514–518.
14. Stevens V, Geiger K, Concannon C, Nelson RE, Brown J, Dumyati G. Inpatient costs, mortality and 30–day re-admission in patients with central-line-associated bloodstream infections. *Clin Microbiol Infect.* 2014;20(5):O318–24.
15. Platt V, Osenkarski S. Improving vascular access outcomes and enhancing practice. *J Infus Nurs.* 2018;41(6):375–382.
16. O’Grady NP, Alexander M, Burns LA, et al. Guidelines for the prevention of intravascular catheter-related infections. *Clin Infect Dis.* 2011;52(9):e162–e193. doi: 10.1093/cid/cir257
17. Sabri A, Szalas J, Holmes KS, Labib L, Mussivand T. Failed attempts and improvement strategies in peripheral intravenous catheterization. *Bio-Medical Materials and Engineering.* 2013;23(1-2):93–108. doi: 10.3233/BME-120735



References for tabs 1–5 (cont'd)

18. Legemaat M, Carr PJ, van Rens RM, van Dijk M, Poslawsky IE, van de Hoogen A. Peripheral intravenous cannulation: complication rates in the neonatal population: a multicenter observational study. *J Vasc Access*. 2016;17(4):360-365. doi: 10.5301/jva.5000558
19. van Loon FHJ, van Hoof LWE, de Boer HD, et al. The modified A-DIVA scale as a predictive tool for prospective identification of adult patients at risk of a difficult intravenous access: a multicenter validation study. *J Clin Med*. 2019;8(2):144. doi: 10.3390/jcm8020144
20. Goff DA, Larsen P, Brinkley J, et al. Resource utilization and cost of inserting peripheral intravenous catheters in hospitalized children. *Hosp Pediatr*. 2013;3(3):185. doi: 10.1542/hpeds.2012-0089
21. Helm R, Klausner JD, Klemperer JD, Flint LM, Huang E. Accepted but unacceptable: peripheral IV catheter failure. *J Infus Nurs*. 2015;38(3):189-203. doi: 10.1097/NAN.0000000000000100
22. Forsberg A, Engström Å. Critical care nurses' experiences of performing successful peripheral intravenous catheterization in difficult situations. *J Vasc Nurs*. 2018;36(2):64-70. doi: 10.1016/j.jvn.2018.02.002
23. Carr PJ, Rippey JCR, Cooke ML, et al. Factors associated with peripheral intravenous cannulation first-time insertion success in the emergency department. A multicentre prospective cohort analysis of patient, clinician and product characteristics. *BMJ Open*. 2019;9(4):e022278. doi: 10.1136/bmjopen-2018-022278
24. van Loon FHJ, Puijin LAPM, Houterman S, Bouwman ARA. Development of the A-DIVA scale: a clinical predictive scale to identify difficult intravenous access in adult patients based on clinical observations. *Medicine*. 2016;95;16. doi: 10.1097/MD.00000000000003428
25. Cooke M, Ullman AJ, Ray-Barruel G, Wallis M, Corley A, Rickard CM. Not “just” an intravenous line: Consumer perspectives on peripheral intravenous cannulation (PIVC). An international cross-sectional survey of 25 countries. *PLoS One*. 2018;13(2):e0193436. doi: 10.1371/journal.pone.0193436
26. Robinson-Reilly M, Paliadelis P, Cruickshank M. Venous access: the patient experience. *Support Care Cancer*. 2016;24(3):1181-1187. doi: 10.1007/s00520-015-2900-9
27. Zempsky WT. Optimizing the management of peripheral venous access pain in children: evidence, impact, and implementation. *Pediatrics*. 2008;122;S121. doi: 10.1542/peds.2008-1055c
28. Walco GA. Needle pain in children: contextual factors. *Pediatrics*. 2008;122 Suppl 3:S125-S129. doi: 10.1542/peds.2008-1055d
29. Kennedy RM, Luhmann J, Zempsky WT. Clinical implications of unmanaged needle-insertion pain and distress in children. *Pediatrics*. 2008;122 Suppl 3:S130-S133. doi: 10.1542/peds.2008-1055e
30. Fields JM, Piela NE, Ku BS. Association between multiple IV attempts and perceived pain levels in the emergency department. *J Vasc Access*. 2014;15(6):514-518. doi: 10.5301/jva.5000282



References for tabs 1–5 (cont'd)

Tab 2: Points to practice

31. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S190 63.B.4.c
32. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S190 63.B.4.c
33. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S208
34. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S213
35. Gorski LA, *Manual of I.V. Therapeutics Evidence-based practice for Infusion Therapy*. 8th ed. Philadelphia, F.A. Davis Company. 2023. Pg 325; 4, cannula size
36. Gorski LA, *Manual of I.V. Therapeutics Evidence-based practice for Infusion Therapy*. 8th ed. Philadelphia, F.A. Davis Company. 2023.
39. Gorski LA, *Manual of I.V. Therapeutics Evidence-based practice for Infusion Therapy*. 8th ed. Philadelphia, F.A. Davis Company. 2023. Pg 325; 2, condition of the vein
40. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S81 27.IA.5
41. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S213
42. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S63 22.C
43. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S63 22.D
44. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S64 22.F
45. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S74 26.2
46. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S75 II.C.1-4

Tab 3: Insertion considerations

47. Gorski LA, *Manual of I.V. Therapeutics Evidence-based practice for Infusion Therapy*. 8th ed. Philadelphia, F.A. Davis Company. 2023. Pg 329; 6, tourniquet
48. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S109 38.1
49. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S109 38.2
50. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S110 38.I



References for tabs 1–5 (cont'd)

Tab 4: PIVC management

51. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S113 41.1, 41.2
52. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S113 41.C
53. MacElory B. Needleless connectors: improving patient care while keeping health care professionals safe. *J Infus Nurs*. 2018; Jan/Feb: 20–22
54. Williams A. Catheter Occlusion in Home Infusion: The Influence of Needleless Connector Design on Central Catheter Occlusion. *J Infus Nurs*. 2018;41(1):52-57.

Tab 5: PIVC complications

55. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S211
56. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S138 46.B
57. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S138 46.C
58. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S139 46.E
59. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S138 46.C
60. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S203
61. Safdar N, Maki DG. The pathogenesis of catheter-related bloodstream infection with noncuffed short-term central venous catheters. *Intensive Care Med*. 2004;30(1):62-67. doi: 10.1007/s00134-003-2045-z
62. Agency for Healthcare Research and Quality (AHRQ). Toolkit for Using the AHRQ Quality Indicators. PSI07 Central Venous Catheter-Related Bloodstream Infections PDF. Published 2011. Updated March 2017. Accessed July 22, 2019. http://www.ahrq.gov/sites/default/files/wysiwyg/professionals/systems/hospital/qitoolkit/combined/d4d_
63. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S154 50.B
64. Gorski LA, *Manual of I.V. Therapeutics Evidence-based practice for Infusion Therapy*. 8th ed. Philadelphia, F.A. Davis Company. 2023.
65. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S210
66. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S149 49.A
67. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S113 41
68. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S150 49.B
69. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S114 41.D.2
70. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S207



References for tabs 1–5 (cont'd)

71. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S206
72. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S142 47.C
73. Alexander M, Corrigan A, Gorski L, Hankins J, Perucca R. *Infusion Nursing An Evidence Based Approach*. 3rd ed. St. Louis, MO: Saunders Elsevier; 2010.
74. Gorski L, Hadaway L, Hagle ME, McGoldrick M, Orr M, Doellman D. (2021). *Infusion therapy: Standards of practice*. INS Digital Press. S81 27.A
75. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S143 47.D.3-5
76. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S143 47.E
77. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S204
78. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S165 54.A
79. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S109 38.A
80. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S109 38.A
81. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S81 27.I.A.5
82. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S165 54.I
83. Gorski LA, Hadaway L, Hagle ME, et al. *Infusion Therapy Standards of Practice*, 8th ed. Infusion Nurses Society; 2021. S110 38.L

