VENTRALEX™ ST Hernia Patch

Self-Expanding Bioresorbable Coated Permanent Mesh Patch
and Strap for Soft Tissue Reconstruction
with *SORBAFLEX™ Memory Technology*

Instructions for Use

Mode d'emploi • Gebrauchsanweisung
Istruzioni per l'uso • Instrucciones de uso
Gebruiksaanwijzing • Instruções de utilização
Οδηγίες χρήσης • Brugervejledning
Bruksanvisning • Käyttöohjeet
Bruksanvisning • Instrukcja użycia
Használati útmutató • Návod k použití
Kullanım Talimatları • 使用说明
사용 설명서 • Инструкции по применению

 Single Use

STERILE | EO

 Do not resterilize

R_x only

 25°C
(77°F)
2°C
(36°F)

BAIRD

DAVOL INC.

VENTRALEX™ ST Hernia Patch

Self-Expanding Bioresorbable Coated Permanent Mesh Patch
and Strap for Soft Tissue Reconstruction
with *SORBAFLEX™ Memory Technology*

DESCRIPTION:

VENTRALEX™ ST Hernia Patch is a self-expanding bioresorbable coated, partially absorbable, sterile prosthesis, containing 2 distinct layers stitched with PTFE monofilament, forming a positioning pocket and strap. The top layer is monofilament polypropylene mesh and the bottom layer is SEPRAMESH™ IP Composite mesh. SEPRAMESH™ IP Composite is co-knitted using polypropylene (PP) and absorbable polyglycolic acid (PGA) fibers to result in a two-sided mesh with a PP surface and a PGA surface. The mesh is coated on the PGA surface with a bioresorbable, chemically modified sodium hyaluronate (HA), carboxymethylcellulose (CMC) and polyethylene glycol (PEG) based hydrogel (see Figure 1).

The fascial side of the mesh allows a prompt fibroblastic response through the interstices of the mesh, allowing for tissue ingrowth into the mesh. The visceral side of the mesh is a bioresorbable coating, separating the mesh from underlying tissue and organ surfaces to minimize tissue attachment to the mesh. Shortly after placement, the biopolymer coating becomes a hydrated gel that is resorbed from the site in less than 30 days.

A depth marker on the positioning strap for the small VENTRALEX™ ST Hernia Patch (**REF 5950007**) is designed to facilitate placement of the small patch through a trocar.


Figure 1

The device contains SORBAFLEX™ Memory Technology which provides memory and stability to the device, facilitating ease of initial insertion, proper placement, and fixation of the device. The SORBAFLEX™ Memory Technology is comprised of an extruded polydioxanone (PDO) absorbable monofilament that is contained within a knitted polypropylene mesh tube. The SORBAFLEX™ PDO monofilament fully degrades in vivo by means of hydrolysis. The PDO monofilament has been found to elicit an inflammatory response during absorption. Absorption is essentially complete in 6-8 months.

INDICATIONS:

The VENTRALEX™ ST Hernia Patch is indicated for use in the reinforcement of soft tissue, where weakness exists, in procedures involving soft tissue repair, including repair of hernias and deficiencies caused by trocars.

CONTRAINDICATIONS:

1. Do not use the VENTRALEX™ ST Hernia Patch in infants or children, whereby future growth will be compromised by use of such mesh material.
2. Do not use the VENTRALEX™ ST Hernia Patch for the reconstruction of cardiovascular defects.
3. Literature reports that there may be a possibility for adhesion formation when the polypropylene is placed in contact with the bowel or viscera.

WARNINGS:

1. The device is supplied sterile. Inspect the packaging to be sure it is intact and undamaged prior to use.
2. This device is for single use only. Do not resterilize. Product should be used once exterior foil pouch has been opened. Do not store for later use. Unused portions of this prosthesis should be discarded.
3.  Do not cut or reshape any portion of the VENTRALEX™ ST Hernia Patch (as this could impact its effectiveness), except for the polypropylene positioning strap. Care should be taken not to cut or nick the SORBAFLEX™ PDO monofilament. If the SORBAFLEX™ PDO monofilament is cut or damaged during insertion or fixation, additional complications may include bowel or skin perforation and infection.
4. Follow proper folding techniques for all patches as described in these Instructions for Use as other folding techniques may potentially compromise the SORBAFLEX™ PDO monofilament.
5. Ensure proper orientation; the bioresorbable coated side of the prosthesis should be oriented against the bowel or sensitive organs. Do not place the polypropylene side against the bowel. There may be a possibility for adhesion formation when the mesh (including strap) is placed in direct contact with the bowel or viscera.
6. To ensure a strong repair, the prosthesis should be secured with tacks or sutures through the polypropylene mesh straps or positioning pocket.
7. Excess positioning strap material above the fixation line must be cut off and discarded to eliminate excess material from remaining in the body.
8. When used to repair deficiencies caused by trocars, the device should be used under endoscopic guidance or direct visualization.
9. The use of any permanent mesh or patch in a contaminated or infected wound could lead to fistula formation and/or extrusion of the prosthesis.
10. If an infection develops, treat the infection aggressively. Consideration should be given regarding the need to remove the prosthesis. An unresolved infection may require removal of the prosthesis.
11. To prevent recurrences when repairing hernias, the prosthesis should be large enough to extend beyond the margins of the defect.

PRECAUTIONS:

1. Please read all instructions prior to use.
2. Only physicians qualified in the appropriate surgical techniques should use this prosthesis.
3.  Care should be taken not to cut or nick the SORBAFLEX™ PDO monofilament during fixation.
4. The safety and effectiveness of VENTRALEX™ ST Hernia Patch has not been evaluated in clinical studies in the presence of malignancies in the abdominopelvic cavity.

ADVERSE REACTIONS:

Possible complications include seroma, adhesions, hematoma, inflammation, extrusion, fistula formation, infection, allergic reaction, and recurrence of the hernia or soft tissue defect. If the SORBAFLEX™ PDO monofilament is cut or damaged during insertion or fixation, additional complications may include bowel or skin perforation and infection.

INSTRUCTIONS FOR USE:

Preparation

It is recommended that VENTRALEX™ ST Hernia Patch be completely immersed in sterile saline for 1-3 seconds immediately prior to placement in order to maximize the flexibility of the prosthesis. The safety and effectiveness of the VENTRALEX™ ST Hernia Patch in combination with solutions other than saline have not been tested.

Surface Orientation

It is extremely important that this product be oriented correctly to function as intended. The visceral side of the VENTRALEX™ ST Hernia Patch is designed to temporarily separate tissue surfaces and minimize tissue attachment to the mesh. Place the bioresorbable coated side of the prosthesis against those surfaces where minimal tissue attachment is desired, i.e. against bowel or other visceral structures. The uncoated polypropylene mesh side should face the surface where tissue ingrowth is desired. The uncoated polypropylene mesh surface should never be placed against the bowel or other visceral structures.

Hernia Repair

Select the patch that is approximately twice the size of the hernia defect. Fold the patch parallel to the opening between the straps with the bioresorbable coated side of the prosthesis facing out for insertion into the defect (see Figure 2). The free or unattached end of the positioning strap remains outside of the body during placement. The positioning strap is manipulated to facilitate the proper positioning of the patch. Gently pulling up on the positioning strap will allow the patch to flatten itself against the underside of the abdominal wall. The mesh positioning strap can be pulled apart to gain access to the inner positioning pocket.


Figure 2

Fixation

To secure the patch, the mesh positioning strap is pulled up and apart and the device is then secured to the margins of the defect through the strap or positioning pocket. Bard permanent or absorbable fixation devices or nonabsorbable monofilament sutures are recommended to properly secure the prosthesis. If other fixation devices are used, they must be indicated for use in hernia repair. Excess positioning strap material above the fixation line must be cut off and discarded. The method of securing the prosthesis should be determined by surgeon preference and the nature of the reconstruction to provide for adequate tissue fixation and to prevent reherniation.

Trocar Defect Repair with the Small VENTRALEX™ ST Hernia Patch

Fold the VENTRALEX™ ST Hernia Patch in half with the bioresorbable coated side of the prosthesis facing out and hold it in the folded position with an endoscopic grasper. Deploy the prosthesis all the way through the trocar and into the intra-abdominal space. For a reference point, the small VENTRALEX™ ST Hernia Patch is constructed with a depth marker (blue line) on the positioning strap 8 inches (203 mm) above the patch (see Figure 3). A minimum 10 mm diameter trocar should be used to introduce the patch. Release the prosthesis from the grasper when the patch is completely through the trocar. The SORBAFLEX™ PDO monofilament allows the patch to open and lay flat in the intra-abdominal space. Pull up on the positioning strap to hoist the prosthesis up to the distal end of the trocar. Pull out the trocar and pass it over the strap. Pull up firmly on the positioning strap to make sure that the prosthesis is completely covering the defect. While pulling up on the positioning strap, suture close the anterior fascia, catching the positioning strap in between the margins of the fascia. Cut off excess positioning strap immediately above the suture and discard.


Figure 3

Trocar Defect Repair with the Medium or Large VENTRALEX™ ST Hernia Patch

Remove the trocar and follow the instructions found in the "Hernia Repair" and "Fixation" sections above.

TRACEABILITY:

A traceability label which identifies the type, size and lot number of the prosthesis is attached to every package. This label should be affixed to the patient's permanent medical record to clearly identify the device which was implanted. If you experience a product failure, please contact Davol, Inc. at 1-800-556-6275 for instructions on returning the product.

Bard, Davol, SorbaFlex and Ventralex are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Sepramesh is a trademark and/or registered trademark of Genzyme Corporation licensed to C. R. Bard, Inc. or an affiliate.

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Copyright © 2011 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2011 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2011 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2011 C. R. Bard, Inc. Reservados todos los derechos.
Copyright © 2011, C. R. Bard, Inc. Alle rechten voorbehouden.
Copyright © 2011 C. R. Bard, Inc. Alle rettigheder forbeholdes.