

BD Product Name: **BD Pyxis™ MedStation™ 3500**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

### Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020

The patches include updates that may adversely affect Windows Server 2008 32 Bit BD Pyxis™ MedStation 3500 Consoles. These patches will not be and should not be installed on BD Pyxis™ MedStation 3500 Consoles at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering

| Patch Name  | Description   | Patch ID                  | Notes |
|---|---|---------------------------|-------|
| 2020-07 Security Only Quality Update for Windows Server 2008 for x86-based Systems                      | This security update includes quality improvements. Key changes include: Security updates to Windows App Platform and Frameworks, Windows Apps, Windows Graphics, Windows Fundamentals, Windows Kernel, and Windows SQL components.                       | <a href="#">KB4565529</a> | None  |
| Security Only Update for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008 SP2                | The security update addresses the vulnerability by correcting how .NET Framework validates the source markup of XML content.  | <a href="#">KB4566469</a> | None  |
| 2020-07 Security and Quality Rollup for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008 SP2 | The security update addresses the vulnerability by correcting how .NET Framework validates the source markup of XML content.  | <a href="#">KB4566520</a> | None  |
| 2020-07 Security Monthly Quality Rollup for Windows Server 2008 for x86-based Systems                   | This security update includes improvements and fixes that were a part of update KB4561670 (released June 9, 2020) and addresses the following issues:<br>Security updates to Windows App Platform and Frameworks, Windows Apps, Windows Graphics, Windows | <a href="#">KB4565536</a> | None  |



|  |   |                           |      |
|--|---|---------------------------|------|
|  | Fundamentals, Windows Kernel, the Microsoft Scripting Engine, and Windows SQL components.   |                           |      |
| 2020-07 Cumulative Security Update for Internet Explorer 9 for Windows Server 2008 for x86-based systems | This security update resolves vulnerabilities in Internet Explorer.   | <a href="#">KB4565479</a> | None |
| 2020-05 Security Only Update for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008 SP2         | The update addresses the vulnerability by correcting how .NET Framework activates COM objects.  | <a href="#">KB4556406</a> | None |
| 2018-10 Security Monthly Quality Rollup for Windows Server 2008 for x86-based Systems                    | This security update includes improvements and fixes that were a part of update KB4458315 (released September 20, 2018) and addresses the following issues:<br>Addresses an issue in which all guest virtual machines running Unicast NLB fail to respond to NLB requests after the virtual machines restart.<br>Security updates to Windows Media Player, Microsoft Office Graphics, Microsoft Graphics Component, Windows Storage and Filesystems, and the Microsoft JET Database Engine. | <a href="#">KB4463097</a> | None |
| 2020-06 Security Only Quality Update for Windows Server 2008 for x86-based Systems                       | This security update includes quality improvements. Key changes include: Security updates to Windows App Platform and Frameworks, Windows Kernel, Microsoft Graphics Component, Windows Input and Composition, Windows Media, Windows Cloud Infrastructure, Windows Peripherals, Windows File Server and Clustering, Windows Hybrid Storage Services, and the Microsoft JET Database Engine.  | <a href="#">KB4561645</a> | None |

|  |  |                                  |             |
|--|--|----------------------------------|-------------|
| <p>2019-03 Security Monthly Quality Rollup for Windows Server 2008 for x86-based Systems</p> | <p>This security update includes improvements and fixes that were a part of update KB4487022 (released February 19, 2019) and addresses the following issues:</p> <p>Addresses an issue in which the abbreviated Japanese Era name did not refer to the typographic ligature.</p> <p>Addresses an issue that fails to recognize the first character of the Japanese era name as an abbreviation and may cause date parsing issues.</p> <p>Security updates to Internet Explorer, Windows App Platform and Frameworks, Windows Server, Windows Hyper-V, Windows Storage and Filesystems, Windows Fundamentals, Windows Kernel, Windows MSXML, and the Microsoft JET Database Engine.</p>  | <p><a href="#">KB4489880</a></p> | <p>None</p> |
| <p>2019-01 Security Monthly Quality Rollup for Windows Server 2008 for x86-based Systems</p> | <p>This security update addresses the following issues:</p> <p>Provides protections against an additional subclass of speculative execution side-channel vulnerability known as Speculative Store Bypass (CVE-2018-3639) for AMD-based computers. These protections aren't enabled by default. For Windows Server guidance, follow the instructions in KB4072698. Use this guidance document to enable mitigations for Speculative Store Bypass (CVE-2018-3639).</p> <p>Additionally, use the mitigations that have already been released for Spectre Variant 2 (CVE-2017-5715) and Meltdown (CVE-2017-5754).</p> <p>Addresses an issue that affects PowerShell remoting loop back using non-administrator accounts. For more details, see Windows Security change affecting PowerShell.</p> | <p><a href="#">KB4480968</a></p> | <p>None</p> |

|   |  |                           |      |
|---|--|---------------------------|------|
|   | Security updates to Windows Kernel, Windows Storage and Filesystems, Windows Wireless Networking, and the Microsoft JET Database Engine.   |                           |      |
| 2019-05 Security Monthly Quality Rollup for Windows Server 2008 for x86-based Systems | <p>This security update includes improvements and fixes that were a part of update KB4493460 (released April 25, 2019) and addresses the following issues:</p> <p>Addresses an issue that prevents the Microsoft Visual Studio Simulator from starting.</p> <p>Addresses an issue that may prevent applications that rely on unconstrained delegation from authenticating after the Kerberos ticket-granting ticket (TGT) expires (the default is 10 hours).</p> <p>Security updates to Windows App Platform and Frameworks, Microsoft Graphics Component, Windows Storage and Filesystems, Windows Cryptography, Windows Kernel, Windows Server, and the Microsoft JET Database Engine.</p> | <a href="#">KB4499149</a> | None |
| 2019-07 Security Monthly Quality Rollup for Windows Server 2008 for x86-based Systems | <p>This security update includes improvements and fixes that were a part of update KB4503271 (released June 20, 2019) and addresses the following issues:</p> <p>Provides protections against a variant (CVE-2019-1125) of the Spectre Variant 1 speculative execution side channel vulnerability.</p> <p>Security updates to Windows Server, Microsoft Graphics Component, Windows Shell, Windows Input and Composition, and Windows Kernel.</p>  | <a href="#">KB4507452</a> | None |

BD Product Name: **BD Pyxis™ Supply**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

### Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020 will not install on Windows server 2008 32-bit BD Pyxis™ Supply v9.4.

The patches include updates that may adversely affect 32-bit BD Pyxis™ Supply v9.4. These patches will not be and should not be installed on 32-bit BD Pyxis™ Supply v9.4 at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

| Patch Name  | Description  | Patch ID                  | Notes |
|---|--|---------------------------|-------|
| 2020-05 Security Only Update for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008 SP2        | The update addresses the vulnerability by correcting how .NET Framework activates COM objects.                               | <a href="#">KB4556406</a> | None  |
| 2020-07 Security and Quality Rollup for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008 SP2 | The security update addresses the vulnerability by correcting how .NET Framework validates the source markup of XML content. | <a href="#">KB4566520</a> | None  |
| 2020-07 Security Only Update for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008 SP2        | The security update addresses the vulnerability by correcting how .NET Framework validates the source markup of XML content. | <a href="#">KB4566469</a> | None  |



BD Product Name: **BD Care Coordination Engine**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

### Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020 will not install on Windows server 2008 64-bit BD Pyxis™ Care Coordination Engine DOD.

The patches include updates that may adversely affect 64-bit BD Pyxis™ Care Coordination Engine DOD. These patches will not be and should not be installed on 64-bit BD Pyxis™ Care Coordination Engine DOD at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

| Patch Name   | Description   | Patch ID                  | Notes |
|--|---|---------------------------|-------|
| 2020-02 Security Monthly Quality Rollup for Windows Server 2008 for x64-based System | This security update includes improvements and fixes that were a part of update KB4534303 (released January 14, 2020) and addresses the following issues:<br>Security updates to Microsoft Graphics Component, Windows Input and Composition, Windows Shell, Windows Fundamentals, Windows Cryptography, Windows Hyper-V, Windows Core Networking, Windows Peripherals, Windows Network Security and Containers, Windows Storage and Filesystems, the Microsoft Scripting Engine, and Windows Server. | <a href="#">KB4537810</a> | None  |
| 2020-02 Security Only Quality Update for Windows Server 2008 for x64-based Systems   | This security update includes quality improvements. Key changes include:<br>Security updates to Microsoft Graphics Component, Windows Input and Composition, Windows Shell, Windows Fundamentals, Windows Cryptography, Windows Hyper-V, Windows Core Networking, Windows Peripherals, Windows  | <a href="#">KB4537822</a> | None  |



|   |   |                           |      |
|---|---|---------------------------|------|
|   | Network Security and Containers, Windows Storage and Filesystems, and Windows Server.   |                           |      |
| 2020-05 Security and Quality Rollup for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008 SP2 for x64 | The update addresses the vulnerability by correcting how .NET Framework activates COM objects.  | <a href="#">KB4556402</a> | None |
| 2020-07 Security Only Update for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008 SP2 for x64        | The security update addresses the vulnerability by correcting how .NET Framework validates the source markup of XML content.  | <a href="#">KB4566469</a> | None |
| 2020-07 Security Only Update for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008 SP2 for x64        | The security update addresses the vulnerability by correcting how .NET Framework validates the source markup of XML content.  | <a href="#">KB4566469</a> | None |
| Windows Malicious Software Removal Tool   | This tool prevalent malicious software (including Blaster, Sasser, and Mydoom) and helps remove any infection that is found.  | <a href="#">KB890830</a>  | None |
| Security Update for Microsoft Silverlight   | This security update to Silverlight includes fixes outlined in KB 4013867. This update is backward compatible with web applications built using previous versions of Silverlight. | <a href="#">KB4013867</a> | None |

BD Product Name: **BD Alaris™ Systems Manager**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

### Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020 will not install on Windows server 2008 R2 64-bit Alaris™ Systems Manager

The patches include updates that may adversely affect 64-bit BD Alaris™ Systems Manager. These patches will not be and should not be installed on 64-bit BD Alaris™ Systems Manager at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

| Patch Name                  | Description  | Patch ID                  | Notes |
|-----------------------------|--|---------------------------|-------|
| Security Update for Windows | This security update includes quality improvements. Key changes include: Updates the 2020 start date for daylight saving time (DST) in the Kingdom of Morocco. For more information, see KB4557900.<br>Security updates to Windows App Platform and Frameworks, Windows Apps, Windows Input and Composition, Windows Kernel, Internet Information Services, Windows Network Security and Containers, Windows Active Directory, the Microsoft JET Database Engine, and Windows Storage and Filesystems. | <a href="#">KB4556843</a> | None  |
| Security Update for Windows | This security update includes quality improvements. Key changes include: Security updates to Windows Kernel, Windows App Platform and Frameworks, Windows Graphics, Windows Media, Windows Cloud Infrastructure, Windows Fundamentals, Windows Core Networking, and the Microsoft JET Database Engine.   | <a href="#">KB4550965</a> | None  |





|   |  |                           |      |
|---|--|---------------------------|------|
| Security Update for Windows   | This security update includes quality improvements. Key changes include: Security updates to Windows App Platform and Frameworks, Microsoft Graphics Component, Windows Fundamentals, Windows Authentication, Windows Cryptography, Windows Core Networking, Windows Storage and Filesystems, Windows Peripherals, and Windows Server .  | <a href="#">KB4541500</a> | None |
| 2020-02 Security Only Quality Update for Windows Server 2008 R2 forx64-based Systems                                  | This security update includes quality improvements. Key changes include: Security updates to Microsoft Graphics Component, Windows Input and Composition, Windows Media, Windows Shell, Windows Fundamentals, Windows Cryptography, Windows Hyper-V, Windows Core Networking, Windows Peripherals, Windows Network Security and Containers, Windows Storage and Filesystems, and Windows Server. | <a href="#">KB4537813</a> | None |
| 2020-03 Security Only Quality Update for Windows Server 2008 R2 forx64-based Systems                                  | This security update includes quality improvements. Key changes include: Security updates to Windows App Platform and Frameworks, Microsoft Graphics Component, Windows Fundamentals, Windows Authentication, Windows Cryptography, Windows Core Networking, Windows Storage and Filesystems, Windows Peripherals, and Windows Server .  | <a href="#">KB4541500</a> | None |
| 2020-05 Extended Security Updates (ESU) Licensing Preparation Package for Windows Server 2008 R2 forx64-based Systems | This update provides the complete set of licensing changes to enable installation of the ESU MAK add-on key, which is one of the steps to prepare for installation of Extended Security Updates. (For the full set of steps, please see KB4522133). A reboot is required after installing this update.   | <a href="#">KB4538483</a> | None |
| 2020-05 Security Only Quality Update for Windows Server 2008 R2 forx64-based Systems                                  | This security update includes quality improvements. Key changes include:   | <a href="#">KB4556843</a> | None |

|   |  |                           |      |
|---|--|---------------------------|------|
|   | <p>Updates the 2020 start date for daylight saving time (DST) in the Kingdom of Morocco. For more information, see KB4557900.</p> <p>Security updates to Windows App Platform and Frameworks, Windows Apps, Windows Input and Composition, Windows Kernel, Internet Information Services, Windows Network Security and Containers, Windows Active Directory, the Microsoft JET Database Engine, and Windows Storage and Filesystems.</p>             |                           |      |
| 2020-06 Security Only Quality Update for Windows Server 2008 R2 forx64-based Systems                        | <p>This security update includes quality improvements. Key changes include: Security updates to Windows App Platform and Frameworks, Windows Kernel, Microsoft Graphics Component, Windows Input and Composition, Windows Media, Windows Shell, Windows Cloud Infrastructure, Windows Fundamentals, Windows Storage and Filesystems, Windows File Server and Clustering, Windows Hybrid Storage Services, and the Microsoft JET Database Engine.</p> | <a href="#">KB4561669</a> | None |
| 2020-07 Cumulative Security Update for Internet Explorer 11 for Windows Server 2008 R2 forx64-based systems | <p>This security update resolves vulnerabilities in Internet Explorer.</p>   | <a href="#">KB4565479</a> | None |
| 2020-07 Security Monthly Quality Rollup for Windows Server 2008 R2 forx64-based Systems                     | <p>This security update includes improvements and fixes that were a part of update KB4561643 (released June 9, 2020) and addresses the following issues:</p> <p>Security updates to Windows App Platform and Frameworks, Windows Apps, Windows Graphics, Windows Input and Composition, Windows Fundamentals, Windows Kernel, Windows Remote Desktop, Internet Explorer, the Microsoft Scripting Engine, and Windows SQL components.</p>             | <a href="#">KB4565524</a> | None |

|   |   |                           |      |
|---|---|---------------------------|------|
| 2020-07 Security Only Quality Update for Windows Server 2008 R2 forx64-based Systems                        | This security update includes quality improvements. Key changes include: Security updates to Windows App Platform and Frameworks, Windows Apps, Windows Graphics, Windows Input and Composition, Windows Fundamentals, Windows Kernel, Windows Remote Desktop, and Windows SQL components.  | <a href="#">KB4565539</a> | None |
| 2020-04 Security Only Quality Update for Windows Server 2008 R2 forx64-based Systems                        | This security update includes quality improvements. Key changes include: Security updates to Windows Kernel, Windows App Platform and Frameworks, Windows Graphics, Windows Media, Windows Cloud Infrastructure, Windows Fundamentals, Windows Core Networking, and the Microsoft JET Database Engine.  | <a href="#">KB4550965</a> | None |
| Security Update for Windows   | The security update addresses the vulnerability by correcting how .NET Framework validates the source markup of XML content.  | <a href="#">KB4565579</a> | None |
| 2020-02 Cumulative Security Update for Internet Explorer 11 for Windows Server 2008 R2 forx64-based systems | This security update resolves vulnerabilities in Internet Explorer.   | <a href="#">KB4537767</a> | None |
| 2020-02 Security Monthly Quality Rollup for Windows Server 2008 R2 forx64-based Systems                     | This security update includes improvements and fixes that were a part of update KB4534310 (released January 14, 2020) and addresses the following issues:<br>Security updates to Internet Explorer, Microsoft Graphics Component, Windows Input and Composition, Windows Media, Windows Shell, Windows Fundamentals, Windows Cryptography, Windows Hyper-V, Windows Core Networking, Windows Peripherals, Windows Network Security and Containers, Windows Storage and Filesystems, the Microsoft Scripting Engine, and Windows Server. | <a href="#">KB4537820</a> | None |
| 2020-09 Security Monthly Quality Rollup for Windows   | This security update includes improvements and fixes that were a part of update KB4571729 (released   | <a href="#">KB4577051</a> |      |

|  |   |                                  |  |
|--|---|----------------------------------|--|
| <p>Server 2008 R2 for x64-based Systems</p>  | <p>August 11, 2020) and addresses the following issues:</p> <p>Updates time zone information for Yukon, Canada.</p> <p>Addresses a security vulnerability issue with user proxies and HTTP-based intranet servers. After you install this update, HTTP-based intranet servers cannot leverage a user proxy to detect updates by default. Scans that use these servers will fail if the clients do not have a configured system proxy. If you must leverage a user proxy, you must configure the behavior by using the Windows Update policy “Allow user proxy to be used as a fallback if detection using system proxy fails.”</p> <p>This change does not affect customers who secure their Windows Server Update Services (WSUS) servers that use the Transport Layer Security (TLS) or Secure Sockets Layer (SSL) protocols. For more information, see Improving security for devices receiving updates via WSUS.</p> <p>Security updates to Windows App Platform and Frameworks, Windows Graphics, Windows Media, Windows Cloud Infrastructure, Windows Authentication, Windows Cryptography, Windows Kernel, Windows Hybrid Cloud Networking, Windows Peripherals, Windows Storage and Filesystems, Windows Network Security and Containers, the Microsoft Scripting Engine, and Windows SQL components.</p> |                                  |  |
| <p>2020-09 Security Only Quality Update for Windows Server 2008 R2 for x64-based Systems</p> | <p>This security update includes quality improvements. Key changes include: Security updates to Windows App Platform and Frameworks, Windows Graphics, Windows Media, Windows Cloud Infrastructure, Windows Authentication, Windows Cryptography, Windows Kernel,</p>   | <p><a href="#">KB4577053</a></p> |  |

|  |   |                           |  |
|--|---|---------------------------|--|
|  | Windows Hybrid Cloud Networking, Windows Peripherals, Windows Storage and Filesystems, Windows Network Security and Containers, and Windows SQL components. |                           |  |
| 2020-09 Security and Quality Rollup for .NET Framework 3.5.1, 4.5. 2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2, 4.8 for Windows Server 2008 R2 for x64 | This patches improves and protect from security issues.   | <a href="#">KB4576628</a> |  |
| 2020-08 Security Monthly Quality Rollup for Windows Server 2008 R2 forx64-based Systems  | This security update includes improvements and fixes that were a part of update KB4565524 (released July 14, 2020)  | <a href="#">KB4571729</a> |  |

BD Product Name: **BD Pyxis™ MedStation™ 4000**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).



BD Product Name: **BD Pyxis™ MedStation™ ES**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).



BD Product Name: **BD Pyxis™ Anesthesia Station ES**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).





BD Product Name: **BD Pyxis™ Anesthesia system 3500**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).



BD Product Name: **BD Pyxis™ Security Module**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).



BD Product Name: **BD Pyxis™ Anesthesia Station 4000**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).



BD Product Name: **BD HealthSight™ Viewer**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).



BD Product Name: **BD Pyxis™ Order Viewer**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).



BD Product Name: **BD Pyxis™ Logistics**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).



BD Product Name: **BD Pyxis™ Data Agent**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).



BD Product Name: **Identity Management**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).





BD Product Name: **BD Pyxis™ CIISafe™**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).



BD Product Name: **BD Pyxis™ CUBIE™ Replenishment System**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).



BD Product Name: **BD Pyxis™ PharmoPack™**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).



BD Product Name: **BD Pyxis™ Barcode Labeling**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).



BD Product Name: **BD Pyxis™ IV Prep**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).



BD Product Name: **BD Health Sight™ Patient Association**

Date of Critical or Security Patches: September 2020

Abstract: Critical or Security Patches – September 2020

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for September 2020. Patches delivered in the month were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

