REF 50-9050**PleurX™ Peritoneal Catheter Mini Kit**

- en** Peritoneal Catheter Mini Kit
For Peritoneal Placement Only
- fr** Mini Kit Cathéter Péritonéal
Pour un placement péritonéal uniquement
- de** Peritonealkatheter-Mini-Kit
Nur zur peritonealen Anlage
- it** Mini kit catetere peritoneale
Esclusivamente per
posizionamento peritoneale
- es** Mini Kit de Catéter Peritoneal
Únicamente para colocar en el peritoneo
- pt** Mini-kit de Cateter Peritoneal
Apenas para colocação peritoneal
- nl** Peritoneale katheter minikit
Uitsluitend voor peritoneale plaatsing
- sv** Peritoneal kateter, mini-kit
Endast för peritoneal placering
- fi** Peritoneaalikatetriminiseti
vain vatsakalvoon asettamiseen
- no** Peritonealt kateter - Minisett
Kun for peritoneal plassering

USA Rx Only© 2010, 2011, CareFusion Corporation or
one of its subsidiaries. All rights reserved.

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For Single Use Only

Do not use if
package is damaged.

For product inquiries or technical assistance:
pleurx-info@carefusion.com
www.carefusion.com/pleurx

Point-Lok is a trademark of
Smiths-Medical, ASD, Inc.
PleurX is a trademark or registered
trademark of CareFusion Corporation,
or one of its subsidiaries.

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










RC081195-2
McGaw Park, IL
Richard Cisneroz
12-07-11
12275967

Proofed by: _____ Date: _____
Dimensions checked: _____ Copy checked: _____

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en

Kit Contents**Placement Components**

- 1 PleurX Peritoneal Catheter, 15.5 Fr.
- 1 Needle, 18 G
- 1 Syringe, 12 ml
- 1 J-Tip Guidewire
- 1 Peel-Away Introducer, 16 Fr.
- 1 Tunneler
- 1 Point Lok Sharps Safety Device
- 1 CSR Wrap

Drainage Components



- 1 Drainage Line with Lockable Access Tip
- 1 Needle, 17 G x 1" (2.5 cm)
- 1 5-in-1 Drainage Line Adapter
- 1 Valve Cap

Dressing Components

- 8 Gauze Pads, 4" x 4" (10.2 cm x 10.2 cm)
- 1 Foam Catheter Pad

PleurX Peritoneal Catheter Description

The PleurX Peritoneal Catheter consists of a fenestrated silicone catheter with a valve mechanism and a polyester cuff. A barium sulfate stripe runs the entire length of the catheter. The valve is designed to prevent the passage of air or fluid in either direction unless it is accessed with the specifically matched drainage line or vacuum bottles provided by CareFusion.

PleurX Peritoneal Catheter**PleurX Lockable Drainage Line****Indications**

The PleurX Peritoneal Catheter Mini Kit (REF 50-9050) is indicated for:

- Intermittent drainage of symptomatic, recurrent, malignant ascites that does not respond to medical management of the underlying disease
- Palliation of symptoms related to recurrent malignant ascites
- For Peritoneal Placement Only

The PleurX Drainage Kits are indicated for use with either the PleurX Peritoneal Catheter Kit or the PleurX Pleural Catheter Kit.

Contraindications

Use of the PleurX Peritoneal Catheter Mini Kit and the PleurX Drainage Kit is contraindicated in the following situations:

1. When the peritoneal cavity is multi-loculated, and the drainage of a single loculation would not be expected to provide relief of dyspnea or other symptoms.
2. When there is a coagulopathy.
3. When the peritoneal cavity is infected.

Warnings

Do not put anything except the access tip of the lockable drainage line or PleurX Vacuum Bottles into the PleurX Catheter valve since this could damage the valve. A damaged valve may allow air into the body or let fluid leak out through the valve when not draining.

Cautions

For single use only. Re-use may result in a non-functional product or contribute to cross contamination.

Sterile technique should be used when placing and draining the catheter.

Individual patient anatomy, such as thin or weak abdominal wall, may require procedural variations to reduce the risk of leakage around the catheter.

Care must be taken when inserting the needle to avoid puncturing or lacerating the liver or bowel.

Damage to the guidewire may result if withdrawn through the needle.

Exercise care when placing the catheter to prevent it from coming into contact with surfaces such as drapes or towels. Silicone rubber is highly electrostatic and attracts airborne particles and surface contaminants.

Use rubber-shod instruments when handling the catheter. Possible cuts or tears can occur if rubber-shod instruments are not used.

Place a thumb over the end of the sheath as the dilator is removed to avoid peritoneal fluid leakage. Care must be taken not to bend or kink the sheath. Damage to the sheath may prevent passage of the catheter.

Do not use forceps on the introducer to break its handle and/or peel the sheath.

Exercise care when placing ligatures to avoid cutting or occluding the catheter.

In malignant ascites patients, paracentesis-related hypotension is uncommon, but has been documented. Use of IV fluid replacement and/or administration of colloidal agents can reduce the risk of hypotension. Additionally, initial drainage should be no more than 6L in the first 24 hours.

The pinch clamp must be fully closed to occlude the drainage line. When not connected to a suction source, make sure the pinch clamp is fully closed. Otherwise, the drainage line may allow air into the body or let fluid leak out.

When connecting to a vacuum bottle, make sure the pinch clamp on the drainage line is fully closed. Otherwise, it is possible for some, or all of the vacuum in the bottle to be lost.

When draining with glass vacuum bottles, do not use a needle larger than 17 G. If wall suction is used, it must be regulated to no greater than (-)60 cm H₂O, or to drain no more than 400 ml of fluid per minute.

Keep the valve on the PleurX Catheter and the lockable access tip on the drainage line clean. Keep them away from other objects to help avoid contamination.

Make sure that the valve and lockable access tip are securely connected when draining. If they are accidentally separated, they may become contaminated. If this occurs, clean the valve with an alcohol pad and use a new drainage set to avoid potential contamination.

Precautions should be taken to ensure the drainage line is not tugged or pulled.

It is normal for the patient to feel some discomfort or pain when draining fluid. If discomfort or pain is experienced when draining, clamp the drainage line to slow or stop the flow of fluid for a few minutes. Pain may be an indication of infection.

Potential complications of access and drainage of the peritoneal cavity include, but may not be limited to, the following: laceration of liver or bowel, hypotension/circulatory collapse, electrolyte imbalance, protein depletion, ascites leakage, peritonitis, wound infection, tumor growth in the catheter tunnel, and loculation of the peritoneal cavity.

Removal of chylous malignant ascites could exacerbate protein depletion or related nutritional complications.

Sterility


This product has been sterilized. It is for single use only and is not to be resterilized. Do not use if package is damaged. CareFusion will not be responsible for any product that is resterilized, nor accept for credit or exchange, any product that has been opened but not used.

General Guidelines

1. Systemic prophylactic antibiotics may be indicated.
2. The procedure for peritoneal placement can be performed using local anesthetic and sedation. However, depending on patient needs, it may be performed using alternative approaches to anesthesia or sedation.

3. The catheter should be placed under image guidance, using all precautions normally used for percutaneous placement of indwelling, tunneled catheters. Care should be taken to identify and avoid contact with vasculature near the guidewire insertion site.
4. Guidewire insertion site selection should be based upon patient anatomy and presentation with consideration given to any possible adhesions or loculated pockets of fluid. The fenestrated section of the catheter should preferentially be placed low in the peritoneal cavity to maximize access to fluid.
5. Consideration should be given to the patient's ease of access in determining the location of the catheter exit site.

Suggested Placement Procedure

Proper medical and surgical procedures are the responsibility of the physician. The appropriateness of any procedure must be based upon the needs of the patient. **Diagram (3)** illustrates the placement of the PleurX Peritoneal Catheter, as described in the following procedure.

Caution: Individual patient anatomy, such as thin or weak abdominal wall, may require procedural variations to reduce the risk of leakage around the catheter.

1. Position the patient appropriately to access the desired guidewire insertion site.
2. Identify the appropriate insertion site through which to place the guidewire. The guidewire insertion site is typically lateral to the midline, 6–10 cm below the costal margin, and above the patient's beltline. Ultrasound can be used to confirm the guidewire insertion site.
3. Identify the location of the catheter exit site, which is usually 5–8 cm superior and medial to the guidewire insertion site.
4. Surgically prep the patient.
5. Drape and anesthetize the planned insertion and tunneling sites.

Note: It is important to use an oblique approach while passing the introducer needle through the peritoneal wall to minimize the chance of ascitic fluid leakage.

Caution: Care must be taken when inserting the needle to avoid puncturing or lacerating the liver or bowel.

6. Insert the 18 G needle, attached to a syringe, through the peritoneum using an oblique angle.
7. Ensure free aspiration of ascitic fluid, then remove the syringe, leaving the needle in place.
8. Insert the guidewire through the needle and advance the distal end into the peritoneal cavity.
9. Remove the needle, leaving the guidewire in place.

Caution: Damage to the guidewire may result if withdrawn through the needle.

10. Make a 1 cm incision at the guidewire insertion site.
11. Make a second 1-2 cm incision 5-8 cm superior and medial to the guidewire insertion site. This incision will be the catheter exit site. Consider the patient's ease of access in determining its location. **See Diagram (3).**

Note: A smaller incision may provide better security of the catheter.

Note: Take care to ensure that the tunnel track has been anesthetized.

12. Attach the fenestrated end of the catheter onto the tunneler.

Caution: Exercise care when placing the catheter to prevent it from coming into contact with surfaces such as drapes or towels. Silicone rubber is highly electrostatic and attracts airborne particles and surface contaminants.

Caution: Use rubber-shod instruments when handling the catheter. Possible cuts or tears can occur if rubber-shod instruments are not used.

13. Pass the tunneler **(A)** and catheter **(B)** subcutaneously from the second incision down to and out through the first incision at the guidewire insertion site. **See Diagram (4).** Continue to draw the catheter through the tunnel until the polyester cuff lies inside the tunnel, about 1 cm **(C)** from the second incision. **See Diagram (5).** Disconnect the tunneler from the catheter.

Note: If the cuff is advanced further into the tunnel, it can make later removal of the catheter difficult.

14. Thread the 16 Fr. peel-away introducer sheath over the guidewire into the peritoneal cavity.
15. Remove the guidewire and dilator as a unit, leaving the 16 Fr. peel-away introducer sheath in place.

Caution: Place a thumb over the end of the sheath as the dilator is removed to avoid peritoneal fluid leakage. Care must be taken not to bend or kink the sheath. Damage to the sheath may prevent passage of the catheter.

16. Insert the fenestrated end of the catheter into the sheath advancing it until all the fenestrations are within the peritoneal cavity. This can be verified under fluoroscopy as fenestrations are located along the barium sulfate stripe.

17. Peel away the sheath while ensuring the catheter remains in place. Adjust the catheter so that it lies flat in the tunnel without any kinks.

Caution: Do not use forceps on the introducer to break the handle and/or peel the sheath.

18. Close the incision at the guidewire insertion site.
19. Close the incision site around the catheter and suture the catheter to the skin taking care not to restrict the diameter of the catheter. This suture is intended to remain in place at least until there is tissue ingrowth around the cuff.

Caution: Exercise care when placing ligatures to avoid cutting or occluding the catheter.

Drainage Procedure

The drainage procedure can be performed using:

- a) PleurX Vacuum Bottle(s)
- b) PleurX Drainage Line with other vacuum bottle(s) or
- c) Wall Suction

If using PleurX Vacuum Bottle(s), refer to PleurX Drainage Kit Instructions for Use.

Caution: In malignant ascites patients, paracentesis-related hypotension is uncommon, but has been documented. Use of IV fluid replacement and/or administration of colloidal agents can reduce the risk of hypotension. Additionally, initial drainage should be no more than 6L in the first 24 hours.

1. Clamp the drainage line completely closed using the pinch clamp found on the tubing.

See Diagram (6)

Caution: The pinch clamp must be fully closed to occlude the drainage line. When not connected to a suction source, make sure the pinch clamp is fully closed. Otherwise the drainage line may allow air into the body or let fluid leak out.

Caution: When connecting to a vacuum bottle make sure the pinch clamp on the drainage line is fully closed, otherwise, it is possible for some, or all of the vacuum in the bottle to be lost.

2. If using wall suction, attach the 5-in-1 adapter to the Luer fitting on the drainage line. If using a vacuum bottle other than PleurX, attach a 17 G needle to the Luer fitting on the drainage line.

Caution: When draining with glass vacuum bottles, do not use a needle larger than 17 G. If wall suction is used, it must be regulated to no greater than (-)60 cm H₂O, or to drain no more than 400 ml of fluid per minute.

3. Connect the drainage line to the vacuum/suction source.
4. Hold the drainage line near the lockable access tip and remove the protective cover by twisting it and gently pulling. Take care to avoid contaminating the lockable access tip.

See Diagram (7)

Caution: Keep the valve on the PleurX Catheter and the lockable access tip on the drainage line clean. Keep them away from other objects to help avoid contamination.

5. Continue holding the catheter near the valve. Carefully insert the lockable access tip into the catheter valve and advance it completely into the valve. You will feel and hear a click when the lockable access tip and valve are securely connected. **See Diagram (8)**

Caution: Make sure that the valve and the lockable access tip are securely connected when draining. If they are accidentally separated, they may become contaminated. If this occurs, clean the valve with an alcohol pad and use a new drainage set to avoid potential contamination.

6. If desired, lock the access tip to the catheter valve by twisting the access tip until you feel and hear a second click. **See Diagram (9)**

Caution: Precautions should be taken to ensure the drainage line is not tugged or pulled.

7. Release the pinch clamp on the drainage line to begin drainage. You can reduce the flow rate by squeezing the clamp partially closed.

See Diagram (10)

Caution: It is normal for the patient to feel some discomfort or pain when draining fluid. If discomfort or pain is experienced when draining, clamp the drainage line to slow or stop the flow of fluid for a few minutes. Pain may be an indication of infection.

Caution: Potential complications of access and drainage of the peritoneal cavity include, but may not be limited to, the following: laceration of liver or bowel, hypotension/circulatory collapse, electrolyte imbalance, protein depletion, ascites leakage, peritonitis, wound infection, tumor growth in the catheter tunnel, and loculation of the peritoneal cavity.

Caution: Removal of chylous malignant ascites could exacerbate protein depletion or related nutritional complications.

8. If you need to change to a new vacuum bottle or suction source for any reason, squeeze the pinch clamp on the drainage line completely closed. Remove the drainage line from the vacuum/suction source and connect to a new vacuum bottle or suction source. Release the pinch clamp to resume draining.

9. When flow stops or the desired amount of fluid has been removed, squeeze the pinch clamp on the drainage line completely closed.

See Diagram (6)

10. If locked, twist the lockable access tip to unlock it from the catheter valve. **See Diagram (11)**

11. Ensure the drainage line has been unlocked. With the drainage line in one hand and the catheter valve in the other hand, pull the lockable access tip out of the valve in a firm, smooth motion.

See Diagram (12)

12. Clean the catheter valve with an alcohol pad. Do not try to push anything through the valve as damage to the valve may occur.

See Diagram (13)

Warning: Do not put anything except the access tip of the lockable drainage line or PleurX Vacuum Bottles into the PleurX Catheter valve since this could damage the valve. A damaged valve may allow air into the body or let fluid leak out through the valve when not draining.

13. Place the valve cap over the catheter valve and twist it clockwise until it clicks into the locked position. **See Diagram (14)**

14. Clean around the catheter site.

15. Place the foam catheter pad around the catheter.

16. Wind the catheter into loops and place it over the foam pad.

17. Cover the catheter with gauze pads and secure with a self-adhesive dressing.

18. Disconnect the drainage line from the glass vacuum bottle or suction source. Dispose of the used drainage line and/or used vacuum bottles in accordance with applicable local, state, and federal regulations. Used product may present a potential biohazard.

Subsequent Drainage Procedures

Subsequent drainage procedures are to be performed using the PleurX Drainage Line, PleurX Vacuum Bottle, or the PleurX Drainage Kits. Each drainage kit contains the necessary drainage line, vacuum bottle, and other necessary items to perform the drainage procedure.

It is vital that patients and/or caregivers are carefully instructed on how to use the kit to drain malignant ascites. The person(s) responsible for drainage must be able to demonstrate they are capable of performing the procedure.

If the patient/caregiver is not able or willing to perform the drainage, a medical professional should perform the drainage.

It is recommended that the patient is periodically contacted or seen by a clinician to evaluate treatment regimen, assess need for albumin supplementation and evaluate catheter's functional status.

Catheter Removal Procedure

It may be appropriate and/or necessary at a later date to remove the PleurX Peritoneal Catheter. Three successive attempts to drain fluid that result in less than 50 ml of fluid removed may indicate one of the following:

- the ascites has resolved
 - the catheter is loculated away from the fluid
 - the catheter is occluded
1. Place the patient appropriately to access the catheter insertion site.
 2. Aseptically clean the patient's abdomen around the catheter insertion site.
 3. Anesthetize the site.
 4. Remove any remaining sutures securing the catheter.
 5. Using forceps, dissect around the cuff to free it from the ingrowth. Ensure that the cuff is completely free within the tunnel.
 6. Grasp the catheter in one hand and pull with a firm, constant pressure.
 7. Cover the site as appropriate.

Contains Phthalates. The benefit of treatment outweighs the remote possibility of exposure to phthalates.

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