



Alaris™ GW 800 Volumetric Pump

Models: 800TIG2GBD1, 800TIG2RWN1


Quick Reference Guide – en








CE
2797





























- This guide is not intended to be comprehensive instructions for the set-up and operation of the Alaris™ GW 800 Volumetric Pump.
- For complete pump information, refer to the Alaris GW 800 Volumetric Pump Directions For Use.
- All persons operating the pump should consult the Directions For Use.




Symbol	Description
	ON/OFF button - Press once to switch the Pump on. Press and hold down for approximately three seconds to switch the Pump off.
	RUN/HOLD button - Press to start the infusion or to put the infusion on hold.
	CLEAR/MUTE button - Press to silence alarm or advisories for one minute. The audio from the alarm or advisories will resound after this time. Resets numeric values to zero. Note: To re-enable the alarm audio press any other button or key to unmute.
	PRIME/BOLUS button - Primes the Infusion set. Administers bolus during the infusion.
	PRIMARY/SECONDARY button - Switches the Pump between Primary and Secondary infusion modes, if enabled.
	ENTER button - Scrolls between rate, time, VTBI and total volume infused (VI). Enters values for selected infusion / configuration parameters. Confirms the rate during an infusion titration.
	CHEVRON keys - Increases or decreases the infusion rate, TIME limit and VTBI. Press and hold to increase the selection speed. Used to adjust user selectable options

Symbol	Description
	AC POWER indicator - When illuminated the Pump is connected to an AC power supply.
	ALARM Indicator - When flashing indicates the Pump is in an alarm condition.
	MUTE Indicator - When flashing indicates the Pump audio alarm is muted.
RATE	When illuminated the Pump is displaying the infusion rate in millilitres per hour (ml/h).
VTBI	When illuminated the Pump is displaying the Volume To Be Infused (VTBI) in millilitres (ml).
VI	When illuminated the Pump is displaying the Volume Infused (VI) in millilitres (ml).
TIME	When illuminated the Pump is displaying the infusion time in hours : mins.
MICRO	When illuminated the Pump is operating in the Micro Mode. When not illuminated the Pump is in the Standard Mode.
SEC	When illuminated the Pump is operating in the Secondary Mode. When not illuminated the Pump is in the Primary Mode.
ml/h	(Millilitres / hour) When ml is illuminated the Pump displays the rate, VTBI or VI. When the h is illuminated the Pump displays the rate or infusion time.


	<ol style="list-style-type: none"> 1. Prepare the fluid container. Squeeze drip chamber until fluid reaches fill line (if applicable). Close in-line clamp.
	<ol style="list-style-type: none"> 2. Press door latch to open door.
	<ol style="list-style-type: none"> 3. Release flow stop mechanism by pushing the lever arm up and to the right.
	<ol style="list-style-type: none"> 4. Insert the infusion set from left to right into the pumping channel. Avoid any slack and avoid stretching or pulling the tubing.
	<ol style="list-style-type: none"> 5. Engage the flow stop mechanism by pushing the lever left and down.
	<ol style="list-style-type: none"> 6. Close door. Press recesses firmly, as shown, to ensure latch <i>clicks</i> into place.

	<ol style="list-style-type: none"> 1. Switch the pump on.
	<ol style="list-style-type: none"> 2. Open the in-line clamp. Ensure no drops are falling in drip chamber. Prime infusion set by pressing . While FILL is displayed, press  again and hold until all visible air is cleared from the line.
	<ol style="list-style-type: none"> 3. Connect the infusion set to the patient access device. Enter rate using  / . Press  once to confirm rate and scroll to VTBI.
	<ol style="list-style-type: none"> 4. Enter VTBI using  / , or, if using a flow sensor, VTBI can be switched off by scrolling until OFF is displayed. Press  once to confirm VTBI and scroll to VI.
	<ol style="list-style-type: none"> 5. If necessary, press  to clear VI.
	<ol style="list-style-type: none"> 6. Press  to start the infusion.

User Selectable Options

The pump must be on HOLD. Press .

Check the Battery Status

Note: Only available when Pump is operating on battery power.
Press and hold  for 2 seconds until Battery Status is displayed:

BAT=

- Battery has greater than 30 minutes run time.

BAT_




- Battery is low and has approximately 30 minutes or less run time.

Press  button to go to next option.

Selecting the Occlusion Pressure Level (if Enabled)


PrES option is displayed:






- Change occlusion pressure level to high (**HI**), normal (**nor**) or low (**lo**) using  / .
- Press  to confirm and go to next option.

Selecting the Alarm Volume (if Enabled)

tonE option is displayed:







- Change alarm volume between 1 and 7* using  / .
- Press  to confirm and go to next option.

* 1=Low, 7=High

Selecting VTBI / Time Infusion (if Enabled)

CLoc option is displayed:



- Select **On** or **OFF** using  /  keys to turn the VTBI/time infusion** setting on or off.
- Press  to confirm and go to next option.

** Ensure **TIME** indicator is lit, if set to **On**.

Selecting Micro Mode (if Enabled)

0.0 option is displayed:



- Select **On** or **OFF** using  /  keys to turn micro mode*** on or off.
- Press  to confirm.

*** Ensure **MICRO** indicator is lit, if set to **On**.

Panel Lock

Panel lock prevents certain key operation:



Rate Titration

To change the rate while infusing (if titration enabled):

1. Press / to adjust rate.
2. Press to confirm change of rate.


Administer BOLUS: Press

To administer a bolus while infusing:

- Press . While **bol** is displayed, press again and hold to administer the desired bolus volume.

An alarm indicates that the infusion has stopped.

Check the display for an alarm symbol.

Press  to silence the sound. (Exceptions are *Err* and *bAt* alarms, which cannot be silenced)

When the cause has been rectified, press  key to resume infusion.



All alarms will stop the infusion except for *Lo bAt* and *Atten*, where the infusion status will remain the same state as prior to the alarm.





Alarm / Condition	Action
<i>Air Occl</i> Air-in-line Upstream occlusion	<ul style="list-style-type: none"> • Disconnect the infusion set from the IV access device and use the prime function. • Remove occlusion/air and restart the infusion.
<i>bAt</i> Internal battery depleted.	<ul style="list-style-type: none"> • Connect pump to AC power. • Ensure pump is switched on to charge internal battery.
<i>door</i> Door open	<ul style="list-style-type: none"> • Close door and restart infusion.
<i>Err</i> System fault	<ul style="list-style-type: none"> • Switch pump off. • Remove pump from service and have the pump inspected by a qualified service engineer.
<i>FLo Err</i> Flow error Flow sensor attached to secondary infusion set.	<ul style="list-style-type: none"> • Close in-line clamp. • Ensure tubing is loaded correctly. • Ensure ample fluid is in container. • Clear blockage. Close pump door. • Restart infusion. • Ensure flow sensor is attached to primary infusion set.
<i>FLo SEnS</i> Flow sensor connection error.	<ul style="list-style-type: none"> • Connect flow sensor to pump and restart infusion. Ensure fluid in drip chamber is not above fill line.
<i>HI PrES</i> Downstream occlusion	<ul style="list-style-type: none"> • Remove pressure in the infusion set to prevent a post occlusion bolus to the patient. • Remove the cause of the blockage. • Restart the infusion.
<i>bad Set</i> Incorrect infusion set, infusion set incorrectly loaded or infusion set worn.	<ul style="list-style-type: none"> • Use infusion set as specified in DFU. Ensure infusion set is loaded correctly into pump.


Alarm / Condition	Action
End Finished pre-set volume to be infused with KVO switched off.	<ul style="list-style-type: none">• Set new VTBI.
Lo bAt Low battery. (At least 30 mins before <i>bAt</i> alarm)	<ul style="list-style-type: none">• Connect pump to AC power source.
Attn The pump has been left unattended for two minutes and infusion has not started.	<ul style="list-style-type: none">• Attend to pump.

Advisories alert the user and are indicated by an audible alarm, a message on the display or both.

Check the display for advisory message.

Press  to silence the sound.

Advisory / Condition	Action
<i>bol</i> Bolus is being administered.	<ul style="list-style-type: none"> Release  key to continue the infusion once prescribed bolus has been administered.
<i>End</i> VTBI limit reached.	<ul style="list-style-type: none"> Pump will infuse at the KVO rate, until  is pressed.
<i>FILL</i> The pump is priming the infusion set.	<ul style="list-style-type: none"> Ensure all air has been primed out of the infusion set, before starting the infusion.
<i>Hold</i> The pump is on hold.	<ul style="list-style-type: none"> Press  to restart / continue infusion or press  to return to setup.
<i>tEst</i> Automatic infusion set check.	<ul style="list-style-type: none"> Allow test to complete before operating the pump further.



BD Switzerland Sàrl, Route de Crassier 17,
Business Park Terre-Bonne, Batiment A4,
1262 Eysins, Switzerland

bd.com

BD, the BD Logo and Alaris are trademarks of Becton, Dickinson
and Company or its affiliates.
© 2019 BD. All rights reserved.

