
Revisions

BALTSO0191 Version 11.0 Template 4
Inserts

Rev from	Rev to	JOB #
02	03	101-17

Notes:

- BD Catalog Number: 220143
- Blank (Sheet) Size: Length: 8.5" Width: 33.0"
- Number of Pages: 12 Number of Sheets: 1
- Page Size: Length: 8.5" Width: 5.5" Final Folded Size: 5.5" x 8.5"
- Ink Colors: No. of Colors: 3 PMS#: 2755 Blue; 032 Red; Standard Black
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 4

- Vendor Printed ☒ Online/In House Printed ☐ Web ☐
- See specification control no. BALT220143 for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Nancy Carlsen at 2:44 pm, Dec 06, 2017	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY By Natalie Morio at 8:39 am, Dec 07, 2017		
Checked By	THIRD EYE BY By Nichole Graham at 11:05 am, Dec 07, 2017		
Part Number:	L008194	Category and Description Package Insert, BD ProbeTec™ ET CT/GC Amplified DNA Assay Collection Kit for Male Urethral Specimens	Sheet: 1 of 13 Scale: N/A A

BD ProbeTec™ ET *Chlamydia trachomatis*/Neisseria gonorrhoeae (CT/GC) Amplified DNA Assay Collection Kit for Male Urethral Specimens

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INTENDED USE

The **BD ProbeTec™ ET *Chlamydia trachomatis*/Neisseria gonorrhoeae (CT/GC) Amplified DNA Assay Collection Kit for Male Urethral Specimens** is used to collect and transport male patient urethral specimens to the laboratory for testing with the **BD ProbeTec ET *Chlamydia trachomatis* and Neisseria gonorrhoeae Amplified DNA Assays**.

SUMMARY AND EXPLANATION

The **BD ProbeTec ET *Chlamydia trachomatis*/Neisseria gonorrhoeae (CT/GC) Amplified DNA Assay Collection Kit for Male Urethral Specimens** consists of a sterile swab and a CT/GC Diluent tube for specimen transport. The rayon-tipped collection swab consists of a flexible aluminum wire in a plastic shaft. After the specimen is collected, the swab is placed into the CT/GC Diluent tube and the shaft is broken at the score mark. The tube is recapped and transported to the testing facility.

REAGENTS

Materials Provided: The **BD ProbeTec ET *Chlamydia trachomatis*/Neisseria gonorrhoeae (CT/GC) Amplified DNA Assay Collection Kit** contains 100 units. Each unit contains one sterile rayon-tipped collection swab and one CT/GC Diluent tube. The CT/GC Diluent tube contains 2 mL of diluent composed of potassium phosphate, DMSO, glycerol, Polysorbate 20 and 0.03% Proclin™ (preservative).

Materials Required But Not Provided: Gloves.

Storage Instructions: Store Collection Kit at 15–30 °C. Do not use beyond the expiration date. Swab sterility guaranteed if primary swab container is intact.

Precautions: For *in vitro* Diagnostic Use.

Warning

H315 Causes skin irritation. **H317** May cause an allergic skin reaction. **H319** Causes serious eye irritation. **H335** May cause respiratory irritation.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P272** Contaminated work clothing should not be allowed out of the workplace. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P321** Specific treatment (see on this label). **P304+P340** IF INHALED: Remove person to fresh air and keep comfortable for breathing. **P312** Call a POISON CENTER or doctor/physician if you feel unwell. **P363** Wash contaminated clothing before reuse. **P333+P313** If skin irritation or rash occurs: Get medical advice/attention. **P337+P313** If eye irritation persists: Get medical advice/attention. **P362+P364** Take off contaminated clothing and wash it before reuse. **P405** Store locked up. **P403+P233** Store in a well-ventilated place. Keep container tightly closed. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Warnings:

1. Do not pre-wet swab with CT/GC Diluent prior to collecting specimen.
2. Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁻⁴ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

- Optimal performance of the **BD ProbeTec ET** Amplified DNA Assays requires proper specimen collection, handling and transport.
- Use only collection and transport systems labeled for use with the **BD ProbeTec ET** *Chlamydia trachomatis* and *Neisseria gonorrhoeae* Amplified DNA Assays.
- Proper labeling should accompany each specimen to the laboratory.
- The CT/GC Diluent tube contains dimethyl sulfoxide (DMSO). DMSO is harmful by inhalation, in contact with skin and if swallowed. Keep container in a well-ventilated place. Avoid contact with eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water. Wear suitable protective clothing and gloves.
- When breaking the shaft of the swab, take care to avoid splashing of contents.
- Dispose of unused CT/GC Diluent tubes according to the Material Safety Data Sheet.
- ⓧ BD ProbeTec ET** *Chlamydia trachomatis/Neisseria gonorrhoeae* (CT/GC) Amplified DNA Assay Collection Kit for Male Urethral Specimens is for single use only; reuse may cause a risk of infection and/or inaccurate results.

SPECIMEN COLLECTION AND TRANSPORT

- Remove the swab from packaging.
- Insert the swab 2–4 cm into the urethra and rotate for 3–5 sec.
- Withdraw the swab.
- Uncap the CT/GC Diluent tube.
- Fully insert the swab into the CT/GC Diluent tube.
- Break the shaft of the swab at the score mark. Use care to avoid splashing of contents.
- Tightly** recap the tube.
- Label the tube with patient information and date/time collected.
- Transport to laboratory.

Swab Storage and Transport

The CT/GC Diluent tube with swab must be stored and transported to the laboratory and/or test site at 2–27 °C within 4–6 days of collection. Storage up to 4 days has been validated with clinical specimens; storage up to 6 days has been demonstrated with seeded swabs. In addition storage up to 30 days at 2–8 °C has been demonstrated with seeded specimens.

NOTE: If specimens cannot be transported directly to the testing laboratory under ambient temperatures (15–27 °C) and must be shipped, use an insulated container with ice and ship either overnight or with 2-day delivery.

For domestic and international shipments, specimens should be labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and etiologic agents/infectious substances. Time and temperature conditions for storage must be maintained during transport.

REFERENCES

- Clinical and Laboratory Standards Institute. Approved Guideline M29. Protection of laboratory workers from occupationally acquired infections. (Refer to latest addition)
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
- U.S. Department of Health and Human Services. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), U.S. Government Printing Office, Washington, D.C. (Refer to latest addition)
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.

Technical Information: In the United States contact BD Technical Service and Support at 800.638.8663 or www.bd.com.

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Rx Only This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

 Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Proclin is a trademark of Rohm and Hass Co.

© 2017 BD. BD, the BD Logo and BD ProbeTec are trademarks of Becton, Dickinson and Company.