
Rev from	Rev to	JOB #
01	02	7412-15

NOTES:

- BD Catalog Number: 220142
- Blank (Sheet) Size: Length: 8.5" Width: 33.0"
- Number of Pages: 12 Number of Sheets: 1
- Page Size: Length: 8.5" Width: 5.5" Final Folded Size: 5.5" x 8.5"
- Ink Colors: No. of Colors: 2 PMS#: 032 Red; Standard Black
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 4

- See specification control no. N/A for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Nancy Carlsen at 3:41 pm, Nov 09, 2015	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY By Helaine Jones at 4:20 pm, Nov 09, 2015		
Checked By	THIRD EYE BY By Terrence Means at 2:40 pm, Nov 11, 2015		
Part Number: L008195	Category and Description Package Insert, BD ProbeTec ET CT/GC Amplified DNA Assay Collection Kit for Endocervical Specimens	Sheet: 1 of 13	A
		Scale: N/A	

 BD ProbeTec™ ET *Chlamydia trachomatis*/Neisseria gonorrhoeae (CT/GC) Amplified DNA Assay Collection Kit for Endocervical Specimens

English: pages 1 – 2 Italiano: pagine 7 – 8
Dansk: side 2 – 4 Español: páginas 8 – 10
Français: pages 4 – 5 Svenska: sidan 10 – 11
Deutsch: seiten 5 – 7

L008195(02)
2015-11

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteys lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o representante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

The **BD ProbeTec™ ET *Chlamydia trachomatis*/Neisseria gonorrhoeae (CT/GC) Amplified DNA Assay Collection Kit** for Endocervical Specimens is used to collect and transport patient endocervical specimens to the laboratory for testing with the **BD ProbeTec ET *Chlamydia trachomatis* and Neisseria gonorrhoeae Amplified DNA Assays**.

SUMMARY AND EXPLANATION

The **BD ProbeTec ET *Chlamydia trachomatis*/Neisseria gonorrhoeae (CT/GC) Amplified DNA Assay Collection Kit** for Endocervical Specimens consists of a sterile cleaning swab, a sterile collection swab and a CT/GC Diluent tube for specimen transport. The large rayon-tipped cleaning swab is provided to clean the cervical os of blood and mucus prior to collection. The smaller polyurethane-tipped swab on a plastic shaft is used to collect the specimen for testing. After the specimen is collected, the swab is placed into the CT/GC Diluent tube and the shaft is broken at the score mark. The tube is recapped and transported to the testing facility.

REAGENTS

Materials Provided: The **BD ProbeTec ET *Chlamydia trachomatis*/Neisseria gonorrhoeae (CT/GC) Amplified DNA Assay Collection Kit** contains 100 units. Each unit contains one large sterile rayon-tipped cleaning swab, one sterile polyurethane-tipped collection swab and one CT/GC Diluent tube. The CT/GC Diluent tube contains 2 mL of diluent composed of potassium phosphate, DMSO, glycerol, Polysorbate 20 and 0.03% Proclin™ (preservative).

Materials Required But Not Provided: Gloves.

Storage Instructions: Store Collection Kit at 15–30 °C. Do not use beyond the expiration date. Swab sterility guaranteed if primary swab container is intact.

Precautions: For *in vitro* Diagnostic Use.

Warning

H315 Causes skin irritation. **H319** Causes serious eye irritation. **H317** May cause an allergic skin reaction. **H335** May cause respiratory irritation.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P304+P340** IF INHALED: Remove person to fresh air and keep comfortable for breathing. **P333+P313** If skin irritation or rash occurs: Get medical advice/attention. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Warnings:

1. Do not pre-wet swab with CT/GC Diluent prior to collecting specimen.
2. Pathogenic microorganisms including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁻⁴ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.
3. Optimal performance of the **BD ProbeTec ET Amplified DNA Assays** requires proper specimen collection, handling and transport.

4. Use only collection and transport systems labeled for use with the **BD ProbeTec ET *Chlamydia trachomatis* and *Neisseria gonorrhoeae* Amplified DNA Assays**.
5. The large swab is for cleaning the cervical os and should not be used for collecting the patient specimen.
6. Proper labeling should accompany each specimen to the laboratory.
7. The CT/GC Diluent tube contains dimethyl sulfoxide (DMSO). DMSO is harmful by inhalation, in contact with skin and if swallowed. Keep container in a well-ventilated place. Avoid contact with eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water. Wear suitable protective clothing and gloves.
8. When breaking the shaft of the swab, take care to avoid splashing of contents.
9. Dispose of unused CT/GC Diluent tubes according to the Material Safety Data Sheet.
10. **ⓧ BD ProbeTec ET *Chlamydia trachomatis*/*Neisseria gonorrhoeae* (CT/GC) Amplified DNA Assay Collection Kit** for Endocervical Specimens is for single use only; reuse may cause a risk of infection and/or inaccurate results.

SPECIMEN COLLECTION AND TRANSPORT

1. Remove the cleaning swab from packaging.
2. Using cleaning swab, remove excess blood and mucus from the cervical os.
3. Discard the used cleaning swab.
4. Remove the collection swab from packaging.
5. Insert the collection swab into the cervical canal and rotate for 15–30 sec.
6. Withdraw the swab carefully. Avoid contact with the vaginal mucosa.
7. Uncap the CT/GC diluent tube.
8. Fully insert the collection swab into the CT/GC Diluent tube.
9. Break the shaft of the swab at the score mark. Use care to avoid splashing of contents.
10. **Tightly** recap the tube.
11. Label the tube with patient information and date/time collected.
12. Transport to laboratory.

Swab Storage and Transport

The CT/GC Diluent tube with swab must be stored and transported to the laboratory and/or test site at 2–27 °C within 4–6 days of collection. Storage up to 4 days has been validated with clinical specimens; storage up to 6 days has been demonstrated with seeded swabs. In addition, storage up to 30 days at 2–8 °C has been demonstrated with seeded specimens.

NOTE: If specimens cannot be transported directly to the testing laboratory under ambient temperatures (15–27 °C) and must be shipped, use an insulated container with ice and ship either overnight or with 2-day delivery.

For domestic and international shipments, specimens should be labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and etiologic agents/infectious substances. Time and temperature conditions for storage must be maintained during transport.

REFERENCES

1. Clinical and Laboratory Standards Institute. Approved Guideline M29. Protection of laboratory workers from occupationally acquired infections. (Refer to latest addition)
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53–80.
3. U.S. Department of Health and Human Services. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C. (Refer to latest addition)
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021–0045.

Technical Information: In the United States, contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

BD ProbeTec ET *Chlamydia trachomatis*/*Neisseria gonorrhoeae* (CT/GC) udvidet DNA-assay prøvetagningskit til endocervikale præparater

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Proclin is a trademark of Rohm and Hass Co.
BD, BD Logo and BD ProbeTec are trademarks of Becton, Dickinson and Company © 2015 BD
EC REP Becton Dickinson France S.A.S