
Revisions

BALTSO0191 Version 8.0 Template 4
Inserts

Rev from	Rev to	JOB #
01	02	7100-14

NOTES:

1. BD Catalog Number: 231743, 231736
2. Blank (Sheet) Size: Length: 9" Width: 20"
3. Number of Pages: 10 Number of Sheets: 1
4. Page Size: Length: 9" Width: 4" Final Folded Size: 2.25" x 4"
5. Ink Colors: No. of Colors: 1 PMS#: Standard Black
6. Printed two sides: Yes ☒ No ☐
7. Style (see illustrations below): # 4

8. See specification control no. VS-L0001520JAA-PR for material information.
9. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Vendor Spec Held by BD Caribe, LTD.

Label Design	REVISED BY By Nancy Carlson at 4:14 pm, May 06, 2015	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY By Helaine Jones at 11:36 am, May 11, 2015		
Checked By	THIRD EYE BY By Sarah Cork at 4:03 pm, May 26, 2015		
Part Number:	L0001520JAA	Category and Description Package Insert, BD BBL Taxo INH	Sheet: 1 of 11 Scale: N/A A

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INTENDED USE

BD BBL™ Taxo™ INH Test Strips are used for detecting isonicotinic acid and its metabolites in urine.

BD BBL Taxo INH Test Control is used as a positive control for INH Test Strips.

SUMMARY AND EXPLANATION

INH test strips were developed by Kilburn et al.¹ to monitor patients on isoniazid therapy on an outpatient basis. The test aids physicians in determining whether their patients are taking the prescribed medication. The reagent-impregnated paper strip detects the presence of isonicotinic acid in the patient's urine. The strip method is a modification of Kilburn and Kubica's strip procedure for detecting niacin produced by *M. tuberculosis*.²

PRINCIPLES OF THE PROCEDURE

The pyridine ring of isonicotinic acid is split by cyanogen chloride to form a glutacetaldehyde derivative. This derivative condenses with barbituric acid to form a blue-purple polymethine dye.

REAGENTS

BD BBL Taxo INH Test Strips are absorbent paper strips impregnated with chloramine T, potassium thiocyanate, citric acid and barbituric acid as described by Kilburn et al. with some modifications.

BD BBL Taxo INH Test Control is an INH-impregnated disc that will yield a positive result in the test procedure.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Storage Instructions: Store **BD BBL Taxo INH Test Strips** at 2 – 8 °C. Store **BD BBL Taxo INH Test Control** at 15 – 30 °C.

After opening, store both products at the indicated temperature in the dark.

The expiration date applies to the product in its intact container when stored as directed.

Product Deterioration: Do not use a product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND HANDLING

Collect freshly voided urine 10 – 12 h after drug isoniazid administration.

PROCEDURE

Materials Provided: **BD BBL Taxo INH Test Strips** and/or **BD BBL Taxo INH Test Control**.

Materials Required But Not Provided: Test tubes, distilled or deionized water and scissors.

Preparation of Reagents:

1. Obtain one sealed-tube **BD BBL Taxo INH Test Strip** from the jar. Cut off one corner of the plastic tube at the arrow end of the strip.
2. Prepare a positive control for the test by placing one **BD BBL Taxo INH Test Control** disc in 2 mL distilled or deionized water in a test tube. Shake 2 – 3 times over a 15 min period to assure extraction of the INH into the water.

Test Procedure:

1. Squeeze 1/2 inch of the plastic tube of a **BD BBL Taxo INH Test Strip** at the end opposite the arrow between thumb and forefinger. While squeezing the tube, insert open end of the tube below the surface of the urine specimen to be tested. Release pressure. Sufficient specimen should rise in the tube to cover the arrow on the strip. The tube may be left to float in the urine container or may be transferred to a test tube.
2. Perform the test using the solution prepared from **BD BBL Taxo INH Test Control** in parallel with the urine specimen.
3. Incubate at room temperature for 15 – 30 min.
4. Observe and record results.

User Quality Control:

INH Test Strips

Identity Specifications

Appearance: White paper strips encased in a plastic tube. Arrows are imprinted on both sides.

Performance Response – Positive results (green, blue or purple color) were obtained with INH solutions in concentrations of 100 µg/mL to 5 µg/mL.

INH Test Control

Identity Specifications

Appearance: White, 1/4" paper disks with "I" over "C".

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

The appearance of a blue, purple or green color on the strip and in the liquid in the plastic tube is a positive test and indicates the presence of INH or its metabolites in the test specimen. Only the strip may appear colored if the reaction is weak.

LIMITATIONS OF THE PROCEDURE

A large proportion of the INH administered may be acetylated rapidly to isonicotinic acid in some patients. Urine to be tested for the metabolites of INH must be collected within 10 – 12 h of drug administration to maximize the sensitivity of the test.

The potential for interfering substances has not been evaluated for this test.

PERFORMANCE CHARACTERISTICS

INH Test Control Disks

Prior to release, all lots of **BD BBL Taxo** INH Test Control Discs are tested for performance characteristics. Representative samples of the lot are placed into tubes containing 2 mL of processed water. Tubes are shaken to extract the INH. A **BD BBL Taxo** INH Test Strip is added to each tube per label directions and to one tube containing water only. A blue-purple or green color develops on the test strips in the tubes containing the **Taxo** INH Control Disc. No color change is observed on the INH test strip placed in water only.

INH TEST STRIPS

Prior to release, all lots of **BD BBL Taxo** INH Test Strips are tested for performance characteristics. Dilutions of isonicotinic acid are prepared in concentrations of 100 µg/mL, 40 µg/mL, 20 µg/mL, 10 µg/mL and 5 µg/mL. Representative samples of the lot are placed per label directions into tubes containing 1 mL of each of the above dilutions and into one tube containing only processed water (negative control). The tubes are gently shaken and incubated at 20 to 25 °C. After 15 to 30 min, **Taxo** INH Test Strips placed in the INH dilutions show a positive reaction (blue, green, or purple color on the strip and/or in the liquid in the tube). The negative control remains colorless.

AVAILABILITY

Cat. No.	Description
231743	BD BBL™ Taxo™ INH Test Strips
231736	BD BBL™ Taxo™ INH Test Control

REFERENCES

1. Kilburn, J. O., R.E. Beam, H. L. David, E. Sanchez, R. F. Corpe and W. Dunn. 1972. Reagent-impregnated paper strip for detection of metabolic products of isoniazid in urine. *Am. Rev. Resp. Dis.* 106:923.
2. Kilburn, J. O. and G. P. Kubica. 1968. Reagent-impregnated paper strips for detection of niacin. *Amer. J. Clin. Path.* 60:530 – 532.

Technical Information: In the United States, contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Benex Limited
Becton Dickinson Pty Ltd.