# Veritor™ System

## For Rapid Detection of Flu A+B

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For *in vitro* diagnostic use only.

### INTENDED USE

The BD Veritor™ System for Rapid Detection of Flu A+B is a rapid chromatographic immunoassay for the direct and qualitative detection of influenza A and B viral nucleoprotein antigens from nasopharyngeal wash, aspirate and swab in transport media samples from symptomatic patients. The BD Veritor System for Rapid Detection of Flu A+B is a differentiated test, such that influenza A viral antigens can be distinguished from influenza B viral antigens from a single processed sample using a single device. The test is to be used as an aid in the diagnosis of influenza A and B viral infections. A negative test is presumptive and it is recommended that these results be confirmed by viral culture or an FDA-cleared influenza A and B molecular assay. Outside the U.S. a negative test is presumptive and it is recommended that these results be confirmed by viral culture or a molecular assay cleared for diagnostic use in the country of use. FDA has not cleared this device outside the U.S. Negative test results do not preclude influenza viral infection and should not be used as the sole basis for treatment or other patient management decisions. The test is not intended to detect influenza C antigens.

Performance characteristics for influenza A and B nasopharyngeal (NP) washes/aspirates were established during January through March of 2011 when influenza viruses A/2009 H1N1, A/H3N2, B/Victoria lineage, and B/Yamagata lineage were the predominant influenza viruses in circulation according to the *Morbidity and Mortality Weekly Report* from the CDC entitled "Update: Influenza Activity—United States, 2010–2011 Season, and Composition of the 2011–2012 Influenza Vaccine." Performance characteristics may vary against other emerging influenza viruses.

Performance characteristics for influenza A and B NP swabs in transport media were established during January through April of 2012 when influenza viruses A/2009 H1N1, A/H3N2, B/Victoria lineage, and B/Yamagata lineage were the predominant influenza viruses in circulation according to the *Morbidity and Mortality Weekly Report* from the CDC entitled "Update: Influenza Activity—United States, 2011–2012 Season, and Composition of the 2012–2013 Influenza Vaccine." Performance characteristics may vary against other emerging influenza viruses.

If infection with a novel influenza virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to the state or local health department for testing. Virus culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.

### SUMMARY AND EXPLANATION

Influenza illness classically presents with sudden onset of fever, chills, headache, myalgias, and a non-productive cough. Epidemics of influenza typically occur during winter months with estimated 114,000 hospitalizations<sup>1</sup> and 36,000 deaths<sup>2</sup> per year in the U.S. Influenza viruses can also cause pandemics, during which rates of illness and death from influenza-related complications can increase dramatically.

Patients who present with suspected influenza may benefit from treatment with an antiviral agent especially if given within the first 48 hours of onset of illness. It is important to rapidly distinguish influenza A from influenza B in order to allow physicians a choice in selective antiviral intervention. Moreover, it is important to determine if influenza A or B is causing symptomatic disease in a particular institution (e.g., nursing home) or community, so that appropriate preventative intervention can be taken for susceptible individuals. It is therefore important to not only rapidly determine whether influenza is present, but also which type of influenza virus is present as severity and treatment can be different.<sup>3</sup>

Diagnostic tests available for influenza include rapid immunoassay, immunofluorescence assay, polymerase chain reaction (PCR), serology, and viral culture.<sup>4-11</sup> Immunofluorescence assays entail staining of specimens immobilized on microscope slides using fluorescent-labeled antibodies for observation by fluorescence microscopy.<sup>6,12,13</sup> Culture methods employ initial viral isolation in cell culture, followed by hemadsorption inhibition, immunofluorescence, or neutralization assays to confirm the presence of the influenza virus.<sup>13-15</sup>

The BD Veritor System for Rapid Detection of Flu A+B (also referred to as the BD Veritor System and BD Veritor System Flu A+B) is a digital immunoassay (DIA) to qualitatively detect influenza A or B nucleoprotein antigens from respiratory specimens of symptomatic patients with a time to result of 10 minutes. The speed and simplified workflow of the BD Veritor System for Rapid Detection of Flu A+B makes it applicable as a "STAT" influenza A and B antigen detection test providing relevant information to assist with the diagnosis of influenza. All BD Veritor System Flu A+B test devices are interpreted by a BD Veritor System Instrument, either a BD Veritor Reader or BD Veritor Plus Analyzer (the "Analyzer"). When using an Analyzer, procedures to evaluate test devices depend on the workflow configuration chosen. In **Analyze Now mode**, the instrument evaluates assay devices after manual timing of their development. In **Walk Away mode**, devices are inserted immediately after application of the specimen, and timing of assay development and analysis is automated. Connection of the Analyzer to a printer is possible if desired. Additional result documentation capabilities are possible with the implementation of the BD Synapsys™ Informatics Solution, and with the addition of the BD Veritor InfoScan module and BD Veritor Plus Connect. Please refer to the Analyzer *Instructions for Use* for details on these features and contact BD technical support for more information.

## PRINCIPLES OF THE PROCEDURE

The BD Veritor System for Rapid Detection of Flu A+B is a qualitative digital immunoassay for the detection of influenza A and B viral antigens in samples processed from respiratory specimens. When specimens are processed and added to the test device, influenza A or B viral antigens bind to anti-influenza antibodies conjugated to detector particles in the A+B test strip. The antigen-conjugate complex migrates across the test strip to the reaction area and is captured by the line of antibody on the membrane. A positive result for influenza A is determined by the BD Veritor System Instrument when antigen-conjugate is deposited at the Test "A" position and the Control "C" position on the BD Veritor System Flu A+B assay device. A positive result for influenza B is determined by the BD Veritor System Instrument when antigen-conjugate is deposited at the Test "B" position and the Control "C" position in the BD Veritor System Flu A+B assay device. The instrument analyzes and corrects for non-specific binding and detects positives not recognized by the unaided eye to provide an objective digital result.

## REAGENTS

The following components are included in the BD Veritor System for Rapid Detection of Flu A+B kit:

|                                   |                              |                                                                                                                                                                                                                  |
|-----------------------------------|------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| BD Veritor System Flu A+B Devices | 30 devices                   | Foil pouched device containing one reactive strip. Each strip has two test lines of monoclonal antibody specific to either Flu A or Flu B influenza viral antigen and murine monoclonal control line antibodies. |
| RV Reagent C                      | 30 tubes with 100 µL reagent | Detergent with <0.1% sodium azide                                                                                                                                                                                |
| 300 µL Pipette                    | 30 each                      | Transfer pipette                                                                                                                                                                                                 |
| Control A+/B- Swab                | 1 each                       | Flu A Positive and Flu B Negative Control Swab, influenza A antigen (inactive recombinant nucleoprotein) with <0.1% sodium azide                                                                                 |
| Control B+/A- Swab                | 1 each                       | Flu A Negative and Flu B Positive Control Swab, influenza B antigen (inactive recombinant nucleoprotein) with <0.1% sodium azide                                                                                 |

**Materials Required But Not Provided:** BD Veritor™ Plus Analyzer (Cat. No. 256066), Timer, Tube Rack for specimen testing.

**Optional Equipment:** BD Veritor™ InfoScan Module (Cat. No. 256068), USB Printer Cable for BD Veritor™ Analyzer (Cat. No. 443907), Epson Printer model TM-T20 II, BD Veritor Plus Connect (contact BD Technical Services for details).

## Warnings and Precautions:

### Warning

**H302** Harmful if swallowed. **H402** Harmful to aquatic life. **H412** Harmful to aquatic life with long lasting effects.

**P273** Avoid release to the environment. **P264** Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P330** Rinse mouth. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

1. For *in vitro* Diagnostic Use.
2. Test results are not meant to be visually determined. **All test results must be determined using the BD Veritor System Instrument.**
3. If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to the state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.
4. Pathogenic microorganisms, including hepatitis viruses, Human Immunodeficiency Virus and novel influenza viruses, may be present in clinical specimens. "Standard Precautions"<sup>16-19</sup> and institutional guidelines should be followed in handling, storing and disposing of all specimens and all items contaminated with blood and other body fluids.
5. Dispose of used BD Veritor System test devices as biohazardous waste in accordance with federal, state and local requirements.
6. Reagents contain sodium azide, which is harmful if inhaled, swallowed or exposed to skin. Contact with acids produces very toxic gas. If there is contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.
7. Do not use kit components beyond the expiration date.
8. Do not reuse the BD Veritor System test device.
9. Do not use the kit if the Control A+/B- swab and Control B+/A- swab do not yield appropriate results.
10. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
11. To avoid erroneous results, specimens must be processed as indicated in the assay procedure section.

12. FluMist® is made from attenuated live flu virus and although the concentration tested (1%) was non-interfering, it is possible when tested with higher concentrations that an influenza A and/or influenza B false positive may occur.
13. Specific training or guidance is recommended if operators are not experienced with specimen collection and handling procedures.

**Storage and Handling:** Kits may be stored at 2–30 °C. DO NOT FREEZE. **Reagents and devices must be at room temperature (15–30 °C) when used for testing.**

## **SPECIMEN COLLECTION AND HANDLING**

**Specimen Collection and Preparation:** Acceptable specimens for testing with the BD Veritor System for Rapid Detection of Flu A+B include nasopharyngeal (NP) washes, aspirates and swab specimens in transport media. It is essential that correct specimen collection and preparation methods be followed. Specimens obtained early in the course of the illness will contain the highest viral titers.

Inadequate specimen collection, improper specimen handling and/or transport may yield a false negative result; therefore, training in specimen collection is highly recommended due to the importance of specimen quality to accurate test results.

**Specimen Transport Media:** The following transport media have been tested and found to be compatible using moderate positive samples with the BD Veritor System for Rapid Detection of Flu A+B:

- Modified Amies Medium (liquid), ESwab, Liquid Stuart Medium, Amies, Bartel ViraTrans™, BD Universal Viral Transport, Hank's Balanced Salt solution, M4, M4-RT, M5, M6, Normal Saline, Phosphate Buffered Saline.

Samples in these transport media can be stored at 2–8 °C for up to 72 hours. Allow specimens that have been stored at 2–8 °C to warm to room temperature prior to use. Failure to do so may result in inconsistent flow.

Other transport media may be utilized if an appropriate validation exercise is performed.

### **Specimen Transport and Storage:**

Freshly collected specimens should be processed within 1 hour. If necessary, specimens may be stored at 2–8 °C for up to 72 hours and then tested at room temperature. It is essential that correct specimen collection and preparation methods be followed. Do not centrifuge specimens prior to use, as the removal of cellular material may adversely affect test sensitivity.

### **Procedure for Nasopharyngeal Washes/Aspirates:**

- For NP washes/aspirates, sample volumes of 1 to 3 mL are recommended. If transport medium is used, minimal dilution of specimens is recommended.
- Excessive wash volumes should be avoided as they may result in decreased test sensitivity.
- Process specimen as described in "Test Procedure."

### **Procedure for Nasopharyngeal Swabs in Transport Media:**

- For NP swabs in transport media, a minimal volume of transport media (1 mL) is recommended to minimize dilution.
- Process specimen as described in "Test Procedure."

TEST PROCEDURE


**NOTES:** Allow specimens that have been stored at 2–8 °C to warm to room temperature prior to use. Thoroughly mix all specimens prior to removal of an aliquot for processing. Do not centrifuge specimens.

Prepare for testing

The following steps assume that users of a BD Veritor Plus Analyzer have chosen and set all configuration options, and that the Analyzer is ready to use. To choose or change these settings, see the BD Veritor Plus Analyzer *Instructions for Use*, section 4.7. A printer is not necessary to display results. However, if your facility has chosen to connect the BD Veritor Plus Analyzer to a printer, check that the printer is plugged into a power source, paper supply is adequate and any necessary network connections are enabled before testing.

For each patient specimen or control swab:


- Step 1:** Remove one **RV Reagent C** tube/tip and one BD Veritor System Flu A+B device from its foil pouch immediately before testing.
- Step 2:** Label one BD Veritor System device and one **RV Reagent C** tube for each specimen and control to be tested.
- Step 3:** Place the labeled **RV Reagent C** tube(s) in the designated area of the tube rack.

All reagents and specimens must be at room temperature prior to processing


- Step 4:** Process the Sample or Control:
- a. For NP washes, aspirates and swab specimens in transport media:**

1. Vortex or thoroughly mix specimen. Do not centrifuge.
2. Remove and discard the cap from the **RV Reagent C** tube corresponding to the sample to be tested.
3. Using the transfer pipette, transfer 300 µL of the specimen into the **RV Reagent C** tube. Discard pipette after use.


**b. For kit swab controls:**

1. Remove and discard the cap from the **RV Reagent C** tube corresponding to the sample to be tested.
2. Using the transfer pipette add 300 µL of distilled or deionized water to the **RV Reagent C** tube.
3. Insert the control swab into the tube and vigorously plunge the swab up and down in the fluid for a minimum of 15 seconds.
4. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.

- Step 5:**
- a. Press the attached tip firmly onto the **RV Reagent C** tube containing the processed specimen or control (threading/twisting not required).
  - b. Vortex or mix thoroughly by swirling or flicking the bottom of the tube.


**NOTE:** Do not use tips from any other product, including other products from BD or other manufacturers.

| After step 5, choose from the model and workflow option below before continuing to step 6: |                                                                                                                                                                                  |                                                                                              |                                                                                                        |                                                                                              |
|--------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|
|                                                                                            | BD Veritor Reader or Analyzer in Analyze Now mode                                                                                                                                | BD Veritor Plus Analyzer in Walk Away mode                                                   | BD Veritor Plus Analyzer with the BD Veritor InfoScan module in Analyze Now mode---or---Walk Away mode |                                                                                              |
| Instructions in section:                                                                   |  <b>A</b>  | <b>B</b>  | <b>C</b>            | <b>D</b>  |

**Step 6A: Adding the specimen**

- Invert the **RV Reagent C** tube and hold the tube vertically (approximately one inch above the labeled BD Veritor System Flu A+B device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well of a labeled BD Veritor System Flu A+B device.


**NOTE:** Squeezing the tube too close to the tip may cause leakage

**Step 7A: Timing development**

- After adding the sample, allow the test to run for **10 minutes** before inserting into the BD Veritor Instrument.
- **NOTE:** If running test under laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.

**Step 8A: Using the BD Veritor Instrument:**

- During incubation time, turn the BD Veritor Instrument on by pressing the power button once.
- Insert the assay device when the 10 minute assay development time is complete. Follow the on screen prompts to complete the procedure.
- The status of the assay analysis process appears in the display window.

**Do not touch the Instrument or remove the test device**

**Step 9A: Record the Result**


- When analysis is complete, the test result appears in the display window.

**ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).**

To use Walk Away mode – connect the AC power adapter to the Analyzer and a power source


## Step 6B: Starting Walk Away mode:

- Turn the Analyzer on by pressing the blue power button once.
- When the display window reads: "INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE"
  - **Double-click** the blue power button.

## Step 7B: Adding the specimen

- When the display window reads "ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY":
  - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System Flu A+B device sample well).
  - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System Flu A+B device.


**NOTE:** Squeezing the tube close to the tip may cause leakage

## Step 8B: Start the development and reading sequence

- Immediately insert the test device into the slot on the right side of the Analyzer.

**The test device must remain horizontal to prevent spilling the specimen out of the sample well.**

- "DO NOT DISTURB TEST IN PROGRESS" appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- A countdown timer in the display window shows the remaining analysis time.

**Do not touch the Analyzer or remove the test device during this time. Doing so will abort the assay analysis.**


## Step 9B: Record the result

- When analysis is complete, the test result appears in the display window.

**ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 60 minutes (when AC power adapter is connected).**

## Step 6C: Adding the specimen

- Invert the tube, holding it vertically (approximately one inch above the BD Veritor System Flu A+B device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well of a labeled BD Veritor System Flu A+B device. **NOTE: Squeezing the tube close to the tip may cause leakage.**

## Step 7C: Timing development

- Allow the test to develop for 10 minutes.
- If running the test in a laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.



## Step 8C: Using the Analyzer

During incubation time, turn the BD Veritor Plus Analyzer on by pressing the blue button once.


The display window briefly shows “SCAN CONFIG BARCODE.” This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration instructions. Ignore this message and postpone this process when an assay is awaiting analysis.

- When assay development time is complete and the Analyzer display window reads: “INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE”:
  - Insert the BD Veritor System Flu A+B device into the BD Veritor Plus Analyzer.

## Step 9C: Using the Bar Code scanner

- Follow the prompts on the display window to complete any required barcode scans of:
  - OPERATOR ID
  - SPECIMEN ID and/or
  - KIT LOT NUMBER
 in accordance with site requirements and Analyzer settings.

- Prompts for each scanning step appear in the display window for only 30 seconds. Failure to complete scans during that time will cause the Analyzer to default to the beginning of step 8C. To restart this step, remove and reinsert the test device to initiate a new sequence.
- Move the barcode slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.
- The Analyzer can record the Kit Lot Number in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user. BD recommends against the use of expired materials.

- After required scans are completed, the Analyzer displays a countdown timer and test analysis begins.
- Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.
- When analysis is complete a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed.

If a printer is not connected, record the result before removing the assay device.

**ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).**

## Step 10C: Remove the test device

- Pull the device out. The display will show “INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE” to indicate the Analyzer is ready to perform another test.

If the Veritor Plus Analyzer is connected to an LIS, a steady ENVELOPE symbol will appear to indicate that results are awaiting transmission. In the event that a network connection is not detected while the ENVELOPE symbol is still displayed, the Analyzer will queue all untransmitted results and attempt to transmit them when reconnected. If it is powered off during this time, it will attempt to transmit as soon as power is restored, and connection is re-established. A flashing envelope indicates that data are in the process of being transmitted.

To use Walk Away mode - connect the AC power adapter to the Analyzer and a power source


## Step 6D: Starting Walk Away mode

- Turn the Analyzer on by pressing the blue power button once.  
The display window will briefly show "SCAN CONFIG BARCODE." This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration instructions. Ignore this message and postpone this process when an assay is awaiting analysis.
- When the display window reads: "INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE"
  - Double-click the blue power button.

## Step 7D: Using the Bar Code scanner

- Follow the prompts on the display window to complete any required barcode scans of:
  - OPERATOR ID
  - SPECIMEN ID and/or
  - KIT LOT NUMBER


According to site requirements and Analyzer settings.

- Prompts for each scanning step appear in the display window for only 30 seconds. Failure to complete scans during that time will cause the Analyzer to default to the beginning of step 6D. To restart this step, double-click the power button.
- Move the barcode slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.
- The Analyzer can record the Kit Lot Number in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user. BD recommends against the use of expired materials.


## Step 8D: Add the specimen to the test device

- When the display window reads: "ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY":
  - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System Flu A+B device sample well).
  - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System Flu A+B device. **NOTE: Squeezing the tube close to the tip may cause leakage.**

## Step 9D: Start the development and reading sequence

- Immediately insert the test device into the slot on the right side of the Analyzer.
- The test device must remain horizontal to prevent spilling the specimen out of the sample well.**
- "DO NOT DISTURB TEST IN PROGRESS" appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
  - A countdown timer in the display window shows the remaining analysis time.
- Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.**
- When analysis is complete, a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed. **If a printer is not connected, record the result before removing the assay device.**

**ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 60 minutes (when AC power adapter is connected).**

## Step 10D: Remove the test device

- Pull the device out. The display will show INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE to indicate the Analyzer is ready to perform another test. Note that the Analyzer returns to Analyze Now mode at the conclusion of each read sequence.

If the Veritor Plus Analyzer is connected to an LIS, a steady ENVELOPE symbol will appear to indicate that results are awaiting transmission. In the event that a network connection is not detected while the ENVELOPE symbol is still displayed, the Analyzer will queue all untransmitted results and attempt to transmit them when reconnected. If it is powered off during this time, it will attempt to transmit as soon as power is restored, and connection is re-established. A flashing envelope indicates that data are in the process of being transmitted.

**OPTIONAL TEST PROCEDURE:** Use this procedure to test for both INFLUENZA A+B and RSV using a single NP wash, aspirate or swab specimen in transport medium.

**Note: The BD Veritor™ System for Rapid Detection of RSV (Cat. No. 256042) is required for this procedure in addition to the BD Veritor™ System for Rapid Detection of Flu A+B (Cat. No. 256041).**

**IMPORTANT NOTE:** This optional procedure allows for use of the remaining processed sample from Step 5 above to test additionally for RSV. The sample to be tested with the RSV kit must be from a patient less than 20 years of age as indicated in the BD Veritor RSV laboratory kit package insert. **THE PROCESSED SAMPLE SHOULD BE TESTED WITHIN 15 MINUTES.**

**Allow specimens that have been stored at 2–8 °C to warm to room temperature prior to use. Thoroughly mix all specimens prior to removal of an aliquot for processing. Do not centrifuge specimens.**

1. Collect specimen from the patient and follow Steps 1–5 of the test procedure above to prepare the sample for testing.
2. Using the sample from Step 5, continue the test procedure using the test device for RSV and the same workflow configuration used to obtain the Flu A+B result.
3. Refer to the product insert for BD Veritor System for Rapid Detection of RSV (Cat. No. 256042) for the test procedure and full description of the BD Veritor RSV test. Follow the Instructions in the insert and the Instrument on screen prompts to complete the test procedure and obtain results. Refer to the product insert for the BD Veritor System RSV Kit for result interpretation.

**INTERPRETATION OF RESULTS**

The BD Veritor System Instrument (purchased separately) must be used for interpretation of all test results. Operators should not attempt to interpret assay results directly from the test strip contained within the BD Veritor System Flu A+B assay device. With some specimens, up to four lines may be visible on the test device. The Instrument will appropriately interpret the result.

| Display                  | Interpretation                                          |
|--------------------------|---------------------------------------------------------|
| FLU A: +<br>FLU B: -     | Positive Test for Flu A (influenza A antigen present)   |
| FLU A: -<br>FLU B: +     | Positive Test for Flu B (influenza B antigen present)   |
| FLU A: -<br>FLU B: -     | Negative Test for Flu A and Flu B (no antigen detected) |
| RESULT INVALID           | Result Invalid                                          |
| POSITIVE CONTROL INVALID | Test Invalid. Repeat the test.                          |
| NEGATIVE CONTROL INVALID | Test Invalid. Repeat the test.                          |

**Invalid Test –** If the test is invalid, the BD Veritor System Instrument will display “RESULT INVALID” or “CONTROL INVALID” and the test or control must then be repeated. Because true dual positives are exceptionally rare, the BD Veritor System Instrument reports dual positive influenza A and influenza B results as “Result Invalid.” Specimens generating a “Result Invalid” should be retested. Upon retesting, if the specimen produces a “Result Invalid” the user may want to consider other methods to determine whether the sample is positive or negative for influenza virus.

**REPORTING OF RESULTS**

- Positive Test** Positive for the presence of influenza A or influenza B antigen. A positive result may occur in the absence of viable virus.
- Negative Test** Negative for the presence of influenza A and influenza B antigen. Infection due to influenza cannot be ruled-out because the antigen present in the sample may be below the detection limit of the test. A negative test is presumptive and it is recommended that these results be confirmed by viral culture or an FDA-cleared influenza A and B molecular assay.
- Invalid Test** Test result is inconclusive. Do not report results. Repeat the test.

## QUALITY CONTROL:

To utilize the Analyzer's QC documentation capability, specimen barcode scanning must be enabled on an Analyzer equipped with the BD Veritor InfoScan module. Please refer to the Analyzer *Instructions for Use*, section 4, to choose or change this configuration.

Each BD Veritor System Flu A+B device contains both positive and negative internal/procedural controls:

1. The internal positive control validates the immunological integrity of the device, proper reagent function, and assures correct test procedure.
2. The membrane area surrounding test lines functions as a background check on the assay device.

**The BD Veritor System Instrument evaluates the positive and negative internal/procedural controls after insertion of each BD Veritor System test device. The BD Veritor System Instrument prompts the operator if a quality issue occurs during assay analysis. Failure of the internal/procedural controls will generate an invalid test result. NOTE: The internal controls do not assess proper sample collection technique.**

### External Positive and Negative Controls:

Flu A positive/B negative and Flu B positive/A negative control swabs are supplied with each kit. These controls provide additional quality control material to assess that the test reagents and the BD Veritor System Instrument perform as expected. Prepare kit control swabs and test using the same procedure (either **Analyze Now** or **Walk Away** mode) as used for patient specimen swabs. When using the barcode scanning feature to document QC procedures, scan the barcode on the control swab packaging when prompted for a Specimen ID.

**Your laboratory's standard Quality Control procedures and applicable local, state and/or federal regulations or accreditation requirements dictate the performance of external quality control procedures.**

BD recommends external controls be run once for:

- each new kit lot,
- each new operator,
- each new shipment of test kits,
- as required by internal quality control procedures and in accordance with local, state and federal regulations or accreditation requirements.

## LIMITATIONS OF THE PROCEDURE

- Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.
- The contents of this kit are to be used for the qualitative detection of influenza type A and B antigens from NP wash, aspirate and swab in transport media specimens.
- The BD Veritor System for Rapid Detection of Flu A+B is capable of detecting both viable and non-viable influenza particles. The BD Veritor System for Rapid Detection of Flu A+B performance depends on antigen load and may not correlate with other diagnostic methods performed on the same specimen.
- Results from the BD Veritor System for Rapid Detection of Flu A+B test should be correlated with the clinical history, epidemiological data and other data available to the clinician evaluating the patient.
- A false-negative test result may occur if the level of viral antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly; therefore, a negative test result does not eliminate the possibility of influenza A or influenza B infection, and should be confirmed by viral culture or an FDA-cleared influenza A and B molecular assay.
- Positive test results do not rule out co-infections with other pathogens.
- Positive test results do not identify specific influenza A virus subtypes.
- Negative test results are not intended to rule out other non-influenza viral or bacterial infections.
- Children tend to shed virus for longer periods of time than adults, which may result in differences in sensitivity between adults and children.
- Positive and negative predictive values are highly dependent on prevalence rates. Positive test results are more likely to represent false positive results during periods of little/no influenza activity when disease prevalence is low. False negative test results are more likely during peak influenza activity when prevalence of disease is high.
- This device has been evaluated for use with human specimen material only.
- Monoclonal antibodies may fail to detect, or detect with less sensitivity, influenza A viruses that have undergone minor amino acid changes in the target epitope region.
- The analytical reactivity of this device has not been established for avian or swine origin influenza strains other than those included in the "Strain Reactivity" tables.
- The BD Veritor System Instrument reports dual positive influenza A and influenza B results as "Result Invalid." True dual positives are exceptionally rare. Specimens generating a "Result Invalid" should be retested. Upon retesting, if the specimen produces a "Result Invalid" the user may want to consider other methods to determine whether the sample is positive or negative for influenza virus.

EXPECTED VALUES

The rate of positivity observed in respiratory testing will vary depending on the method of specimen collection, handling/transport system employed, detection method utilized, the time of year, age of the patient, geographic location and most importantly, local disease prevalence.

The overall prevalence observed with an FDA-cleared influenza A and B molecular assay in the U.S. during the 2010–2011 clinical study was 23.9% for influenza A and 7.5% for influenza B. At the clinical site located in Hong Kong, the prevalence observed with the same FDA-cleared influenza A and B molecular assay was 7.2% for influenza A and 3.4% for influenza B.

The overall prevalence observed with an FDA-cleared influenza A and B molecular assay in the U.S. during the 2011–2012 clinical study was 31.7% for influenza A and 4.5% for influenza B. At the clinical sites located in Japan, the prevalence observed with the same FDA-cleared influenza A and B molecular assay was 0% for influenza A and 89% for influenza B.

PERFORMANCE CHARACTERISTICS

Explanation of Terms:

- PPA: Positive Percent Agreement =  $a / (a+c) \times 100\%$
- NPA: Negative Percent Agreement =  $d / (b+d) \times 100\%$
- P: Positive
- N: Negative
- C.I.: Confidence Interval

| New Test Method | Comparator Method |       |
|-----------------|-------------------|-------|
|                 | P                 | N     |
| P               | a                 | b     |
| N               | c                 | d     |
| Total           | (a+c)             | (b+d) |

Clinical Performance NP Washes/Aspirates 2010–2011:

Performance characteristics for the BD Veritor System for Rapid Detection of Flu A+B test were established using NP wash/aspirate specimens in multi-center clinical studies conducted at two U.S. trial sites and one Hong Kong trial site during the 2010–2011 respiratory season. A total of 1,502 prospective specimens (1,002 in the U.S and 500 in Hong Kong) were evaluated using the BD Veritor System for Rapid Detection of Flu A+B test and PCR. Five specimens were not evaluable because of data reconciliation issues, an additional 13 were excluded because of insufficient sample volume for reference method testing and 13 samples were excluded as “Result Invalid” (for an invalid rate of 0.9% [13/1,484]).

The prospective specimens consisted of NP washes and aspirates from symptomatic patients. 49% of the samples were from females and 51% from males. 56.6% were from patients less than or equal to 5 years of age. 21.9% of the patients tested were in the 6–21 year age group, 5.7% were from 22–59 years of age and 15.8% were obtained from people greater than or equal to age 60 (the patient age was not provided for 0.1% of samples).

The performance of the BD Veritor System for Rapid Detection of Flu A+B test was compared to an FDA-cleared Influenza A and B molecular assay (PCR).

The performance is presented in Table 1 below.

Table 1: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for All NP Wash/Aspirate Specimens – All Sites

| Clinical kit:<br>BD Flu A                                                                     | Reference PCR |       |       |
|-----------------------------------------------------------------------------------------------|---------------|-------|-------|
|                                                                                               | P             | N     | Total |
| P                                                                                             | 224           | 29    | 253   |
| N                                                                                             | 46            | 1,172 | 1,218 |
| Total                                                                                         | 270           | 1,201 | 1,471 |
| Reference Method: PCR<br>PPA: 83.0% (95% C.I. 78.0–87.0%)<br>NPA: 97.6% (95% C.I. 96.6–98.3%) |               |       |       |

| Clinical kit:<br>BD Flu B                                                                     | Reference PCR |       |       |
|-----------------------------------------------------------------------------------------------|---------------|-------|-------|
|                                                                                               | P             | N     | Total |
| P                                                                                             | 74            | 3     | 77    |
| N                                                                                             | 17            | 1,377 | 1,394 |
| Total                                                                                         | 91            | 1,380 | 1,471 |
| Reference Method: PCR<br>PPA: 81.3% (95% C.I. 72.1–88.0%)<br>NPA: 99.8% (95% C.I. 99.4–99.9%) |               |       |       |

An additional 263 frozen retrospective specimens were evaluated with the BD Veritor System for Rapid Detection of Flu A+B test. Twelve samples were excluded because there was insufficient sample volume for reference method testing, one sample was excluded as a PCR “Unresolved” and one sample was excluded as “Result Invalid” [for an invalid rate of 0.4% (1/250)]. The retrospective specimens consisted of NP washes and aspirates from symptomatic patients. 44.9% of the samples were from females and 55.1% from males. 87.5% were from patients less than or equal to 5 years of age.

The performance is presented in Table 2 below.

**Table 2: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for Retrospective NP Wash/Aspirate Specimens**

| Clinical kit:<br>BD Flu A                                                                     | Reference PCR |     |       |
|-----------------------------------------------------------------------------------------------|---------------|-----|-------|
|                                                                                               | P             | N   | Total |
| P                                                                                             | 58            | 2   | 60    |
| N                                                                                             | 5             | 184 | 189   |
| Total                                                                                         | 63            | 186 | 249   |
| Reference Method: PCR<br>PPA: 92.1% (95% C.I. 82.7–96.6%)<br>NPA: 98.9% (95% C.I. 96.2–99.7%) |               |     |       |

| Clinical kit:<br>BD Flu B                                                                     | Reference PCR |     |       |
|-----------------------------------------------------------------------------------------------|---------------|-----|-------|
|                                                                                               | P             | N   | Total |
| P                                                                                             | 29            | 2   | 31    |
| N                                                                                             | 10            | 208 | 218   |
| Total                                                                                         | 39            | 210 | 249   |
| Reference Method: PCR<br>PPA: 74.0% (95% C.I. 58.9–85.4%)<br>NPA: 99.0% (95% C.I. 96.6–99.7%) |               |     |       |

**Clinical Performance NP Swabs in Transport Media 2011–2012; U.S. and Japan Combined**

Performance characteristics for the BD Veritor System for Rapid Detection of Flu A+B test were established using NP swabs in transport media in multi-center studies conducted at six clinical trial sites located in geographically diverse areas within the United States and five clinical sites in Japan using a total of 292 samples.

The combined results are presented in Table 3 below.

**Table 3: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for NP Swabs in Transport Media – U.S. and Japan Combined**

| Clinical kit:<br>BD Flu A                                                                     | Reference PCR |     |       |
|-----------------------------------------------------------------------------------------------|---------------|-----|-------|
|                                                                                               | P             | N   | Total |
| P                                                                                             | 52            | 6   | 58    |
| N                                                                                             | 12            | 222 | 234   |
| Total                                                                                         | 64            | 228 | 292   |
| Reference Method: PCR<br>PPA: 81.3% (95% C.I. 70.0–88.9%)<br>NPA: 97.4% (95% C.I. 94.4–98.8%) |               |     |       |

| Clinical kit:<br>BD Flu B                                                                     | Reference PCR |     |       |
|-----------------------------------------------------------------------------------------------|---------------|-----|-------|
|                                                                                               | P             | N   | Total |
| P                                                                                             | 77            | 2   | 79    |
| N                                                                                             | 13            | 200 | 213   |
| Total                                                                                         | 90            | 202 | 292   |
| Reference Method: PCR<br>PPA: 85.6% (95% C.I. 76.8–91.4%)<br>NPA: 99.0% (95% C.I. 96.5–99.7%) |               |     |       |

**Clinical Performance NP Swabs in Transport Media 2011–2012; U.S**

Performance characteristics for the BD Veritor System for Rapid Detection of Flu A+B test were established using NP swabs in transport media in multi-center studies conducted at six clinical trial sites located in geographically diverse areas within the United States. A total of 217 prospective specimens were evaluated using the BD Veritor System for Rapid Detection of Flu A+B test and PCR. Two specimens could not be evaluated because of data reconciliation issues, one was eliminated because of an invalid control reading and 13 were excluded because the PCR results were unresolved.

The specimens consisted of NP swabs in transport media from symptomatic patients. 55.8% of the samples were from females and 44.2% from males. 16.1% were from patients less than or equal to 5 years of age, 25.3% were from patients 6–21 years of age, 47.5% were from patients 22–59 years of age and 11.1% were obtained from patients greater than or equal to 60 years of age.

The performance of the BD Veritor System for Rapid Detection of Flu A+B test was compared to an FDA-cleared Influenza A and B molecular assay (PCR).

The results are presented in Table 4 below.

**Table 4: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for NP Swabs in Transport Media – U.S.**

| Clinical kit:<br>BD Flu A                                                                     | Reference PCR |     |       |
|-----------------------------------------------------------------------------------------------|---------------|-----|-------|
|                                                                                               | P             | N   | Total |
| P                                                                                             | 52            | 6   | 58    |
| N                                                                                             | 12            | 131 | 143   |
| Total                                                                                         | 64            | 137 | 201   |
| Reference Method: PCR<br>PPA: 81.3% (95% C.I. 70.0–88.9%)<br>NPA: 95.6% (95% C.I. 90.8–98.0%) |               |     |       |

| Clinical kit:<br>BD Flu B                                                                   | Reference PCR |     |       |
|---------------------------------------------------------------------------------------------|---------------|-----|-------|
|                                                                                             | P             | N   | Total |
| P                                                                                           | 7             | 0   | 7     |
| N                                                                                           | 2             | 192 | 194   |
| Total                                                                                       | 9             | 192 | 201   |
| Reference Method: PCR<br>PPA: 77.8% (95% C.I. 45.3–93.7%)<br>NPA: 100% (95% C.I. 98.0–100%) |               |     |       |

**Clinical Performance NP Swabs in Transport Media 2011–2012; Japan**

Performance characteristics for the BD Veritor System for Rapid Detection of Flu A+B test were established using NP swabs in transport media in multi-center studies conducted at five clinical trial sites in Japan. A total of 93 prospective specimens were evaluated using the BD Veritor System for Rapid Detection of Flu A+B test and PCR. Two specimens were excluded as the results were undetermined with the comparator assay.

The specimens consisted of NP swabs in transport media from symptomatic patients. 49.5% of the samples were from females and 50.5% from males. 31.2% were from patients less than or equal to 5 years of age, 63.4% were from patients 6–21 years of age and 5.4% were from patients 22–59 years of age (there were no specimens from patients greater than or equal to 60 years of age).

The results are presented in Table 5 below.

**Table 5: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for NP Swabs in Transport Media – Japan**

| Clinical kit:<br>BD Flu A                                                              | Reference PCR |    |       |
|----------------------------------------------------------------------------------------|---------------|----|-------|
|                                                                                        | P             | N  | Total |
| P                                                                                      | 0             | 0  | 0     |
| N                                                                                      | 0             | 91 | 91    |
| Total                                                                                  | 0             | 91 | 91    |
| Reference Method: PCR<br>No Data for PPA Calculation<br>NPA: 100% (95% C.I. 95.9–100%) |               |    |       |

| Clinical kit:<br>BD Flu B                                                                     | Reference PCR |    |       |
|-----------------------------------------------------------------------------------------------|---------------|----|-------|
|                                                                                               | P             | N  | Total |
| P                                                                                             | 70            | 2  | 72    |
| N                                                                                             | 11            | 8  | 19    |
| Total                                                                                         | 81            | 10 | 91    |
| Reference Method: PCR<br>PPA: 86.4% (95% C.I. 77.3–92.2%)<br>NPA: 80.0% (95% C.I. 49.0–94.3%) |               |    |       |

## Reproducibility

The reproducibility of the BD Veritor System for Rapid Detection of Flu A+B test was evaluated at three clinical laboratory sites in 2010–2011. The reproducibility panel was composed of 30 simulated influenza A or B samples. These included moderate positive samples, low positive samples (near the assay limit of detection), high negative samples (i.e., containing very low concentrations of virus such that positive results occur ~5% of the time) and negative samples. The panel was tested by two operators at each site for five consecutive days. The results are summarized below.

| Reproducibility Results – Percent of Flu A Positives |                                          |                                          |                                          |                                        |
|------------------------------------------------------|------------------------------------------|------------------------------------------|------------------------------------------|----------------------------------------|
| Sample                                               | Site 1                                   | Site 2                                   | Site 3                                   | Total                                  |
| High negative H1N1 A                                 | 3.3% (1/30)<br>(95% C.I. 0.6–16.7%)      | 0.0% (0/30)<br>(95% C.I. 0.0–11.3%)      | 0.0% (0/30)<br>(95% C.I. 0.0–11.3%)      | 1.1% (1/90)<br>(95% C.I. 0.2–6.0%)     |
| Low positive H1N1 A                                  | 93.3% (28/30)<br>(95% C.I. 78.7–98.2%)   | 86.7% (26/30)<br>(95% C.I. 70.3–94.7%)   | 93.3% (28/30)<br>(95% C.I. 78.7–98.2%)   | 91.1% (82/90)<br>(95% C.I. 83.4–95.4%) |
| Moderate positive H1N1 A                             | 100.0% (30/30)<br>(95% C.I. 88.6–100.0%) | 96.7% (29/30)<br>(95% C.I. 83.3–99.4%)   | 100.0% (30/30)<br>(95% C.I. 88.6–100.0%) | 98.9% (89/90)<br>(95% C.I. 94.0–99.8%) |
| High negative H3N2 A                                 | 16.7% (5/30)<br>(95% C.I. 7.3–33.6%)     | 3.3% (1/30)<br>(95% C.I. 0.6–16.7%)      | 0.0% (0/30)<br>(95% C.I. 0.0–11.3%)      | 6.7% (6/90)<br>(95% C.I. 3.1–13.8%)    |
| Low positive H3N2 A                                  | 93.3% (28/30)<br>(95% C.I. 78.7–98.2%)   | 86.7% (26/30)<br>(95% C.I. 70.3–94.7%)   | 93.3% (28/30)<br>(95% C.I. 78.7–98.2%)   | 91.1% (82/90)<br>(95% C.I. 83.4–95.4%) |
| Moderate positive H3N2 A                             | 100.0% (30/30)<br>(95% C.I. 88.6–100.0%) | 100.0% (30/30)<br>(95% C.I. 88.6–100.0%) | 96.7% (29/30)<br>(95% C.I. 83.3–99.4%)   | 98.9% (89/90)<br>(95% C.I. 94.0–99.8%) |
| Negatives                                            | 0.8% (1/120)<br>(95% C.I. 0.1–4.6%)      | 0.0% (0/120)<br>(95% C.I. 0.0–3.1%)      | 0.0% (0/119)<br>(95% C.I. 0.0–3.1%)      | 0.3% (1/359)<br>(95% C.I. 0.0–1.6%)    |

| Reproducibility Results – Percent of Flu B Positives |                                        |                                          |                                          |                                        |
|------------------------------------------------------|----------------------------------------|------------------------------------------|------------------------------------------|----------------------------------------|
| Sample                                               | Site 1                                 | Site 2                                   | Site 3                                   | Total                                  |
| High negative B                                      | 3.3% (1/30)<br>(95% C.I. 0.6–16.7%)    | 0.0% (0/30)<br>(95% C.I. 0.0–11.3%)      | 0.0% (0/30)<br>(95% C.I. 0.0–11.3%)      | 1.1% (1/90)<br>(95% C.I. 0.2–6.0%)     |
| Low positive B                                       | 90.0% (27/30)<br>(95% C.I. 74.4–96.5%) | 63.3% (19/30)<br>(95% C.I. 45.5–78.1%)   | 82.8% (24/29)<br>(95% C.I. 65.5–92.4%)   | 78.7% (70/89)<br>(95% C.I. 69.0–85.9%) |
| Moderate positive B                                  | 96.7% (29/30)<br>(95% C.I. 83.3–99.4%) | 100.0% (30/30)<br>(95% C.I. 88.6–100.0%) | 100.0% (30/30)<br>(95% C.I. 88.6–100.0%) | 98.9% (89/90)<br>(95% C.I. 94.0–99.8%) |
| Negatives                                            | 0.0% (0/210)<br>(95% C.I. 0.0–1.8%)    | 0.0% (0/210)<br>(95% C.I. 0.0–1.8%)      | 0.0% (0/210)<br>(95% C.I. 0.0–1.8%)      | 0.0% (0/630)<br>(95% C.I. 0.0–0.6%)    |

## Analytical Studies

### Analytical Sensitivity (Limit of Detection)

The limit of detection (LOD) for the BD Veritor System for Rapid Detection of Flu A+B test was established for a total of 8 influenza strains: 5 influenza A and 3 influenza B. The LOD for each strain represents the lowest concentration producing a positivity rate of ≥95% based on testing 20 to 60 replicates.

| Type | Influenza Viral Strain   | Calculated LOD (TCID <sub>50</sub> /mL) | Calculated LOD (EID <sub>50</sub> /mL) | No. Positive / Total | % Positive |
|------|--------------------------|-----------------------------------------|----------------------------------------|----------------------|------------|
| A    | A/Brisbane/10/2007 H3N2  | 7.27 x 10 <sup>2</sup>                  | N/A                                    | 57/60                | 95%        |
| A    | A/Brisbane/59/2007 H1N1  | 3.30 x 10 <sup>2</sup>                  | N/A                                    | 57/60                | 95%        |
| A    | A/California/7/2009 H1N1 | 5.00 x 10 <sup>3</sup>                  | N/A                                    | 57/60                | 95%        |
| A    | A/Victoria/3/75 H3N2     | 3.11 x 10 <sup>3</sup>                  | N/A                                    | 59/60                | 98.3%      |
| A    | A/Anhui/1/2013 H7N9      | N/A                                     | 5.42 x 10 <sup>6</sup>                 | 59/60                | 98.3%      |
| B    | B/Brisbane/60/2008       | 7.42 x 10 <sup>3</sup>                  | N/A                                    | 58/60                | 96.7%      |
| B    | B/Florida/4/2006         | 1.30 x 10 <sup>3</sup>                  | N/A                                    | 58/60                | 96.7%      |
| B    | B/Lee/40                 | 4.44 x 10 <sup>4</sup>                  | N/A                                    | 20/20                | 100%       |

TCID<sub>50</sub>/mL = 50% Tissue Culture Infectious Dose

EID<sub>50</sub>/mL = 50% Egg Infectious Dose

## Strain Reactivity with Influenza A and B Viruses

The BD Veritor System for Rapid Detection of Flu A+B test was evaluated using a panel of influenza strains. Each strain was diluted and tested in triplicate until a point where not all of the replicates were positive. The dilution prior to that is provided in the table below as a minimal detected concentration. All influenza A strains showed positive Flu A test results and negative Flu B test results. Conversely, all of the influenza B strains showed positive Flu B test results and negative Flu A test results. Although this test has been shown to detect novel avian influenza A (H7N9) and H3N2v cultured viruses the performance characteristics of this device with clinical specimens that are positive for novel avian influenza A (H7N9) and H3N2v influenza viruses have not been established. The BD Veritor System Flu A+B assay can distinguish between influenza A and B viruses, but it cannot differentiate influenza A subtypes.

| Strain                                   | Subtype | Minimal Detected Concentration             |
|------------------------------------------|---------|--------------------------------------------|
| A/Brisbane/59/2007                       | H1N1    | $3.3 \times 10^2$ TCID <sub>50</sub> /mL*  |
| A/California/7/2009                      | H1N1    | $5.0 \times 10^3$ TCID <sub>50</sub> /mL*  |
| A/Denver/1/57                            | H1N1    | $4.45 \times 10^4$ CEID <sub>50</sub> /mL  |
| A/FM/1/47                                | H1N1    | $7.91 \times 10^4$ CEID <sub>50</sub> /mL  |
| A/Fujian-Gulou/1896/2009                 | H1N1    | $4.5 \times 10^5$ CEID <sub>50</sub> /mL   |
| A/Mal/302/54                             | H1N1    | $2.22 \times 10^5$ CEID <sub>50</sub> /mL  |
| A/New Caledonia/20/1999                  | H1N1    | $2.5 \times 10^3$ TCID <sub>50</sub> /mL   |
| A/New Jersey/8/76                        | H1N1    | $1.58 \times 10^3$ CEID <sub>50</sub> /mL  |
| A/NWS/33                                 | H1N1    | $1.58 \times 10^4$ CEID <sub>50</sub> /mL  |
| A/PR/8/34                                | H1N1    | $6.31 \times 10^2$ TCID <sub>50</sub> /mL  |
| A/Solomon Island/03/2006                 | H1N1    | $2.5 \times 10^3$ TCID <sub>50</sub> /mL   |
| A/Washington/24/2012                     | H1N1    | $3.16 \times 10^4$ EID <sub>50</sub> /mL   |
| A/Weiss/43                               | H1N1    | $7.03 \times 10^6$ CEID <sub>50</sub> /mL  |
| A/WS/33                                  | H1N1    | $7.91 \times 10^2$ CEID <sub>50</sub> /mL  |
| A/Aichi/2/68                             | H3N2    | $7.91 \times 10^3$ CEID <sub>50</sub> /mL  |
| A/Brisbane/10/2007                       | H3N2    | $7.27 \times 10^2$ TCID <sub>50</sub> /mL* |
| A/California/02/2014                     | H3N2    | $1.45 \times 10^2$ TCID <sub>50</sub> /mL  |
| A/Hong Kong/8/68                         | H3N2    | $8.89 \times 10^4$ CEID <sub>50</sub> /mL  |
| A/Moscow/10/99                           | H3N2    | $5.8 \times 10^6$ TCID <sub>50</sub> /mL   |
| A/Perth/16/2009                          | H3N2    | $1.0 \times 10^6$ TCID <sub>50</sub> /mL   |
| A/Port Chalmers/1/73                     | H3N2    | $3.95 \times 10^4$ CEID <sub>50</sub> /mL  |
| A/Switzerland/9715293/2013               | H3N2    | $3.25 \times 10^2$ TCID <sub>50</sub> /mL  |
| A/Texas/50/2012                          | H3N2    | $1.75 \times 10^3$ TCID <sub>50</sub> /mL  |
| A/Wisconsin/67/2005                      | H3N2    | $2.5 \times 10^5$ TCID <sub>50</sub> /mL   |
| A/Victoria/3/75                          | H3N2    | $3.11 \times 10^3$ TCID <sub>50</sub> /mL* |
| A/Indiana/08/2011                        | H3N2v   | $1 \times 10^4$ TCID <sub>50</sub> /mL     |
| A/Indiana/10/2011                        | H3N2v   | $7.9 \times 10^6$ CEID <sub>50</sub> /mL   |
| A/Kansas/13/2009                         | H3N2v   | $1.0 \times 10^3$ TCID <sub>50</sub> /mL   |
| A/Minnesota/11/2010                      | H3N2v   | $7.9 \times 10^5$ CEID <sub>50</sub> /mL   |
| A/Pennsylvania/14/2010                   | H3N2v   | $1.26 \times 10^6$ CEID <sub>50</sub> /mL  |
| A/West Virginia/06/2011                  | H3N2v   | $7.9 \times 10^3$ TCID <sub>50</sub> /mL   |
| A/Anhui/01/2005                          | H5N1    | 0.512 HA                                   |
| A/Vietnam/1203/2004                      | H5N1    | 0.512 HA                                   |
| A/Northern Pintail/Washington/40964/2014 | H5N2    | $6.28 \times 10^5$ EID <sub>50</sub> /mL   |
| A/Pheasant/New Jersey/1355/1998          | H5N2    | 0.256 HA                                   |
| A/Gyrfalcon/Washington/41088-6/2014      | H5N8    | $1.98 \times 10^6$ EID <sub>50</sub> /mL   |
| A/Mallard/Netherlands/12/2000            | H7N7    | 0.256 HA                                   |
| A/Anhui/1/2013                           | H7N9    | $5.42 \times 10^6$ CEID <sub>50</sub> /mL* |
| A/Chicken/Hong Kong/G9/1997              | H9N2    | 1.024 HA                                   |

\*Values taken from preceding Analytical Limit of Detection Table

- a. EID<sub>50</sub> = 50% Egg Infectious Dose
- b. TCID<sub>50</sub> = 50% Tissue Culture Infectious Dose
- c. CEID<sub>50</sub> = 50% Chicken Embryo Infectious Dose
- d. HA = Hemagglutination Assay

| Strain                                         | Minimal Detected Concentration                 |
|------------------------------------------------|------------------------------------------------|
| B/Brazil/178/96                                | 2.32 x 10 <sup>4</sup> TCID <sub>50</sub> /mL  |
| B/Brisbane/33/2008 (Victoria Lineage)          | 2.45 x 10 <sup>5</sup> CEID <sub>50</sub> /mL  |
| B/Brisbane/60/2008                             | 7.42 x 10 <sup>3</sup> TCID <sub>50</sub> /mL* |
| B/Brisbane/72/97                               | 1.00 x 10 <sup>4</sup> TCID <sub>50</sub> /mL  |
| B/Canada/548/99                                | >0.64 HA                                       |
| B/Egypt/393/99                                 | >1.28 HA                                       |
| B/Florida/2/2006                               | 1.08 x 10 <sup>4</sup> TCID <sub>50</sub> /mL  |
| B/Florida/4/2006                               | 1.30 x 10 <sup>3</sup> TCID <sub>50</sub> /mL* |
| B/Fujian/93/97                                 | 3.95 x 10 <sup>5</sup> TCID <sub>50</sub> /mL  |
| B/Fukushima/220/99                             | 9.33 x 10 <sup>2</sup> TCID <sub>50</sub> /mL  |
| B/Guangdong-Liwan/1133/2014 (Yamagata Lineage) | 9.0 x 10 <sup>5</sup> CEID <sub>50</sub> /mL   |
| B/Guangxi/547/98                               | 2.32 x 10 <sup>5</sup> TCID <sub>50</sub> /mL  |
| B/Hawaii/01/97                                 | >6.4 HA                                        |
| B/Hong Kong/5/72                               | 1.11 x 10 <sup>4</sup> CEID <sub>50</sub> /mL  |
| B/Hong Kong/219/98                             | >1.0 HA                                        |
| B/Hong Kong/259/2010 (Victoria Lineage)        | 1.35 x 10 <sup>6</sup> CEID <sub>50</sub> /mL  |
| B/Jiangsu/10/2003                              | 1.16 x 10 <sup>4</sup> TCID <sub>50</sub> /mL  |
| B/Johannesburg/5/99                            | 3.95 x 10 <sup>4</sup> TCID <sub>50</sub> /mL  |
| B/Lee/40                                       | 4.44 x 10 <sup>4</sup> CEID <sub>50</sub> /mL* |
| B/Lisbon/03/96                                 | >0.08 HA                                       |
| B/Malaysia/2506/2004                           | 5.0 x 10 <sup>4</sup> TCID <sub>50</sub> /mL   |

| Strain                                    | Minimal Detected Concentration                |
|-------------------------------------------|-----------------------------------------------|
| B/Maryland/1/59                           | 3.51 x 10 <sup>2</sup> CEID <sub>50</sub> /mL |
| B/Massachusetts/2/2012 (Yamagata Lineage) | 1.25 x 10 <sup>6</sup> CEID <sub>50</sub> /mL |
| B/Mass/3/66                               | 1.58 x 10 <sup>5</sup> CEID <sub>50</sub> /mL |
| B/Montana/5/2012                          | 3.14 x 10 <sup>5</sup> EID <sub>50</sub> /mL  |
| B/Ohio/11/96                              | >0.16 HA                                      |
| B/Ohio/1/05                               | 1.34 x 10 <sup>5</sup> TCID <sub>50</sub> /mL |
| B/Phuket/3073/2013                        | 6.08 x 10 <sup>3</sup> TCID <sub>50</sub> /mL |
| B/Puerto Mont/10427/98                    | 0.02 HA                                       |
| B/Russia/69                               | 3.9 x 10 <sup>2</sup> TCID <sub>50</sub> /mL  |
| B/Shandong/7/97                           | 1.58 x 10 <sup>6</sup> TCID <sub>50</sub> /mL |
| B/Shanghai/04/97                          | 1.58 x 10 <sup>5</sup> TCID <sub>50</sub> /mL |
| B/Shenzhen/135/97                         | 3.16 x 10 <sup>4</sup> TCID <sub>50</sub> /mL |
| B/Sichuan/116/96                          | 0.016 HA                                      |
| B/Taiwan/2/62                             | 2.81 x 10 <sup>2</sup> CEID <sub>50</sub> /mL |
| B/Texas/06/2011 (Yamagata Lineage)        | 6.2 x 10 <sup>5</sup> CEID <sub>50</sub> /mL  |
| B/Texas/02/2013 (Victoria Lineage)        | 2.75 x 10 <sup>4</sup> CEID <sub>50</sub> /mL |
| B/Utah/09/2014 (Yamagata Lineage)         | 6.3 x 10 <sup>3</sup> CEID <sub>50</sub> /mL  |
| B/Victoria/504/00                         | 4.64 x 10 <sup>4</sup> TCID <sub>50</sub> /mL |
| B/Wisconsin/01/2010 (Yamagata Lineage)    | 7.0 x 10 <sup>2</sup> CEID <sub>50</sub> /mL  |
| B/Yamagata/16/88                          | 9.75 x 10 <sup>3</sup> TCID <sub>50</sub> /mL |
| B/Yamanashi/166/98                        | 4.88 x 10 <sup>4</sup> TCID <sub>50</sub> /mL |

\*Values taken from preceding Analytical Limit of Detection Table

a. EID<sub>50</sub> = 50% Egg Infectious Dose

b. TCID<sub>50</sub> = 50% Tissue Culture Infectious Dose

c. CEID<sub>50</sub> = 50% Chicken Embryo Infectious Dose

d. HA = Hemagglutination Assay

On January 12, 2017 the US FDA published an order reclassifying antigen based rapid influenza virus detection systems (RIDTs) intended to detect influenza virus antigen directly from clinical specimens from Class I to Class II with special controls. One of these special controls is a requirement for annual analytical reactivity testing of contemporary influenza strains identified by the Centers for Disease control (CDC) using a standardized dilution protocol. Results obtained using the BD Veritor System are available for viewing at [bd.com/veritorsystem](http://bd.com/veritorsystem).

## Analytical Specificity (Cross Reactivity)

The BD Veritor System for Rapid Detection of Flu A+B test was evaluated with a total of 51 microorganisms. The 37 bacteria and yeast were tested at a target concentration of approximately  $10^7$  CFU/mL (CFU – Colony Forming Units) with the exception of *Staphylococcus aureus*, which was tested at a final concentration of  $10^6$  CFU/mL. The 14 viruses were evaluated at concentrations of  $10^3$  to  $10^{10}$  TCID<sub>50</sub>/mL. Of the 51 microorganisms tested, none showed cross-reactivity in either the Flu A or Flu B tests.

|                                    |                                                    |                                    |
|------------------------------------|----------------------------------------------------|------------------------------------|
| <i>Bacteroides fragilis</i>        | <i>Neisseria meningitidis</i>                      | <i>Streptococcus</i> sp. Group G   |
| <i>Bordetella pertussis</i>        | <i>Neisseria mucosa</i>                            | <i>Streptococcus salivarius</i>    |
| <i>Candida albicans</i>            | <i>Neisseria</i> sp. ( <i>Neisseria perflaus</i> ) | <i>Veillonella parvula</i>         |
| <i>Chlamydia pneumoniae</i>        | <i>Neisseria subflava</i>                          | Adenovirus, type 1                 |
| <i>Corynebacterium diphtherium</i> | <i>Peptostreptococcus anaerobius</i>               | Adenovirus, type 7                 |
| <i>Escherichia coli</i>            | <i>Porphyromonas asaccharolyticus</i>              | Cytomegalovirus                    |
| <i>Fusobacterium nucleatum</i>     | <i>Prevotella oralis</i>                           | Enterovirus                        |
| <i>Haemophilus influenzae</i>      | <i>Propionibacterium acnes</i>                     | Epstein Barr Virus                 |
| <i>Haemophilus parainfluenzae</i>  | <i>Proteus mirabilis</i>                           | HSV Type 1                         |
| <i>Kingella kingae</i>             | <i>Pseudomonas aeruginosa</i>                      | Human Coronavirus OC43             |
| <i>Klebsiella pneumoniae</i>       | <i>Serratia marcescens</i>                         | Human Coronavirus 2229E            |
| <i>Lactobacillus</i> sp.           | <i>Staphylococcus aureus</i>                       | Human metapneumovirus (HMPV-27 A2) |
| <i>Legionella</i> sp.              | <i>Staphylococcus epidermidis</i>                  | Human Parainfluenza                |
| <i>Moraxella catarrhalis</i>       | <i>Streptococcus mutans</i>                        | Measles virus                      |
| <i>Mycobacterium tuberculosis</i>  | <i>Streptococcus pneumoniae</i>                    | Mumps virus                        |
| <i>Mycoplasma pneumoniae</i>       | <i>Streptococcus pyogenes</i>                      | Respiratory syncytial virus        |
| <i>Neisseria gonorrhoeae</i>       | <i>Streptococcus</i> sp. Group C                   | Rhinovirus                         |

## Interfering Substances

Various substances were evaluated with the BD Veritor System for Rapid Detection of Flu A+B test. These substances included whole blood (2%) and various medications. No interference was noted with this assay for any of the substances tested.

| Substance                | Concentration | Substance                    | Concentration |
|--------------------------|---------------|------------------------------|---------------|
| 4-Acetamidophenol        | 10 mg/mL      | Homeopathic Allergy Medicine | 10 mg/mL      |
| Acetylsalicylic acid     | 20 mg/mL      | Ibuprofen                    | 10 mg/mL      |
| Albuterol                | 0.083 mg/mL   | Loratidine                   | 100 ng/mL     |
| Amantadine Hydrochloride | 500 ng/mL     | Menthol Throat Lozenges      | 10 mg/mL      |
| Ayr Saline Nasal Gel     | 10 mg/mL      | Mometasone                   | 500 ng/mL     |
| Beclomethasone           | 500 ng/mL     | Mupirocin                    | 500 ng/mL     |
| Budesonide               | 500 ng/mL     | Oseltamivir                  | 500 ng/mL     |
| Chlorpheniramine maleate | 5 mg/mL       | Oxymetazoline                | 0.05 mg/mL    |
| Dexamethasone            | 10 mg/mL      | Phenylephrine                | 1 mg/mL       |
| Dextromethorphan         | 10 mg/mL      | Pseudoephedrine HCl          | 20 mg/mL      |
| Diphenhydramine HCl      | 5 mg/mL       | Purified Mucin Protein       | 1 mg/mL       |
| Fexofenadine             | 500 ng/mL     | Ribavirin                    | 500 ng/mL     |
| FluMist                  | 1%            | Rimantadine                  | 500 ng/mL     |
| Flunisolide              | 500 ng/mL     | Three OTC mouthwashes        | 5%            |
| Fluticasone              | 500 ng/mL     | Tobramycin                   | 500 ng/mL     |
| Four OTC nasal sprays    | 10%           | Triamcinolone                | 500 ng/mL     |
| Four OTC throat drops    | 25%           | Whole Blood                  | 2%            |
| Guaiacol Glyceryl Ether  | 20 mg/mL      | Zanamivir                    | 1 mg/mL       |

Of the 44 substances tested in this study, none exhibited interfering reactions when tested with influenza A and influenza B positive samples. Based on the data, the substances tested at the indicated concentration levels did not interfere with the BD Veritor System for Rapid Detection of Flu A+B test.

## AVAILABILITY

### Cat. No. Description

|        |                                                                |
|--------|----------------------------------------------------------------|
| 256041 | BD Veritor™ System for Rapid Detection of Flu A+B, 30 tests    |
| 256042 | BD Veritor™ System for Rapid Detection of RSV, 30 tests        |
| 256051 | BD Veritor™ System Flu A+B Control Swab Set, 10 pairs of swabs |
| 256066 | BD Veritor™ Plus Analyzer                                      |

To network a BD Veritor Plus Analyzer to an LIS, contact BD Technical Services for details.

## REFERENCES

1. Simonsen L., Fukuda K, Schonberger LB, Cox NJ. Impact of influenza epidemics on hospitalizations. *J. Infect. Dis.* 2000; **181**:831–7.
2. Thompson WW, Shay DK, Weintraub E, et al. Mortality associated with influenza and respiratory syncytial virus in the United States. *JAMA* 2003; **289**: 179–86.
3. Treanor, J.J., Hayden, F.G., Vrooman, P.S., et al. 2000. Efficacy and safety of the oral neuraminidase inhibitor oseltamivir in treating acute influenza: a randomized controlled trial. *JAMA.* **283**:1016–1024.
4. Kaiser, L., Couch, R.B., Galasso, G.J., Glezen, W.P., Webster, R.G., Wright, P.F., and Hayden, F.G. 1999. First international symposium on influenza and other respiratory viruses: summary and overview Kapalua, Maui, Hawaii, December 4–6, 1998. *Antiviral Res.* **42**:149–176.
5. Cox, N.J., and Bender, C.A. 1995. The molecular epidemiology of influenza viruses. *Virology*, **6**:359–370.
6. Todd, S.J., Minnich, L., and Waner, J.L. 1995. Comparison of rapid immunofluorescence procedure with TestPack RSV and Directigen Flu A for diagnosis of respiratory syncytial virus and influenza A virus. *J. Clin. Microbiol.* **33**:1650–1651.
7. Harris, P.O. 1989. Clinical relevance and efficient detection of seven major respiratory viruses. *ACL.* p. 15–19.
8. McElhaney, J.E., Gravenstein, S., Krause, P., Hooton, J.W., Upshaw, C.M., and Drinka, P. 1998. Assessment of markers of the cell-mediated immune response after influenza virus infection in frail older adults. *Clin. Diag. Lab. Immunol.* **5**:840–844.
9. Fan, J., Henrickson, K.J., and Savatski, L.L. 1998. Rapid simultaneous diagnosis of infections with respiratory syncytial viruses A and B, influenza viruses A and B, and human parainfluenza virus types 1, 2, and 3 by multiplex quantitative reverse transcription-polymerase chain reaction-hybridization assay (hexaplex). *Clin. Infect. Disease* **26**:1397–1402.
10. Wright, K.E., Wilson, G.A.R., Novosad, D., Dimock, C., Tan, D., and Weber, J.M. 1995. Typing and subtyping of influenza viruses in clinical samples by PCR. *J. Clin. Microbiol.* **33**:1180–1184.
11. Kendal, A.P. 1985. Influenza Viruses. p. 341–357. *Laboratory Diagnosis of Viral Infections*, *In* H. Lennette, (ed.) Marcel Dekker, Inc., New York.
12. McQuillen, J., Madeley, C.R., and Kendal, A.P. 1985. Monoclonal antibodies for the rapid diagnosis of influenza A and B virus infections by immunofluorescence. *Lancet*. *ii*: 911–914.
13. Guenther, S.H., and Linnemann, C.C., Jr. 1988. Indirect immunofluorescence assay for rapid diagnosis of influenza virus. *Laboratory Medicine*. **19**: 581–583.
14. Minnick, L.L., and Ray, C.G. 1986. Early testing of cell cultures for detection of hemadsorbing viruses. *J. Clin. Microbiol.* **25**: 421–422.
15. Schmidt, N.J., Ota, M., Gallo, D., and Fox, V.L. 1982. Monoclonal antibodies for rapid, strain specific identification of influenza virus isolates. *J. Clin. Microbiol.* **16**: 763–765.
16. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
17. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* **17**:53–80.
18. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S Government Printing Office, Washington, D.C.
19. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EECP). *Office Journal L262*, 17/10/2000, p.021–0045.

Technical Information: In the United States, contact BD Technical Service and Support at 1.800.638.8663 or bd.com.

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## Change History

| Revision | Date    | Change Summary                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|----------|---------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| (12)     | 2020-04 | <p>Changed name of BD Synapsys™ Microbiology Informatics Solution to BD Synapsys™ Informatics Solution.</p> <p>Additional result documentation capabilities made possible with the implementation of the BD Synapsys Informatics Solution, and with the addition of the BD Veritor InfoScan module and BD Veritor Plus Connect (replacing the BD Veritor System Reader and the BD InfoSync Module).</p> <p>In sections C and D, revised Steps 10C and 10D regarding the use of the BD Veritor Plus Analyzer with the BD Veritor InfoScan.</p> <p>Added access information to obtain the document from <a href="http://bd.com/e-labeling">bd.com/e-labeling</a>.</p> |

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For Use: "For Investigational Use Only"

Lower limit of temperature / Долен лимит на температурата / Dolní hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Κατώτερο όριο θερμοκρασίας / Limite inferior de temperatura / Alumine temperaturipiir / Limite inférieure de température / Najnižja dovoljena temperatura / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Temperatuurancij tēmēnri ruksat šeri / 하한 온도 / Žemiausia laikymo temperatūra / Temperatūros žemākā robeža / Laagste temperatuurlimiet / Nedre temperaturgrænse / Dolna granica temperatury / Limite minimo de temperatura / Limitā minimā de temperaturā / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sicaklık alt sınırı / Минимальная температура / 温度下限

#### CONTROL

Control / Контроль / Kontrola / Kontrol / Kontrolle / Mátruras / Kontroll / Contrôle / Controlla / Бақылау / 컨트롤 / Kontrolé / Kontrolle / Controle / Controlo / Контроль / kontroll / Контроль / 对照

#### CONTROL+

Positive control / Положительный контроль / Pozitívni kontrola / Positiv kontrol / Positive Kontrolle / Θετικός μάρτυρας / Control positivo / Positive kontrol / Contrôle positif / Pozitívna kontrola / Pozitív kontrol / Controllo positivo / Оң бақылау / 양성 컨트롤 / Teigiama kontrolė / Pozitívna kontrola / Positiveve controle / Kontrola dodatnia / Controllo positivo / Control pozitiv / Положительный контроль / Pozitif kontrol / Позитивный контроль / 阳性对照试剂

#### CONTROL-

Negative control / Отрицательный контроль / Negatívni kontrola / Negativ kontrol / Negative Kontrolle / Αρνητικός μάρτυρας / Control negativo / Negativeve kontrol / Contrôle négatif / Negativna kontrola / Negativ kontrol / Controllo negativo / Негативтік бақылау / 음성 컨트롤 / Neigiama kontrolė / Negativá kontrol / Negativeve controle / Kontrola ujemna / Kontrollo negativo / Control negativ / Отрицательный контроль / Negatif kontrol / Негативный контроль / 阴性对照试剂

#### STERILE

Method of sterilization: ethylenoxide / Метод на стерилизация: этиленов оксид / Způsob sterilizace: ethylenoxid / Steriliseringmetode: ethylenoxid / Sterilisationsmethode: Ethylenoxid / Μέθοδος αποστείρωσης: αιθυλενοξείδιο / Método de esterilização: óxido de etileno / Steriliseringmetode: etüleenoksiid / Méthode de stérilisation: oxyde d'éthylène / Metoda sterilizacije: etilen oksid / Sterilizálás módzere: etilén-oxid / Metoda di sterilizzazione: ossido di etilene / Стерилизация әдісі – этилен тотығы / 소독 방법: 에틸렌옥사이드 / Sterilizavimo būdas: etileno oksidas / Sterilizēšanas metode: etilēnoksīds / Gesteriliseerd met behulp van ethyleenoxide / Steriliseringsmetode: ethylenoksid / Metoda sterylizacji: tlenek etylu / Metoda de esterilização: óxido de etileno / Metodă de sterilizare: oxid de etilenă / Метод стерилизации: этиленоксид / Metóda sterilizácie: etylénoxid / Metoda sterilizacije: etilen oksid / Steriliseringsmetod: etenoxid / Sterilizasyon yöntemi: etilen oksit / Метод стерилизації: етиленоксидом / 灭菌方法: 环氧乙烷

#### STERILE R

Method of sterilization: irradiation / Метод на стерилизация: ирадиация / Způsob sterilizace: záření / Steriliseringmetode: bestråling / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Steriliseringmetode: kiirgus / Méthode de stérilisation: irradiation / Metoda sterilizacije: zračenje / Sterilizálás módzere: besugárzás / Metoda di sterilizzazione: irradiazione / Стерилизация әдісі – сәуле түсіру / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizēšanas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringmetode: bestråling / Metoda sterylizacji: napromienianie / Método de esterilização: irradiação / Metoda de sterilizare: iradiere / Метод стерилизации: облучение / Metóda sterilizácie: ožiarovanie / Metoda sterilizacije: ozračevanje / Steriliseringsmetod: strålning / Sterilizasyon yöntemi: iradyasyon / Метод стерилизації: опроміннювання / 灭菌方法: 辐射

Biological Risks / Биологични рискове / Biologická rizika / Biologisk fare / Biogefährdung / Βιολογικοί κίνδυνοι / Riesgos biológicos / Biologised riskid / Risques biologiques / Biološki rizik / Biológiai veszélyes / Rischio biologico / Биологичнык тауекелдер / 생물학적 위험 / Biologinis pavojus / Biologické riziko / Biologisch risico / Zagrożenia biologiczne / Perigo biológico / Riscuri biologice / Биологическая опасность / Biologické riziko / Biološki rizici / Biologisk risk / Biyolojik Riskler / Биологична небезпека / 生物学风险

Caution, consult accompanying documents / Внимание, направте справка в придружаващите документи / Pozor! Prostudujte si příloženou dokumentaci! / Forsigtig, se ledsagende dokumenter / Achtung, Begleitdokumente beachten / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα / Precaución, consultar la documentación adjunta / Ettevaatus! Lugeda kaasnevat dokumentatsiooni / Attention, consulter les documents joints / Upozorenje, koristi prateću dokumentaciju / Figyelem! Olvassa el a mellékelt tájékoztatót / Attenzione: consultare la documentazione allegata / Абайлаңыз, тиісті құжаттармен танысыңыз / 주의, 동봉된 설명서 참조 / Dmesio, zürüklite pridedamus dokumentus / Piesardzība, skatīt pavaddokumentus / Voorzichtig, gaadpleeg bijgevoegde documenten / Forsikligt, se vedlagt dokumentasjon / Należy zapoznać się z dołączonymi dokumentami / Cuidado, consulte a documentação fornecida / Atenție, consultați documentele însoțitoare / Внимание: см. прилагаемое документацию / Výstraha, pozri sprievodné dokumenty / Pažnja! Pogledajte priložena dokumenta / Obs! Se medföljande dokumentation / Dikkat, birlikte verilen belgelere başvurun / Увара: див. супутню документацию / 小心, 请参阅附带文档。

Upper limit of temperature / Горен лимит на температурата / Horní hranice teploty / Øvre temperaturgrænse / Temperaturobergrenze / Ανώτερο όριο θερμοκρασίας / Limite supérieur de température / Ülemine temperatuuripiir / Limite supérieure de température / Gornja dozvoljena temperatura / Felső hőméréskelti határ / Limite superiore di temperatura / Temperatuurkynn rúskat epligen joagary uer / 상한 온도 / Aukščiausia laikymo temperatūra / Augšējā temperatūras robeža / Hoogste temperatuurlimiet / Øvre temperaturgrense / Górna granica temperatury / Limite máximo de temperatura / Limită maximă de temperatură / Верхний предел температуры / Horná hranica teploty / Gornja granica temperature / Øvre temperaturgräns / Sicaklık üst sınırı / Максимальна температура / 温度上限

Keep dry / Пазете сухо / Skladujte v suchém prostředí / Opbevares tørt / Trocklagern / Φυλάξτε το στεγνό / Mantener seco / Hoida kuiv / Conserver au sec / Držati na suhom / Százraz helyen tartandó / Tenere all'asciutto / Kırkraz kıyindeusta / 건조 상태 유지 / Laikykite sausai / Uzglabāt saus / Droog houden / Hoides tørt / Przechowywać w stanie suchym / Manter seco / A se feri de umezeală / He dopuskać padańia wlo / Uchovávať v suchu / Držite na suvom mestu / Förvaras tørt / Kuru bir şekilde muhafaza edin / Берегти від вологи / 请保持干燥

Collection time / Време на събиране / Čas odběru / Opsamlingsstidspunkt / Entnahmehurzeit / Ώρα συλλογής / Hora de recogida / Kogumisaeg / Heure de prélèvement / Sati prikupljanja / Mintavétel időpontja / Ora di raccolta / Жинау уакыты / 수집 시간 / Paėmimo laikas / Savākšanas laiks / Verzameltijd / Tid prøvetaking / Godzina pobrania / Hora de colecta / Ora colectării / Время сбора / Doba odboru / Vreme prikupljanja / Uppsamlingsstid / Toplama zamani / Час забору / 采集时间

Peel / Обенете / Otefievte zde / Abn / Abziehen / Αποκαλλήστε / Desprendre / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Ўстіті қабатын алып таста / 벗기다 / Pléstí čia / Atfīmēt / Schillen / Trek av / Oderwać / Destacar / Se dezlipeste / Отклеить / Odrhňte / Oljuštiti / Dra isår / Ayırma / Відклеїти / 撕下

Perforation / Перфорация / Perforace / Perforering / Διάρθρωση / Perforación / Perforatsioon / Perforacija / Perforálás / Perforazione / Тесик тесу / 절취선 / Perforacia / Perforăcia / Perforatie / Perforacja / Perforação / Perforare / Перфорация / Perforálás / Perforasyon / Перфорация / 穿孔

Do not use if package damaged / Не използвайте, ако опаковката е повредена / Nepoužívejte, je-li obal poškozený / Má ikke anvendes hvis emballagen er beskadiget / Inhal beschädigter Packung nicht verwenden / Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά. / No usar si el paquete está dañado / Mitte kasutada, kui pakend on kahjustatud / Ne pas l'utiliser si l'emballage est endommagé / Ne koristiti ako je oštećeno pakiranje / Ne használni, ha a csomagolás sérült / Non usare se la confezione è danneggiata / Егер пакет бұзылған болса, пайдаланба / पैکیزیگا 손상된 경우 사용 금지 / Jei pakuotė pažeista, nenaudoti / Nelietoti, ja ierakojums bojāts / Niet gebruiken indien de verpakking beschadigd is / Má ikke brukes hvis pakke er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Não usar se a embalagem estiver danificada / A nu se folosi dacă pachetul este deteriorat / He ismolyzovat pri porvrennenii upakovki / Nepoužívejte, ak je obal poškozený / Ne koristite ako je pakovanje oštećeno / Använd ej om förpackningen är skadad / Ambalaj hasar görmüşse kullanmayın / He використувувати за пошкодженої упаковки / 如果包装破損, 請勿使用

Keep away from heat / Пазете от топлина / Nevystavujte přilísnému teplu / Má ikke udsættes for varme / Vor Wärme schützen / Κρατήστε το μακριά από τη θερμότητα / Mantener alejado de fuentes de calor / Hoida eemal valgusest / Protéger de la chaleur / Držati dalje od izvora topline / Övja a melegtől / Tenere lontano dal calore / Салкын жерде сакта / 열을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargat no karstuma / Beschermen tegen warmte / Má ikke utsettes for varme / Przechowywać z dala od źródeł ciepła / Manter ao abrigo do calor / A se feri de căldură / He harrevat / Uchovávať mimo zdroja tepla / Držite dalje od toplote / Får ej utsättas för värme / Isidan uzak tutun / Берегти від дії тепла / 请远离热源

Cut / Срежете / Odstfihňte / Klip / Schneiden / Κόψτε / Cortar / Lögata / Découper / Reži / Vágja ki / Tagliare / Kecişiz / 잘라내기 / Kirpti / Nogriezt / Knippen / Kutt / Odciać / Cortar / Decupati / Отрезать / Odsthנית / Iseći / Klipp / Kesme / Pozpiżati / 剪下

Collection date / Дата на събиране / Datum odběru / Opsamlingsdato / Entnahmdatum / Ημερομηνία συλλογής / Fecha de recogida / Kogumiskuupäev / Date de prélèvement / Dani prikupljanja / Mintavétel dátuma / Data di raccolta / Жинаған тізбекүні / 수집 날짜 / Paėmimo data / Savākšanas datums / Verzameldatum / Dato prøvetaking / Data pobrania / Data de colecta / Data colectării / Дата сбора / Dátum odboru / Datum prikupljanja / Uppsamlingsdatum / Toplama tarihi / Дата забору / 采集日期

$\mu\text{L}$ /test /  $\mu\text{L}$ /тест /  $\mu\text{L}$ /Test /  $\mu\text{L}$ /εξέταση /  $\mu\text{L}$ /prueba /  $\mu\text{L}$ /teszt /  $\mu\text{L}$ /테스트 /  $\mu\text{L}$ /тест /  $\mu\text{L}$ /tyrimas /  $\mu\text{L}$ /pārbaude /  $\mu\text{L}$ /teste /  $\mu\text{L}$ /анализ /  $\mu\text{L}$ /检测

Keep away from light / Пазете от светлина / Nevystavujte světlu / Má ikke udsættes for lys / Vor Licht schützen / Κρατήστε το μακριά από το φως / Mantener alejado de la luz / Hoida eemal valgusest / Conserver à l'abri de la lumière / Držati dalje od svetla / Fény nem érheti / Tenere al riparo dalla luce / Қараңғыланған жерде ұста / 빛을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargat no gaismas / Niet blootstellen aan zonlicht / Má ikke utsettes for lys / Przechowywać z dala od źródeł światła / Manter ao abrigo da luz / Ferij de lumină / Хранити в темноті / Uchovávať mimo dosahu svetla / Držite dalje od svetlosti / Får ej utsättas för ljus / Işıktan uzak tutun / Берегти від дії світла / 请远离光线

Hydrogen gas generated / Образован е водород газ / Možnost úniku plyného vodíku / Frembringer hydrogengas / Wasserstoffgas erzeugt / Δημιουργία αερίου υδρογόνου / Producción de gas de hidrógeno / Vesinikigaasi tekitatud / Produit de l'hydrogène gazeux / Sadržaji hydrogen vodik / Hidrogeén gázt fejleszt / Produzione di gas idrogeno / Газетекес сүрети пайда болды / 수소 가스 생성됨 / Iiskirvia vandenilio dujas / Rodas ūdeņraža / Waterstofgas gegenereerd / Hydrogengass generert / Powoduje powstawanie wodoru / Produção de gás de hidrogénio / Generare gaz de hidrogen / Выделение водорода / Vyrábene použitím vodíka / Oslobada se vodonik / Genererad välgas / Açığa çıkan hidrojen gazı / Реакция з виділенням водню / 会产生氢气

Patient ID number / ID номер на пациента / ID pacienta / Patientens ID-nummer / Patienten-ID / Αριθμός αναγνώρισης ασθενούς / Número de ID del paciente / Patsiendi ID / No d'identification du patient / Identifikacijski broj pacijenta / Beteg azonosító száma / Numero ID paziente / Пациенттин идентификациялык номери / 환자 ID 번호 / Paciento identifikavimo numeris / Pacienta ID numurs / Identificationnummer van de patiënt / Pasientens ID-nummer / Numer ID pacienta / Número da ID do doente / Număr ID pacient / Идентификационный номер пациента / Identifikačné číslo pacienta / ID broj pacijenta / Patientnummer / Hasta kimlik numarası / Идентификатор пациента / 患者标识号

Fragile, Handle with Care / Чупливо, Работете с необходимото внимание. / Křehké. Při manipulaci postupujte opatrně. / Forsigtig, kan gå i stykker. / Zerbrechlich, vorsichtig handhaben. / Ευθραστο. Χειριστείτε το με προσοχή. / Frágil. Manipular con cuidado. / Örn, kásitsege ettevaatlikult. / Fragile. Manipuler avec précaution. / Lomljivo, rukujte pažljivo. / Törékeny! Óvatosan kezelendő. / Fragile, maneggiare con cura. / Сынгыш, абайлап пайдаланыңыз. / 조심 깨지기 쉬운 처리 / Trapu, elkítés atsargiai. / Trausis; rikoties uzmanigi / Breekbaar, voorzichtig behandelen. / Ømtålig, håndter forsigtig. / Krucha zawartość, przenosić ostrożnie. / Frágil, Manuseie com Cuidado. / Frágil, manipulați cu atenție. / Хрупкое! Обращаться с осторожностью. / Křehké, vyžaduje sa opatrná manipulácia. / Lomljivo - rukujte pažljivo. / Bräckligt. Hantera försiktigt. / Kolay Kırılır, Dikkatli Taşın. / Тендітпа, зерпатитса з обережністю / 易碎，小心轻放

R<sub>x</sub> Only

This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

[bd.com/e-labeling](http://bd.com/e-labeling)

|                            |                      |                        |
|----------------------------|----------------------|------------------------|
| <b>Europe, CH, GB, NO:</b> |                      | <b>+800 135 79 135</b> |
| <b>International:</b>      |                      | <b>+31 20 794 7071</b> |
| AR                         | +800 135 79 135      | LT 8800 30728          |
| AU                         | +800 135 79 135      | MT +31 20 796 5693     |
| BR                         | 0800 591 1055        | NZ +800 135 79 135     |
| CA                         | +1 855 805 8539      | RO 0800 895 084        |
| CO                         | +800 135 79 135      | RU +800 135 79 135     |
| EE                         | 0800 0100567         | SG 800 101 3366        |
| GR                         | 00800 161 22015 7799 | SK 0800 606 287        |
| HR                         | 0800 804 804         | TR 00800 142 064 866   |
| IL                         | +800 135 79 135      | US +1 855 236 0910     |
| IS                         | 800 8996             | UY +800 135 79 135     |
| LI                         | +31 20 796 5692      | VN 122 80297           |

Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, Maryland 21152 USA

Benex Limited  
Pottery Road, Dun Laoghaire  
Co. Dublin, Ireland

BD, the BD logo, Synasys, and Veritor are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2020 BD. All rights reserved.

# **Insert ends here**

DO NOT PRINT THIS PAGE

## **CLSI document should not be included in print file to vendor.**

**NOTES:**

The CLSI document is to be updated at the same time as the package insert and should be released as a separate PDF document to the web.

The package insert should also be released as a separate PDF document to the web.

The package insert and CLSI document should be released as one document to SAP.

Final art provided to Buyer/print vendor should not include this page or the CLSI portion of this document.

Proofreader: please proof the CLSI document along with the package insert.

**Procedure\* BD Veritor™ System For Rapid Detection of Flu A+B**

Laboratory kit configured for testing liquid nasopharyngeal wash, aspirate and swab in transport media samples.

| <b>Prepared by</b> | <b>Date Adopted</b> | <b>Supersedes Procedure #</b> |
|--------------------|---------------------|-------------------------------|
|                    |                     |                               |

| <b>Review Date</b> | <b>Revision Date</b> | <b>Signature</b> |
|--------------------|----------------------|------------------|
|                    |                      |                  |
|                    |                      |                  |
|                    |                      |                  |
|                    |                      |                  |

| <b>Distributed to</b> | <b># of Copies</b> | <b>Distributed to</b> | <b># of Copies</b> |
|-----------------------|--------------------|-----------------------|--------------------|
|                       |                    |                       |                    |
|                       |                    |                       |                    |
|                       |                    |                       |                    |
|                       |                    |                       |                    |

---

**\*Any modifications to this document are the sole responsibility of the facility. This “Sample Procedure” is not intended as a substitute for your facility procedure manual, instrument manual, or reagent labeling/package insert. This “Sample Procedure” is intended as a model for use by your facility to meet the needs of your laboratory.**

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Diagnostic tests available for influenza include rapid immunoassay, immunofluorescence assay, polymerase chain reaction (PCR), serology, and viral culture.<sup>4–11</sup> Immunofluorescence assays entail staining of specimens immobilized on microscope slides using fluorescent-labeled antibodies for observation by fluorescence microscopy.<sup>6,12,13</sup> Culture methods employ initial viral isolation in cell culture, followed by hemadsorption inhibition, immunofluorescence, or neutralization assays to confirm the presence of the influenza virus.<sup>13–15</sup>

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

**NOTE: Do not use tips from any other product, including other products from BD or other manufacturers.**

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


## Step 6B: Starting Walk Away mode:

- Turn the Analyzer on by pressing the blue power button once.
- When the display window reads: “INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE”
  - **Double-click** the blue power button.

## Step 7B: Adding the specimen

- When the display window reads “ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY”:
  - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System Flu A+B device sample well).
  - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System Flu A+B device.

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
- “DO NOT DISTURB TEST IN PROGRESS” appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- A countdown timer in the display window shows the remaining analysis time.

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
- Immediately insert the test device into the slot on the right side of the Analyzer.
- The test device must remain horizontal to prevent spilling the specimen out of the sample well.
- "DO NOT DISTURB TEST IN PROGRESS" appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
  - A countdown timer in the display window shows the remaining analysis time.

Do not touch the Analyzer or remove the test device during this process.

Doing so will abort the assay analysis.

- When analysis is complete, a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed. **If a printer is not connected, record the result before removing the assay device.**

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**Note:** The BD Veritor™ System for Rapid Detection of RSV (Cat. No. 256042) is required for this procedure in addition to the BD Veritor™ System for Rapid Detection of Flu A+B (Cat. No. 256041).

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## QUALITY CONTROL:

**To utilize the Analyzer's QC documentation capability, specimen barcode scanning must be enabled on an Analyzer equipped with the BD Veritor InfoScan module. Please refer to the *Analyzer Instructions for Use*, section 4, to choose or change this configuration.**

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**BD** recommends external controls be run once for:

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**Analytical Studies**

**Analytical Sensitivity (Limit of Detection)**

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## REFERENCES

1. Simonsen L., Fukuda K, Schonberger LB, Cox NJ. Impact of influenza epidemics on hospitalizations. *J. Infect. Dis.* 2000; 181:831–7.
2. Thompson WW, Shay DK, Weintraub E, et al. Mortality associated with influenza and respiratory syncytial virus in the United States. *JAMA* 2003; 289: 179–86.
3. Treanor, J.J., Hayden, F.G., Vrooman, P.S., et al. 2000. Efficacy and safety of the oral neuraminidase inhibitor oseltamivir in treating acute influenza: a randomized controlled trial. *JAMA.* 283:1016–1024.
4. Kaiser, L., Couch, R.B., Galasso, G.J., Glezen, W.P., Webster, R.G., Wright, P.F., and Hayden, F.G. 1999. First international symposium on influenza and other respiratory viruses: summary and overview Kapalua, Maui, Hawaii, December 4–6, 1998. *Antiviral Res.*, 42:149–176.
5. Cox, N.J., and Bender, C.A. 1995. The molecular epidemiology of influenza viruses. *Virology*, 6:359–370.
6. Todd, S.J., Minnich, L., and Waner, J.L. 1995. Comparison of rapid immunofluorescence procedure with TestPack RSV and Directigen Flu A for diagnosis of respiratory syncytial virus and influenza A virus. *J. Clin. Microbiol.*33:1650–1651.
7. Harris, P.O. 1989. Clinical relevance and efficient detection of seven major respiratory viruses. *ACL.* p. 15–19.
8. McElhaney, J.E., Gravenstein, S., Krause, P., Hooton, J.W., Upshaw, C.M., and Drinka, P. 1998. Assessment of markers of the cell-mediated immune response after influenza virus infection in frail older adults. *Clin. Diag. Lab. Immunol.* 5:840–844.
9. Fan, J., Henrickson, K.J., and Savatski, L.L. 1998. Rapid simultaneous diagnosis of infections with respiratory syncytial viruses A and B, influenza viruses A and B, and human parainfluenza virus types 1, 2, and 3 by multiplex quantitative reverse transcription-polymerase chain reaction-hybridization assay (hexaplex). *Clin. Infect. Disease* 26:1397–1402.
10. Wright, K.E., Wilson, G.A.R., Novosad, D., Dimock, C., Tan, D., and Weber, J.M. 1995. Typing and subtyping of influenza viruses in clinical samples by PCR. *J. Clin. Microbiol.*33:1180–1184.
11. Kendal, A.P. 1985. Influenza Viruses. p. 341–357. *Laboratory Diagnosis of Viral Infections*, *In* H. Lennette, (ed.) Marcel Dekker, Inc., New York.
12. McQuillen, J., Madeley, C.R., and Kendal, A.P. 1985. Monoclonal antibodies for the rapid diagnosis of influenza A and B virus infections by immunofluorescence. *Lancet*.ii: 911–914.
13. Guenther, S.H., and Linnemann, C.C., Jr. 1988. Indirect immunofluorescence assay for rapid diagnosis of influenza virus. *Laboratory Medicine.*19: 581–583.
14. Minnick, L.L., and Ray, C.G. 1986. Early testing of cell cultures for detection of hemadsorbing viruses. *J. Clin. Microbiol.*25: 421–422.
15. Schmidt, N.J., Ota, M., Gallo, D., and Fox, V.L. 1982. Monoclonal antibodies for rapid, strain specific identification of influenza virus isolates. *J. Clin. Microbiol.*16: 763–765.
16. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
17. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.*17:53–80.
18. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S Government Printing Office, Washington, D.C.
19. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EECP). *Office Journal L262*, 17/10/2000, p.021–0045.

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