

BD Control Set for the BD ProbeTec™ *Chlamydia trachomatis*/Neisseria gonorrhoeae/Trichomonas vaginalis (CT/GC/TV) Q^x Amplified DNA Assays

CE
2797

IVD

R_x Only



8089071(06)

2021-10

English

REF 441925

INTENDED USE

The Control Set for the BD ProbeTec™ *Chlamydia trachomatis* / *Neisseria gonorrhoeae* / *Trichomonas vaginalis* (CT/GC/TV) Q^x Amplified DNA Assays contains Positive and Negative Controls that are intended for the qualitative Quality Control of the automated BD ProbeTec™ *Chlamydia trachomatis* Q^x Amplified DNA Assay, the BD ProbeTec™ *Neisseria gonorrhoeae* Q^x Amplified DNA Assay, and/or the BD ProbeTec™ *Trichomonas vaginalis* Q^x Amplified DNA Assay when tested with the BD Viper™ System in Extracted Mode.

SUMMARY AND EXPLANATION

Quality control must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

The BD ProbeTec™ *Chlamydia trachomatis*/Neisseria gonorrhoeae/Trichomonas vaginalis (CT/GC/TV) Q^x Positive and Negative Controls must be included with each BD Viper™ run. Controls must be positioned according to the BD Viper™ Instrument User's Manual. The CT/GC/TV Q^x Positive Control will monitor for substantial reagent failure only. The CT/GC/TV Q^x Negative Control monitors for reagent and/or environmental contamination.

The CT/GC/TV Q^x Positive Control contains cloned CT/GC/TV target regions. This control may be used for internal quality control or users may develop their own internal quality control material.¹ Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organizations. Refer to CLSI C24-A3 for additional guidance on appropriate internal quality control testing practices.¹ The Positive Control contains approximately 2,400 copies per mL of each pCTB4 and pGCint3 linearized plasmids, and approximately 4,000 copies per mL of TVAP651 linearized plasmid.

REAGENTS

Materials Provided

Each Control Set for the BD ProbeTec™ *Chlamydia trachomatis*/Neisseria gonorrhoeae/Trichomonas vaginalis (CT/GC/TV) Q^x Amplified DNA Assays contains: 24 CT/GC/TV Q^x Positive Control Tubes containing approximately 2,400 copies per mL of each pCTB4 and pGCint3 linearized plasmids, and approximately 4,000 copies of TVAP651 linearized plasmid in carrier nucleic acid, and 24 CT/GC/TV Q^x Negative Control Tubes containing carrier nucleic acid alone. The concentration of the pCTB4, pGCint3, and the TVAP651 plasmids are determined by UV spectrophotometry.

Materials Required But Not Provided

BD ProbeTec™ *Chlamydia trachomatis*/Neisseria gonorrhoeae/Trichomonas vaginalis (CT/GC/TV) Q^x Amplified DNA Assays and accessories, Nitrile gloves.

Storage and Handling Requirements

Reagents may be stored at 2–33 °C. Do not freeze.

WARNINGS AND PRECAUTIONS

1. For in vitro diagnostic use. For Use by Trained Laboratory Personnel.
2. The Control Set for the BD ProbeTec™ *Chlamydia trachomatis*/Neisseria gonorrhoeae/Trichomonas vaginalis (CT/GC/TV) Q^x Amplified DNA Assays is used to evaluate the performance of the BD ProbeTec™ *Chlamydia trachomatis*/Neisseria gonorrhoeae/Trichomonas vaginalis (CT/GC/TV) Q^x Amplified DNA Assays and should not be used with other methods.
3. For additional specific warnings, cautions, and notes specific to the BD Viper™, consult the BD Viper™ Instrument User's Manual.
4. Do not re-hydrate the controls prior to loading in the BD Viper™ Lysing Rack.
5. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.

Dispose of all used reagents and any other contaminated disposable materials following procedures for infectious or potentially infectious waste. It is the responsibility of each laboratory to handle solid and liquid waste according to their nature and degree of hazard and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

QUALITY CONTROL PREPARATION

The CT/GC/TV Q^x Positive and Negative Controls (i.e., assay controls) do not require addition of fluid by the user prior to loading in the BD Viper™ Lysing Rack.

Interpretation of Quality Control Results

The CT/GC/TV Q^x Positive Control and the CT/GC/TV Q^x Negative Control must test as positive and negative, respectively, in order to obtain patient results. If controls do not perform as expected, the run is considered invalid and patient results will not be reported by the instrument. If either of the controls does not provide the expected results, repeat the entire run using a new set of controls, new extraction tubes, new extraction reagent trough, new lysis trough and new microwells. If the repeat QC does not provide the expected results, contact BD Technical Services.

Refer to the BD ProbeTec™ *Chlamydia trachomatis* (CT) Q^x Amplified DNA Assay (BD Catalog Number 441126) / BD ProbeTec™ *Neisseria gonorrhoeae* (GC) Q^x Amplified DNA Assay (BD Catalog Number 441124) / BD ProbeTec™ *Trichomonas vaginalis* Q^x Amplified DNA Assay (BD Catalog Number 441917) package inserts for further interpretation of the quality control results.

LIMITATIONS OF THE PROCEDURE

The CT/GC/TV Q^x Control Set may not adequately control for specimen processing. To prepare and test specimen processing controls, refer to the BD ProbeTec™ *Chlamydia trachomatis* (CT) Q^x Amplified DNA Assay (BD Catalog Number 441126) / BD ProbeTec™ *Neisseria gonorrhoeae* (GC) Q^x Amplified DNA Assay (Catalog Number 441124) / BD ProbeTec™ *Trichomonas vaginalis* Q^x Amplified DNA Assay (Catalog Number 441917) package inserts.

REFERENCE

1. Clinical and Laboratory Standards Institute. 2006. Approved Guideline C24-A3. Statistical quality control for quantitative measurement procedures: principles and definitions, 3rd ed. CLSI, Wayne, PA.

Technical Service and Support: In the United States contact BD at 1.800.638.8663 or bd.com.

For regions outside of the United States, contact your local BD representative or bd.com.

EU Only: Users shall report any serious incident related to the device to the Manufacturer and National Competent Authority.

Outside EU: Contact your local BD representative for any incident or inquiry related to this device.

Refer to the Eudamed website: <https://ec.europa.eu/tools/eudamed> for Summary of Safety and Performance.

Change History

| Revision | Date | Change Summary |
|----------|---------|--|
| 06 | 2021-10 | Added CE Notified Body 2797 for IVDR 2017/746. Added Intended Use, Intended User, Serious Incident statement and Safe Disposal statement. Added IVD, R _x Only, eIFU with URL and Do not Use if Package is Damaged symbol. Updated Australian and New Zealand Sponsor addresses. Updated EC REP address.. Updated Technical Information and Eudamed link. Updated Symbols Glossary. Added CH REP symbol with address. |

SYMBOLS GLOSSARY [L006715(06) 2021-08]

Some symbols listed below may not apply to this product.

US Customers only: For symbol glossary, refer to bd.com/symbols-glossary

| Symbol | Meaning | Symbol | Meaning |
|--------|---|--------|--|
| | Manufacturer | | Patient number |
| | Authorized representative in the European Community | | This way up |
| | Authorised representative in Switzerland | | Do not stack |
| | Date of manufacture | | Single sterile barrier system |
| | Use-by date | | Contains or presence of phthalate: combination of bis(2-ethylhexyl) phthalate (DEHP) and benzyl butyl phthalate (BBP) |
| | Batch code | | Collect separately Indicates separate collection for waste of electrical and electronic equipment required. |
| | Catalogue number | | CE marking: Signifies European technical conformity |
| | Serial number | | Device for near-patient testing |
| | Sterile | | Device for self-testing |
| | Sterilized using aseptic processing techniques | | This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." |
| | Sterilized using ethylene oxide | | Country of manufacture "CC" shall be replaced by either the two letter or the three letter country code. |
| | Sterilized using irradiation | | Collection time |
| | Sterilized using steam or dry heat | | Cut |
| | Do not re-sterilize | | Peel here |
| | Non-sterile | | Collection date |
| | Do not use if package is damaged and consult <i>instructions for use</i> | | Keep away from light |
| | Sterile fluid path | | Hydrogen gas is generated |
| | Sterile fluid path (ethylene oxide) | | Perforation |
| | Sterile fluid path (irradiation) | | Start panel sequence number |
| | Fragile, handle with care | | End panel sequence number |
| | Keep away from sunlight | | Internal sequence number |
| | Keep dry | | Medical device |
| | Lower limit of temperature | | Contains hazardous substances |
| | Upper limit of temperature | | Ukrainian conformity mark |
| | Temperature limit | | Meets FCC requirements per 21 CFR Part 15 |
| | Humidity limitation | | UL product certification for US and Canada |
| | Biological risks | | Unique device identifier |
| | Do not re-use | | |
| | Consult <i>instructions for use</i> or consult <i>electronic instructions for use</i> | | |
| | Caution | | |
| | Contains or presence of natural rubber latex | | |
| | In vitro diagnostic medical device | | |
| | Negative control | | |
| | Positive control | | |
| | Contains sufficient for <n> tests | | |
| | For IVD performance evaluation only | | |
| | Non-pyrogenic | | |



Becton, Dickinson and Company
7 Loveton Circle
Sparks, Maryland 21152 USA



Becton Dickinson Ireland Ltd.
Donore Road, Drogheda
Co. Louth, A92 YW26, Ireland



BD Switzerland Sàrl
Terre Bonne Park – A4
Route de Crassier 17
1262 Eysins, Switzerland

Australian and New Zealand Sponsors:
Becton Dickinson Pty Ltd.
66 Waterloo Road
Macquarie Park NSW 2113, Australia
Becton Dickinson Limited
14B George Bourke Drive
Mt. Wellington Auckland 1060, New Zealand

BD, the BD Logo, ProbeTec, and Viper are trademarks of Becton, Dickinson and company or its affiliates. All other trademarks are the property of their respective owners. © 2021 BD. All rights reserved.