
# Revisions

BALTSO0191 Version 8.0 Template 4  
Inserts

| Rev from | Rev to | JOB #   |
|----------|--------|---------|
| 08       | 09     | 7236-15 |

## NOTES:

- BD Catalog Number: 490500, 490501
- Blank (Sheet) Size: Length: 11 Width: 8.5"
- Number of Pages: 26 Number of Sheets: 13
- Page Size: Length: 11" Width: 8.5" Final Folded Size: 8.5" x 11"
- Ink Colors: No. of Colors: 2 PMS#: 032 Red; Standard Black
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 1 (3-hole punched); Stapled top left corner

- BD Diagnostics inspects per procedure 838-07096-00 upon entry to their plant.
- See specification control no. N/A for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

|              |                                                                    |                                                                                                                                                            |                                                                                                                                                                          |
|--------------|--------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Label Design | REVISED BY<br>By Nancy Carlsen at 12:46 pm, Oct 20, 2015           | COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION. |  <b>BD</b> Becton, Dickinson and Company<br>7 Loveton Circle<br>Sparks, MD 21152 USA |
| Proofer      | PROOFING APPROVED BY<br>By Helaine Jones at 10:19 am, Oct 21, 2015 |                                                                                                                                                            |                                                                                                                                                                          |
| Checked By   | THIRD EYE BY<br>By Sarah Cork at 12:59 pm, Oct 21, 2015            |                                                                                                                                                            |                                                                                                                                                                          |
| Part Number: | 779-05848-10                                                       | Category and Description<br>Package Insert,<br>BD SurePath with ProEx C<br>Immunocytochemical Test Kit                                                     | Sheet: 1 of 27<br>Scale: N/A<br><div style="font-size: 2em; font-weight: bold; text-align: center;">A</div>                                                              |

# SurePath® with ProEx™ C Immunocytochemical Test Kit

|                                                                                                                                                                |                                                                  |    |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------|----|
| REF 490501                                                                                                                                                     | SurePath® with ProEx™ C<br>Immunocytochemical Test Kit           | 75 |
| REF 490500                                                                                                                                                     | SurePath® with ProEx™ C<br>Immunocytochemical Test Counterstains | 75 |
| English: pages 1 – 4    Deutsch: Seiten 9 – 12    Español: páginas 17 – 20<br>Français : pages 5 – 8    Italiano: pagine 13 – 16    Português: páginas 21 – 24 |                                                                  |    |

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyts lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőtől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o representante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Inštrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD

## INTENDED USE

The **SurePath with ProEx C** Immunocytochemical Test is intended for the qualitative evaluation of aberrant S-Phase induction in cervical cytology specimens collected in **SurePath** Preservative Fluid. The test results provide adjunctive information for cervical cytology diagnosis.

## SUMMARY AND EXPLANATION

Minichromosome maintenance (MCM) proteins play an essential role in eukaryotic DNA replication. Each of the MCM proteins has DNA-dependent ATPase motifs in their highly conserved central domain. Levels of MCM proteins generally increase in a variable manner as normal cells progress from G0 into the G1/S phase of the cell cycle. Topoisomerase II alpha (TOP2A) is an essential nuclear enzyme involved in DNA replication and is a target for many anti-cancer drugs used for cancer therapy. Decreased expression of TOP2A is a predominant mechanism of resistance to several chemotherapeutic agents. A significant variation in the range of expression of this protein has been noted in many different tumors. TOP2A is predominant in proliferating cells and is modified in M phase by phosphorylation at specific sites, which is critical for mitotic chromosome condensation and segregation.

**ProEx C** Antibody Cocktail contains mouse monoclonal anti-MCM2 and anti-TOP2A purified from tissue culture supernatant and diluted in buffered saline solution containing protein stabilizers and 0.09% sodium azide.

## PRINCIPLES OF PROCEDURE

The **SurePath with ProEx C** Immunocytochemical Test kit contains an antibody cocktail, detection and counterstain reagents necessary to complete a three-step manual or automated immunocytochemical staining procedure for routinely prepared thinlayer cervical specimens. Following incubation of the sample with a proprietary, mouse antibody cocktail, selective monoclonal antibody binding, indicative of a positive test, is visualized using a unique, ready-to-use, enzyme-linked antibody chromogen system. The enzyme reagent is a secondary goat anti-mouse, horseradish peroxidase conjugate linked to a dextran polymer backbone. Addition of a specific chromogen results in formation of a visible chromogenic product localized at the antigen-antibody binding sites. The specimen is then counterstained with hematoxylin, a bluing agent is applied and the slide is coverslipped. Results are interpreted using a light microscope. A positive result, indicative of high-grade cervical disease, is achieved when the nucleus in cells of interest stain brown.

A reference gallery of potentially positive slides may be created using automated imaging equipment. The gallery then can be reviewed to determine a positive result or negative result.

The **SurePath with ProEx C** Immunocytochemical Test is applicable for both manual and automated staining.

## REAGENTS PROVIDED

The following materials are included, which are sufficient for 75 thinlayer preparations:

| Vial No. | Description                                                                                                                                                               |          |
|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| 1a       | <b>Peroxidase Block:</b> Buffered hydrogen peroxide plus stabilizer and proprietary components                                                                            | 2–8 °C   |
| 1b       | <b>Protein Block:</b> Purified casein plus proprietary combination of proteins in modified PBS with preservative and surfactant                                           | 2–8 °C   |
| 2        | <b>Ab Cocktail:</b> Ready-to-use monoclonal antibody cocktail supplied in Tris buffered solution with Tween 20. Contains stabilizing proteins and anti-microbial reagent. | 2–8 °C   |
| 3a       | <b>Mouse Probe:</b> Binds to mouse monoclonal antibodies                                                                                                                  | 2–8 °C   |
| 3b       | <b>Polymer Reagent:</b> Polymer conjugated with horseradish peroxidase that binds to Mouse Probe Reagent                                                                  | 2–8 °C   |
| 4a       | <b>DAB Substrate Buffer:</b> Substrate buffer used in the preparation of the DAB chromogen                                                                                | 2–8 °C   |
| 4b       | <b>DAB:</b> 3,3'-diaminobenzidine chromogen solution                                                                                                                      | 2–8 °C   |
| 5        | <b>Hematoxylin:</b> Aqueous based Mayer's Hematoxylin                                                                                                                     | 15–30 °C |
| 6        | <b>Bluing Agent:</b> Tris buffered saline with Tween 20 and 0.09% NaN <sub>3</sub>                                                                                        | 15–30 °C |

## MATERIALS AND REAGENTS REQUIRED BUT NOT SUPPLIED

- Absorbent Wipes
- BD SureDetect™** SiHa Cell Line
- Deionized or Distilled Water
- Ethanol (95% and 100%)
- Glass Coverslips
- Gloves
- Humid Chamber
- Light Microscope (10x, 20x [optional], 40x objectives)
- Mounting Media
- Pipettes and Pipette Tips (capable of delivering 20 µL, 200 µL and 1000 µL volumes)
- BD SureDetect** Slide Preparation Buffer 10X (Pretreatment Buffer)
- Staining Jars or Baths
- Timer (capable of 1–60 minute intervals)
- Tris Buffered Saline (TBS)
- Tween 20
- Universal Mouse IgG Negative Control
- Vortexer
- Xylene or Xylene Substitutes
- Steamer/ waterbath

## WARNINGS AND PRECAUTIONS

### 490500 SurePath with ProEx C Immunocytochemical Counterstains

#### Warning

**H302** Harmful if swallowed. **H315** Causes skin irritation. **H319** Causes serious eye irritation.

**P264** Wash thoroughly after handling. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P302+P352** IF ON SKIN: Wash with plenty of soap and water. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

### 490501 SurePath with ProEx C Immunocytochemical Test Kit

#### Danger

**H341** Suspected of causing genetic defects. **H350** May cause cancer.

**P281** Use personal protective equipment as required. **P308+P313** IF exposed or concerned: Get medical advice/attention. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

## PRECAUTIONS

1. For *in vitro* Diagnostic Use.
2. Prepared **SurePath** slides must be placed into the pretreatment buffer as soon as they are prepared. The slides must remain in the pretreatment buffer for at least 1 hour, but not longer than 72 hours prior to immunocytochemistry.
3. Do not allow the slides to dry out at any time during the procedure. Slides that have been allowed to dry out during the procedure may increase background.
4. Specimens and all materials exposed to specimens should be handled as if capable of transmitting infection and disposed of with proper precautions. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents and specimens. If reagents come in contact with sensitive areas, wash with copious amounts of water.
5. Minimize microbial contamination of reagents to avoid nonspecific staining.
6. Incubation times, temperatures or methods other than those specified may give erroneous results.
7. Reagents have been diluted for optimal performance. Further dilution may result in loss of antigen staining.
8. Kit components are lot specific. Do not substitute kit components with different manufactured lot numbers. Lot numbers appear on packaging labels.
9. Do not use the **SurePath** with **ProEx C** Test kit after the expiration date stamped on the package. The user must validate conditions if reagents are stored under any conditions other than those specified in the package insert.
10. There are no obvious signs to indicate degradation of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed, which cannot be explained by variations in laboratory procedures or a problem with the **SurePath** with **ProEx C** Test kit is suspected, contact BD technical support.
11. Wear appropriate Personal Protective Equipment to avoid reagent contact with eyes and skin. Refer to the Material Safety Data Sheet (MSDS) for additional information.
12. The **SurePath** with **ProEx C** Test is intended for use with cervical cytology specimens collected in **SurePath** Preservative Fluid. Compatibility with conventional and monolayer preparations other than **SurePath** has not been evaluated.

## INSTRUCTIONS FOR USE

### Specimen Preparation

Consult the User's Manual for the **BD PrepStain™** Slide Processor for instructions regarding the preparation of slides from residual **SurePath** cervical cytology specimens.

Add 8 mL of **SurePath** Preservative Fluid to the residual sample in the **SurePath** vial (approx. 2 mLs). The diluted sample should be processed on the **BD PrepMate** Instrument using the standard technique and on the **BD PrepStain** Instrument.

Prepared **SurePath** slides must be placed into the **BD SureDetect** Slide Preparation Buffer 10X as soon as they are prepared (Refer to the **BD SureDetect** Slide Preparation Buffer Product Insert for directions on preparation of a working slide solution). The slides must remain in the buffer for at least 1 hour, but not longer than 72 hours prior to immunocytochemistry.

An epitope retrieval procedure must be used for optimal kit performance. This procedure involves soaking prepared slides in a working solution of **BD SureDetect** Slide Preparation Buffer for a minimum of 1 hour at room temperature followed by heating slides in the pretreatment buffer to 95 °C. Slides are held at 95 °C for 15 minutes and are allowed to cool at room temperature for 20 minutes. The use of a calibrated waterbath or steamer capable of maintaining the required temperature is recommended. Laboratories located at higher elevations should determine the best method of maintaining the required temperature. The staining procedure should be initiated immediately following epitope retrieval and cool down. Deviations from the described procedure may adversely affect results.

### Reagent Preparation

Prepare the following reagents prior to staining.

- Tris Buffered Saline with 0.05% Tween 20 (TBST)
  - Prepare TBS according to manufacturer's specifications.
  - If not already present in TBS, add Tween 20 to a final concentration of 0.05%.
- Substrate-Chromogen Solution (DAB) (volume sufficient for 5 slides)
  - Transfer 1 mL of DAB Substrate Buffer (vial 4a) to a test tube.
  - Add one drop (20–30 µL) of DAB Chromogen (vial 4b). Mix thoroughly and apply to slides with a pipette.
  - Prepare Substrate-Chromogen solution fresh daily.
  - Staining quality is not affected by precipitate that may develop in the solution.

## STAINING PROTOCOL

(PERFORMED AT ROOM TEMPERATURE, 20–25 °C)

### Staining Procedural Notes

- Read all these instructions carefully and become familiar with all components prior to use (see Precautions).
- All reagents should be equilibrated to room temperature (20–25 °C) prior to immunostaining.
- All incubations should be performed at room temperature unless noted.
- Do not allow slides to dry out during the staining procedure. Dried cellular preparations may display increased non-specific staining. Protect slides that may be exposed to drafts. Slides should be placed in a humid chamber for prolonged incubations.

### Epitope Retrieval

- Place prepared slides in a coplin jar containing a working solution of **BD SureDetect** Slide Preparation Buffer for a minimum of 1 hour up to a maximum of 72 hours.
- Incubate in a waterbath or steamer for 15 minutes at 95 °C.
- Remove the entire coplin jar with slides from the waterbath or steamer and allow slides to cool in the buffer for 20 minutes.
- Rinse the slides with deionized H<sub>2</sub>O and transfer to a clean coplin jar containing TBST.

### Blocking Reagent

- Tap off excess buffer.
- Load slides into a prepared humidity chamber (filled with water moistened paper towels or gauze).
- Apply 200 µL Peroxidase Block reagent (vial 1a) to cover the cell deposition area.
- Incubate 5 minutes (±1 minute).
- Rinse slides in TBST, 3 changes, 2 minutes each.

### Protein Block

- Tap off excess buffer.
- Load the slides into the prepared humidity chamber (filled with water moistened paper towels or gauze).
- Apply 200 µL of Protein Block (vial 1b) to completely cover cell deposition area.
- Incubate 5 minutes (±1 minute).
- DO NOT RINSE.

### Primary Antibody Cocktail

- Tap off excess Protein Block.
- Load the slides into the prepared humidity chamber (filled with water moistened paper towels or gauze).
- Apply 200 µL **ProEx C Ab Cocktail** (vial 2) to completely cover cell deposition area.
- Incubate 30 minutes at room temperature.
- Rinse each slide individually with TBST using a wash bottle (do not focus the flow directly on the cell deposition area). Load slides into a slide rack.
- Rinse slides in TBST, 3 changes, 2 minutes each.

### Detection Chemistry

- Tap off excess buffer.
- Load slides into prepared humidity chamber (filled with water moistened paper towels or gauze).
- Apply 200 µL Mouse Probe Reagent (vial 3a) to completely cover cell deposition area.
- Incubate 20 minutes (±1 minute).
- Rinse slides in TBST, 3 changes, 2 minutes each.
- Tap off excess buffer.
- Load slides into prepared humidity chamber (filled with water moistened paper towels or gauze).
- Apply 200 µL of Polymer Reagent (vial 3b) to cover cell deposition area.
- Incubate for 20 minutes (±1 minute).
- Rinse slides in TBST bath, 3 changes, 2 minutes each.
- Tap off excess buffer.
- Load the slides into the prepared humidity chamber (filled with water moistened paper towels or gauze).
- Apply 200 µL of Substrate-Chromogen solution (DAB) completely cover cell deposition area.
- Incubate for 5 minutes (±1 minute).
- Rinse slides for 5 minutes in deionized H<sub>2</sub>O.

### Counterstain

- Rinse slides in TBST, 1 change for 2 minutes.
- Load slides into prepared humidity chamber (filled with water moistened paper towels or gauze).
- Apply 200 µL of Hematoxylin counterstain (vial 5) to completely cover cell deposition area.
- Incubate for 1 minute (±10 seconds).
- Rinse slides for 3 minutes in running H<sub>2</sub>O.
- Load slides into prepared humidity chamber (filled with water moistened paper towels or gauze).
- Blue slides by applying 200 µL Bluing Agent (vial 6) for 1 minute (±10 seconds).
- Repeat running water rinse for 1 minute.

### Mounting

- Immerse slides in 95% ethanol, 1 minute or 25 dips.
- Immerse slides in absolute alcohol, 4 changes, 1 minute each or 25 dips.
- Clear with xylene, 3 changes, 1 minute each or 25 dips.
- Coverslip slides with non-aqueous, permanent mounting media using glass coverslips.

### STORAGE

The storage limit of the **SurePath with ProEx C** Immunocytochemical Test Kit (Antibody and Detection reagents) is up to 12 months from date of manufacture at 2–8 °C.

The storage limit of the **SurePath with ProEx C** Immunocytochemical Test Counterstains (Hematoxylin and Bluing Reagent) is up to 12 months from date of manufacture at 15–30 °C.

Once opened, reagents are stable for thirty (30) days when stored at recommended temperatures. Do not use reagents beyond the expiration date printed on the label.

### QUALITY CONTROL

Variability in results is often derived from differences in specimen handling which deviates from recommended test procedures. Consult the proposed quality control guidelines of the Clinical and Laboratory Standards Institute, *Quality Assurance for Immunocytochemistry* for additional information.

**BD SureDetect** SiHa Cell Control is available as a positive control from BD Diagnostics. Each vial contains a cervical cancer cell line, which is processed in a similar manner as the clinical specimens. A universal mouse IgG negative control can be used as a negative control. A positive and negative control slide should be included in each staining procedure. The staining results should be used as an indication of the validity of the staining run.

### INTERPRETATION OF STAINING

**Patient and Control Specimens:** A cytotechnologist or pathologist should evaluate the stained slides using a light microscope. Cells can be reviewed manually or stored in an electronic image gallery derived from a light microscope.

**Control Slides:** The positive and negative control slides should be examined prior to the review of patient specimens to ascertain that all reagents functioned properly. The presence of a brown reaction product (3,3'-diaminobenzidine tetrahydrochloride, DAB) in the nuclei of the SiHa Cell Control slide stained with the **SurePath with ProEx C** Test is indicative of positive reactivity. The universal mouse IgG negative control slide stained with the **SurePath with ProEx C** Test should have no brown staining nuclei and should show staining only from the hematoxylin counterstain.

**Scoring of slides is a 4-step process.**

#### Step 1: Is the specimen adequate?

The Bethesda System (TBS) for Reporting Cervical Cytology (2nd Ed.) states, "An adequate liquid-based preparation should have an estimated minimum of at least 5000 well-visualized/well-preserved squamous cells." The same criteria should be applied to **SurePath with ProEx C** slides. However, as with a routine Pap preparation, any specimen with abnormal cells, which is exhibiting a positive molecular reaction, is, by definition, satisfactory for evaluation. If the answer at this step is "yes," proceed to the next step, if the answer is "no" the result is **Unsatisfactory for Evaluation**.

#### Step 2: Is there moderate to intense brown nuclear staining in epithelial cells?

To answer "yes" to this step look for brown staining that is easily visualized. If just a faint amount, or "blush," of brown is seen this is not enough to warrant a rendering of positive. If no brown nuclear stain is seen, the result is reported as **Negative**. If adequate brown stain is visualized, proceed to the next step.

#### Step 3: Is the cell with nuclear brown staining a squamous or glandular cell?

If the answer is yes, proceed to the next step. If the answer is no, this is a **Negative** test result.

#### Step 4: Is the cell ≥ ASC (atypical squamous cell) or AGC (atypical glandular cell)?

Use the same morphological criteria outlined in The Bethesda System for Reporting Cervical Cytology (2nd Ed.) to determine if the squamous cell containing the brown nucleus is ≥ ASC (atypical squamous cells). If the cell is considered ≥ASC (or ≥AGC) then this is a **Positive** test result. This includes ASC-US, ASC-H, LSIL, HSIL, and cancer. If the cell is glandular in appearance the TBS criteria for determining if a cell is ≥ AGC (atypical glandular cells) applies. This includes endocervical AGC, endometrial AGC, AIS, and adenocarcinoma. If the cell in question is consistent with NILM (negative for intraepithelial lesion or malignancy) this is a **Negative** test result.

## REFERENCES

1. Wright TC, Schiffman M, Solomon D, Cox JT, Garcia F, Goldie S, Hatch K, Noller KL, Roach N, Runowicz C and Saslow D. Interim Guidance for the Use of Human Papillomavirus DNA Testing as an Adjunct to Cervical Cytology for Screening. *Obstet. Gynecol.* 2004. 102:304-309.
2. Ishimi Y, Okayasu I, Kato C, Kwon H, Kimura H, Yamada K and Song S. Enhanced expression of MCM proteins in cancer cells derived from uterine cervix. *Eur. J. Biochem* 2003. 270:1089-1101.
3. Lei M and Tye BK. Initiating DNA synthesis: from recruiting to activating the MCM Complex. *J Cell Science* 2001. Vol. 114 (8): 1447-1454.
4. Solomon D, Nayar R. The Bethesda System for Reporting Cervical Cytology. Definitions, Criteria, and Explanatory Notes. 2nd Edition. 2004.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or [www.bd.com/ds](http://www.bd.com/ds).

## SurePath with ProEx C Immunocytochemical Test Kit

|            |                                                                |    |
|------------|----------------------------------------------------------------|----|
| REF 490501 | SurePath with ProEx C<br>Immunocytochemical Test Kit           | 75 |
| REF 490500 | SurePath with ProEx C<br>Immunocytochemical Test Counterstains | 75 |

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TriPath Imaging, Inc.  
780 Plantation Drive  
Burlington, NC 27215 USA

Benex Limited  
Pottery Road, Dun Laoghaire  
Co. Dublin, Ireland

**Australian Sponsor:**

Becton Dickinson Pty Ltd.  
4 Research Park Drive  
Macquarie University Research Park  
North Ryde, NSW 2113 Australia

Developed with technology from Millennium Pharmaceuticals, Inc.

**MILLENNIUM** and **MILLENNIUM** are trademarks of  
Millennium Pharmaceuticals, Inc.

Millennium Pharmaceuticals, Inc.  
40 Landsdowne Street  
Cambridge, MA 02139 USA  
[www.millennium.com](http://www.millennium.com)

Detection Reagents supplied by

BioCare Medical  
4040 Pike Ln.  
Concord, CA 94520

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2015