
Revisions

BALTSO0191 Version 10.0 Template 4
Inserts

Rev from	Rev to	JOB #
N/A	01	8754-16

NOTES:

- BD Catalog Number: 491266, 491435
- Blank (Sheet) Size: Length: 11" Width: 8.5"
- Number of Pages: 20 Number of Sheets: 10
- Page Size: Length: 11" Width: 8.5" Final Folded Size: 8.5" x 11"
- Ink Colors: No. of Colors: 2 PMS#: 032 Red, Standard Black
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 1 (3-hole punched); Stapled top left corner

- BD Diagnostics inspects per procedure 838-07096-00 upon entry to their plant.
- Vendor Printed ☒ Online/In House Printed ☐ Web ☐
- See specication control no. N/A for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Sonia Thompson at 9:23 am, Oct 13, 2016	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY By Natalie Morio at 2:08 pm, Oct 13, 2016		
Checked By	THIRD EYE BY By Sarah Henderson at 3:27 pm, Oct 16, 2016		
Part Number: 500017024	Category and Description Package Insert, BD SurePath Manual Method	Sheet: 1 of 21	A
		Scale: N/A	

REF 491266	BD SurePath™ Manual Method Kit	Σ	480
REF 491435	BD SurePath™ Manual Method Kit – Japan	Σ	480
English: pages 1 – 3 Deutsch: Seiten 7 – 9 Español: páginas 13 – 15 Français: pages 4 – 6 Italiano: pagine 10 – 12 Português: páginas 16 – 18			

Contact your local BD representative for instructions. / **Свържете се с местния представител на BD за инструкции.** / **Pokyny vám poskytne místní zástupce společnosti BD.** / **Kontakt den lokale BD repræsentant for at få instruktioner.** / **Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες.** / **Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga.** / **Ota yhteyttä lähimpään BD:n edustajaan ohjeiden saamiseksi.** / **Kontaktiraj lokalnog predstavnika BD za upute.** / **A használati utasítást kérje a BD helyi képviselőtől.** / **Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз.** / **Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi.** / **Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo.** / **Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies.** / **Kontakt din lokale BD-representant for mer informasjon.** / **Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD.** / **Contacte o representante local da BD para instruções.** / **Pentru instrucțiuni, contactați reprezentantul local BD.** / **Для получения указаний обратитесь к местному представителю компании BD.** / **Inštrukcie získate u miestneho zástupcu spoločnosti BD.** / **Obratite se svom lokalnom predstavniku kompanije BD za uputstva.** / **Kontaktka narmaste BD-representant för anvisningar.** / **Talimatlar için yerel BD temsilcinize temasa geçin.** / **За інструкціями зверніться до місцевого представника компанії BD.**

INTENDED USE

The **BD SurePath™** Manual Method is a method for producing liquid-based cell preparations (LBPs). The **BD SurePath** Manual Method is intended as a replacement for the conventional Pap smear preparation method for use in cervical cancer screening.

BD SurePath Preservative Fluid is an appropriate collection and transportation medium for gynecologic specimens tested with **BD ProbeTec™ Chlamydia trachomatis** (CT) Qx and **Neisseria gonorrhoeae** (GC) Qx Amplified DNA Assays. Refer to the assay package inserts for instructions on using **BD SurePath** Preservative Fluid to prepare specimens for use with these assays.

SUMMARY AND EXPLANATION

Cervical cytology screening by the Papanicolaou (Pap) method involves the microscopic examination of cell samples that have been taken primarily from the ecto- and endocervix, smeared on glass slides and stained using the Pap procedure.¹⁻³ Cervical cytology screening with the Pap smear has decreased the mortality rates of invasive cervical carcinoma by 50 to 70 percent.⁴ Because cervical cytology is a screening test, abnormal findings must be confirmed histologically.

Specimen collection and preparation are extremely important to accuracy in Pap smears. Randomization or uniform sub-sampling is essential for complete accuracy. The conventional Pap smear technique does not provide for mixing of the sample prior to slide preparation. Due to the entanglement of cells in mucus on the sampling device, the cells actually transferred to the slide may not be representative of the total population collected. The cells are transferred to the slide in relation to where they happen to be on the sampling device. Many cells are left on the device.⁵

The non-homogeneity of a typical cervical sample can make conventional smears difficult to prepare, screen and interpret. Large areas of the conventional slide are often covered with debris, inflammatory cells and sheets of epithelial cells that can obscure valuable diagnostic material. In addition, if the smear is not fixed immediately after preparation, cellular morphology can be distorted as the smear dries (air-drying artifact).

The **BD SurePath** Manual Method is a method for converting a liquid suspension of a cervical sample into a consistently stained, homogeneous **BD SurePath** Liquid-based Pap Test slide while maintaining diagnostic cell clusters.⁶⁻⁹ The process includes cell preservation, randomization, enrichment of diagnostic material, pipetting, and sedimentation to create a cellular preparation. The result of the preparation process is a **BD SurePath** Liquid-based Pap Test slide for use in routine cytology screening and categorization as defined by the Bethesda System.¹⁰

PRINCIPLES OF THE PROCEDURE

The **BD SurePath** Manual Method is a procedure for the preparation of LBPs of cervical cells. Gynecologic specimens are collected by qualified medical personnel using broom-type sampling devices (e.g. Rovers® Cervex-Brush®, Rovers Medical Devices B.V., Oss - The Netherlands) or combination plastic spatula and endocervical brush devices (e.g. Cytobrush Plus® GT and Pap Perfect® Plastic Spatula, CooperSurgical, Inc.) with detachable heads. The head of the brush is removed from the handle and placed into a vial of **BD SurePath** Preservative Fluid. The vial is capped, labeled and sent with appropriate paperwork to the laboratory for processing.

In the laboratory, the preserved sample is mixed by vortexing, and transferred into a tube containing **BD** Density Reagent. An enrichment step, consisting of centrifugal sedimentation through **BD** Density Reagent, partially removes non-diagnostic debris and excess inflammatory cells from the sample. After centrifugation, the tube containing the enriched cellular component is reconstituted with buffered ionized water and the cellular material is re-suspended with a pipettor using an aspirate/dispense sequence. The sample material is then transferred to a **BD** Settling Chamber mounted on a **BD SurePath** PreCoat Slide. Gravity sedimentation occurs during a short incubation. Excess material is decanted. The **BD SurePath** Liquid-based Pap Test slide is stained cleared and coverslipped with the cells presented in a circle, 13 mm in diameter. The **BD SurePath** Liquid-based Pap Test slide is examined by trained cytotechnologists and pathologists in conjunction with other relevant patient information.

LIMITATIONS OF THE PROCEDURE

- Gynecologic specimens for preparation using the **BD SurePath** Manual Method should be collected using a broom-type sampling device or combination plastic spatula and endocervical brush device with detachable head(s) according to the standard collection procedure provided by the manufacturer. Wooden spatulas should not be used. Endocervical brush/plastic spatula combinations that are not detachable should not be used.
- The production and evaluation of **BD SurePath** Liquid-based Pap Test preparations should be performed only by personnel who have been trained by BD Diagnostics or others authorized by BD Diagnostics to provide such training.
- Proper performance of the device requires the use of only those supplies supported by BD Diagnostics, or recommended by BD Diagnostics. Used supplies and product should be disposed of properly in accordance with institutional and governmental regulations.
- All supplies are intended for single use only and cannot be reused.
- A volume of 8.0 ± 0.5 mL of the sample collected in the **BD SurePath** Collection Vial is required for the **BD SurePath** LBC Test process.

WARNINGS AND PRECAUTIONS

Cytologic specimens may contain infectious agents. Wear suitable protective clothing, gloves, and eye/face protection. Follow appropriate biohazard precautions when handling samples.

BD SurePath Preservative Fluid contains an aqueous solution of denatured ethanol. The mixture contains small amounts of methanol and isopropanol. Do not ingest.

Warning

H226 Flammable Liquid and Vapour.

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. **P233** Keep container tightly closed. **P240** Ground/bond container and receiving equipment. **P241** Use explosion-proof electrical/ventilating/lighting/equipment. **P242** Use only non-sparking tools. **P243** Take precautionary measures against static discharge. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303+P361+P353** IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. **P370+P378** In case of fire: Use for extinction: CO₂, powder or water spray. **P403+P235** Store in a well-ventilated place. Keep cool. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

BD Density Reagent contains sodium azide. Do not ingest.

Warning

H302 Harmful if swallowed.

P264 Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P330** Rinse mouth. **P501** Dispose of contents/ container in accordance with local/regional/national/international regulations.

PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. For professional use only.
3. Good laboratory practices should be followed and all procedures for use of the **BD SurePath** Manual Method should be strictly observed.
4. Reagents should be stored at room temperature (15–30 °C) and used prior to their expiration dates to assure proper performance.
5. Microbial contamination of reagents may give incorrect results.
6. Substitution of other than **BD SurePath** PreCoat Slides may result in less than optimal results.
7. Avoid splashing or generating aerosols. Use appropriate hand, eye and clothing protection.
8. **BD SurePath** Preservative Fluid was tested for antimicrobial effectiveness against: *Escherichia coli*, *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Candida albicans*, *Mycobacterium tuberculosis* and *Aspergillus niger* and found to be effective. **BD SurePath** Preservative samples inoculated with 10⁶ CFU/mL of each species yielded no growth after 14 days (28 days for *Mycobacterium tuberculosis*) of incubation under standard conditions. However, universal precautions for safe handling of biological fluids should be practiced at all times.

OPTIONAL ALIQUOT REMOVAL

- Sufficient volume is available in the **BD SurePath** Collection Vial to allow removal of up to 0.5 mL of homogeneous mixture of cells and fluid for ancillary testing, prior to the **BD SurePath** Liquid-based Pap Test while still allowing sufficient volume for Pap testing.
- While there is no evidence that removal of an aliquot from the **BD SurePath** Collection Vial affects the quality of the specimen for cytology testing, rare instances of misallocation of pertinent diagnostic material may occur during this process. Healthcare providers may need to acquire a new specimen if the results do not correlate with the clinical history of the patient. Furthermore, cytology addresses different clinical questions than sexually transmitted disease (STD) testing; therefore, aliquot removal may not be suitable for all clinical situations. If necessary, a separate specimen may be collected for STD testing rather than taking an aliquot from the **BD SurePath** Collection Vial.
- Aliquot removal from low-cellularity specimens may leave insufficient material in the **BD SurePath** Collection Vial for preparation of a satisfactory **BD SurePath** Pap Liquid-based Test.
- Aliquot must be removed prior to processing the **BD SurePath** Liquid-based Pap Test. Only one aliquot may be removed from the **BD SurePath** Collection Vial prior to performing the **BD SurePath** Liquid-based Pap Test, regardless of the volume of the aliquot.

Procedure

1. In order to ensure a homogenous mixture, the **BD SurePath** Collection Vial must be vortexed for 10–20 seconds at 3,000 rpm and the 0.5 mL aliquot must be removed within one minute of vortexing.
2. A polypropylene aerosol barrier pipette tip that is sized appropriately for the volume being withdrawn must be used for aliquot removal. Note: Serological pipettes should not be used. Good laboratory practices must be followed to avoid introducing contaminants into the **BD SurePath** Preservative Fluid Collection Vial or the aliquot. Aliquot removal should be performed in an appropriate location outside an area where amplification is performed.
3. Visually check the aliquot material in the pipette for evidence of large gross particulates or semi-solids. Evidence of such material encountered while withdrawing the aliquot material should prompt return of all the material to the specimen vial and disqualify the specimen for ancillary testing prior to performing the Pap test.
4. For instructions on processing the aliquot using the **BD ProbeTec** CT Q^x and GC Q^x Amplified DNA Assays, refer to the assay Package Inserts provided by the manufacturer.

MATERIALS REQUIRED

Materials Provided

2 x 240 – **BD** Settling Chambers
5 x 96 – **BD SurePath** PreCoat Slides

Materials Required But Not Provided

BD SurePath Collection Vials
BD Density Reagent
BD Centrifuge Tubes
BD Syringing Pipettes
BD Aspirator Tips
BD PrepMate™ Automated Accessory
Centrifuge
Slide Racks
Easy Aspirator (optional)
Broom type sampling device or endocervical brush/plastic spatula with detachable head(s)
Vortex Mixer
Precision Pipettes with Disposable Tips
Deionized Water (pH 7.5 to 8.5)
Isopropanol and Reagent Grade Alcohol
Staining Reagents
Clearing Agent, Mounted Media, Coverslips

STORAGE

- **BD SurePath** Preservative Fluid without cytologic samples may be stored at room temperature (15–30 °C) for up to 36 months from date of manufacture.
- The storage limit for **BD SurePath** Preservative Fluid with cytologic samples is 6 months at refrigerated temperatures (2–10 °C) or 4 weeks at room temperature (15–30 °C).
- **BD SurePath** Preservative Fluid containing cytologic sample intended for use with the **BD ProbeTec** CT Q^x and GC Q^x Amplified DNA Assays can be stored and transported for up to 30 days at 2–30 °C prior to transfer to the Liquid-based Cytology Specimen (LBC) Dilution Tubes for the **BD ProbeTec** Q^x Amplified DNA Assays.

PROCEDURES

1. After the sample is collected using a Rovers Cervex-Brush or equivalent sampling device, the brush head is rinsed directly into the fluid, removed from the handle and dropped into a **BD SurePath** Collection Vial. The vial is then tightly capped, labeled, and sent to the laboratory.
2. When sample vials have been accessioned into the lab, place each vial into a tray with a labeled centrifuge tube pre-filled with 4 mL of **BD** Density Reagent and a labeled **BD SurePath** PreCoat Slide. **BD** Density Reagent must be added to the centrifuge tube before the sample is added or performance will be reduced.
3. Vigorously vortex each sample vial for 10–20 seconds at 3,000 rpm. (Sufficient volume is available in the **BD SurePath** Collection Vial to allow removal of up to 0.5 mL of homogenous mixture of cells and fluid for ancillary testing, while still allowing sufficient volume for Pap testing. The aliquoting may be performed after this vortexing step in the **BD SurePath** LBC Test process.)
4. Use the **BD PrepMate** Automated Accessory and **BD** Syringing Pipettes to transfer 8 mL of the sample into the appropriately labeled centrifuge tube containing **BD** Density Reagent. See the **BD PrepMate** Operator's Manual for instructions.
5. Place the tubes into a centrifuge rack. Arrange tubes according to the placement sequence diagram in the **BD PrepMate** Operator's Manual.
6. Placement sequence is critical and must be balanced. Balance the centrifuge tubes by adding **BD SurePath** Preservative Fluid if necessary.
7. Centrifuge specimens for 2 ± 0.25 minutes at 200 ± 25 rcf.
8. Remove centrifuge tube racks from the centrifuge.
9. Aspirate supernatant using one of the following methods.
 - Use the Easy Aspirator to aspirate supernatant.
 - Turn on the tube vac for the Easy Aspirator system and adjust the pressure to 8–10 in. Hg for a Schuco vacuum pump or 5.5 in. Hg for a KNF vacuum pump. Allow the pump to come to equilibrium vacuum pressure before beginning aspiration.
 - Place the Easy Aspirator block on a rack of disposable **BD** Aspirator Tips (clear) so that the tip adaptors on the block fit down into 12 tips. Apply moderate pressure to attach the tips onto the Easy Aspirator block. All 12 tips are needed on the aspirator block even if the centrifuge tube rack is not full.
 - Hold the Easy Aspirator block with tips over the centrifuge tubes to be aspirated. Slowly lower the aspirator tips into the supernatant, staying

just below the dropping fluid level, until the aspirator head rests evenly across the tops of the centrifuge tubes. At this point you should hear the tips drawing air into the vacuum tubing.

- Carefully withdraw the aspirator head with the tips from the centrifuge tube rack. **BD Aspirator Tips** are single use only and should be discarded after one sample transfer.
- To prevent clogs, run water through the aspirator before shutting off the vacuum pump. Perform this rinse after the last c-tube rack is aspirated and with tips still seated on the aspirator block.
- When aspiration is complete, load the Easy Aspirator block onto the tip ejector. Hold the aspirator head in front of and at the same plane as the white Delrin wedge in the top of the tip ejector. Slide the aspirator head along the top of the wedge so that the posts align into the slots of the tip ejector.
- Sliding the aspirator head into the tip ejector should eject all of the tips into the disposal tray. Withdraw the aspirator head by pulling it up and out of the tip ejector.
- Turn off the vacuum using the on/off switch on the vacuum. This will stop the vacuum and release the tips from the Easy Aspirator block into the waste bucket below the tip ejector.

Or

- Use disposable transfer pipettes to aspirate supernatant.
10. Centrifuge the tubes for 10 ± 1 minutes at 800 ± 50 rcf to concentrate the diagnostic component into a cell pellet at the bottom of the tube.
 11. Remove the tube rack from the centrifuge. In a single, rapid motion, decant the supernatant by inverting each tube rack 180 degrees so as not to disturb the cell pellet. While inverted, blot the tubes carefully on absorbent paper, making sure that the cell pellet stays in the tube. Turn the rack upright after 3 to 5 seconds.
 12. Place the slides into a Slide Rack and lock a **BD Settling Chamber** onto each slide. The position of each numbered **BD SurePath** PreCoat Slide on the Slide Rack must match to the position of the corresponding centrifuge tube.
 13. Add 4 mL of buffered deionized water (pH 7.5–8.0) to each specimen tube.
 14. Working with one sample tube at a time, use a clean disposable transfer pipette tip to mix the sample eight (8) times. Immediately transfer 800 μ L of cell suspension into the correspondingly numbered **BD Settling Chamber/BD SurePath** PreCoat Slide. Repeat for each sample.
 15. Allow 10 minutes for full sedimentation to occur. After the sedimentation, gently invert the Slide Racks to decant the remaining fluid and blot the excess liquid on absorbent paper.
 16. Rinse each **BD Settling Chamber** with 500 μ L of denatured ethanol and decant. Repeat alcohol rinse and decant the remaining fluid and blot the excess liquid on absorbent paper, allowing the **BD Settling Chamber** to remain inverted for at least 1 minute.
 17. Remove the **BD Settling Chamber** from each slide being careful to not disturb the sample deposition area.
 18. Stain and coverslip the **BD SurePath** Liquid-based Pap Test slides.

RESULTS AND INTERPRETATION

- All diagnostic criteria currently utilized in cytology laboratories for conventional Pap smears are applicable to **BD SurePath** Liquid-based Pap Test preparations.
- Any abnormal or questionable screening observations should be referred to a pathologist for review and diagnosis. Any cellular morphologic changes are significant and should be noted.

BIBLIOGRAPHY

1. Papanicolaou GN: A New Procedure for Staining Vaginal Smears. *Science* 1942; 95:438-439.
2. King A, Clay K, Felmar EG, Heustis DG, Karns RM, Krahel P, Tench WD: The Papanicolaou Smear. *West J Med* 1992; 156:202-204.
3. Mandelblatt J, Gopaul I, Wistreich M: Gynecological Care of Elderly Women: Another Look at Papanicolaou Smear Testing. *J Am Med Assoc* 1986; 256:367-371.
4. Koss LG: The Papanicolaou Test for Cervical Cancer Detection: A Triumph and a Tragedy. *J Am Med Assoc* 1989; 261:737-743.
5. Hutchinson ML, Isenstein LM, Goodman A, Hurley AA, Douglass KL, Mui KK, Patten FW, Zahniser DJ: Homogeneous Sampling Accounts for the Increased Diagnostic Accuracy Using the ThinPrep[®] Processor. *Am J Clin Pathol* 1994; 101:215-219.
6. McGoogan E, Reith A: Would Monolayers Provide More Representative Samples and Improved Preparations for Cervical Screening? Overview and Evaluation of Systems Available. *Acta Cytol* 1996; 40:107-119.
7. Bishop JW: Comparison of the CytoRich System with Conventional Cervical Cytology: Preliminary Data on 2,032 Cases from a Clinical Trial Site. *Acta Cytol* 1997; 41:15-23.

8. Geyer JW, Hancock F, Carrico C, Kirkpatrick M: Preliminary Evaluation of CytoRich: An Improved Automated Cytology Preparation. *Diagn Cytopathol* 1993; 9:417-422.
9. Grohs HK, Zahniser DJ, Geyer JW: Standardization of Specimen Preparation Through Mono/Thin-Layer Technology in Automated Cervical Cancer Screening. Edited by HK Grohs, OAN Husain. New York, Igaku-shoin, 1994, pp 176-185.
10. Solomon D, Nayar R (editors): *The Bethesda System for Reporting Cervical Cytology*. New York, Springer Verlag, 2004.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

U.S. Patent Number: 8,617,895.

Rovers Cervex-Brush is a registered trademark of Rovers Medical Devices B.V.

Cytobrush Plus and Pap Perfect are registered trademarks of CooperSurgical, Inc.

© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.