
### Revisions

Rev from	Rev to	ECO #
0508	2010/02	5270-10

**Notes:**

1. BD Cat. Number 270333, 270309
2. Blank (Sheet) Size: Length: 10.375" Width: 20"  
 Number of Pages: 8 Number of Sheets: 1  
 Page Size: Length 10.375" Width 5" Final Folded Size: See VS
3. Style (see illustrations below): # 4

4. VS Controlled by BD Caribe, LTD.
5. Ink Colors: Printed two sides  Yes  No  
 No. of Colors: 1 PMS# 2755 Blue
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	<b>Becton, Dickinson and Company</b> 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 0210761JAA		Category and Description Package Insert, RPR Card Antigen Suspension	Sheet: 1 of 9 Scale: N/A	A

# BD RPR Card Antigen Suspension

For Use in Performing the Macro-Vue™ RPR Card Tests  
(Qualitative and Quantitative) for the Serologic Detection of Syphilis

English: pages 1 – 2 Italiano: pagine 4 – 5  
Français : pages 2 – 3 Português: páginas 5 – 6  
Deutsch: Seiten 3 – 4 Español: páginas 6 – 7

0210761JAA  
2010/02

Pokyny vám poskytnú miestni zástupcov spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзініздің жергілікті BD өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

## INTENDED USE

The BD Macro-Vue™ RPR Card Tests (Qualitative and Quantitative) are for the serologic detection of syphilis. The RPR Card Antigen Suspension is a component of the test.

## SUMMARY AND EXPLANATION

RPR Card Antigen Suspension is a reagent component part used in performing the BD Macro-Vue RPR (Rapid Plasma Reagin) 18 mm Circle Test which is a nontreponemal testing procedure for the serologic detection of syphilis and is described in *A Manual of Tests for Syphilis*.<sup>1</sup> The test techniques described in this manual are the most widely used in the United States and have been thoroughly evaluated and standardized.

## PRINCIPLES OF THE PROCEDURE

The RPR Card Antigen Suspension is a carbon-particle cardiolipin antigen<sup>2,3,4</sup> which detects "reagin," an antibody-like substance present in sera from syphilitic persons, and occasionally from persons with other acute or chronic conditions.

When a specimen contains reagin, flocculation occurs with a coagglutination of the carbon particles of the RPR Card antigen, which appear as black clumps against the white background of the plastic-coated card. This coagglutination can be read macroscopically. By contrast, nonreactive specimens appear to have an even light-gray color.

## REAGENTS

The ingredients\* of the RPR Card Antigen Suspension are<sup>1</sup> 0.003% cardiolipin, 0.020 – 0.022% lecithin, 0.09% cholesterol, 0.0125 M EDTA, 0.01 M Na<sub>2</sub>HPO<sub>4</sub>, 0.01 M KH<sub>2</sub>PO<sub>4</sub>, 0.2% sodium azide (preservative), 0.02% charcoal (specially prepared, Becton Dickinson), 10% choline chloride, w/v, and deionized/distilled water.

\*Adjusted and/or supplemented as required to meet performance criteria.

## Warnings and Precautions:

For *in vitro* Diagnostic Use.

RPR Card Antigen Suspension is recommended only for use in the BD Macro-Vue RPR Card Test procedures. Prior to using the antigen, it is essential that the direction sheet accompanying each BD Macro-Vue RPR Card Test Kit be reviewed.

**Warning:** RPR Card Antigen Suspension contains sodium azide, which is harmful if swallowed. Contact with acids liberates very toxic gas. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

**Storage:** Refrigerate (2 to 8°C) the RPR Card Antigen Suspension. Once the antigen has been placed in the dispensing bottle (Cat. No. 270509) and refrigerated (2 to 8°C), the reactivity remains satisfactory for approximately three months, or until the expiration date, if it occurs sooner. Label the dispensing bottle with the antigen lot number, expiration date and date antigen was placed in the bottle.

Storage of the antigen in bright sunlight or temperatures above 30°C (86°F) should be avoided; such conditions may cause a rough appearance of the antigen when used with nonreactive sera. If the ampule of antigen is frozen during shipment, it can be reconstituted once by warming to room temperature; repeated freezing and thawing should be avoided. Immediate use of a refrigerated antigen may result in decreased sensitivity of the test. Therefore, upon removal from the refrigerator, allow the antigen to warm to room temperature (23 to 29°C) before use.

The antigen suspension should not be used beyond the expiration date.

## PROCEDURE

**Materials Provided:** RPR Card Antigen Suspension, with pre-installed ampule opener, a reagent component part of the BD Macro-Vue RPR Card Test Kits.

### Materials Required But Not Provided:

1. Antigen plastic dispensing bottles (Cat. No. 270509).
2. Dispensing square cut 20 gauge needles (Cat. No. 273509).
3. BD Macro-Vue RPR Card Test cards (Cat. Nos. 271849, 271949 or 272001).
4. Control sera with established patterns of graded reactivity which should be included in each day's testing to confirm optimal reactivity of the antigen suspension.
5. A rotator, 100 ± 2 rpm, circumscribing a circle of 2 cm in diameter, with automatic timer, friction drive and a cover containing a moistened blotter or sponge. See section "Availability" in BD Macro-Vue RPR Card Test Kit direction sheets for RPR Card Test Control Cards and BD Macro-Vue Card Test Rotator with Humidifying Cover.
6. Saline (0.9%) for use in quantitative testing. Prepare by adding 900 mg dry NaCl, ACS to 100 mL with deionized/distilled water.
7. Serum nonreactive to syphilis in 0.9% Saline; required for diluting test specimens giving a Reactive result at the 1:16 dilution.

Also required are the necessary equipment and labware used in serology testing laboratory for preparation, storage and handling of serologic specimens.

## TEST PROCEDURE, QUALITY CONTROL, RESULTS, LIMITATIONS OF THE PROCEDURE, EXPECTED VALUES AND PERFORMANCE CHARACTERISTICS

This package is a reagent component part of the **BD Macro-Vue™** RPR Card Test Kits. Refer to the direction sheet included in all kits.

### AVAILABILITY

#### Cat. No. Description

RPR Card Antigen Suspension (dated and refrigerated product):

270333 3 mL, ampule, for use with the 18 mm Circle Card Test, Box of three.

270309 3 mL, ampule, for use with the 18 mm Circle Card Test, Box of ten.

### REFERENCES

1. Creighton, E.T. 1990. Rapid plasma reagin (RPR) 18-mm circle card test, p. 99-108. *In* S.A. Larsen, E. F. Hunter, and S. J. Kraus, (ed.), *A manual of tests for syphilis*, 8th ed. American Public Health Association, Washington, D.C.
2. Portnoy, J., J. H. Brewer, and A. Harris. 1962. Rapid plasma reagin card test for syphilis and other treponematoses. *Pub. Health Rep.* 77:645-652.
3. Portnoy, J. 1963. Modifications of the rapid plasma reagin (RPR) card test for syphilis, for use in large scale testing. *Am. J. Clin. Pathol.* 40:473-479.
4. Falcone, V.H., G.W. Stout, and M.B. Moore, Jr. 1964. Evaluation of rapid plasma reagin (circle) card test. *Public Health Rep.* 79:491-495.

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Benex Limited