### Revisions

Rev from	Rev to	ECO #
0208	2010/07	5182-09

**Notes:**

1. BD Cat. Number Various
2. Blank (Sheet) Size : Length: 8.5"      Width: 11"  
 Number of Pages: 36      Number of Sheets: 9  
 Page Size: Length 8.5"      Width 5.5"      Final Folded Size: N/A
3. Style (see illustrations below): # 5

4. See Specification Control Number BALT0216647 for Material Information
5. Ink Colors: Printed two sides  Yes       No  
 No. of Colors: 1      PMS# 2755
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	Date	<p style="font-size: small;">COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION</p> <p style="font-weight: bold; font-size: large;">BD</p> <p>Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA</p>	<p>Category and Description</p> <p>Package Insert, Directigen Meningitis Individual Tests</p>		Sheet: 1 of 37	A
Proofer	Date				Scale: N/A	
Checked By	Date		<p>Part Number:</p> <p style="font-size: 1.2em;">0216647JAA</p>			

## Meningitis Individual Tests

For the detection of *Haemophilus influenzae* type b, *Streptococcus pneumoniae*, Group B *Streptococcus* or *Neisseria meningitidis* Group B / *Escherichia coli* K1.

English:	pages	1 – 6	Italiano	pagine	17 – 22
Français:	pages	7 – 12	Español:	páginas	22 – 27
Deutsch	Seiten	12 – 17	Tables:	pages	28 – 30
			Procedural Chart	page	31

0216647JAA  
2010/07

Pokyny vám poskytnú miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcije ziskate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндік жергілікті BD өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

### INTENDED USE

The **Directigen™** Meningitis Individual Tests are presumptive latex slide agglutination tests for the direct qualitative detection of antigens to *H. influenzae* type b, *S. pneumoniae*, *N. meningitidis* Group B and *Escherichia coli* K1 in cerebrospinal fluid (CSF), serum or urine. The test can also be used for the direct qualitative detection of antigens to Group B *Streptococcus* in CSF and serum. In addition, the test kit provides confirmation and serogrouping capabilities from suspected colonies of *H. influenzae* type b, *S. pneumoniae*, Group B *Streptococcus*, and *N. meningitidis* Group B. Visible agglutination occurs when a sample containing any of these bacterial antigens is reacted with its respective antibody-coated latex beads.

### SUMMARY AND EXPLANATION

The diagnosis of bacteremia and meningitis, especially in young children, can be difficult. As many as 55% of children are seen by a physician and started on antibiotics before meningitis is detected.<sup>1</sup> Detection of microbial antigens in the CSF is a rapid and helpful method for diagnostic microbiology. It may be the single most important test in cases of partially treated meningitis since Gram stain and culture may be negative. Detection of a specific antigen is a clinically significant finding and a valuable aid when choosing antimicrobial therapy.<sup>2</sup>

*Haemophilus influenzae* type b, *Neisseria meningitidis*, and *Streptococcus pneumoniae* have been reported to be the three causative agents responsible for approximately 84% of cases of bacterial meningitis.<sup>3</sup>

Group B *Streptococcus* and *Escherichia coli* K1 are major bacterial pathogens in the newborn.<sup>4-7</sup> Strains of Group B streptococci and *E. coli* K1 frequently colonize in the vagina and/or rectum and may be associated with maternal septicemia and neonatal septicemia, pneumonia and meningitis.<sup>8</sup> The *E. coli* K1 polysaccharide antigen has been shown to be structurally and immunologically similar to *Neisseria meningitidis* Group B antigen.<sup>9,10</sup> The **Directigen** *Neisseria meningitidis* Group B Latex Reagent does not differentiate the two antigens, but can be useful in the diagnosis of neonatal *E. coli* K1 meningitis. The high morbidity and mortality associated with Group B streptococci and *E. coli* K1 in newborns make rapid, accurate identification of these organisms extremely important.<sup>4</sup>

Immunological methods for detecting characteristic exoantigens of pathogenic microorganisms in patient fluids (CSF, serum, urine) are typically faster than traditional methods such as culture. These techniques include counterimmunoelectrophoresis (CIE) and latex agglutination.<sup>11-14</sup> The latex agglutination procedure has been found to be more rapid and sensitive than CIE in the detection of purified antigen.<sup>15-17</sup>

### PRINCIPLES OF THE PROCEDURE

Specific antibodies are bound to the surface of latex beads. Latex particle aggregation becomes large enough to allow rapid visualization of positive agglutination in the presence of specific antigens. These specific soluble polysaccharide antigens accumulate in CSF, serum or urine as a result of infection by *H. influenzae* type b, *S. pneumoniae*, *N. meningitidis* Groups A,B,C,Y or W135 and *Escherichia coli* K1; and in CSF or serum as a result of infection by Group B *Streptococcus*. These antigens can be detected with the **Directigen** Meningitis Combo Test Kit or the Meningitis Individual Kits (see "Availability").

## REAGENTS

### Directigen Meningitis Individual Tests:

#### *Haemophilus influenzae* type b:

<b>Reagent 1</b>	(1.0 mL),	Anti- <i>H. influenzae</i> type b Rabbit Polyclonal Antibody-Coated Latex Suspension,
<b>Reagent A</b>	(0.5 mL),	Control Latex, Rabbit Immunoglobulin-Coated Latex Suspension, each with 0.2% sodium azide (preservative).
<b>Control +</b>	(1.5 mL),	Polyvalent Positive Antigen Control, Group B <i>Streptococcus</i> , <i>N. meningitidis</i> Groups A, C, Y and W135, <i>H. influenzae</i> type b and <i>S. pneumoniae</i> antigens,
<b>Control –</b>	(1.5 mL),	Negative Antigen Control, glycine buffered saline, each with 0.2% sodium azide (preservative).

#### *Streptococcus pneumoniae*:

<b>Reagent 1</b>	(1.0 mL),	Anti- <i>S. pneumoniae</i> Rabbit Polyclonal Antibody-Coated Latex Suspension,
<b>Reagent A</b>	(0.5 mL),	Control Latex, Rabbit Immunoglobulin-Coated Latex Suspension, each with 0.2% sodium azide (preservative).
<b>Control +</b>	(1.5 mL),	Polyvalent Positive Antigen Control, <i>N. meningitidis</i> Groups A, C, Y and W135, <i>H. influenzae</i> type b and <i>S. pneumoniae</i> antigens,
<b>Control –</b>	(1.5 mL),	Negative Antigen Control, glycine buffered saline, each with 0.2% sodium azide (preservative).

#### Group B *Streptococcus*:

<b>Reagent 1</b>	(1.0 mL),	Anti-Group B <i>Streptococcus</i> Rabbit Polyclonal Antibody-Coated Latex Suspension,
<b>Reagent A</b>	(0.5 mL),	Control Latex, Rabbit Immunoglobulin-Coated Latex Suspension, each with 0.2% sodium azide (preservative).
<b>Control +</b>	(1.5 mL),	Polyvalent Positive Antigen Control, Group B <i>Streptococcus</i> , <i>N. meningitidis</i> Groups A, C, Y and W135, <i>H. influenzae</i> type b and <i>S. pneumoniae</i> antigens,
<b>Control –</b>	(1.5 mL),	Negative Antigen Control, glycine buffered saline, each with 0.2% sodium azide (preservative).

#### *Neisseria meningitidis* Group B / *E. coli* K1:

<b>Reagent 1</b>	(1.0 mL),	Anti- <i>N. meningitidis</i> Group B / <i>E. coli</i> K1 Mouse Monoclonal Antibody-Coated Latex Suspension,
<b>Reagent A</b>	(0.5 mL),	Control Latex, Mouse Immunoglobulin-Coated Latex Suspension, each with 0.2% sodium azide (preservative).
<b>Control +</b>	(1.5 mL),	Polyvalent Positive Antigen Control, <i>N. meningitidis</i> Groups A, B, C, Y and W135, <i>H. influenzae</i> type b, Group B <i>Streptococcus</i> and <i>S. pneumoniae</i> antigens,
<b>Control –</b>	(1.5 mL),	Negative Antigen Control, glycine buffered saline, each with 0.2% sodium azide (preservative).

### Precautions: For *in vitro* Diagnostic Use.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions" <sup>18–21</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Prior to discarding, sterilize specimen containers and other contaminated materials by autoclaving.

**Reagents:** Do not use beyond expiration date. Upon removal from the refrigerator, allow the reagents to warm to room temperature (15 – 30°C) before use.

To assure proper drop delivery, reagent dispensing bottle must be held vertically dispensing one free-falling drop at a time.

**Warning:** Reagents contain sodium azide, which is very toxic by inhalation, in contact with skin, and if swallowed. Contact with acids liberates very toxic gas. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

**Controls:** Do not use the kit if **Control +** and **Control –** do not yield appropriate results.

**Test Cards:** Cards must be flat for proper reactions. If necessary, flatten cards by bowing back in a direction opposite to that of the curl. Care should be taken not to finger-mark the test areas, since this may result in an oily deposit and improper test results. Use each card once and discard. Store cards in the original package in a dry area at room temperature. When spreading within confines of test areas, avoid scratching the card surface with the plastic stirrers. If the specimen does not spread to the outer perimeter of test area, use another test area of the card.

**Test Slide (Glass):** If the glass slide is used, disinfect the slide after each use, and wash before reuse. A phenolic disinfectant is suggested. Be sure all detergent and/or disinfectants are thoroughly rinsed from the slide before reuse.

**Rotation:** The recommended speed for mechanical rotation is 100 ± 2 rpm, but rotation between 95 – 110 rpm does not significantly affect the results obtained. The rotator should circumscribe a circle approximately 2 cm in diameter in the horizontal plane. A moistened humidifying cover should be used to prevent drying of test specimens during rotation.

**Storage of Reagents:** Upon receipt, remove the carton containing reagents and refrigerate at 2 – 8°C. DO NOT FREEZE. Reagents should be recapped and returned to refrigeration when not in use, taking care not to mix color-coded caps.

## SPECIMEN COLLECTION AND HANDLING

Refer to appropriate texts for details of specimen collection and handling procedures. Specimens should be tested as soon as possible; however, if the sample cannot be tested immediately, it should be stored at 2 – 8°C (for up to 48 h), or at -20°C.

Serum must be separated from whole blood prior to testing or storage.

#### SPECIMEN PRETREATMENT

The use of covered glass test tubes (borosilicate) in a 100°C (boiling) water bath is the most effective means of specimen pretreatment. Although glass test tubes in a heat block may be used, greater variation in the thorough heating of a sample may be noted as a result of variations in tube size and sample volume.

#### PROCEDURES

The testing area, all reagents, test specimens and test components should be at room temperature (15 – 30°C) when used.

**Materials Provided:** All materials as listed under "Reagents," work station and accessories.

**Materials Required But Not Provided:** **Directigen** Meningitis Combo Test Cards, **Directigen** Specimen Buffer, Rotator, 100 ± 2 rpm, Humidifying cover, Micropipettor, 50 µL delivery, Pipette tips } (see "Availability")

Also required are the necessary equipment and labware used for preparation, storage and handling of CSF, serum and urine specimens.

#### Specimen Preparation (CSF):

1. When testing CSF specimens for *H. influenzae* type b, *S. pneumoniae* and Group B *Streptococcus*, heat specimens for 3 min at 100°C (e.g., water bath or heat block) and allow to cool to room temperature before use.  
For optimal sensitivity, *N. meningitidis* Group B / *E. coli* K1 testing should be performed on unheated specimens (see "Limitations of the Procedure").
2. For CSF specimens showing turbidity, centrifuge after heating for 10 min at 1400 x g prior to testing. The supernatant fluid is to be used as the test specimen.
3. Test specimen (heated or unheated, as appropriate) as described in "Test Procedure."

#### Specimen Preparation (Serum):

1. Dilute serum specimens of at least 0.4 mL 1:1 with **Directigen** Specimen Buffer and mix.
2. Heat specimens for 5 min at 100°C (e.g., water bath or heat block) and allow to cool to room temperature before use.
3. Using a wooden applicator stick, break up the protein "clot" formed, and vortex vigorously (approximately 5 sec).
4. Centrifuge at a minimum of 1400 x g for 15 min.
5. Test supernatant fluid as described in "Test Procedure."

#### Specimen Preparation (Unconcentrated Urine):

1. Dilute urine specimens of at least 0.4 mL 1:1 with **Directigen** Specimen Buffer and mix.
2. Heat specimens for 5 min at 100°C (e.g., water bath or heat block) and allow to cool to room temperature before use.
3. Centrifuge at a minimum of 1400 x g for 10 min.
4. Test supernatant fluid as described in "Test Procedure."

#### Specimen Preparation (Concentrated Urine):

1. Urine samples that are turbid or have particulate material should be centrifuged at 1400 x g for 10 min before concentrating.
2. Urine samples may be concentrated 25-fold with a Minicon B-15 concentrator (Amicon Corporation, Danvers, Massachusetts, USA).
3. Dilute at least 200 µL of urine concentrate 1:1 with **Directigen** Specimen Buffer and mix.
4. Heat specimens for 5 min at 100°C (e.g., water bath or heat block) and allow to cool to room temperature before use.
5. Centrifuge at a minimum of 1400 x g for 10 min.
6. Test supernatant fluid as described in "Test Procedure."

#### Preparation for Confirmation of Colonies from Culture:

1. Locate suspected colonies on the agar surface from 18 – 24 h cultures that meet morphological and Gram stain characteristics of organisms that are appropriate for testing with **Directigen** Meningitis latex reagents. **Note:** A Gram stain should be performed prior to testing to ensure that organisms are appropriate for testing with the **Directigen** Meningitis latex reagents.
2. Pipette 0.5 mL (approximately 10 drops) of Control – reagent into a small glass test tube (10 x 75 mm or equivalent).
3. Select several (2 – 3) isolated colonies of similar morphology from the original or subculture plate using a sterile loop and suspend into the above tube to achieve a suspension equal to a McFarland #1 turbidity standard. **Note:** Over inoculating will yield an excessively heavy suspension which may cause erroneous results.
4. Heat the suspension for 3 min at 100°C (e.g., boiling water bath or heat block) and allow to cool to room temperature before testing.
5. Centrifuge at a minimum of 1400 x g for 10 min.
6. Test supernatant as described under "Test Procedure." **Note:** If atypical agglutination patterns are observed, refer to "Limitations of the Procedure."

#### User Quality Control:

Include **Control +** and **Control –** testing with each batch of specimens tested as described in step 1, "Test Procedure."

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

#### **Test Procedure (CSF, serum, urine, concentrated urine, and colony confirmation):**

Remove the reagents from refrigerated storage and slide the reagent tray into the slot in the kit Work Station.

**Refer to Procedural Chart Illustrations, page 31.**

The Meningitis Combo test card or glass slide may be used to test samples. If using the Meningitis Combo card for more than one sample, label circles appropriately. Before use of the glass slide, thoroughly clean with a lint-free tissue.

1. Dispense one drop of **Control +** onto the circle corresponding to Reagent 1 of row "+." Place one drop of **Control -** onto the circle corresponding to Reagent 1 of row "-."
2. Micropipette 50 µL of test sample onto the circle corresponding to Reagent 1 of row "S" and circle A.
3. Holding the dispensing bottle by the cap, vigorously swing (without inverting) to thoroughly mix each reagent. Before uncapping the bottle, gently tap the base on the counter top to assure no latex remains in the tip.
4. Dispense one drop of **Reagent A** (Control Latex Suspension) onto circle A.
5. Dispense one drop of **Reagent 1** onto the circles "+," "-" and "S" in the column corresponding to Reagent 1.
6. Mix the samples and latex reagents in each circle with a plastic stirrer, alternately using first one end, and then the opposite end for the next circle. Discard the stirrer.
7. Place the card or glass slide on a mechanical rotator and rotate at a speed of 100 ± 2 rpm for 10 min. Use a moistened humidifying cover to prevent evaporation.
8. Immediately at the end of 10 min, read the test results macroscopically under a high intensity incandescent light.

#### **INTERPRETATION OF TEST RESULTS**

Record **Control +** and **Control -** test results first.

The **Control +** should yield agglutination with **Reagent 1** within 10 min. The **Control -** should show no agglutination. Agglutination in the circle containing **Reagent A** renders the reaction uninterpretable.

Record patient test results:

**Positive Test** – Should show agglutination with **Reagent 1**. Any degree of agglutination with **Reagent 1** indicates the presence of the corresponding antigen. *Agglutination of Reagent A renders the reaction uninterpretable.*

**Negative Test** – Should show no agglutination with **Reagent 1** or **Reagent A**.

#### **LIMITATIONS OF THE PROCEDURE**

These latex agglutination tests are not intended as a substitute for bacterial culture. Confirmatory diagnosis of bacterial meningitis infection is only possible with appropriate culture procedures and should be routinely performed. Samples with extremely low levels of antigen, for instance early in the course of infection, may yield false negative results. Furthermore, samples with exceedingly high antigen concentrations may exhibit prozone effects producing inappropriately negative results. Although not extensively studied, prozone phenomena have only been observed in specimens seeded with extremely high antigen levels, and not in clinical specimens.

This assay is for the qualitative detection of antigens to *H. influenzae* type b, *S. pneumoniae*, *N. meningitidis* group B, and *E. coli* K1 in CSF, serum or urine; and to Group B *Streptococcus* in CSF and serum. In addition, the test kit can also be used with suspected colonies of *H. influenzae* type b, Group B *Streptococcus*, *S. pneumoniae*, and *N. meningitidis* Group B. Performance with other specimen types has not been assessed.

The **Directigen** Meningitis Individual Tests cannot be used with blood culture media.

A positive or negative Group B streptococcal result is only indicative of the presence or absence of Group B streptococcal antigen and is not diagnostic for the presence or absence of Group B streptococcal disease.

This device should not be used as a substitute for bacterial culture in the diagnosis of Group B streptococcal septicemia and/or meningitis.

A positive or negative result is a presumptive result for Group B streptococcal antigen. The result must be confirmed by culture.

The only *infant* specimens recommended for the direct qualitative detection of Group B streptococcal antigen are serum and cerebrospinal fluid. Testing infant urine with devices for direct qualitative detection of Group B streptococcal antigen cannot be recommended at this time; currently there is insufficient performance data supporting the use of this test on urine as a reliable predictor of Group B streptococcal disease.

Pneumococcal and *H. influenzae* type b strains not possessing a capsular antigen may not be detected by immunological techniques.

Cross reactions of *H. influenzae* reagent are known to occur in the presence of *Escherichia coli* K100. *Other cross reactions, uninterpretable and false positives may occur which are reduced or eliminated by proper specimen preparation (see "Specimen Pretreatment" and "Procedures – Specimen Preparation").*

Studies with CSF specimens seeded with *N. meningitidis* Group B and *E. coli* K1 antigen have suggested that some loss in sensitivity may result after heating as compared to unheated samples.

Untreated serum and urine specimens may yield uninterpretable results with the latex reagents. See "PROCEDURES – Specimen Preparation." CSF specimens or treated urine specimens showing any degree of turbidity should be centrifuged for 10 min at 1400 x g after heating, and prior to testing. The supernatant fluid is to be used as the test specimen.

Urine specimens processed according to "PROCEDURES – Specimen Preparation" yielding uninterpretable or suspected false positive results may be filtered using a MILLIPORE™ MilllexHA 0.45 µm filter (#SLHA-0250S) and retested with the reactive test latex(es) and the corresponding control latex(es).

Colony confirmation samples that yield uninterpretable or atypical agglutination reactions may require a 1:10 dilution of the supernatant with Control – reagent. Care should be taken to avoid bacterial suspensions greater than a McFarland #1 turbidity standard.

Latex agglutination tests of urine have potential usefulness in providing supportive evidence for the presence of pathogenic organism in infants presenting with fever without localizing signs.<sup>22</sup> However, the excretion of urinary antigen may persist in infants for 3 – 14 days, or up to 30 days, following the administration of *Haemophilus influenzae* type b (Hib) vaccines. Frequency and duration of antigen depends on the patient and the type of vaccines given. Patient vaccination history is important for the correct interpretation of positive urinary Hib antigen results in vaccinated infants.<sup>22-24</sup>

Cross reactions have been produced by some strains of the viridans streptococci group (i.e. *S. mitis*, *S. sanguis* II, *S. salivarius*)<sup>25</sup> with the **Directigen** *S. pneumoniae* latex reagent. Additional tests may be needed to differentiate the viridans streptococci group from *S. pneumoniae*.

A Gram stain should be performed prior to use with the colony confirmation procedure. Certain bacterial isolates obtained from blood culture plates which are morphologically similar to but Gram stain inconsistent with the suspect colony may exhibit non-specific agglutination and/or false positive reactions (i.e. beta-hemolytic nonpathogenic *Neisseria* species) with the Group B Strep latex.

#### EXPECTED VALUES

*Haemophilus influenzae* type b, *Neisseria meningitidis*, and *Streptococcus pneumoniae* have been reported to be the three causative agents responsible for approximately 84% of cases of bacterial meningitis.<sup>3</sup>

#### PERFORMANCE CHARACTERISTICS

**Sensitivity:** In a retrospective clinical evaluation, 177 cerebrospinal fluids (heated and unheated), 31 serum, and 40 urine specimens were tested by the **Directigen** Meningitis latex reagents and counterimmunoelectrophoresis (CIE) (Table 1, see page 28). The 248 specimens tested were originally derived from patients with culture or CIE confirmed cases of either *Haemophilus influenzae* type b, *Streptococcus pneumoniae*, *Neisseria meningitidis* Group B or Group B streptococcal infection.

*E. coli* K1 sensitivity was determined by using whole cell *E. coli* K1 seeded into body fluids and compared to commercially available latex agglutination reagent. The **Directigen** *Neisseria meningitidis* Group B Latex (Reagent 1) could detect *E. coli* K1 in CSF (4/4), sera (4/4), urine (3/4) and concentrated urine (4/4). The endpoint sensitivities were equal to or better than those observed with another commercially available latex reagent.<sup>15</sup>

The **Directigen** Meningitis latex reagents were equal to, or more sensitive than CIE in detecting antigen in CSF, serum and urine specimens.

Each **Directigen** latex reagent was also compared with counterimmunoelectrophoresis for its ability to detect partially purified antigen seeded in CSF, serum, and urine. Commercial antiserum specific for each antigen was used in the CIE test. The **Directigen** latex reagents were more sensitive than CIE.<sup>15</sup>

Additionally, a retrospective study was performed on CSF (heated and unheated) and urine specimens from patients with culture confirmed Group B streptococcal infection. **Directigen** detected antigen in seven of eight CSF specimens, 22 of 25 unconcentrated urine specimens, and ten of ten 25-fold concentrated urine specimens.

**Specificity:** The specificity of the **Directigen** Meningitis reagents was determined by testing prospective CSF, serum and urine specimens in clinical studies conducted at five institutions (Table 2, see page 29). The data indicates specificity ranges of 97 – 100% depending on the latex and specimen tested.

**Colony Confirmation:** Latex reagents were tested using suspensions of isolated colonies that met morphological characteristics of the organisms. Performance characteristics are listed in Table 3 (see page 30).

#### AVAILABILITY

##### Cat. No. Description

**Directigen™** Kits:

252260 *H. influenzae* type b Test Kit, (30 patient tests, 60 controls).

251960 *S. pneumoniae* Test Kit, (30 patient tests, 60 controls).

255460 Group B Strep Test Kit, (30 patient tests, 60 controls).

255560 *N. meningitidis* Group B/*E. coli* K1 Test Kit, (30 patient tests, 60 controls).

250160 *N. meningitidis* Groups A, C, Y and W135 Test Kit, (30 patient tests, 60 controls).

252360 Meningitis Combo Test Kit, (30 patient tests, 60 controls).

252480 **Directigen™** Meningitis Combo Test Card (single use), Box 30.

256391 **Directigen™** Specimen Buffer, 8 mL.

252350 **Directigen™** Meningitis Negative Control Reagent, 9 mL, Box 4.

278051 **Macro-Vue™** Card Test Rotator with humidifying cover, 100 ± 2 rpm, automatic timer, friction drive, Model-51-1L.

## REFERENCES

1. Rothrock, S.G., S.M. Green, J. Wren, D. Letai, L. Daniel-Underwood, E. Pillar: Pediatric bacterial meningitis: is prior antibiotic therapy associated with unaltered clinical presentation?, *Annals Emerg. Med.* 27:146–152, 1992.
2. Fasola, E. and P. Ferrieri: Laboratory diagnostic methods for central nervous system infections, *Neuro. Clinics N. Amer.* 3:279–290, 1992.
3. Morbidity and Mortality Weekly Report (CDC), Surveillance summary: bacterial meningitis and meningococemia, United States, 28:277–279, 1978.
4. Finegold, S.M. and Martin, W.J.: *Bailey and Scott's Diagnostic Microbiology*, 6th Ed., C.V. Mosby Co., 1982.
5. Baker, C.J.: Summary of the workshop on perinatal infections due to group B *Streptococcus*, NIH News, *J. Infect. Dis.*, 136:137–152, 1977.
6. Mulder, C.J.J., Van Alphen, L. and Zanen, H.C.: Neonatal meningitis caused by *Escherichia coli* in the Netherlands., *J. Infect. Dis.*, 150:935, 1984
7. Sarff, L.D., McCracken, G.H., Jr., Schiffer, M.S., Gloche, M. P., Robbins, J.B. and Orskov, F.: Epidemiology of *Escherichia coli* K1 in healthy and diseased newborns., *Lancet*, 7:1099–1104, 1975.
8. Blaunstein, H., Tucker, E.B. and Gibson, B.C.: Identification and significance of *Streptococcus agalactiae*, *Am. J. Clin. Path.*, 57:207–213, 1969.
9. Kasper, D.L., Winkelhake, J.L., Zollinger, W.D., Brandt, B.L. and Artenstein, M.S.: Immunochemical similarity between polysaccharide antigens of *Escherichia coli* 07:K1(L): NM and group B *Neisseria meningitidis*, *J. Immunol.*, 110:262–268, 1973.
10. Lively, M.R., A.S. Gilbert and C. Moreno: Sialic acid polysaccharide antigens of *Neisseria meningitidis* and *Escherichia coli*: esterification between adjacent residues., *Carbohydr. Res.*, 94:193–203, 1981.
11. Edwards, E.A.: Immunologic investigations of meningococcal disease, group-specific *Neisseria meningitidis* antigens present in the serum of patients with fulminant meningococemia., *J. Immunol.* 106:314–317, 1971.
12. Whittle, H. C., Egler, J.J., Tugwell, P. and Greenwood, B.M.: Rapid bacteriological diagnosis of pyogenic meningitis by latex agglutination. *Lancet*, 2:619–621, 1974.
13. Ingram, D.L., Anderson, P. and Smith, D.H.: Countercurrent immunoelectrophoresis in the diagnosis of disease caused by *Haemophilus influenzae* type b., *J. Ped.*, 81:1156–1159, 1972.
14. Newman, R.B., Stevens, R.W. and Gaafar, H.H.: Latex agglutination test for the diagnosis of *Haemophilus influenzae meningitis*, *J. Lab. Clin. Med.*, 76:107–113, 1970.
15. Data on file at BD Diagnostic Systems.
16. Tilton, R.C., Dias, F. and Ryan, R.W.: Comparative evaluation of three commercial products and counterimmuno-electrophoresis for the detection of antigens in cerebrospinal fluid, *J. Clin. Microbiol.*, 20:231–234, 1984.
17. Rench, M.A., Metzger, T. G. and Baker, C. J.: Detection of group B streptococcal antigen in body fluids by a latex- coupled monoclonal antibody assay, *J. Clin. Microbiol.*, 20:852–854, 1984.
18. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, PA.
19. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53–80.
20. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC) 5<sup>th</sup> ed. U.S. Government Printing Office, Washington, D.C.
21. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021-0045.
22. Jones, R.G., Bass, J.W., Weisse, M.E., and Vincent, J.M. Antigenuria after immunization with *Haemophilus influenzae* oligosaccharide CRM197 conjugate (HbOC) vaccine. *Pediatr. Inf. Dis. J.* 1991; 10:557–9.
23. Goepf, J.G., Hohenboken, M., Almeida-Hill, J., and Santocham, M. Persistent urinary antigen excretion to infants vaccinated with *Haemophilus influenzae* type b capsular polysaccharide conjugated with outer membrane protein from *Neisseria meningitidis*. *Pediatr. Inf. Dis. J.* 1992; 11:2–5.
24. Rothstein, E.P., Madore, D.V., et al. Comparison of antigenuria after immunization with three *Haemophilus influenzae* type b conjugate vaccines. *Pediatr. Inf. Dis. J.* 1991; 10:331–14.
25. Holmberg, H., Danielsson D, Hardie, J., Krook, A., Whiley, R.: Cross-reactions between  $\alpha$ -streptococci and omniserum, a polyvalent pneumococcal serum, demonstrated by direct immunofluorescence, immunoelectrophoresis, and latex agglutination, *J. Clin Microbiol.*, 21:745–748, 1985.

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MILLIPORE is a trademark of Millipore Corporation.
BD, BD Logo, Directigen and Macro-Vue are trademarks of Becton, Dickinson and Company © 2010 BD.
Benex Limited