
## Revisions

Rev from	Rev to	ECO #
0503	0304	2676-04

**Notes:**

1. BD Cat. Number 494775, 494776, 494480
2. Blank (Sheet) Size : Length: 8.5"    Width: 11"  
 Number of Pages: 24    Number of Sheets: 6  
 Page Size: Length 8.5"    Width 5.5"                      Final Folded Size: 8.5" x 5.5"
3. Style (see illustrations below): #5

4. See Specification Control Number 4400362JAA for Material Information
5. Ink Colors: Printed two sides  Yes     No  
 No. of Colors: 2                      PMS# 279 Blue, PMS# 230 Pink
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	<b>Becton, Dickinson and Company</b> 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 4400362JAA		Category and Description Package Insert, Qtest Strep Test 20, 40 & 80 Tests	Sheet: 1 of 25 <hr/> Scale: 1:1	A

# Qtest™ Strep

4400362JAA  
2004/03

**For the Rapid Detection of Group A  
*Streptococcus***

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Italiano: pagine 15 – 19  
Español: páginas 19 – 23

See symbol glossary at end of insert.

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20/40/80

# BD Qtest™ Strep

## For the Rapid Detection of Group A *Streptococcus*

4400362JAA  
2004/03

### INTENDED USE

Qtest™ Strep test is a liposome immunoassay test for the rapid, qualitative detection of group A streptococcal antigen directly from throat swabs.

### SUMMARY AND EXPLANATION

Group A streptococci are organisms that typically cause illness such as tonsillitis, pharyngitis and scarlet fever. If untreated, these infections can lead to complications such as rheumatic fever. The Qtest Strep test is performed directly on throat-swab extracted antigens and the test is completed in approximately two minutes after collection of the specimen. Antibody-coated liposomes are used to detect the extracted antigen.

### PRINCIPLES OF THE PROCEDURE

To perform the test, a throat swab specimen is collected, then streptococcal antigens are extracted from the specimen with Qtest Strep A, B and C reagents. The resulting solution of extracted antigen is added to the Qtest device. The Qtest device contains a membrane that has antibodies to group A *Streptococcus* coated on its surface. As the specimen passes through the membrane, the extracted antigen attaches to the antibodies on the surface.

Qtest Strep 1 and 2 reagents are then added to the Qtest device. When reagent 1 is added to the Qtest device, antibody-coated liposomes containing pink dye attach to the streptococcal antigen previously bound to the membrane. If group A streptococcal antigen is present in the specimen, a pink triangle will appear, indicating a positive result. In the absence of antigen, a pink negative sign (-) will appear, indicating a negative result.

The Qtest Strep test contains built-in quality control features that include both positive and negative internal procedural controls for the Qtest device and color-coded reagents for the extraction step of the assay. The appearance of a pink negative sign (-) when a specimen is negative is the internal positive procedural control and indicates that reagents 1 and 2 were working properly and added in the correct sequence. The membrane background area is the internal negative procedural control for the Qtest device and should be white-to-light pink in color. If the test is not conducted correctly, or the reagents are not performing properly, no distinct test result will appear.

### REAGENTS

**Qtest™ Strep A:** 2M sodium nitrite. (7.0 mL sufficient for 20 or 40 tests or 14.0 mL sufficient for 80 tests).

**Qtest™ Strep B:** 2N acetic acid. (7.0 mL sufficient for 20 or 40 tests or 14.0 mL sufficient for 80 tests).

**Qtest™ Strep C:** Buffered 0.66N sodium hydroxide. (8.0 mL sufficient for 20 tests, 12.0 mL sufficient for 40 tests or 24.0 mL sufficient for 80 tests).

**Qtest™ Strep 1:** Anti-group A *Streptococcus* (rabbit antibody) coated liposomes, with 0.2% sodium azide (preservative). (4.5 mL sufficient for 20 tests, 9.0 mL sufficient for 40 tests or 18.0 mL sufficient for 80 tests).

**Qtest™ Strep 2:** Buffered wash solution (5.5 mL sufficient for 20 or 40 tests or 11.0 mL sufficient for 80 tests).

**Qtest™ Strep Devices:** Membrane coated with Anti-group A *Streptococcus* (rabbit antibody) and group A streptococcal antigen (one per test).

**Sterile Throat Swabs:** One per test.

**DispensTube™ Specimen Holder:** One per test.

**DispensTube™ Tips and Tubes:** One per test.

**Positive Control:** Heat inactivated group A *Streptococcus* with 0.2% sodium azide (preservative). (2.0 mL).

**Negative Control:** Heat inactivated *Streptococcus*, not group A, with 0.2% sodium azide (preservative). (2.0 mL).

### Warnings and Precautions:

For *in vitro* Diagnostic Use.

**Warning:** Qtest Strep A and Positive Control contain combustible and toxic materials. Contact with combustible material may cause fire. Toxic if swallowed. Keep away from combustible material. When using, do not eat or drink. Wear suitable protective clothing, gloves and eye/face protection. In case of accident or if you feel unwell, seek medical advice immediately.

Qtest Strep B contains acid; Qtest Strep C contains strong alkali and causes burns. If these reagents contact skin or eyes, rinse immediately with plenty of water and seek medical advice.

Qtest Strep 1 and the Negative Control contain sodium azide which is very toxic by inhalation, in contact with skin and if swallowed. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. If disposed into a drain, flush with a large volume of water to prevent azide build-up.

**Reagents:** Do not use beyond the expiration date marked on the kit. *Do not mix reagents from different kit lot numbers.* Do not remove the Qtest device from the pouch until ready to use.

To assure proper drop delivery, reagent bottles must be held vertically while dispensing free-falling drops.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>1-4</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Prior to discarding, sterilize containers and other contaminated materials by autoclaving.

**Swabs:** Use the swabs supplied in this kit to obtain throat specimens.

**DispensTube Tip:** Must contain white filter materials to ensure proper test performance.

**Storage:** Refrigerate reagents at 2 – 8°C. DO NOT FREEZE. Allow reagents to reach room temperature before use. See room temperature storage option on procedure card.

**Extractions:** Optimal results are obtained when the throat swab is squeezed during the extraction process.

**Interpretation:** Interpretation of test results is best made immediately after Qtest Strep 2 is absorbed, as color intensity may change over time.

#### SPECIMEN HANDLING

Collect a throat specimen by standard recommended procedures using the swabs supplied in this kit. Consult standard reference procedures, such as the collection method described by Facklam.<sup>5</sup>

Processing samples for Qtest Strep testing is recommended as soon as possible after collection. Swabs supplied with the kit provide optimal results. Rayon or polyester fiber swabs on a plastic shaft, without transport media, may be used. Transport devices containing charcoal, agar, wooden shafts or calcium alginate material should not be used for rapid antigen tests. If culture is desired, streak the swab onto a culture plate using standard laboratory procedure before conducting the Qtest Strep test.


#### PROCEDURES

**Materials Provided:** See “Reagents” for materials provided.


**Materials Required But Not Provided:** Required are the necessary equipment and labware used for preparation, storage and handling of throat specimens.

#### Test Procedures

##### Sample Extraction


**A** Add 3 drops Qtest Strep A into the **DispensTube** sample tube.

**B** Add 3 drops Qtest Strep B into **DispensTube** sample tube. The solution should turn gold. Immediately place the throat swab specimen in the **DispensTube** sample tube. **Mix by rolling the swab against the side of the DispensTube sample tube, squeezing the swab through the tube several times during the extraction process so that the liquid is expressed from the swab and reabsorbed.** Let stand for 1 min but not longer than 3 min.

**CAUTION:** Assay performance depends upon thorough antigen extraction.


**C** Holding swab shaft to the side, add 5 drops Qtest Strep C directly into the **DispensTube** sample tube. Mix the solution by squeezing the tube and moving the swab. The color of the solution should change to red. If the solution remains purple, continue squeezing the swab. Alternatively, an additional drop of Qtest Strep B may be added. Solution should then turn red.


Remove the liquid from the swab by rolling the swab against the wall of the **DispensTube** sample tube and squeezing the flexible sides. Discard the swab.

Insert the **DispensTube Tip** into the **DispensTube** sample tube. The red solution may be held for up to 1 h at room temperature before testing, without affecting the results.


## Sample Testing

Invert **DispensTube** sample tube; allow the fluid to be absorbed into the filter material. Dispense entire contents of the **DispensTube** sample tube into the **Qtest** well. Allow the fluid to be absorbed. Proceed immediately to Step 1.

- 1** Place 3 drops **Qtest** Step 1 into the **Qtest** well. Allow the fluid to be absorbed. Proceed immediately to Step 2.

- 2** Place 3 drops **Qtest** Step 2 into the **Qtest** well. Allow the fluid to be absorbed. Immediately interpret the test results.

**Quality Control:** Each **Qtest** Step device contains built-in procedural controls that are reagent-activated. The appearance of a pink negative sign (-) and a white-to-light pink background area in the **Qtest** well validate the immunological integrity of the test. It also provides assurance that **Qtest** Step 1 and **Qtest** Step 2 were added in the correct sequence.

**NOTE:** A pink negative sign (-) may not be apparent with a strong positive sample reaction.

A procedural control is also provided for the specimen extraction step of the test. A sequence of color changes occurs with each step of the extraction process: red to gold, and then back to red. This assures extraction reagent integrity and proper extraction procedure.

**Positive and Negative Liquid Controls** are supplied with each kit. Liquid controls may be used in place of a patient sample in order to monitor the integrity of the reagents/device, assure that the test procedure has been performed correctly, and/or to demonstrate a positive or negative reaction. *After shaking or vortexing thoroughly*, add 2 drops of Control into a **DispensTube** sample tube containing 3 drops **Qtest** Step A and 3 drops **Qtest** Step B. Add dry swab. Mix. After 1 min, but not longer than 3 min, add 5 drops **Qtest** Step C. Mix well, discard swab and insert **DispensTube** Tip. Dispense entire contents of **DispensTube** sample tube into the **Qtest** device. Proceed with Step 1 "Sample Testing."

**Quality Control:** Positive and Negative Liquid Controls are provided in order to monitor the integrity of the reagents and assure that the test procedure has been performed correctly. Prior to using a new kit or lot number of the **Qtest** Step test, it is recommended that the Positive and Negative Controls be tested in accordance with the laboratory's established Quality Control procedures. Upon observing the expected results, the kit is ready for use with patient samples.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

If expected control results are not obtained, do not report patient results. Contact your local BD representative or Technical Services for assistance.

## RESULTS

Interpretation of test results is best made immediately after Qtest Strep 2 is absorbed, as color intensity may change over time.

Positive


**Positive Test** - A pink triangle appears in the Qtest well. The background area should be white-to-light pink.

**Note:** The negative sign may be visible; however, the presence of any shade of a pink triangle indicates a positive test.

Negative

**Negative Test** - A pink negative sign (-) appears in the Qtest well. The background area should be white-to-light pink and no triangle should be visible.

Uninterpretable

No distinct test result will appear if reagents are added incorrectly or are not performing properly. For questions, contact your local BD representative or Technical Services for assistance.

## LIMITATIONS OF THE PROCEDURE

Respiratory infections, including pharyngitis, can be caused by streptococci from serogroups other than group A as well as other pathogens. The Qtest Strep test will not differentiate asymptomatic carriers of group A *Streptococcus* from those exhibiting streptococcal infection. In rare cases, test specimens heavily colonized with *Staphylococcus aureus* can yield false positive results. If clinical signs and symptoms are not consistent with clinical test results, follow-up culture and grouping procedures should be performed.

Use of liquid transport systems may result in reduced intensity of test results relative to testing with the kit swab or a comparable dry swab.

## PERFORMANCE CHARACTERISTICS<sup>6</sup>

The ability of the Qtest Strep test to identify group A *Streptococcus* from throat swab specimens was determined by correlation with conventional culture techniques during internal trials. One hundred eighty eight (188) throat swab specimens were collected on swabs from patients with pharyngitis. Swabs in transport medium were shipped overnight from 15 clinical sites to an internal laboratory for processing. Immediately after arrival, each swab was inoculated onto a sheep blood agar plate. Swabs were either refrigerated or frozen while inoculated plates were incubated at 37°C for 48 h. Beta-hemolytic colonies were identified as group A streptococci by serologic method. After culture results were determined, stored swabs were selected on the basis of culture classification result and the Qtest Strep test was performed.

Swabs were selected to yield a population with a prevalence of infection between 25 – 30%, dominated by 3+ and 4+ culture classification (as is usually seen during strep season).

Culture Versus Qtest Strep			
		Culture	
		+	-
Qtest Strep	+	49	5
	-	3	131

Specimen Mix in Positive Sample Population	
Culture Classification	% in Sample Pop.
1+	2
2+	11
3+	56
4+	31

Sensitivity 94.2%  
 Specificity 96.3%  
 Overall Accuracy 95.7%

## AVAILABILITY

Cat. No.	Description
494775	Qtest™ Strep 20 Test Kit
494776	Qtest™ Strep 40 Test Kit
494780	Qtest™ Strep 80 Test Kit

## REFERENCES

1. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, PA.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
3. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
5. Facklam, R.R. U.S. Dept. of Health and Human Services, PHS CDC, Pub. No. CDC 77-13.
6. Data on file, BD Diagnostic Systems.

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BENEX Limited