
## Revisions

SO 0191-5

Rev from	Rev to	ECO #
03	2010/06	5372-10

**Notes:**

1. BD Cat. Number 260683, 260684, 260685
2. Blank (Sheet) Size: Length: 17.0"    Width: 22.0"  
 Number of Pages: 16    Number of Sheets: 1  
 Page Size: Length 8.5"    Width 5.5"    Final Folded Size: 4.25" x 5.5"
3. Style (see illustrations below): # 7

4. See Specification Control Number VS-80-1041-9PR for Material Information
5. Ink Colors: Printed two sides  Yes     No  
 No. of Colors: 1    PMS# 2755
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA		
Proofer	Date			Category and Description Package Insert, GasPak EZ Gas Generating Pouch Systems	Sheet: 1 of 17
Checked By	Date				Scale: N/A
Part Number: 8010419					

# **BD GasPak™ EZ Gas Generating Pouch Systems**

English: pages 1 – 3      Italiano: pagine 7 – 9  
Français : pages 3 – 5      Español: páginas 9 – 11  
Deutsch: Seiten 5 – 7      Japanese: pages 11 – 12

8010419  
2010/06

Pokyny vám poskytnú miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / **Свържете се с местния представител на BD за инструкции.** / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / **Өзіңіздің жергілікті BD өкіліне жүгініп нұсқау алыңыз.** / Kontaktiraj lokalnog predstavnika BD za upute.

## INTENDED USE

The GasPak™ EZ Gas Generating Pouch Systems are single-use systems that produce atmospheres suitable to support the primary isolation and cultivation of anaerobic, microaerophilic, or capnophilic bacteria by use of gas generating sachets inside single-use resealable pouches.

## SUMMARY AND EXPLANATION

In 1965, Brewer and Allgeier introduced a disposable hydrogen generator envelope, which was later modified to include carbon dioxide generation and an internal catalyst.<sup>1,2</sup> The GasPak EZ Gas Generating Pouch Systems offer bio-performance and reliability without the need for a catalyst or an activation procedure. The GasPak EZ Gas Generating Sachets contain all the ingredients needed to create a specific atmosphere for specimen incubation. The sachet and specimens are placed in the GasPak EZ resealable pouch. The pouch is sealed and incubated.

## PRINCIPLES OF THE PROCEDURE

The GasPak EZ Gas Generating Sachet consists of a reagent sachet containing inorganic carbonate, activated carbon, ascorbic acid and water. When the sachet is removed from the outer wrapper, the sachet becomes activated by exposure to air. The activated reagent sachet and specimens are placed in the resealable pouch and sealed. The sachet rapidly reduces the oxygen concentration within the pouch. At the same time, inorganic carbonate produces carbon dioxide.

For the cultivation of anaerobic bacteria, the GasPak EZ Anaerobe Pouch System with Indicator Sachets produce an anaerobic atmosphere within 2.5 h with greater than or equal to 10% carbon dioxide within 24 h.

For the cultivation of microaerophilic bacteria, the GasPak EZ Campy Pouch System Sachets produce an atmosphere with approximately 5 – 15% oxygen.

For the cultivation of carbon dioxide requiring organisms, the GasPak EZ CO<sub>2</sub> Pouch System Sachets produce an atmosphere with greater than 3% carbon dioxide.

## REAGENTS

### Warnings and Precautions:

For *in vitro* Diagnostic Use.

Observe established precautions against microbiological hazards throughout all procedures. Prior to discarding, sterilize specimen containers and other contaminated material by autoclaving.

After use, the sachet may remain warm. Discard sachet after reaction is complete (when cool).

After use, discard pouch in a proper waste container. Do not reuse.

DO NOT use sachet if outer foil is damaged or open in any manner.

Avoid direct sunlight and excessive temperatures.

DO NOT open resealable pouch until incubation is complete.

GasPak EZ Pouch System Sachets are intended for use in GasPak EZ resealable pouches. Their use in other systems may yield misleading results.

GasPak EZ resealable pouches are designed for GasPak EZ Pouch System Sachets. Do not use with other gas generating envelopes (e.g., BBL™ GasPak™/GasPak™ Plus Envelopes).

**Storage Instructions:** On receipt, store sachets at 2 – 25°C in a dry environment.

GasPak EZ Pouch System Sachets are ready to use. The expiration date is for unopened, intact and properly stored sachets. Do not open until ready to use.

GasPak EZ Dry Anaerobic Indicator should appear white before use. Do not use indicator if blue in color before use. Do not touch indicator tablet or open indicator package to remove tablet.

**Product Deterioration:** Do not use reagent sachets if outer package has been damaged or opened. Do not use resealable pouches that show rips or tears or other irregularities.

**Active Ingredients:** Ascorbic acid, activated carbon and water

**Warning:** Activated carbon is irritating to eyes, respiratory system and skin.

## PROCEDURE

**Materials Provided:** GasPak EZ Gas Generating Pouch System. Each system requires the use of a resealable pouch and sachet (See "Availability").

**Materials Required But Not Provided:** Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

### Test Procedure:

1. Place the inoculated plates inside the resealable pouch. The GasPak EZ Campy Pouch System requires at least 2 Petri dishes and will hold up to 4 petri dishes. The GasPak EZ CO<sub>2</sub> Pouch System and the GasPak EZ Anaerobe Pouch System with Indicator can be used for 1-4 Petri dishes.

**Note:** For optimum growth using the GasPak EZ Campy Pouch System, place a paper towel or cotton ball moistened with 5 mL of water inside the pouch.

2. One GasPak EZ Pouch System Sachet is used in each resealable pouch. Remove the GasPak EZ Pouch System Sachet from the carton. Remove the outer foil packaging.
3. Place the activated sachet in the GasPak EZ resealable pouch with the plates. The sachet should be placed between the plates and the pouch. If using the GasPak EZ Anaerobe Pouch System with Indicator, there is no need to add a separate anaerobic indicator.
4. Close the pouch by pressing the zipper part of the pouch together.
5. Incubate the GasPak EZ Pouch System at a temperature appropriate for the organism being cultured, but not to exceed 45°C.
6. After incubation, open the pouch, remove the plates and dispose of the GasPak EZ Pouch System Sachet, indicator and pouch in the appropriate manner.

### User Quality Control:

Each GasPak EZ Gas Generating Pouch System should be tested periodically for its ability to provide adequate conditions for the growth of appropriate bacteria.

System	QC Organism
GasPak™ EZ Anaerobe Pouch System with Indicator	<i>Bacteroides fragilis</i> ATCC™ 25285
GasPak™ EZ Campy Pouch System	<i>Campylobacter jejuni</i> ATCC 29428
GasPak™ EZ CO <sub>2</sub> Pouch System	<i>Neisseria gonorrhoeae</i> ATCC 19424

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

## RESULTS

Anaerobic conditions are achieved rapidly within 2.5 h, with greater than or equal to 10% carbon dioxide within 24 h at 35°C. Visible condensate should occur within 30 min of activation. Blood-containing agar plates appear reduced within 2 – 4 h at 35°C. The anaerobic indicator should appear reduced (white) within 2 – 4 h at 35°C.

Microaerophilic conditions are achieved rapidly with an oxygen concentration of 5 – 15% at 35°C.

Carbon dioxide enriched conditions are achieved rapidly with a carbon dioxide concentration of greater than 3% at 35°C.

## PERFORMANCE CHARACTERISTICS

### GasPak™ EZ Anaerobe Pouch System with Indicator

Prior to release, representative samples of each lot of GasPak EZ Anaerobe Pouch System with Indicator are tested for performance characteristics.

Each sample sachet of the GasPak EZ Anaerobe Pouch System with Indicator is placed in a resealable incubation pouch with 4 media filled Petri dishes. The sachet is activated and the pouch is sealed according to label directions and incubated at 35 ± 2°C. After approximately 2.5 h, a 22-gauge 1.5" needle is inserted into the incubation pouch and a sample of the gas is removed, placed in a gas chromatograph, and analyzed for the amount of oxygen present in the incubation pouch. The percent oxygen for each pouch is less than or equal to 1% and the mean percent for all samples tested for oxygen is less than 0.7%. After approximately 24 h incubation at 35 ± 2 °C, another sample is taken from the pouch and analyzed for the amount of carbon dioxide. The percent carbon dioxide for each pouch is greater than or equal to 10%. All anaerobic dry indicators are reduced.

### GasPak™ EZ Campy Pouch System

Prior to release, representative samples of each lot of GasPak EZ Campy Pouch System are tested for performance characteristics.

Each sample sachet of the GasPak EZ Campy Pouch System is placed in a resealable incubation pouch along with 2 media filled Petri dishes. The sachet is activated and the pouch is sealed according to label directions and incubated at 35 ± 2°C. After approximately 24 h, a 22-gauge 1.5" needle is inserted into the incubation pouch and a sample of the gas is removed, placed in a gas chromatograph, and analyzed for the amount of oxygen present in the incubation pouch. The percent oxygen for each pouch is 5 to 15%.

### GasPak™ EZ CO<sub>2</sub> Pouch System

Prior to release, representative samples of each lot of GasPak EZ CO<sub>2</sub> Pouch System are tested for performance characteristics.

Each sample sachet of the GasPak EZ CO<sub>2</sub> Pouch System is placed in a resealable incubation pouch along with 2 media filled Petri dishes. The sachet is activated and the pouch is sealed according to label directions and incubated at 35 ± 2°C. After approximately 24 h, a 22-gauge 1.5" needle is inserted into the incubation pouch and a sample of the gas is removed, placed in a gas chromatograph, and analyzed for the amount of carbon dioxide present in the incubation pouch. The percent carbon dioxide for each pouch is greater than 3%.

#### AVAILABILITY

Cat. No.	Description
260683	<b>GasPak™</b> EZ Anaerobe Pouch System with Indicator (Containing 20 sachets, 20 resealable pouches and 20 Dry Anaerobic Indicators)
260685	<b>GasPak™</b> EZ Campy Pouch System (Containing 20 sachets and 20 resealable pouches)
260684	<b>GasPak™</b> EZ CO <sub>2</sub> Pouch System (Containing 20 sachets and 20 resealable pouches)

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ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BBL and GasPak are trademarks of Becton, Dickinson and Company. 2006 BD.
BBL, BD, BD Logo and GasPak are trademarks of Becton, Dickinson and Company. ©2010 BD.
Benex Limited