Revisions

SO 0191-5

Rev from	Rev to	ECO #
2010/08	2010/12	5522-10

Notes:

1. BD Cat. Number 441124
2. Blank (Sheet) Size: Length: 8.5" Width: 11"
 Number of Pages: 140 Number of Sheets: 35
 Page Size: Length 8.5" Width 5.5" Final Folded Size: 8.5" x 5.5"
3. Style (see illustrations below): # 5

4. See Specification Control Number BALT 8081409 for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8081409		Category and Description Package Insert, BD ProbeTec GC Qx Amplified DNA Assay		Sheet: 1 of 141
				Scale: N/A
				A

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INTENDED USE

The **BD ProbeTec™ GC Q^x Amplified DNA Assay**, when tested with the **BD Viper™ System** in Extracted Mode, uses Strand Displacement Amplification technology for the direct, qualitative detection of *Neisseria gonorrhoeae* DNA in clinician-collected female endocervical and male urethral swab specimens, patient-collected vaginal swab specimens (in a clinical setting), and male and female urine specimens. The assay is also intended for use with gynecological specimens collected in **BD SurePath™ Preservative Fluid** or **PreservCyt™ Solution** using an aliquot that is removed prior to processing for either the **BD SurePath™** or **ThinPrep™** Pap test. The assay is indicated for use with asymptomatic and symptomatic individuals to aid in the diagnosis of gonococcal urogenital disease.

SUMMARY AND EXPLANATION

The World Health Organization estimates that 62 million new cases of infection due to *Neisseria gonorrhoeae* are diagnosed each year.¹ In the United States, gonorrhea is the second most commonly reported infectious disease with over 358,000 cases in 2006.² During this period, for the first time in 10 years the reported rate of infection was higher among women than men.² Infection of women is often asymptomatic and if left untreated can lead to pelvic inflammatory disease, infertility, ectopic pregnancy and chronic pelvic pain. In men, symptoms of acute urethritis and dysuria usually cause infected individuals to present for treatment before serious sequelae result. Transmission of *N. gonorrhoeae* occurs through sexual contact but can also take place in the birth canal leading to neonatal conjunctivitis.

Because of the high frequency of asymptomatic infections, the US Preventive Services Task Force has published recommendations for screening young, sexually active women and those who are older and considered at increased risk of infection in order to prevent complications and reduce transmission.³ The Advisory Committee on Human Immunodeficiency Virus (HIV) and Sexually Transmitted Disease (STD) Prevention also encourages active control programs that target treatable STDs as a primary intervention in the HIV epidemic.⁴ Nevertheless, a recent rise in fluoroquinolone resistance has reduced the options available to combat *N. gonorrhoeae* infection such that administration of cephalosporins is the only treatment recommended by the Centers for Disease Control and Prevention.⁵

N. gonorrhoeae are gram-negative, oxidase-positive diplococci that can be observed in Gram-stained smears of urethral discharge, usually within neutrophils. Culture of *N. gonorrhoeae* can be difficult because the organism does not survive long outside the host and is highly susceptible to adverse environmental conditions such as lack of humidity and temperature extremes. Although culture of urogenital swabs remains an important tool in the diagnosis of *N. gonorrhoeae* infection due to the continued need for monitoring of antimicrobial susceptibility, use of molecular methods that amplify and detect specific nucleic acid sequences is increasing due to their applicability to both swab specimens and more easily collected urine specimens.^{6,7}

When used with the **BD Viper System**, the **BD ProbeTec GC Q^x Amplified DNA Assay** involves automated ferric oxide-based extraction of DNA from clinical specimens using **BD FOX™** Extraction technology with the chemical lysis of cells, followed by binding of DNA to para-magnetic particles, washing of the bound nucleic acid and elution in an amplification-compatible buffer. When present, *N. gonorrhoeae* DNA is then detected by Strand Displacement Amplification (SDA) of a specific target sequence in the presence of a fluorescently-labeled detector probe.^{8,9}

PRINCIPLES OF THE PROCEDURE

The **BD ProbeTec GC Q^x Amplified DNA Assay** is designed for use with the **BD ProbeTec Chlamydia trachomatis/Neisseria gonorrhoeae (CT/GC) Q^x specimen collection and transport devices**, applicable reagents, the **BD Viper System** and **BD FOX** Extraction. Specimens are collected and transported in their respective transport devices which preserve the integrity of the *N. gonorrhoeae* DNA over the specified ranges of temperature and time.

Urine and swab specimens undergo a pre-warm step in the **BD Viper Lysing Heater** to dissolve mucus and homogenize the specimen. After cooling, the specimens are loaded onto the **BD Viper System** which then performs all the steps involved in extraction and amplification of target DNA, without further user intervention. For gynecological specimens that are collected and transported in **BD SurePath Preservative Fluid** or **PreservCyt Solution**, the pre-warm step is not necessary; i.e., an aliquot is simply transferred to a Liquid-Based Cytology Specimen (LBC) Dilution Tube for the **BD ProbeTec Q^x Amplified DNA Assays** prior to loading on the instrument. The specimen is transferred to an Extraction Tube that contains ferric oxide particles in a dissolvable film and dried Extraction Control. A high pH is used to lyse the bacterial cells and liberate their DNA into solution. Acid is then added to lower the pH and induce a positive charge on the ferric oxide, which in turn binds the negatively charged DNA. The particles and bound DNA are then pulled to the sides of the Extraction Tube by magnets and the treated specimen is aspirated to waste. The particles are washed and a high pH Elution Buffer is added to recover the purified DNA. Finally, a Neutralization Buffer is used to bring the pH of the extracted solution to the optimum for amplification of the target.

The **BD ProbeTec GC Q^x** Amplified DNA Assay is based on the simultaneous amplification and detection of target DNA using amplification primers and a fluorescently-labeled detector probe.^{8,9} The reagents for SDA are dried in two separate disposable microwells: the Priming Microwell contains the amplification primers, fluorescently-labeled detector probe, nucleotides and other reagents necessary for amplification, while the Amplification Microwell contains the two enzymes (a DNA polymerase and a restriction endonuclease) that are required for SDA. The **BD Viper System** pipettes a portion of the purified DNA solution from each Extraction Tube into a Priming Microwell to rehydrate the contents. After a brief incubation, the reaction mixture is transferred to a corresponding, pre-warmed Amplification Microwell which is sealed to prevent contamination and then incubated in one of the two thermally-controlled fluorescent readers. The presence or absence of *N. gonorrhoeae* DNA is determined by calculating the peak fluorescence (Maximum Relative Fluorescent Units [MaxRFU]) over the course of the amplification process and by comparing this measurement to a predetermined threshold value.

In addition to the fluorescent probe used to detect amplified *N. gonorrhoeae* target DNA, a second fluorescently-labeled oligonucleotide is incorporated in each reaction. The Extraction Control (EC) oligonucleotide is labeled with a different dye than that used for detection of the *N. gonorrhoeae*-specific target and is used to confirm the validity of the extraction process. The EC is dried in the Extraction Tubes and is re-hydrated upon addition of the specimen and extraction reagents. At the end of the extraction process, the EC fluorescence is monitored by the **BD Viper** instrument and an automated algorithm is applied to both the EC and *N. gonorrhoeae*-specific signals to report specimen results as positive, negative, or EC failure.

REAGENTS

Each **BD ProbeTec GC Q^x** Reagent Pack contains:

- GC Q^x Amplified DNA Assay Priming Microwells, 12 x 96: each Priming Microwell contains approximately 30 pmol oligonucleotides, 45 pmol fluorescently-labeled detector probe, 100 nmol dNTPs, with stabilizers and buffer components.
- GC Q^x Amplified DNA Assay Amplification Microwells, 12 x 96: each Amplification Microwell contains approximately 140 units DNA polymerase and 500 units restriction enzyme, with stabilizers and buffer components.

NOTE: Each microwell pouch contains one desiccant bag.

Control Set for the **BD ProbeTec CT/GC Q^x** Amplified DNA Assays: 24 CT/GC Q^x Positive Control Tubes containing approximately 2400 copies each of pCTB4 and pGCint3 linearized plasmids in carrier nucleic acid, and 24 CT/GC Q^x Negative Controls Tubes containing carrier nucleic acid alone. The concentrations of the pCTB4 and pGCint3 plasmids are determined by UV spectrophotometry.

Q^x Swab Diluent for the **BD ProbeTec Q^x** Amplified DNA Assays: 48 tubes each containing approximately 2 mL of potassium phosphate/potassium hydroxide buffer with DMSO and preservative.

Liquid-Based Cytology Specimen (LBC) Dilution Tube for the **BD ProbeTec Q^x** Amplified DNA Assays (LBC Specimen Dilution Tube): 400 tubes each containing approximately 1.7 mL of Tris/Sodium Chloride solution and preservative.

BD FOX Extraction Tubes: 48 strips of 8 tubes, each containing approximately 10 mg of iron oxide in a dissolvable film and approximately 240 pmol fluorescently-labeled Extraction Control oligonucleotide.

BD Viper Extraction Reagent and Lysis Trough: each 4-cavity Extraction Reagent trough contains approximately 16.5 mL Binding Acid, 117 mL Wash Buffer, 35 mL Elution Buffer, and 29 mL Neutralization Buffer with preservative; each Lysis Trough contains approximately 11.5 mL Lysis Reagent.

INSTRUMENT, EQUIPMENT AND SUPPLIES

Materials Provided: **BD Viper** Instrument, **BD Viper** Instrument Plates, **BD Viper** Pipette Tips, **BD Viper** Tip Waste Boxes, **BD Viper** Amplification Plate Sealers (Black), **BD Viper** Lysing Heater, **BD Viper** Lysing Rack, **BD Viper** Neutralization Pouches, Specimen Tubes and Caps for use on the **BD Viper** System (Extracted Mode), Urine Preservative Transport for the **BD ProbeTec Q^x** Amplified DNA Assays (Q^x UPT), **BD ProbeTec™** Q^x Collection Kit for Endocervical or Lesion Specimens, Male Urethral Specimen Collection Kit for the **BD ProbeTec CT/GC Q^x** Amplified DNA Assays, Vaginal Specimen Transport for the **BD ProbeTec CT/GC Q^x** Amplified DNA Assays, **BD ProbeTec** Accessories, Liquid-Based Cytology Specimen (LBC) Dilution Tube Caps for the **BD ProbeTec Q^x** Amplified DNA Assays, **BD Viper** Liquid-Based Cytology Specimen Rack.

Materials Required But Not Provided: Nitrile gloves, 1% (v/v) sodium hypochlorite*, DNA AWAY™, 3% (w/v) hydrogen peroxide, *Neisseria gonorrhoeae* ATCC™ 19424 (diluted in phosphate buffered saline) or Bio-Rad AmpliTrol™ CT/GC, *Chlamydia trachomatis* ATCC VR-879 (Serovar H) or VR-902B (LGV II) (diluted in phosphate buffered saline), displacement pipettes, polypropylene aerosol-resistant pipette tips capable of delivering 0.5 ± 0.05 mL, and a vortex mixer.

*Mix 200 mL of bleach with 800 mL of water. Prepare fresh daily.

Storage and Handling Requirements: Reagents may be stored at 2 – 33°C. Unopened Reagent Packs are stable until the expiration date. Once a pouch is opened, the microwells are stable for 6 weeks if properly sealed or until the expiration date, whichever comes first. Do not freeze.

Warnings and Precautions

General:

1. For *in vitro* Diagnostic Use.
2. Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁰⁻¹³ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.
3. For additional specific warnings, cautions and notes specific to the **BD Viper**, consult the **BD Viper** System User's Manual.

Specimen:

4. For collection of endocervical swab specimens, use only the **BD ProbeTec Q^x** Collection Kit for Endocervical or Lesion Specimens.

5. For patient-collection and transport of vaginal swabs, use only the Vaginal Specimen Transport for the **BD ProbeTec CT/GC Q^x Amplified DNA Assays**.
6. For collection of male urethral swab specimens, use only the Male Urethral Specimen Collection Kit for the **BD ProbeTec CT/GC Q^x Amplified DNA Assays**.
7. For urine specimens, use only the Q^x UPT or unpreserved (neat) urine.
8. Under or over filling Specimen Tubes or the Q^x UPT with urine may affect assay performance. Over filling the tubes may also result in liquid overflow on the **BD Viper** deck, and could cause contamination.
9. For male urethral and female endocervical swab specimens, specimens must be collected and tested before the expiration date of the Q^x Swab Diluent tube.
10. For vaginal specimens, specimens must be collected and processed before the expiration date of the Vaginal Specimen Transport. Once expressed, specimens must be tested before the expiration date of the Q^x Swab Diluent tube.
11. For urine specimens, specimens must be tested before the expiration date of the Q^x UPT.
12. For liquid-based cytology specimens, use only the Liquid-Based Cytology Specimen (LBC) Dilution Tube for the **BD ProbeTec Q^x Amplified DNA Assays**.
13. Liquid-based cytology solutions contain flammable substances. Do not place specimens transferred to the LBC Specimen Dilution Tubes in the **BD Viper** Lysing Rack or the Lysing Heater. Specimens transferred to the LBC Specimen Dilution Tubes should be placed in the **BD Viper** LBC Specimen Rack.
14. For testing with the **BD ProbeTec CT/GC Q^x Amplified DNA Assays** on the **BD Viper** System in Extracted Mode, be sure to obtain aliquots of specimens collected in **BD SurePath** Preservative Fluid or PreservCyt Solution prior to processing for either the **BD SurePath** or ThinPrep Pap test. Failure to do so may result in erroneous results.
15. The **BD ProbeTec CT/GC Q^x Amplified DNA Assays** may not be used with **BD SurePath** or PreservCyt residual specimens.
16. Do not run PreservCyt specimens that have been treated with glacial acetic acid on the **BD Viper** System in Extracted Mode. Extraction Control failures or False Negative results may occur.
17. Use only polypropylene aerosol-resistant pipette tips to transfer specimens to the LBC Specimen Dilution Tubes.
18. Liquid-based cytology specimens must be tested before the expiration date of the LBC Specimen Dilution Tube.

Assay/Reagent:

19. This reagent pack is for testing endocervical and patient-collected vaginal swabs (in a clinical setting), male urethral swabs, liquid-based cytology specimens, and male and female urine specimens with the **BD Viper** System in Extracted Mode.
20. The Q^x UPT contains **NAP Guard™** (approximately 742.5 mM K₂EDTA). **NAP Guard** may be irritating to the eyes, skin and respiratory system. In case of contact with eyes, rinse opened eye immediately with plenty of water and seek medical advice if symptoms persist. After contact with skin, wash immediately with plenty of soap and water. If inhaled, seek medical attention in case of problems.
21. Use only sample and control tubes with pierceable caps on the **BD Viper** System in Extracted Mode. Do not remove pierceable caps prior to running the instrument. Be sure to replace any punctured pierceable caps with new pierceable caps prior to running the instrument.
22. Do not interchange or mix kit reagents from kits with different lot numbers.
23. The Q^x Swab Diluent for the **BD ProbeTec Q^x Amplified DNA Assays** contains dimethyl sulfoxide (DMSO). DMSO is harmful by inhalation, in contact with skin and if swallowed. Avoid contact with eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water.
24. Do not test the Q^x Swab Diluent tube from the Endocervical/Lesion or the Male Urethral Specimen Collection Kits if received in the laboratory without the swab present. A false negative test result may occur.
25. Use only the **BD Viper** pipette tips as supplied by BD with the **BD Viper** System.
26. The **BD Viper** Extraction Reagent and Lysis Troughs contain corrosive substances. The use of personal protection equipment, such as nitrile gloves, safety glasses, and lab coats is strongly recommended when handling these reagents. These solutions have a strong caustic effect, and may cause severe burns, to skin and mucous membranes. Avoid contact with the eyes and skin. Avoid breathing fumes, vapors, or spray. Harmful if swallowed. Do not eat or drink in the vicinity of these reagents. In case of contact, immediately remove contaminated clothing. Wash the skin with water and soap and rinse thoroughly. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
27. Use only the **BD Viper** Amplification Plate Sealers (Black) on the Amplification plates with the **BD Viper** System. Using the clear sealers for sealing the Amplification plates may cause erroneous results.
28. Reagent pouches containing unused Priming Microwells and Amplification Microwells **MUST** be carefully resealed after opening. Verify that desiccant is present prior to resealing the reagent pouches.
29. Because the CT/GC Q^x Positive control is used for both CT Q^x and GC Q^x testing, correct positioning of the microwell strips is important for final results reporting.
30. The plate containing the Amplification Microwells **MUST** be properly sealed with the **BD Viper** Amplification Plate Sealer (Black) prior to moving the plate from the **BD Viper** System. Sealing ensures a closed reaction for amplification and detection and is necessary to avoid contamination of the instrument and work area with amplification products. **Do not remove sealing material from microwells at any time.**
31. Priming Microwells with residual fluid (after transfer of liquid from the Priming Microwells to the Amplification Microwells) represent a source of target contamination. Carefully seal Priming Microwells with plate sealer prior to disposal.
32. To prevent contamination of the work environment with amplification products, use the disposal bags provided in the Accessory kit to dispose of tested Amplification Microwells. Make sure the bags are properly closed before disposal.

33. Although dedicated work areas are not required because the **BD Viper** design reduces the possibility of amplicon contamination in the testing environment, other precautions for controlling contamination, particularly to avoid contamination of specimens during manipulation, are necessary.
34. **CHANGE GLOVES** if they come in contact with specimen or appear to be wet, to avoid contaminating other specimens. Change gloves before leaving work area and upon entry into work area.
35. In the event of contamination of the work area or equipment with specimens or controls, thoroughly clean the contaminated area with 1% (v/v) sodium hypochlorite, DNA AWAY, or 3% (w/v) hydrogen peroxide (do not use hydrogen peroxide from a bottle that has remained open for longer than 8 days) and rinse thoroughly with water. Allow surface to dry completely before proceeding.
36. In case of a spill on the **BD Viper** Lysing Rack, immerse the rack in 1% (v/v) sodium hypochlorite for 1 – 2 min. Do not exceed 2 min. Thoroughly rinse the rack with water and allow to air dry.
37. Clean the entire work area – counter tops and instrument surfaces – with 1% (v/v) sodium hypochlorite on a daily basis. Thoroughly rinse with water. Allow surfaces to dry completely before proceeding with additional testing.
38. Contact BD Technical Services in the event of an unusual situation, such as a spill into the **BD Viper** instrument or DNA contamination that cannot be removed by cleaning.
39. Acid and Base spill kits should be on hand in the event of a spill of extraction reagents.

SWAB SPECIMEN COLLECTION, STORAGE AND TRANSPORT

For swab specimens, performance data in this package insert have been established with the **BD ProbeTec** collection kits listed. Performance with collection devices other than those listed has not been evaluated.

- **BD ProbeTec Q^x Collection Kit for Endocervical or Lesion Specimens.**
- Vaginal Specimen Transport for the **BD ProbeTec CT/GC Q^x Amplified DNA Assays**
- Male Urethral Specimen Collection Kit for the **BD ProbeTec CT/GC Q^x Amplified DNA Assays**

Swab Specimen Collection

Endocervical Swab Specimen Collection using BD ProbeTec Q^x Collection Kit for Endocervical or Lesion Specimens

1. Remove the cleaning swab from packaging.
2. Using the polyester fiber-tipped cleaning swab with the white shaft, remove excess blood and mucus from the cervical os.
3. Discard the used cleaning swab.
4. Remove the pink collection swab from packaging.
5. Insert the collection swab into the cervical canal and rotate for 15 – 30 s.
6. Withdraw the swab carefully. Avoid contact with the vaginal mucosa.
7. Uncap the Q^x Swab Diluent tube.
8. Fully insert the collection swab into the Q^x Swab Diluent tube.
9. Break the shaft of the swab at the score mark. Use care to avoid splashing of contents.
10. **Tightly** recap the tube.
11. Label the tube with patient information and date/time collected.
12. Transport to laboratory.

Vaginal Swab Patient Collection Procedure using Vaginal Specimen Transport for the BD ProbeTec CT/GC Q^x Amplified DNA Assays

NOTE: Ensure that patients read the Patient Collection Instructions before providing them with a collection kit.

1. Wash hands with soap and water. Rinse and dry.
2. It is important to maintain a comfortable balance during the collection procedure.
3. Twist the cap to break the seal. Pull the cap with attached swab from the tube. Do not touch the soft tip or lay the swab down. If you touch or drop the swab tip or the swab is laid down, discard the swab and request a new vaginal swab.
4. Hold the swab by the cap with one hand so that the swab tip is pointing toward you.
5. With your other hand, gently spread the skin outside the vagina. Insert the tip of the swab into the vaginal opening. Point the tip toward your lower back and relax your muscles.
6. Gently slide the swab no more than 2 inches into the vagina. If the swab does not slide easily, gently rotate the swab as you push. **If it is still difficult, do not attempt to continue.** Make sure the swab touches the walls of the vagina so that moisture is absorbed by the swab.
7. Rotate the swab for 10 – 15 s.
8. Withdraw the swab without touching the skin. Place the swab in the tube and cap securely.
9. After collection, wash hands with soap and water, rinse, and dry.
10. Return the tube with the swab to the nurse or clinician as instructed.
11. Label with patient information and date/time collected.
12. Transport to laboratory.

Male Urethral Swab Specimen Collection using Male Urethral Specimen Collection Kit for the BD ProbeTec CT/GC Q^x Amplified DNA Assays

1. Remove the swab from packaging.
2. Insert the swab 2 – 4 cm into the urethra and rotate for 3 – 5 s.
3. Withdraw the swab.
4. Uncap the Q^x Swab Diluent tube.
5. Fully insert the collection swab into the Q^x Swab Diluent tube.

6. Break the shaft of the swab at the score mark. Use care to avoid splashing of contents.
7. **Tightly** recap the tube.
8. Label the tube with patient information and date/time collected.
9. Transport to laboratory.

Swab Storage and Transport

Table 1 provides instructions for storage and transport conditions to the laboratory and/or test site for swab specimens. The endocervical and the male urethral swab specimens must be stored and transported to the laboratory and/or test site within 30 days after collection if kept at 2 – 30°C or within 180 days after collection if kept frozen at -20°C. Patient-collected vaginal swab specimens must be stored and transported to the laboratory and/or test site within 14 days after collection if kept at 2 – 30°C or within 180 days after collection if kept frozen at -20°C. Patient-collected vaginal swab specimens that are expressed in Q^x Swab Diluent may be stored and processed within 30 days after expression if kept at 2 – 30°C or within 180 days after the date of expression if kept frozen at -20°C.

Table 1: Swab Specimen Storage and Transport

SWAB SPECIMEN TYPE TO BE PROCESSED	FEMALE ENDOCERVICAL SWAB SPECIMEN/MALE URETHRAL SWAB SPECIMEN		VAGINAL SWAB SPECIMEN			
			DRY VAGINAL SWAB SPECIMEN (COLLECTION SITE)		EXPRESSED VAGINAL SWAB SPECIMEN (TEST SITE)	
Temperature Condition for Transport to Test Site and Storage	2 – 30°C	-20°C	2 – 30°C	-20°C	2 – 30°C	-20°C
Process Specimen According to Instructions	Within 30 days of collection	Within 180 days of collection	Express and process within 14 days of collection	Express and process within 180 days of collection	Within 30 days of expression	Within 180 days of expression

For U.S. and international shipments, specimens should be labeled in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and etiologic agents/infectious substances. Time and temperature conditions for storage must be maintained during transport.

URINE SPECIMEN COLLECTION, STORAGE AND TRANSPORT

For urine specimens, performance has been established with the Q^x UPT and with urine collected in a sterile, plastic, preservative-free, specimen collection cup (i.e., neat urine without preservatives). Performance with other collection methods and collection devices has not been established.

Urine Specimen Collection

1. The patient should not have urinated for at least 1 h prior to specimen collection.
2. Collect the specimen in a sterile, preservative-free specimen collection cup.
3. The patient should collect the first 20 – 60 mL of voided urine (the first part of the stream – NOT midstream) into a urine collection cup.
4. Cap and label with patient identification and date/time collected.

Urine Transfer to Q^x UPT

NOTE: Urine specimens should be transferred from the collection cup to the Q^x UPT within 8 h of collection if the urine specimen has been stored at 2 – 30°C. Urine Specimens stored at 2 – 8°C can be held up to 24 h prior to transfer to the Q^x UPT.

Wear clean gloves when handling the Q^x UPT tube and urine specimen. If gloves come in contact with the specimen, immediately change them to prevent contamination of other specimens.

1. Open the Q^x UPT Collection and Transport Kit and remove the Q^x UPT and transfer pipette from their packaging.
2. Label the Q^x UPT with the patient identification and date/time collected.
3. Hold the Q^x UPT upright and firmly tap the bottom of the tube on a flat surface to dislodge any large drops from inside the cap. Repeat if necessary.
4. Uncap the Q^x UPT and use the transfer pipette to dispense urine into the tube. The correct volume of urine has been added when the fluid level is between the purple lines on the fill window located on the Q^x UPT label. This volume corresponds to approximately 2.0 – 3.0 mL of urine. DO NOT overfill or under fill the tube.
5. Discard the transfer pipette in a biohazard waste container.

NOTE: The transfer pipette is intended for use with a single specimen.

6. Tighten the cap securely on the Q^x UPT.
7. Invert the Q^x UPT 3 – 4 times to ensure that the specimen and reagent are well mixed.

Q^x UPT Urine Storage and Transport

Store and transport Q^x UPT urine specimens at 2 – 30°C and pre-warm them within 30 days of transfer to the Q^x UPT. Specimens may be stored in the Q^x UPT at -20°C for up to 180 days prior to pre-warming.

Neat Urine Storage and Transport

Store and transport neat urine specimens from the collection site to the test site at 2 – 8°C and pre-warm them within 7 days of collection. Neat urine stored at 2 – 30°C must be pre-warmed within 30 h of collection. Neat urine specimens may also be stored frozen at -20°C for up to 180 days prior to pre-warming.

Table 2: Urine Specimen Storage and Transport

Urine Specimen Type to be Processed	Q ^x UPT			NEAT		
Urine handling options prior to transfer to Q ^x UPT	Store urine specimen at 2 – 30°C and transfer to Q ^x UPT within 8 h of collection or Store urine specimen at 2 – 8°C and transfer to Q ^x UPT within 24 h of collection or Transfer to Q ^x UPT immediately					
Temperature Condition for Storage and Transport to Test Site	2 – 8°C	2 – 30°C	-20°C	2 – 8°C	2 – 30°C	-20°C
Process and Test Specimen According to Instructions	Within 30 days after transfer to Q ^x UPT		Within 180 days after transfer to Q ^x UPT	Within 7 days of collection	Within 30 h of collection	Within 180 days of collection

LBC SPECIMEN COLLECTION, STORAGE AND TRANSPORT

BD SurePath or PreservCyt specimens must be collected using either an endocervical broom or a brush/spatula combination as described in the **BD SurePath** or PreservCyt product insert. Once collected, **BD SurePath** or PreservCyt specimens can be stored and transported in their original vials for up to 30 days at 2 – 30°C prior to transfer to LBC Specimen Dilution Tubes.

Specimen Transfer to LBC Specimen Dilution Tubes

A 0.5 mL aliquot of either the **BD SurePath** or PreservCyt specimen must be transferred from the original vial to the LBC Specimen Dilution Tube prior to processing for either the **BD SurePath** or ThinPrep Pap test.

Wear gloves when handling the LBC Specimen Dilution Tube and the **BD SurePath** or PreservCyt specimen vial. If gloves come in contact with the specimen, immediately change them to prevent contamination of other specimens.

BD SurePath Specimen Transfer

NOTE: Refer to the **BD PrepStain™** Slide Processor Product Insert for instructions on removing an aliquot from the **BD SurePath** specimen vial prior to performing the **BD SurePath** liquid-based Pap test.

1. Label an LBC Specimen Dilution Tube with patient identification information.
2. Remove the cap from the LBC Specimen Dilution Tube.
3. Transfer 0.5 mL from the specimen vial to the LBC Specimen Dilution Tube. Avoid pipetting fluid from the bottom of the vial. Discard pipette tip.

NOTE: A separate pipette tip must be used for each specimen.

4. Tighten the cap on the LBC Specimen Dilution Tube securely.
5. Invert the LBC Specimen Dilution Tube 3 – 4 times to ensure that the specimen and diluent are well mixed.

PreservCyt Specimen Transfer

NOTE: Refer to the ThinPrep™ 2000/3000 System Operator's Manual Addendum for instructions on removing an aliquot from the PreservCyt specimen vial prior to performing the ThinPrep Pap test.

1. Label an LBC Specimen Dilution Tube with patient identification information.
2. Remove the cap from the LBC Specimen Dilution Tube.
3. Transfer 0.5 mL from the specimen vial to the LBC Specimen Dilution Tube. Avoid pipetting fluid from the bottom of the vial. Discard pipette tip.

NOTE: A separate pipette tip must be used for each specimen.

4. Tighten the cap on the LBC Specimen Dilution Tube securely.
5. Invert the LBC Specimen Dilution Tube 3 – 4 times to ensure that the specimen and diluent are well mixed.

Storage and Transport of Specimens Transferred to the LBC Specimen Dilution Tubes

After transfer to an LBC Specimen Dilution Tube the diluted specimen can be stored at 2 – 30°C for up to 30 days. Diluted specimens may also be stored at -20°C for up to 90 days.

SWAB SPECIMEN PROCESSING

Processing procedure for the BD ProbeTec Q^x Collection Kit for Endocervical or Lesion Specimens or the Male Urethral Specimen Collection Kit for the BD ProbeTec *Chlamydia trachomatis/Neisseria gonorrhoeae* (CT/GC) Q^x Amplified DNA Assays

NOTE: If specimens are refrigerated or frozen, make sure they are brought to room temperature and mixed by inversion prior to proceeding.

1. Using the tube layout report, place the Q^x Swab Diluent tube with black pierceable cap in order in the BD Viper Lysing Rack and lock into place.
2. Repeat step 1 for additional swab specimens.
3. Specimens are ready to be pre-warmed.
4. **Change gloves** before proceeding to avoid contamination.

Processing procedure for the Vaginal Specimen Transport for the BD ProbeTec *Chlamydia trachomatis/Neisseria gonorrhoeae* (CT/GC) Q^x Amplified DNA Assays

NOTE: Wear clean gloves when handling the vaginal swab specimen. If gloves come in contact with specimen, immediately change them to prevent contamination of other specimens.

NOTE: If specimens are refrigerated or frozen, make sure they are brought to room temperature prior to expression.

1. Label a pre-filled Q^x Swab Diluent tube for each swab specimen to be processed.
2. Remove the cap and insert the swab specimen into the Q^x Swab Diluent. Mix by swirling the swab in the Q^x Swab Diluent for 5 – 10 s.
3. Express the swab along the inside of the tube so that liquid runs back into the bottom of the tube.
4. Remove the swab carefully from the Q^x Swab Diluent tube to avoid splashing.
5. Place the expressed swab back into the transport tube and discard with biohazardous waste.
6. Tightly recap the Q^x Swab Diluent tube with the **black pierceable cap**.
7. Repeat steps 1 – 6 for additional swab specimens.
8. Using the tube layout report, place the tube in order in the **BD Viper** Lysing Rack and lock into place.
9. Specimens are ready to be pre-warmed.
10. **Change gloves** before proceeding to avoid contamination.

URINE SPECIMEN PROCESSING

NOTE: If specimens are refrigerated or frozen, make sure they are brought to room temperature and mixed by inversion prior to proceeding.

Processing procedure for the Q^x UPT

1. Make sure the urine volume in each Q^x UPT tube falls between the lines indicated on the tube label. Under or over filling the tube may affect assay performance. Over filling the tube may also result in liquid overflow on the **BD Viper** deck, and could cause contamination.
2. Make sure the Q^x UPT tube has a **black pierceable cap**.
3. Repeat steps 1 and 2 for additional Q^x UPT tube specimens.
4. Using the tube layout report, place the Q^x UPT tube in order in the **BD Viper** Lysing Rack and lock into place.
5. Specimens are ready to be pre-warmed.
6. **Change gloves** before proceeding to avoid contamination.

Processing procedure for unpreserved (Neat) urine specimens

NOTE: Wear clean gloves when handling the urine specimen. If gloves come in contact with specimen, immediately change them to prevent contamination of other specimens.

1. Label a Specimen Tube for use on the **BD Viper** System (Extracted Mode) with the patient identification and date/time collected.
2. Swirl the urine cup to mix the urine specimen and open carefully.
NOTE: Open carefully to avoid spills which may contaminate gloves or the work area.
3. Uncap the tube and use a pipette to transfer the urine specimen into the tube. The correct volume of urine has been added when the fluid level is between the purple lines on the fill window located on the label. This volume corresponds to approximately 2.0 – 3.0 mL of urine. DO NOT overfill or under fill the tube.
4. Tighten a **black pierceable cap** securely on each tube.
5. Repeat steps 1 through 4 for each urine specimen. Use a new pipette or pipette tip for each sample.
6. Using the tube layout report, place the neat urine specimens in order in the **BD Viper** Lysing Rack and lock into place.
7. Specimens are ready to be pre-warmed.
8. **Change gloves** before proceeding to avoid contamination.

NOTE: The pre-warm step must be started within 30 h of collection if the urine has been stored at 2 – 30°C; within 7 days of collection if stored at 2 – 8°C; or within 180 days if stored frozen at -20°C.

PROCESSING PROCEDURE FOR LBC SPECIMENS TRANSFERRED TO THE LBC SPECIMEN DILUTION TUBES

NOTE: Do not place specimens transferred to the LBC Specimen Dilution Tubes in the **BD Viper** Lysing Rack or the **BD Viper** Lysing Heater. Specimens transferred to the LBC Specimen Dilution Tubes should be placed in the **BD Viper** LBC Specimen Rack.

NOTE: If specimens are frozen, make sure they are thawed completely at room temperature and mixed by inversion prior to proceeding.

1. Make sure the LBC Specimen Dilution Tube has a blue pierceable cap.
2. Using the tube layout report, place the LBC Specimen Dilution Tube containing the specimen in order in the **BD Viper** LBC Specimen Rack and lock into place.
3. Specimens are ready to be tested on the **BD Viper** System in Extracted Mode.
4. Change gloves prior to proceeding to avoid contamination.

QUALITY CONTROL PREPARATION

NOTE: Do not re-hydrate the controls prior to loading in the BD Viper Lysing Rack.

1. Using the tube layout report, place CT/GC Q^x Negative Controls into the appropriate positions in the **BD Viper** Lysing Rack.
2. Using the tube layout report, place CT/GC Q^x Positive Controls into the appropriate positions in the **BD Viper** Lysing Rack.
3. Controls are ready to be pre-warmed with the specimens if desired.

PRE-WARM PROCEDURE FOR SWAB AND URINE SPECIMENS

NOTE: The pre-warm procedure must be applied to all swab and urine specimens to ensure that the specimen matrix is homogeneous prior to loading on the BD Viper System. Failure to pre-warm specimens may have an adverse impact on performance of the BD ProbeTec CT/GC Q^x assays and/or BD Viper System. Swab and urine specimens must be pre-warmed; however, pre-warming of the controls is optional.

NOTE: Refrigerated or frozen specimens must be brought to room temperature prior to pre-warming.

1. Insert the **BD Viper** Lysing Rack into the **BD Viper** Lysing Heater.
2. Pre-warm the specimens for 15 min at 114 +/- 2°C.
3. Remove the Lysing Rack from the Lysing Heater and let specimens cool at room temperature for a minimum of 15 min before loading into the **BD Viper** instrument.
4. Refer to the Test Procedure for testing specimens and controls.
5. After pre-warming, specimens may be stored for 7 days at 2 – 30°C or for 180 days at -20°C without additional pre-warming prior to testing on the **BD Viper** System.

TEST PROCEDURE

Refer to the **BD Viper** Instrument User's Manual (Extracted Mode Operation) for specific instructions for operating and maintaining the components of the system. The optimum environmental conditions for the GC Q^x assay were found to be 18 – 27°C and 20 – 85% Relative Humidity.

QUALITY CONTROL

Quality control must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

The Control Set for the **BD ProbeTec** CT/GC Q^x Amplified DNA Assays is provided separately. One Positive and one Negative Control must be included in each assay run and for each new reagent kit lot number. Controls must be positioned according to the **BD Viper** Instrument User's Manual. The CT/GC Q^x Positive Control will monitor for substantial reagent failure only. The CT/GC Q^x Negative Control monitors for reagent and/or environmental contamination. Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organizations. Refer to CLSI C24-A3 for additional guidance on appropriate internal quality control testing practices.¹⁴ The Positive Control contains approximately 2400 copies per mL of pCTB4 and pGCint3 linearized plasmids.

The Extraction Control (EC) oligonucleotide is used to confirm the validity of the extraction process. The EC is dried in the Extraction Tubes and is re-hydrated by the **BD Viper** System upon addition of the specimen and extraction reagents. At the end of the extraction process, the EC fluorescence is monitored by the instrument and an automated algorithm is applied to both the EC and *N. gonorrhoeae*-specific signals to report specimen results as positive, negative, or EC failure.

General QC Information for the BD Viper System:

The location of the microwells is shown in a color-coded plate layout screen on the LCD Monitor. The plus symbol (+) within the microwell indicates the positive QC sample. The minus symbol (-) within the microwell indicates the negative QC sample.

A QC pair must be logged in for each reagent kit lot number and for each plate to be tested. If QC pairs have not been properly logged in, a message box appears that prevents saving the rack and proceeding with the run until complete.

A maximum of two QC pairs per rack is permitted. Additional control materials may be added provided they are logged in as samples.

NOTE: The BD Viper System will re-hydrate the controls during the assay run. Do not attempt to hydrate the assay controls prior to loading them into the BD Viper Lysing Rack.

Running one plate on a BD Viper System:

The first two positions (A1 and B1) are reserved for the positive (A1) and negative (B1) controls, respectively. The first available position for a patient sample is C1.

Running two plates on a BD Viper System:


For plate one, the first two positions (A1 and B1) are reserved for the positive (A1) and negative (B1) controls, respectively. The first available position for a patient sample is C1. For plate two (full plate) the last two positions (G12 and H12) are reserved for the positive (G12) and negative (H12) controls, respectively. For plate two (partial plate) the last two positions after the last patient sample are automatically assigned as the positive and negative controls, respectively.

Interpretation of Quality Control Results:

The CT/GC Q^x Positive Control and the CT/GC Q^x Negative Control must test as positive and negative, respectively, in order to obtain patient results. If controls do not perform as expected, the run is considered invalid and patient results will not be reported by the instrument. If either of the controls does not provide the expected results, repeat the entire run using a new set of controls, new extraction tubes, new extraction reagent trough, new lysis trough and new microwells. If the repeat QC does not provide the expected results, contact BD Technical Services.

If the *N. gonorrhoeae*-specific signal is greater than or equal to a threshold of 125 Maximum Relative Fluorescent Units (MaxRFU), the EC fluorescence is ignored by the algorithm. If the *N. gonorrhoeae*-specific signal is less than a threshold of 125 MaxRFU, the EC fluorescence is utilized by the algorithm in the interpretation of the result.

Table 3: Interpretation of Quality Control Results

Control Type	Tube Result Report Symbol	GC Q ^x MaxRFU	QC Disposition
GC Q ^x Positive Control	OK	≥125	QC Pass
GC Q ^x Positive Control		<125	QC Failure
GC Q ^x Positive Control	 or  or 	Any value	QC Failure
GC Q ^x Negative Control	OK	<125	QC Pass
GC Q ^x Negative Control		≥125	QC Failure
GC Q ^x Negative Control	 or  or  or 	Any value	QC Failure

Refer to the Interpretation of Test Results for a description of Tube Result Report symbols.

INTERPRETATION OF TEST RESULTS

The **BD ProbeTec GC Q^x Amplified DNA Assay** uses fluorescent energy transfer as the detection method to test for the presence of *N. gonorrhoeae* in clinical specimens. All calculations are performed automatically by the **BD Viper** software.

The presence or absence of *N. gonorrhoeae* DNA is determined by calculating the peak fluorescence (MaxRFU) over the course of the amplification process and by comparing this measurement to a predetermined threshold value. The magnitude of the MaxRFU score is not indicative of the level of organism in the specimen. If the *N. gonorrhoeae*-specific signal is greater than or equal to a threshold of 125 MaxRFU, the EC fluorescence is ignored by the algorithm. If the *N. gonorrhoeae*-specific signal is less than a threshold of 125 MaxRFU, the EC fluorescence is utilized by the algorithm in the interpretation of the result. If assay control results are not as expected, patient results are not reported. See the Quality Control section for expected control values. Reported results are determined as follows.

Table 4: Interpretation of Test Results for the GC Q^x Assay

Tube Report Result	GC Q ^x MaxRFU	Report	Interpretation	Result
	≥125	<i>N. gonorrhoeae</i> DNA detected by SDA.	Positive for <i>N. gonorrhoeae</i> . <i>N. gonorrhoeae</i> organism viability and/or infectivity cannot be inferred since target DNA may persist in the absence of viable organisms.	Positive
	<125	<i>N. gonorrhoeae</i> DNA not detected by SDA.	Presumed negative for <i>N. gonorrhoeae</i> . A negative result does not preclude <i>N. gonorrhoeae</i> infection because results are dependent on adequate specimen collection, absence of inhibitors, and the presence of sufficient DNA to be detected.	Negative
	<125	Extraction control failure. Repeat test from initial specimen tube or obtain another specimen for testing.	<i>N. gonorrhoeae</i> , if present, is not detectable.	Extraction Control Failure
	Any value	Extraction Transfer Failure. Repeat test from initial specimen tube or obtain another specimen for testing.	<i>N. gonorrhoeae</i> , if present, is not detectable.	Extraction Transfer Failure
	Any value	Liquid Level Failure. Repeat test from initial specimen tube or obtain another specimen for testing.	<i>N. gonorrhoeae</i> , if present, is not detectable.	Liquid Level Failure
	Any value	Error. Repeat test from initial specimen tube or obtain another specimen for testing.	<i>N. gonorrhoeae</i> , if present, is not detectable.	Error

SPECIMEN PROCESSING CONTROLS

Specimen Processing Controls may be tested in accordance with the requirements of appropriate accrediting organizations. A positive Specimen Processing Control tests the entire assay system. For this purpose, known positive specimens can serve as controls by being processed and tested in conjunction with unknown specimens. Specimens used as processing controls must be stored, processed, and tested according to the package insert. If a known positive specimen is not available, additional options for Specimen Processing Controls are described below:

A. Preparation of Specimen Processing Controls in BD ProbeTec Q^x Swab Diluent

ATCC *Neisseria gonorrhoeae*:

Assay a stock culture of *N. gonorrhoeae* (ATCC # 19424) prepared as described below:

1. Thaw a vial of *N. gonorrhoeae* stock culture, received from ATCC and immediately inoculate a chocolate agar plate.
2. Incubate at 37°C in 3 – 5% CO₂ for 24 – 48 h.
3. Resuspend colonies from the chocolate agar plate with phosphate buffered saline (PBS).
4. Dilute cells in PBS to a 1.0 McFarland turbidity standard (approximately 3 x 10⁸ cells/mL).
5. Prepare 10-fold serial dilutions to a 10⁻⁵ dilution of the McFarland (at least 4 mL final volume) in PBS.
6. Place 0.1 mL of the 10⁻⁵ dilution in a **BD ProbeTec Q^x Swab Diluent** tube and tightly recap using a **black pierceable cap**.
7. Using the tube layout report, place the Specimen Processing Control(s) in order in the **BD Viper** Lysing Rack and lock into place.
8. Process the controls according to the Pre-warming Procedure and then follow the Test Procedure.

Bio-Rad AmpliTrol - *Chlamydia trachomatis* & *Neisseria gonorrhoeae*:

NOTE: Refer to manufacturer's processing instructions.

1. Add the appropriate volume of Bio-Rad AmpliTrol CT/GC to a **BD ProbeTec Q^x Swab Diluent** tube and tightly recap using a **black pierceable cap**.
2. Mix the solution by vortexing or with inversion.
3. Using the tube layout report, place the Specimen Processing Control(s) in order in the **BD Viper** Lysing Rack and lock into place.
4. Process the controls according to the Pre-warming Procedure and then follow the Test Procedure.

B. Preparation of Specimen Processing Controls in LBC Specimen Dilution Tubes

ATCC *Neisseria gonorrhoeae*

1. Grow *N. gonorrhoeae* culture overnight on chocolate agar plates.
2. Resuspend *N. gonorrhoeae* colonies in phosphate buffered saline (PBS).
3. Prepare a McFarland #1 turbidity standard from the resuspended colonies.
4. Prepare 10-fold serial dilutions of the McFarland #1 suspension to 10⁻⁵.
5. Add 0.1 mL of 10⁻⁵ dilution of *N. gonorrhoeae* to an LBC Specimen Dilution Tube containing 0.5 mL of **BD SurePath** Preservative Fluid or PreservCyt Solution. Tightly recap the LBC Specimen Dilution Tube using the blue pierceable cap.
6. Invert the LBC Specimen Dilution Tube 3 – 4 times to ensure that the contents are well mixed.
7. Using the tube layout report, place the Specimen Processing Control(s) in order in the **BD Viper** LBC Specimen Rack and lock into place.
8. Specimen Processing Controls are ready to be tested on the **BD Viper** System in Extracted Mode.
9. Change gloves prior to proceeding to avoid contamination.

ATCC *Chlamydia trachomatis* and *Neisseria gonorrhoeae*:

1. Thaw vial of *C. trachomatis* LGV II or serovar H cells received from ATCC.
2. Prepare 10-fold serial dilutions to 10⁻⁵ in PBS.
3. Grow *N. gonorrhoeae* culture overnight on chocolate agar plates.
4. Resuspend *N. gonorrhoeae* colonies in PBS.
5. Prepare a McFarland #1 turbidity standard from the resuspended colonies.
6. Prepare 10-fold serial dilutions of the McFarland #1 suspension to 10⁻⁵.
7. Add 0.1 mL of 10⁻⁵ dilution of *C. trachomatis* and 0.1 mL of 10⁻⁵ dilution of *N. gonorrhoeae* to an LBC Specimen Dilution Tube containing 0.5 mL of **BD SurePath** Preservative Fluid or PreservCyt Solution. Tightly recap the LBC Specimen Dilution Tube using the blue pierceable cap.
8. Invert the LBC Specimen Dilution Tube 3 – 4 times to ensure that the contents are well mixed.
9. Using the tube layout report, place the Specimen Processing Control(s) in order in the **BD Viper** LBC Specimen Rack and lock into place.
10. Specimen Processing Controls are ready to be tested on the **BD Viper** System in Extracted Mode.
11. Change gloves prior to proceeding to avoid contamination.

Bio-Rad AmpliTrol *Chlamydia trachomatis* and *Neisseria gonorrhoeae*:

NOTE: Refer to manufacturer's processing instructions.

1. Add the appropriate volume of Bio-Rad AmpliTrol CT/GC to an LBC Specimen Dilution Tube containing 0.5 mL of **BD SurePath** Preservative Fluid or PreservCyt Solution. Tightly recap the LBC Specimen Dilution Tube using the blue pierceable cap.
2. Invert the LBC Specimen Dilution Tube 3 – 4 times to ensure that the contents are well mixed.

- Using the tube layout report, place the Specimen Processing Control(s) in order in the **BD Viper** LBC Specimen Rack and lock into place.
- Specimen Processing Controls are ready to be tested on the **BD Viper** System in Extracted Mode.
- Change gloves prior to proceeding to avoid contamination.

MONITORING FOR THE PRESENCE OF DNA CONTAMINATION

At least monthly, the following test procedure should be performed to monitor the work area and equipment surfaces for the presence of DNA contamination. Environmental monitoring is essential to detect contamination prior to the development of a problem.

- For each area to be tested, use a clean collection swab from the **BD ProbeTec Q^x** Collection Kit for Endocervical or Lesion Specimens.
- Dip the swab into the **BD ProbeTec Q^x** Swab Diluent tube and wipe the first area* using a broad sweeping motion.
- Fully insert the collection swab into the Q^x Swab Diluent tube.
- Break the shaft of the swab at the score mark. Use care to avoid splashing of contents.
- Tightly recap the tube using the **black pierceable cap**.
- Repeat for each desired area.
- After all swabs have been collected, expressed in diluent, process according to the Pre-warming Procedure and then follow the Test Procedure.

*Recommended areas to test include: **Instrument deck:** Pipette Tip Station Covers (2); Tube Processing Station: Tube Alignment Block and Fixed Metal Base; Deck Waste Area, Priming and Warming Heaters/Stage; Extraction Block; Plate Sealing Tool; Tip Exchange Stations (2); **Instrument Exterior:** Upper Door Handle; Lower Door Handle; Waste Liquid Quick Release Valve; LCD Monitor (Touchscreen); Keyboard/Scanner; Staging Area; Locking Plate and Fixed Metal Base; **Accessories:** Tube Lockdown cover, **BD Viper** Lysing Rack/Table Base; **BD Viper** Lysing Heater; Metal Microwell Plates; Timer; Laboratory Bench Surfaces.

If an area gives a positive result or if contamination is suspected, clean the area with fresh 1% (v/v) sodium hypochlorite, DNA AWAY, or 3% (w/v) hydrogen peroxide. (Do not use hydrogen peroxide from a bottle that has remained open for longer than 8 days). Make sure the entire area is wetted with the solution and allowed to remain on the surface for at least 2 min or until dry. If necessary, remove excess cleaning solution with a clean towel. Wipe the area with a clean towel saturated with water and allow the surface to dry. Retest the area. Repeat cleaning process until negative results are obtained. If the contamination does not resolve, contact BD Technical Services for additional information.

LIMITATIONS OF THE PROCEDURE

- This method has been tested only with endocervical, vaginal, male urethral swab specimens, **BD SurePath** or PreservCyt specimens collected with cytobrush/spatula or broom device, and male and female urine specimens. Performance with other specimen types has not been assessed.
- Optimal performance of the test requires adequate specimen collection and handling. Refer to the "Specimen Collection and Transport" sections of this insert.
- Endocervical specimen adequacy can only be assessed by microscopic visualization of columnar epithelial cells in the specimen.
- Collection and testing of urine specimens with the **BD ProbeTec GC Q^x** Amplified DNA Assay is not intended to replace cervical exam and endocervical sampling for diagnosis of urogenital infection. Cervicitis, urethritis, urinary tract infections and vaginal infections may result from other causes or concurrent infections may occur.
- The **BD ProbeTec GC Q^x** Amplified DNA Assay for male and female urine specimen testing should be performed on first catch random urine specimens (defined as the first 20 – 60 mL of the urine stream).
- The effects of other potential variables such as vaginal discharge, use of tampons, douching, and specimen collection variables have not been determined.
- A negative test result does not exclude the possibility of infection because test results may be affected by improper specimen collection, technical error, specimen mix-up, concurrent antibiotic therapy, or the number of organisms in the specimen which may be below the sensitivity of the test.
- As with many diagnostic tests, results from the **BD ProbeTec GC Q^x** Amplified DNA Assay should be interpreted in conjunction with other laboratory and clinical data available to the physician.
- The **BD ProbeTec GC Q^x** Amplified DNA Assay should not be used for the evaluation of suspected sexual abuse or for other medico-legal indications. Additional testing is recommended in any circumstance when false positive or false negative results could lead to adverse medical, social, or psychological consequences.
- The **BD ProbeTec GC Q^x** Amplified DNA Assay cannot be used to assess therapeutic success or failure since nucleic acids from *N. gonorrhoeae* may persist following antimicrobial therapy.
- The **BD ProbeTec GC Q^x** Amplified DNA Assay provides qualitative results. No correlation can be drawn between the magnitude of the positive assay signal (MaxRFU) and the number of cells in an infected sample.
- The predictive value of an assay depends on the prevalence of the disease in any particular population. See Table 5 for hypothetical predictive values when testing varied populations.
- Because the Positive Control for the **BD ProbeTec CT/GC Q^x** Amplified DNA Assays is used in testing for both *C. trachomatis* and *N. gonorrhoeae*, correct positioning of the microwell strips is important for final results reporting.
- Use of the **BD ProbeTec GC Q^x** Amplified DNA Assay is limited to personnel who have been trained in the assay procedure and the **BD Viper** System.
- The reproducibility of the **BD ProbeTec GC Q^x** Amplified DNA Assay was established using seeded simulated swabs and seeded CT/GC Q^x Swab Diluent to simulate urine specimens. These specimens were inoculated with either *N. gonorrhoeae* alone or *N. gonorrhoeae* plus *C. trachomatis*.

16. Performance has not been established for urine specimens in Q^x UPT when fill volumes other than those falling within the purple lines on the fill window (approximately 2.0 mL to 3.0 mL) are used.
17. The **BD ProbeTec Neisseria gonorrhoeae** (GC) Q^x Amplified DNA Assay may cross-react with *N. cinerea* and *N. lactamica*. These organisms have only rarely been isolated from the genital tract.¹⁵⁻¹⁸ Refer to “Performance Characteristics” for further information.
18. The performance of the **BD ProbeTec GC Q^x Amplified DNA Assay** on the **BD Viper System** in extracted mode with swab specimens was evaluated for interference by blood, gynecological lubricants, and spermicides. The performance with urine specimens was evaluated for interference by blood and commonly used over-the-counter pain relievers. No interference was observed with any of the substances at the concentrations tested.
19. The patient-collected vaginal swab specimens are an option for screening women when a pelvic exam is not otherwise indicated.
20. The patient-collected vaginal swab specimen application is limited to healthcare facilities where support/ counseling is available to explain procedures and precautions.
21. The **BD ProbeTec GC Q^x Amplified DNA Assay** has not been validated for vaginal swab specimens collected by patients at home.
22. The performance of vaginal swab specimens has not been evaluated in patients less than 17 years of age.
23. The performance of vaginal swab specimens has not been evaluated in pregnant women.

EXPECTED RESULTS

NOTE: An explanation of symbols and abbreviations used in tables can be found in the Interpretation of Tables section (at end of insert).

A. Prevalence

The prevalence of positive *N. gonorrhoeae* specimens in patient populations depends upon: clinic type, age, risk factors, gender, and test method. The prevalence observed with the GC Q^x Amplified DNA Assay during a multi-center clinical trial for swab and urine specimens ranged from 1.4% to 19.2% for female specimens and 4.8% to 40.5% for male specimens (Table 10A).

The prevalence observed with the GC Q^x Assay during a multi-center clinical trial for **BD SurePath** specimens ranged from 0.0% to 25.9% (Table 10B). The prevalence observed with the GC Q^x Assay during a multi-center clinical trial for PreservCyt specimens ranged from 0.0% to 13.3% (Table 10C).

B. Positive and Negative Predictive Value

Hypothetical positive and negative predictive values (PPV & NPV) for the GC Q^x Assay with swab and urine specimens are shown in Table 5A. Hypothetical positive and negative predictive values (PPV & NPV) for the GC Q^x Assay from the multi-center clinical trial for **BD SurePath** specimens are shown in Table 5B. Hypothetical positive and negative predictive values (PPV & NPV) for the GC Q^x Assay from the multi-center clinical trial for PreservCyt specimens are shown in Table 5C. These calculations are based on hypothetical prevalence and overall sensitivity and specificity (compared to the patient infected status) of 99.3% and 99.3%, for swab and urine specimens, of 100.0% and 99.9% for **BD SurePath** specimens, and of 95.3% and 99.9% for PreservCyt specimens. In addition, PPV and NPV based on actual prevalence, sensitivity and specificity are shown in Tables 8 and 9. PPV was calculated using: (Sensitivity * Prevalence) / (Sensitivity * Prevalence + [1 - Specificity] * [1 - Prevalence]). NPV was calculated using: (Specificity * [1 - Prevalence]) / ([1-Sensitivity] * Prevalence + Specificity * [1-Prevalence]).

Table 5A: GC Hypothetical Positive and Negative Predictive Values (Swabs/Urines) Compared to Patient Infected Status

Prevalence (%)	Sensitivity (%)	Specificity (%)	PPV (%)	NPV (%)
2	99.3	99.3	74.3	100.0
5	99.3	99.3	88.2	100.0
10	99.3	99.3	94.0	99.9
20	99.3	99.3	97.3	99.8
30	99.3	99.3	98.4	99.7
40	99.3	99.3	99.0	99.5
50	99.3	99.3	99.3	99.3

Table 5B: GC Hypothetical Positive and Negative Predictive Values (BD SurePath) Compared to Patient Infected Status

Prevalence (%)	Sensitivity (%)	Specificity (%)	PPV (%)	NPV (%)
2	100.0	99.9	95.3	100.0
5	100.0	99.9	98.1	100.0
10	100.0	99.9	99.1	100.0
20	100.0	99.9	99.6	100.0
30	100.0	99.9	99.8	100.0
40	100.0	99.9	99.9	100.0
50	100.0	99.9	99.9	100.0

Table 5C: GC Hypothetical Positive and Negative Predictive Values (PreservCyt) Compared to Patient Infected Status

Prevalence (%)	Sensitivity (%)	Specificity (%)	PPV (%)	NPV (%)
2	95.3	99.95	97.5	99.9
5	95.3	99.95	99.0	99.8
10	95.3	99.95	99.5	99.5
20	95.3	99.95	99.8	98.8
30	95.3	99.95	99.9	98.0
40	95.3	99.95	99.9	97.0
50	95.3	99.95	99.9	95.5

C. MaxRFU Frequency Distribution

A total of 6284 GC Q^x Assay results from swab and urine specimens was evaluated at seven geographically diverse clinical sites. A frequency distribution of the initial MaxRFU values for the GC Q^x assay is shown in Figure A. The distribution of MaxRFU values from GC Q^x true positive, true negative, false positive and false negative specimens (i.e. from those specimens that yielded results which were discordant with the patient infected status [PIS]) is shown in Table 6A.

A total of 1715 GC Q^x Assay results from **BD SurePath** specimens was evaluated from eleven geographically diverse clinical sites. A frequency distribution of the initial MaxRFU values for the GC Q^x assay is shown in Figure B. The distribution of MaxRFU values from GC Q^x true positive, true negative, false positive and false negative specimens (i.e., from those specimens that yielded results which were discordant with the patient infected status [PIS]) is shown in Table 6B.

A total of 2074 GC Q^x Assay results from PreservCyt specimens was evaluated from eleven geographically diverse clinical sites. A frequency distribution of the initial MaxRFU values for the GC Q^x assay is shown in Figure C. The distribution of MaxRFU values from GC Q^x true positive, true negative, false positive and false negative specimens (i.e., from those specimens that yielded results which were discordant with the patient infected status [PIS]) is shown in Table 6C.

Figure A: Frequency Distribution of MaxRFU for the GC Q^x Assay (Swab and Urine Specimens)


Figure B: Frequency Distribution of MaxRFU for the GC Q^x Assay (BD SurePath Specimens)


Figure C: Frequency Distribution of MaxRFU for the GC Q^x Assay (PreservCyt Specimens)


Table 6A: GC Q^x MaxRFU Ranges for False Negative, False Positive, True Negative and True Positive Results (Swab/Urine Specimens)

MaxRFU Range		0 – 49	50 – 99	100 – 124	125 – 149	150 – 199	200 – 249	250 – 349	350 – 499	500 – 799	≥ 800
Total		5636	16	2	2	0	0	4	3	4	617
FN	FNU	2	0	0							
	FS	1	0	0							
	FUPT	1	0	0							
	Total	4	0	0							
FP	FNU				0	0	0	1	1	0	3
	FS				0	0	0	1	0	0	2
	FUPT				0	0	0	0	1	0	2
	FV				2	0	0	0	0	1	5
	MNU				0	0	0	1	0	1	5
	MS				0	0	0	0	0	0	6
	MUPT				0	0	0	0	1	0	5
	Total				2	0	0	3	3	2	28
TN	FNU	920	3	0							
	FS	918	5	1							
	FUPT	925	0	0							
	FV	913	6	1							
	MNU	655	0	0							
	MS	646	1	0							
	MUPT	655	1	0							
	Total	5632	16	2							
TP	FNU				0	0	0	0	0	0	63
	FS				0	0	0	0	0	0	64
	FUPT				0	0	0	0	0	0	64
	FV				0	0	0	1	0	0	64
	MNU				0	0	0	0	0	0	112
	MS				0	0	0	0	0	2	110
	MUPT				0	0	0	0	0	0	112
	Total				0	0	0	1	0	2	589

Table 6B: GC Q^x MaxRFU Ranges for False Negative, False Positive, True Negative and True Positive Results (BD SurePath Specimens)

MaxRFU Range	0 – 49	50 – 99	100 – 124	125 – 149	150 – 199	200 – 249	250 – 349	350 – 499	500 – 799	≥ 800
FN	0	0	0							
FP				0	0	0	0	0	0	2
TN	1659	2	1							
TP				0	0	0	0	0	0	51
Total	1659	2	1	0	0	0	0	0	0	53

Table 6C: GC Q^x MaxRFU Ranges for False Negative, False Positive, True Negative and True Positive Results (PreservCyt Specimens)

MaxRFU Range	0 – 49	50 – 99	100 – 124	125 – 149	150 – 199	200 – 249	250 – 349	350 – 499	500 – 799	≥ 800
FN	2	0	0							
FP				0	0	0	0	0	0	1
TN	2023	6	1							
TP				0	0	0	0	0	0	41
Total	2025	6	1	0	0	0	0	0	0	42

D. Controls

During the swab/urine clinical evaluation, there were no GC Q^x Positive Control failures from 253 GC Q^x plate runs. For the GC Q^x Negative Control, a failure was observed in 1 of 253 GC Q^x plate runs. During the **BD SurePath** specimen clinical evaluation, there was one GC Q^x Positive Control failure and no GC Q^x Negative Control failures from 120 GC Q^x plates that were run. During the PreservCyt specimen clinical evaluation, there were no GC Q^x Positive Control failures and one GC Q^x Negative Control failure from 142 GC Q^x plates that were run. The CT/GC Q^x Positive and Negative Control MaxRFU values observed in the clinical trials are shown in Table 7.

Table 7: Distribution of MaxRFU Results for the GC Q^x Assay Negative and Positive Controls

Control	Statistic	Swab and Urine Specimen Clinical Study	BD SurePath Specimen Clinical Study	PreservCyt Specimen Clinical Study
GC Q ^x Negative Control	n	252	120	141
MaxRFU	Maximum	17	42	10
	95th Percentile	7	0	0
	Median	0	0	0
	Mean	1	0	0
	5th Percentile	0	0	0
	Minimum	0	0	0
GC Q ^x Positive Control	n	253	120	142
MaxRFU	Maximum	2242	2156	2259
	95th Percentile	2083	1982	2045
	Median	1835	1786	1785
	Mean	1814	1777	1789
	5th Percentile	1502	1478	1555
	Minimum	530	1370	886

PERFORMANCE CHARACTERISTICS

Swab and Urine Specimen Clinical Study

Clinician-collected endocervical and male urethral swab specimens, patient-collected vaginal swab specimens (in a clinical setting), and male and female Q^x UPT and neat urine specimens were collected from 1059 symptomatic and asymptomatic female subjects and 787 symptomatic and asymptomatic male subjects attending OB/GYN, sexually transmitted disease (STD) and family planning clinics at seven geographically diverse clinical sites in North America. Subjects were classified as symptomatic if they reported symptoms such as dysuria, urethral discharge, coital pain/difficulty/bleeding, testicular or scrotum pain/swelling, abnormal vaginal discharge, or pelvic/uterine/adnexal pain. Subjects were classified as asymptomatic if they did not report symptoms. Sixty five female subjects and 13 male subjects were excluded from the data analysis due to age requirement violations, antibiotic treatment in the last 21 days, opting to withdraw from the study after initially consenting, failure to obtain paired swab and urine specimens, urine quantity less than 20 mL, or transport and storage errors related to specimen collection. Therefore, the final data analysis included 994 compliant female subjects and 774 compliant male subjects.

Five specimens were collected from each of the 994 eligible female subjects. A urine specimen was collected and split into Q^x UPT, neat urine and the two reference urine specimen collection devices followed by a vaginal swab specimen and three randomized endocervical swab specimens. Up to four specimens were collected from each of the 774 eligible male subjects. Up to three randomized urethral swab specimens were collected followed by a urine specimen that was split into Q^x UPT, neat urine and the two reference urine specimen collection devices. **BD ProbeTec GC Q^x** assay results were generated from the Q^x UPT and neat urine specimens, the vaginal swab specimen, one endocervical swab specimen and one male urethral swab specimen. The remaining two endocervical swab specimens, up to two male urethral swab specimens, and the two reference urine specimens for each male and female subject were tested using two reference methods: the **BD ProbeTec ET GC/AC** assay and another commercially available NAAT (Nucleic Acid Amplification Test). Specimen testing was conducted either at the site of collection or at a designated **BD Viper** testing site.

All performance calculations were based on the total number of **BD ProbeTec GC Q^x** assays results for endocervical, vaginal and male urethral swab specimens, and male and female Q^x UPT and neat urine specimens compared to a patient infected status (PIS) algorithm for each gender. In the algorithm, the designation of a subject as being infected with GC or not was based on endocervical swab and urine specimen results from the commercially available **BD ProbeTec ET GC/AC** assay and the other commercially available NAAT. Subjects were considered infected with GC if two of the four endocervical swab and urine specimens (or two of the three or four urethral swab and urine specimens) tested positive in the **BD ProbeTec ET GC/AC** assay and the other reference NAAT (one specimen testing positive in each NAAT). Subjects were considered non-infected if less than two reference NAAT results were positive. A total of 6284 **BD ProbeTec GC Q^x** assay results from symptomatic and asymptomatic male and female subjects were used to calculate sensitivity and specificity. Sensitivity and specificity by specimen type and symptomatic status are presented in Table 9A.

Performance of the assay with endocervical swabs, patient collected vaginal swab specimens (in a clinical setting), female UPT and neat urine was assessed in the clinical study. Separate performance was calculated for specimens collected from pregnant females. For the latter, sensitivity compared to patient infected status for FS, FV, FNU, and FUPT was 100% (3/3). In each case, specificity was 100% (24/24) for FS, FV, FNU, and FUPT separately.

Tables 11A and 11B summarize the number of results from symptomatic and asymptomatic subjects designated as infected or non-infected with GC according to the PIS algorithm.

NOTE: An explanation of symbols and abbreviations used in tables can be found in the Interpretation of Tables section (at end of insert).

BD SurePath Specimen Clinical Study

Endocervical swab specimens and **BD SurePath** specimens were collected from 1728 compliant female subjects attending family planning, OB/GYN, and sexually transmitted disease clinics at eleven geographically diverse clinical sites in North America. Subjects were classified as symptomatic if they reported symptoms such as dysuria, coital pain/difficulty/bleeding, abnormal vaginal discharge, or pelvic/uterine/adnexal pain. Subjects were classified as asymptomatic if they did not report symptoms. Thirteen subjects did not have a **BD SurePath** specimen result. Therefore there were 1715 subjects evaluated.

Three randomized endocervical swab specimens and a **BD SurePath** specimen were collected from each female subject. The three reference endocervical swabs were tested with the **BD ProbeTec ET CT/GC/AC** assay, the **BD ProbeTec GC Q^x** assay, and another commercially available NAAT (Nucleic Acid Amplification Test). Sensitivity and specificity for **BD SurePath** specimens were calculated by comparing results to a patient infected status (PIS) algorithm. The designation of positive or negative PIS was based on the endocervical swab specimen results from the three reference methods. At least two positive reference results were required to establish a subject as PIS-positive. At least two negative reference results were required to establish a subject as PIS-negative. The distribution of cervical sampling devices used in the clinical study according to clinical collection site is summarized in Table 8A. Sensitivity and specificity by symptomatic status are presented in Table 9B.

Table 11C summarizes the number of results from symptomatic and asymptomatic subjects designated as infected or non-infected with GC according to the PIS algorithm.

Table 12A summarizes the GC Q^x assay performance for **BD SurePath** specimens compared to PIS by clinic type.

Table 8A: Summary of Cervical Sampling Devices Used in the BD SurePath Specimen Clinical Study

Cervical Sampling Device Used	Clinical Collection Site Number											Total
	1	2	3	4	5	6	7	8	9	10	11	
Broom-Type Device	54	50	511	18	374	0	127	0	0	71	0	1205
Spatula/Cytobrush	0	25	0	0	182	112	32	24	103	8	37	523

PreservCyt Specimen Clinical Study

Endocervical swab specimens and PreservCyt specimens were collected from 2079 compliant female subjects attending family planning, OB/GYN, and sexually transmitted disease clinics at eleven geographically diverse clinical sites in North America. Subjects were classified as symptomatic if they reported symptoms such as dysuria, coital pain/difficulty/bleeding, abnormal vaginal discharge, or pelvic/uterine/adnexal pain. Subjects were classified as asymptomatic if they did not report symptoms. Two subjects were excluded due to an undetermined patient infected status. Three subjects did not have a PreservCyt specimen result. Therefore there were 2074 subjects evaluated.

Three randomized endocervical swab specimens and a PreservCyt specimen were collected from each female subject. The three reference endocervical swabs were tested with the **BD ProbeTec ET CT/GC/AC** assay, the **BD ProbeTec GC Q^x** assay, and another commercially available NAAT (Nucleic Acid Amplification Test). Sensitivity and specificity for PreservCyt specimens were calculated by comparing results to a patient infected status (PIS) algorithm. The designation of positive or negative PIS was based on the endocervical swab specimen results from the three reference methods. At least two positive reference results were required to establish a subject as PIS-positive. At least two negative reference results were required to establish a subject as PIS-negative. The distribution of cervical sampling devices used in the clinical study according to clinical collection site is summarized in Table 8B. Sensitivity and specificity by symptomatic status are presented in Table 9C.

Table 11D summarizes the number of results from symptomatic and asymptomatic subjects designated as infected or non-infected with GC according to the PIS algorithm.

Table 12B summarizes the GC Q^x assay performance for PreservCyt specimens compared to PIS by clinic type.

Table 8B: Summary of Cervical Sampling Devices Used in the PreservCyt Specimen Clinical Study

Cervical Sampling Device Used	Clinical Collection Site Number											Total
	1	2	3	4	5	6	7	8	9	10	11	
Broom-Type Device	89	0	0	45	16	464	272	83	0	99	0	1068
Spatula/Cytobrush	74	154	95	0	0	52	0	209	282	0	145	1011

Table 9A: GC Q^x Assay Performance for Swab and Urine Specimens Compared to Patient Infected Status (by symptomatic status)

Specimen Type	Symptomatic Status	n	Sensitivity	95% C.I.	Specificity	95% C.I.	PPV	NPV	Error Initial/Final
FS	A	450	96.3% (26/27)	(81.0% – 99.9%)	99.5% (421/423)	(98.3% – 99.9%)	92.5%	99.8%	3/0
	S	542	100.0% (38/38)	(90.7% – 100.0%)	99.8% (503/504)	(98.9% – 100.0%)	97.4%	100.0%	2/2
	Total	992	98.5% (64/65)	(91.7% – 100.0%)	99.7% (924/927)	(99.1% – 99.9%)	95.9%	99.9%	5/2
FV ¹	A	449	100.0% (27/27)	(87.2% – 100.0%)	98.6% (416/422)	(96.9% – 99.5%)	82.0%	100.0%	0/0
	S	544	100.0% (38/38)	(90.7% – 100.0%)	99.6% (504/506)	(98.6% – 100.0%)	95.0%	100.0%	0/0
	Total	993	100.0% (65/65)	(94.5% – 100.0%)	99.1% (920/928)	(98.3% – 99.6%)	88.5%	100.0%	0/0
FNU ²	A	450	96.3% (26/27)	(81.0% – 99.9%)	99.3% (420/423)	(97.9% – 99.9%)	89.8%	99.8%	0/0
	S	543	97.4% (37/38)	(86.2% – 99.9%)	99.6% (503/505)	(98.6% – 100.0%)	94.8%	99.8%	0/0
	Total	993	96.9% (63/65)	(89.3% – 99.6%)	99.5% (923/928)	(98.7% – 99.8%)	93.1%	99.8%	0/0
FUPT ³	A	450	100.0% (27/27)	(87.2% – 100.0%)	99.5% (421/423)	(98.3% – 99.9%)	92.7%	100.0%	0/0
	S	543	97.4% (37/38)	(86.2% – 99.9%)	99.8% (504/505)	(98.9% – 100.0%)	97.3%	99.8%	0/0
	Total	993	98.5% (64/65)	(91.7% – 100.0%)	99.7% (925/928)	(99.1% – 99.9%)	95.8%	99.9%	0/0
MS ⁴	A	508	100.0% (12/12)	(73.5% – 100.0%)	99.2% (492/496)	(97.9% – 99.8%)	75.5%	100.0%	0/0
	S	257	100.0% (100/100)	(96.4% – 100.0%)	98.7% (155/157)	(95.5% – 99.8%)	98.0%	100.0%	1/0
	Total	765	100.0% (112/112)	(96.8% – 100.0%)	99.1% (647/653)	(98.0% – 99.7%)	95.0%	100.0%	1/0
MNU ⁴	A	517	100.0% (12/12)	(73.5% – 100.0%)	99.2% (501/505)	(98.0% – 99.8%)	74.6%	100.0%	0/0
	S	257	100.0% (100/100)	(96.4% – 100.0%)	98.1% (154/157)	(94.5% – 99.6%)	97.1%	100.0%	0/0
	Total	774	100.0% (112/112)	(96.8% – 100.0%)	98.9% (655/662)	(97.8% – 99.6%)	93.9%	100.0%	0/0
MUPT ⁴	A	517	100.0% (12/12)	(73.5% – 100.0%)	99.2% (501/505)	(98.0% – 99.8%)	74.6%	100.0%	1/0
	S	257	100.0% (100/100)	(96.4% – 100.0%)	98.7% (155/157)	(95.5% – 99.8%)	98.0%	100.0%	0/0
	Total	774	100.0% (112/112)	(96.8% – 100.0%)	99.1% (656/662)	(98.0% – 99.7%)	95.0%	100.0%	1/0
Total		6284	99.3% (592/596)	(98.3% – 99.8%)	99.3% (5650/5688)	(99.1% – 99.5%)	93.7%	99.9%	7/2 ⁵

¹ Of the 994 female subjects enrolled in the study, one subject did not provide vaginal swab specimens.

² Of the 994 female subjects enrolled in the study, one neat urine specimen was excluded for noncompliant urine specimen storage.

³ Of the 994 female subjects enrolled in the study, one Q^x UPT urine specimen was excluded for noncompliant urine specimen storage.

⁴ Clinical Trial enrollment for asymptomatic male subjects was extended to obtain the total number of clinical positives for this sub-population.

⁵ Three liquid level errors, two extraction control failures, and one extraction transfer error were generated. Two of the three liquid level errors and the two extraction control failures resolved as negative and were included in the sensitivity and specificity calculations. The one liquid level error and one extraction transfer error failed to resolve and were not included in the sensitivity and specificity calculations.

Table 9B: GC Q^x Assay Performance for BD SurePath Specimens Compared to Patient Infected Status (by symptomatic status)

Symptomatic Status	n	Sensitivity	95% C.I.	Specificity	95% C.I.	PPV	NPV	Error Initial/Final
A	1157	100.0% (32/32)	(89.1% – 100.0%)	99.8% (1123/1125)	(99.4% – 100.0%)	93.5%	100.0%	2/0
S	558	100.0% (19/19)	(82.4% – 100.0%)	100.0% (539/539)	(99.3% – 100.0%)	100.0%	100.0%	0/0
Total	1715	100.0% (51/51)	(93.0% – 100.0%)	99.9% (1662/1664)	(99.6% – 100.0%)	96.90%	100.0%	2/0

Table 9C: GC Q^x Assay Performance for PreservCyt Specimens Compared to Patient Infected Status (by symptomatic status)

Symptomatic Status	n	Sensitivity	95% C.I.	Specificity	95% C.I.	PPV	NPV	Error Initial/Final
A	1349	92.3% (24/26)	(74.9% – 99.1%)	100.0% (1323/1323)	(99.7% – 100.0%)	100.0%	99.9%	1/0
S	725	100.0% (17/17)	(80.5% – 100.0%)	99.9% (707/708)	(99.2% – 100.0%)	95.9%	100.0%	0/0
Total	2074	95.3% (41/43)	(84.2% – 99.4%)	99.95% (2030/2031)	(99.7% – 100.0%)	100.0%	99.9%	1/0

Table 10A: GC Q^x Assay Performance for Swab and Urine Specimens Compared to Patient Infected Status (by clinical site)

Specimen Type	Collect Site	Prevalence	n	Sensitivity	95% C.I.	Specificity	95% C.I.	# CT (+) and GC (+)	PPV	NPV
FS ⁶	1	8.4%	155	100.0% (13/13)	(75.3% – 100.0%)	99.3% (141/142)	(96.1% – 100.0%)	5	92.9%	100.0%
	2	10.4%	154	93.8% (15/16)	(69.8% – 99.8%)	99.3% (137/138)	(96.0% – 100.0%)	6	94.0%	99.3%
	3	6.8%	73	100.0% (5/5)	(47.8% – 100.0%)	98.5% (67/68)	(92.1% – 100.0%)	2	82.9%	100.0%
	4	19.0%	105	100.0% (20/20)	(83.2% – 100.0%)	100.0% (85/85)	(95.8% – 100.0%)	6	100.0%	100.0%
	5	1.4%	70	100.0% (1/1)	(2.5% – 100.0%)	100.0% (69/69)	(94.8% – 100.0%)	0	100.0%	100.0%
	6	2.2%	365	100.0% (8/8)	(63.1% – 100.0%)	100.0% (357/357)	(99.0% – 100.0%)	3	100.0%	100.0%
	7	2.9%	70	100.0% (2/2)	(15.8% – 100.0%)	100.0% (68/68)	(94.7% – 100.0%)	0	100.0%	100.0%
FV ⁷	1	8.4%	155	100.0% (13/13)	(75.3% – 100.0%)	99.3% (141/142)	(96.1% – 100.0%)	5	92.9%	100.0%
	2	10.3%	155	100.0% (16/16)	(79.4% – 100.0%)	97.1% (135/139)	(92.8% – 99.2%)	6	79.8%	100.0%
	3	6.8%	73	100.0% (5/5)	(47.8% – 100.0%)	100.0% (68/68)	(94.7% – 100.0%)	2	100.0%	100.0%
	4	19.0%	105	100.0% (20/20)	(83.2% – 100.0%)	97.6% (83/85)	(91.8% – 99.7%)	6	90.7%	100.0%
	5	1.4%	70	100.0% (1/1)	(2.5% – 100.0%)	100.0% (69/69)	(94.8% – 100.0%)	0	100.0%	100.0%
	6	2.2%	365	100.0% (8/8)	(63.1% – 100.0%)	99.7% (356/357)	(98.4% – 100.0%)	3	88.2%	100.0%
	7	2.9%	70	100.0% (2/2)	(15.8% – 100.0%)	100.0% (68/68)	(94.7% – 100.0%)	0	100.0%	100.0%
FNU ⁸	1	8.4%	155	100.0% (13/13)	(75.3% – 100.0%)	98.6% (140/142)	(95.0% – 99.8%)	5	86.8%	100.0%
	2	10.3%	155	93.8% (15/16)	(69.8% – 99.8%)	97.8% (136/139)	(93.8% – 99.6%)	6	83.0%	99.3%
	3	6.8%	73	100.0% (5/5)	(47.8% – 100.0%)	100.0% (68/68)	(94.7% – 100.0%)	2	100.0%	100.0%
	4	19.2%	104	100.0% (20/20)	(83.2% – 100.0%)	100.0% (84/84)	(95.7% – 100.0%)	6	100.0%	100.0%
	5	1.4%	70	100.0% (1/1)	(2.5% – 100.0%)	100.0% (69/69)	(94.8% – 100.0%)	0	100.0%	100.0%
	6	2.2%	366	100.0% (8/8)	(63.1% – 100.0%)	100.0% (358/358)	(99.0% – 100.0%)	3	100.0%	100.0%
	7	2.9%	70	50.0% (1/2)	(1.3% – 98.7%)	100.0% (68/68)	(94.7% – 100.0%)	0	100.0%	98.5%
FUPT ⁹	1	8.4%	155	100.0% (13/13)	(75.3% – 100.0%)	99.3% (141/142)	(96.1% – 100.0%)	5	92.9%	100.0%
	2	10.3%	155	93.8% (15/16)	(69.8% – 99.8%)	99.3% (138/139)	(96.1% – 100.0%)	6	93.9%	99.3%
	3	6.8%	73	100.0% (5/5)	(47.8% – 100.0%)	100.0% (68/68)	(94.7% – 100.0%)	2	100.0%	100.0%
	4	19.2%	104	100.0% (20/20)	(83.2% – 100.0%)	98.8% (83/84)	(93.5% – 100.0%)	6	95.2%	100.0%
	5	1.4%	70	100.0% (1/1)	(2.5% – 100.0%)	100.0% (69/69)	(94.8% – 100.0%)	0	100.0%	100.0%
	6	2.2%	366	100.0% (8/8)	(63.1% – 100.0%)	100.0% (358/358)	(99.0% – 100.0%)	3	100.0%	100.0%
	7	2.9%	70	100.0% (2/2)	(15.8% – 100.0%)	100.0% (68/68)	(94.7% – 100.0%)	0	100.0%	100.0%

continued

Specimen Type	Collect Site	Prevalence	n	Sensitivity	95% C.I.	Specificity	95% C.I.	# CT (+) and GC (+)	PPV	NPV
MS ¹⁰	1	10.5%	313	100.0% (33/33)	(89.4% – 100.0%)	99.6% (279/280)	(98.0% – 100.0%)	11	96.7%	100.0%
	2	40.5%	79	100.0% (32/32)	(89.1% – 100.0%)	95.7% (45/47)	(85.5% – 99.5%)	10	94.1%	100.0%
	4	20.6%	170	100.0% (35/35)	(90.0% – 100.0%)	98.5% (133/135)	(94.8% – 99.8%)	11	94.5%	100.0%
	5	6.0%	182	100.0% (11/11)	(71.5% – 100.0%)	99.4% (170/171)	(96.8% – 100.0%)	5	91.4%	100.0%
	7	4.8%	21	100.0% (1/1)	(2.5% – 100.0%)	100.0% (20/20)	(83.2% – 100.0%)	0	100.0%	100.0%
MNU ¹¹	1	10.5%	313	100.0% (33/33)	(89.4% – 100.0%)	99.3% (278/280)	(94.7% – 99.9%)	11	94.4%	100.0%
	2	40.5%	79	100.0% (32/32)	(89.1% – 100.0%)	95.7% (45/47)	(85.5% – 99.2%)	10	94.1%	100.0%
	4	20.6%	170	100.0% (35/35)	(90.0% – 100.0%)	97.8% (132/135)	(93.6% – 99.5%)	11	92.2%	100.0%
	5	5.8%	191	100.0% (11/11)	(71.5% – 100.0%)	100.0% (180/180)	(98.0% – 100.0%)	5	100.0%	100.0%
	7	4.8%	21	100.0% (1/1)	(2.5% – 100.0%)	100.0% (20/20)	(83.2% – 100.0%)	0	100.0%	100.0%
MUPT ¹²	1	10.5%	313	100.0% (33/33)	(89.4% – 100.0%)	98.9% (277/280)	(96.9% – 99.8%)	11	91.4%	100.0%
	2	40.5%	79	100.0% (32/32)	(89.1% – 100.0%)	97.9% (46/47)	(88.7% – 99.9%)	10	97.0%	100.0%
	4	20.6%	170	100.0% (35/35)	(90.0% – 100.0%)	99.3% (134/135)	(95.9% – 100.0%)	11	97.4%	100.0%
	5	5.8%	191	100.0% (11/11)	(71.5% – 100.0%)	99.4% (179/180)	(96.9% – 100.0%)	5	91.1%	100.0%
	7	4.8%	21	100.0% (1/1)	(2.5% – 100.0%)	100.0% (20/20)	(83.2% – 100.0%)	0	100.0%	100.0%

⁶ 22 of the 65 FS PIS positive subjects were co-infected with CT.

⁷ 22 of the 65 FV PIS positive subjects were co-infected with CT.

⁸ 22 of the 65 FNU PIS positive subjects were co-infected with CT.

⁹ 22 of the 65 FUPT PIS positive subjects were co-infected with CT.

¹⁰ 37 of the 112 MS PIS positive subjects were co-infected with CT.

¹¹ 37 of the 112 MNU PIS positive subjects were co-infected with CT.

¹² 37 of the 112 MUPT PIS positive subjects were co-infected with CT.

Table 10B: GC Q^x Assay Performance for BD SurePath Specimens Compared to Patient Infected Status (by clinical site)

Collection Site	Prevalence	n	Sensitivity	95% C.I.	Specificity	95% C.I.	# CT (+) and GC (+)	PPV	NPV
1	10.8%	74	100.0% (8/8)	(63.1% – 100.0%)	100.0% (66/66)	(94.6% – 100.0%)	7	100.0%	100.0%
2	3.9%	103	100.0% (4/4)	(39.8% – 100.0%)	100.0% (99/99)	(96.3% – 100.0%)	1	100.0%	100.0%
3	0.0%	37	NA	NA	100.0% (37/37)	(90.5% – 100.0%)	0	NA	NA
4	25.9%	54	100.0% (14/14)	(76.8% – 100.0%)	97.5% (39/40)	(86.8% – 99.9%)	4	93.3%	100.0%
5	4.3%	69	100.0% (3/3)	(29.2% – 100.0%)	100.0% (66/66)	(94.6% – 100.0%)	1	100.0%	100.0%
6	1.6%	555	100.0% (9/9)	(66.4% – 100.0%)	99.8% (545/546)	(99.0% – 100.0%)	2	89.0%	100.0%
7	2.0%	511	100.0% (10/10)	(69.2% – 100.0%)	100.0% (501/501)	(99.3% – 100.0%)	5	100.0%	100.0%
8	1.3%	159	100.0% (2/2)	(15.8% – 100.0%)	100.0% (157/157)	(97.7% – 100.0%)	2	100.0%	100.0%
9	0.0%	112	NA	NA	100.0% (112/112)	(96.8% – 100.0%)	0	NA	NA
10	5.6%	18	100.0% (1/1)	(2.5% – 100.0%)	100.0% (17/17)	(80.5% – 100.0%)	0	100.0%	100.0%
11	0.0%	23	NA	NA	100.0% (23/23)	(85.2% – 100.0%)	0	NA	NA

Table 10C: GC Q^x Assay Performance for PreservCyt Specimens Compared to Patient Infected Status (by clinical site)

Collection Site	Prevalence	n	Sensitivity	95% C.I.	Specificity	95% C.I.	# CT (+) and GC (+)	PPV	NPV
1	5.5%	163	88.9% (8/9)	(51.8% – 99.7%)	100.0% (154/154)	(97.6% – 100.0%)	5	100.0%	99.4%
2	5.2%	154	100.0% (8/8)	(63.1% – 100.0%)	99.3% (145/146)	(96.2% – 100.0%)	1	88.7%	100.0%
3	3.2%	95	100.0% (3/3)	(29.2% – 100.0%)	100.0% (92/92)	(96.1% – 100.0%)	2	100.0%	100.0%
4	13.3%	45	100.0% (6/6)	(54.1% – 100.0%)	100.0% (39/39)	(91.0% – 100.0%)	2	100.0%	100.0%
5	0.0%	16	NA	NA	100.0% (16/16)	(79.4% – 100.0%)	0	NA	NA
6	1.6%	516	100.0% (8/8)	(63.1% – 100.0%)	100.0% (508/508)	(99.3% – 100.0%)	2	100.0%	100.0%
7	2.9%	272	87.5% (7/8)	(47.3% – 99.7%)	100.0% (264/264)	(98.6% – 100.0%)	3	100.0%	99.6%
8	0.0%	292	NA	NA	100.0% (292/292)	(98.7% – 100.0%)	0	NA	NA
9	0.0%	282	NA	NA	100.0% (282/282)	(98.7% – 100.0%)	0	NA	NA
10	0.0%	97	NA	NA	100.0% (97/97)	(96.3% – 100.0%)	0	NA	NA
11	0.7%	142	100.0% (1/1)	(2.5% – 100.0%)	100.0% (141/141)	(97.4% – 100.0%)	0	100.0%	100.0%

Table 11A: Analysis of GC Positive/Negative Swab and Urine Specimens from Female Subjects Based on Patient Infected Status

PIS GC	NAAT 1		NAAT 2		BD ProbeTec GC Qx Amplified DNA Assay				Symptomatic Status		
	Endocervical Swab	Urine	Endocervical Swab	Urine	Qx Endocervical Swab	Qx Vaginal Swab	Neat Urine	Qx UPT Urine	A	S	Total
	+	-	+	+	+	-	+	+	+	1	0
+		-	+	-	+	+	-	-	0	1	1
+		-	+	-	+	+	+	+	3	0	3
+		-	+	+	+	+	+	+	1	1	2
+		+	+	-	+	+	+	+	2	1	3
+		+	+	+	+	+	-	+	1	0	1
+		+	+	+	+	+	+	+	19	35	54
Total PIS Positive									27	38	65
-	NA	-	-	-	-	-	-	-	12	2	14
	-	NA	E	-	-	-	NA	NA	0	1	1
	-	NA	-	-	-	-	-	-	1	1	2
	-	I	-	-	-	-	-	-	5	1	6
	-	-	NA	-	-	-	-	-	1	2	3
	-	-	E	-	-	-	-	-	1	0	1
	-	-	-	-	ET	-	-	-	0	1	1
	-	-	-	-	LE	-	-	-	0	1	1
	-	-	-	-	-	NA	-	-	1	0	1
	-	-	-	-	-	-	-	-	390	484	874
	-	-	-	-	-	-	-	+	0	1	1
	-	-	-	-	-	-	+	-	1	1	2
	-	-	-	-	-	-	+	-	4	1	5
	-	-	-	-	-	-	+	+	0	1	1
	-	-	-	-	-	-	+	+	1	0	1
	-	-	-	-	+	-	-	-	0	1	1
	-	-	+	-	-	-	-	-	1	3	4
	-	-	+	-	+	-	-	-	1	0	1
	-	+	-	-	-	-	-	-	1	2	3
+	-	-	-	-	-	-	-	2	3	5	
+	+	-	-	-	+	+	+	1	0	1	
Total PIS Negative									423	506	929

I = Indeterminate

LE = Liquid Level Error

Table 11B: Analysis of GC Positive/Negative Specimens from Male Subjects Based on Patient Infected Status

PIS GC	NAAT 1		NAAT 2		BD ProbeTec GC Q ^x Amplified DNA Assay			Symptomatic Status		
	Urethral Swab	Urine	Urethral Swab	Urine	Q ^x Urethral Swab	Neat Urine	Q ^x UPT Urine	A	S	Total
	+	+	+	+	+	+	+	+	11	81
+		+	NA	+	+	+	+	1	13	14
NA		+	+	+	+	+	+	0	6	6
Total PIS Positive								12	100	112
-	-	I	-	-	-	-	-	4	1	5
	-	I	NA	-	-	-	-	1	0	1
	-	-	E	-	-	-	-	2	0	2
	-	-	-	E	-	-	-	0	1	1
	-	-	-	-	NA	-	-	9	0	9
	-	-	-	-	-	-	-	422	124	546
	-	-	-	-	-	-	+	2	1	3
	-	-	-	-	-	+	-	1	1	2
	-	-	-	-	-	+	+	1	0	1
	-	-	-	-	+	-	-	3	0	3
	-	-	-	+	-	-	-	2	1	3
	-	-	+	-	-	-	-	2	1	3
	-	-	+	+	+	+	-	0	1	1
	-	-	NA	-	-	-	-	29	11	40
	-	+	-	-	-	-	-	1	0	1
	-	NA	-	-	-	-	-	1	0	1
	+	-	-	-	-	-	-	0	1	1
	+	+	NA	-	-	-	-	0	1	1
	NA	-	-	-	-	-	-	22	11	33
	NA	-	-	-	-	-	+	1	0	1
NA	-	+	-	-	-	-	1	0	1	
NA	-	+	+	+	+	+	1	1	2	
NA	+	-	-	-	-	-	0	1	1	
Total PIS Negative								505	157	662

Table 11C: Analysis of GC Positive/Negative BD SurePath Specimens Based on Patient Infected Status

PIS GC	NAAT 1	NAAT 2	NAAT 3	BD ProbeTec GC Q ^x Amplified DNA Assay	Symptomatic Status		
	Swab	Swab	Swab	BD SurePath	A	S	Total
+	-	+	+	+	0	1	1
	+	-	+	+	1	1	2
	+	+	+	+	31	17	48
Total PIS Positive					32	19	51
-	-	-	+	+	1	0	1
	-	+	-	+	1	0	1
	-	I	-	-	2	2	4
	-	-	NA	-	6	1	7
	-	-	-	-	1103	531	1634
	-	-	+	-	6	1	7
	-	+	-	-	5	3	8
	+	-	-	-	1	1	2
Total PIS Negative					1125	539	1664

Table 11D: Analysis of GC Positive/Negative PreservCyt Specimens Based on Patient Infected Status

PIS GC	NAAT 1	NAAT 2	NAAT 3	BD ProbeTec GC Q ^x Amplified DNA Assay	Symptomatic Status		
	Swab	Swab	Swab	PreservCyt	A	S	Total
+	NA	+	+	+	1	3	4
	+	-	+	-	1	0	1
	+	-	+	+	1	0	1
	+	+	NA	+	1	0	1
	+	+	+	-	1	0	1
	+	+	+	+	21	14	35
Total PIS Positive					26	17	43
-	NA	-	-	-	181	79	260
	-	I	-	-	1	0	1
	-	-	NA	-	3	0	3
	-	-	LE	-	2	0	2
	-	-	-	-	1129	624	1753
	-	-	-	+	0	1	1
	-	-	+	-	2	0	2
	-	+	-	-	4	3	7
Total PIS Negative					1323	708	2031

Table 12A: GC Q^x Assay Performance for BD SurePath Specimens Compared to Patient Infected Status (by clinic type)

Clinic Type	Prevalence	n	Sensitivity	95% C.I.	Specificity	95% C.I.	PPV	NPV
Family Planning	1.4%	844	100.0% (12/12)	(73.5% – 100.0%)	99.9% (831/832)	(99.3% – 100.0%)	93.4%	100.0%
OB/GYN	1.8%	548	100.0% (10/10)	(69.2% – 100.0%)	100.0% (538/538)	(99.3% – 100.0%)	100.0%	100.0%
STD	9.0%	323	100.0% (29/29)	(88.1% – 100.0%)	99.7% (293/294)	(98.1% – 100.0%)	97.1%	100.0%

Table 12B: GC Q^x Assay Performance for PreservCyt Specimens Compared to Patient Infected Status (by clinic type)

Clinic Type	Prevalence	n	Sensitivity	95% C.I.	Specificity	95% C.I.	PPV	NPV
Family Planning	0.7%	1187	100.0% (8/8)	(63.1% – 100.0%)	100.0% (1179/1179)	(99.7% – 100.0%)	100.0%	100.0%
OB/GYN	3.0%	367	90.9% (10/11)	(58.7% – 99.8%)	100.0% (356/356)	(99.0% – 100.0%)	100.0%	99.7%
STD	4.6%	520	95.8% (23/24)	(78.9% – 99.9%)	99.8% (495/496)	(98.9% – 100.0%)	95.9%	99.8%

GC Q^x Assay Analytical Sensitivity:

The Limits of Detection (LODs) for the GC Q^x Assay with *Neisseria gonorrhoeae* strain ATCC 19424 in urine and swab specimens when extracted on the **BD Viper** System were determined to be < 50 cells per mL for neat and Q^x UPT urine and < 100 GC cells per mL for expressed vaginal, endocervical swab, **BD SurePath** and PreservCyt specimens.

The GC Q^x Assay on the **BD Viper** System in extracted mode was able to detect 17 GC strains (ATCC 19424, 27628, 27629, 27630, 27632, 27633, 27631, 21823, 51803, 23051, 31407, 31953, 35201, 31397, 31151, 43785, 51804) with ≥ 95% proportion positive at a concentration of 50 cells per mL in CT/GC Q^x Swab Diluent, in **BD SurePath** Preservative Fluid in LBC Specimen Dilution Tubes, and in PreservCyt Solution in LBC Specimen Dilution Tubes.

GC Q^x Assay Analytical Specificity:

DNA from 141 organisms listed in Table 13 was extracted on the **BD Viper** System and tested with the **BD ProbeTec** GC Q^x Amplified DNA Assay. All potential cross-reactive species were tested at ≥ 1x10⁸ cells/mL except where noted. Two *N. cinerea* and two *N. lactamica* strains were shown to cross-react in the GC Q^x assay.

Table 13: Potential Cross-reacting Microorganisms

<i>Acinetobacter calcoaceticus</i>	<i>Enterococcus faecium</i>	<i>Peptostreptococcus asaccharolyticus</i>	<i>Neisseria elongata</i> subsp. <i>glycolytica</i>
<i>Acinetobacter Iwoffii</i>	Epstein Barr Virus***	<i>Peptostreptococcus productus</i>	<i>Neisseria elongata</i> subsp. <i>nitroreducens</i> (2)
<i>Actinomyces israelii</i>	<i>Escherichia coli</i>	<i>Plesiomonas shigelloides</i>	<i>Neisseria elongata</i>
Adenovirus***	<i>Flavobacterium meningosepticum</i>	<i>Propionibacterium acnes</i>	<i>Neisseria flava</i> (4)
<i>Aeromonas hydrophilia</i>	<i>Gardnerella vaginalis</i>	<i>Providencia stuartii</i>	<i>Neisseria flavescens</i> (4)
<i>Alcaligenes faecalis</i> *	<i>Gemella haemolyans</i>	<i>Pseudomonas aeruginosa</i>	<i>Neisseria lactamica</i> (7)
<i>Bacillus subtilis</i> *	<i>Haemophilus influenzae</i>	<i>Salmonella minnesota</i>	<i>Neisseria meningitidis</i> (12)
<i>Bacteroides fragilis</i>	Herpes Simplex Virus **	<i>Salmonella typhimurium</i>	<i>Neisseria mucosa</i> (5)
<i>Candida albicans</i> *	Human papillomavirus (16 and 18)***	<i>Staphylococcus aureus</i>	<i>Neisseria perflava</i> (8)
<i>Candida glabrata</i> *	<i>Kingella kingae</i>	<i>Staphylococcus epidermidis</i>	<i>Neisseria polysaccharea</i> (2)
<i>Candida tropicalis</i> *	<i>Klebsiella pneumoniae</i>	<i>Streptococcus agalactiae</i>	<i>Neisseria sicca</i> (5)
<i>Chlamydia trachomatis</i>	<i>Lactobacillus acidophilus</i> *	<i>Streptococcus mitis</i>	<i>Neisseria subflava</i> (15)
<i>Chlamydia pneumoniae</i> ****	<i>Lactobacillus brevis</i>	<i>Streptococcus mutans</i>	<i>Neisseria weaverii</i> (3)
<i>Chlamydia psittaci</i> *	<i>Lactobacillus jensenii</i> *	<i>Streptococcus pneumoniae</i> *	
<i>Citrobacter freundii</i>	<i>Listeria monocytogenes</i>	<i>Streptococcus pyogenes</i>	
<i>Clostridium perfringens</i>	<i>Mobiluncus mulieris</i>	<i>Streptomyces griseus</i> **	
<i>Corynebacterium renale</i>	<i>Moraxella lacunata</i> *	<i>Trichomonas vaginalis</i> **	
<i>Cryptococcus neoformans</i> *	<i>Moraxella osloensis</i>	<i>Veillonella parvula</i>	
Cytomegalovirus**	<i>Morganella morganii</i>	<i>Vibrio parahaemolyticus</i>	
<i>Edwardsiella tarda</i>	<i>Mycobacterium gordonae</i>	<i>Yersinia enterocolitica</i>	
<i>Enterobacter cloacae</i>	<i>Mycobacterium smegmatis</i>	<i>Branhamella catarrhalis</i> (5)	
<i>Enterococcus faecalis</i>	<i>Peptostreptococcus anaerobius</i>	<i>Neisseria cinerea</i> (2)	

(n) number of strains tested in the **BD ProbeTec** GC Q^x Assay

* Tested at > 1x10⁷ cells or EB/mL; **Tested at > 1x10⁶ cells or viral particles per mL; ***Tested at ≥ 1x10⁶ genomic equivalents per mL;**** tested at ≥ 1x10⁵ TCID₅₀/mL

GC Q^x Interfering Substances

The performance of the **BD ProbeTec GC Q^x Assay** on the **BD Viper System** in extracted mode was evaluated in the presence of potential interfering substances which may be encountered in swab, urine, **BD SurePath** and/or **PreservCyt** specimens. Potential interfering substances were spiked into Q^x UPT urine, vaginal swab, specimen matrices, **BD SurePath** specimens in LBC Specimen Dilution Tubes, and **PreservCyt** specimens in LBC Specimen Dilution Tubes in both the presence and the absence of GC organisms (150 GC cells/mL in urine matrix and 300 GC cells/mL in Swab/LBC Specimen Dilution Tube matrix). Results are summarized in Table 14.

Table 14: GC Q^x Interfering Substances

Interpretation	Swab	Urine	BD SurePath	PreservCyt
No interference observed	Blood (≤ 60%) Seminal fluid Mucus Over the counter vaginal products and contraceptives Hemorrhoidal cream Prescription vaginal treatments Leukocytes (1x10 ⁶ cells/mL) 1x10 ⁶ EB/mL <i>Chlamydia trachomatis</i>	Blood (≤ 1%) Seminal fluid Mucus Antibiotics Analgesics Phenazopyridine Over the counter deodorant sprays and powders Hormones Leukocytes Albumin <1 mg/mL Glucose Acidic urine (pH 4.0) Alkaline urine (pH 9.0) Bilirubin 1x10 ⁶ EB/mL <i>Chlamydia trachomatis</i> Organisms associated with urinary tract infections	Blood (≤ 1%) Seminal fluid Mucus Over the counter vaginal products and contraceptives Hemorrhoidal cream Prescription vaginal treatments Leukocytes (1x10 ⁶ cells/mL) 1x10 ⁶ EB/mL <i>Chlamydia trachomatis</i>	Blood (≤ 1%) Seminal fluid Mucus Over the counter vaginal products and contraceptives Hemorrhoidal cream Prescription vaginal treatments Leukocytes (1x10 ⁶ cells/mL) 1x10 ⁶ EB/mL <i>Chlamydia trachomatis</i>
May cause extraction control (EC) failures	Blood (> 60%)	Not applicable	Not applicable	Glacial Acetic Acid + Blood (≤ 5%/1% V/V)
May cause false negative results	Not applicable	Not applicable	Not applicable	Glacial Acetic Acid + Blood (≤ 5%/1% V/V)

Neat and Q^x UPT Urine Stability

Pools of GC negative male and female urine specimens were used in analytical experiments to support the urine storage and transport stability claims. For neat urine, pools were co-spiked with CT serovar H and GC strain ATCC 19424 at 45 EB per mL and 150 cells per mL, respectively. Neat urine specimens were stored at either 2 – 8°C for 1, 3 or 7 days; or at 30°C for 8, 24 or 30 h; or at -20°C for 180 days. At each time point, samples were removed from storage and tested with the **BD ProbeTec GC Q^x Assay** on the **BD Viper System** in extracted mode. Thirty-two assay replicates were generated for each condition (sample type/temperature/duration). The expected results were obtained with the GC Q^x assay under all conditions tested.

For Q^x UPT urine, pooled specimens were co-spiked with CT serovar H and GC strain ATCC 19424 at 45 EB per mL and 150 cells per mL, respectively. The spiked urine specimen pools were then stored at either 2 – 8°C for 24 h or 30°C for 8 h prior to transfer into Q^x UPT tubes. The Q^x UPT specimens were then stored either at 2 – 8°C for 14, 21 or 30 days; or at 30°C for 14, 21 or 30 days; or at -20°C for 180 days. At each time point Q^x UPT specimens were removed from storage and tested with the **BD ProbeTec GC Q^x Assay** on the **BD Viper System** in extracted mode. Thirty-two assay replicates were generated for each condition (sample type/temperature/duration). The expected results were obtained with the GC Q^x assay under all conditions tested.

Vaginal Dry and Expressed Swab Stability

Pools of GC negative vaginal swab matrix were used in analytical experiments to support the storage and transport stability claims for dry vaginal swab specimens. Pools were co-spiked with CT serovar H and GC strain ATCC 19424 to achieve 90 EB per mL and 300 cells per mL, respectively, when seeded onto swabs and expressed in CT/GC Q^x Swab Diluent. Seeded dry swabs were stored at 2 – 8°C for 3, 7, or 14 days; or at 30°C for 3, 7 or 14 days; or at -20°C for 30, 60 or 180 days. At each time point, dry swabs were removed from storage and expressed into 2 mL of CT/GC Q^x Swab Diluent and evaluated with the **BD ProbeTec GC Q^x Assay** on the **BD Viper System** in extracted mode. Thirty-two assay replicates were generated for each condition (sample type/temperature/duration). The expected results were obtained with the GC Q^x assay under all conditions tested.

Pools of GC negative vaginal swab matrix were used in analytical experiments to support the storage and transport stability claims for expressed vaginal swab specimens. Pools were spiked with CT serovar H and GC strain ATCC 19424 to achieve 90 EB per mL and 300 cells per mL, respectively. The spiked swab matrix was stored at 2 – 8°C for 7, 14 or 30 days; or at 30°C for 7, 14 or 30 days; or at -20°C for 30, 60 or 180 days. At each time point, samples were removed from storage and tested with the **BD ProbeTec GC Q^x Assay** on the **BD Viper System** in extracted mode. Thirty-two assay replicates were generated for each condition (sample type/temperature/duration). The expected results were obtained with the GC Q^x assay under all conditions tested.

Endocervical and Urethral Swab Specimen Stability

Pools of GC negative endocervical swab matrix were used in analytical experiments to support the storage and transport stability claims for endocervical and urethral swab specimens. Pools of swab matrix were spiked with CT serovar H and GC strain ATCC 19424 at 90 EB per mL and 300 cells per mL, respectively. The pools were dispensed in 2 mL volumes into BD sample tubes to simulate “wet” endocervical specimens and stored at either 2 – 8°C for 7, 14 or 30 days; or at 30°C for 7, 14 or 30 days; or at -20°C for 30, 60 or 180 days. At each time point, samples were removed from storage and tested with the **BD ProbeTec GC Q^x Assay** on the **BD Viper System** in extracted mode. Thirty-two assay replicates were generated for each condition (sample type/temperature/duration). The expected results were obtained with the GC Q^x assay under all conditions tested.

Post Pre-warm Specimen Stability

Pools of male and female GC negative neat urine specimens were used in analytical experiments to support the storage stability claims for pre-warmed neat and Q^x UPT urine specimens. Pooled specimens were spiked with CT serovar H and GC strain ATCC 19424 at 45 EB per mL and 150 cells per mL, respectively and either added to Q^x UPT tubes or left untreated as neat urine. Both specimen types were pre-warmed at 114°C for 15 min, and cooled for 15 min. After the pre-warm process, specimen tubes were stored at either 2 – 8°C for 1, 3 or 7 days; or at 30°C for 1, 3 or 7 days; or at -20°C for 30 or 180 days. At each time point samples were removed from storage and tested with the **BD ProbeTec GC Q^x Assay** on the **BD Viper System** in extracted mode. Thirty-two assay replicates were generated for each condition (sample type/temperature/duration). The expected results were obtained with the GC Q^x assay under all conditions tested.

Pools of GC negative vaginal and endocervical swab specimen matrices in CT/GC Q^x Swab Diluent were used in analytical experiments to support the storage stability claims for pre-warmed expressed vaginal, endocervical, and male urethral swab specimens. For both types of matrix, pooled specimens were spiked with CT serovar H and GC strain ATCC 19424 at 90 EB per mL and 300 cells per mL, respectively and aliquotted into 2 mL volumes in BD specimen tubes. The tubes were pre-warmed at 114°C for 15 min and cooled for 15 min. After the pre-warm process, the specimen tubes were stored either at 2 – 8°C for 3 or 7 days; or at 30°C for 3 or 7 days; or at -20°C for 30 or 180 days. At each time point, samples were removed from storage and tested with the **BD ProbeTec GC Q^x Assay** on the **BD Viper System** in extracted mode. Thirty-two assay replicates were generated for each condition (sample type/temperature/duration). The expected results were obtained with the GC Q^x assay under all conditions tested.

BD SurePath Specimen Stability

Pools of CT and GC negative **BD SurePath** clinical specimens were used in analytical experiments to support the storage and stability claims. Pools were co-spiked with CT serovar H and GC strain ATCC 19424 to achieve 90 EB per mL and 300 cells per mL, respectively. The pools were dispensed in 10 mL volumes in **BD SurePath** vials and stored at either 2 – 8°C or 30°C. After 30 days, 0.5 mL from each vial was removed and added to an LBC Specimen Dilution Tube. The specimens in the LBC Specimen Dilution Tube were then stored at 2 – 8°C for 30 days; or at 30°C for 30 days; or at -20°C for 90 days. At each time point, samples were removed from storage and tested with the **BD ProbeTec GC Q^x Assay** on the **BD Viper System** in extracted mode. Twenty-four assay replicates were generated for each condition (temperature/duration). The expected results were obtained with the GC Q^x assay under all conditions tested.

PreservCyt Specimen Stability

Pools of CT and GC negative PreservCyt clinical specimens were used in analytical experiments to support the storage and stability claims. Pools were co-spiked with CT serovar H and GC strain ATCC 19424 to achieve 90 EB per mL and 300 cells per mL, respectively. The pools were dispensed in 20 mL volumes in PreservCyt vials and stored at either 2 – 8°C or 30°C. After 30 days, 0.5 mL from each vial was removed and added to an LBC Specimen Dilution Tube. The specimens in the LBC Specimen Dilution Tube were then stored at 2 – 8°C for 30 days; or at 30°C for 30 days; or at -20°C for 90 days. At each time point, samples were removed from storage and tested with the **BD ProbeTec GC Q^x Assay** on the **BD Viper System** in extracted mode. Twenty-four assay replicates were generated for each condition (temperature/duration). The expected results were obtained with the GC Q^x assay under all conditions tested.

Reproducibility

Reproducibility of the **BD Viper System** using the **BD ProbeTec GC Q^x Assay** was evaluated at three clinical sites on one **BD Viper System** per site. A panel of simulated specimens was tested that comprised CT and GC organisms seeded into swab diluent for the **BD ProbeTec GC Q^x Assay**. Simulated endocervical and urethral specimens contained a clean endocervical swab whereas the simulated urine and vaginal swab specimens did not. Uninoculated swab diluent for the **BD ProbeTec GC Q^x Assay** was used for the GC negative samples. Nine replicates of each panel member were tested every day for five days on each **BD Viper System**. The data are summarized in Table 15A.

Table 15A: Summary of Reproducibility Data for Swab and Urine Specimens on the BD Viper System for the GC Q^x Assay

						Within Run		Between Runs Within Site		Between Site	
Specimen Type	CT EB's/mL	GC Cells/mL	% Correct	95% CI	MaxRFU Mean	SD	%CV	SD	%CV	SD	%CV
Endocervical/ Urethral	0	0	99.3% (134/135)	(95.9%, 100.0%)	13.8	151.3	1096.3	0.0	0.0	0.6	4.3
	30	0	98.5% (133/135)	(94.8%, 99.8%)	28.1	220.7	785.3	0.0	0.0	33.8	120.3
	0	100	100.0% (135/135)	(97.3%, 100.0%)	1859.5	94.1	5.1	0.0	0.0	19.2	1.0
	30	250	100.0% (135/135)	(97.3%, 100.0%)	1847.3	117.6	6.4	0.0	0.0	25.9	1.4
	75	100	100.0% (135/135)	(97.3%, 100.0%)	1855.9	119.4	6.4	0.0	0.0	42.2	2.3
Urine/Vaginal	0	0	99.3% (134/135)	(95.9%, 100.0%)	15.7	162.3	1031.1	0.0	0.0	0.0	0.0
	30	0	100.0% (135/135)	(97.3%, 100.0%)	1.1	3.1	295.8	0.7	69.7	0.5	48.3
	0	100	100.0% (135/135)	(97.3%, 100.0%)	1899.0	86.1	4.5	22.8	1.2	0.0	0.0
	30	250	100.0% (135/135)	(97.3%, 100.0%)	1884.2	94.0	5.0	13.8	0.7	0.0	0.0
	75	100	100.0% (135/135)	(97.3%, 100.0%)	1867.2	87.7	4.7	0.0	0.0	19.2	1.0

A second study was conducted internally to characterize the reproducibility of test results (i.e., proportion positive or negative) at target levels below the analytical Limit of Detection (LOD) of the **BD ProbeTec GC Q^x Assay**. A panel of simulated specimens was tested that comprised GC and CT organisms seeded into Q^x swab diluent at two different levels (1:10, 1:100) each of which was below the analytical LOD for the respective organism. These levels were selected to fall within the dynamic range of the analytical LOD curve of the assay. Fifteen replicates of each panel member were tested every day for five days across three **BD Viper** Systems. The data are summarized in Table 15B.

Table 15B: Characterization of System Reproducibility at Target Levels below the Analytical Limit of Detection for the GC Q^x Assay for Swab and Urine Specimens

Specimen	Dilution of Analytical LOD	% Positive	95% CI (Positive)	MaxRFU Mean (Positive)	% Negative	95% CI (Negative)	MaxRFU Mean (Negative)
Endocervical/ Urethral	1:10	92.9 (209/225)	(88.7, 95.9)	1324.6	7.1 (16/225)	(4.1, 11.3)	41.4
Endocervical/ Urethral	1:100	30.7 (69/225)	(24.7, 37.1)	835.9	69.3 (156/225)	(62.9, 75.3)	7.2
Urine/Vaginal	1:10	90.7 (204/225)	(86.1, 94.1)	1165.9	9.3 (21/225)	(5.9, 13.9)	34.2
Urine/Vaginal	1:100	22.7 (51/225)	(17.4, 28.7)	872.7	77.3 (174/225)	(71.3, 82.6)	7.8

A reproducibility study of the **BD Viper** System using the **BD ProbeTec GC Q^x Assay** was also conducted for Liquid Based Cytology (LBC) specimens at three clinical sites on one **BD Viper** System per site. A panel of simulated specimens comprising CT and GC organisms seeded into LBC Specimen Dilution Tubes containing LBC medium was tested with the **BD ProbeTec GC Q^x Assay**. Uninoculated LBC Specimen Dilution Tubes containing LBC medium were used for the GC negative samples. Nine replicates of each panel member were tested every day for five days on each **BD Viper** System. The data are summarized in Table 15C. Two additional levels were included in the panels to characterize the reproducibility of test results (i.e., proportion positive or negative) at target levels below the analytical Limit of Detection (LOD) of the **BD ProbeTec GC Q^x Assay**. These additional specimens comprised CT and GC organisms seeded into LBC Specimen Dilution Tubes containing LBC medium at dilutions of 1:10 and 1:100 of the respective analytical LODs of each analyte. These levels were selected to fall within the dynamic range of the analytical LOD curves for the **BD ProbeTec CT Q^x** and **GC Q^x** assays. Nine replicates of each panel member were tested every day for five days across the three **BD Viper** Systems. The data are summarized in Table 15D.

Table 15C: Summary of Reproducibility Data for LBC Specimens on the BD Viper System for the GC Q^x Assay

					Within Run		Between Runs Within Site		Between Site	
CT EBs/mL	GC Cells/mL	% Correct	95% CI	Mean MaxRFU	SD	%CV	SD	%CV	SD	%CV
0	0	100.0% (135/135)	(97.3% – 100.0%)	1.21	4.00	330.38	0.00	0.00	0.00	0.00
30	0	100.0% (135/135)	(97.3% – 100.0%)	0.98	7.47	761.30	0.00	0.00	0.17	17.04
0	100	100.0% (135/135)	(97.3% – 100.0%)	1982.77	83.92	4.23	0.00	0.00	0.00	0.00
30	250	100.0% (135/135)	(97.3% – 100.0%)	1983.66	87.76	4.42	0.00	0.00	24.80	1.25
75	100	100.0% (135/135)	(97.3% – 100.0%)	1920.14	81.94	4.27	59.45	3.10	0.00	0.00

Table 15D: Characterization of System Reproducibility at Target Levels Below the Analytical Limit of Detection for the GC Q^x Assay for LBC Specimens

Dilution of Analytical LOD	% Positive	95% CI (Positive)	MaxRFU Mean (Positive)	% Negative	95% CI (Negative)	MaxRFU Mean (Negative)
1:10	74.1 (100/135)	(65.8 - 81.2)	1159.2	25.9 (35/135)	(18.8 - 34.2)	21.2
1:100	8.9 (12/135)	(4.7 - 15.0)	1136.5	91.1 (123/135)	(85.0 - 95.3)	6.6

System Cross Contamination and Carryover

An internal study was conducted to evaluate the risk of producing a false positive result in either the same run on the **BD Viper** System in extracted mode (within run cross-contamination) or in a subsequent run (between run carryover). Testing was conducted using negative and positive samples on three **BD Viper** Systems. Negative samples consisted of CT/GC Q^x Swab Diluent/LBC Specimen Dilution Tube with PreservCyt Solution. Positive samples consisted of a representative analyte (10⁵ CT EB/mL) spiked into CT/GC Q^x Swab Diluent/LBC Specimen Dilution Tube with PreservCyt Solution. The overall rate of cross-contamination (i.e., with alternating columns of positive and negative samples and a prevalence of 50%) was 0.41% (9/2208) for the CT/GC Q^x Swab Diluent and 0.45% (5/1104) for the LBC Specimen Dilution Tube with PreservCyt Solution. The overall rate of carryover contamination (i.e., carryover between successive runs when the prevalence was 50% in the previous run) was 0.36% (8/2208) for the CT/GC Q^x Swab diluent and 0.54% (6/1104) for the LBC Specimen Dilution Tube with PreservCyt Solution. Cross-contamination and carryover rates across the three **BD Viper** Systems are summarized in Tables 16A and 16B.

Table 16A: Cross Contamination and Carryover Contamination (Swab/Urine)

Assay Dispense Mode Selected	BD Viper System	Cross-Contamination			Carryover Contamination		
		n	Positive Results	Percent Positive	n	Positive Results	Percent Positive
Dual Assay	1	736	5	0.68	736	1	0.14
	2	736	0	0.00	736	3	0.41
	3	736	4	0.54	736	4	0.54
	Overall	2208	9	0.41	2208	8	0.36
Single Assay	1	190	0	0.00	186	0	0.00
	2	188	1	0.53	186	1	0.54
	3	188	0	0.00	186	0	0.00
	Overall	566	1	0.18	558	1	0.18

Table 16B: Cross Contamination and Carryover Contamination (LBC Medium)

Medium Type	BD Viper System	Cross-Contamination			Carryover Contamination		
		n	Positive Results	Percent Positive	n	Positive Results	Percent Positive
PreservCyt	1	368	1	0.27	368	1	0.27
	2	368	3	0.82	368	0	0.00
	3	368	1	0.27	368	5	0.45
	Overall	1104	5	0.45	1104	6	0.54

INTERPRETATION OF TABLES

Symbols and Abbreviations

Symbols

(+)	positive
(-)	negative
#	number
%	percentage

Abbreviations

A	Asymptomatic
CI	Confidence Interval
CT	<i>Chlamydia trachomatis</i>
CV	Coefficient of Variation
E	Equivocal
EC	Extraction Control
ET	Extraction Transfer Error
FN	False Negative
FNU	Female Neat Urine
FP	False Positive
FS	Female endocervical swab
FUPT	Female urine in Q ^x UPT
FV	Female vaginal swab
GC	<i>Neisseria gonorrhoeae</i>
HIV	Human Immunodeficiency Virus
I	Indeterminate
IFU	Inclusion Forming Units
LBC	Liquid Based Cytology
LE	Liquid level error
LOD	Limit of Detection
MaxRFU	Maximum relative fluorescent units
MNU	Male Neat Urine
MS	Male urethral swab
MUPT	Male urine in Q ^x UPT
n	number
NA	non-applicable
NAAT	Nucleic Acid Amplification Test
NPA	Negative Percent Agreement
NPV	Negative Predictive Value
OB/GYN	Obstetrics/Gynecology
PA	Percent Agreement
PBS	Phosphate Buffered Saline
PIS	Patient Infected Status
PPA	Positive Percent Agreement
PPV	Positive Predictive Value
QC	Quality Control
S	Symptomatic
SD	Standard Deviation
SDA	Strand Displacement Amplification
STD	Sexually Transmitted Disease
TN	True Negative
TP	True Positive
UPT	Urine Preservative Transport

AVAILABILITY

The following **BD ProbeTec CT/GC Q^x** and **BD Viper** products are also available:

Cat. No.	Description
440724	BD Viper™ Pipette Tips, 960
441392	BD Viper™ Trash Box
441391	BD Viper™ Trash Bags
440818	BD Viper™ Trash Boxes and Bags
440974	BD Viper™ Tube Lockdown Cover
440975	BD Viper™ Lysing Heater (115V)
440976	BD Viper™ Lysing Heater (230V)
440977	BD Viper™ Lysing Rack
440984	Amplification Plate Sealers (Black)
441072	BD Viper™ Liquid Waste Bottle
441074	BD Viper™ Plate Seal Tool
441091	BD Viper™ System
441122	Vaginal Specimen Transport for the BD ProbeTec™ CT/GC Q ^x Amplified DNA Assays, 100 units
441124	BD ProbeTec™ GC Q ^x Amplified DNA Assay Reagent Pack, 1152 tests
441126	BD ProbeTec™ CT Q ^x Amplified DNA Assay Reagent Pack, 1152 tests
441125	Control Set for the BD ProbeTec™ CT/GC Q ^x Amplified DNA Assays, 24 positive and 24 negative
441128	BD Viper™ Extraction Reagent and Lysis Trough, 12 Extraction Reagent Troughs and 12 Lysis Troughs
441129	BD FOX™ Extraction Tubes, 384 tests.
441354	BD Viper™ Neutralization Pouch, 12 pouches
441357	BD ProbeTec™ Q ^x Collection Kit for Endocervical or Lesion Specimens, 100 units
441358	Male Urethral Specimen Collection Kit for the BD ProbeTec™ CT/GC Q ^x Amplified DNA Assays, 100 units
441359	Caps for use on the BD Viper™ (Extracted Mode), 4 x 100
441360	Specimen Tubes and Caps for use on the BD Viper™ (Extracted Mode), 4 x 100
441361	Swab Diluent for the BD ProbeTec™ Q ^x Amplified DNA Assays, 2 mL x 48
441362	BD™ Urine Preservative Transport for the Q ^x Amplified DNA Assays, 100 units
441444	Liquid Based Cytology Specimen (LBC) Dilution Tubes for the BD ProbeTec Q ^x Amplified DNA Assays
441443	Liquid Based Cytology Specimen (LBC) Dilution Tube Caps for the BD ProbeTec Q ^x Amplified DNA Assays

The following strains are available from:

American Type Culture Collection (ATCC)

10801 University Boulevard

Manassas, VA 20110-2209, USA.

ATCC # 19424 *Neisseria gonorrhoeae*

ATCC # VR-879 *Chlamydia trachomatis* (serotype H)

ATCC # VR-902B *Chlamydia trachomatis* LGV II

Bio-Rad AmpliTrol CT/GC is available from:

Bio-Rad Laboratories (Blackhawk Biosystems)

12945 Alcosta Blvd. 2nd Floor

San Ramon, CA 94583

1-800-866-0305

AmpliTrol CT/GC # 00126

REFERENCES

1. World Health Organization. 2001. Global prevalence and incidence of selected curable sexually transmitted infections: overview and estimates. WHO.
2. Centers for Disease Control and Prevention. STD surveillance 2006: national profile. <http://www.cdc.gov/std/stats/gonorrhoea.htm>.
3. US Preventive Services Task Force. 2005. Screening for gonorrhea: recommendation statement. *Ann Fam Med* 3: 262-267.
4. Advisory Committee for HIV and STD Prevention. 1998. HIV prevention through early detection and treatment of other sexually transmitted diseases – United States. *MMWR* 47 (RR-12): 1-24.
5. Centers for Disease Control and Prevention. 2007. Update to CDC's Sexually Transmitted Diseases Treatment Guidelines, 2006: fluoroquinolones no longer recommended for treatment of gonococcal infections. *MMWR* 56: 332-336.
6. Centers for Disease Control and Prevention. 2002. Screening tests to detect *Chlamydia trachomatis* and *Neisseria gonorrhoeae* infections – 2002. *MMWR* 51 (RR-15): 1-40.
7. Centers for Disease Control and Prevention. 2006. Sexually transmitted diseases treatment guidelines, 2006. *MMWR* 55 (RR-11): 1-94.
8. Little, MC, J Andrews, R Moore, et al. 1999. Strand displacement amplification and homogeneous real-time detection incorporated in a second-generation DNA probe system, *BDProbeTec ET*. *Clin Chem* 45: 777-784.
9. Hellyer, TJ and J.G. Nadeau. 2004. Strand displacement amplification: a versatile tool for molecular diagnostics. *Expert Rev Mol Diagn* 4: 251-261.

10. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed., CLSI. Wayne, PA.
11. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hosp. Epidemiol. 17:53-80.
12. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
13. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
14. Clinical and Laboratory Standards Institute. 2006. Approved Guideline C24-A3. Statistical quality control for quantitative measurement procedures: principles and definitions, 3rd ed. CLSI. Wayne, PA.
15. Brunton, WAT, H Young, and DRK Fraser. 1980. Isolation of *Neisseria lactamica* from the female genital tract. Br. J. Vener. Dis. 56 : 325-326.
16. Knapp, JS, and EW Hook. 1988. Prevalence and persistence of *Neisseria cinerea* and other *Neisseria* spp. in adults. J. Clin. Microbiol. 26 : 896-900.
17. Knapp, JS, PA Totten, MH Mulks, and BH Minshew. 1984. Characterization of *Neisseria cinerea*, a nonpathogenic species isolated on Martin-Lewis medium selective for pathogenic *Neisseria* spp. J. Clin. Microbiol. 19 : 63-67.
18. Wilkinson, AE. 1952. Occurrence of *Neisseria* other than the gonococcus in the genital tract. Br. J. Vener. Dis. 28 : 24-27.

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Nº de cat.	Descripción
440724	BD Viper Pipette Tips, 960
441392	BD Viper Trash Box
441391	BD Viper Trash Bags
440818	BD Viper Trash Boxes and Bags
440974	BD Viper Tube Lockdown Cover
440975	BD Viper Lysing Heater (115V)
440976	BD Viper Lysing Heater (230V)
440977	BD Viper Lysing Rack
440984	Amplification Plate Sealers (negro)
441072	BD Viper Liquid Waste Bottle
441074	BD Viper Plate Seal Tool
441091	BD Viper System
441122	Vaginal Specimen Transport for the BD ProbeTec CT/GC Q ^x Amplified DNA Assays, 100 unidades
441124	BD ProbeTec GC Q ^x Amplified DNA Assay Reagent Pack, 1.152 análisis
441126	BD ProbeTec CT Q ^x Amplified DNA Assay Reagent Pack, 1.152 análisis
441125	Control Set for the BD ProbeTec CT/GC Q ^x Amplified DNA Assays, 24 positivos y 24 negativos
441128	BD Viper Extraction Reagent and Lysis Trough, 12 cubetas de reactivo de extracción y 12 cubetas de lisis
441129	BD FOX Extraction Tubes, 384 análisis
441354	BD Viper Neutralization Pouch, 12 bolsas
441357	BD ProbeTec Q ^x Collection Kit for Endocervical or Lesion Specimens (equipo de recogida de muestras endocervicales o de lesiones), 100 unidades
441358	Male Urethral Specimen Collection Kit for the BD ProbeTec CT/GC Q ^x Amplified DNA Assays, 100 unidades
441359	Caps for use on the BD Viper (modo de extracción), 4 x 100
441360	Specimen Tubes and Caps for use on the BD Viper (modo de extracción), 4 x 100
441361	Swab Diluent for the BD ProbeTec Q ^x Amplified DNA Assays, 2 mL x 48
441362	BD Urine Preservative Transport for the Q ^x Amplified DNA Assays, 100 unidades
441444	Liquid Based Cytology Specimen (LBC) Dilution Tubes for the BD ProbeTec Q ^x Amplified DNA Assays
441443	Liquid Based Cytology Specimen (LBC) Dilution Tube Caps for the BD ProbeTec Q ^x Amplified DNA Assays

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA
(800) 638-8663
www.bd.com/ds

Benex Limited
Rineanna House
Shannon Free Zone
Shannon, County Clare, Ireland

AmpliTrol is a trademark of Bio-Rad Laboratories, Inc.
ATCC is a trademark of the American Type Culture Collection.
DNA AWAY is a trademark of Molecular BioProducts.
PreservCyt and ThinPrep are trademarks of Hologic, Inc.
BD, BD Logo, NAP Guard, BD ProbeTec, BD FOX, BD SurePath, BD PrepStain and BD Viper are trademarks of Becton, Dickinson and Company. ©2010 BD.