
### Revisions

Rev from	Rev to	ECO #
0805	0308	4749-08

**Notes:**

1. BD Cat. Number 236500
2. Blank (Sheet) Size : Length: 10"      Width: 22.5"  
 Number of Pages: 6      Number of Sheets: 1  
 Page Size: Length 10"      Width 7.5"      Final Folded Size: 10" x 7.5"
3. Style (see illustrations below): # 3

4. See Specification Control Number N/A for Material Information
5. Ink Colors: Printed two sides  Yes     No  
 No. of Colors: 1      PMS# Standard Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	Date	<p style="font-size: small;">COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION</p> <p style="font-weight: bold; font-size: large;">BD</p> <p>Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA</p>				
Proofer	Date					
Checked By	Date					
Part Number: 8082048		Category and Description Package Insert, BBL Vacutainer Anaerobic Specimen Collector	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%;">Sheet: 1 of 7</td> <td rowspan="2" style="text-align: center; font-size: 2em; font-weight: bold;">A</td> </tr> <tr> <td>Scale: N/A</td> </tr> </table>	Sheet: 1 of 7	A	Scale: N/A
Sheet: 1 of 7	A					
Scale: N/A						

# BD BBL™ Collection and Transport Systems

## BBL™ Vacutainer™ Anaerobic Specimen Collector (Sterile)

English: pages 1 – 2  
Français : pages 2 – 3  
Español: páginas 4 – 5

VDP40042-03  
8082048  
2008/03

### INTENDED USE

The **Vacutainer™** Anaerobic Specimen Collector with swab provides a method to protect anaerobic bacteria from exposure to toxic amounts of oxygen during the time from specimen collection to inoculation of suitable media. The system may also be used for liquid specimens collected with a needle and syringe.

Although primarily designed for the collection of specimens suspected of containing anaerobic organisms, aerobic organisms will survive in this system.

The **Vacutainer** Anaerobic Specimen Collector is particularly suited for instances where volumes are small or dilution of specimens by a transport medium is to be avoided.

### SUMMARY AND EXPLANATION

Key factors in the recovery of anaerobic bacteria are protection from the deleterious effects of oxygen exposure and drying during transport to the laboratory. The **Vacutainer** Anaerobic Specimen Collector provides a moist, anaerobic atmosphere.<sup>1,2</sup> Methods used prior to this system consisted of either placing a swab or specimen in a carbon dioxide-filled container or transport medium, drawing a specimen into a syringe (when sufficient volume is available) and immediately transporting the specimen, or injecting the specimen into a carbon dioxide-flushed vial. Although prolonged storage of any specimen is not recommended, 72 h survival of anaerobes and perhaps longer could be expected in the Anaerobic Specimen Collector.

### PRINCIPLES OF THE PROCEDURE

Each sterile package contains a rayon-tipped swab for the collection of specimens. The outer tube (F) of the **Vacutainer** Anaerobic Specimen Collector contains hydrogen in an oxygen-free atmosphere of gases. After the plunger is depressed, a minimal amount of oxygen is introduced.

The hydrogen combines in the presence of the palladium catalyst (D) with the oxygen to form water.

An indicator disc (E) saturated with a buffered resazurin solution is deep pink in the presence of toxic quantities of atmospheric oxygen. The indicator is colorless (disc white) in the absence of these amounts (resazurin is colorless at Eh-42 mV).

The oxygen present in the inner tube (C) may cause the indicator disc to turn pink in about 10 min. This should revert to colorless within approximately 90 min if the specimen collector is properly used (reduction in the intensity of the color will be noted within 45 min). As with most indicator systems, the indicator response is dependent on the rate of diffusion and gaseous exchange between the solid support and the surrounding atmosphere. The indicator response (color to colorless) to the reduction process is slow compared to the relative speedy lowering of the redox potential (Eh), as evidenced by the formation of moisture on the interior wall of the inner tube.

Moisture from the reaction, moisture present from indicator (E) and resulting anaerobic conditions provide an atmosphere conducive to survival of anaerobic and aerobic bacteria.

### REAGENTS

**Warnings and Precautions:** For *in vitro* Diagnostic Use.

Contents of unopened or undamaged package guaranteed sterile.

**Warning:** Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>3-6</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Prior to discarding, specimen containers, swabs and other contaminated materials must be sterilized by autoclaving.

Whenever possible, specimens should be collected before antimicrobial agents have been administered.

It is important the plunger not be depressed in the immediate vicinity of the patient. Product is glass and breakage could occur.

Collection of an adequate amount of specimen is necessary to permit a complete examination.

Transportation time should be as short as possible and specimens refused if they have been in transport for an excessive period of time (according to institutional guidelines).

All specimens should be handled aseptically.

An uncultured specimen should never be discarded without obtaining a new specimen or consulting with the patient's physician or both.

Proper labelling with the patient's full name, ward number, physician, site of sample, time, date and laboratory requests should accompany each specimen to the laboratory.

**Storage Instructions:** Store at 4 – 25°C; product should be used before expiration date. Do not expose to sunlight.

**Product Deterioration:** DO NOT USE IF PACKAGE IS DAMAGED. DO NOT USE IF INDICATOR DISK IS PINK.

### PROCEDURE

**Material Provided:** Anaerobic Specimen Collector.

**Materials Required But Not Provided:** Sterile alcohol swabs, sterile needle and syringe, appropriate sterile instruments for obtaining tissue specimens, materials required for isolation, identification and for other microbiological procedures for specimen cultivation.

**Instructions:** Only the infected area should be sampled. Surrounding surface areas may be inhabited or contaminated by saprophytic anaerobes.

### Safety Warning Message

- Handle all biologic samples and blood collection "sharps" (lancets, needles and blood collection sets) in accordance with the policies and procedures of your facility. Blood collection "sharps" should be discarded in biohazard containers approved for their disposal.
- Obtain appropriate medical attention in the event of any exposure to biologic samples (for example, through a puncture injury) since the samples may transmit HBV/HCV (HEPATITIS), HIV (AIDS), or other infectious diseases.
- Utilize any built-in used needle protector, if the blood collection device provides one. BD does not recommend reshielding used needles, but the policies and procedures of your facility may differ and should always be followed.


### Specimen Collection

**CAUTION:** Care should be used in sampling the suspected lesion.

\* DO NOT REMOVE STOPPER DURING COLLECTION.

#### A. Collection by Swab:

1. Peel apart package (G) and remove specimen collector.
2. Remove plunger with sterile swab attached (B).
3. Obtain specimen using precaution noted above.\*
4. Replace swab through holes and into inner tube (C).\*
5. Press down on disc portion of plastic plunger (B) with continuous gentle force until plastic disc rests against top of rubber stopper (A) forcing the inner tube (C) into outer tube (F). HOLD TUBE AT 45° ANGLE WHILE DEPRESSING PLUNGER.\*

NOTE: The hole will be plugged by the plastic plunger and the small inner tube will be detached from the stopper.

6. Hold tube at approximately 10 – 30° angle to the floor and rotate with a swirling motion. This is to facilitate mixing of air in inner tube with hydrogen atmosphere in outer tube.\*
7. Transport to laboratory. For collectors containing liquid specimens, maintain in an upright position.\*

The **Vacutainer** Anaerobic Specimen Collector, as is, is not suitable for transportation of specimens through the mail. The user is advised to refer to U.S. Public Health Act (42 USC 27 [A] Section 3, 68 [A] ) or to the U.S. Centers for Disease Control and Prevention, Atlanta, Georgia for requirements concerning the mailing of etiologic agents and specimens. This relates to an outer mailing container of metal, wrapping with absorbent material, and indicating the contents are biohazardous material.

#### B. Collection of Liquid or Purulent Specimens:

1. Collect specimen with sterile syringe and needle. Air trapped in syringe should be expelled by holding syringe and needle upright. Expel air at tip of syringe into alcohol saturated sponge. Set aside with sponge or needle shield on needle.
2. Peel apart package (G) and remove specimen collector.
3. Remove swab plunger unit (B) and expel material into inner tube (C).\*

4. Perform steps 4 – 7 as described above.\*

**User Quality Control:** Check transport device performance by inoculation with pure cultures of stable control organisms. The following quality control cultures are recommended for viability/recovery testing:

*Fusobacterium necrophorum* ATCC™ 27852  
*Peptostreptococcus anaerobius* ATCC 27337

1. Inoculate the inner tube of two **Vacutainer** Anaerobic Specimen Collectors using 1 mL of a 0.5 McFarland suspension. Organism test suspension should be prepared immediately before use by inoculating sterile saline solution using fresh 18 – 24 h plate cultures.
2. Store inoculated collectors upright at room temperature for 18 – 24 h.
3. Remove swab and inoculate appropriate media using a standard four quadrant streak method. Incubate plates at 35 – 37°C in an anaerobic atmosphere.
4. Remove plates after 18 – 24 h of incubation and check for recovery.
5. All plate cultures should demonstrate at least 20 colonies in the first quadrant.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

#### LIMITATIONS OF THE PROCEDURE

Recovery of microorganisms is influenced by numerous factors; e.g., specimen composition, organism type and concentration, transport temperature and time in transit. Temperature extremes during transport should be avoided. Excessive transport time may result in overgrowth and/or loss of viability.

The product is not intended for collection of specimens suspected of containing *N. meningitidis* or *N. gonorrhoeae*.

Although many facultative anaerobic and strict aerobic organisms do survive in this system, it is primarily designed for the collection and holding or transportation of anaerobic specimens.

#### EXPECTED RESULTS AND PERFORMANCE CHARACTERISTICS

The **Vacutainer** Anaerobic Specimen Collector will transport and maintain most anaerobic microorganisms for up to 72 h at room temperature. All organisms will exhibit some loss of viability over time and transport time should be minimized.

Wilkins and Jimenez-Ulate reported survival of twelve strains of anaerobic bacteria in logarithmic phase of growth in this system.<sup>7</sup> They classified the organisms into groups of "very oxygen sensitive," "moderately oxygen sensitive," and "oxygen tolerant." The "very oxygen sensitive" strains were nonviable after 2 h exposure to air but survived the total 48 h holding-period in the **Vacutainer** Anaerobic Specimen Collector.

Hill investigated the **Vacutainer** Anaerobic Specimen Collector with *in vitro* prepared and clinical material for the transportation of polymicrobial specimens containing anaerobes, aerobes and facultative anaerobes.<sup>8</sup> The **Vacutainer** Anaerobic Specimen Collector was more efficient for maintaining the status quo of the microbial population in polymicrobial mixtures than was storage in an aerobic environment or anaerobic glove box. The Anaerobic Specimen Collector was shown to be an excellent transport container for anaerobes.

#### AVAILABILITY

Cat. No.	Description
236500	<b>BBL™ Vacutainer™</b> Anaerobic Specimen Collector (Sterile), Ctn. of 100.

#### REFERENCES

1. Hillier, S.L., and B.J. Moncla. 1995. *Peptostreptococcus*, *Propionibacterium*, *Eubacterium*, and other nonsporeforming anaerobic gram-positive bacteria, p.587-602. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
2. Holden, J. 1992. Collection and transport of clinical specimens for anaerobic culture, p. 2.2.3 (suppl. 1). In H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
3. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
4. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
5. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
6. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
7. Wilkins, T.D., and F. Jimenez-Ulate. 1975. Anaerobic specimen transport device. J. Clin. Microbiol. 2:441-447.
8. Hill, G.B. 1978. Effects of storage in an anaerobic transport system on bacteria in known polymicrobial mixtures and in clinical specimens. J. Clin. Microbiol. 8:680-688.

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