Revisions

Rev from	Rev to	ECO #
1008	2010/07	5440-10

Notes:


1. BD Cat. Number 222321, 229101, 222281, 222291, 222301, 222311, 228801, 228811, 228911, 222521, 222531, 241110, 241111, 241112, 241113
2. Blank (Sheet) Size: Length: 9" Width: 40"
 Number of Pages: 20 Number of Sheets: 1
 Page Size: Length 9" Width 4" Final Folded Size: 2.25" x 4"
3. Style (see illustrations below): # 4


#1

#2


#3


#4

#5

#6

#7

4. See Specification Control Number N/A for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# 2755 Blue
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8085876		Category and Description Package Insert, Neisseria Meningitidis Antisera	Sheet: 1 of 21 Scale: N/A	A

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INTENDED USE

Difco™ Neisseria Meningitidis Antisera are recommended for use in slide agglutination tests for serotyping *Neisseria meningitidis*.

SUMMARY AND EXPLANATION

Neisseria meningitidis, also known as meningococcus, is found in the oropharynx and nasopharynx of humans. Because the organism survives poorly in the environment, humans are the primary reservoir. In asymptomatic persons the carrier state lasts for variable time periods, usually several weeks. The microorganism is transmitted from person-to-person by direct contact with respiratory secretions or airborne droplets.^{1,2}

Typical human specimens for isolating the organism are: cerebrospinal fluid (CSF), blood, skin lesions (in cases where petechiae occur) and nasopharyngeal swabs. *N. meningitidis* occurs in the cervix and vagina of females and can cause serious pelvic disease. Other sites include the anal canal and in males, the urethra.

N. meningitidis is divided into serologic groups based on the presence of either capsular or outer membrane protein antigens. Among the currently recognized groups are A, B, C, D, 29E, H, I, K, L, X, Y, Z, Z' and W135. Groups A, B, C, Y and W135 are most frequently implicated in systemic disease.³ Other serogroups are sporadically isolated from carriers and patients with disease.

N. meningitidis are gram-negative cocci usually occurring in pairs called diplococci. They are strict aerobes and produce the enzyme cytochrome oxidase. The growth of *N. meningitidis* is enhanced by a CO₂-enriched atmosphere.

PRINCIPLES OF THE PROCEDURE

Identification of *N. meningitidis* includes the isolation of the microorganism, biochemical identification and serologic confirmation.

Serologic confirmation involves the reaction in which the microorganism (antigen) reacts with its corresponding antibody. This *in vitro* reaction produces macroscopic clumping called agglutination. The desired homologous reaction is rapid, does not dissociate (high avidity) and binds tightly (high affinity).

Because a microorganism (antigen) may agglutinate with an antibody produced in response to some other species, heterologous reactions are possible. These are characterized as weak in strength or slow in formation. Such unexpected and perhaps unpredictable reactions may lead to some confusion in serologic identification. Therefore, a positive homologous agglutination reaction should support the morphological and biochemical identification of the microorganism.

Homologous reactions occur rapidly and are strong. Heterologous reactions form slowly and are weak.

REAGENTS

Difco Neisseria Meningitidis Antisera are stable, lyophilized, polyclonal rabbit antisera containing approximately 0.1% sodium azide as a preservative.

The Difco Neisseria Meningitidis Antisera detect the following antigenic groups:

Antiserum	Antigenic Group Detected	Antiserum	Antigenic Group Detected
Poly	A, B, C, D	D	D
Poly 2	X, Y, Z*	X	X*
W135	W135*	Y	Y*
A	A*	Z	Z*
B	B*	Z'	Z, Z'*
C	C*		

*Not absorbed for Group D, which is rarely isolated.

When properly rehydrated and used as described (see Procedure) each vial of Difco Neisseria Meningitidis Antiserum is sufficient reagent for 20 slide tests.

Warnings and Precautions

For *in vitro* Diagnostic Use.

The Packaging of This Product Contains Dry Natural Rubber.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, specimens, containers, slides, tubes and other contaminated material must be sterilized by autoclaving. Directions for use should be followed carefully.

WARNING: This product contains sodium azide. Sodium azide is toxic by inhalation, by skin contact, and if swallowed. Contact with acid liberates very toxic gas. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Storage

Store lyophilized and rehydrated **Difco** Neisseria Meningitidis Antisera at 2 – 8°C.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

Expiration date applies to product in its intact container when stored as directed. Do not use if the product is caked, discolored or shows other signs of deterioration.

SPECIMEN COLLECTION AND PREPARATION

N. meningitidis can be recovered on blood agar or chocolate agar. For specific recommendations consult appropriate references.^{4,5} Determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as *N. meningitidis*. After these criteria are met, serologic identification can be performed.

PROCEDURE

Materials Provided: **Difco** Neisseria Meningitidis Antiserum Poly; **Difco** Neisseria Meningitidis Antiserum Poly 2; **Difco** Neisseria Meningitidis Antiserum A; **Difco** Neisseria Meningitidis Antiserum B; **Difco** Neisseria Meningitidis Antiserum C; **Difco** Neisseria Meningitidis Antiserum D; **Difco** Neisseria Meningitidis Antiserum X; **Difco** Neisseria Meningitidis Antiserum Y; **Difco** Neisseria Meningitidis Antiserum Z; **Difco** Neisseria Meningitidis Antiserum Z'; **Difco** Neisseria Meningitidis Antiserum W135.

Materials Required But Not Provided: Agglutination slides; Applicator sticks; Sterile purified water; Sterile 0.85% saline; Bacteriological loop.

Reagent Preparation

To rehydrate **Difco** Neisseria Meningitidis Antisera, add 1 mL sterile purified water and rotate gently to completely dissolve the contents.

Equilibrate all materials to room temperature prior to performing the tests.

Ensure that all glassware and pipettes are clean and free of residues such as detergent.

User Quality Control

At the time of use, test both positive and negative control cultures to check performance of the antisera, techniques and methodology.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

Test isolate for autoagglutination

- From the test culture on chocolate agar, transfer a loopful of growth to a drop of sterile 0.85% saline on a clean slide and emulsify the organism.
- Rotate the slide for one min, then observe for agglutination.
- If agglutination (autoagglutination) occurs, the culture is rough and cannot be tested. Subculture to chocolate agar, incubate and test the organism again as described in steps 1 and 2.
- If no agglutination occurs, proceed with testing the organism.

Choice of antisera

Test the organism first with **Difco** Neisseria Meningitidis Antisera Poly, Poly 2, Group Z' and Group W135. Depending upon the reaction, continue testing according to either a, b or c listed below.

- If agglutination occurs with **Difco** Neisseria Meningitidis Antisera Poly or Poly 2, test the organism with the monospecific antisera indicated:

If positive reaction with this antiserum:	Retest with these antisera:
Poly Poly 2	A, B, C and D X, Y and Z

- If agglutination occurs with **Difco** Neisseria Meningitidis Antiserum Group Z', test the organism with monospecific **Difco** Neisseria Meningitidis Antiserum Group Z. Cultures of *N. meningitidis* group Z' will not agglutinate **Difco** Neisseria Meningitidis Antiserum Group Z. However, cultures of *N. meningitidis* group Z will agglutinate **Difco** Neisseria Meningitidis Antiserum Group Z'. The expected agglutination patterns are as follows:

Organism	Neisseria Meningitidis Antiserum Group Z'	Neisseria Meningitidis Antiserum Group Z
<i>N. meningitidis</i> group Z'	3+	–
<i>N. meningitidis</i> group Z	+	3+

- If agglutination occurs with **Difco** Neisseria Meningitidis Antiserum W135, no further testing is necessary.

Procedural Steps

1. Dispense 1 drop of the **Difco** Neisseria Meningitidis Antiserum to be tested on an agglutination slide.
2. Transfer a loopful of growth of the test organism to the drop of antiserum and mix thoroughly.
3. Rotate the slide for one min and read for agglutination.
4. Repeat steps 1 – 3 for known positive and negative cultures.

RESULTS

1. Read and record results as follows:
 - 4+ 100% agglutination (background clear to slightly hazy)
 - 3+ 75% agglutination (background slightly cloudy)
 - 2+ 50% agglutination (background moderately cloudy)
 - 1+ 25% agglutination (background cloudy)
 - No agglutination
2. The positive control should show 3+ or greater agglutination.
3. The negative control should show no agglutination.
4. For test isolates, a 3+ or greater agglutination within one min is a positive result.

LIMITATIONS OF THE PROCEDURE

1. Correct interpretation of serologic reactions depends upon culture purity and morphological characteristics and biochemical reactions that are consistent with identification of the microorganism as *N. meningitidis*.
2. Serologic methods alone cannot identify the isolate as *N. meningitidis*. Organisms unrelated to *Neisseria*, yet capable of causing meningitis and other species of *Neisseria* can cross-react with meningococcal antisera.
3. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False positive reactions may occur.
4. Rough culture isolates do occur and will agglutinate spontaneously, causing agglutination of negative controls (autoagglutination). Smooth colonies must be selected and tested in serologic procedures.
5. *N. meningitidis* group A and group C may cross-react due to the presence of common capsular polysaccharides.
6. Group Z meningococci may agglutinate group Z' antiserum. Group Z' meningococci will not agglutinate group Z antiserum.
7. **Difco** Neisseria Meningitidis Antisera have been tested using undiluted cultures taken from agar media. These antisera have not been tested using antigen suspensions in saline or other diluents. If the user applies variations in the recommended steps (see **Procedure**) each lot of antiserum must be tested with known control cultures to verify expected reactions under the modified procedure.
8. Rehydrated **Difco** Neisseria Meningitidis Antisera that are cloudy or have a precipitate during the period of use should be discarded.
9. Confirmatory testing of uncommon serogroups X, W135, Y, Z and 29E is suggested.

PERFORMANCE CHARACTERISTICS⁶

The performance of **Difco** Neisseria Meningitidis Antisera, manufactured by Difco Laboratories, was compared to two other commercially available tests for serological grouping of meningococci in a study by van der Ende, Schuurman, Hopman, Fijen and Dankert.⁶ One commercial product used rabbit and mouse antibodies in slide agglutination tests. The other commercial product used rabbit and mouse antibodies in latex agglutination tests. Seventy two (72) reference strains were used in this study.⁶

The table below lists the number of correctly identified isolates obtained for each test system.

Serogroup	No. of Samples	Slide Agglutination Test using Difco™ Neisseria Meningitidis Antisera	Commercial Slide Agglutination Test	Commercial Latex Test
A	10	9	10	9
B	10	9	10	10
C	10	9	9	10
X	6	5	5	*
Y	9	7	9	9†
Z	4	1	4	*
W135	6	5	5	6‡
29E (Z')	2	2	*	*
Other serogroups & nonencapsulated strains which should be negative by the test.	15	14/15	15/17‡	23/27‡
Total	72	61/72	67/72	67/72

* The other 2 commercial kits do not contain antisera to detect *N. meningitidis* groups X, Z and/or 29E.

† Latex test cannot distinguish between *N. meningitidis* serogroups Y and W135.

‡ Serogroup 29E and serogroups X, Z and 29E are also included for the commercial slide and latex tests, respectively.

In this study, the sensitivity of the **Difco** Neisseria Meningitidis Antisera was 92% and the specificity was 67%.⁶ Using only antisera against the common serogroups A, B, C for serogrouping, the sensitivity was 93% and the specificity was 98%.⁶

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Difco is a trademark of Difco Laboratories, Inc., a subsidiary of Becton, Dickinson and Company. BD and BD Logo are trademarks of Becton, Dickinson and Company. © 2010 BD.
Benex Limited