
Revisions

SO 0191-5

Rev from	Rev to	ECO #
2010/11	2011/11	5950-11

Notes:

1. BD Cat. Number: 223591, 223781
2. Blank (Sheet) Size: Length: 9" Width: 36"
 Number of Pages: 18 Number of Sheets: 1
 Page Size: Length: 9" Width: 4" Final Folded Size: 2.25" x 4"
3. Style (see illustrations below): # 4

4. See Specification Control Number N/A for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS # 2755 BLue
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8085878		Category and Description Package Insert, FA Bordetella Pertussis and Parapertussis	Sheet: 1 of 19 <hr/> Scale: N/A	A

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INTENDED USE

Difco™ FA Bordetella Pertussis and Difco FA Bordetella Parapertussis are recommended for use in the direct fluorescent antibody technique for the identification of *Bordetella pertussis* and *Bordetella parapertussis*.

SUMMARY AND EXPLANATION

All members of the genus *Bordetella* are respiratory pathogens of warm-blooded animals. *B. pertussis* and *B. parapertussis* are two species of *Bordetella* that are uniquely human pathogens. These organisms adhere to, multiply among and remain localized in the ciliated epithelial cells of the respiratory tract. *B. pertussis* is the major cause of whooping cough or pertussis. *Bordetella parapertussis* is associated with a milder, less frequently occurring form of the disease.¹ Person to person transmission occurs by the aerosol route.

Pertussis is a highly contagious disease that attacks more than 90% in unimmunized populations.² Toxin production remains the major distinction of *B. pertussis*. Classic pertussis caused by *B. pertussis* occurs in three stages. The first stage or catarrhal stage is characterized by nonspecific symptoms similar to a cold or viral infection. The disease is highly communicable during this stage, and lasts 1 – 2 weeks. In the second stage or paroxysmal stage, the cough increases in intensity and frequency. This stage is marked by sudden attacks of severe, repetitive coughing, often cumulating with the characteristic whoop.³ The whooping sound is caused by the rapid inspiration of air after the clearance of mucus-blocked airways. This stage may last 1 – 4 weeks. The beginning of the convalescent stage is marked by a reduction in frequency and severity of coughing spells. Complete recovery may require weeks or months.

Despite the availability of an effective whole-cell vaccine, pertussis remains a disease of worldwide distribution because many developing nations do not have the resources for vaccinating their populations.⁴ Even in developed nations such as Great Britain and Sweden, major outbreaks have occurred. Pertussis is endemic in the United States, with most disease occurring as isolated cases. A shift in the age group affected by the disease has occurred. In the past, children in the 1 – 5 year age group were more prone to pertussis. Children less than one year of age² have become more susceptible to the disease, because of a decrease in passively transferred maternal antibodies, since adults do not receive booster vaccinations.

Bordetella spp. are tiny gram-negative coccobacilli occurring singly or in pairs and may exhibit a bipolar appearance. They are strict aerobes and some members of the genus are motile. *B. pertussis* and *B. parapertussis* are nonmotile, and produce no acid from carbohydrates. *B. pertussis* will not grow on common blood agar bases or chocolate agar, whereas *B. parapertussis* will grow on blood agar and sometimes chocolate agar. Media for primary isolation consist of potato infusion-based media such as Bordet-Gengou (BG) medium or charcoal-based

media such as Regan-Lowe (RL) medium supplemented with glycerol, peptones and horse or sheep blood.⁵ An antibiotic agent is often added to reduce the growth of normal flora. *B. pertussis* may be recovered from secretions collected from the posterior nasopharynx, bronchoalveolar lavage and transbronchial specimens.

The direct fluorescent antibody test (DFA) has been used for many years with various degrees of success for the rapid, direct detection of *B. pertussis* and *B. parapertussis* in nasopharyngeal specimens.⁵⁻⁷ Eldering, Eveland and Kendrick^{8,9} and Holwerda and Eldering¹⁰ showed the usefulness of the FA procedure, although complete correlation between the agglutination method and FA technique was not obtained. Nonetheless, the FA procedure could detect both smooth and rough cultures of *B. pertussis* and *B. parapertussis* and could also be applied to direct specimens. Further data obtained showed little or no cross reactions between conjugates prepared from *B. pertussis* and *B. parapertussis* cultures. The disadvantage of this procedure is that technical skill and experience of the technicians is required to perform and read the test. The DFA should always be used with and not as a replacement for culture.¹¹ It has been suggested that laboratories proficient in DFA and culture for pertussis should obtain a DFA sensitivity of 60% or better and a specificity of at least 90% over time compared with culture.¹¹

B. pertussis and *B. parapertussis* are slow growing organisms, developing in 3 – 4 days. By employing the fluorescent antibody technique, the time required to detect these organisms can be significantly reduced. The FA procedure may be applied to direct nasopharyngeal smears or may be used to identify young cultures of *B. pertussis* or *B. parapertussis*.

PRINCIPLES OF THE PROCEDURE

The FA technique involves the preparation of a smear from the clinical specimen or culture isolate on a glass slide. Nasopharyngeal swabs are obtained from the patient and inoculated into 0.5 mL of Casamino Acids solution.¹¹ Smears are made from this solution and stained with a specific antibody labeled with a fluorescent marker (fluorescein isothiocyanate or FITC) directed against *B. pertussis* or *B. parapertussis*. After incubation with the antibody preparation, smears are washed in a phosphate-buffered saline, air dried, cover slipped with fluorescent-antibody mounting fluid, and examined under a fluorescent microscope.

REAGENTS

Difco FA Bordetella Pertussis and Difco FA Bordetella Parapertussis are lyophilized, polyclonal, fluorescein-conjugated, chicken antisera. They have been prepared with modifications according to the method of Eldering, Eveland and Kendrick^{8,9} and Holwerda and Eldering.¹⁰ Approximately 0.1% sodium azide is added as a preservative.

Difco FA Buffer, Dried is a phosphate buffer-NaCl mixture, which, upon rehydration, yields a 0.85% saline solution buffered to pH 7.2.

Difco FA Mounting Fluid pH 7.2 is standardized reagent grade glycerin adjusted to pH 7.2.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

The Packaging of This Product Contains Dry Natural Rubber.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures.^{12,13} After use, specimens, containers, slides, tubes and other contaminated material must be sterilized by autoclaving. Directions for use should be followed carefully.

WARNING: This product contains sodium azide. Sodium azide is toxic by inhalation, by skin contact, and if swallowed. Contact with acid liberates very toxic gas. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Storage: Store lyophilized **Difco** FA Bordetella Pertussis and **Difco** FA Bordetella Parapertussis at 2 – 8°C. Aliquots of the titered conjugate may be put into small vials, frozen in the undiluted state and stored below -20°C for optimal stability. The conjugate should not be exposed to repeated freezing and thawing.

Store dehydrated FA Buffer, Dried below 30°C. Store rehydrated FA Buffer, Dried at 2 – 8°C.

Store FA Mounting Fluid pH 7.2 at 15 – 30°C.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

Product Deterioration: Expiration date applies to product in its intact container when stored as directed. Do not use if the product is caked, discolored or shows other signs of deterioration.

SPECIMEN COLLECTION AND PREPARATION

Direct Nasopharyngeal Smears

1. Obtain a nasopharyngeal swab and emulsify in 0.5 mL of sterile 1% **Bacto**™ Casamino Acids.¹¹
2. Hold specimen in Casamino Acids solution for no more than 2 h.
3. Smear emulsified specimen on a clean microscope slide.
4. Allow smear to air dry and fix by gentle heating or by 1 min immersion in 95% ethanol.

Culture Isolates

1. From clinical specimens, isolation of *Bordetella* requires the use of certain media such as Bordet-Gengou Agar or Regan-Lowe Agar. Colonies of *B. pertussis* on Bordet-Gengou Agar are very small, white, opaque, convex and entire. Colonies on Regan-Lowe Agar are round, domed, mercury-silver colored and shiny. For specific recommendations, consult appropriate references.^{3,11} Determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as *Bordetella* sp. After these criteria are met, serological identification can be performed.
2. Pick appropriate colonies and emulsify in approximately 2 mL sterile purified water. Adjust cell density to approximate a McFarland Turbidity Standard No.1.
3. Smear emulsified specimen on a clean microscope slide.
4. Allow smear to air dry and fix by 1 min immersion in 95% ethanol. Remove slides and allow them to air dry.

PROCEDURE

Materials Provided: **Difco** FA Bordetella Pertussis, **Difco** FA Bordetella Parapertussis, **Difco** FA Buffer, Dried, **Difco** FA Mounting Fluid pH 7.2, 1% **Bacto** Casamino Acids

Materials Required But Not Provided: Fluorescent microscope assembly:

Lamps: HBO-50, HBO-100, HBO-200 or Xenon XBO-150; 6X 5A Tungsten. Excitation wavelength: 365 nm. Ocular: 10X. Objective: 10X, 40X (Fluorite). Filters: BG-12 or KP490,

K515 or K530. Condenser: Dark field D1.20-1.40.

Microscope slides, 95% ethanol, Staining jar, Cover slips, McFarland Turbidity Standard No. 1, Staining tray or moisture chamber.

Reagent Preparation: Equilibrate all materials to room temperature prior to performing the tests. Ensure that all glassware and pipettes are clean and free of residues such as detergents.

To rehydrate **Difco** FA Bordetella Pertussis and **Difco** FA Bordetella Parapertussis, add 5 mL sterile purified water and rotate gently to completely dissolve the contents.

The working dilution of the conjugate should be determined shortly after its rehydration. The titer of a conjugate varies with the technique used, the fluorescent microscope used, the filter used and the age of the bulb.

The conjugate should be titrated using a known culture of *B. pertussis* or *B. parapertussis* homologous to the conjugate. Dilutions of the conjugate are made in rehydrated **Difco** FA Buffer. (Prepare FA Buffer, Dried by dissolving one vial or 10 g in 1 L purified water.) The titer is determined as follows:

Dilution of conjugate	Fluorescence
1:5	4+
1:10	4+
1:20	4+
1:40	4+
1:80	2+

In this example, the last 4+ fluorescence is found in a 1:40 dilution of the conjugate. One less dilution is chosen for a margin of safety. Therefore, the working dilution, in this case is 1:20.

To rehydrate FA Buffer, Dried, add 10 g to 1 L of purified water and stir until completely dissolved.

Difco FA Mounting Fluid pH 7.2 is ready for use.

Control Slides: Prepare positive and negative control slides using appropriate homologous and heterologous antigens (i.e., cultures), following the procedure listed under "Specimen Collection and Preparation, Culture Isolates."

Test Procedure

1. Add several drops (one drop equals ~35 μ L) of the appropriate **Difco** FA Bordetella conjugate to the fixed smear.
2. Spread the conjugate over the surface of the smear.
3. Place the slides in a staining tray or moisture chamber.
4. Incubate at room temperature for 30 min.
5. Remove the excess conjugate and place the slide in a staining jar containing **Difco** FA Buffer for 10 min with 2 changes of the buffer followed by 1 rinse in purified water for 2 min.
6. Remove the slide and allow to drain and air dry, or blot with bibulous paper.
7. Add a drop (~35 μ L) of **Difco** FA Mounting Fluid pH 7.2 to the center of the stained area and mount with a cover slip.
8. Examine each smear using a fluorescent microscope with an excitation wavelength of 365 nm under a 40X or 100X objective. Record presence or absence and degree of fluorescence.

User Quality Control

At the time of use, test both positive and negative antigen controls to check performance of the antisera, techniques and methodology.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

1. Read and record results based on intensity of fluorescence as follows:
 - 4+ Maximum fluorescence; brilliant yellow-green peripheral staining
 - 3+ Bright yellow-green peripheral staining

- 2+ Definite, but dull, yellow-green peripheral staining
 - 1+ Barely visible peripheral staining
 - Complete absence of yellow-green peripheral fluorescence
2. The positive control should show a 4+ reaction with the homologous conjugate.
 3. The negative control should not exceed a 1+ reaction with the heterologous conjugate.
 4. For the test smears, a 2+ fluorescence should be considered a positive result.
 5. If the positive control is less than 3+, or the negative control exceeds 1+, the conjugate may have deteriorated or pH may have changed in the FA Buffer or FA Mounting Fluid. Repeat test with new reagents.

LIMITATIONS OF THE PROCEDURE

1. When testing culture isolates, growth from antibiotic-free media should be used to avoid autoagglutination.¹¹
2. Some experience is required to grade the intensity of the fluorescence and to ignore the occasional nonspecific staining of gram-negative diplococci, gram-positive cocci, and diphtheroid-like rods.¹¹
3. When testing culture isolates, the density of the positive control should be adjusted to give 4+ fluorescence with homologous conjugate. The density of culture isolates should be comparable to the positive control in order to standardize the fluorescence.
4. All glassware employed in the preparation, testing and storage of these reagents must be free of detergents or other harmful residues.
5. The fluorescent antibody technique can provide only presumptive identification of *B. pertussis* or *B. parapertussis*. A negative result should not be considered conclusive, as this type of reaction may occur when only a few organisms are present in the specimen. Final identification can only be made after consideration of biochemical, morphological and serological characteristics.
6. Cross reactions of variable staining intensities may occur with a small population of the cells in a *B. bronchiseptica* suspension at any dilution of the conjugate.
7. For both **Difco** FA Bordetella Pertussis and **Difco** FA Bordetella Parapertussis, the reaction should be brilliant and specific with smooth strains of homologous cultures. With some heterologous strains a 1+ reaction may occur. Such minimal reactions should not interfere with the interpretation of the test results.

PERFORMANCE CHARACTERISTICS¹⁴

The performance of **Difco** FA Bordetella Pertussis (DFA), manufactured by Difco Laboratories, a subsidiary of Becton, Dickinson and Company, was compared to culture and an in-house developed PCR assay in a study by Loeffelholz, Thompson, Long and Gilchrist.¹⁴ Three hundred nineteen (319) paired nasopharyngeal swab specimens were tested for the presence of *B. pertussis* in this study. The results of this study are summarized in the table below:

Method	Sensitivity*	Specificity
Culture	15.2%	100%
DFA	52.2%	98.2%
PCR	93.5%	97.1%

*Based on 46 positive specimens that were either (i) culture positive, (ii) PCR and DFA positive, or (iii) PCR or DFA positive, with clinical features indicating pertussis.

AVAILABILITY

Cat. No. Description

- 223591 **Difco**™ FA Bordetella Pertussis, 5 mL
- 223781 **Difco**™ FA Bordetella Parapertussis, 5 mL
- 223143 **Difco**™ FA Buffer, Dried, 6 x 10 mL
- 223142 **Difco**™ FA Buffer, Dried, 100 g
- 223050 **Bacto**™ Casamino Acids, 500 g

REFERENCES

1. Linneman, C.C. and E.B. Pery. 1977. *Bordetella parapertussis*: recent experience and a review of the literature. *Am. J. Dis. Child.* 131:560-563.
2. Bass, J.W. and S.R. Stephenson. 1987. The return of pertussis. *Pediatr. Infect. Dis. J.* 6:141-144.
3. Loeffelholz, M.J. 2003. *Bordetella*, p. 780-788. In Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller and R.H. Tenover, Manual of clinical microbiology, 8th ed., vol. 1. American Society for Microbiology, Washington, D.C.
4. Wright, P.F. 1991. Pertussis in developing countries: definition of the problem and prospects for control. *Rev. Infect. Dis.* 13: S228-S234.
5. Strelbel, P.M., S.L. Cochi, K.M. Farizo, B.J. Payne, S.D. Hanauer and A.L. Baughman. 1993. Pertussis in Missouri: evaluation of nasopharyngeal culture, direct fluorescent antibody testing, and clinical case definitions in the diagnosis of pertussis. *Clin. Infect. Dis.* 16:276-285.
6. Halperin, S.A., R. Bortolussi and A.J. Wort. 1989. Evaluation of culture, immunofluorescence and serology for the diagnosis of pertussis. *J. Clin. Microbiol.* 27:752-757.
7. Onorato, I.M. and S.G.F. Wassilak. 1987. Laboratory diagnosis of pertussis: the state of the art. *Pediatr. Infect. Dis. J.* 6:145-151.
8. Kendrick, P.L., G. Eldering and W.C. Eveland. 1961. Fluorescent antibody techniques. *Am. J. Diseases Children.* 101:149-154.
9. Eldering, G., W.C. Eveland and P.L. Kendrick. 1962. Fluorescent antibody staining and agglutination reactions in *Bordetella pertussis* cultures. *J. Bacteriol.* 83:745-749.
10. Holwerda, J. and G. Eldering. 1963. Culture and fluorescent antibody methods in diagnosis of whooping cough. *J. Bacteriol.* 86:449-451.
11. Pezzlo, M. 1994. Aerobic bacteriology, p. 1.0.1-1.20.47. In Isenberg, H.D. (ed.), *Clinical microbiology procedures handbook*, vol. 1. American Society for Microbiology, Washington, D.C.
12. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC) 5th ed. U.S. Government Printing Office, Washington, D.C.
13. Clinical and Laboratory Standards Institute. 2005. Approved guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
14. Loeffelholz, M.J., C.J. Thompson, K.S. Long and M.J.R. Gilchrist. 1999. Comparison of PCR, culture and direct fluorescent-antibody testing for detection of *Bordetella pertussis*. *J. Clin. Microbiol.* 37:2872.

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Difco and Bacto are trademarks of Difco Laboratories, a subsidiary of Becton, Dickinson and Company.
BD and BD Logo are trademarks of Becton, Dickinson and Company. ©2011 BD.
Benex Limited