## Revisions

Rev from	Rev to	ECO #
0604	2010/06	5369-10

**Notes:**

1. BD Cat. Number 231604,231605
2. Blank (Sheet) Size : Length: See VS    Width: See VS  
 Number of Pages: See VS    Number of Sheets: 1  
 Page Size: Length 9.875"    Width 4.5"                      Final Folded Size: See VS
3. Style (see illustrations below): # 4

4. VS Controlled By BD Caribe, LTD
5. Ink Colors: Printed two sides  Yes     No  
 No. of Colors: 1                      PMS# Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	<b>Becton, Dickinson and Company</b> 7 Loveton Circle Sparks, MD 21152 USA			
Proofer	Date			<b>Category and Description</b>  Package Insert, Sensi-Disc Elution Discs, Anaerobes	Sheet: 1 of 13	A
Checked By	Date				Scale:    N/A	
Part Number: 8800591						

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## INTENDED USE

**BBL™ Sensi-Disc™** Metronidazole Elution Discs are to be used for determining the *in vitro* susceptibility of anaerobic bacteria to metronidazole by the broth-disc method.

## SUMMARY AND EXPLANATION

In 1959 Schneierson<sup>1</sup> described a method for testing the antibiotic susceptibility of bacteria. In this method a known concentration of antibiotic is impregnated on a paper disc and placed in a tube of broth. The antibiotic concentration approximates the concentration of the antibiotic achievable in the blood. This method was named the "broth-disc method." In 1972 Wilkins and Thiel<sup>2</sup> modified this method for testing anaerobes by placing the discs in a pre-reduced medium. The major advantage of the broth-disc method is that most anaerobic isolates can be tested for antibiotic susceptibility. The major disadvantage is that the degree of susceptibility or resistance cannot be determined. In 1990 the NCCLS no longer listed this method as an approved alternative method for susceptibility testing of anaerobes due to its unreliability and poor correlation with standard methods.

## PRINCIPLES OF THE PROCEDURE

A paper disc containing a specific concentration of an antimicrobial is added to a broth medium. The antimicrobial agent is eluted into the broth from the paper disc. The test organism is inoculated into a tube containing the broth medium and the antimicrobial disc and into a control tube containing only the broth medium. After suitable incubation, susceptibility is determined by comparing the amount of growth in the tube with the antimicrobial disc to the amount of growth in the control tube.

## REAGENTS

**Sensi-Disc** brand discs are prepared by impregnating high quality absorbent paper with accurately determined amounts of antimicrobial agents. Discs are clearly marked on both sides with letters and numbers designating the agent and the drug content. **Sensi-Disc** agents are furnished in cartridges containing 50 discs each. The last disc in each cartridge is marked "X" and contains the drug as coded. The cartridges are for use in Single Disc Dispensers.

## Warnings and Precautions:

For *in vitro* Diagnostic Use.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, prepared tubes, specimen containers and other contaminated materials must be sterilized by autoclaving before discarding.

In 1990 the NCCLS no longer listed this method as an approved alternative method for susceptibility testing of anaerobes due to its unreliability and poor correlation with standard methods.

**Storage Instructions:** **Sensi-Disc** Metronidazole Elution Discs should be stored in the dark at -20 to +8°C.

**Product Deterioration:** Discard expired discs. Also discard cartridges from which discs have been frequently removed during a week. Discard containers left out overnight in the laboratory, or else test the discs for performance. The expiration date applies to discs in intact container stored as directed. Do not open until ready to use.

## SPECIMEN COLLECTION AND HANDLING

Specimens suitable for culture may be handled using various techniques. For detailed information, consult appropriate texts.<sup>3,4</sup> Specimens should be obtained before antimicrobial agents have been administered. Provision must be made for prompt delivery to the laboratory.

## PROCEDURE

**Material Provided:** **Sensi-Disc** Metronidazole Elution Discs.

**Materials Required But Not Provided:** Ancillary culture media, quality control organisms and laboratory equipment as required for this procedure.

## Culture Medium

Thioglycollate Medium without Indicator – 135C, supplemented with 5 µg/mL of hemin before autoclaving, and 0.1 µg/mL of vitamin K<sub>1</sub> and 1 mg/mL of sodium bicarbonate just prior to use. Medium should be dispensed in 5 mL quantities in 13 x 100 mm screw-capped test tubes and autoclaved at 121°C for 15 min.

## Inoculum

The inoculum should be prepared by incubating the organism at 35 to 37°C in Thioglycollate Medium without Indicator – 135C with hemin, vitamin K<sub>1</sub> and sodium bicarbonate until visible turbidity has developed. This culture should then be diluted with a sterile clear broth (boiled and cooled just prior to use) to match a barium sulfate turbidity standard (0.5 McFarland standard).

The barium standard may be prepared by adding 0.5 mL of 0.048M BaCl<sub>2</sub> (1.175% w/v BaCl<sub>2</sub> • 2H<sub>2</sub>O) to 99.5 mL of 0.36N H<sub>2</sub>SO<sub>4</sub> (1% v/v). The barium sulfate suspension may be distributed into screw-capped tubes that are sealed tightly and stored in the dark at room temperature for future use. This standard barium sulfate suspension must be mixed vigorously immediately prior to use.

## Test Procedure:

1. Set up and label three tubes of reduced medium.
2. Add one **Sensi-Disc** Metronidazole Elution Disc to each of tubes #2 and #3. Tube #1 serves as a culture control and does not receive a disc. Tube #3 serves as a sterility control for the metronidazole discs.
3. Incubate the tubes at 35 to 37°C for 2 h to permit equilibration of the drug throughout the medium.
4. Inoculate each of tubes #1 and #2 with 0.1 mL of the diluted inoculum as in **Inoculum** section above.
5. Loosen caps and incubate the tubes at 35 to 37°C in an anaerobic atmosphere (e.g., anaerobic glove box, anaerobic jar, etc.).

The length of the incubation period will vary, depending on the rapidity of growth of the test organism and can be determined by the appearance of visible turbidity in tube #1.

**User Quality Control:** Control strains recommended for standardized susceptibility testing. Each time the test is performed, the following strains should be included: *Clostridium perfringens* ATCC™ 13124 (susceptible) and *Propionibacterium acnes* ATCC 6919 (resistant).

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

## RESULTS

The degree of susceptibility of the test organism is determined by recording the presence or absence of growth, as judged by visible turbidity, in the three tubes. For the test to be valid, growth must be present in tube #1 and there must be no growth in tube #3. Growth in tube #3 would indicate disc contamination or faulty technique.

1. If growth is present in tube #1 and not in tube #2 and #3, the organism is inhibited by ≤ 16 µg of metronidazole per milliliter of medium and is considered susceptible.
2. If growth is present in tubes #1 and #2 and is not present in tube #3, the organism is considered to be resistant to metronidazole (>16 µg/mL).
3. If growth is absent from tube #1 or present in tube #3, the test is invalid and should be repeated.

## LIMITATIONS OF THE PROCEDURE

Contaminating organisms may be present in tube #2 even though tube #3 (sterility control) shows no growth; this could be falsely interpreted as resistance of the test organism to metronidazole.

## EXPECTED VALUES

Cumulative Percentage of Inhibition of Anaerobic Bacteria by Metronidazole Chemical.

Organisms	Number of Strains	Percentage Inhibited at				References
		≤ 4 µg/mL	≤ 8 µg/mL	≤ 16 µg/mL	≥ 100 µg/mL only	
<i>Bacteroides fragilis</i> group	679	88	95	98	1	5-15
<i>Bacteroides</i>	476	89	92	93	5	6,11,13-17
<i>Fusobacterium</i>	146	98	99	99	1	11,13-15,18,19
<i>Clostridium perfringens</i>	336	96	99	99	0	11,13-15,20
<i>Clostridium</i>	384	82	92	97	1	11-15,19-21
<i>Peptococcus</i>	253	85	87	87	11	11,14
<i>Peptostreptococcus</i>	177	68	70	72	21	11,14
<i>Veillonella</i>	39	74	79	87	13	11,13
<i>Propionibacterium</i>	203	2	3	3	97	11,13,14
<i>Eubacterium</i>	111	64	68	68	26	11,13,14
<i>Lactobacillus</i>	36	50	56	56	39	11,14
<i>Actinomyces</i>	20	20	20	25	45	11,14
<i>Bifidobacterium</i>	25	36	52	60	32	11,14

## PERFORMANCE CHARACTERISTICS

In clinical trials utilizing discs for determining susceptibility of anaerobes to metronidazole, the incidence of false negatives for metronidazole elution discs was 1.6% compared to the method described by Tally.<sup>22</sup>

## AVAILABILITY

Cat. No.	Description
231604	<b>BBL™ Sensi-Disc™</b> Metronidazole Elution Discs, 80 µg, single cartridge of 50 discs
231605	<b>BBL™ Sensi-Disc™</b> Metronidazole Elution Discs, 80 µg, Pkg. of 10 cartridges

## REFERENCES

- Schneierman, S.S. 1954. A simple rapid disc-tube method for determination of bacterial sensitivity to antibiotics. *Antibiot. Chemother.* 4:125-132.
- Wilkins, T.D. and T. Thiel. 1972. Modified broth-disk method for testing the antibiotic susceptibility of anaerobic bacteria. *Antimicrob. Agents Chemother.* 3:350-356.
- Murray, P.R., E.J. Baron, M.A. Tenover, and R.H. Tenover (ed.). 1999. *Manual of clinical microbiology*, 7th ed. American Society for Microbiology, Washington, D.C.
- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 1998. *Bailey & Scott's diagnostic microbiology*, 10th ed. Mosby, Inc., St. Louis
- Dornbusch, K., C.E. Nord, and T. Wadstrom. 1974. *Scand. J. Infect. Dis.* 6:253-258.
- Hohne, C. 1977. *Gesundheitswes. Disinfekt.* 32:1146-1149.
- Nastro, L.J., and S.M. Finegold. 1972. *J. Infect. Dis.* 126:104-107
- Jokipii, L., and A.M.M. Jokipii. 1977. *J. Antimicrob. Chemother.* 3:571-578.
- Santorio, J., D. Kaye, and M.E. Levison. 1976. *Antimicrob. Agents Chemother.* 10:188-190.
- Whelan, J.P.F., and J.H. Hale. 1973. *J. Clin. Pathol.* 26:393-395.
- Chow, A.W., D. Bednorz, and L.B. Guze. 1977. *Excerpta Medica*, p. 286-292.
- Dubois, J., J.C. Pechere, and P. Turgeon. 1978. *J. Antimicrob. Chemother.* 4:329-334
- Wust, J. 1977. *Antimicrob. Agents Chemother.* 11:631-656.
- Sutter, V.L., and S.M. Finegold. 1977. *Excerpta Medica*, p. 279-285.
- Staneck, J.L., and J.A. Washington II. 1974. *Antimicrob. Agents Chemother.* 6:311-315.
- Murray, P.R., and J.E. Rosenblatt. 1977. *Antimicrob. Agents Chemother.* 11:605-608.
- Wise, R., J.M. Andrews, and K.A. Bedford. 1977. *J. Infect. Dis.* 23:19-24.
- Rodriguez, J.A.C., J.E.G. Sanchez, M.C.S. Gonzales, and J.P. Prieto. 1977. *Pharmatherapeutica* 1:573-582.
- Henderson, D.K., A.W. Chow, and L.B. Guze. 1977. *Antimicrob. Agents Chemother.* 11:679-695.
- Dornbusch, K., and C.E. Nord. 1974. *Med. Microbiol. Immunol.* 160:265-267.
- Tally, F.P., et al. 1974. *Antimicrob. Agents Chemother.* 5:589-593.
- Tally, F.P., N.V. Jacobus, and S.J. Gorbach. 1978. *Antimicrob. Agents Chemother.* 14:436-438.

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**EC REP** Benex Limited