Revisions

SO 0191-5

Rev from	Rev to	ECO #
0604	2010/07	5384-10

Notes:

1. BD Cat. Number 231046, 231047, 231048, 231554
2. Blank (Sheet) Size : Length: N/A Width: N/A
 Number of Pages: N/A Number of Sheets: N/A
 Page Size: Length N/A Width N/A Final Folded Size: N/A
3. Style (see illustrations below): # N/A

4. See Specification Control Number VS-88-0068-1 for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: N/A PMS# N/A
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

VENDOR SPEC CONTROLLED BY BD CARIBE, LTD.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8800681		Category and Description Package Insert, BBL Taxo Discs for Differentiation of Pneumococci	Sheet: 1 of 11 <hr/> Scale: N/A	A

BD BBL™ Taxo™ Discs for Differentiation of Pneumococci

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INTENDED USE

Taxo™ P Discs are recommended for use in the presumptive identification of pneumococci, and are intended for use with pure cultures.

SUMMARY AND EXPLANATION

Taxo P Discs are impregnated with ethylhydrocupreine hydrochloride (optochin), a drug used for pneumonia therapy before sulfonamides became available. The growth of pneumococci, but not of other streptococci, is markedly inhibited by this chemical.^{1,2} Pneumococci may, therefore, be differentiated from other alpha-hemolytic streptococci by the formation of a zone of inhibition around a **Taxo P disc** placed on a blood agar plate heavily inoculated with a pure culture suspected to be *Streptococcus pneumoniae*.²

PRINCIPLES OF THE PROCEDURE

Taxo paper discs are impregnated with various chemicals as an aid in the rapid differentiation and presumptive identification of microorganisms. Evidence may be obtained from the reactions or responses to the impregnated agents to indicate the taxonomy of the organism under test.

The **Taxo P Disc** on a blood agar plate can be used for the rapid presumptive identification of pneumococci by observing zones of inhibition and hemolytic reactions.

REAGENTS

Taxo P Discs are 6 mm discs made from high quality absorbent paper impregnated with approximately 5.0 µg of ethylhydrocupreine hydrochloride per disc. The last disc in each cartridge is marked "X" and contains the agent as coded.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, test plates and other contaminated materials must be sterilized by autoclaving before discarding.

Storage Instructions: On receipt, store at -20 to +8°C. After use, store vial or cartridge at 2 to 8°C to protect product integrity. The expiration date applies to unopened containers stored as directed. Do not open until ready to use.

Use oldest discs first and discard expired discs. Allow containers to come to room temperature before opening. Return unused discs to the refrigerator. Discard containers which have been left out overnight.

SPECIMENS

These discs are not for use directly with clinical specimens or other sources containing mixed flora. The organism to be presumptively identified must first be isolated as separate colonies by streaking the specimen onto appropriate culture media; e.g., **Trypticase™ Soy Agar with 5% Sheep Blood (TSA II)**.

PROCEDURE

Material Provided: **Taxo P Discs**.

Materials Required But Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Test Procedure:

1. With sterile forceps or Single Disc Dispenser, place a **Taxo P disc** onto a **Trypticase Soy Agar with 5% Sheep Blood (TSA II)** subculture plate which has been heavily inoculated with a pure culture of the test organism exhibiting alpha hemolysis on the primary isolation plate.
2. Incubate plate(s) aerobically at 35 ± 2°C for 24 h, or as needed to obtain good growth; incubation in a CO₂-enriched atmosphere will enhance growth but reduce zone size.^{3,4}
3. Measure the diameter of the zone obtained with a millimeter ruler or caliper.

User Quality Control: At the time of use, check performance with pure cultures of stable control organisms producing known desired reactions. The use of *S. pneumoniae* ATCC™

6305 is recommended to demonstrate zone formation (14 mm or larger in diameter) and *S. pyogenes* ATCC 12384 is recommended for use as a negative (no zone) control.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Zones of inhibition of 14 mm or more are formed with pure cultures of *S. pneumoniae*. Other organisms may show zone sizes less than 14 mm in diameter. A diameter between 6 and 14 mm is questionable for pneumococci and the strain should be presumptively identified as a pneumococcus only if it is bile soluble.^{4,5}

LIMITATIONS OF THE PROCEDURE

Taxo P Disc tests are presumptive. Positive results should be followed with more specific tests. Additional tests for pneumococci include bile solubility, the Neufeld (quellung) reaction and biochemical and serological procedures.^{4,5}

Occasional strains of pneumococci not inhibited by optochin have been reported and strains of alpha streptococci have been reported to form zones of approximately 10 to 12 mm when light inocula were used.⁶

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of **Taxo™ P** Discs are tested for performance characteristics. Representative samples of the lot are assayed for ethylhydrocuprein hydrochloride (optochin) content using a chemical assay procedure.

Additionally, plates of **Trypticase** Soy Agar with 5% Sheep Blood are swab-inoculated with cultures diluted 10⁻¹ of *Streptococcus pneumoniae* (ATCC 6305), *S. pneumoniae* (Wilson), *S. pneumoniae* (SPA-18), *S. pneumoniae* (III), *Streptococcus pyogenes* (ATCC 12384) and *Streptococcus faecalis* (ATCC 10541). Representative samples of the lot are placed on the inoculated plates and incubated at 35 ± 2 °C for 24 h. The plates are then read for zone size around the **Taxo P** Disc with a millimeter ruler. The average zone size with *S. pneumoniae* cultures is at least 14 mm in diameter. No zones are observed around the **Taxo P** Disc with the *S. pyogenes* and *S. faecalis* cultures.

AVAILABILITY

Cat. No.	Description
231046	Taxo™ P Discs, Single vial of 50 discs.
231047	Taxo™ P Discs, Pkg. of 6 vials.
231048	Taxo™ P Discs, Single cartridge of 50 discs.
231554	Taxo™ P Discs, Pkg. of 10 cartridges.

REFERENCES

1. Bowers, E.F., and L.R. Jeffries. 1955. Optochin in the identification of *Str. pneumoniae*. J. Clin. Pathol. 8:58-60.
2. Bower, M.K., L.C. Thiele, B.D. Stearman, and I.G. Schaub. 1957. The optochin sensitivity test: a reliable method for identification of pneumococci. J. Lab. Clin. Med. 49:641-642.
3. Ragsdale, A.R., and J.P. Sanford. 1971. Interfering effect of incubation in carbon dioxide on the identification of pneumococci by optochin discs. Appl. Microbiol. 22:854-855.
4. Ruoff, K.L., R.A. Whaley, and D. Beighton. 1999. *Streptococcus*, p. 283-296. In P.R. Murray, E.J. Baron, M.A. Tenover, and R.H. Tenover (ed.). Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
5. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis.
6. Hall, C.T. May 7, 1970. Summary analysis of results for the proficiency testing survey in bacteriology (Jan. 9, 1970). National Communicable Disease Center, Atlanta.

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Benex Limited